

Agilent MassHunter BioConfirm Software

Familiarization Guide

Confirming Protein Molecular Weight

Exercise 1. Interactive Protein Molecular Weight Determination [3](#)

Exercise 2. Viewing Compound Information [8](#)

Sequence Matching Workflow [13](#)

Confirming Protein Identity

Exercise 3. Creating a Protein Sequence File [14](#)

Exercise 4. Interactive Protein Sequence Matching [15](#)

Exercise 5. Automated Protein Sequence Matching [18](#)

Confirming Protein Digest Identity

Exercise 6. Interactive Protein Digest Sequence Matching [20](#)

Exercise 7. Automated Protein Digest Sequence Matching [25](#)

Confirming Synthetic Peptide Identity

Exercise 8. Interactive Synthetic Peptide Sequence Matching [27](#)

Exercise 9. Automated Synthetic Peptide Sequence Matching [30](#)

Confirming Oligonucleotide Identity

Exercise 10. Interactive Oligonucleotide Sequence Matching [32](#)

Exercise 11. Automated Oligonucleotide Sequence Matching [35](#)

Comparing Protein Digest Files

Exercise 12. Comparing Protein Digest Files [37](#)



Agilent Technologies

Where to find more information

- *Agilent MassHunter BioConfirm Software Quick Start Guide*
- Online help provides in-depth information and can be displayed in the following ways:
 - Click **Contents**, **Index**, or **Search** from the Qualitative Analysis software Help menu.
 - Press the **F1** key to get more information about a window or dialog box.

How to use this guide

Try to do these familiarization exercises initially using the steps listed in the first column. Then if you need more information, follow the detailed instructions in the second column.

Exercise 1. Interactive Protein Molecular Weight Determination

This exercise shows you how to open a data file, integrate the chromatogram, extract spectra, deconvolute and view results. Deconvolution software does charge state deconvolution of mass spectra of large molecules with high charge states, such as proteins and large oligonucleotides.



Before you start

Copy the data file used for Exercises 1, 2, 4, and 5 onto your hard disk as follows:

- 1 Copy the **myoglobin.d** data file from the **Data** directory on the Qualitative Analysis setup disk to your computer hard drive.
- 2 Make sure you have both read and write permissions for the folder you just created on your computer. This is required if you want to save results.
 - a In Windows Explorer, right-click the **myoglobin.d** folder and click **Properties** from the shortcut menu.
 - b *Clear* the **Read-only Attributes** check box if it is marked.
 - c In the Confirm Attribute Changes dialog, click **Apply changes to this folder, subfolders, and files**, then click **OK**.

Steps	Detailed Instructions	Comments
1 Open the data file.	<ol style="list-style-type: none"> a Click File > Open Data File. b Locate the myoglobin.d folder. c Click Open. 	The TIC is automatically displayed in the Chromatogram Results window.
2 Integrate and extract peak spectra.	Right-click on the TIC and click Integrate and Extract Peak Spectra from the shortcut menu. See Figure 1 on page 5.	<i>Alternate method:</i> Click Actions > Integrate and Extract Peak Spectra .
3 Open the Deconvolute (MS) Method Editor section.	Select Deconvolute (MS): Maximum Entropy from the BioConfirm Workflow section of the Method Explorer.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.
4 Set the deconvolution range to 16000-18000 Da.	On the Deconvolution tab of the Deconvolute (MS) section of the Method Editor, enter 16000-18000 for Mass range .	See Figure 2 on page 6.

Exercise 1. Interactive Protein Molecular Weight Determination

Steps	Detailed Instructions	Comments
5 Set the mass step to 0.1 Da.	On the Deconvolution tab of the Deconvolute (MS) section of the Method Editor, enter 0.1 for Mass step .	See Figure 2 on page 6.
6 Select the extracted MS peak spectrum.	Click on the spectrum in the MS Spectrum Results window.	
7 Deconvolute the spectrum.	Click  on the Method Editor toolbar to start the deconvolution process.	Tip: Steps 2 and 7 can be combined by clicking the Chromatograms > Integrate and Deconvolute Peak Spectra menu item.
8 Review deconvolution results.	The results appear in the Deconvolution Results window. See Figure 3 on page 6.	For information on changing the display of data in the Deconvolution Results window, see <i>online help</i> .
9 View peak information.	<p>a Click on the spectrum in the Deconvolution Results window to select it. Right-click on the spectrum and click View MS Peak List 1 from the shortcut menu.</p> <p>b Click on the Abund. column heading to sort results by abundance.</p>	Mass (m/z), Abundance, and Fit score are listed for each peak in the spectrum. See Figure 4 on page 7.
10 Close the peak list.	Click  on the main toolbar.	
11 View compound information for the deconvoluted spectrum.	See Exercise 2 on page 8 .	

Exercise 1. Interactive Protein Molecular Weight Determination




Figure 1 Results of integration and spectra extraction for myoglobin.d

Exercise 1. Interactive Protein Molecular Weight Determination




Figure 2 Deconvolution parameters for myoglobin.d




Figure 3 Results of deconvolution for myoglobin.d

Exercise 1. Interactive Protein Molecular Weight Determination

MS Peaks One: + Scan (3.168-3.666 min) Deco... X

Mass	Abund.	Fit
16951.42	812466	10
16936.2	59821	9
16935.84	59549	9
16935.84	55998	9
16973.91	46128	9
16959.41	39141	10
16961.58	38535	9
16962.22	38320	9
16962.81	38026	9
16960.08	37925	10
16960.8	37872	10
16963.41	36696	9
16982.15	22195	9
16894.76	21634	8
16908.74	21288	8
17005.26	20242	8
17005.7	20201	8
16921.09	18846	8
16994.99	15786	8
16981.47	12950	9
17015.55	11679	8
16880.88	10395	8
16877.5	9880	8
17029.43	9653	7


Figure 4 Peak information for the deconvoluted spectrum for myoglobin.d (partial list, sorted by Abundance)

Exercise 2. Viewing Compound Information

This exercise shows you how to view compound information for deconvoluted spectra.

Steps	Detailed Instructions	Comments
1 Deconvolute myoglobin.d spectrum.	See “Exercise 1. Interactive Protein Molecular Weight Determination” on page 3.	You do not need to repeat the deconvolution steps, if you have already done them in Exercise 1.
2 View the compound list.	Click  on the main toolbar to display the Compound List window. Alternate method: Click Window Layouts > Load Layout on the Configuration menu, select BioConfirm-IntactProtein-MaximumEntropy-Default and click Open . This opens and reformats the Compound List to show the appropriate information for a deconvolution operation.	See Figure 5 on page 10.
3 Click on mass 16951.5 in the compound list.		If linked navigation is turned on, associated data in the following windows are automatically displayed and selected: <ul style="list-style-type: none"> • Deconvolution Results window • A compound spectrum that displays all the charge states from the original m/z data for that specific protein mass in the MS Spectrum Results window
4 Select the ion set spectrum in the Spectrum Results window for the mass 16951.5.		
5 View the charge states found for the protein along with their ppm error in the MS Peak List 2 window.	Click  on the main toolbar to open the MS Peak List 2 window after clicking on the spectrum to select it.	The following information is displayed for the ion set spectrum: <ul style="list-style-type: none"> • Mass • Abundance • Charge state • Diff (ppm) See Figure 6 on page 10.

Exercise 2. Viewing Compound Information

Steps	Detailed Instructions	Comments
6 Switch from List mode to Overlay mode in the MS Spectrum Results window.	Click  on the toolbar in the MS Spectrum Results window.	See Figure 7 on page 11.
7 Zoom in on the m/z 848.5 peak in the raw data	Right-drag to expand the area around m/z 848.5 in the MS Spectrum Results window.	See Figure 8 on page 11.
8 Select compound 1 in the compound list.	Click on the first line of the Compound List table.	Notice that the ion set peak for that mass shows a peak label. See Figure 9 on page 12.
9 Select compound 2 in the compound list.	Click on the second line of the Compound List table.	Notice that the ion set peak for another peak is highlighted. See Figure 10 on page 12.
10 Print a compound report.	<ol style="list-style-type: none"> Display the Compound Report section in the Method Editor by selecting Compound Report from the BioConfirm Workflow section in the Method Explorer. Review the options in this section. Verify that the sections that you want included in the report are marked. Display the Common Reporting Options section in the Method Editor by selecting Common Reporting Options from the BioConfirm Workflow section in the Method Explorer. Review the parameters in both the Templates and Options tabs. Click Compound Report from the File > Print menu to print the report. 	

Exercise 2. Viewing Compound Information

Compound List							
Show/Hide	File	RT	Mass	Fit	Min Z	Max Z	Height
<input checked="" type="checkbox"/>	myoglobin.d		16951.4187	10	8	27	471668
<input checked="" type="checkbox"/>	myoglobin.d		16934.9128	9	10	27	38105
<input checked="" type="checkbox"/>	myoglobin.d		16973.7846	9	9	23	27258
<input checked="" type="checkbox"/>	myoglobin.d		16994.468	8	11	26	21525
<input checked="" type="checkbox"/>	myoglobin.d		17006.9	8	16	23	16767
<input checked="" type="checkbox"/>	myoglobin.d		16894.7621	8	10	26	13741
<input checked="" type="checkbox"/>	myoglobin.d		16908.7668	8	10	26	13443
<input checked="" type="checkbox"/>	myoglobin.d		17030.5653	8	9	26	11999

Figure 5 Left half of Compound List window for myoglobin.d

MS Peaks Two: + Scan						
m/z	Abund	Charge	Calc Mass	Cpd Mass	Diff (ppm)	
628.8431	1817	27	16951.58	16951.42	-9.25	
652.9813	8770	26	16951.33	16951.42	5.18	
679.0674	28136	25	16951.51	16951.42	-5.44	
707.3181	61747	24	16951.47	16951.42	-2.86	
738.0282	103618	23	16951.49	16951.42	-4.12	
771.5311	143514	22	16951.53	16951.42	-6.54	
808.2222	176393	21	16951.52	16951.42	-5.91	
848.5842	192965	20	16951.54	16951.42	-7.36	
893.1934	188626	19	16951.54	16951.42	-7.2	
942.7561	164620	18	16951.48	16951.42	-3.8	
998.151	122827	17	16951.45	16951.42	-1.72	
1060.4724	88543	16	16951.45	16951.42	-1.64	
1131.1031	69263	15	16951.44	16951.42	-1.35	
1211.8224	55680	14	16951.42	16951.42	0.21	
1304.9627	37144	13	16951.42	16951.42	-0.29	
1413.6242	22264	12	16951.41	16951.42	0.7	
1542.0526	16663	11	16951.5	16951.42	-4.89	
1696.1567	8049	10	16951.5	16951.42	-4.63	
1884.5227	1445	9	16951.64	16951.42	-13.14	
2119.9559	249	8	16951.59	16951.42	-10.15	

Figure 6 MS Peaks Two window for myoglobin.d

Exercise 2. Viewing Compound Information




Figure 7 MS Spectrum Results window for myoglobin.d (Overlay Mode)




Figure 8 MS Spectrum Results window for myoglobin.d (zoomed)

Exercise 2. Viewing Compound Information




Figure 9 Ion set peak label for Compound 1




Figure 10 Ion set peak label for Compound 2

Sequence Matching Workflow

The steps outlined below show the workflow for sequence matching with Agilent MassHunter BioConfirm Software.

Step 1 - Open the data file of interest.

Step 2 - Open a Qualitative Analysis Method or create a new one.

Step 3 - Find compounds by molecular feature or by integration and deconvolution.

Step 4 - Select the sequences to match.

If the sequence you want to match is not in the method, then:

- Import or create a sequence.

Step 5 - Edit sequences if necessary:

- Set the sequence type: Protein, Synthetic peptide, Protein Digest, Oligonucleotide.
- Add or edit the sequence text.
- Apply or edit modifications (*not available for Oligonucleotide sequences in this version*)
- Apply or edit links (*not available for Oligonucleotide sequences in this version*)
- Assign or edit digest reagents (*Protein Digest sequences only*).
- Select matching rules.

Step 6 - Review/set other match sequence method parameters.

Step 7 - Start the sequence matching process.

Step 8 - Review the results in the Compound Identification Results windows.

Step 9 - For protein digests only:

- View sequence coverage results in the Sequence Coverage Map window.

Step 10 - Print report.


Exercise 3. Creating a Protein Sequence File

This exercise guides you through the creation of a myoglobin sequence file that you will use in “[Exercise 4. Interactive Protein Sequence Matching](#)” on page 15 and “[Exercise 5. Automated Protein Sequence Matching](#)” on page 18.


Steps	Detailed Instructions	Comments
1 Display the Define and Match Sequences section in the Method Editor window.	Click BioConfirm Workflow > Define and Match Sequences in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.
2 Create a new sequence.	a Click Sequence > New Sequence . b The Sequence Editor window opens automatically with a new sequence displayed for editing.	The new sequence is also added to the list of sequences in the Sequences tab of the Define and Match Sequences section of the Method Editor window.
3 Enter a Sequence Name .	In the Sequence Editor window, type in Myoglobin in the Sequence Name box.	
4 Select Protein as the Sequence Type .	In the Sequence Editor window, select Protein as the Sequence Type.	
5 Enter the amino acid sequence shown below into the Sequence Editor box.	Type in individual amino acids one at a time between the N-term and C-term symbols.	Use the single-character (letter) amino acids abbreviations, as shown in the Amino acid list on the left side of the Sequence Editor window.
GLSDGEWQQVLNVWGKVEADIAGHGQEV LIRLFTGHPETLEK FDKFKHLKTEAEMKASEDLKKHGT VVLTALGGILKKKGHHEAE LKPLAQSHATKHKIPIKYLEFISDAIIHVLHSHKHPGDFGADAQG AMTKALELFRNDIAAKYKELGFQG		Tip: If you are reading this document as a PDF file on your computer, you can copy and paste the sequence into the Sequence Editor window.
Note: The myoglobin sequence does not have any links or modifications, but some sequences do. In that case, add links and modifications as described in the <i>Quick Start Guide</i> or <i>online help</i> .		
6 Save the sequence as the name <i>iii_myoglob.psq</i> , where <i>iii</i> represents your initials.	a Click Save as on the Sequences tab of the Define and Match Sequences section of the Method Editor window. b Type <i>iii_myoglob</i> in the File name box. c Click Save .	The sequence is saved as a .psq file that can be imported for use in other methods as described in Exercise 4 or referenced from worklists as described in Exercise 5.

Exercise 4. Interactive Protein Sequence Matching


This exercise shows you how to set method parameters, match an intact protein sequence, and view the results. This exercise uses the *iii_myoglobin.psq* sequence file created in Exercise 3 and the **myoglobin.d** data file copied in Exercise 1.

Steps	Detailed Instructions	Comments
1 Open the method to use as a starting point for the new method.	<ol style="list-style-type: none"> Click Method > Open. Select the BioConfirmIntactProtein-Default.m folder. Click Open. 	
2 If the myoglobin.d data file is not already open, open it.	<ol style="list-style-type: none"> Click File > Open Data File. Locate the myoglobin.d folder. Click Open. 	The TIC is automatically displayed in the Chromatogram Results window.
3 Display the Find Compounds by Molecular Feature section in the Method Editor window.	Select BioConfirm Workflow > Find by Molecular Feature in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.
4 Find compounds.	<ol style="list-style-type: none"> Review the settings and modify them if necessary. Click  on the Method Editor toolbar to start the compound search. Review the results in the Compound List window. 	In this case we are using the default method parameters. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i> .
5 Change the layout to BioConfirm-IntactProtein-LMFE.	<ol style="list-style-type: none"> Click Window Layouts > Load Layout on the Configuration menu. Select BioConfirm-IntactProtein-LMFE.xml. Click Open. 	
6 Display the Define and Match Sequences section in the Method Editor window.	Click BioConfirm Workflow > Define and Match Sequences in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.

Exercise 4. Interactive Protein Sequence Matching


Steps	Detailed Instructions	Comments
7 Import the myoglobin sequence.	<ol style="list-style-type: none"> Click the Sequence tab in the Define and Match Sequences section of the Method Editor window. Click Import. Select iii_myoglobin.psq and click Open. 	<p>The iii_myoglobin.psq sequence file was created in Exercise 3.</p> <p>For this exercise, we will use the sequence as is, but you can add modifications and links to sequences as described in <i>online help</i> and the <i>Quick Start Guide</i>.</p>
8 Select protein matching rules.	<ol style="list-style-type: none"> Right-click in the Sequence Editor window and click Edit Matching Rules from the shortcut menu to open the Rules dialog box. Select the following tests to use for matching the theoretical masses of proteins to those from MS data: <ul style="list-style-type: none"> • Intact protein • Predicted modifications Click OK to close the Rules dialog box. 	<p>Use <Ctrl>+click to select multiple tests from the list.</p> <p>Note that these tests are already selected in the BioConfirmIntactProtein-Default method that you loaded in Step 1.</p>
9 Select Sequence as the match source.	<ol style="list-style-type: none"> Click the Source tab in the Define and Match Sequences section of the Method Editor window. Click Qualitative method and mark the Sequences check box. 	<p>Save the method for use in the Exercise 5 as follows:</p> <ol style="list-style-type: none"> Click Method > Save As. Type the File name iii_myoglobin.m, where <i>iii</i> represents your initials. Click Save.
10 Start the match search.	Right-click in the Compound List window and click Match Sequences from the shortcut menu.	<p>Alternate methods:</p> <ul style="list-style-type: none"> • Click  on the Method Editor toolbar. • Click Sequences > Match Sequences. • Click Match Sequences on the Method Editor shortcut menu. • Click Match Sequences on the Data Navigator Compounds shortcut menu.

Exercise 4. Interactive Protein Sequence Matching

Steps	Detailed Instructions	Comments
11 Review the results.	<p>a Click  on the toolbar to display the Compound Identification Results window.</p> <p>b When you open the window, the window displays the results for the first compound that is highlighted in the Data Navigator (that is marked to show).</p>	<p>Alternate method:</p> <ul style="list-style-type: none">• Click View > Compound Identification Results.

Exercise 5. Automated Protein Sequence Matching

This exercise guides you through the setup of a worklist to automatically confirm the presence of myoglobin in a previously acquired sample. This exercise uses the *iii_myoglob.psq* sequence file created in Exercise 3 and the *myoglobin.d* data file copied in Exercise 1.

Steps	Detailed Instructions	Comments
1 If not already open, open the method <i>iii_myoglobin.m</i> .	a Click Method > Open . b Select the <i>iii_myoglobin.m</i> folder. c Click Open .	This method was created in “Exercise 4. Interactive Protein Sequence Matching” on page 15.
2 Display the Define and Match Sequences section in the Method Editor window.	Click BioConfirm Workflow > Define and Match Sequences in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.
3 Select Worklist as the match source.	a Click the Source tab in the Define and Match Sequences section of the Method Editor window. b Click Worklist .	This selection causes the software to get the sequence from the worklist rather than the method as described in Exercise 4.
4 Display the Worklist Automation > Worklist Actions section in the Method Editor.	Click Worklist Automation > Worklist Actions in the Method Explorer window.	
5 Select the appropriate worklist actions.	Select the following worklist actions in the Available actions list: <ul style="list-style-type: none"> • Find Compounds by Molecular Feature • Match Sequences • Generate Compound Report 	In this case we are using the default method parameters. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i> .
6 Add the selected actions to the Actions to be run list.	Click the  button. The items are added to the end of the list.	Actions will be executed in the order they appear in the list. You can reorder them using the Up and Down arrow buttons to the right of the list
7 Save the method.	Click Method > Save .	

Exercise 5. Automated Protein Sequence Matching

Steps	Detailed Instructions	Comments
8 Create a worklist of one sample in the MassHunter Workstation Data Acquisition software.	a Display the Worklist pane. b Click Worklist > Add Sample . A new sample row is added to the Worklist table.	If you plan to do batch data analysis using the worklist, consider using the DA Reprocessor tool that is installed with Agilent MassHunter Workstation Software - Data Acquisition program. In that case, skip Step 12 below.
9 Specify the myoglobin sequence file, <i>iii_myoglobin.psq</i> , in the worklist.	a Click Worklist > Add Column . b When the Add Column dialog box appears, select Protein as the Column Type. c Enter the Column name as myoglobin . d Select the <i>iii_myoglobin.psq</i> file as the Value. e Click OK .	The <i>iii_myoglobin.psq</i> file was created in Exercise 3.
10 Enter <i>iii_myoglobin.m</i> in the Method column of the worklist.	The <i>iii</i> represents your initials.	This method was saved in Step 9 above.
11 Enter myoglobin.d in the Data File column of the worklist.		
12 Set up to run the worklist for data analysis only.	a Click Worklist > Worklist Run Parameters . b Select DA Only as Part of method to run . c Select the paths for the DA method and data file, then click OK .	
13 Run the worklist.	Click Worklist > Run .	
14 Review the printed Compound reports.		

Exercise 6. Interactive Protein Digest Sequence Matching

This exercise shows you how to confirm protein digests interactively.


Before you start

Copy the files used for Exercises 6 and 7 onto your hard disk as follows:

- 1 Copy the **enolase-Chip-final.d** data file from the **Data** folder on the Qualitative Analysis setup disk to the **MassHunter\Data** folder on your computer hard drive.
- 2 Copy the **EnolaseDigest.psq** sequence file from the **Data** folder on the Qualitative Analysis setup disk to the **MassHunter\ProteinSequences** folder on your computer hard drive.
- 3 Make sure you have both read and write permissions for the folder you just created on your computer. This is required if you want to save results.
 - a In Windows Explorer, right-click the **enolase-Chip-final.d** folder and click **Properties** from the shortcut menu.
 - b *Clear* the **Read-only Attributes** check box if it is marked.
 - c In the Confirm Attribute Changes dialog, click **Apply changes to this folder, subfolders, and files**, then click **OK**.

Steps	Detailed Instructions	Comments
1 Open the method to use as a starting point for the new method.	<ol style="list-style-type: none">a Click Method > Open.b Select the BioConfirmProteinDigest-Default.m folder.c Click Open.	
2 Open the demo data file.	<ol style="list-style-type: none">a Click File > Open Data File.b Locate the Enolase-Chip-final.d folder.c Click Open.	The TIC is automatically displayed in the Chromatogram Results window.





Exercise 6. Interactive Protein Digest Sequence Matching

Steps	Detailed Instructions	Comments
3 Review the parameters in the Find Compounds by Molecular Feature section in the Method Editor window.	<ol style="list-style-type: none"> Select BioConfirm Workflow > Find by Molecular Feature in the Method Explorer window. Review the settings on the various tabs of the Find Compounds by Molecular Feature sections of Method Editor. Click the Extraction tab and set the Mass range to 300-1700 <i>m/z</i>. In the Extraction tab, set the peak height filter to 500 counts. 	<ul style="list-style-type: none"> If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu. Change the default parameters as described in the next steps. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i>. A very low peak height filter can result in greater sequence coverage but requires much more time to process.
4 For MS/MS data, set parameters to extract MS/MS spectra.	<ol style="list-style-type: none"> Click the Results tab of Find Compounds by Molecular Feature in the Method Editor. Mark the Extract MS/MS Spectrum check box. Mark Deisotope MS/MS Spectrum. 	
5 Find compounds.	<ol style="list-style-type: none"> Click  on the Method Editor toolbar to start the compound search. When processing is complete, review the results in the Compound List window. 	
6 Change the layout to BioConfirm-ProteinDigest.	<ol style="list-style-type: none"> Click Window Layouts > Load Layout on the Configuration menu. Select BioConfirm-ProteinDigest.xml. Click Open. 	
7 Import the sequence.	<ol style="list-style-type: none"> Click BioConfirm Workflow > Define and Match Sequences in the Method Explorer window. Click the Sequence tab in the Define and Match Sequences section of the Method Editor window. Click Import. Select EnolaseDigest.psq. Click Open. 	<p>The enolase digest sequence is automatically displayed in the Sequence Editor window.</p> <p>For this exercise, we will use the sequence as is, but you can add modifications and links to sequences as described in <i>online help</i> and the <i>Quick Start Guide</i>.</p>

Exercise 6. Interactive Protein Digest Sequence Matching

Steps	Detailed Instructions	Comments
8 Assign or edit digest reagents.	<ol style="list-style-type: none"> Click the Edit button in the Define and Match Sequence dialog to open the Sequence Editor window. Right-click in the Sequence Editor window and click Edit Digest Reagents from the shortcut menu to open the Digest Reagents dialog box. Notice that the Reagent selected is Trypsin for EnolaseDigest. Set the maximum number of Missed Cleavages to allow to 2. Click OK. 	You can customize the list of available reagents using the Chemical Data Dictionary; see <i>online help</i> for more information.
9 View the digest list.	Right-click on the sequence in the Sequence Editor window and click Digest Current Sequence from the shortcut menu to digest the sequence and display the results in the Digest List window.	<i>Alternate method:</i> Click Sequence > Digest Current Sequence .
10 Select protein digest matching rules	<ol style="list-style-type: none"> To open the Rules dialog box, right-click in the Sequence Editor window and click Edit Matching Rules from the shortcut menu. Click Predicted Modifications, then under Selected list, click Carbamylation to move it to the Available list. Select the following tests: <ul style="list-style-type: none"> Complete Digest Incomplete Digest Predicted Modifications Click OK to close the Rules dialog box. 	<ul style="list-style-type: none"> Use <Ctrl>+click to select multiple tests from the list. The enolase sample was denatured thermally without using urea, so modification by carbamylation is not present. Note that these tests were automatically selected when you opened the BioConfirmProteinDigest-Default method in Step 1.
11 Set Match Sequence parameters.	<ol style="list-style-type: none"> Click the Source tab in the Define and Match Sequences section of the Method Editor window. Click Qualitative method and mark the Sequences check box. 	<p>For MS/MS data, you can adjust the following parameters on the Scoring tab:</p> <ul style="list-style-type: none"> MS/MS score to increase or decrease its contribution to the overall Score (Bio). MS/MS scored peak intensity and MS/MS matched ion score contribute to Score (Bio MS/MS).

Exercise 6. Interactive Protein Digest Sequence Matching


Steps	Detailed Instructions	Comments
12 Save the method for use in the Exercise 7.	<ol style="list-style-type: none"> a Click Method > Save As. b Type the File name <i>iii_Enolase-Chip-Final.m</i>, where <i>iii</i> represents your initials. c Click Save. 	
13 Start the match search.	Right-click in the Compound List window and click Match Sequences from the shortcut menu	<p><i>Alternate methods:</i></p> <ul style="list-style-type: none"> • Click  on the Method Editor toolbar. • Click Sequence > Match Sequences. • Click Match Sequences on the Method Editor shortcut menu. • Click Match Sequences on the Data Navigator Compounds shortcut menu.
14 Review the results.	<ol style="list-style-type: none"> a Click  on the toolbar to display the Compound Identification Results window. b When you open the window, the window displays the results for the first compound that is highlighted in the Data Navigator (that is marked to show). c Select another sequence match result to view by selecting a different compound in the Data Navigator or Compound List windows. 	
15 View sequence coverage results.	<ol style="list-style-type: none"> a Display the Sequence Coverage Map window by clicking  on the toolbar. b Click to select a different sequence coverage result in the Data Navigator window (under Matched Sequences). 	<p>Alternate method: You can also display the Sequence Coverage Map window in these ways:</p> <ul style="list-style-type: none"> • Click Sequence > View Sequence Coverage Map. • Click View > Sequence Coverage Map.
16 Save the results	<ol style="list-style-type: none"> a Click  to save your results to the data file folder. 	

Exercise 6. Interactive Protein Digest Sequence Matching

Steps	Detailed Instructions	Comments
17 Repeat the interactive processing with enolase-oxidized-chip-final.d.	<p>a Load the data file enolase-oxidized-chip-final.d (see step 2).</p> <p>b Select MFE and verify the processing parameters (step 3).</p> <p>c Find compounds (step 5).</p> <p>d Match sequences (step 13).</p> <p>e Save the results to the second data file (step 16).</p>	<p>Most of the processing parameters used for the first data file are the same for the second data file.</p> <p>The data files and results of this exercise will be used in "Exercise 12. Comparing Protein Digest Files" on page 37.</p>
To view more information.	<p>Click the following items on the Sequence Coverage Map window shortcut menu to view more information about the sequence:</p> <ul style="list-style-type: none">• Applied Modifications• Applied Links• Applied Reagents• Applied Matching Rules• Show Sequence Description	

Exercise 7. Automated Protein Digest Sequence Matching

This exercise guides you through the setup of a worklist to automatically confirm the presence of serotransferrin in a previously acquired sample.

Steps	Detailed Instructions	Comments
1 Open the method.	a Click Method > Open . b Select the <i>iii_Enolase-Chip-Final.m</i> folder. c Click Open .	This method was created in Exercise 6 (<i>iii</i> represents your initials).
2 Display the Worklist Automation > Worklist Actions section in the Method Editor.	Click Worklist Automation > Worklist Actions in the Method Explorer window.	
3 Select the appropriate worklist actions.	Select the following worklist actions in the Available actions list: <ul style="list-style-type: none"> • Find Compounds by Molecular Feature • Match Sequences • Generate Compound Report 	In this case we are using the default method parameters. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i> .
4 Add the selected actions to the Actions to be run list.	Click the  button. The items are added to the end of the list.	Actions will be executed in the order they appear in the list. You can reorder them using the Up and Down arrow buttons to the right of the list
5 Select Worklist as the match source.	a Click the Source tab in the Define and Match Sequences section of the Method Editor window. b Click Worklist .	This causes the software to get the sequence from the worklist rather than from the method.
6 Save the method.	Click Method > Save .	

Exercise 7. Automated Protein Digest Sequence Matching

Steps	Detailed Instructions	Comments
7 Create a worklist of one sample in the MassHunter Workstation Data Acquisition software.	<p>a Display the Worklist pane.</p> <p>b Click Worklist > Add Sample. A new sample row is added to the Worklist table.</p> <p>c Click Worklist > Add Column.</p> <p>d When the Add Column dialog box appears, select Protein as the Column Type.</p> <p>e Enter the Column name as EnolaseDigest.</p> <p>f Select EnolaseDigest.psq as the Value.</p> <p>g Click OK.</p>	If you plan to do batch data analysis using the worklist, consider using the DA Reprocessor tool that is installed with Agilent MassHunter Workstation Software - Data Acquisition program. In that case, skip Step 10 below.
8 Enter <i>iii</i> Enolase-Chip-Final.m in the Method column of the worklist.		The <i>iii</i> represents your initials.
9 Enter Enolase-Chip-Final.d in the Data File column of the worklist.		
10 Set up to run the worklist for data analysis only.	<p>a Click Worklist > Worklist Run Parameters.</p> <p>b Select DA Only as Part of method to run.</p> <p>c Select the paths for the DA method and data file, then click OK.</p>	
11 Run the worklist.	Click Worklist > Run .	
12 Review the printed Compound reports.		

Exercise 8. Interactive Synthetic Peptide Sequence Matching

This exercise shows you how to set method parameters, import a sequence, match a synthetic peptide sequence, and view the results.


Before you start

Copy the files used for Exercise 8 and 9 onto your hard disk as follows:



- 1 Copy the **SynPep3.d** data file from the **Data** folder on the Qualitative Analysis setup disk to the **MassHunter\Data** folder on your computer.
- 2 Copy the **SynPep3.psq** sequence file from the **Data** folder on the Qualitative Analysis setup disk to the **MassHunter\ProteinSequences** folder on your computer.
- 3 Make sure you have both read and write permissions for the data folder you just created on your computer. This is required if you want to save results.
 - a In Windows Explorer, right-click the **SynPep3.d** folder and click **Properties** from the shortcut menu.
 - b *Clear* the **Read-only Attributes** check box if it is marked.
 - c In the Confirm Attribute Changes dialog, click **Apply changes to this folder, subfolders, and files**, then click **OK**.

Steps	Detailed Instructions	Comments
1 Open the method BioConfirmSyntheticPeptide-Default to use as a starting point for the new method.	<ol style="list-style-type: none"> a Click Method > Open. b Select the BioConfirmSyntheticPeptide-Default.m folder. c Click Open. 	
2 Open the SynPep3.d data file.	<ol style="list-style-type: none"> a Click File > Open Data File. b Locate the SynPep3.d folder. c Click Open. 	The TIC is automatically displayed in the Chromatogram Results window
3 Display the Find Compounds by Molecular Feature section in the Method Editor window.	Select BioConfirm Workflow > Find by Molecular Feature in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.

Exercise 8. Interactive Synthetic Peptide Sequence Matching


Steps	Detailed Instructions	Comments
4 For MS/MS data, set parameters to extract MS/MS spectra.	<p>a Click the Results tab of Find Compounds by Molecular Feature in the Method Editor.</p> <p>b Mark the Extract MS/MS Spectrum check box.</p> <p>c Mark Deisotope MS/MS Spectrum.</p>	
5 Find compounds.	<p>a Review the settings and modify them if necessary.</p> <p>b Click  on the Method Editor toolbar to start the compound search.</p> <p>c Review the results in the Compound List window.</p>	In this case we are using the default method parameters. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i> .
6 Change the layout.	<p>a Click Window Layouts > Load Layout on the Configuration menu.</p> <p>b Select BioConfirm-SyntheticPeptide.xml.</p> <p>c Click Open.</p>	
7 Display the Define and Match Sequences section in the Method Editor window.	Click BioConfirm Workflow > Define and Match Sequences in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.
8 Import the sequence.	<p>a Click the Sequence tab in the Define and Match Sequences section of the Method Editor window.</p> <p>b Click Import.</p> <p>c Select SynPep3.psq.</p> <p>d Click Open.</p>	For this exercise, we will use the sequence as is, but you can add modifications and links to sequences as described in <i>online help</i> and the <i>Quick Start Guide</i> .
9 Select peptide matching rules.	<p>a Right-click in the Sequence Editor window and click Edit Matching Rules from the shortcut menu to open the Rules dialog box.</p> <p>b Select the following tests to use for matching the theoretical masses of oligonucleotides to those from MS data:</p> <ul style="list-style-type: none"> • Intact Peptide • Extra Amino Acid • Missing Amino Acid • Fmoc blocking groups <p>c Click OK to close the Rules dialog box.</p>	Use <Ctrl>+click to select multiple tests from the list.

Exercise 8. Interactive Synthetic Peptide Sequence Matching

Steps	Detailed Instructions	Comments
10 Select Sequences as the match source.	<ol style="list-style-type: none"> a Click the Source tab in the Define and Match Sequences section of the Method Editor window. b Click Qualitative method and mark the Sequences check box. 	<p>Save the method for use in the Exercise 9 as follows:</p> <ol style="list-style-type: none"> a Click Method > Save As. b Type the File name <i>iii_SynPep3.m</i>, where <i>iii</i> represents your initials. c Click Save.
11 Start the match search.	Right-click in the Compound List window and click Match Sequences from the shortcut menu.	<p>Alternate methods:</p> <ul style="list-style-type: none"> • Click  on the Method Editor toolbar. • Click Sequence > Match Sequences. • Click Match Sequences on the Method Editor shortcut menu. • Click Match Sequences on the Data Navigator Compounds shortcut menu.
12 Review the results.	<ol style="list-style-type: none"> a Click  on the toolbar to display the Compound Identification Results window. b When you open the window, the window displays the results for the first compound that is highlighted in the Data Navigator (that is marked to show). c Select another sequence match result to view by selecting a different compound in the Data Navigator or Compound List windows. 	<p>Alternate method:</p> <ul style="list-style-type: none"> • Click View > Compound Identification Results.

Exercise 9. Automated Synthetic Peptide Sequence Matching

This exercise guides you through the setup of a worklist to automatically confirm the presence of SynPep3 in a previously acquired sample.

Steps	Detailed Instructions	Comments
1 Open the <i>iii_SysPep3.m</i> method.	a Click Method > Open . b Select the <i>iii_SysPep3.m</i> folder. c Click Open .	This method was created in Exercise 8 (<i>iii</i> represents your initials).
2 Display the Worklist Automation > Worklist Actions section in the Method Editor.	Click Worklist Automation > Worklist Actions in the Method Explorer window.	
3 Select the appropriate worklist actions.	Select the following worklist actions in the Available actions list: <ul style="list-style-type: none"> • Find Compounds by Molecular Feature • Match Sequences • Generate Compound Report 	In this case we are using the default method parameters. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i> .
4 Add the selected actions to the Actions to be run list.	Click the  button. The items are added to the end of the list.	Actions will be executed in the order they appear in the list. You can reorder them using the Up and Down arrow buttons to the right of the list
5 Select Worklist as the match source.	a Click the Source tab in the Define and Match Sequences section of the Method Editor window. b Click Worklist .	This causes the software to get the sequence from the worklist rather than from the method.
6 Save the method.	Click Method > Save .	

Exercise 9. Automated Synthetic Peptide Sequence Matching

Steps	Detailed Instructions	Comments
7 Create a worklist of one sample in the MassHunter Workstation Data Acquisition software.	a Display the Worklist pane. b Click Worklist > Add Sample . A new sample row is added to the Worklist table. c Click Worklist > Add Column . d When the Add Column dialog box appears, select Protein as the Column Type. e Enter the Column name as SynPep3 . f Select SynPep3.psq as the Value. g Click OK .	If you plan to do batch data analysis using the worklist, consider using the DA Reprocessor tool that is installed with Agilent MassHunter Workstation Software - Data Acquisition program. In that case, skip Step 10 below.
8 Enter <i>iii</i> SynPep3.m in the Method column of the worklist.	The <i>iii</i> represents your initials.	
9 Enter SynPep3.d in the Data File column of the worklist.		
10 Set up to run the worklist for data analysis only.	a Click Worklist > Worklist Run Parameters . b Select DA Only as Part of method to run . c Select the paths for the DA method and data file, then click OK .	
11 Run the worklist.	Click Worklist > Run .	
12 Review the printed Compound reports.		

Exercise 10. Interactive Oligonucleotide Sequence Matching

This exercise shows you how to set method parameters, import a sequence, match an oligonucleotide sequence, and view the results.


Before you start

Copy the files used for Exercises 10 and 11 onto your hard disk as follows:



- 1 Copy the **DNA-2ug-r001.d** data file from the **Data** folder on the Qualitative Analysis setup disk to the **MassHunter\Data** folder on your computer hard drive.
- 2 Copy the **21mer_oligo.psq** sequence file from the **Data** folder on the Qualitative Analysis setup disk to the **MassHunter\ProteinSequences** folder on your computer hard drive.
- 3 Make sure you have both read and write permissions for the data folder you just created on your computer. This is required if you want to save results.
 - a In Windows Explorer, right-click the **DNA-2ug-r001.d** folder and click **Properties** in the shortcut menu.
 - b *Clear* the **Read-only Attributes** check box if it is marked.
 - c In the Confirm Attribute Changes dialog, click **Apply changes to this folder, subfolders, and files**, then click **OK**.

Steps	Detailed Instructions	Comments
1 Open the method.	<ol style="list-style-type: none"> a Click Method > Open. b Select the BioConfirmOligonucleotideSmall-Default.m folder. c Click Open. 	
2 Open the data file.	<ol style="list-style-type: none"> a Click File > Open Data File. b Locate the DNA-2ug-r001.d folder. c Click Open. 	The TIC is automatically displayed in the Chromatogram Results window
3 Display the Find Compounds by Molecular Feature section in the Method Editor window.	Select BioConfirm Workflow > Find by Molecular Feature in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.

Exercise 10. Interactive Oligonucleotide Sequence Matching


Steps	Detailed Instructions	Comments
4 Find compounds.	<p>a Review the settings and modify them if necessary.</p> <p>b In the Extraction tab, change Restriction retention time to 0.5-9.6 minutes, Restrict m/z to 800-2300 m/z, and Use peaks with height to 400 counts.</p> <p>c In the Charge state tab, change the Isotope model to unbiased and limit the charge state to a maximum of 10.</p> <p>d In the Results tab, select to extract an MFE spectrum and an ECC for each compound.</p> <p>e Click  on the Method Editor toolbar to start the compound search.</p> <p>f When processing is complete, review the results in the Compound List window.</p>	<p>These boundary conditions are used with the MFE procedure to avoid long processing times.</p> <p>For guidance to adjust these parameters, see the <i>Quick Start Guide</i> or <i>online Help</i>.</p>
5 Change the layout to BioConfirm-Oligonucleotide.	<p>a Click Window Layouts > Load Layout on the Configuration menu.</p> <p>b Select BioConfirm-Oligonucleotide.xml.</p> <p>c Click Open.</p>	
6 Display the Define and Match Sequences section in the Method Editor window.	Click BioConfirm Workflow > Define and Match Sequences in the Method Explorer window.	If the BioConfirm workflow is not available in Method Explorer, select it from the Configuration > Configure for Workflow menu.
7 Import the sequence.	<p>a Click the Sequence tab in the Define and Match Sequences section of the Method Editor window.</p> <p>b Click Import.</p> <p>c Select 21mer_oligo.psq.</p> <p>d Click Open.</p>	
8 Select oligonucleotide matching rules.	<p>a Right-click in the Sequence Editor window and click Edit Matching Rules from the shortcut menu to open the Rules dialog box.</p> <p>b Select the following matching rules to use for matching the theoretical masses of peptides to those from MS data:</p> <ul style="list-style-type: none"> • Intact oligonucleotide • Oligonucleotide truncation <p>c Click OK to close the Rules dialog box.</p>	<p>Use <Ctrl>+click to select multiple matching rules from the list.</p> <p>Note that these matching rules are already selected in the BioConfirmOligonucleotide-Small-Default method that you loaded in Step 1.</p>

Exercise 10. Interactive Oligonucleotide Sequence Matching

Steps	Detailed Instructions	Comments
9 Select Sequence as the match source.	<ol style="list-style-type: none"> Click the Source tab in the Define and Match Sequences section of the Method Editor window. Click Qualitative method and mark the Sequences check box. 	<p>Save the method for use in the Exercise 11 as follows:</p> <ol style="list-style-type: none"> Click Method > Save As. Type the File name <i>iii_Oligo.m</i>, where <i>iii</i> represents your initials. Click Save.
10 Start the match search.	Right-click in the Compound List window and click Match Sequences from the shortcut menu.	<p>Alternate methods:</p> <ul style="list-style-type: none"> Click  on the Method Editor toolbar. Click Sequence > Match Sequences. Click Match Sequences on the Method Editor shortcut menu. Click Match Sequences on the Data Navigator Compounds shortcut menu.
11 Review the results.	<ol style="list-style-type: none"> Click  on the toolbar to display the Compound Identification Results window. When you open the window, the window displays the results for the first compound that is highlighted in the Data Navigator (that is marked to show). Select another sequence match result to view by selecting a different compound in the Data Navigator or Compound List windows. 	<p>Alternate methods:</p> <ul style="list-style-type: none"> Click View > Compound Identification Results.

Exercise 11. Automated Oligonucleotide Sequence Matching

This exercise guides you through the setup of a worklist to automatically confirm the presence of 21mer_oligo sequence in a previously acquired sample.

Steps	Detailed Instructions	Comments
1 Open the <i>iii_oligo.m</i> method.	a Click Method > Open . b Select the <i>iii_oligo.m</i> folder. c Click Open .	This method was created in Exercise 10 (<i>iii</i> represents your initials).
2 Display the Worklist Automation > Worklist Actions section in the Method Editor.	Click Worklist Automation > Worklist Actions in the Method Explorer window.	
3 Select the appropriate worklist actions.	Select the following worklist actions in the Available actions list: <ul style="list-style-type: none"> • Find Compounds by Molecular Feature • Match Sequences • Generate Compound Report 	In this case we are using the default method parameters. For some data files, you will need to use different parameters as described in the <i>Quick Start Guide</i> or <i>online help</i> .
4 Add the selected actions to the Actions to be run list.	Click the  button. The items are added to the end of the list.	Actions will be executed in the order they appear in the list. You can reorder them using the Up and Down arrow buttons to the right of the list
5 Select Worklist as the match source.	a Click the Source tab in the Define and Match Sequences section of the Method Editor window. b Click Worklist .	This causes the software to get the sequence from the worklist rather than from the method.
6 Save the method.	Click Method > Save .	

Exercise 11. Automated Oligonucleotide Sequence Matching

Steps	Detailed Instructions	Comments
7 Create a worklist of one sample in the MassHunter Workstation Data Acquisition software.	<p>a Display the Worklist pane.</p> <p>b Click Worklist > Add Sample. A new sample row is added to the Worklist table.</p> <p>c Click Worklist > Add Column.</p> <p>d When the Add Column dialog box appears, select Protein as the Column Type.</p> <p>e Enter the Column name as 21mer_oligo.</p> <p>f Select 21mer_oligo.psq as the Value.</p> <p>g Click OK.</p>	If you plan to do batch data analysis using the worklist, consider using the DA Reprocessor tool that is installed with Agilent MassHunter Workstation Software - Data Acquisition program. In that case, skip Step 10 below.
8 Enter iii_oligo.m in the Method column of the worklist.	The <i>iii</i> represents your initials.	
9 Enter 11mer_oligo.d in the Data File column of the worklist.		
10 Set up to run the worklist for data analysis only.	<p>a Click Worklist > Worklist Run Parameters.</p> <p>b Select DA Only as Part of method to run.</p> <p>c Select the paths for the DA method and data file, then click OK.</p>	
11 Run the worklist.	Click Worklist > Run .	
12 Review the printed Compound reports.		

Exercise 12. Comparing Protein Digest Files

This exercise shows you how to compare compounds in two protein digest files.

Before you start

- Do “[Exercise 6. Interactive Protein Digest Sequence Matching](#)” on page 20 to get method and results files for this exercise.

Steps	Detailed Instructions	Comments
1 Open the data files to compare. In this exercise, we are using: <ul style="list-style-type: none"> • enolase-chip-final.d • enolase-oxidized-chip-final.d 	a Click File > Open Data File . b On the Open Data File dialog box, select the enolase-chip-final.d folder. c Select the following options: <ul style="list-style-type: none"> • Load results method • Load result data d Click Open . e Repeat <i>Steps a - d</i> to open enolase-oxidized-chip-final.d .	The data files used in this exercise have already been processed in Qualitative Analysis to find and identify compounds. These results are loaded when you open the data files.
2 View method parameters.	To view method parameters that were used to find and identify compounds, display the following sections of the Method Editor: <ul style="list-style-type: none"> • BioConfirm Workflow > Find by Molecular Feature • BioConfirm Workflow > Define and Match Sequences 	<ul style="list-style-type: none"> • Some of the method parameters are shown in Figure 11 on page 39, and Figure 12 and Figure 13 on page 40. • Note that 95% coverage is obtained for the Enolase digest sequence; see Figure 14 on page 40.
3 Start the compare protein digest files wizard.	Click Wizard > Compare Protein Digest Files: Compare Existing Results . The Select Reference File and Sample File(s) page is displayed.	For your own data files, you can find and identify compounds in either of the following ways: <ul style="list-style-type: none"> • Manually as described in “Exercise 6. Interactive Protein Digest Sequence Matching” on page 20, <i>or</i> • By clicking Wizard > Compare Protein Digest Files: Find Results, Identify and Compare to have the wizard guide you through the process.

Exercise 12. Comparing Protein Digest Files

Steps	Detailed Instructions	Comments
4 Select the files to use as Reference and Sample files.	<p>a From the list of opened files, select enolase-chip-final.d and click Select reference file.</p> <ul style="list-style-type: none"> The selected file name appears in the Reference file text box. The remaining file, enolase-oxidized-chip-final.d, is automatically moved to the list of Sample files. <p>See Figure 15 on page 41.</p> <p>b Click Next to open the next page of the wizard.</p>	<p>For your own data files:</p> <ul style="list-style-type: none"> To set the order that the sample files will be processed, use the Up and Down arrow buttons to the right of the list of samples. To open additional data files, click the Browse (...) button. To remove a selected Reference or Sample file, double-click on the file or click the Up arrow button. To remove <i>all</i> Sample files, click the double Up arrow button.
5 Set the compound correlation parameters on the Alignment Information page.	<p>a Set the retention time (RT window) and RT window tolerance values.</p> <p>b Set the Mass window and mass window tolerance values.</p> <p>See Figure 16 on page 41.</p> <p>c When finished setting the parameters, click Finish.</p> <p>A progress bar is displayed while the compounds are being compared.</p>	<p>This process may take several minutes to complete, depending on the number of samples and size of the data files.</p>
6 Review the results in the MassHunter Comparative Analysis program.	<p>Results are presented as follows in the Compound List window:</p> <ul style="list-style-type: none"> Column 1 (purple shading): correlated compound information Column 2 (red shading): reference file information Column 3 (blue shading): sample file information <p>See Figure 17 on page 42.</p>	

Exercise 12. Comparing Protein Digest Files

Method Editor: Find Compounds by Molecular Feature

Find Compounds by Molecular Feature | Method Items

Mass Filters | Mass Defect | Peak Filters (MS/MS) | Results

Extraction | Ion Species | Charge State | Compound Filters

Extraction algorithm

Target data type: Small molecules (chromatographic)

Input data range

☐ Restrict retention time to: minutes

☒ Restrict m/z to: 300.0000-1700.0000 m/z

Peak filters

☐ Use peaks with signal-to-noise >= 5.0 (Profile spectra only)

☒ Use peaks with height >= 500 counts (Profile and centroid spectra)

Allowed ion species

Positive ions

- ☒ +H
- ☒ +Na
- ☒ +K
- ☐ +NH4

Negative ions

- ☒ -H
- ☒ -Cl
- ☐ -Br
- ☐ -HCOO
- ☐ -CH3COO
- ☐ -CF3COO

Neutral losses

- ☐ H2O
- ☐ H3PO4

☐ Salt dominated positive ions (M+H may be weak or missing)

Figure 11 Method Parameters - Find Compounds by Molecular Feature, part 1

Method Editor: Find Compounds by Molecular Feature

Find Compounds by Molecular Feature | Method Items

Mass Filters | Mass Defect | Peak Filters (MS/MS) | Results

Extraction | Ion Species | Charge State | Compound Filters

Isotope grouping

Peak spacing tolerance: 0.0025 m/z, plus 7.0 ppm

Isotope model: Peptides

Charge state

☒ Limit assigned charge states to a maximum of: 15

☐ Treat ions with unassigned charge as singly-charged

Previous results

☒ Delete previous compounds

New results

☒ Highlight first compound

☐ Highlight all compounds

Chromatograms and spectra

☐ Extract MFE spectrum

☐ Extract raw spectrum

☐ Extract ECC

☐ Extract EIC

☐ Prefer profile for raw spectrum, if available

☐ Clip extracted raw spectrum

Asymmetric (m/z) - 5.0000 + 10.0000

☒ Extract MS/MS Spectrum

☐ Separate MS/MS spectrum per CE

☒ Average MS/MS spectrum for all CEs

Precursor tolerance: +/- 20.00 ppm

☒ Deisotope MS/MS spectrum

Display limits

☐ Display only the largest 100 compounds

Figure 12 Method Parameters - Find Compounds by Molecular Feature, part 2

Exercise 12. Comparing Protein Digest Files




Figure 13 Method Parameters - Define and Match Sequences




Figure 14 Sequence Coverage List Map for Enolase Digest Sequence

Exercise 12. Comparing Protein Digest Files




Figure 15 Reference and Sample files selected




Figure 16 Alignment Information Page with default parameters

Exercise 12. Comparing Protein Digest Files




Figure 17 MassHunter Qualitative Compare Program Window

www.agilent.com

In this book

This guide has familiarization exercises to demonstrate the use of Agilent MassHunter BioConfirm Software.

This guide is valid for the B.04.00 revision or higher of the G3298AA Agilent MassHunter BioConfirm Software, until superseded.

Microsoft® and Windows® are registered trademarks of Microsoft Corporation in the US and/or other countries.

© Agilent Technologies, Inc. 2011

Printed in USA
Revision A, February 2011



G6829-90007



Agilent Technologies