

# Efficient Analysis of Condensation Resins

## Agilent PLgel 5 $\mu$ m MIXED-D Columns

### Technical Overview

#### Introduction


The PLgel 5  $\mu$ m MIXED-D column is specifically designed for the analysis of polymers, paints and resin systems where material above 400,000 MW is unlikely to be present. High pore volume, combined with the 5  $\mu$ m efficiency, provides excellent resolution for low MW polymers and oligomers. Two, or even three, PLgel 5  $\mu$ m MIXED-D columns are the perfect replacement for the popular  $10^4/500$  Å or  $10^4/10^3/500$  Å/ $100$  Å column combinations.



**Agilent Technologies**

A variety of resins, produced by condensation polymerization, can be analyzed using three columns.

Column: 3 x PLgel 5  $\mu$ m MIXED-D, 300 x 7.5 mm  
(part number PL1110-6504)  
Eluent: THF  
Flow Rate: 1.0 mL/min  
Detection: RI



**Figure 1. Analysis of resins using PLgel 5  $\mu$ m MIXED-D columns**

[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 21, 2010

SI-01828



**Agilent Technologies**