



Princeton Applied Research

# WARRANTY

PerkinElmer Instruments, Inc. warrants each instrument of its own manufacture to be free of defects in material and workmanship. Obligations under this Warranty shall be limited to replacing, repairing or giving credit for the purchase price, at our option, of any instrument returned, shipment prepaid, to our Service Department for that purpose within ONE year of delivery to the original purchaser, provided prior authorization for such return has been given by an authorized representative of PerkinElmer Instruments, Inc.

This Warranty shall not apply to any instrument, which our inspection shall disclose to our satisfaction, to have become defective or unworkable due to abuse, mishandling, misuse, accident, alteration, negligence, improper installation, or other causes beyond our control. This Warranty shall not apply to any instrument or component not manufactured by PerkinElmer Instruments, Inc. When products manufactured by others are included in PerkinElmer Instruments, Inc. equipment, the original manufacturer's warranty is extended to PerkinElmer Instruments, Inc. customers.

PerkinElmer Instruments, Inc. reserves the right to make changes in design at any time without incurring any obligation to install same on units previously purchased.

THERE ARE NO WARRANTIES THAT EXTEND BEYOND THE DESCRIPTION ON THE FACE HEREOF. THIS WARRANTY IS IN LIEU OF, AND EXCLUDES ANY AND ALL OTHER WARRANTIES OR REPRESENTATIONS, EXPRESSED, IMPLIED OR STATUTORY, INCLUDING MERCHANTABILITY AND FITNESS, AS WELL AS ANY AND ALL OTHER OBLIGATIONS OR LIABILITIES OF PERKINELMER INSTRUMENTS, INC., INCLUDING, BUT NOT LIMITED TO, SPECIAL OR CONSEQUENTIAL DAMAGES. NO PERSON, FIRM OR CORPORATION IS AUTHORIZED TO ASSUME FOR PERKINELMER INSTRUMENTS, INC. ANY ADDITIONAL OBLIGATION OR LIABILITY NOT EXPRESSLY PROVIDED FOR HEREIN EXCEPT IN WRITING DULY EXECUTED BY AN OFFICER OF PERKINELMER INSTRUMENTS, INC.

## SHOULD YOUR EQUIPMENT REQUIRE SERVICE

A. Contact the Customer Service Department (865-482-4411) or your local representative to discuss the problem. In many cases it will be possible to expedite servicing by localizing the problem to a particular plug-in circuit board.

B. If it is necessary to send any equipment back for service, we need the following information.

- |  |  |
|--|--|
| 1. Model number and serial number.                     | 5. Your telephone number and extension.  |
| 2. Your name (instrument user).                        | 6. Symptoms (in detail, including control settings).   |
| 3. Your address.                                       | 7. Your purchase order number for repair charges (does not apply to repairs in warranty).                            |
| 4. Address to which the instrument should be returned. | 8. Shipping instructions (if you wish to authorize shipment by any method other than normal surface transportation.) |

C. U.S. CUSTOMERS — Ship the equipment being returned to:

PerkinElmer Instruments  
801 S. Illinois Avenue  
Oak Ridge, TN 37831  
ATTN: Customer Service

PHONE: 865-482-4411  
FAX: 865-483-0396

D. CUSTOMERS OUTSIDE OF U.S.A. — To avoid delay in customs clearance of equipment being returned, please contact the factory or the nearest factory distributor for complete shipping information.

---

All trademarks used herein are the property of their respective owners.

# TABLE OF CONTENTS

1. GENERAL .....	1
2. SYSTEM SPECIFICATIONS .....	3
3. INTERCONNECTIONS .....	5
3.1. Connecting the Potentiostat to the Power Booster .....	5
3.2. 263A Cell Cable Connections .....	7
3.3. KEPCO Front-Panel Setup .....	7
3.4. External Cell Connection and Initial Check .....	7



# 1. GENERAL

The Model 263A Power Booster is designed for special electrochemical applications such as battery research, etching, plating, electrosynthesis, and corrosion studies. Three models are available:

- The 263A/10A provides  $\pm 10$  amps (nominal<sup>1</sup>) at  $\pm 20$  volts (nominal) compliance
- The 263A/20A provides  $\pm 20$  amps (nominal) at  $\pm 20$  volts (nominal) compliance
- The 263A/8A provides  $\pm 8$  amps (nominal) at  $\pm 50$  volts (nominal) compliance

This add-on system consists of the 263A/94 2A option; a specially configured 20-10A, 20-20A, or 50-8A bipolar power supply from KEPCO<sup>®</sup>, Inc.; an internal PerkinElmer Instruments Power Booster Interface (see the block diagram in Fig. 1); and cabling.

The Power Booster interfaces the 263A and the external cell, acts as the cell switch, and provides the drive signal for the KEPCO power supply. Experiments can be controlled manually from the front panel of the 263A or from a PC using the appropriate PerkinElmer Instruments software.

The Power Booster Interface also provides cell current (scaled) required by the 263A I/E converter. With the 263A current range set to 10  $\mu$ A, 1  $\mu$ A reading on 263A = 1 A actual cell current; and 10  $\mu$ A reading on 263A = 10 A actual cell current.

The system frequency response has been modified to provide stability for the complete closed-loop response for electrochemical cells with a very high capacitance and very low resistance.

The Power Booster Interface operates in power-booster or normal mode — a simple cable connection and flip of a switch convert between modes. This manual discusses boosted mode only. For normal operation, see the potentiostat's user manual.

---

<sup>1</sup>Determined by the KEPCO.




Fig. 1 Block Diagram for the 263A Power Booster System.

Connections and settings are identical among the three Power Booster Interface options. However, for simplicity, we discuss and illustrate only the 10-A unit in this manual.

**CAUTION** *Do not use auto-ranging current in Power Booster mode!*  
 The 263A should be set to the 10- $\mu$ A current range on the front panel or via PC/software control.

## 2. SYSTEM SPECIFICATIONS

**Bandwidth** 5 KHz at -3 db with 1 V applied at counter with 1- $\Omega$  cell (1-W load)

**Voltage Rise Time** <70  $\mu$ s (10–90%, measured at Counter,  $\pm 10$  V/10-Hz square wave into a 1-W load, high-speed mode).

**Potential Noise** <500  $\mu$ V rms at E Monitor Out.

**Applied Potential Accuracy**  $\pm 0.5\%$  of reading,  $\pm 2$  mV (nominal; determined by the KEPCO).

**Current Measurement Rise Time** <75 $\mu$ s (10–90%, measured at I Monitor Out,  $\pm 10$  A/10-Hz square wave).

**Current Measurement Accuracy**  $\pm 1\%$  of reading,  $\pm 20$  mA.

**Rated Output Current** 20 A (nominal) for the 263A/20A system; 10 A for the 263A/10A system; and 8 A (nominal) for the 263A/8A system.

**Rated Output Voltage** 20 V at 20 A (nominal) for the 263A/20A and 263A/10A; 50 V at 8 A (nominal) for the 263A/8A.

**Overload Shutdown**  $\pm 20$  V or  $\pm 20$  A (nominal) for the 263A/20A;  $\pm 20$  V or  $\pm 10$  A (nominal) for the 263A/10A; and  $\pm 50$  V or  $\pm 8$  A (nominal) for the 263A/8A.<sup>2</sup>

---

<sup>2</sup>The KEPCO power supply allows you to limit the overload shutdown (that is, to create shutdown set points at lower than the maximum value. Refer to the KEPCO manual for more information.)





## 3. INTERCONNECTIONS

### 3.1. Connecting the Potentiostat to the Power Booster

The 263A and KEPCO are connected by 2 marked cables provided with the interface unit (see also the diagram in Fig. 2).

**CAUTION** Make all connections, as well as changes between Normal and Boost modes, with the power *off*.

**NOTE** Do not use the **C0345** cable (the label is on a white tab) in power-boosted mode. C0345 is for normal operation only (for instructions on using it, see the *263A User's Guide*).

**NOTE** In 4-terminal connections where the sense lead is used, the 263A must be in high-speed mode.

#### Boost Mode


On the 263A rear panel, flip the **BOOST/NORMAL** toggle switch to **BOOST**.

#### S-COM Connection

Connect the S-COM brown wire to **S COM** (or **COM**; they are shorted together) on the KEPCO rear panel. Connect the opposite end of the wire to screw-clamp connector **3** on the rear panel of the Power Booster Interface.

#### S-OUT Connection

Connect the S-OUT red wire to **S OUT** on the KEPCO rear panel. Connect the opposite end of the wire to screw-clamp connector **4** on the rear panel of the Power Booster Interface.



**Fig. 2** Connecting the 263A, KEPCO Power Supply, and Power Booster Interface.

### KEPCO Cardedge/Double BNC Cable

Connect the cardedge (female) connector to the **card edge (male)** projecting from the KEPCO rear panel (both connectors are keyed to ensure correct orientation).

On the double-BNC end of the cable, connect the **E/IN** connector to the rear panel of the 263A at the BNC connector marked **E IN**. Connect the **I/OUT** BNC portion of the cable to the rear panel of the potentiostat at the BNC connector marked **I OUT**.

## 3.2. 263A Cell Cable Connections

### CO378 Screw Clamp/Spring Clamp Cable

Connect the red and brown leads to the 263A rear panel at screw clamps **1** and **2**, respectively (see Fig. 2). The opposite end of this cable terminates in 2 spring-clamp connectors for the working (brown wire/green spring clamp) and counter (red wire/red spring clamp) leads.

### CO377 LEMO/Spring Clamp Cable (replaces CO345 during boost mode)


Connect this cable to the 263A front panel at the female LEMO connector. The opposite end of this cable terminates in white and green spring-clamp connectors for the reference and sense leads, respectively.

## 3.3. KEPCO Front-Panel Setup

- The **MODE** switch on the front panel of the KEPCO should be set to Voltage mode (left position).
- The **Voltage CONTROL** switch (red toggle on left side) should be in the OFF (down) position.
- The **Current CONTROL** switch (red toggle on right side) should be in the OFF (down) position.
- The **SENSE/COMMON** (black banana jacks at bottom-center) should be connected with a shorting bracket.
- The **OUTPUT/SENSE** (red banana jacks at bottom-center) should be connected with a shorting bracket.

## 3.4. External Cell Connection and Initial Check

Use the provided cell cables (Section 3.2) to connect the electrochemical cell being tested (see Fig. 3).



**Fig. 3 Connecting the Power Booster Interface to the 263A Differential Electrometer.**


**DANGER** Dangerous and possibly lethal voltages may be present on the electrometer's red lead. Never touch the red lead while mains power is supplied to the electrometer. For more information, see the electrometer manual.

Once the interconnections are complete, perform the following check to ensure the system is functioning properly:

1. Connect a 1- $\Omega$  dummy cell resistor as shown in Fig. 4.

**WARNING** The dummy cell resistor must have a power rating of at least 9 W. Three 3- $\Omega$ /3-W resistors linked in parallel are acceptable.

2. Connect the **Reference** and **Counter** to one side of the resistor (Reference connection closest to cell), and **Sense** and **Working** to the other side (Sense connection closest to the cell).
3. Turn on the 263A power.



**Fig. 4 Dummy Cell Connections.**

4. Turn on the KEPCO power supply.
5. From the front panel of the 263A, select the 10- $\mu$ A current range.
6. Set the 263A to high-speed mode.
7. Apply a  $\pm 3$ -V potential via the front panel (see the 263A *User's Guide*) and turn the cell on.
8. The current reading on the front-panel display should be  $\pm 3 \mu\text{A}$  ( $+3 \mu\text{A}$  if  $-3 \text{ V}$  are applied), and the E reading should be  $\pm 3 \text{ V}$ .
9. Turn the cell off.
10. Turn the 263A and KEPCO power off.
11. Disconnect the cell.

**Part No. 223119**  
**Manual Revision A**  
Printed in USA

**4352 1199**