# amoeba Dual Purpose Digital Microscope

44325

44326

### INSTRUCTION MANUAL MODEL #44325 / #44326

### INTRODUCTION

Congratulations on your Celestron microscope purchase. Your new Amoeba Dual Purpose Microscope is a precision optical instrument made of high quality materials to ensure durability and long life and is designed to give you a lifetime of enjoyment with minimal maintenance.

Amoeba can be used for virtually anything that traditional microscopes are used for: Hobbyists, Educators, Industrial Inspection, Teachers, Students, Science Applications, Law Enforcement Investigation and general use by consumers.

Amoeba is a dual purpose microscope – it works well as a stereoscopic (low power) microscope for solid objects and as a biological (high power) microscope for specimen slides.

Before attempting to use your Amoeba, please read through these instructions thoroughly to familiarize yourself with its functions and operations. Refer to the microscope diagram to locate parts discussed in this manual.

Amoeba provides powers of 10x, 60x and 200x (powers with a 14" monitor). It is well suited for examining solid objects such as coins, stamps, rocks, relics, insects, plants, skin, gems, circuit boards, various materials and many other objects. Amoeba is equally well suited to examine specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria and more at low, medium and high powers.

You can also explore the exciting microscopic world by creating your own specimen slides. An excellent book to read that describes various activities, projects and techniques is "The World of the Microscope" by Usborne publishing.

With the software included, you can observe magnified images, capture video or take snapshots with Microsoft Windows 7, XP, Vista and some 2000 operating systems. You can also use Amoeba with some other image capture software (Amcap, etc.) that you may have on your Windows operating system. With most Apple MAC Operating Systems, you can take videos and snapshots but will require image/photo capture software for MAC (e.g., combination of iChat with Photobooth, etc.).

**Note:** A CD/DVD drive is required and an open USB port (Version 2.0 interface works best although version 1.1 will work at a lower resolution).

MODEL #44325 AND #44326	SPECIFICATIONS
Power (Magnification)	10x, 60x and 200x using a 14" Monitor
USB Cable	2.0 power from a computer (1.1 is compatible, however video resolution will be lower)
Illuminator	Top and Bottom White LEDs
Digital Camera	1.3 MP CMOS Sensor – 1280x1024 pixel array
Plain Stage with Clips	82 mm (3.2") Diameter
Focuser	Coarse Focus with Dual Knobs
Weight /Dimensions	18 oz (510 g), LxWxH 5.3"x6.5"x10.0" (135x165x254 mm)

### **SPECIFICATIONS**

# GETTING TO KNOW YOUR AMOEBA MICROSCOPE

> Microscope Module	> Microscope Stand with Base
> USB Cable – 2.0	> CD-ROM - ULead Photo Management Software
> Instruction Manual	> Tweezers — Plastic
> Eyedropper – Plastic	> Needle Probe - Plastic
> 8 Pieces Prepared Slides	> 4 Pieces Blank Slides

### SETTING UP AMOEBA

- Carefully remove the microscope tube (module) and base with stand (silver part) from box and set on a table, desk or other flat surface.
- Carefully slide the microscope tube into the arms of the microscope base and stand until the wide section around the power ring rests against the arms.
- 3. The metal terminals at the back of the microscope tube should make contact with the terminals at the back of the arm (on base). These terminals supply power from the computer's USB port to the microscope base for the bottom illuminator.

1. Microscope Tube and Base

2. Microscope Base Arms

3. Terminals on the Back

### USING AMOEBA FOR VIEWING AND/OR IMAGING

Before using Amoeba for viewing or imaging on your computer screen, ensure the driver is recognized by your computer. Amoeba utilizes a "Plug and Play" system where the computer should automatically recognize the driver type on Microsoft Windows 7, XP, Vista and some 2000 operating systems as well as most Apple MAC versions 10.4.9 and later operating systems.

**USB Cable Socket** 

Plug the USB cable (larger socket end) into a USB port on your computer and plug the small socket of the USB cable into the microscope tube. Each computer system will have slightly different messages. For example, on some Windows operating systems, a message may read "Installing Device Driver Software" and shortly thereafter you may see a message "Your Devices Are Ready to Use." Follow the instructions shown on the screen.

You will need to install ULead Photo Management Software (on Windows operating systems) to actually view and/or image. All imaging should be done with the software program as the photo button on the microscope is not usable. Prior to installing the software, unplug the USB cable until the installation is complete. You can also use some different image/photo capture software that you may have on your Windows operating system (Amcap, etc.) if you desire. **Note:** To use Amoeba with most MAC operating systems, you will need MAC Version 10.4.9 or higher and must use a MAC compatible image/photo capture program (iChat with PhotoBooth, etc.). To find PhotoBooth on most Macs, go to your hard drive and click Applications. Choose PhotoBooth from the list of applications. If a camera is already attached to the computer, PhotoBooth will show up as your default. You will need to change this to "Camera" or "USB 2.0 Camera." You can now take photos with Amoeba using the capture button. Please note: Many keyboards on Macs have a USB port. Amoeba will not work through this USB port because it does not have enough power. You must connect the USB directly to your computer tower.

### INSTALLING ULEAD PHOTO MANAGEMENT SOFTWARE

Amoeba is packaged with a ULead Photo Management Software CD-ROM, which allows you to use your microscope for visual purposes and/ or for imaging with your computer when using Microsoft Windows operating systems.

Before installing the software, remove the microscope's USB cable from your computer. Do not connect the cable until the software installation is complete. Insert the CD into your CD-ROM drive. The software should automatically begin.

Open Windows Explorer and choose the drive which the CD is inserted and start the selection menu. Double click Run AutoRun.exe. or on some operating systems "Auto Play" opens automatically. If the software (for Windows operating systems) does not begin automatically. double click "My Computer" and right click the CD-ROM icon. The software should then start. Double click on the menu item "Install ULead VideoStudio 7" or you may see a list of files and select "Setup" - (sailboat icon) to begin installation.

- InstallShield Wizard will appear. Select a language from the following choices; Chinese (simplified or traditional),
- English, French, German, Italian, Japanese, Korean or Spanish. Next, select the default choices ("Next", etc.) until "Finish" is seen and click it.
- > Enter the software serial number found on the CD-ROM cover. Keep the cover in case you need to reinstall the software for any reason.

Messages will appear (each operating system may show a different message) to guide you through the installation process.

**Note:** During installation you may see a pop up screen about installing "Real Player" or "Quick Time" and just click "Cancel". From the initial set up screen files, you can also choose to download "Install USB 2.0 PC Camera" which is AmCap for viewing or imaging.

After the software is properly installed, you need to restart your computer. You will have a desktop icon called "ULead Video Studio". If this icon does not appear, you need to select the program from your "Installed Programs" list when you use your Amoeba.

If you want to upgrade to a more sophisticated ULead software program, please visit www.ulead.com.

### AMOEBA OPERATION

PLUG IN THE USB CABLE TO THE COMPUTER AND AMOEBA TO BEGIN VIEWING OBJECTS AND SLIDES – You can view or image coins, stamps, currency, gems, plants, foods, circuit boards, specimen slides and many other objects. Read the information for focusing, power and illumination (see below) before proceeding.

**Note:** For solid objects you can view up to 10 mm thickness at all powers or up to 15 mm at 10x. However, for large solid objects, read the instructions (#5 below).

#### FOCUSING AND CHANGING POWER

(MAGNIFICATION) – When Amoeba is placed directly over the object or slide, use the focus knob to obtain a sharp focus. Rotate past the focus and back and you will see how to obtain the best focus position. With practice, it will become easier to find.

Slide Placed on the Stage

60x Power

- 1. Always start with lowest power (10x) by turning the power ring to the 10x position.
- For viewing a specimen slide, position the slide under the stage clips as close to the centerof the stage as possible. For viewing a solid object, position the object as close as possible to the center of the stage.
- 3. To start the ULead Photo Management Software program, a green LED will light up to indicate power is ON. Your specimen slide or object should appear in the computer screen's preview area. If not, click on Capture/Capture Settings. Turn the focus knob in either direction until you get a sharp image. You may have to move the specimen slide or object to center in the preview screen. Centering objects will become easier with practice.

4. Now you can choose a higher power (60x) or the maximum power (200x) by rotating the power ring until it clicks into position. You will need to refocus your slide or object after changing power. At high power, your views will be greatly magnified but somewhat darker and sometimes more challenging to find. The most enjoyable views are seen at lower powers which have a wider field of view and display brighter images, especially with solid objects.

5. Viewing or imaging large solid objects (over 10 mm thick or those that will not fit on the stage area) is possible using the microscope tube by itself (not connected to the base and stand). With the USB cable attached, remove the microscope tube from the base and stand. The top illuminator is used for this type of viewing or imaging.

Low power (10x) can be used directly on the object and up to 10 mm away. Medium power (60x) can be used on objects 5 mm to 10 mm away. High power (200x) cannot be used for this type of viewing or imaging. Remember to refocus each time you move the microscope.

NOTE: For imaging or during set up when using some operating systems; you may see windows or ULead error message such as "cannot start the capture graph", "cannot pre-allocate capture file space", "window color format has changed",etc. you should disregard any of these messages and click ok.

#### IMAGING

On some operating systems you may see error messages pop up when trying to image and if so click ok or ignore them. On some operating systems, the LED lights may turn off when a video is started but immediately turn the LED's back on as the video will be recording.

#### ILLUMINATION

Illumination is powered by the USB cable plugged into your computer and microscope tube. The USB cable also provides power for viewing and imaging using ULead Photo Management Software.

- 1. The top illuminator (light shining down onto solid objects) turns on when pressing the light button once.
- The bottom illuminator (light shining up from the base through specimen slides) turns on when pressing the light button twice.
- Both illuminators (top and bottom) can be turned on by pressing the light button (again) for viewing semi-transparent objects or slides.
- 4. Press the light button (again) to turn off all illumination.

Note: When power is on, a white LED will light up inside the top of the microscope tube which has no function other than to provide a visual effect.

### CARE, MAINTENANCE AND WARRANTY

Your Amoeba is a precision optical instrument and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- > Store Amoeba in a clean and dry place
- > Be very careful if using your Amoeba in direct sunlight to prevent damage to the microscope or your eyes
- > Never point Amoeba towards the Sun or its camera can be damaged and cease to operate
- > Clean the outside surfaces with a moist cloth
- > Never disassemble or clean internal parts of the microscope. This should be done by qualified technicians at the factory or other authorized repair facilities
- > Do not disassemble or take apart the

microscope or damage can occur

- > When handling glass specimen slides, use care as the edges can be sharp
- > Remove the USB cable from the computer's USB port when not in use

#### WARRANTY

Your Amoeba has a two year limited warranty. Please visit the Celestron website for detailed information on all Celestron microscopes at www.celestron.com.

#### FCC Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- · Reorient or relocate the receiving antenna.
- · Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

2835 Columbia St. Torrance, California 90503 U.S.A. ©2012 Celestron

All rights reserved. / Printed in China 04-12 Designed and intended for those 13 years of age and older

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