# **PAJ Series**

**NTEP / MC  
Certified Models**

**Instruction Manual  
Manuel d'instruction**

1. INTRODUCTION

1.1 Safety Precautions

Please follow these safety precautions:

- Verify that the AC Adapter input voltage matches the local AC power supply.
- Use the balance only in dry locations.
- Do not operate the balance in hostile environments.
- Do not drop loads on the platform.
- Do not place the balance upside down on the platform or platform mounting cone.
- Service should be performed only by authorized personnel.

1.2 Controls

| Button:                                      | Functions:          |                                        |
|----------------------------------------------|---------------------|----------------------------------------|
| <b>Zero - On</b><br><b>Off</b><br><b>Yes</b> | Short Press:        | Turns balance on, sets display to zero |
|                                              | Long Press:         | Turns balance off                      |
|                                              | Short press (Menu): | Selects or accepts setting             |
| <b>Unit</b><br><b>Menu</b><br><b>No</b>      | Short Press:        | Steps through active units and modes   |
|                                              | Long Press:         | Enters Menu                            |
|                                              | Short press (Menu): | Steps through available settings       |
| <b>Tare</b>                                  | Long press (Menu):  | Exit menu or abort out of menu item    |
|                                              | Short Press:        | Performs Tare function                 |
|                                              | Long Press:         | no function                            |
| <b>Print</b><br><b>Cal</b>                   | Short Press:        | Sends data                             |
|                                              | Long Press:         | Initiates Span Calibration             |

2. INSTALLATION

2.1 Package Contents

Gold Models


- Instruction Manual
- Power Adapter
- Balance
- Pan
- Pan Support
- Wind-Ring
- Warranty Card

Carat Models

- Instruction Manual
- Power Adapter
- Balance
- Pan
- Glass Doors and Panels
- Gem Scoop (Carat models only)
- Warranty Card

## 2.2 Install Components

### Gold Models

1) Install Wind-Ring

2) Install Pan Support


3) Install Pan

### Carat Models

1) Install Side Doors - Insert fully into Top Frame then down over retainer.

2) Install Front and Rear Panels - Insert bottom edge in groove then press until locked.


3) Install Top Door


4) Install Pan

## 2.3 Level Balance

Level the balance on a firm, steady surface. Avoid locations with excessive air current, vibrations, heat sources or rapid temperature changes.

## 2.4 Connect Power

For use with CSA Certified (or equivalent approved) power source, which must have a limited circuit output.

## 2.5 Initial Calibration

**Without InCal** - Power on the balance by pressing **0/T**. Press and hold **Print/Cal** until [**CAL**] is displayed. The display flashes the calibration mass needed. To select the alternate calibration weight press **No**. Put the calibration mass on the pan. The display flashes [**busy**], then [**CLEAR PRN**]. Remove the mass. When calibration is complete, [**done**] is displayed.

**InCal** - Press and hold **Print/Cal** until [**CAL**] is displayed. The display flashes [**busy**], when calibration is complete, [**done**] is displayed.

Note: Calibrations should be performed after a warm up time of 60 minutes.

## 3. OPERATION

Count, Percentage or specific units of measure must be activated in the MODE or UNIT menus if they are not initially available.

### 3.1 Weighing Mode

Repeatedly press **Unit** until the desired unit icon is displayed.

Press **Zero** to zero the balance and then place objects to be weighed on the pan.

Press **Tare** to zero the weight of a weighing container.

### 3.2 Count Mode

- Use the Count mode to count parts of uniform weight.

To access Count Mode, press **Unit** until the display shows [**Count**].

Establish an Average Piece Weight (APW) - Each time a new type of part is to be counted, the nominal weight of one piece (APW) must be established using a small quantity of pieces.

With [**CLR APW**] displayed, press **No** to use the previously saved APW, or press **Yes** to establish a new APW.

The display indicates the number of pieces to be used to establish the new APW. If a different sample size is preferred, press **No** until the desired sample size is displayed (10, 20, 50 or 100). Put the specified number of pieces on the pan. Press **Yes** to accept new APW or **No** to abort.

Count - Place the quantity to be counted on the pan.


APW Optimization - Since the weight of each piece varies slightly, APW Optimization may be used to increase

the accuracy of the count. The balance automatically recalculates the Average Piece Weight when the number of parts on the pan is less than three times the original sample size. The display shows [**APW.OP1**] each time the APW is optimized.


### 3.4 Weigh Below Feature

Remove Weigh Below Cover


Attach wire or string to Hook

Suspend sample

## 4. SETTINGS

### 4.1 Menu Navigation

**ENTER MENU** - When the balance is on, press and hold **Unit/Menu** until [**ENTER**] appears. Release the button and the Calibrate [**.CAL.**] menu will display.

**MENU NAVIGATION** - Select menus, menu items and settings through use of the **Yes** and **No** buttons. Solid arrows point to the content displayed when **Yes** is pressed, Dashed lines when **No** is pressed.

**CHANGE SETTINGS** - To select the displayed Setting, press **Yes**. To move to the next Setting, press **No**.

**EXIT MENU** - When [**.END.**] is displayed, press **Yes** to exit the menu function, or press **No** to return to the Cal menu. Note: Press and hold **No** at any time to exit quickly.

## 4.2 Calibration Menu [**.CAL.**]

InCal or Span calibration should be performed daily and when the room temperature changes.

InCal [**inCAL**] calibrates the balance using an internal mass. (InCal models only)

Span Calibration [**SPAN**] uses two weight values: zero and a weight between 50% and 100% of the capacity of the balance.

Linearity Calibration [**L IN**] uses three weight values: zero, 50% of the capacity and full capacity. Generally this calibration is not required unless testing shows that the linearity error exceeds the linearity tolerance in the specification table. (Not available in InCal models.)

Calibration Adjust [**CAL Adj**] allows adjustment to the result of the internal calibration by  $\pm 99$  divisions. (InCal models only)

## 4.3 Setup Menu [**.SETUP.**]

Filter [**FILTER**] - Use the low setting (SET LO) when environmental disturbances are not present. Use the medium setting (SET MED) for normal environments. Use the high setting (SET HI) when vibrations or air currents are present.

Auto Tare [**A-TARE**] - The initial item placed on the balance is assumed to be a container so it is zeroed out. The next item is then weighed. When the pan is cleared the balance resets, waiting for a container. (SET OFF, SET ON). Local weights and measures authorities may require that Auto Tare be disabled.

Global Reset [**GBL RST**] - Resets all settings to factory default values. (RESET)

## 4.4 Mode Menu [**.MODE.**]

Count Mode [**COUNT**] - [SET ON, SET OFF]

Average Piece Weight (APW) Optimization [**APW OPT**] - (SET ON, SET OFF)

## 4.5 Unit Menu [**.UNIT.**]

The Unit menu is used to enable or disable a specific unit. (SET ON, SET OFF) The unit is indicated by a small character next to Unit in the display (g = grams).

## 4.6 Print Menu [**.PRINT.**]

Auto Print [**APRINT**] - Data will be continuously sent when [**CONTINU**] is set. Interval [**INTER**] sends data every 1 to 3600 seconds. When Stable [**STABLE**] will send data when the balance detects a new stable reading. This can be a weight value only [**LOAD**] or it can also include a stable zero [**ZERO**]. [**OFF**] disables automatic printing.

Content [**CONTENT**] - The content in the data transmission can be modified. Each of the following settings can be set on or off. Number Only [**NUMBERS**] will only send the numeric result. Balance ID [**BAL. ID**] will add the Balance serial number for traceability purposes. Reference [**REFERENCE**] will add reference information relevant to the current mode. GLP [**GLP**] will send additional items to allow proper documentation of laboratory results.

Line Format [**LINE Fmt**] - Single line format [**SINGLE**] will put all the data in one line separating each output with a comma (.). Multi line format [**MULTI**] will put each data output on a new line. Multi +4 [**MULTI+4**] will add 4 line spaces between each output.

```

^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^
User ID: .....
Bal ID: 123456789
Proj ID: .....
Time: .....:.....:.....
Date: ...../...../.....
120.01 g
^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^^

```

```

_____ GLP (ON)
_____ Balance ID (ON)
_____ GLP (ON)
_____ GLP (ON)
_____ GLP (ON)
_____ Result
_____ Multi Line Format with 4 line feeds
_____ Line Feed-2
_____ Line Feed-3
_____ Line Feed-4

```

## 4.7 RS232 Menu [RS232.]

**Baud** [BAUD] - The RS232 baud rate can be set to 600, 1200, 2400, 4800, 9600 and 19200.

**Parity** [PARITY] - Parity can be sent to 7 bits-even parity [7 EVEN], 7bits-odd parity [7 odd], 7bits-no parity [7 No] or 8bit-no parity [8 No].

**Handshake** [HANDSH] - Handshake can be set to off [SET OFF], X on-X off [on-off], or hardware [HANDSH].

## 4.8 Lock Menu [LOC.]

When a Lock Menu item is set on the indicated menu cannot be changed.

[LOC CAL] - Calibration menu lock

[LOC SET] - Setup menu lock

[LOC Mod] - Mode menu lock

[LOC Unit] - Unit menu lock

[LOC Print] - Print menu lock

[LOC 232] - RS232 menu lock

## 5.0 LEGAL FOR TRADE

When the balance is used in trade or legally controlled applications it must be set up, verified and sealed in accordance with local weights and measures regulations. It is the responsibility of the purchaser to ensure that all pertinent legal requirements are met.

### 5.1 Settings

Before the verification and sealing perform the following steps:

1. Verify the menu settings meet the local weights and measures regulations.
2. Perform a calibration as explained in Section 2.5.
3. Set the position of the security switch as shown in Section 5.3.

**Note:** When the security switch is set on the following menu settings cannot be changed:

CAL, SETUP, MODE, UNIT and LOC.

### 5.2 Verification

The local weights and measures official or authorized service agent must perform the verification procedure.

### 5.3 Sealing Access to the Balance Settings

The local weights and measures official or authorized service agent must apply a security seal to prevent tampering with the settings. Refer to the illustrations below for the sealing methods.

Un-Locked

Locked with Paper Seal

Locked with Wire Seal

## 6.0 MAINTENANCE

### 6.1 Troubleshooting

| Symptom            | Possible Cause                                                              | Remedy                                                                  |
|--------------------|-----------------------------------------------------------------------------|-------------------------------------------------------------------------|
| Cannot turn on     | No power to balance                                                         | Verify connections and voltage                                          |
| Poor accuracy      | Improper calibration<br>Unstable environment                                | Perform calibration<br>Move balance to suitable location                |
| Cannot calibrate   | Unstable environment<br>Incorrect calibration masses                        | Move the balance to suitable location<br>Use correct calibration masses |
| Cannot access mode | Mode not enabled                                                            | Enter menu and enable mode                                              |
| Cannot access unit | Units not enabled                                                           | Enter menu and enable units                                             |
| Err 5              | Average Piece Weight too small                                              | Add additional samples                                                  |
| Err 7.0            | Time out                                                                    |                                                                         |
| Err 8.1            | Pan has load during power on                                                | Remove weight from pan and re-zero                                      |
| Err 8.2            | Pan was removed prior to power on                                           | Install pan and re-zero                                                 |
| Err 8.3            | Weight on pan exceeds capacity                                              | Remove weight from the pan                                              |
| Err 8.4            | Pan was removed during weighing                                             | Re-install pan                                                          |
| Err 9.5            | Factory calibration data corrupted                                          | Contact the authorized dealer                                           |
| Err 9.8            | Factory calibration data corrupted                                          | Perform calibration                                                     |
| Error 53           | EEPROM Checksum error                                                       | Contact the authorized dealer                                           |
| REF Err            | Reference Weight is too small                                               | Add additional samples                                                  |
| LOWrEF             | Reference Weight is too low for accurate parts counting or percent weighing | Add additional samples or continue to weigh with less accurate results  |

### 6.2 Service Information

If the troubleshooting section does not resolve or describe your problem, contact your authorized Ohaus service agent. Please visit our web site, [www.ohaus.com](http://www.ohaus.com) to locate the Ohaus office nearest you.

### 6.3 Accessories

|                          |                   |
|--------------------------|-------------------|
| Security device          | 76288-01          |
| Auxiliary Display        | PAD7              |
| Printer – Thermal        | STP-103           |
| Printer - Impact         | CBM-910 or SRP275 |
| Cable Kit - STP-103      | 80252581          |
| Cable Kit - CBM-910      | 80252571          |
| Data collection Software | SW12              |
| Ohaus Collect            | 80500746          |

## 7. TECHNICAL DATA

Ambient conditions - The technical data is valid under the following ambient conditions:

Ambient temperature: 10°C to 30°C

Relative humidity: 15 % to 80 % at 31°C non-condensing, decreasing linearly to 50% at 40°C

Height above sea level: Up to 2000 m



Operability is assured at ambient temperatures between 5°C and 40°C

Power - AC adapter, Balance power input 8-14.5 VAC, 50/60Hz 4VA or 8-20 VDC, 4W



Protection - Protected against dust and water, Pollution degree: 2, Installation category: Class II

### 7.1 Drawings

#### Gold Models

#### Carat Models

### 7.2 Specifications

#### External Calibration Models

| MODEL           | PAJ2003N                   | PAJ812N                                               | PAJ2102N       | PAJ2101N | PAJ3102N       | PAJ3101N | PAJ4102N       | PAJ4101N |
|-----------------|----------------------------|-------------------------------------------------------|----------------|----------|----------------|----------|----------------|----------|
| Max             | 2000ct                     | 810g                                                  | 2100g          |          | 3100g          |          | 4100g          |          |
| Min             | 1ct                        | 5g                                                    |                |          |                |          |                |          |
| d=              | 0.005ct                    | 0.01g                                                 |                | 0.1g     | 0.01g          | 0.1g     | 0.01g          | 0.1g     |
| e=              | 0.05ct                     | 0.1g                                                  |                |          |                |          |                |          |
| Approval        | Class II                   |                                                       |                |          |                |          |                |          |
| Tare Range      | To capacity by subtraction |                                                       |                |          |                |          |                |          |
| Stabilization   | 3 Seconds                  |                                                       |                |          |                |          |                |          |
| Span cal mass   | 200g or 400g               | 500g or 800g                                          | 1000g or 2000g |          | 2000g or 3000g |          | 2000g or 4000g |          |
| Lin. cal masses | 200g, 400g                 | 500g, 800g                                            | 1000g, 2000g   |          | 2000g, 3000g   |          | 2000g, 4000g   |          |
| Pan Size        | Ø4.7in/Ø12cm               | Ø7.1in / Ø18cm (PAJ models 6.1 x 7.1in / 16.8 x 18cm) |                |          |                |          |                |          |
| Net Weight      | 10.1 lb / 4.6kg            | 10 lb / 4.5kg                                         |                |          |                |          |                |          |

InCal Models

| MODEL           | PAJ603CN                   | PAJ1003CN    |
|-----------------|----------------------------|--------------|
| Max             | 600ct                      | 1000ct       |
| Min             | 0.2ct                      |              |
| d=              | 0.001ct                    |              |
| e=              | 0.01ct                     |              |
| Approval        | Class II                   |              |
| Tare Range      | To capacity by subtraction |              |
| Stabilization   | 3 Seconds                  |              |
| Span cal mass   | 50g or 100g                | 100g or 200g |
| Lin. cal masses | 50g, 100g                  | 200g, 500g   |
| Pan dia.        | Ø3.5in / Ø9cm              |              |
| Net Weight      | 10.1 lb / 4.6kg            |              |

7.3 Communication


7.3.1 Commands

The RS232 Interface allows a computer to control the balance as well as receiving data such as the displayed weight. The balance will return "ES" for invalid commands.

| Command | Function                                                                      |
|---------|-------------------------------------------------------------------------------|
| IP      | Immediate Print of displayed weight.                                          |
| P       | Print displayed weight (uses stable on/off menu setting).                     |
| CP      | Continuous Print.                                                             |
| xP      | Interval Print x = Print Interval (1 -3600 sec)                               |
| T       | Same as pressing Zero Key.                                                    |
| ON      | Turns balance ON.                                                             |
| OFF     | Turns balance OFF.                                                            |
| PSN     | Show Serial Number.                                                           |
| PV      | Version: Print product name, software revision and LFT ON (if LFT is set ON). |
| PU      | Print current mode/unit                                                       |
| x#      | Set PC ref wt (x) in grams                                                    |
| P#      | Print PC ref wt                                                               |
| x%      | Set % ref wt (x) in grams                                                     |
| P%      | Print % ref wt                                                                |

7.3.2 RS232 (DB9) Pin Connections

- Pin 2: Balance transmit line (Tx/D)
- Pin 3: Balance receive line (Rx/D)
- Pin 5: Ground signal (GND)
- Pin 7: Clear to send (hardware handshake) (CTS)
- Pin 8: Request to send (hardware handshake) (RTS)

## 7.4 Compliance

| Compliance to the following standards is indicated by the corresponding mark on the product. |                                                    |
|----------------------------------------------------------------------------------------------|----------------------------------------------------|
| Mark                                                                                         | Standard                                           |
|             | CAN/CSA-C22.2 No. 61010-1-04; UL Std. No. 61010A-1 |

### FCC Note

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

### Industry Canada Note

This Class A digital apparatus complies with Canadian ICES-003.

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### ISO 9001 Registration

In 1994, Ohaus Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritas Quality International (BVQI), confirming that the Ohaus quality management system is compliant with the ISO 9001 standard's requirements. On May 15, 2003, Ohaus Corporation, USA, was re-registered to the ISO 9001:2000 standard.

## LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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PN 80252972 © Ohaus Corporation 2009, all rights reserved / tous droits réservés.