


# Phenolmager™ HT

## USER GUIDE



An Automated  
Quantitative Pathology  
Imaging System



## **Notice**

This manual is published by Akoya Biosciences, Inc. The information in this document is subject to change without notice and should not be construed as a commitment by Akoya Biosciences, Inc. Akoya assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Akoya Biosciences be liable for incidental or consequential damages in connection with or arising from the use of this manual.

This manual describes how to use Phenolmager HT (formerly Vectra Polaris) software.

For more information contact:

Akoya Biosciences, Inc.  
100 Campus Drive, 6th Floor, Marlborough, MA, 01752, USA

Phone: +1 855-896-8401

Fax: +1 855-404-0061

Email: [CustomerCare@akoyabio.com](mailto:CustomerCare@akoyabio.com)

Email (outside the US): [CustomerCareEMEA@akoyabio.com](mailto:CustomerCareEMEA@akoyabio.com)

Web site: <http://www.akoyabio.com>

This software covered by US Patent 7,555,155; 7,655,898; and patents pending.

---

# Table of Contents

<b>Chapter 1 Preface</b>	<b>5</b>
<b>Chapter 2 Instrument Safety</b>	<b>11</b>
1 Required Training	11
2 Electrical Safety	11
Power Cord Selection	12
Fuses	12
Cables and Adapters	13
3 Mechanical Safety	13
4 Weight Warning	13
5 Bright Light	13
<b>Chapter 3 Introduction</b>	<b>15</b>
1 Intended Use	15
2 Principles of Operation	15
Example Applications	17
3 Theory of Imaging	17
Light	18
Human Perception of Light Intensity and Color	18
Light Absorbance and Reflection	19
Fluorescence	19
Photobleaching	20
Filter Sets for Conventional vs. Multispectral Imaging	21
Multispectral Analysis	22
4 Whole Slide Scanning	24
5 Multispectral Field Imaging	25
Configuring the Multispectral Fields	26
Processing the Multispectral Fields	26
<b>Chapter 4 Hardware Reference</b>	<b>27</b>
1 Front View	27
2 Top View	29
3 Right-Side View	30
4 Left-Side Connectors	31
5 Slide Carrier Hotel	32
6 Slide Carrier	33
7 Barcode Reader	35
8 System Computer and Monitor	35

9 Specifications.....	36
<b>Chapter 5 Operation.....</b>	<b>38</b>
1 System Startup.....	38
Turn on the Phenolmager HT Instrument.....	38
Launch the Phenolmager HT Software.....	39
2 System Shutdown.....	39
3 Inspecting Slides and Slide Carriers.....	39
4 Loading Slides into the Slide Carriers.....	40
5 Loading Slide Carriers into the Slide Carrier Hotel.....	41
6 Removing Slide Carriers from the Slide Carrier Hotel.....	43
7 Removing Slides from the Slide Carriers.....	43
<b>Chapter 6 Software Operation.....</b>	<b>44</b>
1 Software Overview.....	44
Gear Menu .....	46
2 System Dashboard.....	47
Disk Space .....	48
Brightfield References .....	48
Fluorescence References .....	50
Compensation Information .....	52
3 Creating and Editing Protocols.....	53
Studies .....	53
Creating Protocols .....	54
Editing Protocols .....	56
Setting Exposures .....	69
4 Scanning Slides.....	75
The Carrier .....	76
Carrier Status .....	76
Slide Status .....	78
Setting Up Scan Rules .....	78
<b>Chapter 7 Maintenance.....</b>	<b>87</b>
1 Cleaning the Instrument Exterior.....	87
2 Cleaning the Monitor.....	87
3 Cleaning the Power and Communication Ports.....	88
4 Cleaning the Slide Carriers.....	88
5 Replacing the Fuses.....	88
<b>Chapter 8 Appendix A: Akoya TIFF Specification.....</b>	<b>89</b>
<b>Chapter 9 Appendix B: EULA.....</b>	<b>96</b>
<b>Index .....</b>	<b>112</b>



# 1 Preface

## Copyright

---

This manual is published by Akoya Biosciences, Inc., 100 Campus Drive, 6th Floor, Marlborough, MA 01752 USA. Copyright 2020, Akoya Biosciences, Inc. All rights reserved, including but not limited to those rights to reproduce this publication or parts thereof. Reproduction of this publication or parts thereof, or the products it describes, by any means or in any form is expressly prohibited without the written permission of Akoya Biosciences.

## Content

---

Any errors or omission which may have occurred in this publication despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection. Akoya provides this publication “As Is” without warranty of any kind, either express or implied, including but not limited to the implied warranties of merchantability or fitness for a particular purpose. Some states or jurisdictions do not allow disclaimer of express or implied warranties in certain transactions; therefore, this statement may not apply to you. Akoya reserves the right to revise this publication and to make changes from time to time in the content hereof without obligation of Akoya to notify any person of such revision or changes. Further, Akoya may make modifications to the product described in this manual at any time without any obligation to notify any person of such modifications.

## Proper Equipment Operation

---

### WARNINGS:



- To reduce the risk of electric shock, do not remove the instrument panels. No user serviceable parts are inside. Refer to qualified service personnel if help is required.
- Use this product only in the manner described in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

### Avertissements:



- Pour réduire le risque de choc électrique, ne pas retirer le couvercle. Ce produit ne contient aucune pièce pouvant être réparée par l'utilisateur. Au besoin, confier l'appareil à un réparateur qualifié.
- Ce produit ne doit être utilisé que comme décrit dans ce manuel. Si cet appareil est utilisé d'une manière autre que celle spécifiée par le fabricant, la protection fournie par l'appareil peut être entravée.

## Contact Us

---

If you have a question about a product that is not answered in this manual or online help, or if you need assistance regarding this product, or for information about training for this product, please contact Akoya Customer Care:

Phone: +1 855-896-8401

Fax: +1 855-404-0061

Email: [CustomerCare@akoyabio.com](mailto:CustomerCare@akoyabio.com)

Email (outside the US): [CustomerCareEMEA@akoyabio.com](mailto:CustomerCareEMEA@akoyabio.com)

Website: [www.akoyabio.com](http://www.akoyabio.com)

Before you call, have the following information available for the technical representative:

- Product serial number
- Software version (found by choosing About from the main Help menu)

## Product Service and Customer Support Plans

---

Akoya offers a full range of services to ensure your success. From our original factory warranty through a comprehensive line of customer support plans, Akoya offers you Field Service Engineers and in-house Specialists who are dedicated to supporting your hardware, software, and application development needs.

Our programs can include such useful services as:

- Preventive maintenance
- Diagnostic servicing performed on-site by Akoya field service engineers or remotely via Technical Support
- Validation performed on-site by Akoya field service engineers
- Extended use of the Akoya Technical Support Center
- Software updates
- Parts, labor, and travel expense coverage
- Other customized services upon request

## CE

---

This device complies with all CE rules and requirements.



**NOTE:** Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**REMARQUE:** Tout changement ou modification apporté à cet instrument non expressément approuvé par l'entité responsable de la conformité peut annuler l'autorisation d'opérer l'appareil accordée à l'utilisateur.

## FCC

### RADIO AND TELEVISION INTERFERENCE

**NOTE:** This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

In order to maintain compliance with FCC regulations shielded cables must be used with this equipment. Operation with non-approved equipment or unshielded cables is likely to result in interference to radio & television reception.

## KC

This device complies with MSIP (Ministry of Science, ICT, and Future Planning) EMC Registration requirements. This instrument is registered as a Class A instrument for business use only. Product seller and user should notice that this equipment is not for household use.




### A급 기기 (업무용 정보통신기기)











이 기기는 업무용으로 전자파적합등록을 한 기기이오니  
판매자 또는 사용자는 이 점을 주의하시기 바라며, 만약  
잘못판매 또는 구입하였을 때에는 가정용으로 교환하시기  
바랍니다.










## Table of Symbols

[Table 1](#)<sup>7)</sup> contains symbols that identify particularly important information and alert you to the presence of hazards. These symbols may appear in this manual and/or on the product it describes.

Table 1. Important Symbols

Symbol Symbole	Description
	<b>DANGER:</b> An imminently hazardous situation, which, if not avoided, will result in death or serious injury.  <b>DANGER:</b> Situation présentant un danger imminent qui, s'il n'est pas éliminé, peut entraîner des blessures graves, voire la mort.
	<b>WARNING:</b> Caution. Refer to the User's documentation. (ISO 7000-0434B)  <b>AVERTISSEMENT:</b> Attention. Se reporter à la documentation de l'utilisateur.
	<b>NOTE:</b> A cautionary statement; an operating tip or maintenance suggestion; may result in instrument

	<p>damage if not followed.</p> <p><b>REMARQUE:</b> Énoncé indiquant une précaution à prendre, un conseil de fonctionnement ou une suggestion d'entretien; son non-respect peut provoquer des dommages à l'instrument.</p>
	<p>Hazardous voltage; risk of electric shock. (IEC 60417-6042)</p> <p>Tension dangereuse; risque de blessure par électrocution.</p>
	<p>Crush hazard. Risk of body parts, hair, jewelry, or clothing getting caught in a moving part. (ISO 3864)</p> <p>Danger d'écrasement. Faire attention que les parties corporelles, les cheveux, les bijoux ou les vêtements ne soient pas pris dans une pièce mobile.</p>
	<p>Risk of fire. (ISO 3864)</p> <p>Risque d'incendie.</p>
	<p>Risk of explosion. (ISO 3864)</p> <p>Risque d'explosion.</p>
	<p>Lifting hazard. May result in injury. (ISO 3864)</p> <p>Levage dangereux. Peut entraîner des blessures.</p>
	<p>Protective ground symbol. (IEC 60417-5019)</p> <p>Symbole de terre de protection.</p>
	<p>Ground symbol. (IEC 60417-5017)</p> <p>Symbole de terre.</p>
	<p>Fuse. (IEC 60417-5016)</p> <p>Fusible.</p>
	<p>Alternating current. (IEC 60417-5032)</p> <p>Courant alternatif.</p>
	<p>On (power). (IEC 60417-5007)</p> <p>Marche (alimentation).</p>

	Off (power). (IEC 60417-5008) Arrêt (alimentation).
	CE compliance mark. Marque de conformité CE.
	Serial Number (ISO 7000-2498) Numéro de série.
	Date of Manufacture. (ISO 7000-2497) Date de fabrication.
	WEEE symbol (EN50419:2005). Do not dispose of as unsorted municipal waste.
	Unified Korea Certification Mark (KC Mark)
	Signifies that the unit has passed safety tests for grounding, power line transience, and current leakage.  Signifie que l'appareil a réussi les tests de sécurité pour la mise à la terre, le courant transitoire de ligne d'alimentation et la perte de courant.
	Input. Entrée.
	Output. Sortie.
Equipment labels are color coded:  Les étiquettes de l'appareil sont codées couleur:	Yellow Caution, risk of danger Red Stop Blue Mandatory action Green Safe condition or information  Jaune Attention, danger potentiel Rouge Arrêter Bleu Intervention obligatoire Vert Condition sûre ou informations de sécurité



## 2 Instrument Safety

The following safety information about the Phenolmager HT is included in this documentation. Read and review all safety information before operating the Phenolmager HT.

- [Required Training](#) <sup>11</sup>
- [Electrical Safety](#) <sup>11</sup>
- [Mechanical Safety](#) <sup>13</sup>
- [Weight Warning](#) <sup>13</sup>
- [Bright Light](#) <sup>13</sup>

### 2.1 Required Training

Ensure that all personnel involved with the operation of the instrument have:

- Received instruction in general safety practices for laboratories.
- Received instruction in specific safety practices for the instrument.

#### WARNING



*Use this product only in the manner described in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.*

### 2.2 Electrical Safety

The Phenolmager HT is powered by a 100-120VAC/200-240VAC, 50-60Hz (±10%) input power supply.

The wall outlet or the power cable connector on the left side of the instrument should be accessible after the system's installation, to enable trained service personnel to safely disconnect power from the system during servicing.

#### WARNING



*Do not operate the system in an environment with explosive or flammable gases.*

#### WARNING





- *DO NOT remove instrument covers. There are no user serviceable parts inside. The covers are intended to be removed only by qualified Akoya Biosciences or PerkinElmer service personnel; they are not intended to be removed during operation or for maintenance by users. Contact Akoya Biosciences technical support if help is required (see [Contact Us](#)<sup>[6]</sup>).*
- *Do not operate the system if there has been a malfunction of the system door or slide loading components. Contact Akoya Biosciences technical support if help is required (see [Contact Us](#)<sup>[6]</sup>).*
- *Do not operate the system in places where it may be splashed with liquid.*

For further electrical safety information, refer to the following sections:

- [Power Cord Selection](#)<sup>[12]</sup>
- [Fuses](#)<sup>[12]</sup>
- [Cables and Adapters](#)<sup>[13]</sup>

### 2.2.1 Power Cord Selection

Contact Akoya Biosciences Technical Support (see [Contact Us](#)<sup>[6]</sup>) to order replacement power cords.

#### WARNING



- *Use only the power supply cord set provided with the Phenolmager HT system. If the correct cord set for the location was not provided, contact Akoya Biosciences Technical Support (see [Contact Us](#)<sup>[6]</sup>) for a replacement. Do not use power supply cords with inadequate ratings.*
- *Use only a properly grounded power outlet when connecting the system to power.*

### 2.2.2 Fuses

Contact Akoya Biosciences Technical Support (see [Contact Us](#)<sup>[6]</sup>) to order replacement fuses.

#### WARNING



*The fuses in this instrument are only replaceable by trained Akoya Biosciences or PerkinElmer personnel.*

### 2.2.3 Cables and Adapters

Some cables and adapters supplied with the system have proprietary specifications.

#### WARNING



*Do not connect components supplied by Akoya Biosciences using unqualified cables or adapters. Contact Akoya Biosciences technical support (see [Contact Us](#)<sup>[6]</sup>) to order replacement cables and adapters.*

## 2.3 Mechanical Safety

#### WARNING



*Instrument components may move during operation. Always keep body parts, hair, jewelry, and clothing away from the instrument during operation.*

#### WARNING



*Procedures which could result in injury may be performed only by operators who have been warned of the potential hazards and have received adequate training in performing the procedures in the safest possible manner.*

## 2.4 Weight Warning

#### WARNING



*LIFTING HAZARD. The Phenolmager HT instrument weighs 185 lbs. (84 kg). Do not move the Phenolmager HT instrument. Installing, servicing, and moving the Phenolmager HT instrument should be performed only by qualified Akoya Biosciences or PerkinElmer service personnel. Contact Akoya Biosciences technical support if [help is required](#) (see [Contact Us](#)<sup>[6]</sup>).*

## 2.5 Bright Light

#### WARNING



*BRIGHT LIGHT HAZARD. The interior of the Phenolmager HT system includes a barcode reader with a Class 2 LED Light. Do not look into the bright light to avoid eye injury.*

CAUTION - CLASS 2 LASER RADIATION WHEN OPEN  
DO NOT STARE DIRECTLY INTO THE BEAM

## 3 Introduction

This manual describes the use and functionality of the Phenolmager HT Automated Quantitative Pathology Imaging System. It includes operating instructions, functional descriptions, troubleshooting, illustrations, and other relevant information.

This section of the manual contains the following topics:

- [Intended Use](#) <sup>[15]</sup>
- [Principles of Operation](#) <sup>[15]</sup>
- [Theory of Imaging](#) <sup>[17]</sup>
- [Whole Slide Scanning](#) <sup>[24]</sup>
- [Multispectral Field Imaging](#) <sup>[25]</sup>

### 3.1 Intended Use


The Phenolmager HT is a multimodal digital pathology instrument that integrates both multispectral analysis and automated slide scanning that allows researchers to visualize, analyze, quantify and phenotype immune cells in situ in FFPE tissue sections and TMAs.

**NOTE:** *Akoya's Phenolmager HT Imaging System is for research use only.  
Not for use in diagnostic procedures.*

### 3.2 Principles of Operation

Akoya's Phenolmager HT is an automated imaging system for performing whole slide scans of tissue sections and microarrays (TMAs), and for acquiring multispectral (MSI) regions of interest. The system has been optimized to image samples stained with Akoya's Opal™ multiplexed fluorescent immunohistochemistry reagent kits and is also compatible with typical brightfield staining reagents.

The Phenolmager HT is configured to store and inventory up to 20 slide carriers, each holding up to 4 tissue slides, for a total of 80 slides.



*Figure 1. The Phenolmager HT System*

The Phenolmager HT has been designed to expand and support a high-throughput workflow and includes:

- True whole slide scanning of slides at 1.0  $\mu\text{m}/\text{pix}$ , 0.5  $\mu\text{m}/\text{pix}$ , and 0.25 $\mu\text{m}/\text{pix}$
- Review and annotation of whole slides scans for multispectral field acquisition
- Acquisition of multispectral fields
- Analysis of multispectral whole-slide regions and fields including protein expression and phenotyping

To implement the full capabilities of the instrument and workflow, the Phenolmager HT system includes the following Akoya software:

- Phenolmager HT: Operator-centric software for performing whole slide scans and acquiring multispectral fields of interest. The Phenolmager HT software runs on the workstation connected to the Phenolmager HT instrument.
- Phenochart™: Whole-slide viewer and annotator of fluorescent and brightfield scans acquired by the Phenolmager HT. Phenochart allows the user to view the whole slide (zoom, pan, etc.), and make decisions (annotations) on next steps for the sample. Annotations in Phenochart are also used to record the workflow actions for each slide scan. Annotations include reviewer requested multispectral fields; automated (inForm) field requests; and reviewer edits, approvals, and rejections. The annotation file is a fully auditable transaction log. Phenochart is freely distributed and can be used by multiple users who wish to view or review slide scans taken by the Phenolmager HT.
- inForm® Tissue Finder: Software typically used for the analysis of MSI images. inForm supports features such as tissue classification and training, cell phenotyping, protein expression measurements, and data export. It can be run on the Phenolmager HT computer and other Microsoft® Windows® 10 computers. Additional inForm software seats beyond those that come with each Phenolmager HT system are available for purchase.

Phenolmager HT workflows range from simple two-step procedures (e.g., acquire whole slide scan and review) to automated acquisition of regions of interest selected by the user or the Phenolmager HT itself. An example fluorescence workflow might include the following steps:

1. Stain tissue with Akoya Biosciences Opal™ fluorescent IHC reagents.
2. Acquire whole slide fluorescent imagery using the Phenolmager HT.
3. Review the whole slide imagery with Phenochart and annotate regions of interest for MSI analysis.
4. Acquire the MSI regions with Phenolmager HT.
5. Use inForm to phenotype cells and measure protein expression levels in the acquired MSI regions.

For example applications for the Phenolmager HT, see [Example Applications](#).

### 3.2.1 Example Applications

Some examples of Phenolmager HT applications include:

- Whole slide scanning and multispectral interrogation of tissue samples and microarrays stained with Opal™ reagent kits
- Whole slide scanning of tissue samples stained with H&E and conventional IHC stains
- Phenotypic analysis and protein expression of immune and cancer cells in the context of the tumor microenvironment.

## 3.3 Theory of Imaging

This section introduces some important concepts used by Akoya's Phenolmager HT imaging systems, including:

- [Light](#)<sup>[18]</sup>
- [Human Perception of Light Intensity and Color](#)<sup>[18]</sup>
- [Light Absorbance and Reflection](#)<sup>[19]</sup>
- [Fluorescence](#)<sup>[19]</sup>
- [Photobleaching](#)<sup>[20]</sup>
- [Filter Sets for Conventional vs. Multispectral Imaging](#)<sup>[21]</sup>
- [Multispectral Analysis](#)<sup>[22]</sup>

### 3.3.1 Light

For purposes of this discussion, light refers to the part of the electromagnetic spectrum that can be seen by the human eye and the nearby ultraviolet and infrared wavelengths. While the physical description of light can be highly complex, we will focus on these wavelengths of light and how they interact with physical and biological materials.




Figure 2. The Electromagnetic Spectrum

### 3.3.2 Human Perception of Light Intensity and Color

#### Response to Illumination

The human eye is a highly adaptive light detector. It is significantly more sensitive in low light than in bright light. When light levels change, it takes some time for the eye to fully adjust. This is the reason people need to “dark adapt” in a darkened room before they are able to observe weak fluorescence through a microscope.

Humans can see both in very dark and very bright settings. Because the eye is so adaptable to various lighting conditions, humans are unable to quantify absolute levels of light. In any given situation, the eye has a limited ability to discriminate levels of illumination. US Department of Defense research indicates that most humans can only distinguish approximately 30 to 35 levels of gray, ranging from black to white.

The eye's response to illumination is not linear. It more closely approximates a logarithmic function. The result is that the human eye cannot see small proportional changes in brightness.

Contrasting the eye with conventional microscope imagery, any sensor that has 8-bit resolution can detect 256 levels of gray. As the number of bits of resolution increase, the number of gray levels also increases. A 12-bit sensor yields 4096 levels of gray. Digital electronic sensors are linear in response to light levels.



## Ability to Distinguish Colors

While the eye is relatively poor at discriminating intensity, it is very good at distinguishing colors. Most individuals are able to discriminate thousands of colors. However, no two individuals see a given color in exactly the same way. The eye contains three different types of color sensors, also known as cone cells. While the arrangement of cone cells is generally standard from person to person, the ratio of each type of cone cell varies, as does their actual physical arrangement within the eye. These minor variations (along with the brain's interpretation of the color) lead to the differences in perceived color between individuals.

### 3.3.3 Light Absorbance and Reflection

We perceive objects based on the way they transmit, absorb, and/or reflect light.

**Absorbance** and **reflection** work in tandem. **Absorbance** refers to the wavelengths of light that are 'taken in by' the objects. This means that an object that we perceive as red has absorbed all visible wavelengths of light except red. The red wavelengths are **reflected** back to the eye of the observer.

**Transmittance** refers to light emitting objects such as light sources, and fluorescing or phosphorescing objects. An object we would perceive as red in transmission is one that transmits primarily red wavelengths, while absorbing or reflecting other wavelengths.

In brightfield light microscopy we observe light that passes through a specimen. Except for a few pigments and inclusions, biological specimens are essentially invisible. To impart contrast, we employ some absorbing dye, or specific optical arrangement. It is this need for contrast that led to the initial development of biological stains and stain protocols and subsequently to phase contrast and other optical contrast enhancing techniques.

**Optical Density (OD)** is used to measure the interaction of light with absorbing materials. The science of absorbing spectroscopy is based on the Beer-Lambert law. When absorbing images collected in brightfield are converted to OD images, the information contained in each pixel is quantitative, as to the amount of absorbing material present. Akoya's brightfield multispectral imagery is automatically converted to optical density at acquisition time, enabling quantitative analysis.

### 3.3.4 Fluorescence

Many biological and natural materials give off light of a particular color when exposed to light of another color. This property is a type of luminescence. There are two types of **luminescence**:

- **Fluorescence** refers to luminescence that occurs when the light is emitted rapidly after illumination (around one-millionth of a second).
- If the light emission takes longer than one-millionth of a second, the luminescence is called **phosphorescence**.

Materials that exhibit fluorescence have proven extremely useful as labels or indicators in many biological systems.

Fluorescence light emission is different than light absorption. Each fluorescent molecule generates light. Fluorescent light can be measured quantitatively because it does not interact with other

materials. While it would seem that fluorescence could be measured more accurately than absorbed light, there are a number of factors that complicate such measurements. For example, light scatters, it is affected by the local environment (such as pH), and the measurements can be affected by surrounding molecules.

### Stokes Shift

When you excite a specimen with a particular (shorter) wavelength (such as blue), the specimen then shines in a different (longer) wavelength (such as red, orange, or yellow). The difference between the wavelength of the (shorter) exciting light and the wavelength of the (longer) emitted light is called the **Stokes Shift**, which is based on **Stokes Law**.

The wavelength or color you use to excite the specimen (i.e. the 'excitation light') and the color the specimen glows (i.e. the 'emission light') depend on the dye involved. For any given fluorescence dye, there will be a range of excitation wavelengths that will excite fluorescence. This range of excitation wavelengths is known as the absorption spectrum. Each dye also emits across a range of wavelengths, known as the emission spectrum. Figure 3 contains an example of excitation and emission spectra, showing Stokes Shift and the overlap of the spectra.




Figure 3. Stokes Shift

In addition, many biological materials are naturally fluorescent. This is known as autofluorescence. In particular, many vitamins, some hormones, and a variety of biological enzymes and structural proteins are naturally fluorescent. These materials often fluoresce strongly enough to interfere with specific fluorescence labeling studies.

### 3.3.5 Photobleaching

Because dyes can be damaged by intense light, reducing the emission signal ('photobleaching'), it is important to limit the time they are exposed to excitation light or to bright light during routine handling. Usually, blue or UV light is the most damaging. The Phenolmager HT uses an electronically-gated excitation source synchronized with its camera so the sample is only exposed to light while the camera is taking an image. Also, the Phenolmager HT front door is made of a translucent plastic that absorbs harsh blue and UV light.

When using Phenolmager HT, avoid spending a long time in the Protocol Exposures Editor while fixed on any one spot of the sample, since it takes a live image stream.


These steps enable repeated measurements with minimal effect on the sample.

### 3.3.6 Filter Sets for Conventional vs. Multispectral Imaging


Filters used for conventional fluorescence imaging are often designed so they only transmit a very narrow range of wavelengths of light. Limiting the measurement to wavelength bands where the dye is inherently most responsive helps distinguish the desired dye from other dyes or background signals in the sample. In this way, it is possible to image several dyes, provided that their spectra are separated rather than overlapping. Based on the properties of common dyes, this puts an upper limit of ~4 on how many dyes can be imaged in any one sample.

An alternative approach is to image the sample multispectrally. In this case, a broad emission filter is used, and a tunable filter is engaged in the imaging path. The instrument takes pictures at several wavelengths within the emission band, so maps out the full shape of the dye response. This enables analysis software to identify what dye(s) are present, in what amounts, in each pixel, by spectral decomposition ("unmixing"). It also enables identifying, and removing, contributions from autofluorescence.

These two approaches are illustrated in Figure 4 and Figure 5



*Figure 4. Narrowband excitation (solid line) and emission (dashed line) filter for conventional imaging.*



*Figure 5. Narrowband excitation (solid line) and long-pass emission filter (dashed line) for multispectral imaging.*

Ideally, the excitation filter would match the excitation maximum of the fluorescence label being used, and the emission filter would include the emission maximum. In practical terms, the filter maxima may be slightly different from the ideal case, due to limitations of filter manufacturing and because for many dyes the Stokes shift is small, so the maxima are quite close to one another.

To image multiple dyes conventionally, one selects dyes that have very distinct excitation and/or emission response, and selects filters that are narrow enough to predominantly transmit the signal of only one dye at a time (Figure 6 and Figure 7).




Figure 6. Imaging four dyes by conventional fluorescence methods.





Figure 7. Seven dyes is too many to image by conventional methods.

This shows the emission spectra for four dyes (DAPI, FITC, Cy3, and Cy5.5) that are spectrally fairly distinct, and can be separated by conventional imaging. This works because there is only a little overlap between adjacent dyes. Excitation filters and dye response are omitted for clarity, but are similarly separable.

This approach breaks down when more dyes are present, or when it is important to account for the effects of autofluorescence. Figure 7 shows the spectra of 7 dyes, and with this many dyes there is no way to isolate their signals using conventional imaging techniques, or to account for autofluorescence.


### 3.3.7 Multispectral Analysis

The Phenolmager HT imaging system offers a unique solution to the problem of separating the signals from highly multiplexed samples. Multispectral analysis is based on the fact that all fluorescent materials produce a unique spectral emission. If you excite a material and examine the emitted fluorescence over a range of wavelengths, the resulting emission intensities can generate an “emission spectrum”. This spectrum is different for each specific fluorescent material. For many fluorescent labels of biological interest, the emission spectra overlap, and may be further obscured by autofluorescence from the specimen. Multispectral imaging provides a way to distinguish between many overlapping emission spectra within the same area, overcoming the limitations of conventional filter-based imaging. With the additional information provided by the LCTF during multispectral field imaging the system can distinguish between dyes with fully overlapping spectra within a single channel (Figure 8).



*Figure 8. Unmixed multispectral image of human breast cancer tissue stained against CD4 (Green), CD20 (Red), CD8 (Yellow), FoxP3 (Orange), CD68 (Purple), Cytokeratin (Light Blue), and DAPI (Dark Blue) using Opal reagents.*

In general, multispectral analysis generates the spectral curves for the various fluorescent dyes or materials in a specimen. In addition, it generates a spectral curve for the autofluorescence that almost always is present to some degree. Using spectral analysis algorithms, the contribution of the individual fluorescence spectra are separated. The result is a set of images representing each spectrum that contributes to the final image (Figure 9).



**Figure 9. Removing auto-fluorescence with multispectral imaging.** A) Conventional fluorescence image of tonsil tissue stained against CD4 with Opal 520 dye. B) Overlay of all 35 layers of a multispectral image acquired using five filter cubes. C) Emission spectra of pure autofluorescence (black line), pure Opal 520 (green line), and the mixture observed in 'B' (gray dashed line). The spectral references are used to 'unmix' the contribution of autofluorescence and Opal 520 at each pixel. D) Image of unmixed Opal 520 signal extracted from 'B' with > 10-fold higher signal-to-background compared to 'A' because the autofluorescence contribution has been separated into the component image in 'E'.

## 3.4 Whole Slide Scanning

Phenolmager HT scans slides using the following process:

### 1. Color Overview:

Phenolmager HT takes a low power color overview of all four slides in each carrier, including the label for each slide. This initial step is performed regardless of whether you are using a fluorescence or brightfield protocol and is used to identify the presence of slides in the carrier and capture their labels.

### 2. Coverslip Finding:

For each slide, Phenolmager HT will then find the coverslip using this overview scan. The coverslip defines the potential scan area.

For fluorescence protocols, you can further restrict this area by making a closed loop with a red, green, or blue Sharpie® marker. If closed-loop markup is present, the system will only scan within the loop. This is useful if your tissue is faintly stained, punctate, or if you have highly fluorescing non-tissue material (PAP pen, for example.)

For brightfield protocols, the coverslip will define the potential scan area. Closed-loop markup is not available for brightfield scans.

### 3. Slide Height Finding:

Phenolmager HT engages specialized height-sensing optics to measure the top of the coverslip at up to 9 locations. This gives an initial focus estimate based on the expected coverslip thickness.

#### 4. Fluorescence Overview (Fluorescence protocols only):

If you are using a fluorescence protocol for this slide, Phenolmager HT will take another overview within the coverslip (or closed-loop markup), this time in fluorescence.

#### 5. Tissue Finding:

Using the corresponding overview, Phenolmager HT will automatically detect the sample on the slide. The resulting area will be scanned. If requested within the protocol, Phenolmager HT will scan the entire area within the coverslip (or closed-loop markup).

#### 6. Focus Finding:

Phenolmager HT will measure focus at multiple points on the tissue to determine the best focus. It uses the sample map from the previous step to choose the measurement grid, and continues until the grid is fully measured.

If the measured tissue height is irregular, Phenolmager HT increases the grid density and takes more readings until it finds the readings are regular at the newly finer scale.

All focus measurements include a dust-rejection algorithm, and an overall consistency check is applied as well, to further reduce the likelihood of dust-induced focus errors.

#### 7. Scanning:

Phenolmager HT then scans the slide.

Brightfield scans are conventional color scans that have been color and background corrected.

Fluorescence scans are multi-layered, with one layer for each filter you chose. To avoid photobleaching, the system uses a pulsed LED so the sample is only exposed to light during the time that the camera is taking a picture.

When scans are complete, they can be opened in Phenochart, the whole slide viewing application.

**NOTE:** *Some dyes narrowly express in a single filter. Other dyes may express in multiple filters and may appear in more than one layer in a Phenolmager HT fluorescence scan. For example, Opal 570 will have signal in both Cy3 and TexasRed filters. If your sample is highly multiplexed, multiple dyes may appear in the same channel.*

## 3.5 Multispectral Field Imaging

After the Phenolmager HT has completed a whole slide scan of the tissue, individual fields can be selected for multispectral imaging.

This section contains:

- [Configuring the Multispectral Fields](#) <sup>[26]</sup>
- [Processing the Multispectral Fields](#) <sup>[26]</sup>



### 3.5.1 Configuring the Multispectral Fields

Multispectral imagery is acquired using the following process:

#### 1. Selection of Multispectral Fields:

Fields for multispectral imaging are selected on a previously scanned slide. Using Phenochart, you can select individual fields or regions of interest. If desired, you can also select fields using a trained inForm algorithm. See the Phenochart documentation for more information.

### 3.5.2 Processing the Multispectral Fields

After you configure the slide for field acquisition, Phenolmager HT will perform the following actions:

#### 1. Color Overview:

Phenolmager HT takes a low power color overview of all four slides in each carrier, including the label for each slide. This initial step is performed regardless of whether you are using a fluorescence or brightfield protocol and is used to identify the presence of slides in the carrier.

#### 2. Slide Registration:

Using the color overview along with the slide's original overview, Phenolmager HT accounts for any shift or rotation of the slide to ensure that the multispectral region locations are accurate. The slide edges are used to account for any rotation or horizontal shift. The coverslip edges to account for vertical shifting.

#### 3. Slide Height Finding:

Phenolmager HT engages specialized height-sensing optics to measure the top of the coverslip at up to 9 locations. This gives an initial focus estimate based on the expected coverslip thickness.

#### 4. Acquisition of Multispectral Fields:

Phenolmager HT will then travel to each multispectral region site, autofocus, correct for exposure if requested, and acquire the multispectral image.

Multispectral imagery can then be viewed and analyzed in inForm.

## 4 Hardware Reference

This section identifies and describes the Phenolmager HT system hardware. It also lists the Phenolmager HT technical specifications.

### WARNING:



*Lifting Hazard. Do not move the Phenolmager HT instrument. Installing, servicing, and moving the Phenolmager HT instrument should be performed only by qualified Akoya Biosciences or PerkinElmer service personnel. Contact Akoya Biosciences technical support if help is required (see [Contact Us](#)).*

This section contains the following information:

- [Front View](#) <sup>27</sup>
- [Top View](#) <sup>29</sup>
- [Right-Side View](#) <sup>30</sup>
- [Left-Side Connectors](#) <sup>31</sup>
- [Slide Carrier Hotel](#) <sup>32</sup>
- [Slide Carrier](#) <sup>33</sup>
- [Barcode Reader](#) <sup>35</sup>
- [System Computer and Monitor](#) <sup>35</sup>
- [Specifications](#) <sup>36</sup>

### 4.1 Front View

Figure 10 shows the front view of the Phenolmager HT instrument with the system door *open*. The parts identified are described in Table 1.




Figure 10. Phenolmager HT Front View - System Door Open

Table 1. Front View - Hardware Components	
Part	Description
<b>Status Progress Display</b>	Each light represents an individual slide and indicates the progress of the slide currently being scanned.
<b>Slide Carrier Hotel</b>	Houses up to 20 slide carriers. See <a href="#">Slide Carrier Hotel</a> <sup>[32]</sup> for more detailed information.
<b>Door Sensor</b>	Opens or closes the <a href="#">System Door</a> <sup>[29]</sup> when a hand is placed in front of the sensor.

Figure 11 shows the front view of the Phenolmager HT instrument with the system door *closed*. The parts identified are described in Table 2.




Figure 11. Phenolmager HT Front View - System Door Closed

Table 2. Front View - System Door Closed- Hardware Components	
Part	Description
Slide Carrier Status Lights	Each <a href="#">Slide Carrier Status Light</a> <sup>[32]</sup> represents an individual slide carrier and indicates the processing status of each <a href="#">Slide Carrier</a> <sup>[33]</sup> . See <a href="#">Table 7</a> <sup>[32]</sup> for the color codes of the slide carrier status lights.
System Door	When closed, covers the <a href="#">Slide Carrier Hotel</a> <sup>[32]</sup> and slide processing can occur. When open, reveals the <a href="#">Slide Carrier Hotel</a> <sup>[32]</sup> .

**WARNING:**



Do not operate the system if there has been a malfunction of the system door or slide loading components. Contact Akoya Biosciences technical support if help is required (see [Contact Us](#)<sup>[6]</sup>).

## 4.2 Top View

Figure 12 shows the top view of the Phenolmager HT instrument. The part identified is described in Table 3




Figure 12. Phenolmager HT Top View

Table 3. Top View - Hardware Components	
Part	Description
Cover Handle	Used by only service personnel to remove instrument covers during service and installation.

**WARNING:** *DO NOT remove instrument covers. There are no user serviceable parts inside. The covers are intended to be removed by qualified service personnel only; they are not intended to be removed during operation or for maintenance by users. Contact Akoya Biosciences Technical Support if help is required (see [Contact Us](#)<sup>67</sup>).*

### 4.3 Right-Side View

Figure 13 shows the right-side of the Phenolmager HT instrument. The part identified is described in Table 4





Figure 13. Phenolmager HT Right-Side View

Table 4. Right-Side View- Hardware Components	
Part	Description
<b>System Power Switch</b>	Turns the Phenolmager HT instrument ON   or OFF 

## 4.4 Left-Side Connectors

Figure 14 identifies the connectors on the bottom left-side of the Phenolmager HT instrument. The connectors identified are described in Table 5




Figure 14. Phenolmager HT Left-Side Connectors

### WARNING:

- Appliance inlet is a disconnecting device. Place device or equipment in a manner so that disconnecting device is accessible at all times.
- Use only the power supply cord set provided with the Phenolmager HT system. If the correct cord set for the location was not provided, contact Akoya Biosciences Technical Support (see [Contact Us](#)<sup>[36]</sup>) for a replacement.
- Use only a properly grounded power outlet when connecting the system to power.

Table 5. Left-Side View - Connection Components	
Part	Description
<b>Power Connector</b>	Connects to the power cord to provide power to the Phenolmager HT instrument.
<b>USB 2.0 Connector</b>	Connects a USB 2.0 cable to a USB 2.0 port on the <a href="#">System Computer Connectors</a> <sup>[36]</sup> .
<b>USB 3.0 Hub Connector</b>	Connects a USB 3.0 HUB to a USB 3.0 port on the <a href="#">System Computer Connectors</a> <sup>[36]</sup> .

# 4.5 Slide Carrier Hotel

The Slide Carrier Hotel houses up to 20 [Slide Carriers](#)<sup>[33]</sup> and is visible when the [System Door](#)<sup>[29]</sup> opens. Figure 15 shows a closeup of the Slide Carrier Hotel and identifies its components. The components are described in Table 6.

**WARNING:** Do not operate the system if there has been a malfunction of the system door or slide loading components. Contact Akoya Biosciences technical support if help is required (see [Contact Us](#)<sup>[6]</sup>).




Figure 15. Slide Carrier Hotel

Table 6. Slide Carrier Hotel - Hardware Components	
Part	Description
Slide Carrier Status Lights	Indicates the processing status of each slide carrier. <a href="#">Table 7</a> <sup>[32]</sup> lists the color codes associated with the slide carrier status lights.
Hotel Slot	Each hotel slot holds one <a href="#">Slide Carrier</a> <sup>[33]</sup> . The hotel slots are numbered 1-20, starting at the top of the hotel.
Slide Carrier	Houses up to four slides. See <a href="#">Slide Carrier</a> <sup>[33]</sup> for more detailed information.

Table 7. Slide Carrier Status Lights - Color Codes	
Color	Status



None	Slide carrier hotel is empty.
White	Initial state of slot after slide carrier is inserted and no rules have been applied via software.
Solid Red	Hotel slot is malfunctioning.
Blinking Red	Slide carrier is not completely inserted.
Blue	Awaiting Processing.
Blinking Blue	Prioritized and awaiting processing.
Yellow	Processing instructions are either incomplete or invalid.
Blinking Green	Processing.
Solid Green	Processing Complete.
Orange	Processing complete but a processing error occurred.

## 4.6 Slide Carrier

Slide Carriers are an integral component of Phenolmager HT. The Slide Carrier holds up to four microscope slides and helps protect slides from damage. The microscope slide positions are labeled 1 to 4 and correspond to slide identification in the Phenolmager HT software.

16 shows a closeup of a Slide Carrier and identifies its components. The identified components are described in Table 8



Figure 16. Phenolmager HT Slide Carrier


Table 8. Phenolmager HT Slide Carrier - Hardware Components	
Part	Description
Slide Position	Each slide position is populated with one slide. Slide position 1 is on the far left; slide position 4 is on the far right.
Slide	Each microscopic slide is manually loaded into the slide positions. <a href="#">Table 9</a> <sup>[34]</sup> identifies the details of the slide format.
Insert Indicator	An icon  that indicates the side of the slide carrier to be inserted into a <a href="#">Hotel Slot</a> <sup>[32]</sup> .
Tab Cover	Holds the spring-loaded tabs in place.
Spring-Loaded Tab	When a slide is inserted into a slide slot, it is gently placed up against a spring-loaded tab.
Carrier Handle	The side of the slide carrier to hold when inserting and removing a slide carrier from a <a href="#">Hotel Slot</a> <sup>[32]</sup> . A unique number is printed on top of the handle for slide carrier identification purposes.

Table 9. Details of Slide Formats			
Slide Type	Width	Height	Thickness
Metric	25.0 ± 1.0mm	75.0 ± 1.0mm	1.00 ± 0.10mm

English	25.4mm	76.2mm	1.00 ± 0.10mm
---------	--------	--------	---------------

## 4.7 Barcode Reader

The Phenolmager HT is equipped with an internal Barcode Reader. The Barcode Reader reads the barcode on the slide labels.

### WARNING



**BRIGHT LIGHT HAZARD.** The interior of the Phenolmager HT system includes a barcode reader with a Class 2 LED Light. Do not look into the bright light to avoid eye injury.

CAUTION - CLASS 2 LASER RADIATION WHEN OPEN  
DO NOT STARE DIRECTLY INTO THE BEAM

## 4.8 System Computer and Monitor

The Phenolmager HT system includes a widescreen monitor and a computer pre-installed with Phenolmager HT, Phenochart, and inForm software. A wireless keyboard and mouse are also included. See below for a description of the computer connectors.

Figure 17 shows the Phenolmager HT computer connectors. The connectors identified are described in Table 20




Figure 17. System Computer Connectors

Table 10. System Computer Connections	
Part	Description
USB 2.0 Port	Connects a USB 2.0 cable to the USB 2.0 connector on the instrument's <a href="#">Left-Side Connectors</a> <sup>31</sup> .
USB 3.0 Port	Connects a USB 3.0 HUB to the USB 3.0 HUB connector on the instrument's <a href="#">Left-Side Connectors</a> <sup>31</sup> .

## 4.9 Specifications

This section lists the technical specifications for the Phenolmager HT instrument. Technical specifications are subject to change without notice.

**NOTE:** *Phenolmager HT is for research use only. Not for use in diagnostic procedures.* **System**

<b>Size (H x W x L)</b>	28" (72 cm) x 30" (77 cm) x 27.25" (69.2 5cm)
<b>Weight</b>	185 lbs. (84 kg)
<b>Spectral Range</b>	440 nm to 780 nm
<b>Objectives</b>	4x, 10x, and 20x
<b>Pixel Resolution</b>	0.25 um, 0.5 um, or 1.0 um
<b>File Format</b>	Akoya Biosciences .qptiff format for whole slide scans, proprietary .im3 file format for multispectral data; 24-bit Windows-compatible bitmap for RGB/Mono imagery
<b>Operating System</b>	Microsoft® Windows 10, 64-bit
<b>RAM</b>	16GB

### Environmental

<b>Operating Temperature</b>	59°F to 83°F (18°C to 28°C)
<b>Operating Humidity</b>	30% - 50% non-condensing
<b>Storage Temperature</b>	59°F to 86°F (15°C to 30°C)
<b>Storage Humidity</b>	0% - 80% relative humidity, non-condensing
<b>Altitude</b>	Up to 2000m (6560ft.)

<b>Shipping Temperature (up to 72 hours max)</b>	14°F to 113°F (-10°C to +45°C)
<b>Pollution Degree</b>	2
<b>Indoor Use Only</b>	

### Electrical

One properly grounded AC power outlet for the computer, monitor, and instrument must be located within 6 feet (1.8 m) of the location.

<b>Input Voltage</b>	100 to 230 VAC ( $\pm 10\%$ ), 500 W, 50/60 Hz  System does not have transient overvoltage protection.
<b>Computer Interface</b>	USB 2.0, USB 3.0
<b>Fuse</b>	4A Littelfuse® 250V, 5mm x20mm

### Barcode Reader Light Source (White laser diode)

<b>Wavelength</b>	400 nm to 750
<b>Power Output</b>	70mAmps

## 5 Operation

This section describes how to operate the Phenolmager HT hardware. It is important to read and understand [Instrument Safety](#) <sup>11</sup> before using the system. If you are not familiar with the Phenolmager HT system hardware, [Hardware Reference](#) <sup>27</sup> provides a hardware component description of each [Reference](#) in the Phenolmager HT system.

This section includes the following procedures:

- [System Startup](#) <sup>38</sup>
- [System Shutdown](#) <sup>39</sup>
- [Inspecting Slides and Slide Carriers](#) <sup>39</sup>
- [Loading Slides into the Slide Carriers](#) <sup>40</sup>
- [Loading Slide Carriers into the Slide Carrier Hotel](#) <sup>41</sup>
- [Removing Slide Carriers from the Slide Carrier Hotel](#) <sup>43</sup>
- [Removing Slides from the Slide Carriers](#) <sup>43</sup>

### 5.1 System Startup

This section describes the procedures needed to start the Phenolmager HT system:

- [Turn on the Phenolmager HT Instrument](#) <sup>38</sup>
- [Launch the Phenolmager HT Software](#) <sup>39</sup>

#### 5.1.1 Turn on the Phenolmager HT Instrument

To start the Phenolmager HT instrument:

1. If necessary, plug the Phenolmager HT power cord into an appropriate power outlet.

**WARNING:**

- Use only the power supply cord set provided with the Phenolmager HT system. If the correct cord set for the location was not provided, contact Akoya Biosciences Technical Support (see [Contact Us](#) <sup>6</sup>) for a replacement. Do not use power supply cords with inadequate ratings.
- Use only a properly grounded power outlet when connecting the system to power.
- Appliance inlet is a disconnecting device. Place device or equipment in a manner so that disconnecting device is accessible at all times.

2. Turn on the computer and allow Windows<sup>®</sup> 10 to start.
3. Switch the [System Power Switch](#)<sup>[31]</sup> to the ON (I) position. The [Status Progress Display](#)<sup>[28]</sup> flash green to indicate a quiescent state. The [Door Sensor](#)<sup>[28]</sup> is responsive and functional.

### 5.1.2 Launch the Phenolmager HT Software

To launch the Phenolmager HT software:

1. Double-click the Phenolmager HT icon on the Windows<sup>®</sup> 10 desktop.
2. The Phenolmager HT home page opens (see [Software Overview](#)<sup>[24]</sup>).

## 5.2 System Shutdown

This section describes how to shut down the Phenolmager HT

system. To shutdown the Phenolmager HT system:

1. Exit the Phenolmager HT software. If open, the [System Display](#)<sup>[29]</sup> closes.
2. Select Shut down from the Windows<sup>®</sup> Start Menu to power down the computer.
3. Switch the [System Power Switch](#)<sup>[31]</sup> on the right side of the instrument to the OFF (O) position. The status lights turn off and the system shuts down.

## 5.3 Inspecting Slides and Slide Carriers

Before inserting slides into slide carriers, both should be inspected for potential defects.

This section describes the steps needed to properly inspect the slides and slide carriers.

### Inspecting Slides:

- Verify the slides meet the required formats and dimensions (see [Table 9](#)<sup>[34]</sup>).
- Do not use broken or damaged slides, or slides with broken or damaged coverslips.
- Use only slides that are free of debris, fingerprints, and dust.

### Inspecting Carriers:

- Verify the slide Carrier Tab Cover is secure.
- Do not use slide carriers that are warped or bent.


- If any sticky residue is on the carrier handle or outer surfaces, clean before use (see [Cleaning the Slide Carriers](#)<sup>[88]</sup>).

## 5.4 Loading Slides into the Slide Carriers

After the Slides and Slide Carriers have been successfully inspected, slides can be loaded into the Slide Carriers.

To load a slide into a slide carrier:

1. Place the slide carrier onto a flat surface.
2. Hold the microscope slide by the label end with the coverslip side up.
3. Gently push the opposite end of the slide into the desired [Slide Position](#)<sup>[34]</sup> until the slide is up against the [Spring-Loaded Tab](#)<sup>[34]</sup>.



*Figure 18. Insert Slide into Slide Carrier*

4. Gently lay the label end of the slide completely into the slide position. The tab should push the slide against the opposite wall of the slide position.
5. Load up to four slides into each slide carrier. Slides need not be contiguous.
6. Click the **Enter Slide IDs** button (see [Slide IDs](#)<sup>[78]</sup>) in the Phenolmager HT software. The Enter Slide IDs dialog opens.
7. Type the applicable slide IDs into the text boxes and click the **OK** button.



8. To load additional slides into another slide carrier, repeat steps 1<sup>[40]</sup> to 7<sup>[40]</sup> above. Slide carriers may be stacked on top of one another for easy handling and storage (see Figure 19).



Figure 19. Slide Carrier Stack


## 5.5 Loading Slide Carriers into the Slide Carrier Hotel

The Slide Carrier Hotel stores Slide Carriers before and after microscope slide scanning. The slide carrier hotel can store up to 20 slide carriers for a total of 80 slides. The Hotel Slots are numbered 1-20, starting at the top of the hotel.

**NOTE:** Slide carriers can be loaded into the slide carrier hotel before or after launching the software.

To load a slide carrier into the slide carrier hotel:

1. Inspect the slide carrier (see [Inspecting Slide Carriers](#)<sup>[39]</sup>).
2. Verify the Phenolmager HT System is on (see [System StartUp](#)<sup>[40]</sup>).
3. If the System Door is closed, place your hand in front of the Door Sensor to open it. The system doors slides to the left to reveal the slide carrier hotel.
4. Hold the slide carrier by the Carrier Handle with the slide labels facing upwards.

5. Gently insert the end of the slide carrier with the Insert Indicator  into a Hotel Slot.



*Figure 20. Ends of a Slide Carrier*

6. When the Slide Carrier is halfway into the Hotel Slot, the carrier engages a roller and creates a small amount of friction. Push the carrier further into the slot until it is fully seated.



*Figure 21. Push the Slide Carrier until Fully Seated*

- If the carrier is inserted correctly, a click noise indicates that the carrier has been seated correctly into the slot. The Slide Carrier Status Light ([Slide Carrier Status Lights](#)<sup>32)</sup>) next to the associated Hotel Slot turns white.
  - If the carrier is not inserted correctly, the Slide Carrier Status Light next to the Hotel Slot blinks red until the slide carrier is removed from the slot or is positioned correctly.
7. Continue to load up to 20 slide carriers into the slots.
  8. When done, place your hand in front of the Door Sensor to close the System Door.

## 5.6 Removing Slide Carriers from the Slide Carrier Hotel

When slide scanning is complete, the [Slide Carriers](#)<sup>33)</sup> can be removed from the [Slide Carrier Hotel](#)<sup>32)</sup>.

To remove a slide carrier from the hotel:

1. Pull the [Carrier Handle](#)<sup>34)</sup> until the slide carrier is free from its [Hotel Slot](#)<sup>32)</sup>. The [Slide Carrier Status Lights](#)<sup>32)</sup> next to the Hotel Slot turns off, indicating no slide carrier is inside the Hotel Slot.
2. Lay the slide carrier onto a flat surface.

## 5.7 Removing Slides from the Slide Carriers


After [Slide Carriers](#)<sup>33)</sup> are removed from the [Slide Carrier Hotel](#)<sup>32)</sup>, the microscopic slides can be removed from the Slide Carriers.

To remove a slide from a slide carrier:

1. Gently push against the slide label end of the microscope slide using your index finger, compressing the tab on the far end of the slide.
2. Using the same finger, lift the label end of the microscope slide from the slot.
3. Grab the label end of the slide with your thumb and index finger.
4. Remove the remaining end of the slide from the slide carrier.

## 6 Software Operation

This section describes how to operate the Phenolmager HT software.



*Figure 22. Phenolmager HT Software*

This section includes the following topics:

- [Software Overview](#)<sup>[44]</sup>
- [System Dashboard](#)<sup>[47]</sup>
- [Creating and Editing Protocols](#)<sup>[53]</sup>
- [Scanning Slides](#)<sup>[75]</sup>

### 6.1 Software Overview

The Home Page links to the pages needed to maintain and run the Phenolmager HT.




Figure 23. Phenolmager HT Home Page

The Home Page buttons link to the following pages:

[Check Dashboard](#)<sup>[46]</sup> - Check your remaining disk space and acquire references.

[Edit Protocol](#)<sup>[46]</sup> - Create or edit brightfield and fluorescent protocols and studies.

[Scan Slides](#)<sup>[46]</sup> - Select the rules to scan your slides and scan your slides based on those rules.

[Launch Phenochart](#)<sup>[46]</sup> - Launch the Phenochart viewing software.

You can also get to those pages by using the [Gear Menu](#)<sup>[46]</sup> (in the upper right):




Figure 24. Gear Menu



Additionally, you can use the [Gear Menu](#)<sup>[46]</sup> to check your version number or perform other functions, and launching the online help.

Below is an overview of each page.

### Dashboard:

Use the Dashboard to check the status of your system prior to imaging. It includes the following:

- Disk space - checks the available space on the disk drive where images will be stored.
- Acquire References - views and acquires brightfield references and fluorescent references.

For detailed information on how to use the System Dashboard, see [System Dashboard](#)<sup>[47]</sup>.

### Edit Protocol:

Use the Edit Protocol page to define protocols and take snapshots.

A protocol describes how a sample is to be imaged, including the imaging mode (brightfield or fluorescence), and the spatial resolution (magnification) for the whole slide scan and for multispectral regions of interest. For fluorescent imaging, it also describes the exposure times and what bands to use for focusing and imaging.

For detailed information on how to create and edit protocols and studies (see [Creating and Editing Protocols](#)<sup>[53]</sup>).

### Scan Slides:

Use the Scan Slides page declare how slides should be imaged. Slide scanning can be started and stopped from this page.

For detailed information on how to perform whole slide scans (see [Scanning Slides](#)<sup>[75]</sup>).

### Launch Phenochart:

Use this button to launch the Phenochart program.

For detailed information on how to use Phenochart, see the Phenochart User's Manual or online Help from within the Phenochart program.

## 6.1.1 Gear Menu

Use the Gear menu (in the upper right hand corner of the Home Page) to:

- Link to the following pages:
  - [Check Dashboard](#)<sup>[47]</sup>

- [Edit Protocol](#)<sup>[53]</sup>
- [Scan Slides](#)<sup>[75]</sup>
- [Launch Phenochart](#)<sup>[46]</sup>
- Find out What's New in the current installed version
- Launch the online Help system
- Open the Preferences dialog
- Open the Phenomager HT software Log for reference or troubleshooting
- Open the About window to view the current software version
- View contact information for Akoya Biosciences technical support

## 6.2 System Dashboard

Use the System Dashboard to:

- See the available Disk Space (see [Disk Space](#)<sup>[48]</sup>)
- View and acquire Brightfield References (see [Brightfield References](#)<sup>[48]</sup>)
- View and acquire Fluorescence References (see [Fluorescence References](#)<sup>[50]</sup>)




Figure 25. System Dashboard Window

### 6.2.1 Disk Space

The disk space bar shows where your data is stored and how much space is available on that drive.

### 6.2.2 Brightfield References

#### Viewing Brightfield References

---

Click **View** to see reference imagery and information for the label image, brightfield overview, color image and multispectral image.

- Export For Diagnostics will save this image to aid in technical support.
- Show raw images will display the images without scaling, which is useful for technical support.
- [Compensation Information](#)<sup>52</sup> describes how these reference images are applied.






Figure 26. View Brightfield References Window

## Acquiring Brightfield References

---

Click **Acquire** to take new brightfield references. Use the stage control to:

- Load the Reference Carrier onto the slide.
- Move around to a clean area in the live view (no coverslip lines, tissue, or label) using the stage navigation tool.
- Click within the slide to change positions. You can refine the position by using arrow keys for small movements, and control+arrow for slightly larger movements.
- Click **OK** to take references.




Figure 27. Prepare for Brightfield Reference Acquisition

### 6.2.3 Fluorescence References

#### Viewing Fluorescence References

---

Click **View** to see reference imagery and information for the overview, 20x and 40x resolution references for each filter.

- Export For Diagnostics will save this image to aid in technical support.
- [Compensation Information](#)<sup>52)</sup> describes how these reference images are applied.




Figure 28. View Fluorescence References Window

## Acquiring Fluorescence References

Click **Acquire** to take new fluorescent references. We recommend you take references for all filters at once, but you can take references for a single filter if necessary. Acquiring references for all filters can take over 45 minutes, and the system can be left unattended during this time.

Use the stage control to:

- Load the Reference Carrier onto the stage.
- Move to a tissue-filled location on the fluorescent reference slide using the stage navigation tool.
- Click within the slide to change positions. You can refine the position by using arrow keys for small movements, and control+arrow for slightly larger movements.
- Click **OK** to take references.




Figure 29. Acquire Fluorescence References Window

## 6.2.4 Compensation Information

### Brightfield

Acquired images of a sample are normalized (divided) by the reference images on a pixel-by-pixel, wavelength-by-wavelength basis. This yields the sample transmission  $T$ , which is in the range 0 - 1. This transmission ( $T$ ) is then mapped:

- For multispectral images, the Optical Density is calculated - this is based on  $\text{Log}_{10}(T)$ .
- For simple color images, like a whole slide scan, the transmission is mapped from 0 - 255.

### Fluorescence Protocols

Acquired images are normalized by a shading pattern derived from the reference images. The shading pattern, which is the reference image divided by the mean intensity in the center, is applied on a pixel-by-pixel basis, with one pattern per epi-filter. Here, center means the central 1/16th of the image area. While the exact shape of the shading pattern varies per instrument, the overall effect is to increase the signal near the image edges, and to do little or nothing to the signal from the center of the image.

## 6.3 Creating and Editing Protocols

This section contains the following topics:

- [Studies](#) <sup>[53]</sup>
- [Creating Protocols](#) <sup>[54]</sup>
- [Editing Protocols](#) <sup>[56]</sup>

### 6.3.1 Studies

Phenolmager HT stores scan data by Study.

A study is a group of slides that belong together. This could be an experimental study (e.g. Ki67 markers in breast cancer tissue), all slides from one source, or other groupings. Each study contains one or more slides. Each slide may be scanned more than once, if needed.

- The default location for a study is D:\Data\VectraPolaris\[Study] (where [Study] is the name of the study).
- Whole-slide scans and supporting imagery acquired from specific slides are saved to slide-specific subfolders in the main study folder. See the section on Scanning Slides for more details about imagery.
- Multispectral Fields acquired from a particular slide scan will be stored in an MSI folder within the scan folder.

A protocol defines the set of rules to be used during whole slide and multispectral region acquisition, including imaging mode, pixel resolution, filter cubes, exposure times, and other parameters.

- Protocols have the file extension “.ppr” and are saved in D:\Data\VectraPolaris\[Study].




Figure 30. Study structure

### 6.3.2 Creating Protocols

Before you create a protocol, it is helpful to know how the slide was stained.

- Brightfield protocols are used to acquire imagery from slides stained with H&E or conventional chromogenic IHC methods.
- Fluorescence protocols are used with Akoya Biosciences Opal and other fluorescent dyes.

To create a new protocol:

- Select Edit Protocol from the PhenolMager HT Home Page.
- Click the **New** button. The Create New Protocol window opens.




Figure 31. Create New Protocol Window

- Enter a Protocol name.
- Select Brightfield or Fluorescence under Imaging Mode.
- If you have upgraded your system for multispectral whole slide scans, choose the type of staining kit used for this protocol, or choose to select your own filters for imaging.
- Select a previously created Study or create a new Study.
  - To select a previously created study, click on the study in the Available Studies list. This will highlight the study.
  - To create a new study, enter the **Study Name** in the text box and click the **Create Study** button. The new study will be added to the Available Studies list where you can select the study.
- Click the **Create Protocol** button or the Select Scan Bands button to create the protocol in the selected study.
- If you chose to select your own filters, the Select Scan Bands window opens. Select the bands to use for imaging and click the **OK** button.

To load an existing protocol, click **Load** and select the protocol from the study folder.

### 6.3.3 Editing Protocols

After you have created your brightfield or fluorescence protocol and assigned it to a study, use the Edit Protocol screen to add specific details to the protocol. The next sections are organized by two different types of protocols: brightfield and fluorescence. The sections appear as follows:

- [Brightfield Protocols](#) <sup>56</sup>
- [Brightfield Snapshots](#) <sup>59</sup>
- [Fluorescence Protocols](#) <sup>62</sup>
- [Fluorescence Exposures](#) <sup>69</sup>

#### Brightfield Protocols

After you have created your brightfield protocol and assigned it to a study, the Edit Protocol window (for brightfield protocols) opens.




Figure 32. Edit Protocol Window

Under Slide Scan Settings, select the **Pixel Resolution** in the drop-down list to select the resolution to apply to the whole slide. You can choose 1 um per pixel (nominally 10x), 0.5 um per pixel (nominally 20x), or 0.25 um per pixel (nominally 40x).




Figure 33. Slide Scan Settings- Pixel Resolution Drop-down List

If you plan to take Multispectral Fields, select the Pixel Resolution in the drop-down list to use for imaging. Again, available resolutions are 1 um per pixel (nominally 10x), 0.5 um per pixel (nominally 20x), or 0.25 um per pixel (nominally 40x).






Figure 34. Multispectral Field Settings - Pixel Resolution Drop-down List

## Advanced Settings

There may be situations when you need to fine-tune how Phenolmager HT scans the slides. Click the **Advanced Settings** button to adjust the following functions:

- **Coverslip Thickness** - Phenolmager HT is designed for use with standard #1.5 coverslips, where the distance between the tissue and the top of the sample is approximately 170 microns. While this is the optimal setup, Phenolmager HT can handle other scenarios.
  - If your lab uses thinner coverslips (#0), select “100 um (#0 or similar)” for a Coverslip Thickness.
  - If your lab uses thicker coverslips, or has a significant amount of mounting media between the sample and coverslip, select the option that best matches your slides.

**NOTE:** Phenolmager HT is capable of handling slides with up to 420 microns of combined coverslip or mounting medium above the sample, though imagery at these extremes may be sub-optimal.

- **Scan within the entire coverslip region** - If Phenolmager HT is having difficulty finding your entire tissue sample, check the Scan within the entire coverslip region checkbox. This will increase scan time and file size, but will enable you to complete scanning of difficult samples.
  - **Sample is a TMA** - Phenolmager HT has a specialized algorithm for finding TMA samples. If your samples are TMAs, select Sample is a TMA.
  - **Image Compression** - By default, brightfield whole slide scans in Phenolmager HT are .jpg compressed in order to save disk space. You can adjust the image quality as needed; higher quality will result in larger files. If you would like use lossless compression instead, select LZW rather than JPEG.
  - **Label Barcode Reading** - If you would like Phenolmager HT to decode the barcode on your slide label and save it as part of the scan file, select the type of barcode you are using. By default, Phenolmager HT will not attempt to decode a barcode unless you specifically select this option.
    - If you have only 1D barcodes, select All 1D barcodes (autodetect type).
    - If your lab uses a variety of 2D barcodes, select All 2D barcodes (autodetect type).
- NOTE:** If your lab only uses a specific type of 2D barcodes, you can decrease the decoding time by selecting that specific type. Phenolmager HT supports decoding of Data Matrix, PDF 417, QR Code, Aztec, and Maxicode.
- If your lab is using a combination of 1D and 2D barcodes, select **All barcodes (autodetect type)**. This option will increase scan time and file size, and should only be used if necessary.




Figure 35. Advanced Settings Window

Click **Save** to save the protocol. The study you previously chose will be automatically selected. You can change the protocol name or study if needed.




Figure 36. Save Protocol Window

## Brightfield Snapshots

Once you have selected your scan and multispectral field resolutions, you can take example snapshots of your slides.

Click **Take Snapshots** to load the Brightfield Snapshots editor.




Figure 37. Brightfield Snapshots

If a carrier is not already on stage, click **Load Carrier** to select a carrier that contains the slides you wish to image.




Figure 38. Select Carrier

Use the navigation tool to select an area on the slide. Click on the slide to change positions.




Figure 39. Navigation Tool

For easier navigation, click **Take Overview** to acquire imagery of all slides in your carrier. You can also select that option in the Preferences dialog in the Gear Menu. Phenolmager HT will then automatically acquire overview imagery whenever a carrier is loaded.

Click on the tissue to change positions.




Figure 40. Navigation tool with overview

You can refine the position by using arrow keys (for small movements) or control+arrow (for slightly larger movements).

Click **Autofocus** or use the Stage Height slider to bring the live view into focus.

You can click on a feature within the live view to center on that location or click and drag within the live view to change locations.

To take a snapshot of your current live view, select either the whole slide or the MSI regions radio button, pick a base file name, and click **Snap**.




Figure 41. Whole Slide and MSI Regions Radio Buttons

You can navigate to new places, focus as needed, and take as many snapshots as you want. Snapshots will be numbered incrementally.

Click **Back** when you are ready to return to the Edit Protocol window.

## Fluorescence Protocols

When creating a fluorescence protocol, you will need slides that include positive expression in all markers of interest in order to set suitable exposure times.

After you have created your fluorescence protocol and assigned it to a study, the Edit Protocol window (for fluorescence protocols) opens.




Figure 42. Fluorescence Protocol Editor

Under Overview Scan Rules, select the desired filter that will be used to help find tissue on the slide. This will typically be your DAPI counterstain. If there is no counterstain, then choose a filter that aligns with the expression of your most common stain or auto-fluorescence.




Figure 43. Overview Scan Rules

Under Slide Scan Settings, choose the Pixel Resolution to be used while imaging the slide (see Figure 43). Options include:

- 1.0 um per pixel (nominally 10x)
- 0.5 um per pixel, (nominally 20x)
- 0.5 um per pixel, (nominally 20x, binned from a 40x acquisition)
- 0.25 um per pixel (nominally 40x).




Figure 44. Slide Scan Settings

**NOTE:** Each time you change the pixel resolution, you will need to reset the exposure times.

Button	Description
<b>Auto-Update</b>	You can auto-update your exposures. This is recommended as a first step to give a close approximation of your exposure times. We then recommend you use the exposures editor to fine-tune your exposures.
<b>Manual Update</b>	You can manually update your exposures. In this case, the current exposures remain unchanged and you must use the exposures editor to obtain valid exposures.
<b>Cancel</b>	You can cancel. This will leave the pixel resolution unchanged.




Figure 45. Exposure Times Require Update Window

If you plan to take Multispectral Fields, choose the Pixel Resolution that you want to use for fields imaging. Available resolutions are 10x, 20x, 20x binned, or 40x.






Figure 46. Multispectral Field Settings

The Select Scan Bands button allows you to choose which bands you use to take your slide scan. Imagery for each band will be taken in the order it is shown. By default, Phenolmager HT will take a slide scan using the best selection available for your staining kit, or the bands you selected for a custom scan.

If you would like to remove bands where fluorescent signals are not expressing, uncheck the band. Drag and drop to change the order in which they are acquired. You can change the colors associated with each band; these are the colors that Phenochart will use to display the scan.

**NOTE:** If you remove the Overview or Focus band, Phenolmager HT will select new bands for you. You may want to visit Advanced Settings to confirm these selections.





Figure 47. Select Scan Bands

You can limit the filters displayed in the Select Scan Bands dialog.

- Show Akoya Scan Filters to limit the available filters to installed Akoya filters that are appropriate for slide scanning.
- Show All Akoya Filters to display all installed Akoya filters, both for scanning and taking multispectral fields.

- Show All Filters to display all installed filters, including custom filters that you may have installed.



*Figure 48. Select Scan Bands - Show Drop-down List*

The Select Field Bands button allows you to choose which bands you use to take your multispectral fields. We strongly recommend that you maintain the default band configuration. Non-signal bands are integral to multispectral stain extraction and unmixing.

If you would like to add or remove bands, select or de-select the filter. If there is more than one band available for a filter, you will need to select the band from the drop-down list. Drag and drop to change the order in which they are acquired.

**NOTE:** If you remove the Focus band, Phenolmager HT will select a new band for you. You may want to visit Advanced Settings to confirm this selection.




Figure 49. Select Field Bands - Band Selection

You can limit the filters displayed in the Select Field Bands dialog.

- Show Akoya MSI Filters to limit the available filters to installed Akoya filters that are appropriate for multispectral fields.
- Show All Akoya Filters to display all installed Akoya filters, both for scanning and taking multispectral fields.
- Show All Filters to display all installed filters, including custom filters that you may have installed.




Figure 50. Select Field Bands - MSI Filters

## Advanced Settings

There may be situations when you need to fine-tune how Phenolmager HT scans the tissue. Click the **Advanced Settings** button to adjust the following functions:

- Phenolmager HT is designed for use with standard #1.5 coverslips, where the distance between the tissue and the top of the sample is approximately 170 microns. While this is the optimal setup, Phenolmager HT can handle other scenarios.
  - If your lab uses thinner coverslips (#0), select “100 um (#0 or similar)” for a Coverslip Thickness.
  - If your lab uses thicker coverslips, or has a significant amount of mounting media between the sample and coverslip, select the option that best matches your slides.

**NOTE:** *Phenolmager HT is capable of handling slides with up to 420 microns of combined coverslip or mounting medium above the sample, though imagery at these extremes may be sub-optimal.*

- You can specify which filter will be used for setting focus during whole slide scans and for multispectral field imaging. Choose the filters where your counterstain expresses. This is normally DAPI.
- Saturation protection prevents you from overexposure when acquiring multispectral field imagery. We highly discourage deactivating the Use Saturation Protection setting.
- If Phenolmager HT is having difficulty finding your entire tissue sample, check the 'Scan within the entire coverslip region' checkbox. This will increase scan time and file size, but will enable you to complete scanning of difficult samples.

**NOTE:** *If there is Sharpie<sup>®</sup> markup on your fluorescent slide, Phenolmager HT will use that as the region instead of the coverslip.*

- Phenolmager HT has a specialized algorithm for finding TMA samples. If your samples are TMAs, select Sample is a TMA.
- **Label Barcode Reading** If you would like Phenolmager HT to decode the barcode on your slide label and save it as part of the scan file, select the type of barcode you are using. By default, Phenolmager HT will not attempt to decode a barcode unless you specifically select this option.
  - If you have only 1D barcodes, select All 1D barcodes (autodetect type).
  - If your lab uses a variety of 2D barcodes, select All 2D barcodes (autodetect type).

**NOTE:** If your lab only uses a specific type of 2D barcodes, you can decrease the decoding time by selecting that specific type. Phenolmager HT supports decoding of Data Matrix, PDF 417, QR Code, Aztec, and Maxicode.

- If your lab is using a combination of 1D and 2D barcodes, select All barcodes (autodetect type). This option will increase scan time and file size, and should only be used if necessary.

Advanced Settings

Coverlip Thickness  
170 µm (#1 or #1.5) ▾

Autofocus  
DAPI ▾ For Slide Scans  
DAPI MSI ▾ For Fields

Saturation Protection  
☒ Use Saturation Protection For Fields

Scan Entire Coverslip Region  
☐ Scan within the entire coverslip region

Sample Finding Options  
☐ Sample is a TMA

Label Barcode Reading  
No barcodes ▾

OK Cancel

Figure 51. Fluorescence - Advanced Settings

After you have selected your settings, use the Exposures Editor to [set exposures](#)<sup>[69]</sup> for your protocol.

### 6.3.4 Setting Exposures

After you have selected your scan and multispectral field resolutions, filters, and bands, you need to set exposures. Click **Set Scan Exposures** or **Set Field Exposures** in the Fluorescence Protocol Editor to set exposure times for scanning and field acquisition.




Figure 52. Edit Exposures

If a carrier is not already on stage, click **Load Carrier** to select a carrier that contains representative slides. Select your carrier from the list and load it to the stage.




Figure 53. Select Carrier

Use the navigation tool to select a slide. You will see a live view of your first filter for your imaging rules. In most cases, this will be DAPI. The highlighted entry in the table identifies which filter or band is currently shown.




Figure 54. Navigation tool

For easier navigation, click **Take Overview** to acquire digitally enhanced imagery of all slides in your carrier. You can also select that option in the Preferences dialog in the Gear Menu. Phenolmager HT will then automatically acquire overview imagery whenever a carrier is loaded.




Figure 55. Navigation tool with overview

Click on the tissue to change positions.

You can refine the position by using arrow keys (for small movements) or Ctrl+arrow (for slightly larger movements).

You can click on a feature within the live view to center on that location, or click and drag within the live view to change locations.

Click **Autofocus** or use the **Stage Height slider** to bring the live view into focus. Click **Autoexpose** to have the system find the best exposure for that filter/band.

- After auto-exposing, you may want to auto-focus and auto-expose again to refine your focus and exposure estimates.
- You can override the auto-exposure value by typing a value in the highlighted cell. Values must be between 0.1 and 2000 ms.
- If you would like to turn off the fluorescence illumination and live view, click **Close Shutter**. You will need to re-open the shutter to see a live view and set exposures.

If you plan to use this protocol for multispectral field acquisition, you should set the multispectral field exposures at this time. In most cases, you can use the same area of the slide to set your multispectral field exposures. Select the **Set Field Exposures** radio button to set exposures for multispectral fields. Focus and click **Autoexpose**.

- As with whole slide scan exposures, after auto-exposing, you may want to auto-focus and auto-expose again to refine your focus and exposure estimates.
- You can override the auto-exposure value by typing a value in the highlighted cell. Values must be between 0.1 and 2000 ms.
- If you would like to turn off the fluorescence illumination and live view, click **Close Shutter**. You will need to re-open the shutter to see a live view and set exposures.

Repeat the steps above for all filters and bands in your protocol. You may need to change locations and/or slides to find the best signals for setting your exposures.



## Filter Cube Recommendations when using Akoya Biosciences Opal™ Reagents

If you are using your Phenolmager HT system with Akoya's Opal™ multiplex IHC detection kits or with individual Opal fluorophores, please refer to Table 11 and Table 12 for exposure time recommendations for acquiring spectral library images.

When acquiring multispectral slide scans with dyes in Opal 4-color / 3-plex or Opal 7-color / 6-plex kits:

Table 11. Multispectral Slide Scan Band								
	DAPI	Opal 480	Opal 520	Opal 570	Opal 620	Opal 690	Opal 780	Sample AF
DAPI								
Opal 480								
Opal 520								
Opal 570								
Opal 620								
Opal 690								
Opal 780								

When acquiring multispectral fields:

Table 12. Multispectral Field Band							
	DAPI	Opal 480	FITC	Cy3	Texas Red	Cy5	Opal 780
DAPI							
Opal 480							
Opal 520							
Opal 540							
Opal 570							
Opal 620							
Opal 650							

Opal 670							
Opal 690							
Opal 780							

**Key:**

- Colors are signal bands: autoexpose or manually set an appropriate exposure per sample.
- Dark Gray is an autofluorescence band: autoexpose or manually set an appropriate exposure per sample.
- Light Gray is a non-signal band: manually set exposure according to Table 13.

Table 13. Exposure Times for Non-signal Bands				
	Resolution			
	10x	20x	20x (40x binned)	40x
Multispectral Slide Scans	6 ms	25 ms	2.5 ms	8 ms
Multispectral Fields	35ms	150 ms	15 ms	50 ms

## Snapshots

When you have set your exposures, you can take snapshots to see sample imagery of how your exposures perform on a given area of the slide. Select either the **Set Scan Exposures** or the **Set Field Exposures** radio button, pick a base file name, and click the **Snap** button. You can navigate to new places and take as many snapshots as you want. They will be numbered incrementally.

Figure 56. FL Exposures Panel

Click the **Back** button when you are ready to return to the Edit Protocol window.

**NOTE:** The exposures you just set are only valid for the resolutions you picked. If you change your scan or multispectral field resolution, you will need to revisit the exposures editor to update your exposures.

## 6.4 Scanning Slides

The section describes the procedures involved when scanning slides:

- [The Carrier](#)<sup>76</sup>
- [Carrier Status](#)<sup>76</sup>
- [Slide Status](#)<sup>78</sup>

- [Setting Up Scan Rules](#) <sup>78</sup>

### 6.4.1 The Carrier

Slides are loaded into Phenolmager HT via carriers, which are stacked into the Slide Carrier Hotel.



*Figure 57. Slide Carrier*

The hotel can hold up to 20 carriers, and each carrier can hold up to 4 slides. This means Phenolmager HT can be loaded with up to 80 slides at any given time. Those 80 slides can all be scanned based on the same set of rules, or you can tailor your scan rules to each slide and/or each carrier as needed.

### 6.4.2 Carrier Status

Within the software, each hotel slot is identified as 'Slot N,' where N is 1-20 running top to bottom. When a slide carrier is loaded into the hotel, an LED light changes color to identify the status of that carrier.




Figure 58. Scan Slides - Carrier Status

The User Interface will reflect the current state of each slot in the hotel. If there is no carrier present, you will see:







Each carrier present in the hotel will be represented as follows:







Each carrier is represented by a carrier status icon, four slide status icons, and four slide IDs.

The carrier status icon corresponds to the hotel LED for that particular carrier. Each possible status is represented by a different color.

When you insert the carrier into the hotel, a white icon  indicates that the carrier is present but has no rules to use for slide scanning. If the icon is blinking red , the carrier has not been completely inserted. Try pushing it in further.

Once the carrier has been inserted in the hotel, you can define rules for your scan. The icon will change to blue  when there is at least one slide ready to be scanned. Carriers that have been prioritized will be blinking blue .

If no slides on the carrier can be scanned with the current rules, the icon will change to yellow . This carrier will be ignored until the problem is resolved.

When scanning has started, the icon for the carrier on stage will be blinking green . When the carrier is returned to the hotel, the icon will either be green  (for a successful scan) or orange  (if there was a problem with the scan).

If the icon is solid red , that hotel slot's sensor is malfunctioning. Do not use this slot.

### 6.4.3 Slide Status

Each carrier also has a status icon for the 4 slides it contains.

➡ means the slide has complete rules and is ready to be scanned.

🟢 (blinking) means the slide is currently being scanned.

✅ means the slide was scanned successfully.

Other options include ⚠️, which means the slide has rules, but they cannot be used. (They may be incomplete, in which case you will need to edit the rules for that slide). 🚫 indicates that the slide failed to scan correctly.

Finally, ⚪ will indicate that the slide is being ignored. This is useful if the carrier has fewer than four slides or if you don't want to process that slide at this time.

### 6.4.4 Setting Up Scan Rules

Scan rules must be defined for each slide you would like to process. Scan rules are complete when a study, protocol, task, and slide ID have been assigned.

This section contains:

- [Slide IDs](#) <sup>78</sup>
- [Editing a Single Slide](#) <sup>79</sup>
- [Editing All Slides Within a Single Carrier](#) <sup>80</sup>
- [Editing Multiple Carriers](#) <sup>80</sup>
- [Saving and Restoring Scan Setup](#) <sup>83</sup>
- [Scanning](#) <sup>85</sup>

#### Slide IDs

You must manually enter the ID for each slide in the carrier. You can enter IDs before the carrier is placed into the hotel, or you can enter IDs after the carrier has been placed into the hotel.

If you would like to enter the IDs while the carrier is still in your hand, click **Enter Slide IDs**.





Figure 59. Enter Slide IDs Button

It will bring up this dialog:



Enter Slide IDs

Please enter Slide IDs for the incoming carrier:

Slide 1:

Slide 2:

Slide 3:

Slide 4:

OK Cancel


Figure 60. Enter Slide IDs Window

You can then enter between 1-4 slide IDs for your slides and click **OK**. Immediately put the carrier in the hotel, and PhenolMager HT will assign the slide IDs you entered to that carrier.

If your carrier is already present in the hotel, you can edit the carrier to add the slide IDs along with the other scanning rules.

## Editing a Single Slide

To edit the scanning rules for a single slide, click the **status icon** for that slide. It will bring up this editor:



Edit Slides

Task:

Study:

Protocol:

Slide ID:

☐ Prioritize the containing carrier

OK Cancel

Figure 61. Edit Slides

For this slide, enter:

- **Task** - Choose **Scan Slides** or **Acquire Fields**.

**NOTE:** Select Ignore if you do not wish to process this particular slide.

- **Study** - Select a study.
- **Protocol** - Select any protocol you have saved within the study.
- **Slide ID** - Add your slide ID, if necessary.

Check Prioritize the containing carrier if you want this slide's carrier to be scanned at the front of the line.

## Editing All Slides Within a Single Carrier

To edit the scanning rules for a single carrier, click the **status icon** for that carrier. It will bring up this editor:




Figure 62. Edit Slides

For each slide, you can enter:

- **Task** - Choose **Scan Slides**, **Acquire Fields** or **Ignore**. Select **Ignore** if there is no slide in the slot.
- **Study** - Select a study.
- **Protocol** - Select any protocol you have saved within the study.
- **Slide ID** - Add your slide ID, if necessary.

Check Prioritize this carrier if you want this carrier to be scanned at the front of the line.

## Editing Multiple Carriers

To edit multiple carriers, Click **Configure Tasks**.






Figure 63. Configure Tasks Button

Select the carriers you want to edit. Multiple selection is available.




Figure 64. Configure Carriers

You will also need to select how you want to process your slides.

To process all slides using the same rules, choose Using the same rules for all slides and click **Next...** The Edit Slides window opens.




Figure 65. Edit Slides

- Select your task, study, and protocol.
- Check **Prioritize the containing carriers** if you want the selected carriers to be scanned at the front of the line.

**NOTE:** You will not be able to edit the slide ID, because you are applying rules across multiple slides. If you need to enter slide IDs, you can edit the individual carriers.

To use different processing rules for each carrier position, choose **Using different rules for each carrier slot** and click **Next...** For each position in the carrier, select the task, study, and protocol. Prioritize if desired. Again, you will not be able to edit the slide ID, because you are applying rules across multiple slides.




Figure 66. Edit Slides

## Saving and Restoring Scan Setup

After setting up all the details for your slides, you may want to save your work for later use. This may be useful if the scanning process is interrupted before the task is finished or if you need to return later to re-scan certain slides.

You can save your slide IDs, studies, protocols, and tasks by clicking the **Save Setup** button. Phenolmager HT will also prompt you to save before scanning (when you click **Scan**). At this point, you can enter a name and some notes about this scan setup.




Figure 67. Scan Slides

Saving a setup is not required, but it can help you locate the information later if needed. By default, this information will be deleted after 30 days, but you can override the default and opt to keep the information until you manually delete it.




Figure 68. Save Scan Setup

To restore an earlier scan setup, click **Load Setup**.

All available setups will be shown. You can filter by date to show setups only from a specific time period. Click the headers to sort by date, name, notes, studies, etc.

**NOTE:** You can also delete entries that are no longer needed by selecting those items and clicking **Delete**.

To use a setup, select it from the dialog and click **Load**.

Load Scan Setup

Display for the last: ☐ Day ☒ Week ☐ Month ☐ All

Date	Name	Notes	Studies	#Carriers	Permanent
6/15/2018 5:28:38 PM	PD-L1	520 titrations	180614 PDL1 520 titration and 180614 CD8 480 titrations	1	No
6/15/2018 5:27:12 PM	PDL1, CD8 titrations	CC, Opal titrations, tonel	180614 PDL1 520 titration and 180614 CD8 480 titrations	2	No

Delete Load Cancel

Figure 69. Load Scan Setup

After you have selected a setup, you need to configure the carriers in Phenolmager HT to match the information in the setup.

If the carriers and slides are the same location as when you saved the setup, you only need to click **OK**, because everything is in the right place.

In other cases (for example: when you need to put your slides back in), this dialog will help you figure out where your slides belong.

1. When scan information is available and a carrier is present in the hotel, the system will show the slide IDs (slot 1 in Figure 70).
2. When scan information is available but there is no carrier present in the hotel slot, the slot will be labeled as Missing in red (see slot 5 in Figure 71). The slide IDs for that carrier are listed, so you can locate the correct slides, place them in the carrier, and insert the carrier into the hotel.
3. When no information is available for a carrier that is present in the hotel, the system will show italicized Slide N text for the slide IDs (slots 7, 9, etc). Slides in these carriers were not present in the original setup.

For every item that the scan setup had - meaning, anywhere the system has information on a carrier that should be present (items 1 and 2) - you need to specify what to do before the carrier can be loaded. There are 3 options:

1. **Use as saved:** The system will load all slide information including the scan task.
2. **Ignore task:** The system will load all slide information except for the scan task. The scan task will be set to ignore. This is useful if you are loading the setup so you can rescan a few select items, or if you know that you're going to be changing the task from how it was originally saved.
3. **Don't restore:** No information for this carrier will be reloaded. This is useful if you do not have or do not need the carrier.




Figure 70. Slide Setup


The red text at the bottom of the screen will list any issues that need to be fixed before scanning.

Click **OK** when you have finalized your selections. The scan slides editor will be cleared and replaced with your selected scan information.

## Scanning

The Scan button is enabled when at least one carrier is ready to be scanned. Click **Scan** to start scanning. If any carriers have been prioritized, these carriers will be processed first.

Phenolmager HT will report scanning progress for each slide on the progress dialog and on the front panel LEDs.



*Figure 71. Processing Slides Window*

Click **Pause Scanning** to remove completed carriers and add new carriers to be scanned. When you pause a scan, Phenolmager HT will return the current carrier to the hotel. It will then be safe to open the door. You can remove carriers that are completed and replace them with unscanned carriers. Edit those carriers to set them up for scanning, and reprioritize carriers if needed. When you're ready, click **Scan** to resume scanning.

# 7 Maintenance

This section includes procedures that are to be performed periodically, either to make the Phenolmager HT run better, or to protect its components from damage.

- [Cleaning the Instrument Exterior](#)<sup>[87]</sup>
- [Cleaning the Monitor](#)<sup>[87]</sup>
- [Cleaning the Power and Communication Ports](#)<sup>[88]</sup>
- [Cleaning the Slide Carriers](#)<sup>[88]</sup>
- [Replacing the Fuses](#)<sup>[88]</sup>

## 7.1 Cleaning the Instrument Exterior

Clean the Phenolmager HT instrument exterior as necessary.

The non-electrical exterior parts of the Phenolmager HT can be wiped down with a soft cloth using standard laboratory grade cleaning solutions including:

- 70% ethanol
- 10% bleach
- Clidox®
- Sporidicin®

**WARNING:**



DO NOT spray cleaning solutions directly onto the Phenolmager HT instrument. Sprays and liquids that come into contact with the Phenolmager HT instrument may result in damage to the system or electric shock. Always spray the cleaning solution onto a cloth and then wipe the instrument exterior with the cloth.

## 7.2 Cleaning the Monitor

Clean the monitor as necessary with a soft, lint-free cloth. If needed, dampen the cloth with water or an eyeglass cleaner.

**WARNING:**



DO NOT spray cleaning solutions directly onto the monitor screen. Sprays and liquids that come into contact with the screen may result in damage to the system or electric shock. Always spray the cleaning solution onto a cloth and then wipe the screens with the cloth.

## 7.3 Cleaning the Power and Communication Ports

Dust can be removed from the power and communication ports with a gentle air stream from compressed air cans. DO NOT shake the can before use. Hold can at a distance from the ports to prevent condensation.

### WARNING:



Turn off the electrical power to the Phenolmager HT system by shutting down the system (see [System Shutdown](#)<sup>[39]</sup>) before cleaning any part of the instrument where electrical or fiber optic cables make connections.

## 7.4 Cleaning the Slide Carriers

It is important to keep carriers free from sticky debris caused by slide labels, tape, or mounting media.

Slide carriers can be cleaned using soap and water with a soft cloth, or using an ultrasonic bath.

## 7.5 Replacing the Fuses

For blown fuses and to order replacement fuses, contact Akoya Biosciences Technical Support (see [Contact Us](#)<sup>[64]</sup>).

### WARNING:



Electric Shock Hazard. DO NOT attempt to replace the fuses. Only qualified Akoya Biosciences or PerkinElmer service personnel can replace the fuses.



## 8 Appendix A: Akoya TIFF Specification

### Background

This describes a TIFF format that Akoya Biosciences uses for its tissue images. The imagery may be a simple RGB image, a set of components extracted from a field or region, or a whole-slide scan. In the latter case, it may be a brightfield (BF) color RGB image or a multiband fluorescence (FL) image.

The goal is to use the same syntax and metadata for all these kinds of images, and minimize the semantic distinctions where possible. Specifically, an extracted component represents signal attributed to a stain or fluorescent dye in a sample. The signal values incorporate image preprocessing such as normalization for exposure time and spectral unmixing when that technique is used.

### Data format

The files are TIFF or BigTIFF images, depending on image size, with multiple images per file.

For images larger than about 2K x 2K pixels, tiled format is used, and the image is provided in several resolutions (pyramidal tiled images). Tile size is 512 x 512 pixels. Images smaller than 2K x 2K use stripped format.

The highest resolution (baseline) image(s) appear first in the file. For each resolution there are N baseline images where N depends on the contents. For BF images, N=1 and each image is an RGB image. For FL images or unmixed component images, N = number of bands, which is usually > 1, and each image is a grayscale image.

A thumbnail RGB image is provided, and this is a good image to use as an icon in graphical image lists. This comes after the baseline images, meaning it is the second image in BF (RGB) images, and the (N+1)st image for FL images or unmixed component images.

Next come the reduced-resolution images (if present). The pyramid contains enough levels that the image size is no larger than 2K x 2K at the coarsest resolution.

For whole-slide scans, there are two more non-tiled images after these: an optional RGB image of the label, and a macro (low-resolution) RGB image of the whole slide.

Overall, the arrangement is:

**Table 1. Images**

Description	RGB/mono	Title/Strip	Resolution	Notes
Baseline image	Varies	Varies	Full	Tiled if > 2K x 2K RGB for BF, else mono
More full-resolution images	Mono	Varies	Full	If N > 1
Thumbnail	RGB	Stripped	~500 x 500	
Half-resolution images	Varies	Varies	Half	Only if baseline is tiled

Quarter, eighth, etc.	Varies	Varies	Quarter, eighth, etc.	Continues until 2K x 2K or smaller
Macro (overview) image of whole slide	RGB	Stripped	~2000 x 4000	Required for whole-slide scans Optional for simple RGB images and extracted components
Label image	RGB	Stripped	~500 x 500	Optional, whole-slide scans

## Detection

Readers can recognize Akoya tissue images via the contents of the “Software” TIFF tag (see below). The file suffix is .qptiff for whole slide scans. Some TIFF tags contain “PerkinElmer” for historical reasons; they remain to retain 3rd party reader compatibility.

## Metadata

Metadata is contained in two locations: within standard TIFF tags as listed in the table below, and within the ImageDescription string, using a set of XML tags described below. These are provided for each image (IFD) in the file, and describe that image rather than the baseline image or the scan as a whole. The ScanProfile tag is only provided on the first, baseline image as it may be large.

**Table 2. TIFF tags**

TIFF Tag	Optional?	Description of contents
Software		Starts with “PerkinElmer-QPI”
ImageDescription		Further metadata in XML format (see next section)
ImageWidth		Width of the image in pixels
ImageLength		Height of the image in pixels
ResolutionUnit		Unit used for resolution and position (see below)
XResolution		Pixel X resolution (see below)
YResolution		Pixel Y resolution (see below)
XPosition	Y	Sample X location in ResolutionUnits. This is ULHC location except for Macro image which reports its image center.
YPosition	Y	Sample Y location in ResolutionUnits. This is ULHC location except for Macro image which reports its image center.
SampleFormat		Integer (1) for BF, FL; or float (3) for extracted components
SMinSampleValue		Minimum signal value in the image
SMaxSampleValue		Maximum signal value in the image
BitsPerSample		8 (FL); 8, 8, 8 (RGB); or 32 (unmixed component)
SamplesPerPixel		1 (FL or unmixed component) or 3 (RGB)
NewSubfileType		0 for full-resolution images, 1 for reduced res images
TileWidth	Y	Tile width (512) if tiled format is used

TileLength	Y	Tile height (512) if tiled format is used
TileOffsets	Y	List of tile offsets, if tiled format is used
TileByteCounts	Y	Size of each (compressed) tile, if tiled format is used
StripOffsets	Y	List of strip offsets, if tiled format is not used
RowsPerStrip	Y	Number of rows per strip, if tiled format is not used
StripByteCounts	Y	Size of each (compressed) strip, if tiled format is not used
PlanarConfiguration		1 (chunky) for RGB images, 2 (planar) otherwise
PhotometricInterpretation		2 (RGB) for RGB images, 1 (BlackIsZero) otherwise
DateTime		Acquisition time
Compression		May be None, CCITT Group 3, PackBits, LZW, or JPEG
JPEG fields	Y	JPEG fields are defined when JPEG compression is used

ResolutionUnit, XResolution and YResolution are required fields in a valid TIFF file. When the true resolution of the image is known, ResolutionUnit will be 3 (cm) and XResolution and YResolution will be pixels/cm. When the true resolution is not known, ResolutionUnit will be 2 (inch) and XResolution and YResolution will be 96 (pixels/inch). Pixels from Akoya Biosciences instruments are always square so XResolution and YResolution will always have the same value.

The TIFF spec is not explicit about the data type and value for SMinSampleValue and SMaxSampleValue; the writer uses the same data type as the image pixels (byte or float).

### Image Description contents

The ImageDescription tag contains a string in XML format. The string contains a top-level <PerkinElmer-QPI-ImageDescription> element. Nested within this element are child elements with the tag names and values as listed in the table below. Elements appear in the order listed. Values are stored as text content of the element. Elements are required unless otherwise specified. See the example below.

**Table 3. Image Description tags**

Tag	Optional?	Contents
DescriptionVersion		Version of the image description field itself, a single number. This document describes version 2 of the field.
AcquisitionSoftware		Software used to acquire the image
Identifier		GUID in string format. This is an identifier for the image file itself.
SlideID	Y	ID of the slide that this image was taken from.
Barcode	Y	Barcode text of the slide this image was taken from.
ComputerName	Y	Name of the computer on which the slide was scanned.
ImageType		A string identifying the type of image within the file (Table 1), with the following values: <ul style="list-style-type: none"> <li>FullResolution</li> <li>ReducedResolution</li> <li>Thumbnail</li> </ul>

		<ul style="list-style-type: none"> <li>• Overview</li> <li>• Label</li> </ul>
IsUnmixedComponent		"True" for unmixed multispectral images, otherwise "False".
ExposureTime		Exposure time as an integer number of microseconds. For unmixed images, this is the exposure time for the dominant wavelength band for the component (FL); or the brightest wavelength in the cube (BF).
SignalUnits		<p>A byte <i>www tttt</i> where the <i>tttt</i> nibble indicates the signal unit type from the following:</p> <ul style="list-style-type: none"> <li>• – raw counts</li> <li>• – normalized (counts/second/gain/full-scale/binning)</li> <li>• – OD (optical density)</li> <li>• – dark-corrected counts</li> </ul> <p>and the <i>www</i> nibble indicates how the signal is weighted across the spectral bands (or colors):</p> <ul style="list-style-type: none"> <li>• – average across all bands</li> <li>• – total summed signal across all bands</li> <li>• – peak signal in highest-valued band</li> </ul> <p>Thus, for example, a value of 68 (hex 44) encodes OD units with peak-signal weighting.</p>
Name	Y	Band name for FL whole slide scans, component name for extracted components. Not present for RGB images
Color	Y	Color to use when rendering this band, as decimal r,g,b byte triplet. Present for FL whole slide scans or extracted components. Not present for RGB images.
Responsivity	Y	Instrument responsivity, if available, for FL whole-slide and unmixed images. See below for details.
Objective	Y	Objective name, if known, otherwise not present.
ScanProfile	Y	Element containing scan and/or and unmix parameters. It is valid XML whose contents are opaque to most readers. It is only provided on the first (baseline) image, and is omitted from all other IFDs.
ProtocolName	Y	Name of the acquisition protocol.
ValidationCode		Used for internal data integrity checks – readers can ignore this.
ImagingCycle	Y	Zero-based acquisition cycle for CODEX imagery
Biomarker	Y	String indicating target such as CD8, Ki67, etc.
ScaleFactor	Y	After allowing for bit-depth conversion, pixel values in this dataset have been scaled by this factor. For example, if a 12-bit image is stored in an 8-bit image, and counts were scaled from 4095 to 255, this would be 1; if counts were scaled from 3102 to 255, ScaleFactor would be 1.3201. Default is 1.
AutofluorescenceSubtracted	Y	String containing 'true' or 'false', indicating whether this dataset has been subtraction-corrected by a background (autofluorescence) image. Treated as 'false' by default.
LampType	Y	<p>String with the FL lamp type</p> <ul style="list-style-type: none"> <li>• Unknown (default)</li> <li>• PriorLumen200</li> <li>• PriorLumen220</li> <li>• XCiteMultiBandLed</li> </ul>

		<ul style="list-style-type: none"> <li>• PriorMimicLed</li> <li>• SolaLed</li> <li>• Colibri7</li> </ul>
InstrumentType	Y	String with the instrument type, one of the following: <ul style="list-style-type: none"> <li>• Unknown (default)</li> <li>• Keyence700</li> <li>• Keyence800</li> <li>• AxioObserver</li> <li>• Dmi8</li> <li>• Polaris</li> <li>• Vectra3</li> </ul>
CameraType	Y	String with the camera sensor type, one of the following: <ul style="list-style-type: none"> <li>• SonyICX285AL</li> <li>• HamamatsuC11440</li> <li>• IMX265</li> <li>• IMX421</li> </ul>
ExcitationFilter	Y	Description of excitation filter with sub-elements listed below If present, all its sub-elements must be provided.
Name		String with the filter public name
Manufacturer		String with the manufacturer or OEM supplier
PartNo		String with the manufacturer part number
Bands		Collection of bands each as described by sub-elements below
Band		Description of a band with sub-elements noted below
Name		String with the band name
Cuton		Cuton wavelength of this band in nm
Cutoff		Cutoff wavelength of this band in nm
Active		True if this band is excited, False if not
EmissionFilter	Y	Description of emission filter with sub-elements listed below. If present, all its sub-elements must be provided.
Name		String with the filter public name
Manufacturer		String with the manufacturer or OEM supplier
PartNo		String with the manufacturer part number
Bands		Collection of bands each as described by sub-elements below
Band		Description of a band with sub-elements noted below
Name		String with the band name
Cuton		Cuton wavelength of this band in nm
Cutoff		Cutoff wavelength of this band in nm
CameraSettings	Y	Description of the camera settings used for this image with sub-elements listed below. If present, all its sub-elements must be provided
Gain		Gain factor, such as 1, 3, etc.
Binning		Binning extent, such as 1 for full-resolution, 2 for 2x2 binning
OffsetCounts		Offset count setting for the camera
BitDepth		Bit resolution used to read out the camera. This may differ from the resolution used to store this image (for example, it

		may be acquired with 12-bit resolution but saved in a QTIFF file with 8-bit resolution)
Orientation		String containing one of the following values: <ul style="list-style-type: none"> <li>• 'Normal' (default)</li> <li>• 'Rotate180'</li> <li>• 'MirrorV'</li> <li>• 'MirrorH'</li> </ul>
ROI	Y	Description of the sensor ROI contributing to the image. Coordinates are sensor pixels before binning. Even if the CameraSettings are provided, this field is optional
X		X origin of ROI. Left-most pixel in sensor array is 0.
Y		Y origin of ROI. Top-most pixel in sensor array is 0.
Width		Width of ROI
Height		Height of ROI

For whole slide images (BF and FL), SignalUnits will be 64 (hex 40) (raw counts, peak signal). For unmixed images, SignalUnits will reflect the unmix settings.

## Instrument Responsivity

The <Responsivity> tag is a container for a list of normalized instrument response values. This tag is present for whole-slide FL images from Phenolmager HT and unmixed FL multispectral images originating from Vectra 3 and Phenolmager HT.

For whole-slide images, the <Responsivity> tag will contain one <Filter> tag. The <Filter> tag contents will be different for each image within the TIFF file, reflecting the filter used to take the image.

For unmixed component images, the <Responsivity> tag will contain one <Band> tag for each band in the original image file. The <Band> tags are repeated for each unmixed component image.

The overview, thumbnail and label images do not have <Responsivity> tags.

Each <Filter> or <Band> tag describes the instrument responsivity for acquisitions using that filter or band. The contents of the <Filter> and <Band> tags are described below.

**Table 4: Contents of Filter and Band tags**

Tag	Contents
Name	The name of the filter (whole-slide image) or band (component image)
Response	The instrument response to the reference artifact, normalized for exposure. This is $\text{raw counts} / (2^{\text{bit depth}} \times \text{exposure time} \times \text{gain} \times \text{binning area})$ , where bit depth is the bit depth of the imagery, exposure time is in seconds, gain is the gain setting of the camera, and binning area is 1 for 1x1, 4 for 2x2, etc.
Date	The date and time of the reference image in UTC, ISO 8601 format.
FilterID	Detailed description of the acquisition filter, as {ExcitationFilter band name}_ {ExcitationFilter manufacturer}:{ExcitationFilter part number} / {EmissionFilter manufacturer}:{EmissionFilter part number}

	The {ExcitationFilter band name}_ prefix may be omitted if there is only a single excitation band.
--	--

## Sample ImageDescription

Sample ImageDescription for the DAPI band of a FL whole-slide image, containing a single <Filter> tag:

```
<?xml version="1.0" encoding="utf-8"?>
< PerkinElmer-QPI-ImageDescription >
  <DescriptionVersion>1</DescriptionVersion>
  <AcquisitionSoftware>VectraScan 1.0.0</AcquisitionSoftware>
  <ImageType>FullResolution</ImageType>
  <Identifier>AABED946-BB58-44FB-95B3-48E177E3BB83</Identifier>
  <IsUnmixedComponent>False</IsUnmixedComponent>
  <ExposureTime>50</ExposureTime>
  <SignalUnits>64</SignalUnits>
  <Name>DAPI</Name>
  <Color>0,0,255</Color>
  <Responsivity>
    <Filter>
      <Name>DAPI</Name>
      <Response>30.7</Response>
      <Date>2015-10-22T13:10:18.0618849Z</Date>
      <FilterID>Semrock:FF02-409/LP-25 Emission / Semrock:FF01-377/50-25 Excitation</FilterID>
    </Filter>
  </Responsivity>
  <Objective>4x</Objective>
  <ScanProfile><!-- this will be a serialized scan protocol. It is valid XML but otherwise opaque -
  </ScanProfile>
  <ValidationCode>4281ff86778db65892c05151d5de738d</ValidationCode>
</ PerkinElmer-QPI-ImageDescription >
```

## 9 Appendix B: EULA

The following is an agreement (the "Agreement") between you and Akoya Biosciences, Inc., 100 Campus Drive, 6th Floor, Marlborough, MA 01752 ("Akoya") for software known as Research Use Only Phenolmager HT (formerly Vectra Polaris), including any updates, upgrades or enhancement provided by Akoya and any and all accompanying documentation (collectively, the "Software").

**AKOYA IS WILLING TO LICENSE THE SOFTWARE TO YOU ONLY UPON THE CONDITION THAT YOU ACCEPT ALL OF THE TERMS CONTAINED IN THIS AGREEMENT. PLEASE READ THE TERMS CAREFULLY BEFORE USING THE SOFTWARE, AS INSTALLING AND/OR USING THE SOFTWARE WILL INDICATE YOUR AGREEMENT TO THEM. IF YOU DO NOT AGREE TO ALL OF THE TERMS AND CONDITIONS IN THIS AGREEMENT, AKOYA IS NOT WILLING TO LICENSE THE SOFTWARE TO YOU, AND YOU (1) MAY NOT USE THE SOFTWARE AND (2) MUST CONTACT AKOYA WITHIN THREE (3) DAYS OF RECEIPT TO ARRANGE FOR THE RETURN OF THE SOFTWARE. BY INSTALLING AND/OR USING THE SOFTWARE, YOU ALSO REPRESENT AND WARRANT THAT YOU HAVE THE AUTHORITY TO AGREE TO AND BE BOUND BY THE TERMS OF THE AGREEMENT.**

- 1. Single Use License.** The Software is licensed to you and not sold, and you acknowledge that there is no transfer of title or ownership of the Software to you. Subject to the terms and conditions of this Agreement and your payment of the license fees for the Software, Akoya hereby grants to you a restricted, non-exclusive, non-transferable, non-assignable, non-sublicensable and revocable single-use license to use, for your internal research purposes only, the executable code version of the Software and the accompanying documentation in hard copy or electronic format. Akoya grants no other licenses hereunder, except for the Research Use Only Phenolmager HT (formerly Vectra Polaris) Software set forth above. Nothing in this Agreement shall be construed as a general internal use license. For the avoidance of doubt, no right to use Phenolmager HT (formerly Vectra Polaris) in Diagnostic Applications is granted under this Agreement. Akoya RESERVES ALL RIGHTS NOT EXPRESSLY GRANTED BY THIS AGREEMENT.
- 2. Specific Restrictions.** You may use the Software only on a single computer. You may make only one (1) copy of the Software, solely for backup purposes, and such copy must reproduce and include all copyright notices indicating Akoya's ownership of the Software. Except as expressly permitted by applicable law or Akoya, you agree that neither you nor a third party acting on your behalf, directly or indirectly, shall: (i) decompile, disassemble or reverse engineer the Software or otherwise attempt to derive the source code of the Software or underlying ideas, techniques or algorithms, including the review of data structures or similar materials produced by the Software; (ii) modify, alter, translate, or create derivative works of the Software; (iii) lease, rent, loan, commercially share or otherwise use the Software in any manner to provide service bureau, commercial time-sharing or other computer services to third parties; (iv) transmit the Software or provide its functionality, in whole or in part, over the Internet or other network; (v) use the Software for any purpose other than solely for your internal research; (vi) use the Software to provide any analytics or diagnostics or otherwise for the benefit of any third party; (vii) sell, distribute, publish, rent, lease, sublicense, assign possession, allow use of or otherwise transfer the Software or any copy of the Software to a third party; (viii) create any computer software program or documentation which is similar to or competitive with the Software; (ix) encumber, disclose or otherwise make available the Software to any third-party; or (x) remove, modify or obscure any copyright, trademark, or other proprietary rights notice in or on the Software. This Agreement may not be assigned in whole or in part, except that either party may assign this Agreement without the other party's prior written consent in the event of a merger with an affiliate or other reorganization with an affiliate. Any attempted assignment contrary to the terms of this Agreement shall be void.
- 3. Ownership of Software.** As between the parties, Akoya owns all right, title and interest, including all existing or future copyrights, trademarks, tools, know-how, and processes, in and to the Software and any copies, modifications, alterations or derivative works thereof. You shall have no title to or



ownership of the Software. As between the parties, Akoya shall retain all rights to market, license and/or sell the Software. The Software (including accompanying documentation) constitutes the unpublished confidential information and trade secrets of Akoya. You (i) acknowledge, and agree not to contest, Akoya's rights in the Software; and (ii) agree not to disclose and to take reasonable steps to protect any and all confidential information of Akoya regarding the Software or that is otherwise disclosed to you in connection with this Agreement. You agree that you are responsible for ensuring your employees and agents are prohibited from taking any action otherwise prohibited to you under this Agreement and to notify Akoya promptly and in writing of the circumstances surrounding any possession, use or knowledge of the Software or any part thereof known to you which is prohibited hereunder. You recognize that the covenants contained in this Agreement are reasonable and necessary to protect the legitimate interests of Akoya, that Akoya would not have entered into this Agreement in the absence of such covenants, and that your breach or threatened breach of such covenants shall cause Akoya irreparable harm and significant injury, the amount of which shall be extremely difficult to estimate and ascertain, thus making any remedy at law or in damages inadequate. Therefore, you agree that Akoya shall be entitled, without the necessity of posting any bond or security, to the issuance of injunctive relief by any court of competent jurisdiction enjoining any breach or threatened breach of such covenants and for any other relief such court deems appropriate. This right shall be in addition to any other remedy available to Akoya at law or in equity.

4. **Termination.** The license granted in Section 1 above is effective until terminated. This Agreement is conditioned upon your continued compliance with the terms and conditions hereof and will terminate automatically without notice from Akoya if you fail to comply with any term or condition of this Agreement. Furthermore, Akoya may terminate this Agreement at any time upon thirty (30) days' notice. Upon termination of this Agreement, your license to use the Software or any part thereof shall immediately terminate, and you shall immediately destroy all copies of the Software (including all accompanying documentation) and any other confidential and proprietary information you have received during or in connection with this Agreement. Termination of this Agreement by Akoya shall not entitle you to a refund of any fees paid for the license of the Software.
5. **Limited Warranty.** The sole and exclusive warranty applicable to software and firmware products provided by Akoya for use with a processor is as follows: Akoya warrants that such software and firmware will substantially conform to Akoya's documentation and program manuals current at the time of shipment to you when properly installed on the processor. The limited warranty set forth above shall not apply if the claimed problem or issue is caused by: the malfunction of computer hardware or software not provided by Akoya; your negligence or fault; your failure to follow the instructions set forth in the documentation; changes in the hardware or other software not authorized by Akoya; modifications or changes to the Software not made or suggested by Akoya; or your failure to implement and maintain a proper backup and recovery system. Your exclusive remedy under this Section 5 shall be, at Akoya's option, the replacement or correction of the problem or issue causing the non-conformance with the above warranty as promptly as possible following notice of such non-conformity, or a refund of the price paid for the Software.
6. **No Other Warranties.** EXCEPT FOR THE LIMITED WARRANTY STATED IMMEDIATELY ABOVE AND TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, THE SOFTWARE IS PROVIDED "AS IS" WITHOUT WARRANTY OF ANY KIND, AND AKOYA EXPRESSLY DISCLAIMS ANY AND ALL IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT OF THIRD PARTY RIGHTS. AKOYA DOES NOT REPRESENT OR WARRANT THAT THE RESULTS OR THE USE OF THE SOFTWARE WILL BE CORRECT, ACCURATE OR RELIABLE, WILL MEET YOUR REQUIREMENTS, OR THAT THE SOFTWARE WILL OPERATE UNINTERRUPTED OR ERROR-FREE, OR THAT DEFECTS IN THE SOFTWARE WILL BE CORRECTED. Some jurisdictions do not allow the exclusion of implied warranties, so the above disclaimer may not apply to you. YOU ASSUME ALL RISK ASSOCIATED WITH THE USE, RESULTS, AND PERFORMANCE OF THE SOFTWARE.

7. **Limitation of Liability.** IN NO EVENT SHALL AKOYA, ITS AFFILIATES OR SUPPLIERS, OR THEIR RESPECTIVE EMPLOYEES, OFFICERS OR AGENTS, BE LIABLE FOR ANY DAMAGES ARISING OUT OF THE DELIVERY, INSTALLATION, PERFORMANCE, OR USE OR INABILITY TO USE THE SOFTWARE, INCLUDING WITHOUT LIMITATION INCIDENTAL, SPECIAL, CONSEQUENTIAL, PUNITIVE, EXEMPLARY, INDIRECT OR DIRECT DAMAGES (INCLUDING WITHOUT LIMITATION DAMAGES FOR LOSS OF PROFITS, LOSS OF DATA, RE-RUN TIME, INACCURATE INPUT, WORK DELAYS, BUSINESS INTERRUPTION OR ANY OTHER COMMERCIAL DAMAGES OR LOSSES), WHETHER IN AN ACTION IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND EVEN IF AKOYA HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, AND REGARDLESS OF WHETHER ANY REMEDY FAILS OF ITS ESSENTIAL PURPOSE. AKOYA'S AGGREGATE LIABILITY ARISING OUT OF OR RELATING TO THIS AGREEMENT, FOR ANY AND ALL DAMAGES THAT YOU MIGHT INCUR (REGARDLESS OF THE FORM OF ACTION GIVING RISE TO SUCH LIABILITY, WHETHER IN CONTRACT, TORT OR OTHERWISE) SHALL NOT EXCEED THE PRICE PAID BY YOU FOR THE LICENSE OF THE SOFTWARE. Some jurisdictions do not allow the exclusion or limitation of liability for consequential or incidental damages. In such jurisdictions, our liability is limited to the greatest extent permitted by law, or the amount you paid for your purchase, whichever is less.
8. **U.S. Government End Users.** The Software qualifies as commercial computer software for purposes of the Federal Acquisition Regulations (FAR) §52.227-19 and the Department of Defense Supplement to the FAR (DFARS) §52.227-7013. If the Software is acquired by a civilian government agency, it is furnished with only the minimum Restricted Rights provided by FAR 52.227-19. If the Software is acquired by a military agency, it is furnished with only the minimum Restricted Rights provided by DFARS 52.227-7013(c)(1)(ii).
9. **Miscellaneous.** This Agreement contains the entire agreement of the parties with respect to the subject matter hereof and supersedes any proposal or prior agreement, written or oral, and any other communications between the parties relating to the subject matter hereof. No modification or waiver of any provision of this Agreement shall be effective unless in writing and signed by the parties. No delay or failure on the part of any party in exercising any right under this Agreement shall impair any such right or any remedy of such party, nor shall it be construed to be a waiver of any continuing breach or default under this Agreement. No action, regardless of form, arising out of the transactions contemplated herein may be brought by either party more than one (1) year after the cause of action has accrued. In the event any provision of this Agreement is held to be unenforceable, the remaining provisions of this Agreement will remain in full force and effect. This Agreement shall be governed by the laws of the State of Delaware without regard to principles of conflicts of laws. Any disputes relating hereto shall be adjudicated only in any court of competent jurisdiction in the State of Delaware, and you hereby consent to the exclusive jurisdiction of those courts. This Agreement shall not be governed by the United Nations Convention on Contracts for the International Sale of Goods, the application of which is expressly excluded. You may not assign or otherwise transfer this Agreement or any of your rights or obligations therein without the prior written consent of Akoya. You may not use the Software for any unlawful purpose nor export or re-export the Software except as authorized by law.

10. **TIFF module: This software uses the LibTIFF 4.0.6 library (<http://remotesensing.org/libtiff/>)**  
Copyright (c) 1988-1997 Sam Leffler Copyright (c) 1991-1997 Silicon Graphics, Inc.

Permission to use, copy, modify, distribute, and sell this software and its documentation for any purpose is hereby granted without fee, provided that (i) the above copyright notices and this permission notice appear in all copies of the software and related documentation, and (ii) the names of Sam Leffler and Silicon Graphics may not be used in any advertising or publicity relating to the software without the specific, prior written permission of Sam Leffler and Silicon Graphics.

THE SOFTWARE IS PROVIDED "AS-IS" AND WITHOUT WARRANTY OF ANY KIND, EXPRESS, IMPLIED OR OTHERWISE, INCLUDING WITHOUT LIMITATION, ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

IN NO EVENT SHALL SAM LEFFLER OR SILICON GRAPHICS BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR ANY DAMAGES WHATSOEVER RESULTING FROM LOSS OF USE, DATA OR PROFITS, WHETHER OR NOT ADVISED OF THE POSSIBILITY OF DAMAGE, AND ON ANY THEORY OF LIABILITY, ARISING OUT OF OR IN CONNECTION WITH THE USE OR PERFORMANCE OF THIS SOFTWARE.

**11. This software uses the Vigna 1.6.0 library (<http://hci.iwr.uni-heidelberg.de/vigna/>)**

This software includes machine-executable object code generated by a source language processor from the Vigna libraries covered by the VIGRA license. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose.

The Vigna library is Copyright © 2008, Ulrich Kothe, Heidelberg Collaboratory for Image Processing

The VIGRA License

=====

(identical to the MIT X11 License)

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

**12. This software uses the Boost 1.53.0 library (<http://www.boost.org/> )**

This software includes machine-executable object code generated by a source language processor from the Boost libraries covered by the Boost license. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose.

**13. This software uses the Quan 1.0.0 library (<https://lists.sourceforge.net/lists/listinfo/quantum-matters>)**

This software includes machine-executable object code generated by a source language processor from the Quan libraries covered by the Boost license. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose.

### **The Boost license is reproduced below:**

Boost Software License - Version 1.0 - August 17th, 2003

Permission is hereby granted, free of charge, to any person or organization obtaining a copy of the software and accompanying documentation covered by this license (the "Software") to use, reproduce, display, distribute, execute, and transmit the Software, and to prepare derivative works of the Software, and to permit third-parties to whom the Software is furnished to do so, all subject to the following:

The copyright notices in the Software and this entire statement, including the above license grant, this restriction and the following disclaimer, must be included in all copies of the Software, in whole or in part, and all derivative works of the Software, unless such copies or derivative works are solely in the form of machine-executable object code generated by a source language processor.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE AND NON-INFRINGEMENT. IN NO EVENT SHALL THE COPYRIGHT HOLDERS OR ANYONE DISTRIBUTING THE SOFTWARE BE LIABLE FOR ANY DAMAGES OR OTHER LIABILITY, WHETHER IN CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

#### **14. This software uses the libjpeg-turbo 1.5.0 library (<http://libjpeg-turbo.virtualgl.org>)**

This software makes use of the libjpeg-turbo library. This software is based in part on the work of the Independent JPEG Group. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by IJG License. Additionally libjpeg-turbo incorporates parts of libjpeg which is covered by the BSD license.

### **The IJG license is reproduced below:**

The authors make NO WARRANTY or representation, either express or implied, with respect to this software, its quality, accuracy, merchantability, or fitness for a particular purpose. This software is provided "AS IS", and you, its user, assume the entire risk as to its quality and accuracy.

This software is copyright (C) 1991-1998, Thomas G. Lane. All Rights Reserved except as specified below.

Permission is hereby granted to use, copy, modify, and distribute this software (or portions thereof) for any purpose, without fee, subject to these conditions:

- (1) If any part of the source code for this software is distributed, then this README file must be included, with this copyright and no-warranty notice unaltered; and any additions, deletions, or changes to the original files must be clearly indicated in accompanying documentation.
- (2) If only executable code is distributed, then the accompanying documentation must state that "this software is based in part on the work of the Independent JPEG Group".
- (3) Permission for use of this software is granted only if the user accepts full responsibility for any undesirable consequences; the authors accept NO LIABILITY for damages of any kind.

These conditions apply to any software derived from or based on the IJG code, not just to the unmodified library. If you use our work, you ought to acknowledge us.

Permission is NOT granted for the use of any IJG author's name or company name in advertising or publicity relating to this software or products derived from it. This software may be referred to only as "the Independent JPEG Group's software".

We specifically permit and encourage the use of this software as the basis of commercial products, provided that all warranty or liability claims are assumed by the product vendor.

ansi2knr.c is included in this distribution by permission of L. Peter Deutsch, sole proprietor of its copyright holder, Aladdin Enterprises of Menlo Park, CA. ansi2knr.c is NOT covered by the above copyright and conditions, but instead by the usual distribution terms of the Free Software Foundation; principally, that you must include source code if you redistribute it. (See the file ansi2knr.c for full details.) However, since ansi2knr.c is not needed as part of any program generated from the IJG code, this does not limit you more than the foregoing paragraphs do.

The Unix configuration script "configure" was produced with GNU Autoconf. It is copyright by the Free Software Foundation but is freely distributable. The same holds for its supporting scripts (config.guess, config.sub, ltconfig, ltmain.sh). Another support script, install-sh, is copyright by M.I.T. but is also freely distributable.

It appears that the arithmetic coding option of the JPEG spec is covered by patents owned by IBM, AT&T, and Mitsubishi. Hence arithmetic coding cannot legally be used without obtaining one or more licenses. For this reason, support for arithmetic coding has been removed from the free JPEG software. (Since arithmetic coding provides only a marginal gain over the unpatented Huffman mode, it is unlikely that very many implementations will support it.) So far as we are aware, there are no patent restrictions on the remaining code.

The IJG distribution formerly included code to read and write GIF files. To avoid entanglement with the Unisys LZW patent, GIF reading support has been removed altogether, and the GIF writer has been simplified to produce "uncompressed GIFs". This technique does not use the LZW algorithm; the resulting GIF files are larger than usual, but are readable by all standard GIF decoders.

We are required to state that "The Graphics Interchange Format(c) is the Copyright property of CompuServe Incorporated. GIF(sm) is a Service Mark property of CompuServe Incorporated."

#### **The BSD for libjpeg license is reproduced below:**

The authors make NO WARRANTY or representation, either express or implied, with respect to this software, its quality, accuracy, merchantability, or fitness for a particular purpose. This software is provided "AS IS", and you, its user, assume the entire risk as to its quality and accuracy.

This software is copyright (C) 1991-2010, Thomas G. Lane, Guido Vollbeding.  
All Rights Reserved except as specified below.

Permission is hereby granted to use, copy, modify, and distribute this software (or portions thereof) for any purpose, without fee, subject to these conditions:

(1) If any part of the source code for this software is distributed, then this README file must be included, with this copyright and no-warranty notice unaltered; and any additions, deletions, or changes to the original files must be clearly indicated in accompanying documentation.

(2) If only executable code is distributed, then the accompanying documentation must state that "this software is based in part on the work of the Independent JPEG Group".

(3) Permission for use of this software is granted only if the user accepts full responsibility for any undesirable consequences; the authors accept NO LIABILITY for damages of any kind.

These conditions apply to any software derived from or based on the IJG code, not just to the unmodified library. If you use our work, you ought to acknowledge us.

Permission is NOT granted for the use of any IJG author's name or company name in advertising or publicity relating to this software or products derived from it. This software may be referred to only as "the Independent JPEG Group's software".

We specifically permit and encourage the use of this software as the basis of commercial products, provided that all warranty or liability claims are assumed by the product vendor.

**15. This software uses the zlib1 1.2.8 library (<http://www.zlib.net>)**

This software makes use of the zlib1 library. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by the zlib1 License.

**The zlib 1 license is reproduced below:**

(C) 1995-2013 Jean-loup Gailly and Mark Adler

This software is provided 'as-is', without any express or implied warranty. In no event will the authors be held liable for any damages arising from the use of this software.

Permission is granted to anyone to use this software for any purpose, including commercial applications, and to alter it and redistribute it freely, subject to the following restrictions:

1. The origin of this software must not be misrepresented; you must not claim that you wrote the original software. If you use this software in a product, an acknowledgment in the product documentation would be appreciated but is not required.

2. Altered source versions must be plainly marked as such, and must not be misrepresented as being the original software.

3. This notice may not be removed or altered from any source distribution.

Jean-loup Gailly	Mark Adler
<a href="mailto:jloup@gzip.org">jloup@gzip.org</a>	<a href="mailto:madler@alumni.caltech.edu">madler@alumni.caltech.edu</a>

**16. This software uses the Fugue icons, version 3.5.6 (<http://p.yusukekamiyamane.com/>)**

This software makes use of the Fugue icons. This content has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The icons are covered by the Creative Commons 3.0 License.

(C) 2013 Yusuke Kamiyamane. All rights reserved.

The Creative Commons 3.0 License is reproduced below:

THE WORK (AS DEFINED BELOW) IS PROVIDED UNDER THE TERMS OF THIS CREATIVE COMMONS PUBLIC LICENSE ("CCPL" OR "LICENSE"). THE WORK IS PROTECTED BY COPYRIGHT AND/OR OTHER APPLICABLE LAW. ANY USE OF THE WORK OTHER THAN AS AUTHORIZED UNDER THIS LICENSE OR COPYRIGHT LAW IS PROHIBITED.

BY EXERCISING ANY RIGHTS TO THE WORK PROVIDED HERE, YOU ACCEPT AND AGREE TO BE BOUND BY THE TERMS OF THIS LICENSE. TO THE EXTENT THIS LICENSE MAY BE CONSIDERED TO BE A CONTRACT, THE LICENSOR GRANTS YOU THE RIGHTS CONTAINED HERE IN CONSIDERATION OF YOUR ACCEPTANCE OF SUCH TERMS AND CONDITIONS.

## 1. Definitions

- a. **"Adaptation"** means a work based upon the Work, or upon the Work and other pre-existing works, such as a translation, adaptation, derivative work, arrangement of music or other alterations of a literary or artistic work, or phonogram or performance and includes cinematographic adaptations or any other form in which the Work may be recast, transformed, or adapted including in any form recognizably derived from the original, except that a work that constitutes a Collection will not be considered an Adaptation for the purpose of this License. For the avoidance of doubt, where the Work is a musical work, performance or phonogram, the synchronization of the Work in timed-relation with a moving image ("synching") will be considered an Adaptation for the purpose of this License.
- b. **"Collection"** means a collection of literary or artistic works, such as encyclopedias and anthologies, or performances, phonograms or broadcasts, or other works or subject matter other than works listed in Section 1(f) below, which, by reason of the selection and arrangement of their contents, constitute intellectual creations, in which the Work is included in its entirety in unmodified form along with one or more other contributions, each constituting separate and independent works in themselves, which together are assembled into a collective whole. A work that constitutes a Collection will not be considered an Adaptation (as defined above) for the purposes of this License.
- c. **"Distribute"** means to make available to the public the original and copies of the Work or Adaptation, as appropriate, through sale or other transfer of ownership.
- d. **"Licensor"** means the individual, individuals, entity or entities that offer(s) the Work under the terms of this License.
- e. **"Original Author"** means, in the case of a literary or artistic work, the individual, individuals, entity or entities who created the Work or if no individual or entity can be identified, the publisher; and in addition (i) in the case of a performance the actors, singers, musicians, dancers, and other persons who act, sing, deliver, declaim, play in, interpret or otherwise perform literary or artistic works or expressions of folklore; (ii) in the case of a phonogram the producer being the person or legal entity who first fixes the sounds of a performance or other sounds; and, (iii) in the case of broadcasts, the organization that transmits the broadcast.
- f. **"Work"** means the literary and/or artistic work offered under the terms of this License including without limitation any production in the literary, scientific and artistic domain, whatever may be the mode or form of its expression including digital form, such as a book, pamphlet and other writing; a lecture, address, sermon or other work of the same nature; a dramatic or dramatico-musical work; a choreographic work or entertainment in dumb show; a musical composition with or without words; a cinematographic work to which are assimilated works expressed by a process analogous to cinematography; a work of drawing, painting, architecture, sculpture, engraving or lithography; a photographic work to which are assimilated works expressed by a

process analogous to photography; a work of applied art; an illustration, map, plan, sketch or three-dimensional work relative to geography, topography, architecture or science; a performance; a broadcast; a phonogram; a compilation of data to the extent it is protected as a copyrightable work; or a work performed by a variety or circus performer to the extent it is not otherwise considered a literary or artistic work.

- g. **"You"** means an individual or entity exercising rights under this License who has not previously violated the terms of this License with respect to the Work, or who has received express permission from the Licensor to exercise rights under this License despite a previous violation.
- h. **"Publicly Perform"** means to perform public recitations of the Work and to communicate to the public those public recitations, by any means or process, including by wire or wireless means or public digital performances; to make available to the public Works in such a way that members of the public may access these Works from a place and at a place individually chosen by them; to perform the Work to the public by any means or process and the communication to the public of the performances of the Work, including by public digital performance; to broadcast and rebroadcast the Work by any means including signs, sounds or images.
- i. **"Reproduce"** means to make copies of the Work by any means including without limitation by sound or visual recordings and the right of fixation and reproducing fixations of the Work, including storage of a protected performance or phonogram in digital form or other electronic medium.

**2. Fair Dealing Rights.** Nothing in this License is intended to reduce, limit, or restrict any uses free from copyright or rights arising from limitations or exceptions that are provided for in connection with the copyright protection under copyright law or other applicable laws.

**3. License Grant.** Subject to the terms and conditions of this License, Licensor hereby grants You a worldwide, royalty-free, non-exclusive, perpetual (for the duration of the applicable copyright) license to exercise the rights in the Work as stated below:

- a. to Reproduce the Work, to incorporate the Work into one or more Collections, and to Reproduce the Work as incorporated in the Collections;
- b. to create and Reproduce Adaptations provided that any such Adaptation, including any translation in any medium, takes reasonable steps to clearly label, demarcate or otherwise identify that changes were made to the original Work. For example, a translation could be marked "The original work was translated from English to Spanish," or a modification could indicate "The original work has been modified.";
- c. to Distribute and Publicly Perform the Work including as incorporated in Collections; and,
- d. to Distribute and Publicly Perform Adaptations.
- e. For the avoidance of doubt:
  - i. **Non-waivable Compulsory License Schemes.** In those jurisdictions in which the right to collect royalties through any statutory or compulsory licensing scheme cannot be waived, the Licensor reserves the exclusive right to collect such royalties for any exercise by You of the rights granted under this License;
  - ii. **Waivable Compulsory License Schemes.** In those jurisdictions in which the right to collect royalties through any statutory or compulsory licensing scheme can be waived, the



Licensor waives the exclusive right to collect such royalties for any exercise by You of the rights granted under this License; and,

- iii. **Voluntary License Schemes.** The Licensor waives the right to collect royalties, whether individually or, in the event that the Licensor is a member of a collecting society that administers voluntary licensing schemes, via that society, from any exercise by You of the rights granted under this License.

The above rights may be exercised in all media and formats whether now known or hereafter devised. The above rights include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. Subject to Section 8(f), all rights not expressly granted by Licensor are hereby reserved.

**4. Restrictions.** The license granted in Section 3 above is expressly made subject to and limited by the following restrictions:

- j. You may Distribute or Publicly Perform the Work only under the terms of this License. You must include a copy of, or the Uniform Resource Identifier (URI) for, this License with every copy of the Work You Distribute or Publicly Perform. You may not offer or impose any terms on the Work that restrict the terms of this License or the ability of the recipient of the Work to exercise the rights granted to that recipient under the terms of the License. You may not sublicense the Work. You must keep intact all notices that refer to this License and to the disclaimer of warranties with every copy of the Work You Distribute or Publicly Perform. When You Distribute or Publicly Perform the Work, You may not impose any effective technological measures on the Work that restrict the ability of a recipient of the Work from You to exercise the rights granted to that recipient under the terms of the License. This Section 4(a) applies to the Work as incorporated in a Collection, but this does not require the Collection apart from the Work itself to be made subject to the terms of this License. If You create a Collection, upon notice from any Licensor You must, to the extent practicable, remove from the Collection any credit as required by Section 4(b), as requested. If You create an Adaptation, upon notice from any Licensor You must, to the extent practicable, remove from the Adaptation any credit as required by Section 4(b), as requested.
- k. If You Distribute, or Publicly Perform the Work or any Adaptations or Collections, You must, unless a request has been made pursuant to Section 4(a), keep intact all copyright notices for the Work and provide, reasonable to the medium or means You are utilizing: (i) the name of the Original Author (or pseudonym, if applicable) if supplied, and/or if the Original Author and/or Licensor designate another party or parties (e.g., a sponsor institute, publishing entity, journal) for attribution ("Attribution Parties") in Licensor's copyright notice, terms of service or by other reasonable means, the name of such party or parties; (ii) the title of the Work if supplied; (iii) to the extent reasonably practicable, the URI, if any, that Licensor specifies to be associated with the Work, unless such URI does not refer to the copyright notice or licensing information for the Work; and (iv) , consistent with Section 3(b), in the case of an Adaptation, a credit identifying the use of the Work in the Adaptation (e.g., "French translation of the Work by Original Author," or "Screenplay based on original Work by Original Author"). The credit required by this Section 4 (b) may be implemented in any reasonable manner; provided, however, that in the case of a Adaptation or Collection, at a minimum such credit will appear, if a credit for all contributing authors of the Adaptation or Collection appears, then as part of these credits and in a manner at least as prominent as the credits for the other contributing authors. For the avoidance of doubt, You may only use the credit required by this Section for the purpose of attribution in the manner set out above and, by exercising Your rights under this License, You may not implicitly or explicitly assert or imply any connection with, sponsorship or endorsement by the Original Author, Licensor and/or Attribution Parties, as appropriate, of You or Your use of the Work,

without the separate, express prior written permission of the Original Author, Licensor and/or Attribution Parties.

- I. Except as otherwise agreed in writing by the Licensor or as may be otherwise permitted by applicable law, if You Reproduce, Distribute or Publicly Perform the Work either by itself or as part of any Adaptations or Collections, You must not distort, mutilate, modify or take other derogatory action in relation to the Work which would be prejudicial to the Original Author's honor or reputation. Licensor agrees that in those jurisdictions (e.g. Japan), in which any exercise of the right granted in Section 3(b) of this License (the right to make Adaptations) would be deemed to be a distortion, mutilation, modification or other derogatory action prejudicial to the Original Author's honor and reputation, the Licensor will waive or not assert, as appropriate, this Section, to the fullest extent permitted by the applicable national law, to enable You to reasonably exercise Your right under Section 3(b) of this License (right to make Adaptations) but not otherwise.

## **5. Representations, Warranties and Disclaimer**

UNLESS OTHERWISE MUTUALLY AGREED TO BY THE PARTIES IN WRITING, LICENSOR OFFERS THE WORK AS-IS AND MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND CONCERNING THE WORK, EXPRESS, IMPLIED, STATUTORY OR OTHERWISE, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF TITLE, MERCHANTIBILITY, FITNESS FOR A PARTICULAR PURPOSE, NONINFRINGEMENT, OR THE ABSENCE OF LATENT OR OTHER DEFECTS, ACCURACY, OR THE PRESENCE OF ABSENCE OF ERRORS, WHETHER OR NOT DISCOVERABLE. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES, SO SUCH EXCLUSION MAY NOT APPLY TO YOU.

**6. Limitation on Liability.** EXCEPT TO THE EXTENT REQUIRED BY APPLICABLE LAW, IN NO EVENT WILL LICENSOR BE LIABLE TO YOU ON ANY LEGAL THEORY FOR ANY SPECIAL, INCIDENTAL, CONSEQUENTIAL, PUNITIVE OR EXEMPLARY DAMAGES ARISING OUT OF THIS LICENSE OR THE USE OF THE WORK, EVEN IF LICENSOR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

## **7. Termination**

- m. This License and the rights granted hereunder will terminate automatically upon any breach by You of the terms of this License. Individuals or entities who have received Adaptations or Collections from You under this License, however, will not have their licenses terminated provided such individuals or entities remain in full compliance with those licenses. Sections 1, 2, 5, 6, 7, and 8 will survive any termination of this License.
- n. Subject to the above terms and conditions, the license granted here is perpetual (for the duration of the applicable copyright in the Work). Notwithstanding the above, Licensor reserves the right to release the Work under different license terms or to stop distributing the Work at any time; provided, however that any such election will not serve to withdraw this License (or any other license that has been, or is required to be, granted under the terms of this License), and this License will continue in full force and effect unless terminated as stated above.

## **8. Miscellaneous**

- o. Each time You Distribute or Publicly Perform the Work or a Collection, the Licensor offers to the recipient a license to the Work on the same terms and conditions as the license granted to You under this License.

- p. Each time You Distribute or Publicly Perform an Adaptation, Licensor offers to the recipient a license to the original Work on the same terms and conditions as the license granted to You under this License.
- q. If any provision of this License is invalid or unenforceable under applicable law, it shall not affect the validity or enforceability of the remainder of the terms of this License, and without further action by the parties to this agreement, such provision shall be reformed to the minimum extent necessary to make such provision valid and enforceable.
- r. No term or provision of this License shall be deemed waived and no breach consented to unless such waiver or consent shall be in writing and signed by the party to be charged with such waiver or consent.
- s. This License constitutes the entire agreement between the parties with respect to the Work licensed here. There are no understandings, agreements or representations with respect to the Work not specified here. Licensor shall not be bound by any additional provisions that may appear in any communication from You. This License may not be modified without the mutual written agreement of the Licensor and You.
- t. The rights granted under, and the subject matter referenced, in this License were drafted utilizing the terminology of the Berne Convention for the Protection of Literary and Artistic Works (as amended on September 28, 1979), the Rome Convention of 1961, the WIPO Copyright Treaty of 1996, the WIPO Performances and Phonograms Treaty of 1996 and the Universal Copyright Convention (as revised on July 24, 1971). These rights and subject matter take effect in the relevant jurisdiction in which the License terms are sought to be enforced according to the corresponding provisions of the implementation of those treaty provisions in the applicable national law. If the standard suite of rights granted under applicable copyright law includes additional rights not granted under this License, such additional rights are deemed to be included in the License; this License is not intended to restrict the license of any rights under applicable law.

**17. This software uses the log4net libraries (<http://logging.apache.org>)**

This software makes use of the log4net libraries. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The libraries are covered by the Apache License.

**18. This software uses the zxing.NET libraries (<https://zxingnet.codeplex.com/>)**

This software makes use of the zxing.NET libraries. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The libraries are covered by the Apache License

The Apache License is reproduced below:

**Apache License**

Version 2.0, January 2004

<http://www.apache.org/licenses/>

## TERMS AND CONDITIONS FOR USE, REPRODUCTION, AND DISTRIBUTION

### 1. Definitions.

"License" shall mean the terms and conditions for use, reproduction, and distribution as defined by Sections 1 through 9 of this document.

"Licensor" shall mean the copyright owner or entity authorized by the copyright owner that is granting the License.

"Legal Entity" shall mean the union of the acting entity and all other entities that control, are controlled by, or are under common control with that entity. For the purposes of this definition, "control" means (i) the power, direct or indirect, to cause the direction or management of such entity, whether by contract or otherwise, or (ii) ownership of fifty percent (50%) or more of the outstanding shares, or (iii) beneficial ownership of such entity.

"You" (or "Your") shall mean an individual or Legal Entity exercising permissions granted by this License.

"Source" form shall mean the preferred form for making modifications, including but not limited to software source code, documentation source, and configuration files.

"Object" form shall mean any form resulting from mechanical transformation or translation of a Source form, including but not limited to compiled object code, generated documentation, and conversions to other media types.

"Work" shall mean the work of authorship, whether in Source or Object form, made available under the License, as indicated by a copyright notice that is included in or attached to the work (an example is provided in the Appendix below).

"Derivative Works" shall mean any work, whether in Source or Object form, that is based on (or derived from) the Work and for which the editorial revisions, annotations, elaborations, or other modifications represent, as a whole, an original work of authorship. For the purposes of this License, Derivative Works shall not include works that remain separable from, or merely link (or bind by name) to the interfaces of, the Work and Derivative Works thereof.

"Contribution" shall mean any work of authorship, including the original version of the Work and any modifications or additions to that Work or Derivative Works thereof, that is intentionally submitted to Licensor for inclusion in the Work by the copyright owner or by an individual or Legal Entity authorized to submit on behalf of the copyright owner. For the purposes of this definition, "submitted" means any form of electronic, verbal, or written communication sent to the Licensor or its representatives, including but not limited to communication on electronic mailing lists, source code control systems, and issue tracking systems that are managed by, or on behalf of, the Licensor for the purpose of discussing and improving the Work, but excluding communication that is conspicuously marked or otherwise designated in writing by the copyright owner as "Not a Contribution."

"Contributor" shall mean Licensor and any individual or Legal Entity on behalf of whom a Contribution has been received by Licensor and subsequently incorporated within the Work.

2. Grant of Copyright License. Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable copyright license to reproduce, prepare Derivative Works of, publicly display, publicly perform, sublicense, and distribute the Work and such Derivative Works in Source or Object form.

3. Grant of Patent License. Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in this section) patent license to make, have made, use, offer to sell, sell, import, and otherwise transfer the Work, where such license applies only to those patent claims licensable by such Contributor that are necessarily infringed by their Contribution(s) alone or by combination of their Contribution(s) with the Work to which such Contribution(s) was submitted. If You institute patent litigation against any entity (including a cross-claim or counterclaim in a lawsuit) alleging that the Work or a Contribution incorporated within the Work constitutes direct or contributory patent infringement, then any patent licenses granted to You under this License for that Work shall terminate as of the date such litigation is filed.

4. Redistribution. You may reproduce and distribute copies of the Work or Derivative Works thereof in any medium, with or without modifications, and in Source or Object form, provided that You meet the following conditions:

(a) You must give any other recipients of the Work or Derivative Works a copy of this License; and

(b) You must cause any modified files to carry prominent notices stating that You changed the files; and

(c) You must retain, in the Source form of any Derivative Works that You distribute, all copyright, patent, trademark, and attribution notices from the Source form of the Work, excluding those notices that do not pertain to any part of the Derivative Works; and

(d) If the Work includes a "NOTICE" text file as part of its distribution, then any Derivative Works that You distribute must include a readable copy of the attribution notices contained within such NOTICE file, excluding those notices that do not pertain to any part of the Derivative Works, in at least one of the following places: within a NOTICE text file distributed as part of the Derivative Works; within the Source form or documentation, if provided along with the Derivative Works; or, within a display generated by the Derivative Works, if and wherever such third-party notices normally appear. The contents of the NOTICE file are for informational purposes only and do not modify the License. You may add Your own attribution notices within Derivative Works that You distribute, alongside or as an addendum to the NOTICE text from the Work, provided that such additional attribution notices cannot be construed as modifying the License.

You may add Your own copyright statement to Your modifications and may provide additional or different license terms and conditions for use, reproduction, or distribution of Your modifications, or for any such Derivative Works as a whole, provided Your use, reproduction, and distribution of the Work otherwise complies with the conditions stated in this License.

5. Submission of Contributions. Unless You explicitly state otherwise, any Contribution intentionally submitted for inclusion in the Work by You to the Licensor shall be under the terms and conditions of this License, without any additional terms or conditions. Notwithstanding the above, nothing herein shall supersede or modify the terms of any separate license agreement you may have executed with Licensor regarding such Contributions.

6. Trademarks. This License does not grant permission to use the trade names, trademarks, service marks, or product names of the Licensor, except as required for reasonable and customary use in describing the origin of the Work and reproducing the content of the NOTICE file.

7. Disclaimer of Warranty. Unless required by applicable law or agreed to in writing, Licensor provides the Work (and each Contributor provides its Contributions) on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied, including, without limitation, any warranties or conditions of TITLE, NON-INFRINGEMENT, MERCHANTABILITY, or FITNESS FOR A PARTICULAR PURPOSE. You are solely responsible for determining the appropriateness of using or redistributing the Work and assume any risks associated with Your exercise of permissions under this License.

8. Limitation of Liability. In no event and under no legal theory, whether in tort (including negligence), contract, or otherwise, unless required by applicable law (such as deliberate and grossly negligent acts) or agreed to in writing, shall any Contributor be liable to You for damages, including any direct, indirect, special, incidental, or consequential damages of any character arising as a result of this License or out of the use or inability to use the Work (including but not limited to damages for loss of goodwill, work stoppage, computer failure or malfunction, or any and all other commercial damages or losses), even if such Contributor has been advised of the possibility of such damages.

9. Accepting Warranty or Additional Liability. While redistributing the Work or Derivative Works thereof, You may choose to offer, and charge a fee for, acceptance of support, warranty, indemnity, or other liability obligations and/or rights consistent with this License. However, in accepting such obligations, You may act only on Your own behalf and on Your sole responsibility, not on behalf of any other Contributor, and only if You agree to indemnify, defend, and hold each Contributor harmless for any liability incurred by, or claims asserted against, such Contributor by reason of your accepting any such warranty or additional liability.

**19. This software uses FTDI drivers (<http://www.ftdichip.com>)**

This software makes use of drivers provided by Future Technology Devices International Limited, which are "as is" and any express or implied warranties, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose are disclaimed. In no event shall future technology devices international limited be liable for any direct, indirect, incidental, special, exemplary, or consequential damages (including, but not limited to, procurement of substitute goods or services; loss of use, data, or profits; or business interruption) however caused and on any theory of liability, whether in contract, strict liability, or tort (including negligence or otherwise) arising in any way out of the use of this software, even if advised of the possibility of such damage. FTDI drivers may be used only in conjunction with products based on FTDI parts. FTDI drivers may be distributed in any form as long as license information is not modified.

Copyright Future Technology Devices International Limited of 2 Seaward Place, Centurion Business Park, Glasgow, Scotland, G41 1HH

**Should you have any question concerning this Agreement, you may contact Akoya by writing to Akoya Biosciences, 100 Campus Drive, 6th Floor, Marlborough MA 01752. You may also call +1 855.896.8401**



# Index

## - B -

Barcode reading  
    brightfield protocol options 57  
    fluorescence protocol options 68  
Brightfield  
    compensation information 52

## - C -

Contact 6

## - D -

Door sensor 28

## - F -

Fluorescence  
    compensation information 52  
Fuses  
    specifications 37

## - H -

Hardware  
    cover handle 30  
    door sensor 28  
    power connector 31  
    slide carrier 32  
    slide carrier hotel 28  
    slide carrier hotel slot 32  
    status lights 29, 32  
    status progress display 28  
    system door 29  
    system power switch 31  
    usb 31, 36

## - I -

inForm 16

## - P -

PerkinElmer TIFF Specification 89  
Phenochart 16  
Power connector 31  
Preface 5  
Protocols  
    brightfield 56  
    brightfield advanced settings 57  
    brightfield barcode reading 57  
    brightfield snapshots 59  
    definition 53  
    fluorescence 62  
    fluorescence advanced settings 68  
    fluorescence barcode reading 68  
    fluorescent snapshots 74

## - S -

Service 6  
Service plans 6  
Slide Carrier 32  
    handle 34  
    insert indicator 34  
    slide positions 34  
    slides 34  
    spring-loaded tab 34  
    status lights 29  
    tab cover 34  
Slides 34  
    format 34  
    status progress display 28  
Snapshots  
    brightfield 59  
    fluorescence 74  
Software  
    saving and restoring setup 83  
    slide IDs 78  
Specifications  
    barcode reader 37  
    electrical 37  
    environmental 36  
    fuses 37  
    system 36  
Status lights 29  
Status Progress Display 28



Stokes Shift 20  
Support 6  
Support Plans 6  
System door 29  
System Power Switch 31

## - T -

Theory of Imaging  
    ability to distinguish colors 19  
    Stokes Shift 20  
TIFF Specification (PerkinElmer) 89  
Training 6

## - U -

USB 31, 36

To learn more visit [AKOYABIO.COM](https://www.akoyabio.com) or email us at [INFO@AKOYABIO.COM](mailto:INFO@AKOYABIO.COM)

For Research Use Only. Not for use in diagnostic procedures.

© 2022 Akoya Biosciences, Inc. All rights reserved. All trademarks are the property of Akoya Biosciences unless otherwise specified.