TaKaRa

TaKaRa PCR Thermal Cycler Dice™

Version

Model TP600 / TP650
Gradient / Standard

Operation Manual 取扱説明書

LR108182

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NOTICE

< Only for US >

Purchase of this instrument conveys a limited non-transferable immunity from suit for the purchaser's own internal research and development and applied fields other than human in vitro diagnostics under one or more of U.S. Patents No. 5,038,852, 5,656,493, 5,333,675, 5,475,610, and 6,703,236 (claims 1-6 only), or corresponding claims in their non-U.S. counterparts, owned by Applera Corporation. No right is conveyed expressly, by implication or by estoppel under any other patent claim, including under U.S. Patent No. 6,814,934, which claims thermal cycler apparatus and systems capable of real-time detection. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Applied Biosystems does not guarantee the performance of this instrument.

< Outside US >

Purchase of this instrument conveys a limited non-transferable immunity from suit for the purchaser's own internal research and development and applied fields other than human in vitro diagnostics under one or more of U.S. Patents No. 5,038,852, 5,656,493, 5,333,675, 5,475,610, and 6,703,236 (claims 1-6 only), or corresponding claims in their non-U.S. counterparts, owned by Applera Corporation. No right is conveyed expressly, by implication or by estoppel under any other patent claim. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Applied Biosystems does not guarantee the performance of this instrument.

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Tel. 077-543-7231

Precautions for Using TaKaRa Scientific Instruments

- Users are recommended to carefully read the content of this manual before operating the TaKaRa PCR Thermal Cycler Dice system.
- After familiarizing themselves with the information in this manual, users are recommended to store it in secure and convenient location from which it is always available for reference

Introduction

All Takara scientific instruments are sold strictly for use in research applications and should not be used by those who are not experts in this field.

Due to ongoing improvements, product specifications and the external appearance are subject to change without prior notification.

Warranty and Responsibilities

Refer to the Warranty card included with this product for detailed description of the Warranty conditions and period.

This product is designed to ultimate safety standards. Nevertheless, should any safety-related issues arise, users are requested to immediately contact the nearest sales agent.

Precautions (safety and accident prevention) for System Use

(1) Danger Warnings

The various messages and labels provided on the product itself and in the various manuals are designed to ensure the system is operated safely to prevent personal injury or damage to property. Users are strongly recommended to read and observe the warnings included with the symbols shown below.

Warnings indicate incorrect system operation and handling procedures that can lead to serious accidents.

Cautions indicate incorrect system operation and handling procedures that can lead to accidents.

(2) General Precautions Regarding Takara Scientific Instruments

It is important to observe the various precautions described below to ensure optimum system performance and prevent accidents before they happen. Ensure safe and correct system operation by carefully reading the precautions.

[Precautions for Installation Location]

Observe the following precautions in order to prevent fires and electrical shocks.

Do not install the system on an unstable or sloped foundation or any other location subject unstable conditions.
If the system tips over or falls, immediately turn off power to the system and all connected devices and disconnect the power supply cables.
Avoid blocking the paths to any air vents included on the product.
Blocked air vents can lead to internal overheating that may result in physical burns or fire.

Avoid installing the system in locations subject to dust or strong vibrations.
Avoid installing in locations subject to direct sunlight and wide temperature fluctuations as this may result in internal overheating and may cause physical burns or fire.
Do not place objects on top of the system.

Do not install the system in locations other than research institutions or related facilities.

[Electrical Precautions]

Observe the following precautions when using electrically powered scientific instruments in order to prevent electric shocks and fires.

Avoid operating system with power sources that do not conform to system power supply specifications.

Avoid sharply twisting, altering, marking or otherwise damaging the power supply cable. Avoid applying stress to, heating or placing objects on the power cable.
Always ground the system to prevent electric shocks due to leakage current.

WARNING

Be sure to compare system power supply requirements with the rated power capacity of the power source to ensure system current requirements do not exceed power socket capacity. Confirm that the total power requirements of all connected equipment do not exceed the rated power supply capacity of the common power source.

The system body and power supply cable are subject to high voltage levels(The maximum voltage is a line voltage AC100V ~ 240V). Avoid inserting metal or flammable objects inside the system or handling it with wet hands, as this can result in electric shock and/or other dangers.

CAUTION

When moving the system, be sure to first disconnect the power cable as well as any peripheral devices.

Be sure to disconnect the power cable prior to cleaning or servicing the system.

[Other precautions]

Observe the following precautions in order to prevent electric shocks and fires.

WARNING

Avoid inserting metal objects, flammable materials or other foreign objects into the various openings on the system case.
Never remove the system cover or attempt to alter it in anyway.

Avoid disassembling or making any direct contact with internal parts of the system.

Should noxious fumes smoke or other abnormal signs appear, immediately terminate operation, turn off system power and disconnect the power cable.

Wait until conditions stabilize before contacting the nearest sales agent or the company technical support line.

Should water or other foreign matter enter the system, immediately terminate operation, turn off system power and disconnect the power cable.

Wait until conditions stabilize before contacting the nearest sales agent or the company technical support line.

CAUTION

Do not sit on the system.
Do not disassemble the system

(3) Warning Safety Labels Applied to System Case

Users are advised to avoid accidents by carefully reading the warnings and cautions contained on warning stickers at key locations on the interior and exterior of the system.

Possible Danger	Warning Type Location of Danger <i>Type d'Avertissement</i> <i>Localisation du Danger</i>	Warning Label	Description of Danger
Electric shocks and fires	Disassembly Prohibited System Interior Démontage Interdit <i>Intérieur de l'Appareil</i>		Never disassemble system, as this may lead to electrical short-circuits, leakage current and/or fire.
Electric shocks	Warning! High Voltage Power cable Attention! Haute tension <i>Cordon d'alimentation</i>		The case and power cable are subject to high voltage levels, so take care to avoid electric shocks.
Electrical shocks and fires	Grounding Mandatory Power cable Mise à la Terre Obligatoire <i>Cordon d'alimentation</i>		Always ground the system to prevent electric shock from leakage current.
Burns	High Temperatures Heat cover & Heat block Haute Température <i>Couvercle Chauffant</i>		Avoid touching the heat cover, which reaches very high temperatures and may cause burns.
Burns and fires	Installation Precaution Air vents Précautions d'Installation <i>Prises d'air</i>		Be sure not to block air vents when installing the system, as this can cause internal overheating.
Personal injury	Do not Open When in Operation System cover Ne Pas Ouvrir en cours de fonctionnement <i>Couvercle</i>		Be sure to close the cover while the system is in operation. Failure to do so can lead to burns and/or damage to the eyes.

This symbol means caution or warning.

This symbol means an action is prohibited.

This symbol means an instruction must be followed.

(4) Electrical Symbols

The symbols shown on this page are adapted to the TaKaRa PCR Thermal Cycler Dice. Read the explanations and make sure you understand what the symbols mean before you interact with the instrument in any way.

This symbol indicates the on position of the main power switch.

This symbol indicates the off position of the main power switch.

Precautions against Transferring or Copying Content of This Document

The content of this document cannot be transferred or copied without the express permission of the copyright holder.

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

Result of Rate Test		MENU
Heat 3.1°C/s	Cool 2.1°C/s	OK
Over 0.1°C	Under 0.4°C	
Started at 15:47:21 03/01/15		
Finished at 15:50:52 03/01/15		
Length of Run 00:03:31		
Result: Normal		

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
Result of Cycle Test		MENU
Heat 3.0°C/s	Cool 2.1°C/s	OK
Over 0.2°C	Under 0.4°C	
Started at 15:39:09 03/01/15		
Finished at 15:42:39 03/01/15		
Length of Run 00:03:30		
Result: Normal		

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
Result of Up Down Test		MENU
Heat	3.2°C/s	Cool 1.4°C/s
Over	0.2°C	Under 0.3°C
Started at	15:59:23 03/01/15	
Finished at	16:05:49 03/01/15	
Length of Run	00:06:26	
Result:	Normal	

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Result of Grad. Test		MENU
Heat	3.1°C/s	Cool 2.0°C/s
Over	0.2°C	Under 0.5°C
Started at	16:21:42 03/01/15	
Finished at	16:31:05 03/01/15	
Length of Run	00:09:23	
Result:	Normal	

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TABLE OF CONTENTS

	page
	Precautions for Using TaKaRa Scientific Instruments
	Table of Contents --- 59
Chapter 1:	Overview --- 60
Chapter 2:	Installation Requirements & Maintenance --- 61
	1. Installation Conditions --- 61
	2. Installation --- 62
	3. Moving the System --- 62
	4. Maintenance --- 63
Chapter 3:	Names of Parts and Their Functions --- 64
Chapter 4:	System Operation --- 67
	1. Introduction --- 67
	2. Turning on System Power --- 68
	3. Starting a System Operation Program --- 69
	4. Terminating or Aborting Operation --- 73
	5. Creating New Program --- 74
	6. Editing and Saving Existing Program --- 80
	7. User Selection and Registration --- 86
	8. Link Function --- 92
Chapter 5:	Diagnosis Functions and Other Tools --- 96
	1. History Function: History of Events --- 96
	2. Settings Function: Time, Date and System Default --- 97
	3. PC Communications Function: Program Send and Receive --- 98
	4. Diagnosis Function: Running Operation Testing Programs --- 101
	5. Delete User Data Function --- 105
Chapter 6:	Troubleshooting --- 107
Chapter 7:	Specifications --- 110
Appendix 1:	Accessories list and packing --- 111
	List of Supplies (PCR Enzymes) --- 112
Appendix 2:	Application 1-4 --- 113
Appendix 3:	Servicing Request and Safety Check --- 119

Chapter 1: Overview

TaKaRa PCR Thermal Cycler Dice™

TaKaRa Code: TP600 Gradient model

TP650 Standard model

Thank you for purchasing the TaKaRa PCR Thermal Cycler Dice™ system.

Until now, conducting polymerase chain reactions (PCRs) involved repeating experiments several times in order to achieve optimum reaction conditions for amplifying the test target.

TP600 Main Feature

The newly developed PCR Thermal Cycler Dice™ (TP600) system is equipped with a gradient function that enables the user to set a maximum 20°C temperature gradient between the left and right rows in the sample block at each temperature condition setting step. This has made it possible to optimize the reaction conditions for a single test and thereby design a new high-performance, functionally enhanced system at a lower cost.

< Enhanced Functionality and Ease-of-Operation >

- Gradient function enables consideration of annealing temperatures.

TP600/TP650 Features

The control panel and display units are integrated in a variable-angle wing panel configuration, enabling the creation of a more compact and stylish system that can be installed almost anywhere.

< Enhanced Functionality and Ease-of-Operation >

- Easy-to-operate slide and lever type lid.
- New variable-angle control panel enhances visibility and ease-of-operation.
- Graphic display makes programming the system much easier.
- Simpler operation through combination of pop-up menus and mobile phone-type alphanumeric character input pad.
- Extended time functions ideally suited to long PCR.
- Extended Temp function that makes touchdown PCR possible.
- Slop function (Dice version 2 and up)
- Link function that chain any User Programs and run. (Dice version 3)
- 8 types of preset protocol that make starting PCR easier.
- Modes that ensure compatibility (temperature rise characteristics) with earlier Takara Thermal Cycler models.
- Air intake vents on the front and sides and exhaust vents on the rear that reduce heat-based interference with other equipment.

< Ease-of-Maintenance >

- System firmware updates and data backup can be performed easily through connection to personal computer.
- Automatic creation of operating history and error logs.

Chapter 2: Installation Requirements and Maintenance

1. Installation Conditions

To ensure optimum performance and long and stable product life, observe the following precautions when selecting an installation location.

Note that installing the system in a location that does not meet the requirements listed below can result in malfunctions and/or system breakdown. Replacement of components and system repairs for any deterioration in performance or damage caused by failure to adhere to the following conditions will not be covered by warranty.

1. Install system in a comparatively dust-free location.
2. Avoid locations subject to corrosive or flammable gases.
3. Avoid locations subject to either strong shocks or continuing low levels of vibration.
4. Install the system away from direct sunlight.
5. Install the system as far away as possible from electrical equipment that generates strong magnetic fields, electrical fields and high frequencies. Avoid locations near equipment that generate noise caused by electrical motors and transformers.
6. Choose a flat location capable of withstanding the design weight of the system.
7. The installation environment must be within 15 to 30°C temperature range and a 20 to 80% relative humidity range.
8. Avoid locations subject to sharp temperature fluctuations, exhaust from air-conditioning systems and condensation caused by sharp humidity fluctuations.
9. Avoid blocking air vents located on the right and left sides, front and rear of the system.
10. Choose an operating location with good ventilation.
11. Do not operate the system in a cold room.
12. Select a room that can provide sufficient power (frequency, voltage and current) for system operation.
13. Do not attempt to move the system while holding it by the lid handle or the control panel, as this can damage to key components. Doing so can also result in personal injury due to part breakage resulting in the system falling down.

<Environmental Conditions>

This equipment is designed to be safe under the following conditions (based on the IEC 61010-1):

1. Indoor use
2. Altitude up to 2000 m
3. Ambient temperature (Operating environment): **15°C to 35°C**
4. Humidity range: **20% to 80%** relative humidity
5. Supply voltage fluctuations not to exceed $\pm 10\%$ of the nominal voltage.
6. Other supply voltage fluctuations as stated by the manufacturer

7. Transient over-voltages according to Installation Categories (Over-voltage Categories) II; For mains supply the minimum and normal category is II;
8. Pollution degree 2 in accordance with IEC 60664-3.

2. Installation

Use the following procedures when installing system.


1. Confirm that the system power switch is set to off.
2. Connect the system side of the power cable to the power connector on the rear of the system.
3. Plug the AC power plug of the cable into a power outlet.
The system is designed to operate normally with a 100 to 240 V AC, 50/60Hz power line.
Make sure that the power outlet is designed for a 3-pronged plug with ground line.

CAUTION

In the interests of safety, and avoiding electric shocks, be sure to connect the system to a properly grounded power outlet.

4. Turn on the system power.
5. The following input prompt screen appears as soon as power is applied.

3. Moving the System

Follow the procedures described below when moving the system to a new location after it has been installed.

- New Location Within the Same Site

Follow the precautions listed below when moving the system to a new location within the same site.

1. Thoroughly clean the system and make sure that there are no tubes or plates remaining in the sample block.
2. Avoid any shock during movement that can affect the operation of this highly sophisticated system.
3. When lifting and carrying the system, be sure to place one hand on either side of the system and lift and carry in a balanced way.
4. Make sure that the new location conforms to the installation requirements

described in section “1. Installation Conditions” in Chapter 2.

5. Follow the installation procedures outlined in section “2. Installation” in Chapter 2.

- Moving to a New Site

When moving the system to a new site, clean it carefully and repack it in the packaging in which it was delivered.

Note that damage sustained by the system during transport due to improper packaging will not be covered by Warranty even if the Warranty period is still valid.

Follow the instructions provided in this chapter when reinstalling the system at the new location.

4. Maintenance

Be careful especially of the following point before beginning maintenance check.

1. Especially the maintenance work without directions of operation should surely turn off a power supply switch for safety.
2. Warn to spill neither water nor a reagent to mechanism parts and electric parts. Moreover, please place neither detergent nor medicine on equipment.
3. Please do not touch with a hand other than the appointed place.

- Routine maintenance work


1. Please clean a main part periodically and suitably.
2. Wipe all outside surfaces of the equipment with a diluted neutral detergent solution.
3. When a liquid is spilt, please wipe off a sample block circumference part as early as possible.
4. Since dust is easy to collect, please clean a sample block frequently. If dust and garbage adhere to a sample block, a temperature control performance will fall and temperature homogeneity will also become bad simultaneously.
5. Please check it that the air hole of the back of equipment, the left and right-hand side is not closed. Please wipe off, when dust has adhered.

Neither strong alkali, high-concentration alcohol nor use of the organic solvent can be performed. It becomes the cause of breakage, such as discoloration, deterioration, and a crack.

Chapter 3 : Names of Parts and Their Functions

Front View


General Description

- (1) Lid Handle: Be sure to perform sliding operations with the lid open.
<Caution> Do not lift the unit using the lid handle.
- (2) Lid: The lid is designed to protect against evaporation resulting from heating of the reaction tube and reaction plate.
<Caution> Avoid inserting fingers as the interior of the unit reaches very high temperatures (up to 110°C).
- (3) Sample Block: Section of the unit in which sample tubes and plates are placed.
<Caution> Do not touch the block as it reaches temperatures up to 110 °C.
- (4) Display: Displays current status of the various system features and functions.

- (5) Key-Pad: Used to enter the various programs.
- (6) Control Panel: Enables adjustment of various angles (gradients).
- (7) Air Vents (Front): Openings used to cool the interior of the unit.
- (8) Control Panel Release Lever: Used for changing the angle of the control panel. Pressing upward releases lock.
- (9) Power Switch: Power switch is equipped with circuit breaker.
- (10) Air Vents, R and L : Openings used to cool the interior of the unit.
<Caution> Do not block any of the air vents.


Rear View


General Description

- (11) Power Cable Connector: Socket for power cable.
<Caution> Make sure to use the power cable provided with the unit.
- (12) Air Vent (rear): Openings used to cool the interior of the unit.
<Caution> Do not to block any of the air vents.
- (13) Brightness Adjust Knob: Used to adjust display brightness.
- (14) Communications Port Connector: RS232C connector(9Pin) used to link unit to computer.
<NOTE> The RS232C cable not included as standard accessory.

Key-Pad

Lid Section

Chapter 4: System Operation

1. Introduction

* This system operation explains both TP600 and TP650. Gradient Function is only for TP600.

* Operation Based on Pop-up Menu Configuration

Operation begins by pressing the **MENU** key, which displays the pop-up menu screens.

Menu items are highlighted using the arrow keys and selected by pressing the **ENTER** key.

Pop-up menus can be closed by simply pressing the **MENU** key.

A part of operation to which the pop-up menu is jumped over with the **ENTER** key and the **CE** key is possible.

* System Operational Flow

The basic flow of system operation is shown below.

1. Settings function
 - Selection of date/time and display format. (Section 2, Chapter 5)
2. Program design
 - User selection and registration (Section 7, Chapter 4)
 - Program creating or editing (Sections 5 and 6, Chapter 4)
 - Program link function (Section 8, Chapter 4)
3. Launching program
 - Operating program procedure. (Sections 2 and 3, Chapter 4)
4. Terminating or Aborting system operation
 - (Section 4, Chapter 4)
5. Maintenance and System tools
 - For description of inspection, diagnosis, backup and deletion procedures. (Sections 1, 4, 3 and 5, Chapter 5)

* Key Entry System

Notes that alphanumeric keys are assigned multiple characters, selectable by repeatedly pressing the key to cycle through the options.

- 0 key: &, /, (,), 0
- 1 key: space, -, ", #, 1
- 2 key: A, B, C, a, b, c, 2
- 3 key: D, E, F, d, e, f, 3
- 4 key: G, H, I, g, h, i, 4
- 5 key: J, K, L, j, k, l, 5
- 6 key: M, N, O, m, n, o, 6
- 7 key: P, Q, R, S, p, q, r, s, 7
- 8 key: T, U, V, t, u, v, 8
- 9 key: W, X, Y, Z, w, x, y, z, 9

2. Turning On the Power

The TP600/TP650 uses a 100V AC to 240V AC (non-switched), 50/60 Hz power line. Electric power is 490 VA and rated current is 4.9A when operating at 100 V.


[Note] Be sure to check local power line specifications prior to connecting the system.

[Note] Make sure that the power plug is inserted snugly into the power outlet.

1. Turn on the power switch.

The beeper sounds twice and the select-User-Name screen appear.

Using arrow keys, select the User Name, and push the **MENU** key.

[Note] When the pop up menu is active, the arrow key does not function.

2. The following password prompt screen appears if a password has been registered for the selected user.

Note that "Common" and "TaKaRa" have been pre-registered (default) user names that have not been assigned passwords.


[Note]

The above screen appears only when the username has been assigned a password. If not, the screen is skipped and the main system screen appears.

[Note]

When forget password, select "Cancel" from the MENU and . It is possible to go to the main system screen of 2-3 as "common" user.

3. After input the correct password, displays the main system screen shown below. All system operation begins with this screen.

[Note]

* Users who do not possess a password or have forgotten theirs may simply go to MENU, select Cancel and press **ENTER**. This logs the user onto the system as a Common user and displays the main system screen (described in step 3).

* Users logged on as TaKaRa are not able to perform program save operations or change username.


* There are eight standard protocols registered for the TaKaRa user name. Their program names are "< 1Kbp", "< 5Kbp", "< 10Kbp" and ">10Kbp", "Mutagenesis", "siRNA(double)", "Differential(DD)", "BcaBEST RNA".

(See "Chapter 4: [Note] Preset Protocols" on page 90)


3. Starting a System Operation Program

If a system operation program does not yet exist, begin by going to section "5. Creating New Programs" or "6. Editing and Saving Existing Programs."

1. With the main screen displayed, press the **MENU** key, select Run from the pop-up menu and press **ENTER**.

2. This displays the names of programs registered to the user chosen from the Select user screen.
Use the arrow keys to select the desired program, press the **MENU** key, then go to pop-up menu, select Start and press **ENTER**.

[Note]


When the pop up menu is active, the arrow key does not function. Press **MENU** key to disappear the pop up menu and select program by Using arrow key.

[Note]

Other functions in pop up menu


(1) Checking Program Information

To confirm the program information, select View from the pop-up menu and press **ENTER**.

(2) Sorting Program Names

To change program name order, select Sort and press **ENTER**.


Select from the following sort options.


- Program name: sorts alphabetically by name (in ASCII code)
- Date saved: sorts in order of date saved
- Date last used: sorts in order of date used

Then go to pop-up menu, select OK and **ENTER** key.

3. Select or directly Reaction mode.

When the program that adds the link function is selected, it becomes the following displays.

Link: Select 1 or 2 for Link function.

[1.Yes]: When operating the program by using the link function.

[2.No]: When operating only one program.

Please refer to the paragraph of the link function for details.

- * Mode 1: Enables system operation with maximum temperature control.
Recommended for normal operation.
- 2: Enables system operation with high temperature control.
Compatible with TaKaRa PCR Thermal Cycler SP operation.
- 3: Enables temperature control compatible with linked compressor unit.
Compatible with TaKaRa PCR Thermal Cycler MP operation.
- 4: Fast speed turning for disregarding overshoot temperature.
- 5: Enables system operation with high temperature control.
Semi-Compatible with TaKaRa PCR Thermal Cycler PP normal mode operation.

Mode performance outlook:


Ramp Speed: Hi 5 > 4 > 1 > 2 > 3 Low

- * Container (1.Tube 2.Plate) Select desired container.
TaKaRa PCR Tube of 0.2ml Hi-Tube with Dome Cap (Code NJ200) is a thin-wall type tube. Select "1" when using this tube.
- * Volume 5 -150 micro litter: Select reaction volume (Usually, below 100micro litter)


[Note] The above information settings are saved separately for each program. When the program is started, each setting condition is memorized.

- 4 . Select OK from the pop-up menu and press **ENTER**.
- 5 . This displays the program pattern screens shown below.
Note that the patterns show a program that has not yet started.

Extension Data Not Selected

Extension Data Selected

Extension data shown here

[Note]

The Extension mode allows the user to temperature, time increments/decrements, pause and gradient settings.


[Note]

Program names are registered using up to 16 characters, however, the pattern screens can only display 12 characters.

6. Go to the pop-up menu, select Start and press **ENTER**.

This starts the program.

Program-Operating Screen

Changing-temperature: diagonal line section of graph flashes.

Stabilized-temperature: horizontal line flashes.

Temperature indication: "HOT!" is flashing when over 50 °C.

Time remaining for current step: counts down


Definition of ∞: continuous operation (99:99)

[Note]


The "Lid is warming up xx.x °C" message appears if it has not reached 108°C over temperature. The System starts as soon as the lid is warm.

(1) Pausing When System Operation

Go to pop-up menu, select Pause and press **ENTER**.


To resume system operation, select resume from the pop-up menu and press **ENTER**.

(2) Viewing Program Information During Operation


While program is running, go to pop-up menu, select Info and press **ENTER**.


To close Info screen, go to pop-up menu, select OK and press **ENTER**.

4. Terminating or Aborting Operation

1. Press **Stop** key.
2. This displays the following confirmation screen.

3. Go to the pop-up menu, select OK and press **ENTER**.
4. This displays the following termination screen.

Explanation of content of display

Status Program Name

User Name

Program beginning time (hour: minute: second) and the date.

Program end time (hour: minute: second) and the date.

Operation time(hour: minute: second)

The maximum overshoot and maximum undershoot program's operating


Result

5. Going to the pop-up menu and selecting OK again displays the main program screen.


5. Creating New Program

New programs are created by going to the top screen selecting and editing a basic program pattern (format), then assigning it a name. When creating a new program from an existing user file or simply editing the file, refer to the instructions provided in section 7, "User Selection and Registration."

1. From the main program screen, go to the pop-up menu, select New and press **ENTER**.

2. This displays the New Program screen. Programs are created by editing segments (Temp count) and numerical values.

Definitions of Terms

Segment: Represents the smallest unit which is equivalent to 1 Temp unit.

Pattern: Represents a combination of segments for which numerical values are continuously set.

Procedure for setting Values

Use arrow keys to move cursor.

Press **ENTER** or an arrow key after input numerical values of temperature (Temp), time (time) and cycle count (Cyc).

Temperature using numerals: For example, for 94.0°C, "940".

Temperature setting range: 4.0 to 99.9°C


Time using numerals: For example, for 1 minute and 30 seconds, "0130".

Time setting range: 00:01 (s) to 99:59 (99 minutes, 59 seconds);
when operating in continuous (∞), "9999".


Cycle count: 1 to 99

Procedure for Inserting/Deleting Patterns and Segments

· Patterns can be inserted or deleted by placing the cursor on the pattern start position (Example: placing cursor on “1”, “94.0” or “10:00” of “1 Temp” or “3”, “94.0” or “0:30” of “3 Temp”) and selecting Insert/Delete from the pop-up menu.

* Segments can be inserted or deleted by placing cursor anywhere except the pattern start position (Example: placing cursor on “55.0/0:30” or “72.0/1:00” for “3 Temp 25 Cyc” on the screen below) and selecting Insert/Delete from the pop-up menu.

* Segment (Temp) count in pattern is modified by inputting a new value (1 to 9) in the first Temp position in pattern and pressing .

Max. segment (Temp) count per program: 12

[Note] Using a segment (Temp) count of 7 or more requires the double memory capacity as two 6-segment (Temp) programs.)

[Note] When the pop-up menu is "MENU2", it is possible to switch to "MENU1" by pushing a right and left arrow key.

Using Extension Mode

Extension mode provides the following functions. But model TP650 does not have the Gradient Function.

- **Extension Settings**
 Temperature and time settings must be increased/decreased each time cycle count is increased.
Temperature extension setting is changed by enabling Extension mode, moving the cursor to the TempEX position and input a value. Go to MENU, select “+/-” and press **ENTER**.
Time extension setting is changed by moving the cursor to the TimeEX position and input a value. Go to MENU, select “+/-” and press **ENTER**.
- **Pause Setting**

With Pause activated, the program pauses when operation moves to the selected pause segment. Program operation is resumed by going to MENU, selecting Resume and pressing.

[Note]

Pause is not designed to activate after the selected segment has been executed.

With the Extension mode enabled, go to MENU, select **Pause** and press **ENTER** to activate **Pause**. To cancel pause, go to MENU, select **Pause** and press again.


- **Slope Setting**

This is used to set the desired temperature heating/cooling speed setting. Do this by setting the desired target temperature and ramp time settings. The resulting temperature change slope will then be used to perform linear temperature control. Move cursor to **TempEX**, go to MENU, select **Slope** and press **ENTER**.


- **Gradient** (TP600 only)

With the cursor positioned on **Temp** or **Time**, go to MENU, select **Grad** and press Enter to select a setting.

- (1) Enable Extension mode by selecting More from pop-up menu and pressing .

- (2) Extension and Pause Settings

Position cursor at the “Ex” line below the Temp or Time lines, desired value and press **ENTER**. Then select increase (+), decrease (-) or Pause from pop-up menu and press again.

Max. temp. Ex value: ±9.9°C

Max. time Ex value: ±9:59 (9 minutes 59 seconds)

Pause setting by moving cursor to TimeEx of the target segment, select Pause in pop-up-up menu and press **ENTER**.


The figure below shows an example of a program setting procedure using 55°C and 30 s as the temperature and time settings, with Time setting decremented by 5 s each cycle.

The Pause is set for activation when the 55°C segment is completed and operation moves

to 72°C.

Procedure

- 1) Use the arrow keys to activate (highlight) the selected position.
- 2) Use numerical pad keys to "5".
- 3) Press **MENU** key to display pop-up menu then use arrow key to select (-).
- 4) Press **ENTER** key.
- 5) Use the arrow keys to activate (highlight) the selected position.
- 6) Press **MENU** key to display pop-up menu then use arrow key to select (Pause).
- 7) Press **ENTER** key.

[Note]

The extension setting can be canceled by pushing **CE** key or entering "0".

The Pause setting can be canceled by pushing **CE** key or going to MENU and selecting Pause again.


[Note] Decreases every five seconds but It doesn't become less than one second.

[Note] Pop-up menu content varies according to cursor position.

-The pop-up displays Save, Insert and Delete when the cursor is positioned at Temp or Time.


- Save, +, -, **Slope** and Cancel are displayed when positioned at TempEx.

- Save, +, -, **Pause** and Cancel when positioned at TimeEx.

(3) Slope Function Settings

With the cursor placed on TempEx, go to MENU, select Slope and press **ENTER**. This displays the following Slope screen.

Slope Setting Range

Temperature: 30.0 to 99.9°C

Time: 1:30 (1 minute 30 seconds) to 99:59 (99 minutes 59 seconds)

[Note] The Minimum Time Setting for Slope is 1 min and 30 s. Any setting shorter than the minimum setting will automatically be adjusted to this value.

[Note] The Slope setting is canceled by moving the cursor to the designated position and pressing the **CE** key or going to MENU, selecting Slope and pressing again. These steps reset the current slope setting to its default value (0.0).

(4) Gradient Function Settings (only for TP600)

With the cursor positioned on Temp or Time, go to MENU, select Grad and press **ENTER** to select a setting. It becomes Gradient temperature setting screen.

	1 Temp	3 Temp	25 cyc	2 Te	MENU1
					Save
					Insert
Temp	94.0	94.0	55.0	72.0	72.0
Ex	0.0	0.5	0.0	0.0	0.0
Time	10:00	0:30	0:30	0:30	5:00
Ex	0:00	0:00	0:00	0:00	0:00
					Grad.
					Cancel

It becomes Gradient temperature setting screen.

Gradient Setting			MENU1
(01) 45.0	(02) 46.0	(03) 47.0	OK
(04) 49.0	(05) 51.4	(06) 53.8	Clear
(07) 56.2	(08) 58.6	(09) 61.0	
(10) 63.0	(11) 64.0	(12) 65.0	

the low-temperature setting at position (01: The left end of a block) and the high-temperature setting at position (12: The right end of a block).

Temperature setting range: 40.0 to 75.0 °C

Temperature difference: 6.0 to 20.0 °C

[Note]

- * Exceeding the temperature settings or the temperature setting range may destabilize controlled temperature.
- * The temperature values for (02) to (11) are calculated values.
- * Temperatures cannot be automatically incremented/decremented for each step.

Go to pop-up menu, select OK and press **ENTER**.

When displaying the temperature gradient, the left temperature is shown on the Temp line and the right temperature (beginning with G) on the Temp Ex line.

Go to pop-up menu, select Save and press **ENTER**.

	1 Temp	3 Temp	25 cyc	2 Tem	
Temp	94.0	94.0	G45.0	72.0	72
Ex	0.0	0.5	G65.0	0.0	0
Time	10:00	0:30	0:30	0:30	5:
Ex	0:00	0:00	0:00	0:00	0:

MENU1

Save

Insert

Delete

Grad.

Cancel

3. Save the newly created program.

Entering program name.

Save Program as

Program **TaKaRa Taq**

User (Common)

MENU

Save As

User

Cancel

Go to pop-up menu, select SaveAs and press **ENTER**.

[Note]

User-created programs: maximum memory capacity for programs (6-Temp size) is **200**

Program name character: Maximum **16**

Max. number of registered users **22** excluding pre-registered <Common> and <TaKaRa> usernames.

[Note]

New programs cannot be registered to the TaKaRa folder, as the Save and SaveAs functions are disabled (blacked out as ---) in the pop-up menu. Select User in the pop-up menu and save the program under a different username.

Save Program

Program **TaKaRa Taq**

User(TaKaRa)

MENU


User

Cancel


6. Editing and Saving Existing Programs

Existing programs can be saved under a new name (newly created program) or under the same name (edited program).

1. Go to the top screen, select Edit and press **ENTER**.

2. This displays a list of programs registered to the User. Choose the program to be edited, select Edit in the pop-up menu and press **ENTER**.


This displays the existing program and enables editing of the segment (Temp count) and numerical values.

Use the arrow keys to move the cursor.


- * Press or an arrow key after setting temperature (Temp), cycle count (Cyc) and numerical values.

Temperature setting: 4.0 to 99.9 °C
 Time settings: 00:01(1 second) to 99:59(99minutes 59seconds)
 99:99 (when operating in continuous (∞))
 Cycle count: maximum 99

- * Blocks can be inserted or deleted by going to the start position (Example: placing cursor on "1" of "1 Temp" block) and selecting Insert/Delete from the pop-up menu.

- * Segments can be inserted or deleted by placing cursor anywhere except the pattern start position (Example: placing cursor on "55.0/0:30" or "72.0/1:00" for "3 Temp 25 Cyc" on the screen below) and selecting Insert/Delete from the pop-up menu.

- * Pattern (Temp) count can be modified by inputting a new value (1 to 9) at the first Temp position and pressing.
 Max. Temp count per program: 12
 (Using a Temp count of 7 or more requires as much memory as two 6-Temp programs.)
- * Segments (Temp count) can be added or deleted by selecting Insert/Delete from the pop-up menu while the cursor is positioned anywhere except at the start of a block.


Using Extension Mode

Extension mode provides the following functions but model TP650 does not have the Gradient Function.


- **Extension Settings**
 Temperature and time settings must be increased/decreased each time cycle count is increased.
Temperature extension setting is changed by enabling Extension mode, moving the cursor to the TempEX position and input a value. Go to MENU, select “+/-” and press **ENTER**.
Time extension setting is changed by moving the cursor to the TimeEX position and input a value. Go to MENU, select “+/-” and press **ENTER**.
- **Pause setting**
 With Pause activated, the program pauses when operation moves to the selected pause segment. Program operation is resumed by going to MENU, selecting Resume and pressing
[Note]
Pause is not designed to activate after the selected segment has been executed.
 With the Extension mode enabled, go to MENU, select **Pause** and press **ENTER** to activate **Pause**. To cancel pause, go to MENU, select **Pause** and press again.
- **Slope Setting**
 This is used to set the desired temperature heating/cooling speed setting. Do this by setting the desired target temperature and ramp time settings. The resulting temperature change slope will then be used to perform linear temperature control. Move cursor to TempEX, go to MENU, select **Slope** and press **ENTER**.
- **Gradient** (TP600 only)
 With the cursor positioned on Temp or Time, go to MENU, select **Grad** and press

Enter to select a setting.

- (1) Enable Extension mode by selecting More from pop-up menu and pressing .

- (2) Extension and Pause Settings

Position cursor at the “Ex” line below the Temp or Time lines, desired value and press **ENTER**. Then select increase (+), decrease (-) or Pause from pop-up menu and press again.

Max. temp. Ex value: $\pm 9.9^{\circ}\text{C}$


Max. time Ex value: $\pm 9:59$ (9 minutes 59 seconds)

Pause setting by moving cursor to TimeEx of the target segment, select Pause in pop-up-up menu and press **ENTER**.

The figure below shows an example of a program setting procedure using 55°C and 30 s as the temperature and time settings, with Time setting decremented by 5 s each cycle. The Pause is set for activation when the 55°C segment is completed and operation moves to 72°C.

Procedure

- 1) Use the arrow keys to activate (highlight) the selected position.
- 2) Use numerical pad keys to “5”.
- 3) Press **MENU** key to display pop-up menu then use arrow key to select (-).
- 4) Press **ENTER** key.
- 5) Use the arrow keys to activate (highlight) the selected position.
- 6) Press **MENU** key to display pop-up menu then use arrow key to select (Pause).
- 7) Press **ENTER** key.


[Note] The extension setting can be canceled by pushing **CE** key or entering "0", and the Pause setting is cancelled by going to MENU and selecting Pause again.

[Note] Pop-up menu content varies according to cursor position.

-The pop-up displays Save, Insert and Delete when the cursor is positioned at Temp or Time.


- Save, +, -, **Slope** and Cancel are displayed when positioned at TempEx.

- Save, +, -, **Pause** and Cancel when positioned at TimeEx.

(3) Slope Function Settings

With the cursor placed on TempEx, go to MENU, select Slope and press **ENTER**. This displays the following Slope screen.

Slope Setting Range

Temperature: 30.0 to 99.9°C


Time: 1:30 (1 minute 30 seconds) to 99:59 (99 minutes 59 seconds)

[Note] The Minimum Time Setting for Slope is 1 min and 30 s. Any setting shorter than the minimum setting will automatically be adjusted to this value.

[Note] The Slope setting is canceled by moving the cursor to the designated position and pressing the **CE** key or going to MENU, selecting Slope and pressing again. These steps reset the current slope setting to its default value (0.0).

(4) Gradient Function Settings (only for TP600)

¥ With the cursor positioned on Temp or Time, go to MENU, select Grad and press **ENTER** to select a setting. It becomes Gradient temperature setting screen.

It becomes Gradient temperature setting screen.

Gradient Setting			MENU 1
(01) 45.0	(02) 46.0	(03) 47.0	OK
(04) 49.0	(05) 51.4	(06) 53.8	Clear
(07) 56.2	(08) 58.6	(09) 61.0	
(10) 63.0	(11) 64.0	(12) 65.0	

the low-temperature setting at position (01: The left end of a block) and the high-temperature setting at position (12: The right end of a block).

Temperature setting range: 40.0 to 75.0 °C
 Temperature difference: 6.0 to 20.0 °C

[Note]

- * Exceeding the temperature settings or the temperature setting range may destabilize controlled temperature.
- * The temperature values for (02) to (11) are calculated values.
- * Temperatures cannot be automatically incremented/decremented for each step.

Go to pop-up menu, select OK and press **ENTER**.

When displaying the temperature gradient, the left temperature is shown on the Temp line and the right temperature (beginning with G) on the Temp Ex line.

Go to pop-up menu, select Save and press **ENTER**.

	1 Temp	3 Temp	2 5 Cyc	2 Te	MENU 1
Temp	94.0	94.0	G45.0	72.0	Save
Ex	0.0	0.5	G65.0	0.0	Insert
Time	10:00	0:30	0:30	0:30	Delete
Ex	0:00	0:00	-0:05	Pause	Grad.
					Cancel


3. Saving an Editing Program.

There are two procedures for saving an existing program.

- 1) Saving the edited program under the existing name. (Save)*1
- 2) Saving the edited program under a new program name. (Save)*2
- 3) Saving the edited program under a new program name. (SaveAs)*3

Go to pop-up menu, select Save/SaveAs and press **ENTER**.

	1 Temp	3 Temp	2 5 cyc	2 Te	MENU 1
Temp	94.0	94.0	55.0	72.0	Save
Ex	0.0	0.5	0.0	0.0	+
Time	5:00	0:30	0:30	1:00	-
Ex	0:00	0:00	-5:00	Pause	Pause
					Cancel

- *1) Press **ENTER** under the existing name. Or select Save from the pop-up menu and press **ENTER**. :Overwrite (Save)
- *2) Press **ENTER** under the new name. Or select Save from the pop-up menu and press **ENTER**. :Saved but the former program is erased.
- *3) Press **ENTER** under the new name. Or select SaveAs from the pop-up menu and press **ENTER**. :The former program remain.

[Note] The message of "Change Program Name" appears when SaveAs is selected without changing the name and is done, and the program name change is demanded. Please input a new name.

[Note] Programs cannot be saved to the TaKaRa username using the Save or SaveAs functions. Select User from the pop-up menu and save the program under a different username.

[Note] User-created programs: maximum memory capacity for programs (6-Temp size) is **200**.

Program name character: Maximum **16**.


Maximum number of registered users **22** excluding pre-registered <Common> and <TaKaRa> usernames.

7. User Selection and Registration


Options include 1. Select user-name (User), 2. Register new user (New) and 3. Change user name (Editor).

1. Select user-name

- (1) Go to main system screen pop-up menu, select User and press ENTER.

- (2) This displays the Select Username screen


<NOTE>

- * This screen appears at system startup when the power is turned on.
- * The defaulted usernames <Common> and <TaKaRa> are pre-registered.
- * Maximum number of user-names: 22
- * Maximum user-name length: 10 characters

Use the arrow keys to select the desired username.
Select OK on the pop-up menu and press ENTER.

- (3) Enter password when the password entry screen appears.


Password length: 0 to 9 characters that appear as "*" when entered
Use the CE key to clear password:

If the password is mistyped, use the CE key to clear and re-enter password.


- (4) Select OK on the pop-up menu, press ENTER and return to the top screen.

2. Registering New User

- (1) Go to main system screen pop-up menu, select User and press ENTER.

- (2) Select New on the pop-up menu and press ENTER.

- (3) Enter the new user name, go to pop-up menu, select Save As and press ENTER.

- (4) Enter the password if necessary.

Enter password.

Password characters and length: numbers 0 to 9 (1 to 6 digits)


This displays the entered numbers on the screen.

If the numbers are mistyped, use the CE key to clear and re-enter password.

To avoid registering a password, simply leave the field blank, select OK from the pop-up menu and press ENTER.

3. Changing Existing Username

- (1) Go to pop-up menu, select User and press ENTER.

- (2) Select the user name to be changed, go to pop-up menu, select OK and press ENTER.


When the password is set, the password is demanded. Enter password when the password entry screen appears and return to the top screen.


- (3) Confirm the user name you desired. Go to pop-up menu, select User and press ENTER.

- (4) confirm the user name that wants to change reverses by the User selection screen. Go to pop-up menu, select Edit and press ENTER.

- (5) Enter a new use name. Go to pop-up menu, select Save and press ENTER.

- (6) Enter password when the password entry screen appears. Select OK on the pop-up menu, press ENTER and return to the top screen.

Password characters and length: numbers 0 to 9 (1 to 6 digits)

This displays the entered numbers on the screen.

Use the CE key to clear password:

If numbers are mistyped, use the CE key to clear and re-enter password.

To avoid registering a password, simply leave the field blank, select OK from the pop-up menu and press ENTER.

To confirm the User name was changed, main system screen is displayed.

<NOTE>

Username Edit Functions

- Up arrow key : insert space
- Down arrow key : backspace
- Left arrow key : cursor left
- Right arrow key : cursor right
- Character input keys
 - 0 key: &,/,(),0
 - 1 key: Space,-, ",#,1
 - 2 key: A,B,C,a,b,c,2
 - 3 key: D,E,F,d,e,f,3
 - 4 key: G,H,I,g,h,i,4
 - 5 key: J,K,L,j,k,l,5
 - 6 key: M,N,O,m,n,o,6
 - 7 key: P,Q,R,S,p,q,r,s,7
 - 8 key: T,U,V,t,u,v,8
 - 9 key: W,X,Y,Z,w,x,y,z,9

Enter only numbers for number fields.

<NOTE>

* Pre-registered user-names <Common> and <TaKaRa> may not be changed.
* The following protocols are registered to the <TaKaRa> folder: "< 1Kbp", "< 5Kbp", "< 10Kbp", ">10Kbp", "Mutagenesis", "siRNA(double)", "Differential(DD)" and "BcaBEST RNA".

* Preset Protocols: These consist of programs (preset protocols) programmed for use according to the length of the target.

< 1 Kbp: (up to 1 Kbp) 94°C for 30 s, 55°C for 30 s and 72°C for 1min at 30 cycles, 4°C/

< 5 Kbp: (up to 5 Kbp) 94°C for 30 s, 55°C for 30 s, 72°C for 5min at 30 cycles, 4°C/

< 10 Kbp: (up to 10 Kbp) 98°C for 10 s and 68°C for 10 min at 30 cycles, 4°C/

> 10 Kbp: (more than 10 Kbp) 98°C 10 s and 68°C 15 min at 30 cycles, 4°C/

* Preset Protocols(with Slope function: See next page.):

Mutagenesis: 94°C/10min, 94°C 37°C/60min and 37°C/15min< 1 Kbp

siRNA(double): 95°C/1min, 70°C/1min, 70°C 37°C/60min and 37°C/


Differential(DD):94°C/30s, 40°C/1min, 40°C 72°C/2min and 72°C/30s

BcaBEST RNA: 30°C/1min, 30°C 65°C/15min and 65°C/20min

Mutagenesis

94°C/10min, 94°C 37°C/60min and 37°C/15min


Mutagenesis

siRNA(double)

95°C/1min 70°C/1min, 70°C 37°C/60min and 37°C/


siRNA(double)

Differential(DD)

94°C/30s, 40°C/1min, 40°C 72°C/2min and 72°C/30s


Differential (DD)

BcaBEST RNA

30°C/1min, 30°C 65°C/15min and 65°C/20min

BcaBEST RNA

8. Link Function

The Link function is used to execute multiple programs consecutively. The Cycler Dice is capable of inputting a maximum of 12 segments with one program. Users wishing to execute even larger programs have the option of running multiple programs consecutively.

Program names are assigned the “#” symbol followed by a single digit number. The Link function executes programs consecutively in the order of the assigned numbers. The Cycler Dice is not optimally suited to operating in a modular configuration in which the required modules are linked together.

Creating Link Programs

A single program name is assigned the “#” symbol followed by single digit numbers at the end.

Example: PLANT_TEST#1
 PLANT_TEST#2
 PLANT_TEST#3
 PLANT_TEST#4

↓ Programs numbered #1 to#4
will be executed in that order.

[Note] The Link function ignores all numbers except the first digit following the “#” symbol.

(Example) LINK#1
 LINK#2 Long
 LINK#3 Short
 LINK#4


↓ Execution commences with
Program #1 and continues to #4.

[Note] Avoid Program Names with Same First Digit

Although possible to save programs to the system that have the same first digit after the “#” symbol if they have differing subsequent digits, they will be regarded as identical by the Link function. They will be executed in the order of their system memory addresses, rather than when created or how named, resulting in confusion.

Example: Program Names


LINK#1
LINK#2
LINK#2 Long
 LINK#3 Short
 LINK#4

Link Function Program Execution Procedure

1. Go from program username screen to main screen and press **MENU**. Then select Run from pop-up menu and press **ENTER**.

Example: Program Created for Username <Common>

In MENU select OK and press **ENTER**.


2. Select first program.

Example: PLANT_TEST#1
 PLANT_TEST#2
 PLANT_TEST#3
 PLANT_TEST#4

Programs numbered #1 to #4 will be executed in that order.

Select "PLANT_TEST#1", go to OK and press **ENTER**.

3. Go to pop-up menu, select Start and press **ENTER** to display Select Reaction Mode screen.


4. Enter desired reaction settings, confirm "1" is selected for Link field (line 4) then choose OK in MENU and press **ENTER**.
5. Next, go to MENU, select Start and press **ENTER** to start execution of program PLANT_TEST#1.
 Execution of PLANT_TEST#2 begins automatically as PLANT_TEST#1 completes.
 Execution continues with PLANT_TEST#3 and PLANT_TEST#4 and terminates

[Note] Execution can skip the beginning programs and start with PLANT_TEST#2, for example, and continue through PLANT_TEST#3 and PLANT_TEST#4.

[Note] If a gap occurs in the order, execution terminates with the last program before the gap.


Example: Execution that begins with TEST#1 terminates with TEST#2.

Executing Programs Intended for Link Function Individually

Programs created specifically for Link execution can also be executed separately.

1. Select program name, go to MENU, select Start and Press **ENTER** to display the Select Reaction Mode screen.
2. Enter "2" for Link then select OK in MENU and press **ENTER**.

Selecting "2. No" for Link cancels the Link function.

3. Going to MENU, selecting Start and pressing **ENTER** executes the program until it is finished.

[Note] Display on running with Link function.

- 1) Program Execution Screen
 Program Title displays name of program in progress. Remaining Time indicates time remaining for execution of current program and Info displays the starting time and expected finishing time of current program in progress.
- 2) Finished Program Screen
 The Finished screen appears program name, program execution starting and completion time and length of run-time of the only last program.
- 3) Program History Screen
 Finished program name is comprised of the name before the "#" and [#starting program number] – [#finishing program number].

Example: LINK#1
 LINK#2 Long
 LINK#3 Short
 LINK#4

When the Link programs finished History screen registered a program name as "LINK#1 - #4".

History 2/2		
User:	<Common>	
Program:	LINK#1-#4	
Mode :	1 Thin 50 µl	
Started at	11:50:10	05/01/21
Finished at	15:28:15	05/01/21

[Note] To view the program history, go to the main system screen, select Tools from the pop-up MENU and press **ENTER**. Next, select History then choose OK in the MENU and press **ENTER**.
 (Refer to section "1. Saving History")

"Started at" indicates the time that the first program commenced and "Finished at" the time the last program was completed.


Chapter 5: Diagnosis Functions and Other Tools

The following functions are provided in the Tools menu.


- 1) History: records history of activity
- 2) Settings: enters time and date and system default settings.
- 3) PC communications: sends and receives programs
- 4) Diagnosis: runs operation testing programs
- 5) Delete user data: deletes saved usernames and programs

Selecting Tool Menu

Go to main system screen pop-up menu, select Tools and press ENTER.

Choose desired item from the displayed Select Tools menu, select OK from the pop-up menu and press ENTER.


<NOTE>

The descriptions of the various functions provided below assume that they have already been selected from the Select Tools menu (by selecting OK in pop-up menu and pressing ENTER).

1. History Function: History of Events

The History function records a log of activity for the most recently executed program. Executing a new program updates the history log by overwriting previous data. Select History, choose OK from the pop-up menu and press ENTER.

The number of recorded events is shown to the right of "History" on the top line. The example shown on screen indicates that information is displayed for 3 of the past 4 events (and does not indicate page count). More events can be viewed by pressing the Down arrow key.

Maximum recorded events: 40

<NOTE>

To display the history of continuation, the pop-up menu is erased and Push down arrow key or pushing the MENU key twice.

When the denominator of the fraction on the History side is four or more on the history screen, the history is displayed on the following screen. Please do not forget and look.


<NOTE>

History display when link function is used

(example)

- LINK_#1
- LINK_#2Long
- LINK_#3Short
- LINK_#4

The history record is displayed as LINK_#1-#4 when running from LINK_#1 to LINK_#4.

- The start time of the first program is displayed as start time.
- The finish time of the last program is displayed as a finish time.

2. Settings Function: Time, Date and System Default


The Settings function includes time, date, system standby temperature, display screen color and highlighting, system sleep time and sound settings.

2-1 Time and Date Settings

1. Choose Setting, then go to pop-up menu, select OK and press ENTER.
2. Choose Date & Time, then go to pop-up menu, select OK and press ENTER.


3. Enter display format and date/time settings, then go to pop-up menu, select OK and press ENTER.

For example, the time (14:15:00) and date (January 15, 2003) are entered as follows.
 Date Style = [1], Date = [030115], Time = [141500]
 Then go to pop-up menu, select OK and press Enter key.

2-2 Default Settings

1. Choose Setting, then go to pop-up menu, select OK and press ENTER.
2. Choose Default Setting, then go to pop-up menu, select OK and press ENTER.

3. Enter the desired settings, go to pop-up menu, select OK and press ENTER.
 - * Stand-by Temperature: block temperature setting during standby mode (default is 25 °C)
 - * LCD Back Color: screen background settings (1: blue, 2: white), character color complements background color
 - * LCD Sleep: enables system sleep mode (shuts off display backlight) when there has been no input while program is idle for a specified interval. (after 1:30 or 2:00 minutes:seconds)
 The LCD is designed to wake up by pressing MENU.
 Will not wake up with other key operation.
 - * Buzzer: Key Pad: pressing key toggles options (1: beep on, 2: beep off).
 - * Buzzer : Finished: selects between options (1: beep on, 2: beep off)
 (Beep count: 3 at program pause, 6 at program finish)

3. PC Communications Function: Program Send and Receive

This function is used to save program data to a computer or send it from a computer to other devices. It is also used to perform system backup and transfer operations.

3-1 PC Communication Settings (Windows 2000 and XP Version)

1. From Start Menu
 Go to Programs, Accessories, Communications and open HyperTerminal.
 If not included in the Start menu, open HyperTerminal by going to the following directory:
 C:/Program Files/Windows NT/hypertrm.exe
2. HyperTerminal Screen
 New Connection Name: Example: TaKaRa
 Connect using: COM1
 COM1Properties → Port Settings: Bits per second: 9600, Data bits: 8, Parity: Stop bit,
 Flow control : Xon/Xoff
 ASCII I Setup : go to File, Properties and select the Settings tab.
 Click ASCII Setup button and enter the following settings:
 Line delay: 2 milliseconds, Character delay: 0 milli seconds
 Uncheck all other options.

When Receiving: go to Transfer and select Capture text, then specify file name and click the Start button.

When sending: go to Transfer , select Send Text File and select the send file name.

<Important Notes>


The RS232C cable with 9-pin D-sub connector required for this feature is not included as a standard accessory.

3-2 Transferring a Program to Computer (**Dice PC**)


1. Choose PC Communication, go to pop-up menu, select OK and press ENTER.


2. Choose Send a program to PC, go to pop-up menu, select OK and press ENTER.

3. This displays the program list.
Choose program to send, go to pop-up menu, select OK and press ENTER.

4. On the PC, use HyperTerminal.
In the HyperTerminal screen, to go to Transfer and select Capture, then save the file by clicking the Start button.
For a description of the PC settings, refer to “PC Communications Settings.”
5. Transfer the desired program.

Go to pop-up menu, select Start and press ENTER.


The “Finished” message appears when transfer is complete.

Select Cancel from the pop-up menu, then press ENTER to be returned to the main


menu screen.

3-3 Transferring Programs From a Computer (PC Dice)


1. Choose PC Communication, go to pop-up menu, select OK and press ENTER.

2. Choose Receive a Program from PC, go to the pop-up menu, select OK and press ENTER.


3. This displays the Receive a Program screen.
Go to pop-up menu, select Start and press ENTER.
Confirm that the "Receiving ..." message appears.

4. On the PC, use HyperTerminal.
Go to Transfer, choose Send Text File and select the send file name.
Send the program from the PC to Dice.


<NOTE> For the correct PC transmission settings, refer to "PC Communication Settings" in section 3-1, "Transferring a Program to a Computer."

5. Information for the received program is displayed when the process is completed.

Go to pop-up menu, select OK and press ENTER.

- Save program under new file name.

<NOTE> For a description of the file saving procedure, refer to “3. Saving New Programs” in section 5 of Chapter 4.

4. Diagnosis Function: Running Operation Testing Programs

This function performs self-diagnosis tests.


These tests are used to check whether the system is operating normally when errors occur in PCR reactions.

The Diagnosis function performs the following four tests.


- (1) Rate Test
Tests heating and cooling rates.
- (2) Cycle Test
Tests cycle times and stability of control operations.
- (3) Up-Down Test
Measures time required to cool down to 4 °C.
- (4) Gradient Test (model TP600 only)
Measures gradient temperatures.

Confirm that standby temperature is set at 25°C prior to commencing tests.

- Choose Diagnosis, go to pop-up menu, select OK and press ENTER.

- This displays a list of self-diagnosis tests.
Choose the desired test item, then go to pop-up menu, select OK and press ENTER.

<NOTE>

Choosing the target test item, then selecting Result from the pop-up menu and pressing the ENTER. Displays the last results for each test item.

- When the diagnosis test screen appears, select Start from pop-up menu and press ENTER.

Test results are displayed when the tests are completed.

Reference values are for ambient temperature in the 20 to 23 °C range.


(1) Rate Test program execute screen

Program details:

25 °C for 30 seconds, 94 °C for 30s and 40 °C for 30 seconds.

Results Displayed


Reference Values: Max. Heating Rate >3.0°C/second


Max. cooling rate >2.0°C/second

(2) Cycle Test Program execute screen

Program details: 94°C for 30 seconds, 55°C for 30 seconds and 72°C for 30 seconds.


Results display


Reference Value: overshoot < 0.5°C
undershoot < 0.5 °C
cycle time < 4 minutes

- (3) Up-Down Test program executes cream
Program Details: 98 °C for 1 minute and 4 °C for 1 minute.


Results displayed


Reference Value: cycle finish time < 10 min

- (4) Gradient program test (model TP600 only)
Program Details: 94 °C for 30 seconds, G50/G60 for 5 min and 72 °C for 1 minute.

Gradient items are displayed in the menu during gradient control, enabling monitoring of the temperatures in each row.
After the timer has counted for 2 or 3 min during the gradient temperature program, selecting OK for Grad from the pop-up menu displays the following monitoring screen.
(This will not be saved.)

Results are displayed at the center of the sample block.

Result of Grad. Test		MENU
Heat 3.1°C/s	Cool 2.0°C/s	OK
Over 0.2°C	Under 0.5°C	
Started at 16:21:42 03/01/15		
Finished at 16:31:05 03/01/15		
Length of Run 00:09:23		
Result: Normal		

Reference value: < +/- 1 °C by 2 to 3 minutes after Timer begins.


Select OK from pop-up screen and press ENTER key to return to the main system screen.

5. Delete User Data Function

This function deletes a username, a user program and all programs in a Use folder. Users can now delete all programs in the User folder with one stroke, instead of having to go to the folder and delete each program separately. Note that the User folder (user name) is not deleted. A User name can be deleted using the "Delete User Name" function.

Procedure


1. Go to the User folder to be deleted.
 From the main screen, go to MENU, select User and press Enter.
 Select the user name to be deleted, and then go to MENU, choose OK and press Enter to move to the desired folder. The User folder in this example is named <Common>.
2. Use  to select Tools from MENU and press Enter to display the Select Tools screen.

3. In Select Tools, select Delete User Data, choose OK from MENU and press Enter to move to the Delete User Data screen.

4. Select the desired function from the Delete User Data screen.

Delete User Data provides the following functions.

- 1) Delete User Program
 Deletes an individual program in the user folder.
- 2) Delete User Name
 Deletes a user name.
 Note that the "Programs exist!" message appears when there are programs remaining in the user folder. All programs in the user folder must be deleted before deleting a user name.


3) Delete All Program in User

Deletes all programs in the user folder without deleting the user folder (user name) itself.

Use the "Delete User Name" function to delete a user name.

5-1) Using "Delete User Program"

With folder programs displayed, use to select desired program, then go to MENU, choose Delete and press Enter.

When the following confirmation message appears, select OK from MENU and press Enter.

5-2) Using "Delete User Name"

When the confirmation message appears, select OK from MENU and press Enter.


<NOTE>

The "Programs are exist!" message appears when programs are still present in the username file and the delete operation is terminated.

Making sure to delete the programs first.

5-3) Using "Delete All Program in User"

Selecting OK from MENU and pressing Enter displays the flashing "Seek and Erasing" message. Deletion is complete in about 15 seconds and the main screen is displayed.

Chapter 6: Troubleshooting

Any errors that occur during system operation are logged in the system History.

The “See History message” is displayed under Result in the system-closing screen.

(The “See History message” is also displayed when the system is paused.)

Error messages are displayed on the main system screen when errors are detected during standby mode.

1. Events and Error Code Table

Event #	Message	Description
0	Started at	Indicates program start time.
1	Finished at	Indicates program finish time.
2	Stopped at	Indicates time program stopped by user.
3	Paused at	Indicates time of system pause.
4	Resumed at	Indicates time system operation resumed.
5	Power failed at	Indicates time of system power failure.
6	Recovered at	Indicates time of recovery from system power failure.
Errors #		
7	Heat rate	Indicates abnormal heating rate.
8	Cool rate	Indicates abnormal cooling rate.
9	Over shoot	Indicates overshooting error.
10	Under shoot	Indicates undershooting error.
11	Center Sensor Open	Indicates disconnected center sensor.
12	Center Sensor Short	Indicates short in center sensor.
13	Left Sensor Open	Indicates disconnected left sensor.
14	Left Sensor Short	Indicates short in left sensor.
15	Right Sensor Open	Indicates disconnected right sensor.
16	Right Sensor Short	Indicates short in right sensor.
17	Lid Sensor Open	Indicates disconnected Lid sensor.
18	Lid Sensor Short	Indicates short in Lid sensor.
19	HS Sensor Open	Indicates disconnected HS sensor.
20	HS Sensor Short	Indicates short in HS sensor.
21	Fan Motor Open	Indicates disconnected fan motor.
22	Fan Motor Short	Indicates short in fan motor.
23	Lid control	Indicates lid temperature control error.
24	Block control	Indicates block temperature control error.

* * Note that only one error is displayed at a time, even if several occur simultaneously.

2. Troubleshooting

Symptom	Cause	Action
System fails to operate and display is off even when power is applied.	Power is not coming to the system. * Power supply cable is not properly plugged into the socket. * Power is not coming to the socket.	* Check power cable connections on system and power outlet sides. * Check installation facility power supply.
Power switch stays off even when turned on.	Fault in internal system circuitry. (Note that power switch includes an over current breaker.)	Immediately disconnect the power cable and have the system serviced.
Heating and cooling rates are abnormally low.	* Operating environment temperature may be unsuitable. * Electronic cooling element may be damaged.	Operate the system in a temperature environment of 2 to 23°C and run self-diagnosis tests. If outside rated performance, have the system serviced.
Cycle time is too long.	* Operating environment temperature may be unsuitable. * Electronic cooling element may be damaged. * Faulty temperature sensor.	Operate the system in a temperature environment of 2 to 23°C and run self-diagnosis tests. If outside rated performance, have the system serviced.
The "See operation manual" message is displayed as the self-diagnosis test result.	This is due to running the test before lid-heating time is up. Therefore, the time required to complete the test exceeded the allowable time.	Allow system to warm up for at least 10 minutes, then rerun the test.
The normal key entry sound is absent.	* Beeper setting may currently be "no". * Faulty keys or keypad wiring.	Go to Default Setting in Tool Menu and check the key entry sound setting.
The display goes off.	* LCD sleep setting may be set to "30 min". * Faulty backlight. * Faulty LCD.	Go to Default Setting in Tool Menu and check LCD sleep setting.
Display is too dark (too light).	Backlight brightness not adjusted properly.	Go to rear of display and adjust Display Brightness Knob.
Lid will not close.	Foreign matter, between lid heater and sample block. * Faulty lid lock mechanism.	Remove the foreign object or matter.
Error number appears.	Refer to list of error codes.	Confirm the nature of the error and request system servicing as required.

Power failure

Operation will be resumed from the stopped state and a program will be ended.

When a heated lid becomes by power failure below in a predetermined temperature, resume operation after warming a heated lid to an appropriate temperature

When you operate in your absence at night and an experiment is completed, Please check the following method whether abnormalities have occurred.

1. Operation End Screen Display

Finished ABC Test		MENU
User: <Common>		OK
Started at	09:47:12 05/01/20	
Finished at	09:54:32 05/01/20	
Length of Run	00:07:20	
Over 0.0 °C	Under 0.2 °C	
Result	Power failed	

Power Failed: It is shown during operation that there was a power failure.

At the time of a normal end: It is stopped when Normal and operation are stopped on the way,

2. History Display (/MENU/Tools => History: See page of 96)

History 3/6		Three events of total six events.
User: <Common>		
Program: ABC test		
Mode:1	Thin	50 µl
Started at	09:47:12	05/01/20
Power Failed	09:50:34	05/01/20
Recovered at	09:50:37	05/01/20

The first time of stop.

The first time of restart.

History 6/6		Last events of total six events.
Power Failed	10:25:10	05/01/20
Recovered at	10:28:32	05/01/20
Stopped at	11:47:12	05/01/20

The second time of power stopped.

The second time of power recovered.

When the electric current is cut off two or more times, and others have an error, it indicates by a history at the order of time. (A maximum of 40 events)

【Attention】 If a program is newly performed, History at the time of the past operation will be eliminated.

【Note】 Please refer to the page 107 of an Operation manual about 「Events and Error Code Table」.

Chapter 7: Specifications

Product name	TaKaRa PCR Thermal Cycler Dice™
Model No.	TP600 Gradient model TP650 Standard model
Operating environment	Ambient temperature: 15 °C to 30 °C Humidity: 20 to 80%RH
Temperature accuracy (*)	+/- 0.5°C (30.0 °C to 99.9 °C) +/- 2.0°C (4.0 °C to 15.0 °C)
Temperature uniformity (*)	+/- 0.5° (30.0 °C to 99.9 °C)
Temperature setting range	4.0 °C to 99.9 °C
Time setting range	1 second to 99minutes 59seconds, (Input "9999")
Display/entry panel unit angle	4 positions from perpendicular to 70 °C
Tubes and plates	0.2 mL X 96 tubes 96-well micro-plate: (12 X 8 rows)
PCR reaction volume setting	5 micro liter to 150 micro liter
Program capacity	User-names: max. 22 (excluding "Common" and "TaKaRa") Number of programs: max. 200
Lid temperature	Approx. 110 °C for the 20.0 °C to 99.9 °C range Approx. 40 °C for the 4.0 °C to 19.9 °C range
Added functions	Gradient function (only TP600) History function Diagnosis function PC communications (RS-232) (PC transmission software not provided)
Gradient temperature setting	40.0 °C to 75.0 °C (temp. difference: 6.0 °C to 20.0 °C)
Safety devices	Power switch with 10A over-current breaker Lid Heater temperature fuse Thermal guard for reaction block heater plate
Power source	100 to 240V AC, 50/60 Hz, 490VA
Dimensions	260(W) x 345(D) x 260(H) mm
Weight	11.5 kg
Battery**	CR2032
Conformity standard LR108182  	CAN/CSA-C22.2 No.10101.1-92, No.1010.1B-97 UL Std No.61010A-1 EN61000-6-3-1:2001, EN61000-6-1-1:2001, EN61010-1:2001 CE: The unit with CE mark complies with EC directives 89/336/EEC, 93/68/EEC and 73/23/EEC.
RoHS	TP600 model S/N D-6051 and over, TP650 model S/N E-1301 and over DIRECTIVE 2002/95/EC the restriction of the use of certain hazardous substances in electrical and electronic equipment.

(*) Performance specifications apply when system is operated in the 20°C to 23°C temperature range.


(**) Please request exchange of a lithium battery from a service engineer.

Appendix 1

Confirm that the following items are included in the system package when delivered

Accessories list and packing

Standard Accessories		Quantity
1.	Power cable	1
2.	Operation Manual (booklet)	1
3.	Warranty sheet (Included Operation manual)	1
4.	0.2 ml Hi-Tube Dome Cap (NJ200)	100
5.	0.2 ml Hi-8-Tube (NJ300)	12
6.	0.2 ml Hi-8-Dome Cap (NJ301)	12

TaKaRa Supplied PCR tubes (Option)

Cat.#	Product name	Size
NJ200	0.2 ml Hi-Tube Dome Cap	1,000 tubes
NJ201	0.2 ml Hi-Tube Dome Cap Recovery	1,000 tubes
NJ202	0.2 ml Hi-Tube Flat Cap	1,000 tubes
NJ203	0.2 ml Hi-Tube Flat Cap Recovery	1,000 tubes
NJ300	0.2 ml Hi-8-Tube	125 strips
NJ301	0.2 ml Hi-8-Dome Cap	125 strips
NJ302	0.2 ml Hi-8-Flat Cap	125 strips
NJ111	TaKaRa 96 well PCR Hi-Plate	50plates (12x8 well/sheet)
NJ811	TaKaRa PCR Hi-Caps (8chains)	8x120 tubes

List of Supplies (PCR Enzymes)

Cat.#		Product Name	Size
RR001A		<i>TaKaRa Ex Taq™</i>	250 U
RR001B			1,000 U (250 U × 4)
RR001C	(B × 3)		3,000 U
RR01AM		<i>TaKaRa Ex Taq™</i> (Mg ²⁺ free Buffer)	250 U
RR01BM	(AM × 4)		1,000 U
RR01CM	(AM × 12)		3,000 U
RR006A		<i>TaKaRa Ex Taq™</i> Hot Start Version	250 U
RR006B	(A × 4)		1,000 U
RR002A		<i>TaKaRa LA Taq™</i> (Mg ²⁺ free Buffer)	125 U
RR002M		<i>TaKaRa LA Taq™</i> (Mg ²⁺ plus Buffer)	250 U
RR02AG		<i>TaKaRa LA Taq™</i> with GC Buffer	125 U
RR042A		<i>TaKaRa LA Taq™</i> Hot Start Version	125 U
RR042B	(A × 4)		500 U
R001A		<i>TaKaRa Taq™</i>	250 U
R001B			1,000 U (250 U × 4)
R001C	(B × 3)		3,000 U
R001AM		<i>TaKaRa Taq™</i> (with Mg ²⁺ free Buffer)	250 U
R001BM	(AM × 4)		1,000 U
R001CM	(AM × 12)		3,000 U
R007A		<i>TaKaRa Taq™</i> Hot Start Version	250 U
R007B	(A × 4)		1,000 U
RR070A		SpeedSTAR™ HS DNA Polymerase	250 U
RR070B	(A × 4)		1,000 U
RR003A		<i>Premix Taq™</i> (<i>Ex Taq™</i> Version) (500 μl × 6)	120 PCR reaction (50 μl PCR)
R004A		<i>Premix Taq™</i> (<i>TaKaRa Taq™</i> Version) (500 μl × 6)	120 PCR reaction (50 μl PCR)
RR005A		PerfectShot™ Ex Taq (Loading dye mix)	25 μl × 48 tubes
9002A		Taq Antibody	250 U
9002B	(A × 4)		1,000 U
R010A		PrimeSTAR® HS DNA Polymerase	250 U
R010B	(A × 4)		1,000 U
R040A		PrimeSTAR® HS (Premix)	500 μl × 5

For more detailed information, refer to the Takara Bio general catalog or check for online catalog at our website (<http://www.takara-bio.com>) or (<http://www.takara-bio.co.jp>)

Appendix 2

Application 1

Gradient PCR Using the TaKaRa PCR Thermal Cycler Dice™ (TP600) System

8-kbp amplification performed in 2 or 3 steps using DNA template

(Reaction Mix)

TaKaRa Ex Taq (5 units/ μ l)	0.25 μ l
10 \times Ex Taq buffer	5 μ l
dNTP Mixture (2.5 mM each)	4 μ l
Primer 1 (20 μ M)	0.5 μ l
Primer 2 (20 μ M)	0.5 μ l
DNA (0.5 ng/ μ l)	1 μ l
Sterilized water	38.75 μ l
Total volume	50 μl

(PCR Conditions)

- 2 step

98°C	10 sec] 30 cycles
52 to 72°C	5 min	

(Gradient temperature of 20°C)

Gradient



- 3 step

98°C	10 sec] 30 cycles
45 to 65°C	30 sec	
72°C	5 min	

(Gradient temperature: 20°C)

(Results)



· Results for 2-step PCR

M : *Hind* III marker

Higher annealing/extension temperatures result in reduced low-molecular nonspecific amplification and increased amplification of the target amplicon.

· Results for 3-step PCR

M : *Hind* III marker

Higher annealing/extension temperatures give reduced low-molecular nonspecific amplification and increased amplification of the target amplicon.

Application 2

One-Step RT-PCR

The performance of TaKaRa PCR Thermal Cycler Dice™ (Cat.# TP600) was evaluated through One Step RNA PCR by using TaKaRa One Step RNA PCR Kit (Cat.# RR024A) in comparison with TaKaRa PCR Thermal Cycler SP. RT-PCR was performed in two conditions :

- In the presence of AMV RTase XL in the reaction mixture.
- Without AMV RTase XL in the reaction mixture.

(Reaction Mix)

10 × One Step RNA PCR Buffer	5 μl
25mM MgCl ₂	10 μl
dNTP Mixture (10 mM each)	5 μl
RNase inhibitor(40 U/μl)	1 μl
AMV RTase XL (5 U/μl) or DEPC treated water	1 μl
AMV Optimized Taq (5 U/μl)	1 μl
Control F-1primer 1 (20 μM)	1 μl
Control R-1primer 2 (20 μM)	1 μl
Positive Control RNA	1 μl
RNase Free H ₂ O	24 μl
<hr/>	
Total volume	50 μl

(PCR Conditions)

50°C 15 min (RT reaction)
 94°C 2 min (inactivation of RTase)
 94°C 30 sec
 60°C 30 sec
 72°C 1.5 min

30 cycles
 PCR

M 1 2 M 1 2

(Result)

Lane.	Sample
1.	RTase(-)
2.	RTase(+)
M.	100 bp DNA Ladder

(reaction)

Dice Lane 1 : E6
 Dice Lane 2 : E7

Dice SP
 (TP600) (TP400)

After the reaction was completed, a portion of reaction was applied for agarose gel electrophoresis. In both cases using TaKaRa PCR Thermal Cycler Dice™ (Cat.# TP600) and TaKaRa PCR Thermal Cycler SP, control reaction product 462 bp was detected only under the existence of RTase (+).

Application 3

Comparison of 8-kbp DNA Region Amplification for TaKaRa PCR Thermal Cycler Dice™ vs SP

- Evaluation target: Compare differences in 8-kbp DNA region amplification between TaKaRa PCR Thermal Cycler Dice™ (TP600) and TaKaRa PCR Thermal Cycler SP.

- Reaction tube: TaKaRa Micro PCR Tube (Cat.# 9047)
- Reaction mode: 2 (Selected at PCR Thermal Cycler Dice™ reaction time)

- Reaction volume: 50 µl

< Reaction Mix > < PCR Conditions >

1 x PCR Buffer 94°C 30 seconds, 65°C 10 min at
30 cycles

200 µM dNTP mixture

10 µM each primer

1 ng DNA

1.25 U TaKaRa Taq (Cat.# R001A)


up to 50 µl

All 96 wells were used for the Dice (TP600) reaction. Random samples were placed 3 per row for the SP reaction.

After the PCR reaction, 1% Agarose gel electrophoresis was performed using 5 µl of reaction mixture.

<Result>



Dice (TP600) samples amplified nicely in all wells with hardly any inconsistencies, indicating reactivity equal to or better than the SP (TP400).

Application 4

Comparison of 35-kbp DNA Region Amplification for TaKaRa PCR Thermal Cycler Dice™ vs SP

- Evaluation target : **Compare differences in 35-kbp DNA region amplification between** TaKaRa PCR Thermal Cycler Dice™ (TP600) and TaKaRa PCR Thermal Cycler SP. TaKaRa LA Taq (Cat.# RR002A) was used for the reaction.
- Reaction tube: TaKaRa Micro PCR Tube (Cat.# 9047)
- Reaction mode 2: Selected at TaKaRa PCR Thermal Cycler Dice™ MODE 2.
- Reaction volume: 50 µl (with addition of 25 pg template/50 µl)
- PCR conditions: 94°C 1 min, (98°C 10 s, 68°C 15 min) 30 cycles, 72°C 10 min
- After reaction: 1% Agarose gel electrophoresis was performed using 5 µl of reaction mixture.

Appendix 3 (Please copy and use this page)

Servicing Request and Safety Check

Attention: Distributor of Takara Bio Inc.

Product Name: TaKaRa PCR Thermal Cycler Dice™ Model: TP600 / TP650
Serial No _____ Date of request: _____

Requested by
Name _____ Contacts: (Telephone, Email) _____
Company/Institute name _____
Address _____

Details of request: _____
(Current status) _____

Customers are requested to provide the following information on the type of the function and any safety considerations. Safety check covers all personnel involved in using the system, as well as safety of the location in which the system was operated. If contaminated, the customer is requested to decontaminate the system prior to servicing.

1. Material in use ()
·Pathogenicity: { none, possibly, yes } ·Toxicity: { none, possibly, yes }
·Radioactive substances: { not used, used: nuclide() }

2. System contamination status
no contamination, decontaminated, possibly contaminated
This is to confirm that the system is free of contamination.
This is to confirm that there is risk of contamination and that specific safety measures and decontamination procedures were implemented.

Date (year/month/day): / /

Company/Institution _____ Witness by (signature) _____

This Safety Confirmation Form shall be presented as a statement of the safety of the equipment to be serviced by customers requesting servicing and maintenance work. We request your cooperation as it is designed to ensure the safety of all those who come into contact with the system.
* Note that we may decline to perform servicing and maintenance work if we do not receive the Safety Confirmation Form.

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Thank you for purchasing the TaKaRa PCR Thermal Cycler Dice™ Model TP600/TP650.

In case of malfunctions that occur despite correct operation on the part of the customer, Takara Bio Inc. shall repair the system in accordance with the conditions listed in this Warranty.

Customers may have their device repaired by presenting this Warranty to a Takara Bio Inc. sales agent or at a Takara Bio Inc. service location within the period listed below.

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1 The following repairs shall be subject to a fee, even during the warranty period.

- (a) Malfunctions resulting from incorrect operating procedures, customer repairs or renovations.
- (b) Repairs resulting from moving the system to a different location on site, moving it to a new site, dropping or transporting the system.
- (c) Malfunction or damage resulting from an act of God (including earthquakes or lightning), pollution or abnormally high voltage input.
- (d) In cases of failure to produce this warranty card.
- (e) In cases where the customer has failed to fill in the date of purchase or has altered the information on the warranty card.
- (f) In cases where the cause of the malfunction does not originate in the product itself and improvements are required.
- (g) In cases where causes the malfunction originates in system supplies that need to be replaced.

2 Takara Bio Inc. shall not be responsible for malfunctions for accidents that occur when the product is connected and used with equipment not recommended by the company.

3 Be sure to consult with Takara Bio Inc. or the place of purchase prior to moving the product.

4 Store this Warranty card in a safe location, as it cannot be reissued.

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* This Warranty covers repairs that shall be made at no cost to the customer for the period listed above.

For information on the cost of repairs after the Warranty period, please to consult with a Takara Bio Inc. branch or local sales agent.

* If there is no the TaKaRa's seal in this warranty card, this warranty is invalid.

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