Operating Manual

BlueInOne

CELL

BlueSens

Content

1. ABOUT THIS DOCUMENT	4
1.1 Function.....	4
1.2 Target group.....	4
1.3 Symbols used.....	4
2. FOR YOUR SAFETY	5
2.1 General information.....	5
2.2 Authorized personnel	5
2.3 Proper use.....	5
2.4 Misuse warning	5
2.5 General safety instructions.....	5
2.6 CE conformity.....	5
3. PRODUCT DESCRIPTION	6
3.1 BlueInOne _{CELL}	6
3.2 Measuring principle	6
3.2.1 Pressure measurement principle and moisture measuring principle	7
3.2.2 CO ₂ measuring principle	8
3.2.3 O ₂ measuring principle	8
4. INSTALLATION	10
4.1 General instructions	10
4.1.1 Included in scope of supply.....	11
4.1.2 Unboxing and assembling.....	12
4.2 Mechanical connection	15
4.2.1 Installing the sensor on pipes	16
4.3 Electrical connection	17
4.3.1 General information.....	17
4.3.2 Electrical connection of the sensor	18
4.3.3 Connection via RS232/USB	18
4.3.4 Connection RS485/USB: General Information about the RS485 connection.....	20
4.3.5 Connection RS485/USB: connect one BlueInOne.....	21
4.3.6 Connect two or more <i>BlueInOne</i> via RS485/USB	24
4.3.7 Pin assignment for an individual serial connection plug A: RS232/RS485.....	25
4.3.8 Pin assignment for an individual analogue connection Socket B	26
4.3.9 Set-up after manual analogue and serial pin assignment.....	27
4.3.10 LED status displays of the sensor.....	28
4.3.11 Software BlueVis.....	28
5 MAINTENANCE.....	29
5.1 1-point calibration	29
5.2 Recalibration.....	29
5.3 Filter change – coarse filter	31
5.3.1 Removing the sensor from the flow adapter	31
5.3.2 Changing the filter	31
5.4 Wall-mounting.....	33
6 APPENDIX.....	34
Calibration table.....	34
Technical data	35

1. About this document

1.1 Function

This operating manual provides you with all of the necessary information for quick start-up and safe operation of the **BlueInOne_{CELL}**. Therefore, please read this operation manual before starting operation.

1.2 Target group

This operating manual is intended for use by trained specialist personnel. The content of this manual must be made available to the corresponding personnel and followed by them.

1.3 Symbols used

Hazard!

This symbol indicates a possibly hazardous situation. Failure to observe this safety instruction can result in personal injury.

Caution!

This symbol indicates the possibility of damage to property.

Note!

This symbol indicates helpful, additional information.

List

The bullet point indicates a list in which the order of items has no particular significance.

1 Sequence of actions

Numbers indicate steps to be performed in a certain sequence.

2. For your safety

2.1 General information

The **BlueInOne**_{CELL} was tested and ready for operation when it left our plant.

Please read this operating manual carefully before installing the device and starting operation. The operating manual contains safety instructions that must be observed to ensure safe operation.

The device may never be operated in conditions that do not correspond to the specifications on the type plate.

Maintenance and servicing may only be performed by specially trained personnel who are familiar with the hazards inherent to the work as well as the guarantee terms.

2.2 Authorized personnel

All actions described in this operating manual may only be conducted by specialist personnel who have been authorized by the plant operator. Work on the device other than that described in this manual may only be performed by personnel of the BlueSens gas sensor GmbH company for safety reasons and to ensure compliance with the terms of the guarantee.

2.3 Proper use

The **BlueInOne**_{CELL} is a gas sensor for measuring carbon dioxide and oxygen concentrations in the specified concentration area under the conditions described in the technical data. It is used to monitor metabolism in biological processes such as fermentation.

2.4 Misuse warning

The **BlueInOne**_{CELL} may not be used as a safety component for monitoring gases in systems or be used as a gas warning device.

2.5 General safety instructions

If the device is not used correctly, or used for its intended purpose, application-specific hazards may arise.

Danger!

If the device is incorrectly installed or set, there is a danger of explosions and poisoning.

After installation, check all connections for leaks.

Note! Please note the safety instructions in all the following chapters, especially in the chapters 3 and 4.

2.6 CE conformity

See the last page of this manual for the CE certificate.

3. Product description

3.1 BlueInOne_{CELL}

The BlueInOne_{CELL} is a combined CO₂ and O₂ sensor with automatic humidity and pressure compensation. The gas to be examined is fed into the three measuring chambers via the integrated flow adapter and analyzed during this process for humidity, pressure and the CO₂ and O₂ content. Flow adapters are available for all pipe diameters from 4 mm to 1 ¼". Please remember that the dimensions and the weight of the analyzer can vary depending on the flow adapters in use.


Fig. 1: Schematic diagram of the overall construction

3.2 Measuring principle

The analyzer contains 4 different measuring components:

1. Pressure measuring
2. Humidity measuring
3. CO₂ measuring
4. O₂ measuring

3.2.1 Pressure measurement principle and moisture measuring principle

A piezoresistive, silicon pressure sensor is used as the pressure sensor.

The moisture measurement is made using a capacitive, polymer sensor.

Hazard!

- **Be careful if you want to use the BlueInOne in extreme environment!**
- **Do not use the device in temperatures outside the specifications.**
- **Do not expose the BlueInOne to volatile organic compounds.**
- **The temperature and humidity sensor may be damaged if you expose them to high concentrated acids or bases or to high concentrations of etching substances such as NH_3 , H_2O_2 etc. Corrosive substances may damage the sensor.**

3.2.2 CO₂ measuring principle

The CO₂ sensor consists of an IR light source, a detector and the reflection measuring cell (fig. 2).

The infrared light source is reflected by a gas-filled measuring adapter and the light weakened by the analyte gas is measured by the detector. The light-permeable sapphire disc prevents the sample atmosphere from escaping and contaminating the optical components. The entire sensor is heated so that no moisture can condensate. The heating-up time takes approx. 1 hour. This is the case after initial start-up and also after each time the device is disconnected from the power supply. The sensor does not output any measuring values during the heating-up time.


Fig. 2: Measuring design of the CO₂ sensor

3.2.3 O₂ measuring principle

The oxygen sensor in the BlueInOne_{CELL} is based on a lead-oxygen battery which incorporates a lead anode, an oxygen cathode made of gold, and a weak acid electrolyte. Oxygen molecules enter the electrochemical cell through a non-porous fluorocarbon membrane, diffuse in the acid electrolyte, and are reduced at the gold electrolyte. The current which flows between the electrodes is proportional to the oxygen concentration in the gas mixture being measured. The terminal voltages across the thermistor (for temperature compensation) and resistor are read as a signal, with the change in output voltages representing the change in oxygen concentration.

**Note!**

High concentration of ammonia (NH_3) or ozone (O_3) could minimize the lifetime of the O_2 sensor element.

4. Installation

4.1 General instructions

The BlueInOne_{CELL} was protected on its way to the installation location by packaging. This secures it against the usual transport strains. However, before installation, check if the device has suffered damaged due to improper transport or incorrect storage. If the device is damaged in any way, operation without hazards is not possible and the device may not be installed and taken into operation.

Check whether the enclosed materials, such as gaskets, are suitable for your process conditions (pressure, temperature, etc.).

The installation should only be performed under supervision by a specialist and in compliance with all applicable work safety rules.

4.1.1 Included in scope of supply


Fig. 3: Scope of supply

	Article No.	Description
A	Z-NT-00009	24 V, 1A power supply international
B	Z-KA-00019	Adapter for M12 (8pol) to USB (RS232)
C	Z-KA-00025	Connection cable for BlueInOne, incl. power connector, 0,2m
D		Flow adapter and connection pipe/hose as requested
E	Z-KA-00026	Connector M12 (male)
F	Z-KA-00027	Connector M12 (female)
G		Mounting plate
	Z-XX-00075	Screw-set BlueInOne*
	Z-XX-00053	Filter set BlueInOne*
	Z-XX-00074	Screw driver TORX T20*
	CC-BlueInOne ferm/cell*	Calibration certificate BlueInOne
	BD-BlueInOne ferm/cell*	Manual BlueInOne Ferm/Cell
	BACCheck*	Diagnostics software via RS232, (via download)

*not shown on picture

4.1.2 Unboxing and assembling

The flow adapter of the BlueInOne is not pre-assembled and will be delivered in a small box together with the BlueInOne (fig 3a). Before the installation, it must be connected to the BlueInOne.

For the assembling you will need (fig. 3b):

1. The BlueInOne (A)
2. The TORX T20 screwdriver (art. no.: Z-XX-00074) (B)
3. The mounting plate (C)
4. The flow adapter (D)
5. Bag with 2 flat head screws M4x6 TORX for the wall mounting and 4 screws M4x16 TORX for the flow adapter (E) (art.no. Z-XX-00075).
6. Bag with filter set BlueInOne (art. no.: Z-XX-00053) (F).

To assemble the BlueInOne proceed in the following steps:

1. First place the filters and gaskets on the flow adapter. All needed filters and gaskets are in the filter set for the BlueInOne (art. no.: Z-XX-00053). For each recess there are two different filters. There is a flat filter (fig. 3c: A) and a thicker filter (fig. 3c: B). The flat filter is hydrophobic and the thicker filter is to hold back pollutions.

Fig. 3a: Unboxing

Fig. 3b: Preparing the assembling


Fig. 3c: Filters

2. For each recess there is a gasket (fig. 3d).

Fig. 3d: Filters and gaskets

3. Insert the flat filter into the recesses (fig. 3e).


Fig. 3e: Insert flat filter

4. Place the gaskets in the recesses (fig. 3f)


Fig. 3f: Place gaskets

5. Insert the thicker filters into the recesses (fig. 3g). If the filters don't fit in the recesses, it might be necessary to push them gently into place (fig. 3h).

Fig. 3g: Insert second filter

Fig. 3h: Fit in filter

6. Now place the 4 TORX screws M4x16 (fig. 4b:E) in the bores of the BlueInOne without screwing. This will make it easier to place the BlueInOne on the right spot on the flow adapter (fig. 4i).
7. Now place the 4 TORX screws M4x16 (fig. 4b: E) in the bores of the BlueInOne without screwing. This will make it easier to place the BlueInOne on the right spot on the flow adapter (fig. 4i).
8. Place the BlueInOne on the flow adapter. Tighten the 4 screws with the screw driver TORX.
9. Put the wall mounting plate on the flow adapter and use the two flat head screws M4x6 TORX (fig. 4b: E) to fix the wall mounting (fig. 4k).


Fig. 4i: Put screws in

Fig. 4j: Tighten screws

Fig. 4k: Fit wall mounting

4.2 Mechanical connection

Caution!

The integrated filter is not intended to repel fluids (see chapter 5.3). Never install the sensor so that fluid can run into it.

If water penetrates the flow adapter, it must be rinsed with a sufficient amount of dry air.

After installation, check that the pipe connection is gas-tight.

Do not use solvent-containing sealant like Locktide™ or similar. The analyzer could be damaged by the solvent. Please use Teflon or hemp tape for sealant if needed.

4.2.1 Installing the sensor on pipes

The sensor has a flow adapter (fig. 5: A), through which the gas to be measured flows. It can be adjusted in sizes from 1/8" to 1 1/4".

Various connection pieces for pipes and hoses are available for the sizes 1/8" to 1/4". This allows the sensor to be connected directly to the exhaust line of fermenters.

There are filters for retaining particles placed between flow adapter and the actual measuring cell. These filters are not water-permeable.

The sensor should always be attached so that **NO** fluid can run into the measuring cells. Never install the BlueInOne at the lowest point of the line. Make sure that the gas can flow in a vertical direction (from up to down) through the flow adapter (fig. 5).


Fig. 4: Gas flow BlueInOne

Fig. 5: The various connection options for hoses. The drill hole in the flow adapter (A) has a 1/4" thread. The flow adapter can be exchanged allowing for threads of up to 1 1/4" to be used. Other constructions can be manufactured on request.

4.3 Electrical connection

4.3.1 General information

Caution!

Read the installation instructions carefully to avoid damage to the device.

Proceed step-by-step.

Only use the original plugs, cables and power adapters.

Never connect or disconnect plugs when the device is connected to the power supply.

The device does not have an on/off switch; it starts operation as soon as it is connected to the power supply.

Improper operation can result in damage to the device.

4.3.2 Electrical connection of the sensor

For the electrical connection of the sensor please note the following chapter 4.3. et seqq..

In chapter 4.3.3 to 4.3.6 the different installations with pre-assembled connectors and adapters for the serial to USB connection will be explained.

In the chapters 4.3.7 to 4.3.9 the installation with individual pin assignment for analog signal as well as digital signals will be specified for experienced operators.

In general the connection plug A is for the serial data transfer and connection socket B is for analogue signals (fig. 6).

4.3.3 Connection via RS232/USB

The following Items will be needed for the electrical connection RS232 to USB (fig. 7):

A: Connection cable for BlueInOne, incl. power connector, 0,2m

B: A BlueInOne

C: 24 V, 1A power supply international (ser. no.: Z-NT-00009). Use original equipment only.

D: Adapter for M12 (8pol) to USB (RS232) (ser. no.: Z-KA-00019)

To establish the connection:

1. Please check the USB connector for the "RS232" sticker to assure that you are using the correct USB adaptor (fig. 8).
2. Connect the adapter for M12 (8pol) to USB (RS232) (ser. no.: Z-KA-00019) with the M12 socket of the connection cable for BlueInOne (fig. 9). Tighten the screwed ring.


Fig. 6: Connection plug A and connection socket B of the BlueInOne


Fig. 7: Connection RS232

Fig. 8: RS232-sticker

3. The spare M12 socket with the attached power cable (fig. 9: A) must be connected with the M12 connection port (port A) of the BlueInOne (fig. 10). Put the power adaptor into the mains socket (fig. 10).
4. Connect the power connector to the power plug, tighten the screwed ring and put the USB connector to a spare USB port of the computer (fig. 11). If you are using Windows XP® or Vista® it might be necessary that you first install all service packs from Microsoft and the adequate USB drivers from the BlueInOne data carrier. Contact your responsible administrator if the drivers can't be found or installed. If you are using Windows7™ all needed drivers should be included already.


Fig. 9: connect cables


Fig.10: Connect the BlueInOne


Fig. 11

Note that the BlueInOne has no separate power switch.

The BlueInOne warms up for approx. 60 minutes. During this time, the right-hand pushbutton will light up in yellow. The left-hand pushbutton then lights up in green while the right-hand one remains yellow. This is how the analyzer indicates that a 1-point calibration needs to be carried out. Continue with the calibration (see chap. 5.1).


Fig. 12: Warm- up time

4.3.4 Connection RS485/USB: General Information about the RS485 connection

RS485 is a data bus that can be used for the transfer of measuring data. RS485 together with an adapter to USB is used as a standard for the application in combination with the BlueVis software. To start the measurement with the *BlueVis* software, at least one *BlueInOne* gas analyzer needs to be connected to the personal computer or a process control system. Up to 12 analyzers can be connected by a data bus via the RS485 data bus to a single interface. That means that only one COM-port is needed for all data. Each analyzer needs an individual MODBUS-ID for this data bus so that the data from each sensor can be addressed in the right way. The Modbus-ID is a number from 1 to 247. The MODBUS-ID is pre-installed at the factory. It is possible to change this ID with BlueVis (see the corresponding manual). **The RS485 data bus must be terminated with a terminating resistor (article no.: Z-KA-00012) at the end of the data bus (fig. 13).**


Fig. 13: Installation scheme RS485

4.3.5 Connection RS485/USB: connect one BlueInOne

You may connect between 1 to 12 BlueInOne to run in one single RS485 data bus. See this chapter for the connection of a single BlueInOne and chap. 4.3.6 for the connection of more than one BlueInOne.

To run a *BlueInOne* with a RS485 to USB connection, the following items (fig. 14) are required:

A: Adaptor USB-RS485-M12 (art. no.: Z-KA-00015)

B: The *BlueInOne*

C: A power supply 24V, 1A. Use original equipment only (art. no.: Z-NT-00009).

D: Connection cable BlueInOne (art. no.: Z-KA-00025).

E: Terminating resistor (art. no.: Z-KA-00012).


Fig.14: Accessory for the installation of the BlueInOne via RS485/USB

Fig.15: Terminal resistor

1. Connect the terminating resistor to the M12 female connector of the connection cable *BlueInOne* that is **not** connected to the power cable (fig. 15). Tighten the screwed ring.

2. Connect the adaptor USB-**RS485**-M12 (art. no.: Z-KA-00015) to the M12 male connector of the connection cable *BlueInOne* (fig. 16). Note, that you need the adaptor for **RS485** for this connection and **NOT** the adaptor for **RS232**.

Fig.16

The adaptor can be USB-RS485-M12 identified by the sticker on the USB-connector (fig. 17). Check that it is an **RS485** sticker.

Fig.17: USB connector RS485

3. Connect the spare female connector (with the power supply cable) of the connection cable *BlueInOne* to the male connector of the *BlueInOne* (**Port A**, fig. 18) and plug in the power supply (fig. 18).


Fig.18: Connect *BlueInOne*

4. Connect the USB connector to a USB-port (fig. 20). The USB-connection is for the data transfer only and not for the power supply of the *BlueInOne*. If you are using Windows XP® or Vista® it might be necessary that you first install all service packs from Microsoft and the adequate USB drivers from the *BlueInOne* data carrier. Contact your responsible administrator if the drivers can't be found or installed. If you are using Windows7™ all needed drivers should be already included already.


Fig. 19

5. Connect the power connection to the connection cable *BlueInOne* (fig. 20). The *BlueInOne* is now connected to the power supply and it will start to heat up and the red light will light up for some seconds. Note, that the *BlueInOne* has no separate power switch.

Fig. 20: Warm- up time

The *BlueInOne* warms up for approx. 60 minutes. During this time, the right-hand pushbutton will light up in yellow. The left-hand pushbutton then lights up in green while the righthand one remains yellow. This is how the analyzer indicates that a 1-point calibration needs to be carried out. Please note the

individual specifications of your *BlueInOne* in your data sheet and observe the manual of your *BlueInOne*. If you perform the 1-point calibration with the wrong gas, wrong measuring results will occur. Continue with the 1-point-calibration now (see chap. 5.1).

4.3.6 Connect two or more *BlueInOne* via RS485/USB

Assure the *BlueInOne* have got **different MODBUS IDs** before you start the installation and assign new IDs if needed (see manual *BlueVis*).

To connect two or more *BlueInOne* to a single RS485 data bus the following items are required (fig. 21):

A: One adaptor USB-RS485-M12 (art. no.: Z-KA-00015).

B: Two or more *BlueInOne*.

C: One connection cable *BlueInOne* (art. no.: Z-KA-00025) for each *BlueInOne*.

D: One terminating resistor (art. no.: Z-KA-00012).

E: A power supply for each *BlueInOne* (not shown on fig. 21, art. no.: Z-NT-00009).

F: An optional M12 extension cable (not shown on fig.21, for up to 5 meters length) could be used to extend the distance between the analyzers.

To establish the connection:

1. Connect all the parts as shown (fig. 22). Always use the Port A of every *BlueInOne* for this type of connection. If you have more than two *BlueInOne* to integrate, install the other analyzers accordingly. At the end of the bus a terminating resistor is always necessary and also an adaptor USB-RS485-M12 at the beginning of the bus.
2. Continue with the installation analog to the installation of a single *BlueInOne* (chap. 4.3.5).


Fig 21: Installation of two or more *BlueInOne*

Fig.22: Data bus

4.3.7 Pin assignment for an individual serial connection plug A: RS232/RS485

The following advices in the chapter are for operators that want to manually access the serial signals from the BlueInOne and do not want to use preassembled connectors. The references are for specialized personnel.

Note!

The numbering of the pins and their assignment refer to the socket when seen from behind (fig. 23).

Remove the insulation of the cables a little as possible to avoid short circuits in the plug housing.

Plug A connection assignment


Fig. 23: Assignment of the supplied socket A for plug A

4.3.8 Pin assignment for an individual analogue connection Socket B

The following advices in this chapter are for operators that want to manually access the serial signals from the BlueInOne and do not want to use preassembled connectors. The references are for specialized personnel.

Note!

The numbering of the pins and their assignment refer to the socket when seen from behind (fig. 24).

Remove the insulation of the cables a little as possible to avoid short circuits in the plug housing.

Connection assignment of socket **B** of the sensor


Fig.24: Assignment of the supplied plug **B** for socket **B**

4.3.9 Set-up after manual analogue and serial pin assignment

1. Connect the 24 V DC voltage supply to pin 1 of socket **A**.
2. Connect GND (ground) to pin 2 of socket **A**.
3. The serial query can be performed by connecting the corresponding cables in accordance with fig. 11 to socket **A**.
4. To display the **CO₂** channel, connect the analog measuring device to pin 1 ($R_L < 500 \text{ Ohm}$) and pin 2 (return conductor) of plug **B**.
5. To display the **O₂** channel, connect the analog measuring device to pin 3 ($R_L < 500 \text{ Ohm}$) and pin 4 (return conductor) of plug **B**.
6. Insert plug **B** in socket **B**.
7. Connect socket **A** to plug **A**.

After approx. 1 hour of heating-up time, the sensor still requires adjusting. During the heating-up time the sensor displays approx. 2.3 mA. To make the adjustment, expose the sensor for at least 30 minutes to the gas stated in the individual data sheet of the sensor. The flow rate must be 200ml per minute or higher. If the BlueInOne was exposed to high concentrations of CO₂ flushing for several hours is necessary. To save calibrating gas in this special case, normal fresh air could be used in the first time and the specified gas could be used only in the last 30 minutes. Note that not all BlueInOne analyzers are suitable for a 1-point-calibration with fresh air. Any changes in the gas mixture for 1-point calibration will cause errors in the measurement. See table 1 on the right side to find the standard gas concentrations for the **BlueInOne_{CELL}**. Please also note

Caution!

Connect all cables carefully. Remember that the device requires a heating-up time of 1 hour before it functions according to specifications.

During this time, the analog outputs only output 2.3 mA. The LED rings by the switches also glow and show the corresponding status of the sensor.

Standard concentrations for 1-point calibration always see data sheet to find out for which gas your sensor is calibrated.

	fresh air 0,04 Vol.% CO ₂ , 20,97 Vol.% O ₂ , rest N ₂ *	20,97 Vol.% O ₂ , rest N ₂ *
BlueInOne Cell 10	x	---
BlueInOne Cell 25	x	---
BlueInOne Cell 50	---	x

* Changes in concentration will cause errors!

Table 1: 1-point calibration for **BlueInOne_{CELL}**

the individual data sheet of your BlueInOne. Advices on the data sheet have absolute priority.

8. Afterwards, connect pin 7 to pin 8 (GND) of plug **B** for 5 seconds or press the two integrated switches down for five seconds at the same time.

9. Screw on the strain relief.

The sensor has now been adjusted.

4.3.10 LED status displays of the sensor

LED condition	Description	Possible causes
Green	Everything OK, sensor ready for operation	
Yellow	Heating-up phase, can last up to 60 minutes depending on various factors	Sensors are being heated and must first attain the correct temperature as otherwise the measuring values could be incorrect.
Green/yellow	1-point adjustment required	1st operation, , and after approx. every month of operation
Green, flashing	1-point adjustment running	Adjustment has started
Green/yellow, flashing	Calibration at plant required	Recommended operating hours exceeded
Red	Sensor not working	Sensor is currently starting up, signal too weak, sensor defective


Fig. 25: Switch with LEDs

4.3.11 Software BlueVis

Sensor of BlueInOne series could be connected via RS485-Modbus to a computer. The software BlueVis could display not only the measured gas concentration but also OUR, CER and RQ. In addition the OPC connectivity to other systems is available.

To run BlueVis a dongle and the corresponding license keys are necessary.

5 Maintenance

We recommend you send the sensors to BlueSens for annual maintenance, inspection and calibration.

5.1 1-point calibration

Once a month, the sensor needs a 1-point calibration to gain exact measurements. To make the adjustment, expose the sensor for at least 30 minutes to the gas stated in the individual data sheet of the sensor. The flow rate must be 200ml per minute or higher. If the BlueInOne was exposed to high concentrations of CO₂ flushing for several hours is necessary. To save calibrating gas in this special case, normal fresh air could be used in the first time and the specified gas could be used only in the last 30 minutes. Note that not all BlueInOne analyzers are suitable for a 1-point-calibration with fresh air. Any changes in the gas mixture for 1-point calibration will cause errors in the measurement. See table 1 on the right side to find the standard gas concentrations for the **BlueInOne_{CELL}**. Please also note the individual data sheet of your BlueInOne. Advices on the data sheet have absolute priority.

Afterwards connect pin 5 to pin 8 (GND) on the connection cable or press the two buttons down simultaneously for five seconds (fig. 8).

The 1-point calibration can also be executed with the **BlueVis**-software. See the respective manual for details.

With the serial version, the adjustment can be performed with the **FermVis** or **BACVis** software.

5.2 Recalibration

The sensor should be sent back to the manufacturer or an authorized dealer for annual recalibration.

Standard concentrations for 1-point calibration always see data sheet to find out for which gas your sensor is calibrated.

	fresh air 0,04 Vol.% CO ₂ , 20,97 Vol.% O ₂ , rest N ₂ *	20,97 Vol.% O ₂ , rest N ₂ *
BlueInOne Cell 10	x	---
BlueInOne Cell 25	x	---
BlueInOne Cell 50	---	x

* Changes in concentration will cause errors!

Table 1: 1-point calibration for **BlueInOne_{CELL}**

Fig. 26: Starting the 1-point calibration

You can get further information for our annual maintenance service Blue4Care incl. extension of the warranty up to 6 years on:

www.bluesens.com → Service/Downloads→Blue4Care

5.3 Filter change – coarse filter

5.3.1 Removing the sensor from the flow adapter

1. Unscrew the four screws visible from the top (see fig. 27) and take the sensor off the flow adapter.

5.3.2 Changing the filter

Caution!

The integrated filter is not intended to repel fluids. Never install the sensor so that fluid can run into it.

If water penetrates the flow adapter, it must be rinsed with a sufficient amount of dry air.

2. Remove the filter from the recess (fig. 28).
3. All needed filters and gaskets are in the filter set for the BlueInOne (art. no.: Z-XX-00053). For each recess there are two different filters. There is a flat filter (fig. 29: A) and a thicker filter (fig. 29: B). The flat filter is hydrophobic and the thicker filter is to hold back pollutions. For each recess there is a gasket, a flat filter and a thicker filter in the filter set (fig. 29).

Fig. 27: Screws for fastening the sensor to the flow adapter


Fig. 28: Filter removal


Fig. 29: New filters

4. Insert the flat filter into the recesses (fig. 30).
5. Check the gaskets for damage and replace them as required (fig.31).
6. Insert the thicker filters into the recesses (fig. 32). If the filters don't fit in the recesses, it might be necessary to push them gently into place (fig. 33).
7. Fasten the sensor to the flow adapter again with the 4 screws.

Fig. 30: New filters and gaskets


Fig. 31: Insert flat filter


Fig. 32: Replace gaskets

Fig. 33: Insert second filter

Fig. 34: Fit in filter

5.4 Wall-mounting

The BlueInOne_{CELL} can be mounted nearly anywhere with the aid of a mounting panel.

Fig. 34: Sensor with mounting panel

6 Appendix

Calibration table

Complete calibration can only be conducted by BlueSens. Monthly 1-point calibration can be performed as described in chapter 5.1.

Fill out the table below when this is performed.

[illegible]

Technical data

See enclosed datasheet.

Blue4Care - The maintenance service by BlueSens

BlueSens' gas sensors are high quality measuring devices. To ensure accurate operation and to extend the guarantee for one year we recommend annual maintenance in our factory.

Blue4Care is a full service package for your BlueSens gas sensors. It's not necessary to order and pay for the service directly with your sensor order. Within one year after the sale of the sensor we will inform you and offer the service to you.

If you decide for **Blue4Care**, we will book a date for the maintenance and calibration of your sensors. The typical timeframe to complete the service will be one week exclusive shipping time.

If you service your sensor annually using **Blue4Care**, the devices will remain in our extended guarantee scheme. This means that, excluding user damage, your sensor will keep its original guarantee for years with significantly reduced maintenance costs.

You may also request at any time "**maintenance on demand**". BlueSens will repair or replace broken parts and calibrate the sensor as required by the user. Maintenance on demand does not extend the guarantee

Advantages:

- Maintenance administration by BlueSens
- Payment and Order one year after purchase
- Reduce Rate compare to "Maintenance on demand"
- Fixed cost per year
- Annual extension of guarantee for one year after maintenance (up to 6 years)

Cost: Please visit the service-area on our homepage for prices:

www.bluesens.com → Service/Downloads

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We hereby declare that our product, type:

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complies with the following relevant provisions:

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Applied harmonized standards, in particular:

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