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|                                                                                                                     |       |
|---------------------------------------------------------------------------------------------------------------------|-------|
| 4-Methylmorpholinooxide                                                                                             | L 336 |
| 4-Methylmorpholinooxide x H <sub>2</sub> O                                                                          | L 471 |
| 4-Morpholine carboxaldehyde                                                                                         | L 392 |
| 4-Pregnene-11β,21-diol-3,20-dione 21-acetate                                                                        | L 591 |
| 5'-O-dimethoxytrityl-2'-deoxyadenosine-(N6-benzoyl)-3'-N,N-diisopropylamino-O-(2-cyanoethyl)phosphoramidite         | L 474 |
| 5-Amino-[2,6-dichloro-4(trifluoromethyl)phenyl]-4-[(1R,S)-4-(trifluoromethyl) sulphonyl]-1H-pyrazole-3-carbonitrile | L 612 |
| 5-Amino-1-cyclohexyl-7-(cis-3,5-dimethyl piperazino)-6,8-difluoro-1,4-dihydro-4-oxo-3-quinolinecarboxylic acid      | L 598 |
| 5-Hydroxy-1-methylpyrazole                                                                                          | L 540 |
| 6-Decyloxy-7-ethoxy-4-hydroxy-3-quinoline carboxylic acid ethyl ester                                               | L 627 |
| 6-Fluoro-1-(4-fluorophenyl)-1,4-dihydro-4-oxo-7-(1-piperazinyl)-3-quinoline carboxylic acid hydrochloride           | L 593 |
| 7-(D-α-amino-phenylacetamido)-3-methyl-3-cepheme-4-carboxylic acid                                                  | L 609 |
| 8-Hydroxy-5,7-dinitro-2-naphthalene sulfonic acid                                                                   | L 576 |
| ABS (Acrylonitrile-Butadiene-Styrene-copolymer)                                                                     | L 313 |
| Acesulfame K                                                                                                        | L 343 |
| Aciclovir                                                                                                           | L 657 |
| Acid-etching mixture(Nitric acid, Hydrofluoric acid, Nitrus acid and Hexafluorosilicic acid)                        | L 620 |
| Acrylamide copolymer                                                                                                | L 389 |
| Acrylate mixture                                                                                                    | L 195 |
| Acrylic acid ethyl ester                                                                                            | L 147 |
| Acrylic dispersion lacquer, water-soluble                                                                           | L 622 |
| Acrylonitrile-Butadiene-Styrene-Copolymere (ABS)                                                                    | L 313 |
| Additive for oils                                                                                                   | L 187 |
| Adenosine-5-monophosphate disodiumsalt x 6H <sub>2</sub> O                                                          | L 031 |
| Adenosine-5-triphosphoric acid disodiumsalt                                                                         | L 249 |
| Adhesive component (containing diethylenetriamine DETA)                                                             | L 569 |
| Adhesive, cyanoacrylate                                                                                             | L 118 |
| Aeroshell Fluid 41 (Mineral hydraulic oil with additives)                                                           | L 619 |
| Aerosol spray cans                                                                                                  | L 165 |
| Aldrin (Insecticide)                                                                                                | L 652 |
| Allyl alcohol                                                                                                       | L 454 |
| Almonds, chopped                                                                                                    | L 086 |
| Almonds, flaked using a homogenizer                                                                                 | L 295 |
| Alphametrine                                                                                                        | L 583 |
| Aluminum chloride hexahydrate                                                                                       | L 634 |
| Aluminum triflate                                                                                                   | L 635 |
| Aluminum trifluorosulfonate                                                                                         | L 635 |
| Aminolevulinic acid-HCl (Lyophylisate)                                                                              | L 506 |
| Aminophylline                                                                                                       | L 009 |
| Aminophylline                                                                                                       | L 468 |
| Ammonia                                                                                                             | L 431 |
| Ammonium citrate dibasic                                                                                            | L 018 |
| Ammonium glycyrrhizate                                                                                              | L 501 |
| Ammonium hydrogen fluoride                                                                                          | L 047 |
| Ammonium nitrate                                                                                                    | L 516 |
| Ammonium sulphate                                                                                                   | L 048 |
| Ammonium sulphate                                                                                                   | L 157 |
| Ammonium thiocyanate                                                                                                | L 482 |
| Amoxycillin x 3 H <sub>2</sub> O                                                                                    | L 352 |
| Ampicillin                                                                                                          | L 422 |

|                                                                       |       |
|-----------------------------------------------------------------------|-------|
| Amylase                                                               | L 149 |
| Anhydrotetracycline hydrochloride                                     | L 595 |
| Aniline                                                               | L 030 |
| Animal feed                                                           | L 438 |
| Animal feed dried                                                     | L 445 |
| Animal vaccine CHPPI                                                  | L 421 |
| Antifreeze 4%                                                         | L 186 |
| Antimony triacetate                                                   | L 289 |
| Antimony trioxide                                                     | L 056 |
| Antimony(III) chloride                                                | L 056 |
| Arginin aspartate                                                     | L 298 |
| Ascorbic acid                                                         | L 568 |
| Aspartame                                                             | L 348 |
| Aspartic acid, magnesium salt                                         | L 010 |
| Atropine sulphate x H <sub>2</sub> O                                  | L 500 |
| Baby food "ISOMIL"                                                    | L 223 |
| Barium perchlorate                                                    | L 426 |
| Benserazidhydrochloride 2-[(2,3,4-Trihydroxy)-phenyl]hydrazide-HCl    | L 416 |
| Bentonite / Sand-bentonite                                            | L 146 |
| Benzyl penicillin procaine x H <sub>2</sub> O                         | L 307 |
| Benzylamine                                                           | L 350 |
| Berberine chloride                                                    | L 465 |
| Biodiesel                                                             | L 546 |
| Bis(4-isocyanatophenyl)methan                                         | L 603 |
| Bis(dimethylthiocarbamoyl) disulphide                                 | L 615 |
| Biscuits, whole meal                                                  | L 294 |
| Bitumen emulsion                                                      | L 263 |
| Blood meal                                                            | L 212 |
| Blueberry extract                                                     | L 444 |
| Borago oil, refined                                                   | L 544 |
| Boric acid                                                            | L 494 |
| Borresperse CA (Calcium lignin sulphonate)                            | L 515 |
| Brake fluid                                                           | L 417 |
| Brake fluid                                                           | L 537 |
| Butane liquid gas                                                     | L 338 |
| Butter                                                                | L 104 |
| Butter Fat                                                            | L 458 |
| Butyl putty (Window sealant)                                          | L 399 |
| Butyrolactone                                                         | L 353 |
| Cabbage flakes, white                                                 | L 375 |
| Cacodylic acid sodium salt x-hydrate                                  | L 572 |
| Caesium fluoride                                                      | L 046 |
| Caffeine (Coffeine)                                                   | L 024 |
| Calcium acetate                                                       | L 481 |
| Calcium ammonium nitrate                                              | L 058 |
| Calcium bis-(dihydrogenphosphate)-monohydrate monobasic               | L 487 |
| Calcium carbonate                                                     | L 045 |
| Calcium dobelisate (2,5-di hydroxylbenzo sulphonic acid-calcium salt) | L 330 |
| Calcium folinate                                                      | L 366 |
| Calcium glycerophosphate                                              | L 251 |
| Calcium glyconate-1-hydrate                                           | L 259 |
| Calcium hydrogen phosphate                                            | L 041 |
| Calcium lignin sulphonate (Borresperse CA)                            | L 515 |

|                                                              |       |
|--------------------------------------------------------------|-------|
| Calcium nitrate tetrahydrate                                 | L 610 |
| Calcium oxide and Calcium peroxide                           | L 312 |
| Calcium phosphate, dehydrated dibasic                        | L 488 |
| Calcium sorbate                                              | L 286 |
| Calcium sulphate 0.5 hydrate                                 | L 042 |
| Calcium sulphate-2-hydrate                                   | L 043 |
| Calcium-bis-(2-aminoethyl)hydrogenphosphate                  | L 530 |
| Calcium-L-aspartate-HCl x 2H <sub>2</sub> O                  | L 290 |
| Calciumsulphonate overbasic                                  | L 183 |
| Candy, soft candy                                            | L 091 |
| Caprolactam                                                  | L 238 |
| Caprolactone                                                 | L 410 |
| Caramel mass                                                 | L 413 |
| Carbamic acid, ethylenebis{dithio-},manganese zinc complex   | L 614 |
| Carbazineic acid methyl ester                                | L 281 |
| Carbon tetrachloride                                         | L 001 |
| Carbonates using a KF oven                                   | L 536 |
| Carnallite, potassium magnesium chloride x 6H <sub>2</sub> O | L 059 |
| Cationic bitumen emulsion                                    | L 263 |
| CDP-choline-Na                                               | L 642 |
| Cefalexin                                                    | L 609 |
| Cefalonium                                                   | L 586 |
| Cement                                                       | L 044 |
| Cephalexin monohydrate                                       | L 473 |
| Cerium oxide CeO <sub>2</sub> x H <sub>2</sub> O             | L 484 |
| Cheese                                                       | L 095 |
| Chewing gum                                                  | L 308 |
| Chili extract                                                | L 480 |
| Chloranil (Fungicide)                                        | L 617 |
| Chloroformic acid methyl esters                              | L 017 |
| Chlor-toluol-sulphoryl chloride                              | L 394 |
| Chocolate                                                    | L 071 |
| Chocolate, milk chocolate                                    | L 079 |
| Chocolate, soft                                              | L 028 |
| Chromium dioxide (CrO <sub>2</sub> )                         | L 054 |
| Cigarette tobacco                                            | L 296 |
| Cipofloxacin-hydrochloride                                   | L 365 |
| Cis-butenedioic acid                                         | L 624 |
| Cis-ethylene carboxylic acid                                 | L 624 |
| Citicolin sodium                                             | L 642 |
| Citric acid anhydrous                                        | L 016 |
| Citric acid anhydrous                                        | L 169 |
| Citric acid diammonium salt                                  | L 018 |
| Citric acid monohydrate                                      | L 015 |
| Citric acid monohydrate                                      | L 169 |
| Cloxacillin sodium salt monohydrate                          | L 637 |
| Coal dust                                                    | L 114 |
| Coal dust (Power station coal)                               | L 467 |
| Coal dust, brown coal                                        | L 115 |
| Cobalt oxalate-2-hydrate                                     | L 357 |
| Cocoa beans                                                  | L 093 |
| Cocoa powder                                                 | L 094 |
| Cocofatamine                                                 | L 131 |

|                                                                                                                 |       |
|-----------------------------------------------------------------------------------------------------------------|-------|
| Coffee cream                                                                                                    | L 222 |
| Coffee extract                                                                                                  | L 092 |
| Coffee roasted, ground                                                                                          | L 178 |
| Coffee, green                                                                                                   | L 179 |
| Coffee, ground                                                                                                  | L 061 |
| Coffee, ground                                                                                                  | L 062 |
| Cold seal cover                                                                                                 | L 423 |
| Collodium wool                                                                                                  | L 234 |
| Colour dye                                                                                                      | L 117 |
| Copper salts                                                                                                    | L 276 |
| Copper(II)-nitrate-3-hydrate                                                                                    | L 418 |
| Corn starch                                                                                                     | L 532 |
| Corticosterone 21-acetate                                                                                       | L 591 |
| Cosmetics                                                                                                       | L 172 |
| Cream wafers                                                                                                    | L 284 |
| Crospovidone                                                                                                    | L 472 |
| Cross-linked homopolymer of 1-vinylpyrrolidone-2-on                                                             | L 472 |
| Crude oil                                                                                                       | L 108 |
| Crude oil, coulometric                                                                                          | L 148 |
| Curd cheese                                                                                                     | L 074 |
| Cyanoacrylate adhesives                                                                                         | L 118 |
| Cyanur chloride                                                                                                 | L 139 |
| Cyclophosphamide                                                                                                | L 463 |
| Cypermethrine (Alphamethrine)                                                                                   | L 583 |
| Cysteine, S-carboxymethyl-L-cysteine                                                                            | L 318 |
| Cytidine 5'-diphosphocholine sodium salt dihydrate                                                              | L 642 |
| D-(-)-α-Phenylglycine dane salt(potassium,ethyl) and D-(-)-p-Hydroxyphenylglycine dane salt (potassium, methyl) | L 364 |
| D(-)Pantolactone and DL Pantolactone                                                                            | L 553 |
| D-(-)-p-Hydroxy phenyl glycin                                                                                   | L 430 |
| D-(+)Glucose (dextrose), anhydrous                                                                              | L 597 |
| D-(+)Glucose monohydrate (Dextrose)                                                                             | L 102 |
| Dazomet                                                                                                         | L 575 |
| Decoquate                                                                                                       | L 627 |
| Dehydrated Iron(II)-Sulphate                                                                                    | L 214 |
| Demeclocycline hydrochloride hemihydrate                                                                        | L 655 |
| Dexpanthenol                                                                                                    | L 507 |
| Dextromethorphan hydrobromide                                                                                   | L 505 |
| Dextromethorphan hydrobromide x H <sub>2</sub> O (+)-3-Methoxy-17-methylmorphinan hydrobromide                  | L 512 |
| Dextrose (D-+)-Glucose monohydrate)                                                                             | L 102 |
| Dextrose anhydride                                                                                              | L 237 |
| Diacetyl (2,3-Butanedione)                                                                                      | L 496 |
| Diammonium citrate                                                                                              | L 018 |
| Diatrizoic acid dehydrate                                                                                       | L 529 |
| Dichlorodifluoromethane                                                                                         | L 345 |
| Dicyclopentadiene                                                                                               | L 158 |
| Diesel                                                                                                          | L 427 |
| Diet candy on maltitol basis                                                                                    | L 309 |
| Diethylene oxide                                                                                                | L 034 |
| Dimethyl sulfoxide                                                                                              | L 141 |
| Dimethyl terephthalate (DMT)                                                                                    | L 121 |
| Dimethyl-1-octadecylphosphonate                                                                                 | L 618 |

|                                                                                    |       |
|------------------------------------------------------------------------------------|-------|
| Dimethylformocarbthaldine                                                          | L 575 |
| Diphenyl mercury                                                                   | L 638 |
| Diphenylphosphinic acid                                                            | L 626 |
| Dipyron hydrate (Metamizol)                                                        | L 588 |
| Disinfectant                                                                       | L 198 |
| Disodium ethylene-bis-dithiocarbamate                                              | L 630 |
| Disodium hydrogen phosphate-12-hydrate                                             | L 489 |
| Disodium hydrogen phosphate-12-hydrate                                             | L 325 |
| DMT-dGuanosine(isobutyryl) β cyanoethyl phosphor amidite                           | L 547 |
| DNA-phosphoramidite                                                                | L 474 |
| Dobutamine hydrochloride                                                           | L 230 |
| Dodecylbenzene sulphonic acid                                                      | L 337 |
| Doxorubicinhydrochloride                                                           | L 562 |
| Dyes, extenders for dyes                                                           | L 221 |
| Egg white, dried                                                                   | L 163 |
| Egg yolk, liquid                                                                   | L 103 |
| Electropolishing solution (Hydrochloric acid-phosphoric acid-sulphuric acid 1:1:6) | L 600 |
| Engine oil                                                                         | L 201 |
| Engine oil, used                                                                   | L 109 |
| Engine oil, used                                                                   | L 110 |
| EPN, O-Ethyl-O-(4-nitrophenyl)phosphonothioate                                     | L 592 |
| Epoxy powder                                                                       | L 174 |
| Epoxy resin                                                                        | L 466 |
| Epoxy Styrol-Varnish                                                               | L 323 |
| Erythromycin                                                                       | L 242 |
| Ethoprophos                                                                        | L 579 |
| Ethosuximide (3-Ethyl-3-methyl-2,5pyrrolidindion)                                  | L 510 |
| Ethyl acetate                                                                      | L 013 |
| Ethyl thioglycolate                                                                | L 280 |
| Ethylene carbonate                                                                 | L 361 |
| Ethylene diamine (1,2-Diaminoethane)                                               | L 027 |
| Ethylene diamine and derivatives                                                   | L 026 |
| Ethylene dichloride                                                                | L 347 |
| Ethylene oxide                                                                     | L 203 |
| Ethylenediamine-tetraacetic acid                                                   | L 012 |
| Ethyleneimine homopolymer (Polyamine)                                              | L 279 |
| Fat, emulsifier, beta-carotene, aroma, vitamins                                    | L 457 |
| Fat, emulsifier, lecithin                                                          | L 457 |
| Fat, solid vegetable fat                                                           | L 098 |
| Fatamine, cocofatamine                                                             | L 131 |
| Fatamine, tallowfatamine                                                           | L 020 |
| Fat-emulsion                                                                       | L 218 |
| Fat-emulsion, whey fat                                                             | L 077 |
| Fatty alcohol sulphate sodium salt                                                 | L 434 |
| Feed, dried animal feed                                                            | L 445 |
| Fennel powder                                                                      | L 429 |
| Fertiliser                                                                         | L 152 |
| Fertiliser granulate                                                               | L 395 |
| Fibre glass fleece                                                                 | L 140 |
| Fibre glass insulating board                                                       | L 208 |
| Fibre glass wool                                                                   | L 206 |
| Films, photographic films                                                          | L 196 |

|                                                                        |       |
|------------------------------------------------------------------------|-------|
| Fipronil (Insecticide)                                                 | L 612 |
| Fish extract                                                           | L 479 |
| Fish food                                                              | L 235 |
| Flour and raising agent for rye wholemeal bred                         | L 374 |
| Fluorescein di-sodium salt                                             | L 535 |
| Fluoroaniline                                                          | L 019 |
| Fluorochloridone                                                       | L 581 |
| Fluticasone propionate                                                 | L 524 |
| Lithium hexafluorophosphate                                            | L 520 |
| Fluvastatin sodium                                                     | L 596 |
| Folic acid                                                             | L 383 |
| Formaldehyde                                                           | L 006 |
| Formic acid                                                            | L 171 |
| Fosfomycin sodium                                                      | L 571 |
| Fruit gum with sugar crust                                             | L 228 |
| Fruit gums based on a concentrate with non-standardized sugar compos.  | L 099 |
| Fruit gums based on saccharose                                         | L 100 |
| Fuel oil, heavy                                                        | L 111 |
| Gallic acid                                                            | L 275 |
| Gelatine                                                               | L 162 |
| Gel-emulsions of the cosmetic industry                                 | L 411 |
| Gelly sweets                                                           | L 402 |
| Gentamicin sulfate                                                     | L 448 |
| Ginger, ground                                                         | L 341 |
| Glass powder                                                           | L 415 |
| Gluconic acid (PLGA) mixture of Poly-lactate, Gluconic acid and Starch | L 475 |
| Glucose + saccharose syrup (Fondant)                                   | L 322 |
| Glucose molasses, syrup                                                | L 105 |
| Glucose polymer                                                        | L 483 |
| Glutamic acid, sodium salt                                             | L 008 |
| Glutardialdehyde 50 %                                                  | L 273 |
| Glycerine monostearate                                                 | L 297 |
| Glycerine monothioglycolate                                            | L 173 |
| Glycerol dibehenate                                                    | L 558 |
| Glycerol monooleate                                                    | L 497 |
| Glyoxal solution 40%                                                   | L 267 |
| Glyoxylic acid methylester methyl hemiacetal                           | L 261 |
| Glyphosate (Herbicide) N-(phosphonomethyl) glycine                     | L 327 |
| Grain                                                                  | L 438 |
| Grain, maize and winter wheat                                          | L 101 |
| Grease, polyglycol-based                                               | L 477 |
| Greases                                                                | L 412 |
| Guanidine hydrochloride                                                | L 035 |
| Gum arabic                                                             | L 021 |
| Gypsum, plaster of Paris                                               | L 042 |
| Hand balm                                                              | L 492 |
| Hand cream                                                             | L 492 |
| Hard coal distillate                                                   | L 277 |
| Hard fat (Solid vegetable fat)                                         | L 098 |
| Hardness oil, used                                                     | L 381 |
| Hazelnuts, ground                                                      | L 349 |
| Herbicide (Glyphosate) N-(phosphonomethyl) glycine                     | L 327 |
| Herbs, liquid extract                                                  | L 138 |

|                                                            |       |
|------------------------------------------------------------|-------|
| Hexa cyanoferrate(II)-3-hydrate                            | L 334 |
| Hexafluoroacetone sesquihydrate                            | L 154 |
| Hexahydrotriazine                                          | L 005 |
| Hexamethylene tetramine                                    | L 342 |
| Hexamethylenediamine                                       | L 025 |
| Hexetidine                                                 | L 217 |
| High-additive power transmission oil                       | L 161 |
| Homoveratryl amine                                         | L 300 |
| Honey                                                      | L 233 |
| Hydraulic fluid, used                                      | L 585 |
| Hydraulic oil                                              | L 107 |
| Hydrazine hydrate solution 80%                             | L 435 |
| Hydrazine monohydrate                                      | L 436 |
| Hydrofluoric acid                                          | L 051 |
| Hydrosilicon insulating board                              | L 209 |
| Hydrosilicon wool                                          | L 207 |
| Hydroxyl ammonium sulphate                                 | L 403 |
| Hydroxypropyl methylcellulose                              | L 382 |
| Hydroxyzine pamoate                                        | L 514 |
| Ice cream                                                  | L 523 |
| I-Leucovorin                                               | L 355 |
| Imidacloprid (Insecticide)                                 | L 611 |
| Sulating oil                                               | L 155 |
| Sulating oil                                               | L 462 |
| Iodine                                                     | L 441 |
| Iodomethane                                                | L 552 |
| Ipratropium bromide                                        | L 559 |
| Iron sulphate adherent water                               | L 216 |
| Iron(II)-chloride-4-hydrate                                | L 040 |
| Iron(II)-lactate hydrate                                   | L 469 |
| Iron(III) salts                                            | L 023 |
| Iron(III)-pyrophosphate                                    | L 252 |
| Iron(III)-sulphate, dehydrated                             | L 368 |
| Isomalt (mixture of Glucosylsorbitol and Glucosylmannitol) | L 372 |
| Jam, whole fruit containing sucrose                        | L 089 |
| Jellied fruits containing sucrose                          | L 100 |
| Kerosene                                                   | L 112 |
| Kieserite, magnesiumsulphate x H <sub>2</sub> O            | L 059 |
| L-2',3'-Dideoxy-3'-thiacytidine                            | L 628 |
| Lacquer                                                    | L 192 |
| Lacquer                                                    | L 452 |
| Lacquer PL, aluminised                                     | L 556 |
| Lactic acid 90-92 %                                        | L 272 |
| Lamivudine                                                 | L 628 |
| Lanolin                                                    | L 211 |
| Lauryl mercaptane                                          | L 299 |
| L-Cysteine-hydrochloride, dehydrated                       | L 317 |
| L-Cysteine-hydrochloride-monohydrate                       | L 145 |
| L-Cysteine-hydrochloride-monohydrate                       | L 316 |
| Lead acetate trihydrate                                    | L 653 |
| Leather                                                    | L 409 |
| Lecithin                                                   | L 393 |
| Lecithin, viscous                                          | L 333 |

|                                                                                                      |       |
|------------------------------------------------------------------------------------------------------|-------|
| Letostein                                                                                            | L 262 |
| L-Hydroxytryptophan                                                                                  | L 502 |
| Linseed oil                                                                                          | L 486 |
| Liquorice, hard                                                                                      | L 232 |
| Liquorice, liquid                                                                                    | L 088 |
| Lithium fluoride                                                                                     | L 641 |
| Lithium hexafluorophosphate                                                                          | L 520 |
| Lithium oxide using KF oven                                                                          | L 656 |
| L-Lysine monohydrate                                                                                 | L 528 |
| Lubricating oil                                                                                      | L 545 |
| Lyo-Corks                                                                                            | L 378 |
| Lyophilized material                                                                                 | L 135 |
| Magnesium aspartate-4-hydrate                                                                        | L 498 |
| Magnesium di-aspartate                                                                               | L 270 |
| Magnesium oxide                                                                                      | L 590 |
| Magnesium stearate                                                                                   | L 326 |
| Magnesium sulphate monohydrate                                                                       | L 059 |
| Magnesiumacetate-2-hydrate                                                                           | L 225 |
| Magnesiumcarbonate 90                                                                                | L 573 |
| Magrocol 8000                                                                                        | L 633 |
| Maleic acid                                                                                          | L 624 |
| Maleic acid/Acrylic acid copolymer                                                                   | L 391 |
| Malic acid                                                                                           | L 120 |
| Malonic acid                                                                                         | L 011 |
| Maltitol, aqueous solution                                                                           | L 560 |
| Maltodextrin                                                                                         | L 388 |
| Maltodextrin with flavours                                                                           | L 555 |
| Mancozeb (Fungicide)                                                                                 | L 614 |
| Manganese dioxide                                                                                    | L 159 |
| Mannitol                                                                                             | L 210 |
| Margarine                                                                                            | L 083 |
| Marzipan                                                                                             | L 080 |
| Marzipan using a homogenizer                                                                         | L 293 |
| Mayonnaise                                                                                           | L 084 |
| MCPA sodium monohydrate                                                                              | L 589 |
| Mercaptopurine                                                                                       | L 282 |
| Mercury(II)diphenyl                                                                                  | L 638 |
| Meropenem trihydrate                                                                                 | L 654 |
| meta-Aminoparacresol                                                                                 | L 321 |
| Metamizol                                                                                            | L 588 |
| Metformin hydrochloride                                                                              | L 631 |
| Methampyrone                                                                                         | L 588 |
| Methanesulfonyl chloride                                                                             | L 373 |
| Methotrexate N-[4-[[[(2,4-Diamino-6-pteridiny)]methyl]methylamino] benzoyl]-L-gluta mic acid         | L 460 |
| Methyl 2-[[[N-(4-methoxy-6-methyl-1,3,5-triazin-2-yl)methylamino] carbonyl] amino]sulfonyl] benzoate | L 570 |
| Magnesium aspartate-4-hydrate                                                                        | L 498 |
| Magnesium di-aspartate                                                                               | L 270 |
| Magnesium oxide                                                                                      | L 590 |
| Magnesium stearate                                                                                   | L 326 |
| Magnesium sulphate monohydrate                                                                       | L 059 |
| Magnesiumacetate-2-hydrate                                                                           | L 225 |

|                                                                                                      |       |
|------------------------------------------------------------------------------------------------------|-------|
| Magnesiumcarbonate 90                                                                                | L 573 |
| Magrocol 8000                                                                                        | L 633 |
| Maleic acid                                                                                          | L 624 |
| Maleic acid/Acrylic acid copolymer                                                                   | L 391 |
| Malic acid                                                                                           | L 120 |
| Malonic acid                                                                                         | L 011 |
| Maltitol, aqueous solution                                                                           | L 560 |
| Maltodextrin                                                                                         | L 388 |
| Maltodextrin with flavours                                                                           | L 555 |
| Mancozeb (Fungicide)                                                                                 | L 614 |
| Manganese dioxide                                                                                    | L 159 |
| Mannitol                                                                                             | L 210 |
| Margarine                                                                                            | L 083 |
| Marzipan                                                                                             | L 080 |
| Marzipan using a homogenizer                                                                         | L 293 |
| Mayonnaise                                                                                           | L 084 |
| MCPA sodium monohydrate                                                                              | L 589 |
| Mercaptopurine                                                                                       | L 282 |
| Mercury(II)diphenyl                                                                                  | L 638 |
| Meropenem trihydrate                                                                                 | L 654 |
| meta-Aminoparacresol                                                                                 | L 321 |
| Metamizol                                                                                            | L 588 |
| Metformin hydrochloride                                                                              | L 631 |
| Methampyrone                                                                                         | L 588 |
| Methanesulfonyl chloride                                                                             | L 373 |
| Methotrexate N-[4-[[[(2,4-Diamino-6-pteridiny)]methyl]methylamino] benzoyl]-L-gluta mic acid         | L 460 |
| Methyl 2-[[[N-(4-methoxy-6-methyl-1,3,5-triazin-2-yl)methylamino] carbonyl] amino]sulfonyl] benzoate | L 570 |
| Methylcyano acrylate                                                                                 | L 420 |
| Metiram (Pesticide)                                                                                  | L 648 |
| Milk                                                                                                 | L 085 |
| Milk, chocolate                                                                                      | L 079 |
| Milk powder, dry whole milk                                                                          | L 081 |
| Milk powder, skimmed                                                                                 | L 082 |
| Mineral oil with additives                                                                           | L 190 |
| Mineral oil without additives                                                                        | L 184 |
| Mixture of Magnesium di-aspartate and Potassium aspartate                                            | L 268 |
| Molasses                                                                                             | L 227 |
| Monosultap (Pesticide)                                                                               | L 643 |
| Multivitamin capsules (contents of Multivitamin capsules)                                            | L 439 |
| Multivitamin tablets and powder                                                                      | L 404 |
| N-(phosphonomethyl) glycine Glyphosate (Herbicide)                                                   | L 327 |
| N,N-Dimethyl acetamide                                                                               | L 014 |
| N,N-Dimethyl formamide                                                                               | L 424 |
| N4-phthalylsulfathiazole                                                                             | L 574 |
| Nabam (algicide, fungicide)                                                                          | L 630 |
| N-Acetyl-L-cysteine                                                                                  | L 315 |
| Naphthol yellow S                                                                                    | L 576 |
| n-Butylnitrite                                                                                       | L 407 |
| N-Butyraldehyde and iso-butyraldehyde                                                                | L 248 |
| N-dodecylmercaptane                                                                                  | L 299 |
| Neomycin B trisulfate hydrate                                                                        | L 644 |

|                                                                           |       |
|---------------------------------------------------------------------------|-------|
| N-Formylmorpholine (4-Morpholine carboxaldehyde)                          | L 392 |
| Nickel sulphamate 4-hydrate                                               | L 651 |
| Niclosamide                                                               | L 369 |
| Nitarstone                                                                | L 602 |
| Nitrogen bases                                                            | L 288 |
| Nitromethane                                                              | L 599 |
| N-Methylmorpholine/N-Oxide cellulose mixture                              | L 390 |
| N-Methylpyrrolidone                                                       | L 549 |
| Noodles                                                                   | L 076 |
| N-Vinyl-2-pyrrolidone (stabilised)                                        | L 371 |
| O-(2,6-dichloro-4-methylphenyl) O,O-diethyl phosphorothioate              | L 616 |
| O-2,6-dichloro-p-tolyl O,O-dimethyl phosphorothioate                      | L 616 |
| Octamethylcyclotetrasiloxane                                              | L 519 |
| O-Ethyl O-(4-nitrophenyl) phenylphosphonothioate                          | L 592 |
| O-Ethyl S,S-dipropyl phosphorodithioate                                   | L 579 |
| Oil of oenothera tetragonal                                               | L 447 |
| Oil, cold compressor oil with freon                                       | L 534 |
| Oil, crude oil                                                            | L 108 |
| Oil, engine oil used                                                      | L 109 |
| Oil, engine oil used                                                      | L 110 |
| Oil, fish oil silage                                                      | L 470 |
| Oil, hardening oil for steel                                              | L 160 |
| Oil, hardness oil used                                                    | L 381 |
| Oil, heavy fuel oil                                                       | L 111 |
| Oil, hydraulic oil                                                        | L 107 |
| Oil, insulating oil                                                       | L 462 |
| Oil, linseed oil                                                          | L 486 |
| Oil, mineral oil with additives                                           | L 190 |
| Oil, mineral oil without additives                                        | L 184 |
| Oil, peanut oil                                                           | L 164 |
| Oil, rape seed oil                                                        | L 319 |
| Oil, silicone oil                                                         | L 113 |
| Oil, silicone oil                                                         | L 156 |
| Oil, transformer                                                          | L 462 |
| Oil, vegetable                                                            | L 063 |
| Oil, vegetable                                                            | L 075 |
| Oils using a KF oven                                                      | L 106 |
| Oils with additives                                                       | L 531 |
| Ointment, oral                                                            | L 476 |
| Ointments (Suppository compounds)                                         | L 136 |
| Olefin copolymer in mineral oil                                           | L 185 |
| Organic plastic                                                           | L 401 |
| O-Tolidine                                                                | L 236 |
| Oxalic acid                                                               | L 385 |
| Oxytetracyclinhydrochloride                                               | L 563 |
| Paints, dispersion paints and binding agents                              | L 191 |
| PAPP 2,2-bis[(r-aminophenoxy)phenyl]propane                               | L 451 |
| Parafin liquid                                                            | L 419 |
| Paraformaldehyde                                                          | L 386 |
| Paraquat dichloride x hydrate (1,1'-Dimethyl-4,4'-bipyridiniumdichloride) | L 504 |
| Peanut oil                                                                | L 164 |
| Penicillin                                                                | L 137 |
| Penicillin-G sodium and Penicillin-G-potassium                            | L 166 |

|                                                                       |       |
|-----------------------------------------------------------------------|-------|
| Pepper, powdered                                                      | L 432 |
| Pepper, white                                                         | L 433 |
| Peppermint oil and spearmint oil                                      | L 456 |
| Peppermint oil and spearmint oil                                      | L 406 |
| Perchloric acid in glacial acetic acid                                | L 517 |
| Peroxide derivates                                                    | L 279 |
| Petrol, unleaded                                                      | L 428 |
| Phenol-urea-formaldehyde-condensate                                   | L 130 |
| Phenyl glyoxal hydrate 97 %                                           | L 478 |
| Phenyl hydrazine                                                      | L 151 |
| Phosphine carboxylacid, homopolymer                                   | L 541 |
| Phosphoserinic acid, anhydrous                                        | L 269 |
| Plastic film                                                          | L 526 |
| Plastics with the Karl Fischer oven                                   | L 328 |
| PLEX 5372 E (Polyelectrolyte, high molecular I)                       | L 257 |
| PLEX 5477 E (Polyelectrolyte, high molecular II)                      | L 356 |
| PLGA, mixture of Poly-lactate,Gluconic acid and Starch                | L 475 |
| Poloxamers                                                            | L 521 |
| Polyacrylamide-copolymerisates                                        | L 246 |
| Polyacrylate (liquid) with sodium silicate                            | L 253 |
| Polyacrylic acid                                                      | L 266 |
| Polyacrylic acid                                                      | L 314 |
| Polyamide                                                             | L 126 |
| Polyamide 6                                                           | L 174 |
| Polyamide 66                                                          | L 167 |
| Polyamide 66                                                          | L 174 |
| Polyamide dissolved in NMP, xylene and cyclohexanol                   | L 226 |
| Polyamine                                                             | L 379 |
| Polybutadiene terephthalate                                           | L 174 |
| Polybutene                                                            | L 188 |
| Polycarbonate                                                         | L 127 |
| Polycarbonate                                                         | L 129 |
| Polycarbonate diol                                                    | L 533 |
| Polyelectrolyte, high molecular I (PLEX 5372 E)                       | L 257 |
| Polyelectrolyte, high molecular I I (PLEX 5372 E)                     | L 256 |
| Polyester foil                                                        | L 125 |
| Polyester Resin LP 24/23 for the lime industry                        | L 239 |
| Polyetherpolyol                                                       | L 377 |
| Polyethylene                                                          | L 193 |
| Polyethylene cross-linked                                             | L 128 |
| Polyethylene cross-linked                                             | L 640 |
| Polyethylene glycol 8000                                              | L 633 |
| Polyethylene imine                                                    | L 522 |
| Polyglycol-based grease with additives                                | L 477 |
| Polyisobutylene 900                                                   | L 002 |
| Poly-lactate (PLGA) mixture of Poly-lactate, Gluconic acid and Starch | L 475 |
| Poly-L-lactate PLLA (PLA)                                             | L 577 |
| Polymer comprising maleic and acrylic acid                            | L 199 |
| Polymers, various                                                     | L 174 |
| Polymethyl methacrylate                                               | L 215 |
| Polyol                                                                | L 122 |
| Polyol                                                                | L 182 |
| Polyol with 20% polyurea                                              | L 551 |

|                                                        |       |
|--------------------------------------------------------|-------|
| Polyphenol                                             | L 123 |
| Polypropylene                                          | L 194 |
| Polypropylene filled                                   | L 181 |
| Polysorbate 80                                         | L 608 |
| Polyurethane - chlorinated                             | L 124 |
| Polyurethane component A for adhesives                 | L 119 |
| Polyurethane granules, thermoplastic (TPU)             | L 303 |
| Polyvinyl alcohol                                      | L 490 |
| Polyvinylacetate-based dispersion                      | L 278 |
| Popcorn                                                | L 292 |
| Potassium aspartate and Magnesium di-aspartate mixture | L 268 |
| Potassium chloride                                     | L 625 |
| Potassium clavulanate                                  | L 437 |
| Potassium cyanide                                      | L 258 |
| Potassium hexacyanoferrate(II)-3-hydrate               | L 334 |
| Potassium hexafluorophosphate                          | L 247 |
| Potassium hydrogen aspartate                           | L 499 |
| Potassium hydrogen carbonate                           | L 052 |
| Potassium hydrogen fluoride                            | L 038 |
| Potassium hydrogen fluoride                            | L 060 |
| Potassium magnesium chloride x 6H <sub>2</sub> O       | L 059 |
| Potassium tetraoxalate-2-hydrate                       | L 543 |
| Potato chips                                           | L 368 |
| Potato fibre                                           | L 440 |
| Potato starch                                          | L 170 |
| Potato-raw material for instant                        | L 090 |
| Povidon                                                | L 604 |
| Prilocaine                                             | L 561 |
| Printing ink                                           | L 202 |
| Procymidone (Fungicide)                                | L 613 |
| Proflavine hemisulphate                                | L 354 |
| Promethazine hydrochloride                             | L 582 |
| Propylene carbonate                                    | L 362 |
| Propylene oxide                                        | L 203 |
| Pteric acid                                            | L 594 |
| PVC                                                    | L 204 |
| PYCA (4-Pyridiniummethyl-7-amino cephalosporin-acid)   | L 414 |
| Pyridoxal-5-phosphate                                  | L 143 |
| Pyrogallol                                             | L 004 |
| Pyromellitic dianhydride PMDA                          | L 449 |
| Rape seed oil                                          | L 319 |
| Raw coffee                                             | L 335 |
| Riboflavin phosphate sodium                            | L 509 |
| Roll-on deodorant                                      | L 453 |
| Rubber, vulcanized                                     | L 176 |
| Rusk with the Karl Fischer oven                        | L 332 |
| Rusks                                                  | L 064 |
| Rusks using a KF oven                                  | L 332 |
| Saccharose                                             | L 073 |
| Saccharose + glucose syrup (Fondant)                   | L 332 |
| S-Adenosyl-L-methioninsulfate sylate                   | L 527 |
| Salbutamol hemisulphate                                | L 565 |
| Salicyl-anthranilic acid                               | L 629 |

|                                                           |       |
|-----------------------------------------------------------|-------|
| Sarafloxacin hydrochloride                                | L 593 |
| Sausage, veal                                             | L 065 |
| S-Carboxymethyl-L-cysteine                                | L 318 |
| Seal cover, cold seal                                     | L 432 |
| Sealants                                                  | L 646 |
| Semolina, soft wheat                                      | L 331 |
| Shampoo                                                   | L 175 |
| Shower gel                                                | L 175 |
| SI enamel                                                 | L 116 |
| Silage - fish oil                                         | L 470 |
| Silanol and siloxanes                                     | L 057 |
| Silica sand                                               | L 142 |
| Silicone oil                                              | L 113 |
| Silicone oil                                              | L 156 |
| Silicone oil, hydride functional                          | L 623 |
| Silicone sealing agent (single-component) interim product | L 621 |
| Siloxanes                                                 | L 180 |
| Silver nitrate                                            | L 645 |
| Sodium (4-chloro-2-methylphenoxy) acetate                 | L 589 |
| Sodium acetate                                            | L 550 |
| Sodium azide                                              | L 037 |
| Sodium benzoate                                           | L 339 |
| Sodium chloride                                           | L 304 |
| Sodium dichromate-2-hydrate                               | L 425 |
| Sodium dodecyl sulphate                                   | L 485 |
| Sodium formate                                            | L 538 |
| Sodium glutamate-1-hydrate                                | L 274 |
| Sodium glycerophosphate hydrated                          | L 632 |
| Sodium hydrogencarbonate                                  | L 053 |
| Sodium hypophosphite monohydrate                          | L 036 |
| Sodium methylate                                          | L 287 |
| Sodium molybdate                                          | L 231 |
| Sodium perchlorate                                        | L 376 |
| Sodium polyphosphate 85 %                                 | L 370 |
| Sodium selenite pentahydrate                              | L 605 |
| Sodium tartrate-2-hydrate                                 | L 200 |
| Sodiumcacodylate x-hydrate                                | L 572 |
| Sodium-L-glutamate monohydrate                            | L 070 |
| Sorbitan monostearate                                     | L 495 |
| Sorbitol                                                  | L 320 |
| Sorbitol solution, 70 %                                   | L 607 |
| Sour dough extract                                        | L 358 |
| Soya grist                                                | L 153 |
| Sparfloxacin                                              | L 598 |
| Spearmint oil and peppermint oil                          | L 456 |
| Sponge fingers                                            | L 087 |
| β-Amidinocrotonic acid amide                              | L 283 |
| Starch                                                    | L 072 |
| Starch syrup                                              | L 213 |
| Strawberries (powdered and pieces)                        | L 658 |
| Strontium hydroxide octahydrate                           | L 647 |
| Succinic acid                                             | L 580 |
| Sugars, various                                           | L 168 |

|                                                                               |       |
|-------------------------------------------------------------------------------|-------|
| Sulfacetamide sodium                                                          | L 461 |
| Sulfonyl urea                                                                 | L 408 |
| Sulfuric acid mixtures with hydrofluoric acid                                 | L 050 |
| Sulfuric acid, concentrated                                                   | L 049 |
| Sulphur dioxide, liquid                                                       | L 340 |
| Sulphurous oxychloride (Thionyl chloride)                                     | L 503 |
| Sultanas                                                                      | L 291 |
| Suppository compounds (Ointments)                                             | L 136 |
| Surface preservative, decking protector                                       | L 539 |
| Suxamethonium chloride                                                        | L 511 |
| Sweet mixture, candies                                                        | L 220 |
| Sweets, gelly sweets                                                          | L 402 |
| Talcum powder                                                                 | L 229 |
| Tallowfat amine                                                               | L 020 |
| Tannic acid                                                                   | L 446 |
| Terephthalic acid                                                             | L 007 |
| Tergitol NPEO                                                                 | L 567 |
| Tert.-Butyl-2-bromo-iso-butyrate                                              | L 302 |
| Tetrachloroisophthalonitrile                                                  | L 351 |
| Tetrachloromethane containing chlorine                                        | L 001 |
| Tetrachloro-p-benzoquinone                                                    | L 617 |
| Tetrahydro-3,5-dimethyl-1,3,5-thiadiazine-2-thione                            | L 575 |
| Tetrahydrofurane (Diethylene oxide)                                           | L 034 |
| Tetramethylthioperoxydicarbonic diamide                                       | L 615 |
| Thiacyclopentylidene-3-cyanomethylacetate                                     | L 285 |
| Thiamine hydrochloride (Thiamine chloride hydrochloride)                      | L 513 |
| Thiocol                                                                       | L 245 |
| Thioglycollic acid                                                            | L 265 |
| Thionyl chloride (Sulphurous oxychloride)                                     | L 503 |
| Thiophanate-methyl (Fungicide)                                                | L 650 |
| Thiosulfuric acid S,S'-[2-(dimethylamino) trimethylene] ester monosodium salt | L 643 |
| Thiosultap monosodium                                                         | L 643 |
| Thiourea                                                                      | L 639 |
| Thiram (Fungicide)                                                            | L 615 |
| Tilmicosin                                                                    | L 601 |
| Tin(II)chloride-2-hydrate and Tin(IV) chloride-5-hydrate                      | L 387 |
| Tin(II)-methanesulfonate                                                      | L 578 |
| Tobacco                                                                       | L 134 |
| Tobacco extract                                                               | L 443 |
| Tobacco ribs                                                                  | L 177 |
| Tobramycine sulphate                                                          | L 398 |
| Toffee                                                                        | L 069 |
| Toffee creme                                                                  | L 224 |
| Toffee methyl (Fungicide)                                                     | L 616 |
| Toluylen-2,4-diisocyanate                                                     | L 305 |
| Toothpaste                                                                    | L 029 |
| TPU (Polyurethane granules, thermoplastic)                                    | L 303 |
| Transformer oil                                                               | L 462 |

|                                                                                                      |       |
|------------------------------------------------------------------------------------------------------|-------|
| Transmission oil, synthetic                                                                          | L 564 |
| Tri-basic calcium phosphate                                                                          | L 442 |
| Tribenuron-methyl                                                                                    | L 570 |
| Triethyl citrate                                                                                     | L 525 |
| Trifluoroacetic acid and Trichloroacetic acid                                                        | L 380 |
| Trifluoromethanesulfonic acid Aluminum salt                                                          | L 635 |
| Tri-magnesium phosphate                                                                              | L 250 |
| Trimethoprim                                                                                         | L 587 |
| Trioxane                                                                                             | L 033 |
| Triphenyl phosphine                                                                                  | L 491 |
| Tri-Potassium citrate-1-hydrate                                                                      | L 241 |
| Tris(aminzinc-ethylene-bis-(dithiocarbamate)-tetrahydro-1,2,4,7-dithiadiazocine-3,8-dithione-polymer | L 648 |
| Tris-(hydroxymethyl) aminomethane                                                                    | L 022 |
| Tri-Sodium citrate and Mono sodium citrate                                                           | L 240 |
| Tri-Sodium citrate-2-hydrate                                                                         | L 169 |
| Triton X105                                                                                          | L 566 |
| Truffles and Chocolates, soft                                                                        | L 028 |
| Tungsten metal powder                                                                                | L 554 |
| Turbine oil, used                                                                                    | L 584 |
| Uranine                                                                                              | L 535 |
| Urea                                                                                                 | L 157 |
| Urea-formaldehyde-condensates                                                                        | L 219 |
| Vaccine, animal CHPI                                                                                 | L 421 |
| Vinyl acetate                                                                                        | L 264 |
| Vinyl acetate                                                                                        | L 271 |
| Vinyl chloride                                                                                       | L 344 |
| Vinyl methyl acetamide                                                                               | L 367 |
| Warfarin sodium clathrate                                                                            | L 400 |
| Wax for cosmetic use (1)                                                                             | L 132 |
| Wax for cosmetic use (2)                                                                             | L 133 |
| Wheat flour                                                                                          | L 068 |
| Wheat flour with the Karl Fischer oven                                                               | L 329 |
| Wheat grits                                                                                          | L 067 |
| Wheat semolina, soft                                                                                 | L 331 |
| Whey fat emulsion                                                                                    | L 077 |
| Whey, powdered                                                                                       | L 493 |
| Whole fruit jam containing saccharose                                                                | L 079 |
| Window sealant                                                                                       | L 399 |
| Wine gums                                                                                            | L 066 |
| Xylite                                                                                               | L 324 |
| Yeast                                                                                                | L 384 |
| Yoghurt                                                                                              | L 097 |
| Yttrium trifluoroacetate                                                                             | L 542 |
| Zinc diamyl dithiocarbamate                                                                          | L 405 |
| Zinc diethyl dithiophosphate                                                                         | L 189 |
| Zinc oxide                                                                                           | L 310 |
| Zinc peroxide                                                                                        | L 311 |
| Zinc powder                                                                                          | L 078 |

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