



# **Preva Dental X-Ray System**

## User Manual

**00-02-1576 Rev. D**

**ECN: P1244**

Progeny Dental  
1407 Barclay Boulevard  
Buffalo Grove, Illinois 60089 U.S.A.  
Phone: (888) 924-3800 Fax: (847) 459-5175  
[WWW.progenydental.com](http://WWW.progenydental.com)



## Table of Contents


General Information.....	1
Product Description.....	1
Compliance with Applicable Standards.....	3
Certified Components .....	3
EC Declaration of Conformity.....	4
Authorized Representatives .....	5
Safety .....	5
Explanation of Symbols on Technical Labels.....	6
Obtaining Technical Support.....	6
Operating the Preva Dental X-Ray System .....	7
Using the Operator Panel .....	7
Taking an X-Ray .....	9
Using the 12 in Cone (30-A2033).....	9
Recommended Maintenance .....	10
Regular Maintenance.....	10
Cleaning and Disinfecting .....	10
Checking System Functions.....	11
System Function Checklist.....	12
New Tube Seasoning Procedure.....	13
Solving Performance Issues.....	14
Performance Issues .....	14
Obtaining Technical Support.....	14
Pre-programmed Exposure Times .....	15
System Configuration .....	17
System Configuration Mode.....	17
Adjusting the Display .....	18
Changing Pre-programmed Exposure Settings .....	19
Preprogramming to Digital Sensors.....	20
Showing Current System Configuration .....	22
Changing the Cone Size .....	23
Diagnostic Mode .....	24
Specifications .....	25
Preva Dental X-Ray System .....	25

## General Information

### Product Description

	<p>The Preva Dental X-Ray System is a state-of-the-art, high-frequency intra-oral x-ray machine. The Preva consists of five components, as shown in Figure 1: the Control Unit, the Tubehead, the Articulating Arm, the Horizontal Arm, the Cone, and the Remote Control option.</p>
<b>Control Unit</b>	<p>The Control Unit provides for the input power connection and control of the Tubehead and Operator Panel. It provides automatic line voltage compensation, kVp control and exposure time control. The Control Unit consists of the mounting base and Operator Panel.</p>
<b>Tubehead</b>	<p>The Tubehead contains the x-ray tube, high voltage circuit, and Cone. The Tubehead is shipped already assembled to the Articulating Arm.</p> <p><b>Note:</b> There is a small hole in the plastic handle covering the back of the Tubehead. Under no circumstances should this hole be blocked as it provides an air vent to allow the Tubehead oil to expand and contract as the unit is operated.</p>
<b>Articulating Arm</b>	<p>The Articulating Arm provides the articulation support for the Tubehead and the reach and coverage of the Tubehead to the patient. The Articulating Arm allows smooth movement for precise positioning and does not drift or vibrate when left in position.</p>
<b>Horizontal Arm</b>	<p>The Horizontal Arm helps provide the necessary reach for the Preva. The Horizontal Arm pivots smoothly around a shaft inserted in the top of the Control Unit. The Horizontal Arm contains an access cover to connect the cable from the Horizontal Arm to the Control Unit. The Horizontal Arm is available in three lengths, providing reaches of 56, 66 and 76 inches.</p>
<b>Cone</b>	<p>The Cone establishes the distance from the x-ray tube to the patient's skin. It provides positioning assistance and collimates the x-ray beam to within a defined circle at its end. The Preva is shipped with the standard 8 inch Cone attached to the Tubehead. A 12 inch Cone (30-A2033) can be ordered as an option.</p>
<b>Remote Control</b>	<p>An optional component, the remote control switch is used to make exposures in addition to or replacing the use of the exposure button.</p>
<b>Installation and Service</b>	<p>The Preva Dental X-Ray System should only be installed and serviced by approved Progeny dealer personnel. Contact Progeny at (888) 924-3800 if you need assistance locating an approved dealer.</p>

**Figure 1**  
**Component**  
**Diagram**



## Compliance with Applicable Standards

**Radiation Protection**

The certified components of the Preva Dental X-Ray System comply with Radiation Performance Standards 21 CFR, Subchapter J, at the time of manufacture.

The certified components of the Preva Dental X-Ray System comply with IEC 60601-1-3 Radiation protection/x-ray equipment.

**UL 2601-1 File Number: E181750**

Classified by Underwriters Laboratories Inc. with respect to electrical shock, fire and mechanical hazards only in accordance with UL 2601-1, and CAN/CSA C22.2 NO, 601.1-M90, and to the following particular standards, IEC60601-2-7, IEC60601-2-28.

**EMI/EMC**

IEC60601-1-2

## Certified Components

Component	Reference Number
Tubehead	30-A1027
Control Unit	30-A0010
Cone 8 in.	30-A2016
Cone 12 in.	30-A2033
Cone 8 in. Rectangular	30-A2041
Cone 6mm	30-A2101

## EC Declaration of Conformity

### Name and Description of Product

#### Progeny Preva

Catalog Model	P7017, 76 inch reach 30-A0010, Control 30-A2071, Extension Arm, Long
Catalog Model	P7016, 66 inch reach 30-A0010, Control 30-A2073, Extension Arm, Short
Catalog Model	P7015, 56 inch reach 30-A0010, Control 30-A2074, Extension Arm, Compact
Class: IIb	

### Reference Numbers to which Conformity is Declared

The following regulatory documents apply:  
 UL 2601-1  
 IEC 60601-1-2  
 IEC 60601-1-3  
 IEC 60601-2-7  
 IEC 60601-2-28  
 IEC 60601-2-32  
 Medical Device directive  
 ISO 13485  
 EN46001

### Declaration

Progeny, Inc. declares that the products described herein meet all the applicable Essential Requirements of the EC Medical Device Directive 93/42/EEC in Annex I. For Class IIb products described herein, the product is manufactured, inspected, tested and released in accordance with the approved quality assurance system established in accordance with ISO 13485 and Annex II of the EC Medical Device Directive under the Supervision of the SGS United Kingdom Ltd., a Notified Body.

### Contact

Alan Crema  
 Product Development  
 Progeny, Inc.

## Authorized Representatives

### North America

PROGENY DENTAL  
1407 Barclay Blvd.  
Buffalo Grove, IL 60089  
Phone: 888-924-3800  
Fax: 847-459-5175  
techsupport@progenydental.com

### Europe

CE Partner 4U  
Esdoornlaan 13  
3951DB Maarn  
The Netherlands  
Phone: +31.343.442.524  
Fax: +31.343.442.162

## Safety

### Radiation Safety

Only qualified and authorized personnel may operate this equipment observing all laws and regulations concerning radiation protection.

- The operator at all times must remain 6ft. (2m) from the focal spot and the X-ray beam for operator protection.
- Full use must be made of all radiation safety features on the equipment.
- Full use must be made of all radiation protection devices, accessories and procedures available to protect the patient and operator from x-ray radiation.

### Electrical Safety







- Only qualified and authorized service personnel should remove covers on the equipment.
- This equipment must only be used in rooms or areas that comply with all applicable laws and recommendations concerning electrical safety in rooms used for medical purposes, e.g., IEC, US National Electrical code, or VDE standards concerning provisions of an additional protective earth (ground) terminal for power supply connection.
- Before cleaning or disinfecting, this equipment must always be disconnected from the main electrical supply.
- The Preva Dental X-Ray System is ordinary type medical equipment without protection against ingress of liquids. To protect against short-circuit and corrosion, no water or any other liquid should be allowed to leak inside the equipment.

### Explosion Safety

This equipment must not be used in the presence of flammable or potentially explosive gases or vapors, which could ignite, causing personal injury and/or damage to the equipment. If such disinfectants are used, the vapor must be allowed to disperse before using the equipment.



## Explanation of Symbols on Technical Labels

	Type B: Protection against electric shock (IEC 60601.1-1988)
	Consult written instructions in User's Manual.
	ATTENTION RAYONS-X: OPERATION SEULEMENT PAR DU PERSONNEL AUTORISE. VOIR MANUEL DE L'OPERATEUR.
	WARNING X-RAY THIS X-RAY UNIT MAY BE DANGEROUS TO PATIENT AND OPERATOR UNLESS SAFE EXPOSURE FACTORS AND OPERATING INSTRUCTIONS ARE OBSERVED.
	<i>X-RAY EMISSION</i>
<b>L</b>	Mains HOT WIRE
<b>N</b>	Mains NEUTRAL WIRE
	Earth Ground

## Obtaining Technical Support

### Contact


PROGENY DENTAL  
1407 Barclay Blvd.  
Buffalo Grove, IL 60089  
Phone: 888-924-3800  
Fax: 847-459-5175  
techsupport@progenydental.com

# Operating the Preva Dental X-Ray System

## Using the Operator Panel

### Power On Settings

When the Preva Dental X-Ray System is powered on, the Operator Panel selections are those that were in use when the system was last powered off.



**Figure 2**  
**Preva Operator**  
**Panel**

**Exposure Settings**

When the system is powered on, the Operator Panel, Figure 2, displays the exposure settings (kV, mA, and seconds) for the currently selected tooth, image receptor type, and patient size. Use the Tooth Selection, Image Receptor Type, and Patient Size buttons to select other exposure settings.

- For a table of the factory-programmed exposure settings, refer to the Pre-programmed Exposure Settings tables later in this manual.

**Adjusting Exposure Settings**

Preset exposure settings can be adjusted prior to making an exposure. Use the right arrow to select the exposure setting to adjust. Then use the up and down arrow buttons to adjust the value.

- To save new presets, use System Configuration mode described later in this manual.

**Exposure Button and Ready Indicator**

The Exposure button is used to initiate an x-ray exposure. For a complete exposure, the button must be pressed and held until the Radiation Indicator no longer illuminates and the audible signal is no longer heard. Releasing the Exposure button immediately terminates the x-ray exposure.

**CAUTION!**

**Releasing the Exposure button prior to the completion of the x-ray exposure will result in an incomplete exposure of the image. This may require the operator to re-take the radiograph. When a premature release of the Exposure button occurs, the system will notify the operator momentarily, then return to operating mode.**

**Ready Indicator**

The Ready Indicator illuminates when the system is ready to make an exposure. Immediately after an exposure, the Ready Indicator flashes until the x-ray tube cools down sufficiently to make the next exposure. When the Ready Indicator is flashing, no exposure can be made.

**Radiation Indicators**

The Preva has a visible and an audible Radiation Indicator. When an exposure is in progress, the Radiation Indicator on the Operator Panel is illuminated and an audible tone is heard. The exposure is complete when the Radiation Indicator is extinguished and the audible tone is no longer heard.

## Taking an X-Ray

1. Turn the power switch, located at the upper right of the Control Unit, to the "On" position. The Ready Indicator on the front of the Operator Panel, Figure 2, will light.
2. Verify that the unit is set for the correct Image Receptor Type. The icon for the currently selected Image Receptor Type is illuminated. To change the Image Receptor type, press the Image Receptor Type button until the correct Image Receptor Type is selected.
3. Verify that the system is set for the appropriate Patient Size. The icon for the currently selected Patient Size is illuminated. To change the Patient Size, press the Patient Size button until the correct Patient Size is selected.
4. Verify that the unit is set for the Tooth to be imaged. The icon for the currently selected Tooth is illuminated. To change the Tooth Selection, press the Tooth Selection button until the correct Tooth is selected.
5. If desired, preset exposure settings for the combination of Image Receptor Type, Tooth Selection, and Patient size, selected in steps 2-4, can be adjusted prior to making an exposure. Use the right arrow to select the exposure setting to adjust. Then use the up and down arrow buttons to adjust the value. Skip this step if you are using pre-programmed exposure settings.  
**Note:** When exposure settings are being adjusted, the Tooth Selection, Image Receptor Type, and Patient Size buttons are turned off.
6. Position the Tubehead to the patient using standard accepted positioning procedures.
7. Press and hold the Exposure button until the audible signal is no longer heard and the Radiation Indicator is no longer illuminated. Releasing the Exposure button or coil-cord hand switch at any time will immediately terminate the exposure.  
**Note:** When using the coil-cord hand switch, it is recommended that the operator exit the operatory if possible.  
**Note:** In order to comply with regulations and good safety practices, the technique factors must be visible to the operator from the remote location.
8. Return the Tubehead to the storage position.  
**Note:** Be careful not to strike the Tubehead on the wall when returning it to the storage position.

## Using the 12 in Cone (30-A2033)

The Preva Dental X-Ray System is factory set for use with the standard supplied 8 inch (20 cm) Cone. The 12 inch (30 cm) Cone (30-A2033) is available. Using the longer cone requires longer exposure times. See the System Configuration section of this manual for instructions to set the system to use the longer cone.

## Recommended Maintenance

### Regular Maintenance

In the interest of equipment safety, a regular maintenance program must be established. This maintenance program should consist of cleaning and disinfecting as well as annual system function checking. It is the owner's responsibility to arrange for this service and to assure that the personnel performing this are fully qualified to service Progeny Dental x-ray equipment.

### Cleaning and Disinfecting

#### Cleaning Compounds

The Preva Dental X-Ray System requires disinfection. The cleaning and disinfecting methods described here protect operators and patients in a manner that is safe for the equipment.

Progeny Dental recommends the use of parachlorometaxylenol-based disinfectant products, such as EnviroSystems "EcoTru Professional", or the equivalent.

#### Cleaning Methods

Between each patient, perform the following cleaning and disinfecting steps.

1. Remove gross bio-burden from the cone, handles and structure with a disposable towel moistened with water.
2. Dry the cone, handles and structure with disposable towels.
3. Wipe the cone, handles and structure with the parachlorometaxylenol-based disinfectant product following the disinfectant manufacturer's instructions.
4. Clean any remaining parachlorometaxylenol from the component with water. This additional step prevents possible product discoloration or corrosion.
5. Dry the cone, handles and structure with disposable towels.

**Caution!**

**The Preva Dental X-Ray System is not waterproof. Use only moistened, not saturated towels.**

## Checking System Functions

The following checks must be performed to complete the installation of the Preva Dental X-Ray System and as part of the recommended maintenance as indicated in the User Manual. Failure to perform these checks may result in an installation that does not comply with U.S. Radiation Performance Standards 21 CFR Subchapter J.

**CAUTION!**

**If the Preva Dental X-Ray System does not perform the functions below, advise the owner that the system is not to be used. See the Troubleshooting section of this manual or contact Progeny's Technical Support.**

<b>System Function Checklist</b>		✓
<b>Wall Mounting</b>	Ensure that the wall support is adequate and that the system is properly mounted to the wall.	
<b>Labels</b>	Ensure that all certified components bear labels that include the model and serial number, date of manufacture and a statement of certification as noted elsewhere in this manual.	
<b>Tubehead</b>	Check for oil leaks or other evidence that could indicate internal damage. Replace the Tubehead, if necessary.	
<b>Tubehead Rotation</b>	Ensure that the Tubehead maintains its position around the horizontal axis while remaining easy to rotate and position. Also check the vertical pivot of the Tubehead for easy movement while remaining in position after moving.	
<b>Suspension</b>	Check that all movements are smooth and quiet. Verify that the Tubehead is properly counterbalanced for vertical drift and that the Horizontal and Articulating Arms do not drift horizontally.	
<b>Power Switch</b>	Verify that the switch is working properly and that the Ready Indicator is illuminated when the power switch is in the ON position.	
<b>Operator Panel Controls</b>	With the power switch, located at the upper right of the Control Unit, in the ON position, verify that technique factors appear on the Operator Panel. Also, check the function of the selection buttons for Tooth Selection, Image Receptor Type and Patient Size. Pressing a selection button should cause indicator lamps to indicate the selected item.	
<b>Exposure Button</b>	Verify that the Exposure button on the Operator Panel is functioning properly. To make an exposure, press and hold the Exposure button until the Radiation Indicator is extinguished and the audible signal is no longer heard.	
<b>Exposure Indicators</b>	Make several exposures and verify that the Radiation Indicator illuminates and the audible signal is heard.	
<b>Premature Termination</b>	Select the longest exposure time possible using the up and down arrows. Initiate an exposure but release the Exposure button after a brief period of time before the timer terminates the exposure. Verify that the display indicates "Pretermination Error" and returns to normal operating mode.	
<b>Coil-cord Hand Switch Option</b>	If a coil-cord hand switch is used, inspect the switch housing and coil cord for damage or wear. Replace if evidence of damage is present.	
<b>User Information</b>	Make certain that the user of the system has received the User Manual.	

<b>New Tube Seasoning Procedure</b>	<p>X-ray tubes which sit dormant for several months can become electrically unstable. To remedy this condition it is recommended to perform a “new tube seasoning procedure”. This process will establish stable high voltage operation and, will ultimately extend the life of the tube. Repeat this procedure before returning to normal operation any time the system has been unused for more than two months.</p> <ol style="list-style-type: none"><li>1. Verify system operation.</li><li>2. Energize the system.</li><li>3. Select 60 kilovolts, 7 milliamperes and, the exposure time of one second.</li><li>4. Make five exposures at this level, observing the normal cooling time.</li><li>5. Select 65 kilovolts, 7 milliamperes and the exposure time of one second.</li><li>6. Make five exposures at this level, observing the normal cooling time.</li><li>7. Select 70 kilovolts, 6 milliamperes, and an exposure time of one second.</li><li>8. Make five exposures at this level, observing the normal cooling time.</li><li>9. Proceed with the remainder of the installation.</li></ol>	
-------------------------------------	---	--



## Solving Performance Issues

### Performance Issues

#### Light or Dark X-Ray Images

1. Adjust the selected exposure time, kilovoltage or tube current to produce an acceptable image. If necessary, reprogram the techniques factors, as explained in the System Configuration section of this manual.
2. Verify the kilovoltage and tube current during an exposure using the diagnostic mode, as explained in the System Configuration section of this manual. Alternatively, you may employ a non-invasive meter to evaluate kilovoltage and exposure time.
3. Inspect the condition of the remaining imaging chain components such as the film, chemistry and processor, or the condition of the x-ray sensor and computer.

#### No X-Ray

- If no x-ray is produced, check the following:
1. Verify that the line cord (if one is in use) is properly connected.
  2. Verify that the power switch is in the ON position.












### Obtaining Technical Support












#### Contact

PROGENY DENTAL  
1407 Barclay Blvd.  
Buffalo Grove, IL 60089  
Phone: 888-924-3800  
Fax: 847-459-5175  
techsupport@progenydental.com

## Pre-programmed Exposure Times

Using the choice of digital receptor, there are 3 choices of default exposure times shown below. The tables below show the factory default exposure settings for each combination of Tooth, Image Receptor Type, and Patient Size on the Operator Panel. These exposure settings can be modified using the System Configuration mode. See the System Configuration section for details.

8 inch (20 cm) Cone		Schick Receptor		Progeny Receptor		Dentrix	
Tooth Selection	Setting	Adult 	Child 	Adult 	Child 	Adult 	Child 
Incisor 	kV	60	60	65	65	70	70
	mA	7	7	7	7	6	6
	seconds	0.08	0.04	0.100	0.05	0.08	0.04
Bicuspid 	kV	60	60	65	65	70	70
	mA	7	7	7	7	6	6
	seconds	0.100	0.05	0.125	0.64	0.08	0.04
Bitewing 	kV	60	60	65	65	70	70
	mA	7	7	7	7	6	6
	seconds	0.200	0.100	0.200	0.100	0.125	0.064
Lower Molar 	kV	60	60	65	65	70	70
	mA	7	7	7	7	6	6
	seconds	0.100	0.05	0.125	0.064	0.08	0.04
Upper Molar 	kV	60	60	65	65	70	70
	mA	7	7	7	7	6	6
	seconds	0.125	0.064	0.160	0.08	0.100	0.050

12 inch (30 cm) Cone		Schick Receptor		Progeny Receptor		Dentrix	
Tooth Selection	Setting	Adult 	Child 	Adult 	Child 	Adult 	Child 
Incisor 	kV	60	60	60	60	70	70
	mA	7	7	7	7	6	6
	seconds	0.160	0.08	0.200	0.100	0.160	0.080
Bicuspid 	kV	60	60	60	60	70	70
	mA	7	7	7	7	6	6
	seconds	0.200	0.1	0.250	0.125	0.160	0.080
Bitewing 	kV	60	60	60	60	70	70
	mA	7	7	7	7	6	6
	seconds	0.400	0.200	0.400	0.200	0.250	0.125
Lower Molar 	kV	60	60	60	60	70	70
	mA	7	7	7	7	6	6
	seconds	0.200	0.1	0.250	0.125	0.160	0.080
Upper Molar 	kV	60	60	60	60	70	70
	mA	7	7	7	7	6	6
	seconds	0.25	0.125	0.320	0.160	0.200	0.050

# System Configuration

## System Configuration Mode


### About System Configuration Mode

The Preva Dental X-Ray System has a software-driven system configuration mode. When the Preva is in system configuration mode, you can perform the following procedures:

- Adjusting the Display
- Changing Pre-programmed Exposure Settings
- Changing the Cone Size
- Showing Current System Configuration
- Displaying Diagnostic Data

### Using System Configuration Mode

1. To enter system configuration mode, depress the Tooth Selection and Patient Size Selection buttons on the Operator Panel simultaneously for 5 seconds. The display shows the Main System Configuration menu, as shown in Figure 3, and the Ready Indicator blinks.
2. To select menu items while in system configuration mode, use the up and down arrows to highlight a menu option. Then use the right arrow button as an “Enter” button to select the highlighted option. When changing presets, the right arrow button is also used to select the technique factor.
3. After selecting a menu option, use the up and down arrows to increase or decrease values.



**Figure 3**  
Main System  
Configuration Menu

## Adjusting the Display

### Adjusting Contrast

### Reversing the Image

The Preva Dental X-Ray System allows the operator to adjust the display image.

1. From the system configuration main menu, shown in Figure 3, select ADJUST DISPLAY. You will see the Display Options menu shown in Figure 4.
  2. Selecting EXIT returns the display to the Main System Configuration menu shown in Figure 3.
1. Select ADJUST CONTRAST from the menu. You will see the Progeny logo.
  2. Use the up and down arrows to increase or decrease the contrast between the menu text and the display background.
  3. Press the right arrow to save your settings.
1. Select REVERSE IMAGE from the menu. The text and display background colors will be swapped.
  2. Press the right arrow to save your settings.



Figure 4  
Display Options  
Menu

## Changing Pre-programmed Exposure Settings

The Preva Dental X-Ray System allows the operator to increase or decrease image density for all presets for a receptor simultaneously or to change each of the technique factors for a preset individually. You can also restore factory default settings. For charts of the factory default settings, refer to Factory Default Exposure Settings later in this manual.

**Note:** If the 12 inch cone is going to be used, configure the Preva for use with the 12 inch cone before changing pre-programmed exposure settings.

Configuring the Preva for use with the 12 inch cone will reset exposure settings to the default settings used with the 12 inch cone.

**Note:** Before changing presets, use the table below to write down the presets you are programming.

### Displaying the Preset Options Menu

1. From the Main System Configuration menu, shown in Figure 3, select CHANGE PRESETS. You will see the Preset Options menu shown in Figure 5.
2. Selecting EXIT returns the display to the Main System Configuration menu shown in Figure 3.




Figure 5  
Preset Options  
Menu

### Changing All Receptor Settings Globally

1. Select ALTER DENSITIES from the Preset Options menu. The first Image Receptor Type illuminates. The display shows the selected Image Receptor Type and current density.
2. Using the Image Receptor Type button, select the image receptor to adjust.
3. Use the up and down arrow buttons to specify a percentage by which densities will be increased or decreased for the selected receptor. Densities can be increased in steps of 25% and decreased in steps of 20%.
4. Press the right arrow to save your settings.

## Preprogramming to Digital Sensors

1. Energize the system.
2. Press and hold the "Tooth Selector" switch and, the "Patient Size Selector" switch for five full seconds.
3. When this is completed, the "Menu Options" screen will be displayed.
4. Press the "Down" switch once to reach the "Change Presets" location.
5. Press the "Enter" key.
6. This causes the display to reach the "Preset Options" screen.
7. Press the "Down" switch twice to reach the "Select Receptor" location.
8. Press the "Enter" key.
9. This causes the display to reach the "Select Receptor" menu.
10. Press the "Down" (or "Up") switch until the desired sensor or phosphor plate is highlighted.
11. Press the "Enter" key.
12. This causes the display to reach a verification screen. Select "Yes" or "No" by pressing the "Up" or "Down" switches.
13. This selection returns the display to the "Preset Options" menu.
14. Press the "Down" switch four times to reach the "Exit" location.
15. Press the "Enter" key.
16. This causes the display to return to the "Menu Options" screen.
17. Press the "Down" switch three times to reach the "Exit" location.
18. Press the "Enter" key.
19. A brief indication of "Saving Settings" will be observed, then the system will return to the normal operational mode.

## Changing Presets Individually

1. Select EDIT PRESETS from the Preset Options menu. The display notifies you that you are entering Edit Preset Mode, and Tooth Size, Image Receptor Type, and Patient Size are illuminated.
2. Use the Tooth Selection, Image Receptor Type, and Patient Size Selection buttons to select the preset to change. The display shows the current values for the preset.
3. Use the right arrow button to display the technique factor to change.
4. Use the up and down arrow buttons to set the value for the selected technique factor and preset.
5. Repeat steps 2-4 to change additional presets.
6. When you have completed all changes, depress the Tooth Selection and Patient Size Selection buttons simultaneously for 5 seconds to record the change.

## Select Receptor















This menu item is reserved for future use.

## Recall Presets

1. To return all presets to factory defaults, select RECALL PRESETS from the Preset Options menu. The menu will ask you to confirm your choice.
2. Select YES using the up arrow button and return all presets to factory default settings. Selecting YES will erase any custom presets that have been set up.
3. Select NO using the down arrow button and retain current presets.

## Record Your Pre-programmed Exposure Settings

If the pre-programmed exposure settings do not produce the density desired, adjust the settings using System Configuration mode. Record your settings in the table below.


8 inch (20 cm) Cone		Digital Receptor 		D-speed Film 		E/F Speed Film 	
Tooth Selection	Setting	Adult 	Child 	Adult 	Child 	Adult 	Child 
Incisor 	kV						
	mA						
	seconds						
Bicuspid 	kV						
	mA						
	seconds						
Bitewing 	kV						
	mA						
	seconds						
Lower Molar 	kV						
	mA						
	seconds						
Upper Molar 	kV						
	mA						
	seconds						



## Showing Current System Configuration

The Preva Dental X-Ray System displays the current system configuration. This display is informational only.

1. From the Main System Configuration menu, shown in Figure 3, select CONFIGURE UNIT. You will see the Configuration menu shown in Figure 6.
2. Select SHOW CONFIG. The display will show:
  - Current software version
  - Cone size
  - Diagnostic mode on or off
3. Press any button on the Operator Panel to return to the Configuration menu.




**Figure 6**  
**Configuration Menu**

## Changing the Cone Size

### Using a 12 inch Cone

Selecting SET CONFIG. from the Configuration menu, shown in Figure 6, displays the Set Configuration menu, Figure 7, with options to change the cone size. The Preva Dental X-Ray System is factory set for use with the standard supplied 8 inch (20 cm) Cone. The 12 inch (30 cm) Cone (30-A2033) is available. Using the longer Cone requires longer exposure times, which the Preva automatically selects when you change the Cone size in the Set Configuration menu.

1. From the Main System Configuration menu, shown in Figure 3, select CONFIGURE UNIT. You will see the Configuration menu shown in Figure 6.
2. Select SET CONFIG. You will see the Set Configuration menu, shown in Figure 7.
3. From the Set Configuration menu, use the up and down arrows to highlight 12" CONE SIZE.
4. Press the right arrow button to select the 12" CONE. The display warns you that selecting the 12 inch Cone will override custom presets with the default factory settings for the 12 inch Cone.
5. Using the up arrow, select YES to install presets for the 12 inch Cone.



**Figure 7**  
**Set Configuration**  
**Menu**

## Diagnostic Mode

### About Diagnostic Mode

The Preva Dental X-Ray System has a diagnostic mode in which you can display a summary of maintenance data or display feedback values after each exposure.

### Showing the Maintenance Summary

1. From the Main System Configuration menu, shown in Figure 3, select CONFIGURE UNIT. You will see the Configuration menu shown in Figure 6.
2. Select SET CONFIG. You will see the Set Configuration menu, shown in Figure 7.
3. To display a summary of maintenance data, highlight select SHOW MAINT. The following maintenance data are displayed:
  - Total KJ (kilojoules—total system heat on x-ray tube)
  - Exposure Count
  - Reboots (power up cycles)
  - OT Counts (over-threshold counts)
4. Press any button on the Operator Panel to return to the Configuration menu.

### Showing Feedback Values After an Exposure

If you take an x-ray while in diagnostic mode, the display shows feedback values for that exposure. Until you exit diagnostic mode, the display will continue to show feedback values after each exposure.

1. From the Main System Configuration menu, shown in Figure 3, select CONFIGURE UNIT. You will see the Configuration menu shown in Figure 6.
2. Select SET CONFIG. You will see the Set Configuration menu, shown in Figure 7.
3. From the Set Configuration menu, use the up and down arrows to highlight DIAG MODE ON. Press the right arrow button to turn on diagnostic mode.
4. Exit System Configuration mode by highlighting and selecting EXIT in the Configuration and Main menus.
5. Make an exposure. The display will show the following feedback values:
  - kV
  - mA
  - Filament current
6. Press any button on the Operator Panel to clear the feedback values from the display.
7. To exit diagnostic mode, depress the Tooth Selection and Patient Size Selection buttons simultaneously for 5 seconds to display the Main System Configuration menu. From the Main System Configuration menu, highlight and select CONFIGURE UNIT. Then highlight and select SET CONFIG. On the Set Configuration menu, highlight and select DIAG MODE OFF.

## Specifications

### Preva Dental X-Ray System

The following specifications contain information required to be provided to the user per Federal Regulation 21 CFR.

<b>Line Voltage</b>	110-230 VAC +/- 10% 50/60 Hz
<b>Max Momentary Rating</b>	10 amps
<b>Long Term Rating</b>	1 amp
<b>Maximum Rated Tube Potential</b>	70 KVp
<b>kVp Accuracy</b>	+/- 5% selectable
<b>Tube Current</b>	4-7 mA +/- 1 mA
<b>Exposure Time</b>	10 ms through 2.00 seconds
<b>Timer Accuracy</b>	5% +/- 1 ms
<b>Source to Skin Distance</b>	8 inch (20 cm) 12 inch (30 cm)
<b>Minimum Half Value Layer</b>	1.7 mm Aluminum equivalent at 70 kVp
<b>Minimum Inherent Filtration</b>	2 mm Aluminum equivalent @ 70 kVp
<b>Focal Spot</b>	0.4 mm (IEC 336)
<b>Automatic Cooling Time</b>	15 times the exposure time wait before the next exposure can begin
<b>Leakage Technique Factors</b>	1.5 mA at 70 kVp
<b>Target Angle</b>	12.5 degrees
<b>Operating Temperature</b>	+50 F/+95 F (+10 C/+35 C)
<b>Storage Temperature</b>	-31F/+150 F (-35 C/+66 C)
<b>Maximum Altitude</b>	12,000 ft
<b>Cone Focal Length</b>	8 inch(20 cm) 12 inch(30 cm)
<b>Diameter of X-Ray</b>	2.7 inch (6.9 cm) at the end of the Cone





Figure  
Tube Rating Charts

