

# FD-PHANTOM ACTA Technical Manual and Spare Parts Catalogue

# **Table of Contents**

1	INT	TRODUCTION	3
	1.1	Symbols and Abbreviations	3
	1.2	Warnings and Safety	
	1.3	Technical Specification	6
	1.4	Disposal of waste	
	1.5	Main Parts	
	1.6	Dimensions	11
2	QU	JICK REFERENCE GUIDE	12
3	MA	AINTENANCE AND CLEANING	14
4	UN	IIT TEST PROTOCOL	15
	4.1	Pre-Test Inspection	
	4.2	Electrical	
	4.3	Main Air and Suction Supply	
	4.4	Water Supply	
	4.5	Functional Test	
5	INS	STRUMENT SERVICE	19
_	5.1	Syringe	
	5.2	Curing Light	
	5.3	Scaler	
	5.4	Micromotor	
	5.5	Turbine	
	5.6	Operating Lamp	
	5.7	Water Bottle	
	5.8	Foot Control	
6	SEE	RVICE INSTRUCTIONS	22
•	6.1	Instrument Selection and Whip Arms	
	6.2	Check the Instrument Hoses for Wear	
	6.3	Replacing an Instrument Hose	
	6.4	Check Instrument LED Lights	
	6.5	Instrument O-rings	
	6.6	Check Instrument Tips for Wear	
	6.7	Clean the Syringe Tip	
	6.8	Lubricate the Contra-Angles	25
	6.9	Lubricate Turbine (Generic Instructions Only)	26
	6.10	After Lubrication (Generic Instructions Only)	
	6.11	Check Mechanical Functions of Instrument Bridge	
	6.12	Check for Leaks Inside	
	6.13	Change Arm Brake Rubber of Horizontal Arm	
	6.14	Calibrate the Horizontal Arm Motion	
	6.15	Lubricate Arm Joints and Pneumatics	
	6.16	Horizontal Arm Stop Adjustment	
	6.17 6.18	Main Air Pressure Valve Filter	
	6.19	Replace the Selection Valve Membrane	
_		·	
7		B DIAGRAMS	
	7.1	Q900 Power Supply PCB	
	7.2	Q905 Phantom LED PCB	
	7.3	Q915 Phantom Instrument Control Board	
8	WI	RING AND PNEUMATIC DIAGRAMS	38
	8.1	Pneumatic Diagram	38
	8.2	Wiring Diagram	40
9	SPA	ARE AND CONSUMABLE PARTS	42
	9.1	Consumable Parts	
	9.2	Spare Parts	
	9.3	Assembly Drawings with Bills of Materials	

#### 1 INTRODUCTION

This manual describes how to maintain and repair Finndent Phantom dental training units. The Phantom unit consists of a table surface, instrument bridge, suction arm, Phantom torso and operating light. Installation may only be done under Finndent supervision.

Maintenance and repair may only be done by a certified Finndent Technician.

This manual is valid for the instrument bridge Q915 PCB firmware version "FD915\_v1\_00\_02".

#### **ACCOMPANYING DOCUMENT - USER MANUAL**

For students, dental professionals and teaching staff. Contains use and regular care instructions for Phantom units. Also contains basic information on parts, trouble shooting and parts replacement.

#### ACCOMPANYING DOCUMENT - TECHNICAL AND SPARE PARTS MANUAL

For Finndent authorized dealers. Contains adjustment, repair, calibration, and replacement part information. Wiring and pneumatic diagrams are included.

# 1.1 Symbols and Abbreviations

V Voltage

A Ampere

W Watt

Hz Hertz

Alternating current

FI Approved by Fimko

| / O On / Off (IEC 60417-5007/5008)

IPX1 Protected against dripping water (IEC 60529)

PCB Printed Circuit Board

RPM Revolutions per minute

T5AL250V Fuse description: T = Slow blow / F = Fast blow

5A = number of amperes, for example, 5 amp

L = low breaking capacity glass fuse

250V = rated for 250 volts 35A BC = 35 A breaking capacity

/ L / N Protective Earth (IEC 60417-5019) / Line In / Neutral conductor connection.



Warning: Dangerous voltage (IEC60878).



Do not dispose in normal household waste (Directive 2002/96/EC)/
Date of manufacture (IEC 60529).



Type B equipment (IEC 60878).



Consult operators manual (ISO 7010).

1 INTRODUCTION 1.2 Warnings and Safety

#### 1.2 Warnings and Safety

#### **ELECTRICAL WARNINGS**



This unit contains live mains voltage parts. Before attempting maintenance do all of the following:

- Turn off all units using the power switches on the unit towers.
- Turn OFF the table power supply switches in the table's end supply cabinet.
- Disconnect the live mains voltage by turning off the building circuit breaker that supplies the power. The building circuit breaker must be equipped with a locking device to prevent accidentally turning the power on.
- If the building circuit breaker is located outside the clinic, Finndent also recommends leaving a note of, "DO NOT SWITCH ON, WORKING ON LINE".
- Test that the unit is electrically isolated using an electrical testing device.

Turning the unit off from the power switches is not a guarantee that the mains power is off in all internal parts.

Always use the correct sized fuses.

Do not modify this equipment.

#### CAUTION - CIRCUIT BOARD HANDLING

Circuit boards may be damaged by static electricity. Finndent suggests the following handling guidelines:

It is recommended to wear wrist grounding straps and store the PCB in an anti-static bag.

If not using wrist straps, touch an exposed metal grounded part of the unit or chair before touching the PCB.

Always put PCBs in anti-static bags for handling and shipping.

PCBs only qualify for warranty repair of shipped in an anti-static bag.

Never place PCBs on a conductive surface or give it to a person without first touching the surface or person.

#### **CAUTION – AIR AND WATER SUPPLY**

A shut-off valve is mounted on each unit. Always close the supply valves and relieve system pressure before doing maintenance.

#### CAUTION - APPROVED ACCESSORIES

Only Finndent approved parts and accessories may be used with this equipment. The use of other accessories is not covered under warranty and may interfere with the electrical or mechanical properties of the equipment.

# 1.3 Technical Specification

#### **Unit Manufacturer**

Finndent Oy, Niittylänpolku 16, 00620 Helsinki, Finland.

Telephone: +358 20 743 5115 E-mail: info@finndent.com

#### Device Class According to Medical Devices Directive 93/42/EEC

Class I equipment. This simulation unit is not for clinical use, or use with a patient.

#### **Electrical Protection Class / Grade**

Class I Type B

#### **Operational Safety Grade**

Not suitable for use in an oxygen rich environment in which the concentration of oxygen is greater than 25%, including mixtures of inflammable anaesthesia, oxygen or laughing gas and air.

**Protection Against Fluids** (Standard EN 60529 +A1 Degrees of protection provided by enclosures)

IPXI for foot control

**Backflow Prevention** (Standard EN 1717 Protection against pollution of potable water in water installations and general requirements of devices to prevent pollution by backflow)

The unit is isolated from the mains water supply and takes water from a bottle system.

**Allowable Temperature, Humidity and Pressure Ranges** 

	Temperature	Relative Humidity	Atmospheric
		(no condensation)	Pressure
Operating Conditions:	15 – 40ºC	20 - 75%	80 kPa-110 kPa
Operating Conditions.	13 – 40-0	20-73/6	(0.8bar-1.1bar)
Storage Conditions:	0 - 50ºC	20 - 75% + Keep Dry	70 kPa-110 kPa
Storage Conditions.		20 - 75% + Reep Dry	(0.7bar–1.1bar)
Transportation Conditions:	-20 - +60ºC	5 - 95% + Keep Dry	70 kPa-110 kPa
Transportation Conditions.	-20 - <del>+0</del> 0=C	3 - 95% + Keep Dry	(0.7bar–1.1bar)

If the unit has been stored outside accepted operating conditions for more than a few hours, allow the unit to reach operating conditions before attempting to install or connect any part.

#### **Electrical Features**

Supply Voltage / Frequency:	230V / 50Hz or 60Hz	
Unit power consumption:	500W (approx. maximum)	
Volt-amperes of unit:	460VA	
F1 (PCB Q900 mains fuse):	T6,3AL250V (T6,3A slow blow fuse, Schurter	
	0034.3125 / Littlefuse 021806.3MXP).	
	Breaking capacity (BC) 35A	
F2 (PCB Q900 mains fuse):	T1AL250V (T1A slow blow fuse, Schurter	
	0034.3117 / Littlefuse 0218001.MXP).	
	Breaking capacity (BC) 35A	
F3 (PCB Q900 mains	T4A250V (T4A slow blow board mount radial	
microfuse):	microfuse, Schurter 0034.6621 / Littlefuse	
	37214000411). Breaking capacity (BC) 40A	
F4 (PCB Q900 mains	T1,25A250V (T1,25A slow blow board mount	
microfuse):	radial microfuse, Schurter 0034.6616 / Littlefuse	
	37211250411). Breaking capacity (BC) 40A	
F1 (PCB Q915 Instrument	F5A125V (F5A fast blow board mount square end	
Bridge Microfuse):	block microfuse, Schurter 3-101-056).	
	Breaking capacity (BC) 50A	
F2 (PCB Q915 Instrument	F5A125V (F5A fast blow board mount square end	
Bridge Microfuse):	block microfuse, Schurter 3-101-056).	
	Breaking capacity (BC) 50A	

The power supply cord is non-detachable and may only be replaced by a technician. The equipment is permanently installed.

#### **Mode of Operation**

Unit: Continuous operation with intermittent loading depending on feature use. Typical intermittent instrument use of 10 minutes ON / 5 minutes OFF.

#### Torso Lift Motor:

Non-Continuous operation.

Duty cycle time in normal use: 10% of max. 2min motor ON/18min OFF.

Speed of motor: 33mm/sec (0.033m/sec) maximum speed.

Length of stroke: 250mm (0.25m).

Time to travel in one direction: 7.5 seconds.

#### **Software**

This unit contains no software and has no programmable parts.

#### **Bottled Water Connection**

Pressure range: 2,5 bar (250 kPa) maximum

Flow Rate: ≤ 4 l/minute maximum consumption

Quality of water: particle size < 10µm

Bottle Connector 28PCO

Size: 1.5L maximum

Water lines are always GREEN or CLEAR/WHITE. Water spray adjustment knob has a ring on the base.



#### **Compressed Air Supply**

Pressure on unit 3 bar maximum (300kPa)

valve:

Flow Rate on unit: ≥ 55 l/minute maximum consumption

Quality: dry, oil-free and hygienically clean (medical grade)

Hose Connector 10mm

Hygienic air lines are always blue.

Air spray adjustment knob is smooth on the base.



The turbine air supply and return lines are red in the instrument bridge and black in the unit.

#### Cooling / Power Fluid

The hoses for instrument air or water cooling or for air turbine rotation are white/clear.

#### **Suction Connections**

Vacuum range per unit:	150 – 170 mbar (15-17 kPa)
Required Flow Rate per unit:	≥ 500-800 I/minute consumption
Universal cannula on large suction	250 – 350 l/min
hose:	
Surgical suction cannula on small	50 – 80 l/min
suction hose:	
Saliva cannula on small suction hose:	30-60 l/min
Mains Suction Supply to Table:	50mm
2 x Central Table Suction Hoses:	32mm
Large Suction Hose	17.5mm
Small Suction Hose	11.5mm

#### **Load Capacity**

Foot Control Tray:	5kg
Instrument Tray:	1kg

#### Colour

White painted parts and plastics are RAL 9003 Signal White. Grey plastic parts are RAL 9006.

#### **Installation and Disassembly**

The dental unit must be installed by a Finndent technician. Contact Finndent if the unit requires relocation or disassembly other than listed in this manual.

1 INTRODUCTION 1.4 Disposal of waste

# 1.4 Disposal of waste



Finndent products have been designed and manufactured to be as safe as possible. Any waste materials must be recycled or disposed of in an environmentally friendly manner according to national regulations.

Hazardous materials requiring special waste collection must be disposed of in accordance with national waste and environmental regulations.

When handling waste materials, all precautions must be taken to reduce the associated risks.

When the unit reaches the end of its service life, return it to your dealer for disposal.


Contact your local dealer or national governmental agency for more specific information.

Description	Main Materials	Recycle	Waste Disposal	Special Waste Collection
	Aluminium	Х	•	
Unit frames	Galvanized steel	Х		
	Particle board table top	Х		
	Copper	Х		
Electronics, PCBs,	Steel	Х		
motors	PCBs, Motors, other			X - EC directive
	rebs, Motors, other			2002/96
Rubber			Х	
Glass		Х		
Cleaning products	Chemicals			Х
Oil collector	Chemicals			Х
	Wood	Х		
	Cardboard	Х		
Packaging	Paper	Х		
	Foam and SealedAir padding		Χ	
	Other plastics	Х		
Other accessories	See OEM Product			
Other accessories	Documentation	-	_	-

Table: Guideline for waste disposal

1 INTRODUCTION 1.5 Main Parts

#### 1.5 Main Parts




- 1. Dental lamp.
- 2. Tower.
- 3. ON / OFF switch.
- 4. Power socket.
- 5. Whip arms.
- 6. Instrument bridge.
- 7. Instrument tray.
- 8. Phantom head.
- 9. Torso.
- 10. Foot pedal.
- 11. Suction supply hoses.
- 12. Foot control plug.
- 13. Suction arm.
- 14. Water bottle for all water supply.
- 15. Suction holder (manifold).
- 16. Torso lift control (up-down).
- 17. Power socket.
- 18. Supply access for maintenance (air, water, suction, power).

1 INTRODUCTION 1.6 Dimensions

# 1.6 Dimensions

Phantom dimensions depend on the number of units and lamp arm length:




#### **Arm Motion Dimensions**

Radius of instrument bridge: ~1050mm adjustable. Radius of Faro dental lamp: ~1600mm maximum.

Radius of suction arm: ~650mm.


Arm motion may vary depending on the motion blocks used.

#### **Torso and Table Heights**




# 2 QUICK REFERENCE GUIDE

See User Manual for complete instructions.




#### **Instrument use and Basic Functions**




**Table: Quick Reference instructions** 

# Air and water spray adjustment



# Adjusting dental lamp intensity



Turn knobs left to increase spray and right to decrease spray.


The AIR spray knob has a smooth knob.




The WATER spray has a ring on the base.




#### Raise or lower torso




#### **Position torso**



#### **Position the Phantom Head**









#### Position the headrest

- 1. Push the support bar inward.
- 2. At the same time, rotate the bar.
- 3. Release the bar to re-lock.
- 4. Tilt the head rest up and down.





# 3 MAINTENANCE AND CLEANING

See User Manual for complete preventative and maintenance instructions.

For instruments, lamps and suction manifolds, follow the OEM manufacturer's cleaning instructions.

The following are the basic maintenance required by Finndent. The Phantom Clinic Technicians may add to this list for clinic-specific maintenance and cleaning.

The dental unit must be serviced by a Finndent approved technician once per year. **Annual service and testing is mandatory.** 

		Recommended Frequency			
Part	Maintenance Description		Weekly	Monthly	Yearly
Instrument bridge handles	Remove handles (P/N: 6203230) and clean any dirt from inside so the touch sensing stays sensitive.		Х		
Syringe	Remove metal sleeve and clean inside.		Χ		
Scaler Micromotor Turbine	Measure tips to check for wear regularly and replace as needed. Follow the manufacturer's instructions!		Х		
Micromotor	Clean and lubricate motors regularly.		Х		
Turbine Syringe Micromotor Turbine	Follow the manufacturer's instructions!  Check and replace O-rings regularly as needed. Follow the manufacturer's instructions!			X	Х
Turbine	Check the turbine oil return cup regularly. When full, empty and replace the sponge.		Х	Х	
Faro light handles	Remove handles and thoroughly clean.		Х		
Dental light parabola	Perflex spray or similar to polish surface as needed. Follow the manufacturer's instructions!	х	Х		
Unit Exterior	Wipe, vacuum and clean the exterior of the unit.	Х			
Phantom Head	Clean inside the head. Replace face mask as needed.	Х	Х		
Dürr Comfort suction manifold filter	Check and rinse coarse suction filter (P/N: 8900511) regularly. Replace filter as needed.	х			
Suction handpieces	Remove, disassemble and clean dirt from the inside so the volume control is clean.		Х		
Suction hoses	Replace worn hoses.				Х
Suction hoses	Clean and rinse hoses to prevent bio and dirt build-up.		Х		
Suction hoses	Replace O-rings regularly, if suction handpieces are leaking.			Х	Х
Suction selection valve	Replace membrane (P/N: 5500310) as needed.				Х
Water bottle	Replace water bottle when cracked or dirty.		X	Х	
Water line cleaning	Before the weekend and before long holidays fill the water lines with BioFinndent.		Х	Х	
Unit foot control	Check cable for wear and breaks in the grey insulation. Replace as needed.			Х	Х
Main air pressure regulator filter	Replace air filter.				х

**Table: Maintenance and cleaning schedule** 

4 UNIT TEST PROTOCOL 4.2 Electrical

#### 4 UNIT TEST PROTOCOL

If a problem is suspected during inspection:

- Turn the mains power supply off at the mains fuse panel.
- Turn the clinic air supply off.
- Turn the clinic suction supply off.

It is now safe to investigate and repair.

#### 4.1 Pre-Test Inspection

Check that all parts are installed on the unit (instruments, water bottle, and torso).

Visually check that the green/yellow grounding wires are attached to the frame.

Visually check that the cables and hoses are in place.

Visually check for leaks in the unit or on the floor.

Listen for the sound of leaking air.

Use a current/voltage detector to check that the power sockets have no voltage when the unit is OFF.

Turn the unit ON only if all inspections are OK.

Visually check that the power switch light is ON.

#### 4.2 Electrical

Electrical safety tests are required for Class I equipment. Follow the guidance in IEC 62353:2014 ed2.0b Medical electrical equipment – Recurrent test and test after repair of medical electrical equipment.

#### Protective earth resistance test (IEC 62353:2014 Ch. 5.3.2)

The protective earth resistance shall not exceed  $300m\Omega$ .

Measure the resistance between the protective earth connector of the mains plug and the protectively earthed instrument bridge and torso backrest.

#### Touch current test (IEC 62353:2014 Ch. 5.3.4.4)

The touch current for accessible conductive parts shall not exceed 100µA.

Measure the touch current at the mains plug and at the instrument bridge and the torso backrest.

#### Power socket test

Use a current/voltage detector to check that the power sockets have 230V when the unit is ON.

4 UNIT TEST PROTOCOL 4.5 Functional Test

# 4.3 Main Air and Suction Supply

When the unit and air supply are ON, the main air pressure regulator shall read between 6-8 bar.

Measure the suction of each hose:

- Universal cannula on large suction hose: 250 350 l/min.
- Surgical suction cannula on small suction hose: 50 80 l/min.
- Saliva cannula on small suction hose: 30-60 l/min.

# 4.4 Water Supply

The water bottle is pressurized with air when the unit is on.

Check the air pressure by gently squeezing the bottle. It should feel full / hard.

The water bottle pressure should be 2 bar.

#### 4.5 Functional Test

A full functional test should be performed after installation and after annual maintenance.

A partial functional test can be performed after minor repairs to an individual system.

4 UNIT TEST PROTOCOL 4.5 Functional Test

FINNDENT OY	NNDENT OY FUNCTIONAL INSTALLATION TEST		PHANTOM UNIT
Table number:		Test date: Tested	
number:		by:	
Serial number:		Page:	1/2
	Test Description	Verified	Observation, reading or description of repair done?
V	isual inspection of electrical systems.		
Visual inspect	ion of air, water and suction systems.		
Protective eart	th resistance test (300 mΩ maximum)		
-	Touch current test (100μA maximum)		
Pow	er sockets tested with unit off and on		
Main ai	r pressure regulator reading (6-8 bar)		
Large suction hose v	vith universal cannula (250-350I/min)		
Surgical suction	cannula on small suction hose (50-80 l/min)		
Saliva cannul	a on small suction hose (30-60 I/min)		
FARO LAMP TEST			
Le	Turn the Faro lamp on. Test the half-intensity function. Test the dimming function. eave the lamp on during the unit test.		
SUCTION ARM FUNC	TION		
Chec	rm left and right. Put the suction arm in a usable position. k that the water bottle is pressurized. cion hose and check for suction. Put it back.		
Lift the sma	Ill suction hose and check for suction.  Put it back.		
TORSO FUNCTION			
the Check the position	I suction arm fully up and down using motor (switch is on the suction arm).  Move the torso to a usable height.  Lower the torso to 30°.  brake by pushing gently on the chest.		
PHANTOM FUNCTIO	r the torso to the horizontal position.  N		
	Open the Phantom head. ew the neck lock and move the head. Re-lock.		
INSTRUMENT BRIDG			
	ment bridge by holding the handle to unlock the arm air brake. brake locks the arm motion when the handle is released.		
	se and place in the third hand holder. that the third hand is firmly attached.		

**Table : Functional test** 

4 UNIT TEST PROTOCOL 4.5 Functional Test

FINNDENT OY FUNCTIONAL INST		TALLATION TEST	PHANTOM UNIT
Table number:		Test date: Tested	
rosition number.		by:	
Serial number:		Page:	2/2
Test Desc	cription	Verified	Observation, reading or description of repair done?
SYRINGE FUNCTION			
Spray mixed air-	Lift the 3-function syringe. r only in to the suction hose. Spray air only. water in to the suction hose. ift the micromotor to select.		
	nge and micromotor sprays.		
CURING LIGHT FUNCTION	<u> </u>		
DO NOT Ru	he LED curing light to select. LOOK DIRECTLY AT THE LED. In the A 10 second fast cycle.		
SCALER FUNCTION	ne 20 second slow-rise cycle.		
Slide the pedal left and	Lift the scaler to select. I pedal to turn the scaler on. I right to change the power. press the foot control spray ay. The display LED for spray should turn ON. Test spray.		
Turn the spra	y off and replace the scaler.		
MICROMOTOR FUNCTION			
	ift the micromotor to select. edal to turn the micromotor on.		
With the micromotor sel- spray button to activate mixe spray With the micrmotor selected on the foot control to reven	should turn ON. Test spray.		
Turn the spray off. Turn the re	everse direction off. Replace the micromotor.		
TURBINE FUNCTION			
Slide the pedal left an With the turbine selected, button to activate mixed spr	Lift the turbine to select. Deedal to turn the turbine on. I dright to change the speed. Press the foot control spray ay. The display LED for spray should turn ON. Test spray. Tray off. Replace the turbine.		
WATER LINE RINSE PROGRAM			
Replac Put the instruments in a Press the foot control o	The spray LED should flash.  e the curing light handpiece.  bucket to collect the spray.  entre button to start rinsing.  by pressing the foot control  again.		
Or, allow the program to fin	_		

Table : Functional test

5 INSTRUMENT SERVICE 5.3 Scaler

#### 5 INSTRUMENT SERVICE

Follow the OEM manufacturer's instructions for service.

The supply requirements for air, water, suction and power are listed here.

The expected ranges for instrument outputs are given as generally acceptable values. For OEM specific information, see the OEM manufacturer's instructions.

#### 5.1 Syringe



Inspect the tip and handpiece for calcium and biofilm build up.

Check and clean tip for blockages.

Check the air and water buttons for sensitivity (press harder to open valve and have more flow)

Check hose and handpiece for leaks.

Check the air and water inputs:

Max- water input pressure: 2.5 bars Max. air input pressure: 4.5 bars

Air flow rate: 10 NI/min Water flow rate: 110 cc/min

# 5.2 Curing Light



Measure the curing light intensity level annually with a radiometer.

Light radiation: 440 - 480 nm

Check voltage supply to handpiece: 24 V AC ± 10 % or 33 V DC ± 10 %

#### 5.3 Scaler



Check voltage to scaler: 24 VDC.

Control voltage: 5V

Output frequency: 28 kHz minimum.

Water flow at the end of the handpiece: 5 – 40 mL/min

Check and replace O-rings twice yearly.

Faulty o-rings can lead to leaks

Check tip threads on tip and handpiece – damaged threads can ruin any handpiece they are used on!

Check that tips are OEM parts from the same manufacturer as the scaler.

Non-OEM parts can damage the tip and handpiece, and are not covered under Finndent or OEM Manufacturer warranties.

5 INSTRUMENT SERVICE 5.7 Water Bottle

#### 5.4 Micromotor



Check voltage to micromotor: 36 VDC maximum

Check voltage to LED: 2.5-4 VDC Check LED intensity: 15-40lux Check cooling spray: 6-8 NI/min

Check air spray supply: 0-3 bar, adjustable Check water spray supply: 0-3 bar, adjustable

Check and replace O-rings yearly.

Check rotations per minute with a digital tachometer.

Run the motor at low speed and high speed.

Check for vibration or noise (subjective).

If bearing or mechanical wear is detected, send to Finndent for service.

#### 5.5 Turbine



Spray water pressure: 0-3 bar, adjustable Spray air pressure: 0-3 bar, adjustable Turbine air drive pressure: 2.5 - 3 bar

Return air pressure: 0.25bar Check LED voltage: 3V

Check and replace O-rings yearly.

Lubricate.

Check and replace turbine oil return cup as needed. Check rotations per minute with a digital tachometer.

#### 5.6 Operating Lamp



Gently wipe clean using a lint-free cloth and an alcohol free cleaner.

Check the light intensity: from 3.000 lux to 50.000lux at 700 mm of distance.

Check the light focus: 205 x 100 mm focal point at 700 mm of distance.

#### 5.7 Water Bottle

Check water bottle pressure: 2,5 bar (250kPa) maximum Check the bottle out flow rate: ≤ 4 l/minute maximum 5 INSTRUMENT SERVICE 5.8 Foot Control

# 5.8 Foot Control




Check cable for wear.

Check function and responsiveness.

If the internal sliding plates seem hard to slide, wipe them with paper towel and lubricate them with lithium grease.

Potentiometer plate



Foot control cables are consumable goods and should be replaced if any wear is visible.



A step-by-step photo tutorial on how to re-build a Finndent foot control (repair kit 9801932) is available as a supplemental instruction: 2014.02.27 Technical Bulletin 8000 Foot Control V2.

6 SERVICE INSTRUCTIONS 5.8 Foot Control

#### 6 SERVICE INSTRUCTIONS

The routine service procedures will be explained in this section (where an explanation is needed). For other service instructions and non-routine repairs contact Finndent.

Some procedures may be covered on the Finndent website - technical bulletins.



Always turn the Phantom unit power off before attempting any electrical service! Check that the unit and chair are OFF with a multimeter before attempting any electrical service!

In case of disconnecting the live mains voltage by turning off the building circuit breaker, a locking circuit breaker must be used. Finndent recommends also leaving a note, "DO NOT SWITCH ON, WORKING ON LINE' to prevent unexpected switching power on.



This technical manual include the annual and regular service instructions and details only. If more detail is required to repair an electrical component, please contact Finndent Support (support@finndent.com).



#### RELIEVE THE PRESSURE IN THE WATER LINES BEFORE STARTING SERVICE!

Close the main water valve.

Use the syringe to spray out the remaining water in the system. The system is now de-pressurized.



# TURN OFF THE BUILDING POWER CIRCUIT BREAKER BEFORE WORKING ON ELECTRONICS, WIRING, OPENING THE TRANSFORMER BOX OR WORKING ON MOVING PARTS!

Turn off the building circuit breaker to disconnect the live mains voltage from all components. Also make sure the Phantom units main power switches are off. The building circuit breaker should be equipped with a locking device to prevent accidentally turning the power on. In case of disconnecting the live mains voltage by turning off the building circuit breaker, use safety notes 'DO NOT SWITCH ON, WORKING ON LINE' of unexpected switching power on.



#### ASK BEFORE TURNING THE UNIT POWER ON DURING SERVICE!

Before activating the unit, make sure all other people are clear and ask if it is safe to turn the power on. There is no physical lock-out on the Phantom unit power switch.



#### READ ALL THIRD PARTY INSTRUCTION MANUALS FOR SPECIFIC SERVICE INSTRUCTIONS!

Finndent may provide general tips for servicing accessories, however, these are not the official service instructions and should only be used as reference.



# STEP BY STEP INSTRUCTIONS MAY BE AVAILABLE FROM THE FINNDENT DOCUMENTS WEBSITE!

This technical manual assumes that service personnel have been trained by Finndent, are capable of all basic technical skills and can make direct replacements of all parts. For more detailed and step-by-step photographs of procedures, check the Finndent documents page.

#### 6.1 Instrument Selection and Whip Arms

Normal grey whip arm P/N: 9801715 Soft whip arm for scaler P/N: 9801714

For whip arms, pull instruments forward.

If the instrument selection is not working properly, clean around the holder and selection switch and try again.

If cleaning does not work, the selection microswitch P/N: 3200138 may need replacing.

To adjust the tension in the whip arm (make it lie back or bend forward), twist the brass nut at the base of the whip arm inside the instrument bridge.

The cable in the centre of the whip arm will tighten or loosen. Be careful not to pinch your fingers in the bottom opening.




#### 6.2 Check the Instrument Hoses for Wear

Lift one instrument and press the foot pedal to start the instrument's function.

While the instrument is running, rotate and pull the handpiece in all directions.

Flex and rotate the instrument hose to check for internal wire wear.

Listen and watch the handpiece to see if the instrument continues to function normally in all positions.




If the handpiece is loud or vibrates, the rotary parts may need replacing. Contact Finndent.

Repeat the test for all instruments.

#### 6.3 Replacing an Instrument Hose


Each instrument has a unique hose with different supply lines inside. When ordering a replacement, specify the make and model of the instrument.

Replacement hoses are supplied with no electrical connectors on the end. These must be added before installation. Check the old hose to see where and what type of plugs to install, or ask Finndent for support.



6 SERVICE INSTRUCTIONS 6.5 Instrument O-rings

# Examples of replacement hose connectors and internal wires/tubes:



**Table: Replacement hoses for Phantom Table** 

The Luzzani 3F syringe is attached to the instrument hose. The handpice must be replaced at the same time as the hose.

Luzzani 3F Straight Syringe P/N: 8900199.

# 6.4 Check Instrument LED Lights

While checking the instrument hoses for wear, watch to see if the LED lights work. The LED may need to be activated with a separate button on the handpiece.

If the light never turns on, the LED bulb may be in need of replacing.

If the light stops and starts when the hose is moved, there is possibly internal damage.

Remove the bulb and check the bulb voltage with a multimeter. Bulb voltages will vary, read the manufacturer's instructions for correct voltage.

If the voltage is 0, change the bulb and check again.

#### 6.5 Instrument O-rings

Change all instrument O-rings twice yearly.

To check the O-rings, open the instrument handpiece.


Visually inspect the O-rings (turbine for example locations).

If an O-ring needs replacing, use a blunt, plastic tool to remove.

O-rings may need lubrication. Use the OEM recommended lubricant.

Incorrect lubricant can deteriorate or cause the rubber to swell.

For example,




# 6.6 Check Instrument Tips for Wear



Instruments that use removable tips normally come with a card or instructions on how to check for wear. Follow the manufacturer's instructions for when to replace the tips.

Also check the instruments for signs of damage from improper care, such as:

- Rust or corrosion (cause: mixing carbon and stainless steel instruments)
- Chemical pitting (cause: improper cleaning products)
- Water spots or scum (cause: sterilizer not clean, ultrasonic bath not clean, insufficient rinsing or drying)



Red: Unacceptable Wear - Tip is badly worn. Tip replacement is necessary.

# 6.7 Clean the Syringe Tip



Follow the tip manufacturer's instructions for how to clean tips. Suggested methods may include an ultrasonic bath using a mild disinfection solution, steam autoclave and hand scrubbing.

While cleaning, inspect the tip for rust, chemical pitting and other signs of damage.

#### 6.8 Lubricate the Contra-Angles



This is a generic instruction only. Follow the manufacturer's instructions for lubricating instruments. Use only approved lubricants.

Incorrect lubrication and cleaning can cause liquid to enter the bearings. This can cause corrosion, bearing damage and handpiece overheating.

Remove the handpiece from the instrument hose.

Insert the lubricant spray nozzle in to the open end of the handpiece. Wrap a paper towel around the tip of the handpiece.

Spray for 2-3 seconds. Check the towel. If the oil is not clear, repeat spray until all debris is gone.

Using a smaller spray nozzle, spray in to the tip end of the handpiece.



Grease left on the instrument can affect the dental treatment. Make sure no grease remains on the outside of the handpiece.





# 6.9 Lubricate Turbine (Generic Instructions Only)



Follow the manufacturer's instructions for lubricating instruments. Use only approved lubricants.

Incorrect lubrication and cleaning can cause liquid to enter the bearings. This can cause corrosion, bearing damage and handpiece overheating.

Remove the turbine handpiece from the instrument hose.

Insert the lubricant spray nozzle in to the AIR intake of the handpiece only. The air intake is normally the smaller of the two intakes.

Wrap a paper towel around the tip of the handpiece, making sure to cover the air outlet.

Spray for 2-3 seconds. Check the towel. If the oil is not clear, repeat spray until all debris is gone.

Using a smaller spray nozzle, spray in to the tip end of the handpiece.



Air intake

## 6.10 After Lubrication (Generic Instructions Only)

Stand handpieces up on a towel so the extra grease can drain away. Wipe away any grease.

Connect the handpiece and run for approximately 30 seconds over a towel.

Wipe away any grease that comes out.

Run the handpiece again until no grease comes out.



Grease left on the instrument can affect the dental treatment. Make sure no grease remains on the outside of the handpiece.




Stand upright

# 6.11 Check Mechanical Functions of Instrument Bridge

Hold the instrument bridge handle and move the arm up. Release the handle and check that the pneumatic brake has activated to hold the arm in place. The pneumatic brake should make a quick air noise when activating and deactivating.

Rotate the instrument bridge head. Rotate the instrument tray. The joints should be easily moved and silent.

If the tray joint is worn, it may cause the tray to be not level. The plastic bearing may be replaced.



If the instrument bridge and tray are not level, they can be adjusted with screws in the instrument arm (see instrument arm section).

#### 6.12 Check for Leaks Inside

Open the instrument bridge cover:

The hanging hose cover is held in place by screws on the left side. Use a 4mm hex key to remove.

The whip arm model is locked with a screw in the back centre of the cover.

Check inside the instrument bridge for air and water leaks.

# 6.13 Change Arm Brake Rubber of Horizontal Arm

Brake rubber P/N: 6809730

For detailed, step-by-step instructions on how to disassemble the horizontal arm, see the arm brake technical bulletin (20140130 Technical Bulletin Arm Brake Repair v2). If the arm brake is not locking, and slowly moving out of place, the brake rubber may need changing. It is suggested to change the rubber yearly.

- 1. To remove the brake, the back axle, spring and brake screws must be removed.
- 2. Remove the brake end caps, pull out the plastic bushing and discard the old rubber.
  - 3. Install the new rubber.
  - 4. Re-assemble the horizontal arm.



Centre rubber in the brake body.



Fold rubber over the edge.



Slide in the plastic sleeve.




Brake / grease nipple


#### 6.14 Calibrate the Horizontal Arm Motion

To adjust the tilt of the instrument bridge or make the tray level, use a 19mm wrench on the two front nuts to adjust the tilt alignment.

To adjust the up-down vertical resistance, use a 19mm wrench on the spring nut.




To adjust the horizontal motion resistance, use the Finndent arm wrench on the horizontal arm joint:



# 6.15 Lubricate Arm Joints and Pneumatics

Every year, lubricate the horizontal arm axels and brake:



Every few years, lift the accessory post collar and grease the bearings:



# 6.16 Horizontal Arm Stop Adjustment

The adjustment is done by removing the side panels of the frame to expose the horizontal arm lower support. End stoppers are secured by two M6 screws (DIN912).


At the same time, the dental lamp mechanics may be adjusted.


#### Tools required are:

- 3 mm Hexagon key
- 5 mm Hexagon key
- Phantom cabinet door key


#### Process is done by:

- 1. Open and remove the cabinet door from the right side of the torso
- 2. remove the foot controller shelf and the side panel from the left side of the torso
- 3. remove the back panel to expose the horizontal arm lower support
- 4. Adjust the rear stopper so that the centerline of the horizontal arm does not pass the center line of the table. Adjust the front stopper so that the vertical arm and instrument bridge rotating joint does not go past the torso left side.







1. Open and remove cabinet door




2. Remove foot controller shelf and side panel




3. Remove back panel




4. Horizontal arm lower support



5. Horizontal arm limiters





# 6.17 Water Bottle Cleaning

Rinse all bottles for the bottle once per week using hot water and a disinfectant solution. A bottle brush may be used to clean residue off the inside of the bottles.

#### 6.18 Main Air Pressure Valve Filter

P/N: 9805531

Turn the main air supply off.

Turn the unit off.

- 1. Make sure the gauge reads 0 before opening.
- 2. Press the front tab in, and twist the lower half.
- 3. Pull the lower half down.
- 4. Pull the spindle, filter and O-ring up.
- 5. Replace the O-ring and filter once per year.
- 6. Replace the spindle if needed.
- 7. Reassembly the valve.

Turn the main air pressure on, and turn the unit on to check for air leaks.



# 6.19 Replace the Selection Valve Membrane

P/N: 8900255

Wear protective gloves when handling used parts.

- 1. \*\*Carefully pull out the side tabs.
- 2. Open the valve.
- 3. Replace the rubber membrane.
- 4. Make sure the membrane is pressed in to the valve's edge ring.
- 5. Close the valve and push the tabs in.



\*\*If a tab is pulled too hard, it may half-break, or pull off completely.

A broken valve will not close tightly over the membrane and must be replaced.



# 7 PCB DIAGRAMS

The PCB diagrams for Q900, Q905 and Q915 are provided for troubleshooting only. Do not attempt to rebuild a PCB. Send damaged PCBs to Finndent for repair.

# 7.1 Q900 Power Supply PCB

Version 2 AB for use with 40VAC rectifier.

	1
J1 – Input 230 VAC	J2 – Display out 230 VAC
Connect power wires	1. P GND
from the routing boxes in	2. N
the supply cabinet.	3. L
1. P GND	
2. N	
3. L	
J3 – Chair out 230 VAC	J4 – Display control (power
1. PGND	switch)
2. N	1. L (Mains)
3. L	2. L (Display)
3. L	
	· · ·
15 1: 61:	
J5 – Line filter out	J6 – Rectifier inputs
1. N	1. 40 VAC in (AC1)
2. P GND	2. 40 VAC in (AC1)
3. P GND	3. AC2 in
4. L	4. AC2 in
J7 – 36V DC out	J8 – 36V DC out
1. GND	1. GND
2. + DC	2. + DC
J9 – 36V DC out	J10 – 24 VDC out
3. GND	1. GND
4. + DC	2. 24 VDC
J11 – PGND	J12 – Motor output
1. P GND (PE)	Connect to lift
]	motor's own cable
	1. Motor out 1
	2. Motor out 2
J13 – Switch input	J14 – Suction valve power
1. +5 VDC	Connect to selection
2. Motor safety	valve with cable 4801500
switch	1. +24 VDC
3. Motor safety	2. GND via TR2
switch	
4. N.C.	
5. +5 VDC	
6. TR1 gate	
7. +5 VDC	
8. TR2 gate	
J15 – 24 VDC out	J16 – PGND
1. +24 VDC	1. P GND (PE)
2. N.C.	, ,
3. GND via TR1	
	ı




Table: Q900 pinout information

Connect the DC outs to:

4801495 – Main air pressure filter router power

4801480 - To Q915 power in, J32


# 7.2 Q905 Phantom LED PCB

Version 1.


J2 – I2C		J5 – I2C	
1.	+3.3 VDC	5.	+3.3 VDC
2.	I2C SDA	6.	I2C SDA
3.	I2C SCL	7.	I2C SCL
4.	GND	8.	GND

#### Table: Q905 pinout information


Connect display signal cable 4801490 from the Q905 J2 / J5 to Q915 J 38.




# 7.3 Q915 Phantom Instrument Control Board Version 6AD



J1 – Scale	er connector	J2 – CAN		
1.	+32 VDC	1.	CAN L	
2.	GND	2.	CAN H	
3.	N.C.	3.	N.C.	
4.	Scaler on, floating/GND via TR1	4.	Up safety, floating/GND via TR2	
5.	GND	5.	GND	
6.	Scaler control voltage	6.	+32 VDC	
7.	GND			
J3 – Arm	air brake sensor	J4 – Foo	t control signals	
1.	Bridge handle wiring	Cor	nnect foot control cable 4801409 to the Neutrik	
2.	Bridge handle wiring	pass-through plug. Connect signal cable 4801465 from		
		Neutrik	plug to Q915 J4.	
		1.	Button D (foot pedal pressed)	
		2.	Button B (centre button)	
		3.	Button A (left button, spray on-off)	
		4.	Button E (potentiometer)	
		5.	Button C (right button, motor reverse direction)	
		6.	Control voltage input (from foot control)	
		7.	Control voltage output (to foot control)	
J5 – Up s	afety switch input	J6 – Inst	rument selection 1	
1.	3.3 VDC	Cor	nnect to switch with cable 4801475	
2.	TR2 gate (see J2)	1.	GND	
		2.	Instrument selection signal input	



J7 – UART1	J8 – JTAG (For development, Not connected)
1. UART1 RX	
2. GND	
3. UART1 TX	
4. CTS	
5. DSR	
J9 – Arm air brake valve	J10 – Motor/scaler hand piece 1
1. +24 VDC	1. Motor phase 3
2. GND	2. Motor phase 2/scaler AC1
	3. Motor phase 3/scaler AC2
J11 – Air / Cooling / Water Spray – Solenoid valves 1	J12 – Instrument 1
Connect to solenoid valves with cable 4801470	From instrument hose or controller
1. +24 VDC (air valve)	1. Fiberlight +
2. GND (air valve)	2. Fiberlight –
<ol><li>+24 VDC (cooling air valve)</li></ol>	3. GND
4. GND (cooling air valve)	4. Instrument 1 ID input
5. +24 VDC (water valve)	
6. GND (water valve)	
J13 – Instrument selection 2	J14 – Motor/scaler hand piece 2
Connect to switch with cable 4801475	1. Motor phase 3
1. GND	2. Motor phase 2/scaler AC1
2. Instrument selection signal input	3. Motor phase 3/scaler AC2
J15 – Air / Cooling / Water Spray – Solenoid valves 2	J16 – Instrument 2
Connect to solenoid valves with cable 4801470	From instrument hose or controller
1. +24 VDC (air valve)	Fiber optic light +
2. GND (air valve)	2. Fiber optic light –
3. +24 VDC (cooling air valve)	3. GND
4. GND (cooling air valve)	4. Instrument 2 ID input
5. +24 VDC (water valve)	
6. GND (water valve)	
J17 – Scaler handpiece AC input	J18 – Instrument selection 3
1. AC1	Connect to switch with cable 4801475
2. N.C.	1. GND
3. N.C.	<ol><li>Instrument selection signal input</li></ol>
4. AC2	
J19 – Micromotor controller input	J20 – Motor/scaler hand piece
1. Motor phase 1	1. Motor phase 3
2. Motor phase 2	2. Motor phase 2/scaler AC1
3. Motor phase 3	3. Motor phase 3/scaler AC2
J21 – Air / Cooling / Water Spray – Solenoid valves 3	J22 – Instrument 3
Connect to solenoid valves with cable 4801470	From instrument hose or controller
1. +24 VDC (air valve – blue wires)	1. Fiberlight +
2. GND (air valve – blue wires)	2. Fiberlight –
3. +24 VDC (cooling air valve – white wires)	3. GND
4. GND (cooling air valve – white wires)	4. Instrument 3 ID input
5. +24 VDC (water valve – green wires)	
6. GND (water valve – green wires)	
J23 – Instrument selection 4	J24 – Motor/scaler hand piece
Connect to switch with cable 4801475	1. Motor phase 3
1. GND	2. Motor phase 2/scaler AC1
2. Instrument selection signal input	3. Motor phase 3/scaler AC2




J25 – Air / Cooling / Water Spray – Solenoid valves 4	J26 – Instrument 4
Connect to solenoid valves with cable 4801470	From instrument hose or controller
1. +24 VDC (air valve – blue wires)	1. Fiberlight +
2. GND (air valve – blue wires)	2. Fiberlight –
<ol><li>+24 VDC (cooling air valve – white wires)</li></ol>	3. GND
<ol><li>GND (cooling air valve – white wires)</li></ol>	4. Instrument 4 ID input
5. +24 VDC (water valve - green wires)	
6. GND (water valve - green wires)	
J27 – Instrument selection 5	J28 – Motor/scaler hand piece
Connect to switch with cable 4801475	1. Motor phase 3
1. GND	2. Motor phase 2/scaler AC1
<ol><li>Instrument selection signal input</li></ol>	3. Motor phase 3/scaler AC2
J29 – Air / Cooling / Water Spray – Solenoid valves 5	J30 – Instrument 5
Connect to solenoid valves with cable 4801470	From instrument hose
<ol> <li>+24 VDC (air valve – blue wires)</li> </ol>	1. Fiberlight +
2. GND (air valve – blue wires)	2. Fiberlight –
<ol><li>+24 VDC (cooling air valve – white wires)</li></ol>	3. GND
<ol><li>GND (cooling air valve – white wires)</li></ol>	4. Instrument 5 ID input
5. +24 VDC (water valve – green wires)	
6. GND (water valve - green wires)	
J31 – Instrument power OUT	J32 – Q915 power IN
1. GND	Connect power cable 4801480 from Q900
2. +32 VDC	1. GND
	2. +32VDC
J33 – UART5	J34 – SPI
1. UART5 TX	1. +3.3 VDC
2. GND	2. N.C.
3. UART5 RX	3. TDO
	4. SPI MISO
	5. SPI MOSI
	6. GND
J35 – Proportional valve output	J38 – I2C Display
1. GND	Connect cable 4801490 from Q915 to Q905 J2 / J5.
2. Variable voltage	1. +3.3 VDC
	2. I2C SDA
	3. I2C SCL
	4. GND

Table : Q915 pinout information

# 8 WIRING AND PNEUMATIC DIAGRAMS

8.1 Pneumatic Diagram



- 1. Clinic power, suction and air supply
- 2. 9805531 Air pressure filter-regulator with emergency off Camozzi 325869-01



- 3. Unit power supply, air and water regulation Q900 J1 - unit power IN Table Supply 9805531 Camozzi 325869-01 to Unit Supply SMC VT307-5DI-01F-Q
- 4. Suction to Dürr Comfort manifold via selection valve
- 5. Dürr Comfort Manifold high and low volume suction
- 6. Air to water bottle / water from water bottle to instrument bridge

Unit air SMC VT307-5DI-01F-Q to instrument bridge valve block

Bottle water outlet to 9801580 instrument bridge valve block

7. Water bottle air supply to pressurize the water

Unit Air 5500404 SMC ARJ1020F-M5-04-1 to **Bottle System Inlet** 



8. Instrument bridge air and water solenoid valve block:



9. Turbine oil ("dirty air") return to collection cup in unit. Red/ black hoses.


# 8.2 Wiring Diagram

The power panel and circuit breaker must have a lock-out system.

Q900 J11 and J16 Protective Earth grounds must be connected before use.



Warning: Dangerous voltage



- Building power supply 230V
   1.5mm² minimum core size, insulated power cable, per national regulation
- Main power safety switch (3800405 Katko KSM 316U)
   1.5mm² minimum core size, insulated power cable, per national regulation



3. Electrical routing (570006 FiBox ABS) to Q900 J1
1.5mm² minimum core size, insulated power cable, per national regulation

4. Power Supply System: Toroidal Transformer and Rectifier

Transformer orange to Q900 J5 1 and J5 4
Transformer yellow to Q900 J6 1
Transformer blue to Q900 J6 2
Transformer green to rectifier
Transformer brown to rectifier AC
Rectifier + red to J6 4
Rectifier black to J6 3

## 5. Power Out to Unit Systems:

4801495 Q900 J10 power to unit air solenoid valve


4801500 Q900 J14 power to selection valve 4801480 Q900 J7/8/9 power to Q915 J32 4801460 Q900 J13 to Durr Comfort Manifold Lift motor cable to Q900 J12

- 6. Air, water and power to Durr Comfort Manifold
- 7. 4801409 Foot control cable to Neutrik plug 4801465 From Neutrik plug to Q915 J4
- 8. Faro lamp power to Q900 24VDC
- 9. Instrument bridge 4801470 Instrument bridge solenoid valves to instrument air spray (blue), cooling spray (white) and water spray (green)

4801475 Instrument selection cables from whip arm microswitch to Q915 J6, J13, J18, J23 and J27

4801485 Q915 to micromotor controller 4801505 Q915 to Dassym micromotor controller 4801490 Q915 to Q905 display power and signal

10. Power switch and power socketsQ900 J4 to power switch on towerQ900 to power socket, wiring and socket per national regulation



# 9 SPARE AND CONSUMABLE PARTS

# 9.1 Consumable Parts

Some parts may need replacing before the next scheduled service. The following parts can be ordered from a Finndent dealer:

Туре	Description	Part Number	Details
	Silicone instrument rest - 5 place, whip arm	6809960	
	Instrument arm brake rubber	6809730	Gasket for the brake.
	Unit foot control cable	4801409	Only the cable.
	High volume suction handpiece 7600A010-00	8900582	
	High volume suction hose 7600A010-50	7700020	Replacement hose for large diameter suction
	Saliva ejector suction handpiece 7600A020-00	8900581	
Consumable Goods	Saliva ejector hose 7600A020-50	7700015	Replacement hose for small diameter suction
Goods	Coarse suction filter, Dürr 0725-041-11	8900511	
	Water bottle 1,5L	6000050	
	Turbine oil return cup	8900320	
	Linear slide bearings	5400075	The slide bearings on the Phantom torso lift slides
	Torso cover	6800100	Torso cover
	Headrest cover	6800105	Headrest cover
	P-G face mask with cheek, nose, mouth	6000201	
	1260M magnetic mounting plates for teeth	6000202	

This list includes the consumables that can be used with the FINNDENT Phantom unit. For proper use and maintenance please read carefully all the instructions supplied by the manufacturer of each specific item.

WARNING: Do not modify this equipment without authorization of the manufacturer.

# 9.2 Spare Parts


This list includes the spare parts that are used in the FINNDENT Phantom unit. For proper substitution and maintenance please read carefully all the instructions supplied by the manufacturer.

WARNING: Do not modify this equipment without authorization of the manufacturer.

Description	Part Number	Details
Instrument bridge handle – 70mm	6203230	
Unit foot control refurbish kit	9801932	When the foot control begins to wear out, you can fix it rather than buying new.
Unit foot control	9801930- Student_Model	To buy an all-new foot control.
Dürr Selection valve	8900255	
Lift Motor FD20	3500075	The Phantom torso lift motor.
Tilt gas spring	5600025	Gas spring to position torso.
Tilt gas spring pneumatic cable	5600026	Pneumatic cable from gas spring to lock button.
Tilt release lock button Nr. BUCHSE39x25xM16x1.5	5600022	Lock button on chest of Phantom
Q900 v2ab PCB for power supply	4802900	Power PCB in the unit electrical box.
Q905 v1PCB used together with Q915	4802905	LED Display
Q915 v6ad PCB for instrument control	4802915	Instrument control in the instrument bridge.
Handpiece and hose	8900199	Luzzani 3-way syringe, straight assembly SMMSSG, Type B instrument.
Curing light handpiece and hose	8900321	Starlight S X5 05100003, Type BF instrument.
Curing light fibre optic	8900321-1	Fibre Optic for Starlight S X5 02750009.
Curing light optical protection	8900259	Optical Protection Starlight S / SLER 05100003.
Satelec Scaler Xinetic F55094	8900116	Satelec Scaler Xinetic with Newtron Handpiece, Type BF instrument.
Satelec Scaler Xinetic Hose F00187	8900210	
Micromotor control board 251.30.001	8900371	
Finndent micromotor hose 120.35.001	8900373	
Finndent brushless micromotor FD-33-LFA 106.05.000	8900374	
Turbine coupling assembly	9801266	
Midwest Lux turbine attachment 73001	8900092	

**Table: Phantom spare parts** 

# 9.3 Assembly Drawings with Bills of Materials




3 12 7 1 9 4 MIDDLE SUPPORT BAR FOR TABLE, PHANTOM TABLE SUPPORT BAR ON FRAME, PHANTOM TABLE SUPPORT MAIN FRAME, PHANTOM 6313040 INSIDE CONNECTING U-BAR, PHANTOM TABLE SUPPORT MAIN FRAME TOP BAR TABLE MIDDLE BAR, PHANTOM TABLE LAMINATE N:O 3 TABLE END SHORT BAR 6313045 6313035 6511235 6313060 6001150 6907065 6313030 Part Code 9

FD Phantom - ACTA

PHANTOM TABLE OF 6 UNITS 9890510 TABLE UNIT FRAMES


FINNDEN





Part	Code	Name	P
1		ALLEN SCREW 912/8.8ZN M4*10 AND M4*6	'
2	0900009	6000060 TABLE FOOT ACKURAT, PHANTOM	∞
က	6313075	6313075 SHELF FOR FOOT PEDAL, PHANTOM	<u> </u>
4	6511250	6511250 RIGHT END SIDE PLATE, PHANTOM	7
2	6511260	6511260 CONNECTION PLATE FOR SIDE DOOR, PHANTOM	
9	6511270	6511270 END PLATE OF SUPPLY CABINETS	` '
7	6511275	6511275 BOTTOM FLOOR PLATE OF END SUPPLY CABINETS	
∞	6511280	6511280 FLOOR COVER PLATE FOR UNIT CENTRES	-
6	6511285	6511285 CENTRE TABLE COVERS	
10	9805560	9805560 DOOR TO TRANSFORMER BOX ASSEMBLED	
11	9805555	9805555 BACK DOOR ASSEMBLED, PHANTOM	
12	9805570	9805570 MAINS BOX ASSEMBLED	
13	9805580	9805580 UNIT WALL WITH CONNECTIONS	_






Part	Code	Name	Qty
1	3800100	NEUTRIK RJ45 SOCKET	1
2	3800400	POWER SOCKET WALL BOX	1
3	5500276	SECURING RING DURR	2
4	5500277	HOSE SOCKET15 MM DURR	1
5	5500278	HOSE SOCKET20 MM DURR	1
6	5500287	ELBOW 90 DEGREE DURR	1
7	5500305	DURR DOUBLE CONNECTOR	2
8	5500306	RETAINING NUT DURR CONNECTOR	2
9	5901001	O-RING 20 X 2	6
10	6001005	TUBE CLAMP ABA 15-24	2
11	6511255	PHANTOM DESK COVER PLATE	1
12	8900255	SUCTION PLACE SELECTION VALVE	1




# PHANTOM TABLE OF 6 UNITS 9890510 SUPPLY CABINET VIEW



Part	Code	Name	QTY
1	5000331	ALLEN SCREW 912/8.8ZN M3*10	4
2	6000060	TABLE FOOT ACKURAT, PHANTOM	8
3	5100004	M4 NUT	_
4	6511250	RIGHT END SIDE PLATE, PHANTOM	4
5	6511260	CONNECTION PLATE FOR SIDE DOOR, PHANTOM	2
6	6511270	END PLATE OF SUPPLY CABINETS	2
7	6511275	BOTTOM FLOOR PLATE OF END SUPPLY CABINETS	2
8	9805560	DOOR TO TRANSFORMER BOX ASSEMBLED	6
9	9805555	BACK DOOR ASSEMBLED, PHANTOM	2
10	9805570	MAINS BOX ASSEMBLED	1




# 9805570 SUPPLY CABINET MAINS BOX ASSEMBLY



Part	Code	Name	QTY
1	3800405	3-POLE SAFETY SWITCH KATKO IP67	1
2	5700006	TPE GASKET FIBOX ABS150/75HG	2
3	6511230	PHANTOM MAINS CONNECTION PLATE	1
4	9805531	FILTER CONTROLLER CAMOZZI ASSEMBLED	1






Part	Code	Name	QTY
1	5500265	MANOMETER	1
2	5500383	FILTER REGULATOR CAMOZZI MD1-FR0000	1
3	5500384	MOUNTING BRACKET C114-ST CAMOZZI MD1-FR0000	1
4	5500385	THREADED CARTRIDGE MD1-A-1/4 CAMOZZI MD1-FR0000	1
5	5500386	TRIPLE DOUBLE BANJO CAMOZZI CL-7642-03-6-1/4	1
6	5500387	INTERMEDIATE JOINING CARTRIDGE CAMOZZI MD1-C	1
7	5500388	LOCKABLE ISOLATION VALVE 3/2 WAY CAMOZZI MD1-V01	1
8	5500389	SILENCER 3/2 VALVE CAMOZZI MD1-V01	1
9	5500390	PIPE COUPLER 34-BO11002 10MM-1/4UK	1
10	5500217	FILTER ELEMENT CAMOZZI	1







# 9805555 UNIT SUPPLY CABINET ACCESS DOOR WITH LOCK (38CM WIDE DOOR)




Part	Code	Name	QTY
1	5000396	COUNTERSUNK HEAD SCREW M3X6 DIN965	7
2	6000025	HINGE	1
3	6000055	DOOR LOCK ESSENTRA 491479	1
4	6511265	SIDE DOOR	1
5	6000065	KEY FOR LOCK ESSENTRA 491487	1






Part	Code	Name	QTY
1	6000055	DOOR LOCK ESSENTRA 491479	1
2	6511290	SIDE DOOR, 47CM WIDE	1
-	6000065	KEY FOR LOCK ESSENTRA 491487	1






Part	Code	Name	QTY
1	8900307-20	FARO ALYA LAMP, ARM = 91CM	1
2	6203190	SUPPORT COLLAR FOR LIGHT COLUMN	1
3	6204030	BOTTOM BEARING TO LAMP SUPPORT COLLAR	1
4	6907030	HORIZINTAL ARM, PHANTOM	1
5	9802111	PARALLELOGRAM (BRAKE) ARM, ASSEMBLED	1
6	9805545	TOWER FOR 2 UNITS ASSEMBLED, PHANTOM	1
7	6907065	PHANTOM TABLE UNIT FRAME FOR 2 UNITS	1
8	6001150	TABLE TOP, 15MM THICK	1
9	9801992	Standard Tray assembly FD-8000	1
10	9890500	INSTUMENT BRIDGE ASSEMBLED, PHANTOM	1





Part	Code	Name	QTY
1	5000028	ALLEN SCREW M6X16 DIN912	2
2	5000055	RETAINING SCREW 3X4 DIN916	3
3	5001071	ALLEN SCREW M8X10 DIN912 ZN	1
4	5100020	AXIS NUT KM 8	1
5	5400008	PRESSURE BEARING FOR BRIDGE ROTATION	1
6	5400011	PRESSURE PLATE FOR BRIDGE ROTATION	2
7	5400105	BALL BEARING FOR HORIZONTAL ARM	1
8	6809900	BOX BEARING FD-8000	1
9	6201350	PROTECTIVE BEARING COVER FD-8000	1
10	6201360	LOWER BEARING TO HORIZONTAL ARM FD-8000	1
11	6201370	SUPPORT COLLAR FOR HORIZONTAL ARM, FD-8000	1
12	6907075	ARM COLUMN BEARING SUPPORT PLATE	1
13	6410340	"Q" PIN COLLAR TO INSTRUMENT BRIDGE FD-8000	1
14	6420270	COUNTER PLATE FOR STOPPER	1
15	6420280	STOPPER TO HORIZONTAL ARM, PHANTOM	1






Part	Code	Name	Qty
1		POWER SOCKET	1
2	3200129	MAIN SWITCH MARQUARDT 1805.7109	1
3	3800400	ASSEMBLY BOX ABB HW51-F	1
4	5000063	RETAINING SCREW HEX 4 * 6 DIN 916	3
5	6201380	BEARING FRAME FD-8000	1
6	6511240	OVAL TOWER BASE PLATE, PER 2 UNITS	1
7	6600000	PHANTOM LOGO STICKER	1
8	6907070	TOWER OVAL HOUSING, 2 UNITS	1




# BRIDGE AND ATTACHMENT BOX ASSEMBLY




QTY	7	_	4	4	G WIRES 2
Name	PHANTOM BRIDGE ENCLOSURE	ATTCHMENT BOX	M6 DIN125 WASHER FRAME	M6 X 10 DIN 912 SCREW	1.5MM YELLOW-GREEN GROUNDING WIRES 2
Part Code	6904250	6903085	5200034	5000027	3901045
Part	1	2	3	4	







Part	Code	Name	QTY
1	6904250	PHANTOM BRIDGE ENCLOSURE	1
2	6810450	WHITE PLASTIC BUSHING	4
3	6102680	HANDLE FRAME	2
4	6400050	BRIDGE REINFORCEMENT PLATE	2
5	6203230	SHORT BRIDGE HANDLE	2
6	5001060	HEX SCREW M3X8 DIN913	2
7	5100062	NYLOCK NUT M8 DIN985	2
8	5100007	M8 DIN934 NUT	2
	3901002	BLACK EQUIPMENT WIRE, 1.5MM	2






Part	Code	Name	QTY
1	9801715	DORIO ARMS, NORMAL, GREY	4
2	6809915	INSTRUMENT HOSE BUSHING	5
3	6904250	PHANTOM BRIDGE ENCLOSURE	1
4	5001094	HEX SCREW STANDOFF M3X20	2
5	9802200	WHIP ARM HOLDER ASSEMBLY	1
6	5000405	HEX SCREW M3X20	2
7	9801714	DORIO ARM, SOFT SILICONE, GREY	1





Part	Code	Name	QTY
1	3200138	MICROSWITCH	5
2	5000405	ALLEN SCREW M3X20 DIN912	10
3	6801000	PLASTIC HOLDER FOR DORIO ARM FD-8000+	10
4	6904255	PHANTOM DORIO HOLDER SUPPORT PLATE	1
	4801401	SELECTION MICROSWITCH CABLE	5






INSTRUMENT PCB MOUNTING PLATE

6518210

8900371


MICROMOTOR CONTROL PCB MX-30







Part	Code	Name	QTY
1	3700099	VALVE SMC FOR AIR BRAKE	1
2	5000100	PLASCTIC SCREW PLASTOFAST 1,8x16	2
3	5500216	SEALING RING M5 AVS 2661-M5	3
4	5500303	SILENCER M5	1
5	6000020	STEEL BALL 3mm, RST	1
6	6307945	NIPPLE 1,2MM, M5	2
7	6800560	SIRAI MAIN BLOCK	1
8	6307940	NIPPLE 2,7MM, M5	1











Part	Code	Name	QTY
1	4802905	Q905 LED DISPLAY PCB	2
2	5000395	COUNTERSUNK HEAD SCREW M3X4 DIN965	6
3	5000396	COUNTERSUNK HEAD SCREW M3X6 DIN965	6
4	5001113	ALLEN SCREW M3X4 DIN912	4
5	6301050	PHANTOM COVER FOR INSTRUMENT BOX	1
6	6400150	PHANTOM COVER SUPPORT, L-SHAPE	1
7	6600610	PHANTOM DISPLAY STICKER SET	2
8	6809960	INSTRUMENT PAD FD-8000, 5P	1
9	7700071	PIANO HINGE 600 x 32	1
10	4200015	LIGHT COLUMN FOR DISPLAY LED	4
	4801490	CABLE FROM Q915 TO Q905 DISPLAYS	2
















Part	Code	Name	QTY
1	5000040	ALLEN SCREW ISO 7390 AISI 304	2
2	5001070	ALLEN SCREWI M8X25 DIN912	8
3	5200137	SPRING PLATE M8 DIN 127	8
4	6001410	SLEDGE WJ200UM-01-16 DRYLIN	4
5	6420220	SLEDGE SUPPORT	1
6	6518230	SLEDGE BODY	1
7	_	SLEDGE CaRRIAGE SPACER	4











Qty 1 3200020 TORSO AND SUCTION MOTOR LIFT SWITCH 1 2 5000010 ALLEN SCREW DIN912 M5X10 2 2 3 5000020 HEX SCREW CYLINDRICAL HEAD M5 X 16 DIN912 2 4 5000408 ALLEN SCREW M3 X 6 DIN 912 5 5100020 AXIS NUT KM 8 1 2 6 5100060 NYLOCK NUT M5 DIN985 1 7 5400008 PRESSURE BEARING AXK 4060 2 8 5400011 PRESSURE PLATE AS4060 1 9 6204035 COVER PLATE FOR DÜRR COMFORT HOLDER 10 6420150 SUCTION ARM UPPER STOPPER 1 1 11 6511011 DÜRR HOLDER PLATE 1 12 6900360 SUCTION ARM, TUBE 13 8900294 Suction hose holder comfort 1 14 9805595 PHANTOM BOTTLE ON SUCTION ARM ASSEMBLED 1 1 15 6809900 BOX BEARING FD-8000 1 4801460 Q900 J13 to Durr Comfort Power Suction to Durr Comfort Water hose from bottle to instrument bridge Air hose from to water bottle





Part	Code	Name	QTY
1	5201028	RETAINING RING 38X1.5	1
2	5500216	SEALING RING M5 AVS 2661-M5	4
3	5901015	FLANGE SEALING 30 x 19 x 2	1
4	6000050	WATER BOTTLE 1,5 LTR CLEAR	1
5	6304370	AIR NOZZLE TO BOTTLE WATER FD-7000	1
6	6307940	NIPPLE 2,7/M5	2
7	6307945	NIPPLE 1,2/M5	1
8	6802010	WATER PRESSURE HEAD FD-7000/8000	1
9	6901865	BOTTLE HOLDER FD-7000	1
10	7700058	TUBE 1/4 CLEAR FOREST 0017-690	1







1	9805590	PHANTOM TORSO HEAD AND FIXING BAR ASSEMBLED	1
2	9805610	PHANTOM NECK REST ASSEMBLED	1
3	5000042	ALLEN SCREW M5x8 ISO 7380	3
4	5100060	NYLOCK NUT M5 DIN985	4
5	6518215	GAS SRPRING FASTENER PLATE	2
6	6907045	PHANTOM BACK REST	1
7	5600025	GAS SPRING WOE2K-0-175-556-0XX	1
8	6518225	BUTTON SUPPORT	1
9	6907040	PHANTOM HINGE SUPPORT	2
10	5400098	SLIDING BEARING GFM-2023-7	2
11	6313005	PHANTOM BACKREST HINGE, SIDE	2
12	6313010	PHANTOM BACKREST HINGE, MID	1
13	5000019	ALLEN SCREW 912/8.8ZN M5 X 12	16
14	6001132	SNAP IN FASTENER RICHCO CS-113-3	1
15	6510210	PHANTOM BACK REST SUPPORT BASE	1
16	5200033	WASHER PLATE M5 DIN 125	4
17	5100005	HEX NUT M5 DIN934	4
18	6313095	BUTTON NUT TO PHANTOM TORSO	8
19	6800100	CHEST FRAME	1
20	9805600	PHANTOM SUCTION HOLDER ASSEMBLED	2
21	9805615	FASTENER U-SHAPE BOLT SET 26mm	1


QTY

# > FINNDENT 9805600 THIRD HAND SUCTION HOLDER



QTY	_	_	_	7	က	_	_	_	_	~
Name	ALLEN SCREW M6X20 DIN912	FATENING SCREW M8x8 DIN913	WASHER M8 DIN125	PLASTIC WASHER M8 DIN125 PA	DISC SPRING 16x8,2x0,4	COMPRESSION SPRING 0,4x5x10	FLEXIBLE TUBE FL-S 11/5,0 x 450 mm with M6 female	PHANTOM RACK SUCTION TUBE	PHANTOM SUCTION HOLDER BUSHING	PHANTOM MICRO PIN TO SUCTION HOLDER
Code	5000023	5000125	5200035	5201030	5601038	5601039	6103105	6313085	6313090	6804025
Part Code	1	2	3	4	2	9	7	∞	6	10





Part	Code	Name	QTY
1	6000200	PHANTOM HEAD PK-2 HGB 99-001	1
2	6313080	PHANTOM HEAD FIXING BAR	1
3	9805585	PHANTOM TORSO NECK NOINT	1
4	6000201	PHANTOM P-G FACE MASK	1
5	6000202	PHANTOM 1260M MAGNETIC TEETH MOUNTING PLATE	2