



# VETEQUIP

I N C O R P O R A T E D



Research Inhalation Anesthesia  
Honest Answers from Knowledgeable People  
Products and Information

Volume VI



VetEquip, Inc.  
PO Box 10785  
Pleasanton, CA 94588

800-466-6463	Telephone
925-463-1943	Fax
www.vetequip.com	Web Site
info@vetequip.com	Email

## To Place an Order

925-463-1828	Telephone
800-466-6463	Toll Free (US & Canada)
925-463-1943	Fax
info@VetEquip.com	Email
www.vetequip.com	Online

Visa & MasterCard Accepted

Purchase Order Numbers  
Accepted from Established Clients

Purchase Order Hard Copies  
Accepted from New Clients

## Technical Information

Any information given or written by VetEquip concerning any use or application of any product is intended to be informative in nature and is believed to be true and reliable. It is not meant as a substitute for due diligence on the part of the user. Please consult proper personnel, journals, or references to ensure current validity of the information and that all applicable protocols are being met in your facility.

# Our goal is to design and build safe, efficient and user-friendly anesthesia systems for every research application.

Over the past 25 years, we have spent countless hours in literally hundreds of your facilities. We are confident that there is no other veterinary inhalation anesthesia equipment company with our knowledge and/or experience.

We ask the questions necessary to recommend the best possible anesthesia protocol for your specific needs.

What type of procedures are being done and how long are they?

What animal species are being used?

What physiologic effects of the anesthesia might interfere with the study?

How many procedures are usually done in one day? At one time?

How is the procedure area(s) set up?

Will there be MRI applications or level 3 or 4 containment studies?

Is it cost effective to have more than one person doing the procedures?

Is your facility or department under budgetary constraints?

## We Believe in Our Anesthesia Systems!

SEVEN YEAR, 100%, Total Satisfaction Guarantee!

Try for ONE YEAR. Don't like it? Get a full refund!

Use it for SIX MORE YEARS. Problems? We'll repair or replace it.

## Starting Out

Anesthesia Systems Buyer's Guide .....	2
Glossary of Anesthesia and Service Terms .....	3
VetEquip Anesthesia Book .....	4
Lab Animal Inhalation Anesthesia 101 DVD .....	5

## Complete Anesthesia Systems

COMPAC <sup>5</sup> Lab Animal Anesthesia System .....	6
Mobile COMPAC <sup>5</sup> Lab Animal Anesthesia System .....	7
COMPAC <sup>5</sup> Accessories .....	7
Rodent Circuit Controller (RC <sup>2</sup> ) .....	8
Mobile Rodent Circuit Controller (RC <sup>2</sup> ) .....	8
IMPAC <sup>5</sup> System .....	9
Basic Lab Animal Anesthesia Systems .....	10

## Anesthesia Machines and Machine Components

Anesthesia Machines .....	12
Vaporizers, Fill Devices, Inlets and Outlets .....	13

## Carrier Gas Information and Supplies

Cylinders, Descriptions, Regulators, and Accessories .....	14
Quick Disconnects, Wall Plate Connections, High Pressure Hose .....	15

## Breathing Circuits, Nosecones, Facemasks and Induction Chambers

Non-rebreathing Circuits .....	16
Nosecones, Finger Cot Kit, Facemasks, Diaphragms, Stereotaxic Nosecones .....	17
Induction Chambers .....	18
Rebreathing Circuits, Reservoir Bags .....	19
Rebreathing Machine Components .....	20

## Euthanasia Supplies

3 Stage Solution Systems .....	21
CO <sub>2</sub> Flowmeters, Regulators and Chambers .....	21
Universal Lids, Rat and Mouse Housing Cage Size .....	21

## Waste Gas Information and Supplies

Waste Gas Information, NIOSH & OSHA .....	22
Waste Gas Monitoring and Badges .....	23
Waste Control Valve, VaporGuard Activated Carbon Filters .....	24
DSx Vented, Warmed Surgery Table .....	25
Self Contained Evacuation Systems .....	26
Vented Induction Chambers .....	26
Downdraft Anesthesia Work Station .....	27

## Ventilators

Ventilators .....	28
Interfacing Ventilators with Anesthesia Machines .....	29

## Rodent Intubation Aids

Rodent Workstand, Rat and Mouse Intubation Kits .....	30
---	----


## Anesthesia System Modification and Replacement Components

Miscellaneous Machine and Warm Water Pump Components .....	31
Dual Procedure Circuit and Miscellaneous Replacement Components .....	32

## "VetEquip Gives Me Gas!" The Famous Tee Shirt Inside Back Cover

## First Time Buyer's Guide

If this is the first inhalation anesthesia system you've purchased for your facility, we invite you to begin/enhance/complete your anesthesia equipment research with this "What you need and don't need" information.



1. If you are only working with rats, mice, cats, ferrets and other like sized creatures, you do NOT need a CO<sub>2</sub> Absorber (a.k.a. rebreathing head, SodaSorb container) on your anesthesia machine.
2. You need a "carrier gas" supply. This gas, usually oxygen or medical air, carries the anesthetic vapor into the patient's lungs as it breaths. Its pressure, as it is stored in a cylinder, powers the anesthesia machine.
3. You need a carrier gas flowmeter that will allow you to deliver the carrier gas as slowly as ½ liter per minute to a nosecone or as fast as 1 to 4 liters per minute to an induction chamber.
4. You need a good, reliable and valid vaporizer. A vaporizer is to an anesthesia machine as an engine is to a car. It needs to be "agent specific," originally manufactured for the agent you want to use. The only time that fact does not hold true is when a halothane vaporizer is dedicated to isoflurane or visa versa.
5. You need an oxygen flush assembly. If you are unable to "flush out" your induction chamber you will have no choice but to always work in a non-recirculating fume hood or a chemical fume hood with an activated carbon filter. If any vendor tells you otherwise, ask for independent test results showing their chamber can be precharged and opened *or* unflushed and opened without everyone in the procedure area being exposed above OSHA's Permissible Exposure Limits (PELs).
6. You need an induction chamber with a positive pressure seal. If you don't have one, see #5 above.
7. You need activated carbon filters that are reliable. The PEL for isoflurane is 2 parts per million (ppm), the human nose can't detect isoflurane until it reaches 50 ppm. It's important to be confident that your carbon filter is working even when you can't smell the agent.
8. You need a non-rebreathing circuit that delivers the agent saturated gas from the machine to the patient, and then collects and contains the waste gas from the patient to an adequate disposal system.
9. In the United States, isoflurane Permissible Exposure Limits are an OSHA compliance issue. You must take whatever measures are necessary to keep your procedure area below the PEL of 2 ppm. Neither your procedures, the species you use, the things you're studying, the facility you're working in, or the company, agency or educational institute you work for have any affect on this fact; you must comply with the 2 ppm PEL.

**"I have 12 IMPAC's at our facility.**

**I have no reservation about training new users on VetEquip machines, even for those that have never used gas anesthesia before."**

– Kathy

*Associate Scientist  
Animal Pharmacology*

**"I've been using VetEquip for 18 years - they have a top-notch maintenance program.**

**VetEquip sets the standard that others try to achieve and I want that level of quality."**

– Suzette, RLATG, RVT

*Senior Manager  
Lab Animal Resources*



## Anesthesia System Components Glossary

### Anesthesia Machine

This term includes the frame, yoke(s), flowmeter(s), vaporizer(s), oxygen flush assembly and all hoses and connections.

### Regulator(s)

The regulator(s) are responsible for reducing the supply tank pressure to specified levels.

### Oxygen Flush Assembly

The flush assembly is responsible for purging all waste gas from the induction chamber.

### Vaporizer

The vaporizer is responsible for changing the liquid agent to a vapor which the patient can inhale. It then delivers the vapor into the moving stream of oxygen in controllable amounts.

### Reconditioned Vaporizer

In human medicine, a vaporizer that has been properly cleaned and calibrated for its specific agent and is ready for clinical use. In veterinary use, it could refer to simply a "used" vaporizer or a vaporizer which has been re-labeled for an agent other than the agent for which it was designed. It is recommended that proper questions be asked to determine each seller's definition of the term.

### Remanufactured Vaporizer

May be used in lieu of "reconditioned", frequently means some internal parts are new and some used. Often the new and used parts are from different agent-specific vaporizers, e.g. enflurane rotor valve inside an isoflurane casing. May be used to describe a vaporizer whose *only* new part is an isoflurane dial while the rest of the vaporizer is entirely used enflurane. It is recommended proper questions be asked to determine each seller's definition of this term.

## Service Definitions

### Service

A term reserved for any unscheduled work and any repairs directly related to a specific problem. "Service" and "Preventative Maintenance" are occasionally used interchangeably. "Service" is done on an "as needed" basis. "Preventative Maintenance" is regularly scheduled.

### Preventative Maintenance Service (P/M)

Annual inspection including checks of all components for proper operation and potential problems. Identification and prevention of possible future breakdowns. Replacement of rubber or other disposable goods.

### Calibration


Calibration indicates a vaporizer has been set to standard measuring graduations, under a variety of conditions. Internal components are cleaned of foreign matter and contaminants, new gaskets, seals and wicks are installed. The calibration process takes place while operating under different flowrates, temperatures and time frames.

### Efficacy Test

A meter reading of vaporizer output, determining whether or not the output matches the dial reading at both high and low oxygen flows.

### Certification

Documentation that all components of the machine have been inspected and perform to manufacturer's specifications. Guarantees that the components are working properly and should continue to perform as expected for another year.



## VetEquip Anesthesia Book

The "Need To Know" book. Have you wondered what a "MAC" value is? How is it determined and what factors affect it? Is your anesthesia machine going through more isoflurane than you think it should? How much should it be using? What is the OSHA recommendation for anesthesia waste gas? Do you need an MSDS for isoflurane on file? How about one for the activated carbon in that disposable canister? Did you know you can keep that little nosecone in position with a piece of dental floss? Are you keeping track of how much waste gas you're exposed to? How do you do that anyway? Need a sign out sheet for your departmental anesthesia machine? Wondering if your machine is due for service? When is a machine due; once a year? Every six months? All these answers and more can be found in the VetEquip Anesthesia Book.

The book is a 3 ring binder to which you can add any information you may need. The Book is included with every Lab Animal Anesthesia System or can be purchased separately. When The Book comes with your Lab Animal System, it will include manuals and protocols for that specific system.




999999 VetEquip Anesthesia Book



## Inhalation Anesthesia 101 for the laboratory animal environment

This 45 minute DVD is an excellent training film. Created by VetEquip in cooperation with Columbia University Center for Biomedical Communication with oversight by a representative of the AALAS On Line Learning Center. It qualifies the viewer for one unit of AALAS CE. Every component of a traditional anesthesia machine is described with clear instructions on how and why it is necessary. A basic rat anesthesia protocol is followed step by step. Special attention is given to waste gas control throughout. The DVD is included with every VetEquip Anesthesia Book, or can be purchased separately.




999998 Anesthesia 101 DVD



## The VetEquip COMPAC<sup>5</sup>

A revolutionary anesthesia system that combines the best features of every VetEquip system! The COMPAC<sup>5</sup> has two independent breathing circuits and one to three independent vented induction chambers.

- All circuits, chamber or breathing, are independent of each other
- All circuits have a "flip-switch" choice of ½ LPM, 1 LPM or OFF
- Chambers are heated with standard warm water pump
- All chambers and circuits are vented; no pollution, no chamber flushing
- Set it on a counter top, carry it to the next room, tuck it in a cupboard, it's 13" x 13" x 15"



901812 VetEquip COMPAC<sup>5</sup> Complete Anesthesia System  
Includes Warm Water Pump (# 983105) and Extract All (# 933105)

**Ask About Available Trade In Opportunities!**



## The VetEquip Mobile COMPAC<sup>5</sup>

Combine the safety of the new COMPAC<sup>5</sup> with the convenience of complete portability. Exhaust cabinet and 4-wheeled, counter-balanced leg assembly, combined to create a mobile platform for your COMPAC<sup>5</sup>. Smooth rolling double-wheeled castors allow you to move your unit to any procedure room in your facility. Available as a standard cabinet base or with fold-down shelf surfaces on each side of unit for the ultimate workstation. Add a single, dual, or triple O<sub>2</sub> yoke assembly (see pg. 14) to make full use of the mobile units' capabilities. Cabinet also doubles as storage unit for frequently needed items.

901814 **Mobile COMPAC<sup>5</sup> Multi Patient Anesthesia System**

As Shown

901816 **Mobile COMPAC<sup>5</sup> Multi Patient Anesthesia System**

With dual 12" X 12" fold down shelves



## The VetEquip COMPAC<sup>5</sup> Accessories

The VetEquip COMPAC<sup>5</sup> is designed to be a complete system when used with your existing warm water pump and in-house or freestanding exhaust systems. If you do not currently have one or both, the following equipment is available to complete your COMPAC<sup>5</sup> System:

### Extract-All Cabinet



933105 **Extract-All Cabinet**  
Cabinet has double filters  
18" x 42 1/2" x 13 not incl.  
rear filter

or

### Extract-All Bench Top



933101 **Extract-All Bench Top**  
12" x 10 1/2" x 7" not incl. filter

or

### Filter Mate Portable Extractor



935101 **Filter Mate Portable Extractor**  
14" x 23" x 22 3/4"  
Up to 250 CFM  
Houses a carbon based filter cell,  
surrounded by a negative pressure  
plenum to capture and re-filter any leaks

## Warm Water Pump



T-Pump provides a large water reservoir and temperature settings from 85° to 107° F with a high flow rate to assure uniform temperature across the pad. Has two overheat safety devices as well as a "tip over" shut off. Illuminated on / off switch and folding handle.

983105 **Gaymar Warm Water Pump**

## RC<sup>2</sup> - Rodent Circuit Controller

### All oxygen flowrates are independently preset:


- Any combination of circuits can be used without readjusting oxygen flows
- All circuits are unaffected by the resistance occurring in split lines
- Any combination of non-rebreathing circuit styles can be used
- Any circuit length or combination of lengths can be used



#### 922100 RC<sup>2</sup> - Rodent Circuit Controller System

12" wide, 6" deep, and 12" high

Includes: vaporizer; 2 induction chambers; 4 breathing circuits; 4 nose cones; 6' 22mm evacuation tubing; 10' O<sub>2</sub> hose and 6 VaporGuard filters



#### 922120 Mobile RC<sup>2</sup> - Rodent Circuit Controller System

Mobile RC<sup>2</sup> takes advantage of our time-tested anesthesia machine frame. Mount your RC<sup>2</sup> to the extruded aluminum, counter-balanced, body for the ultimate in portability. Allows the potential for having single, dual or triple O<sub>2</sub> tanks directly on the machine (see pg. 14) RC<sup>2</sup> unit can be easily removed from frame for situations where space is limited.

### Circuits 1, 2, 3, and 4 preset to ½ lpm:

- Sufficient oxygen without excessive waste
- Oxygen flows cannot be set too high causing leaks around nosecones or facemasks

### Circuits 5 and 6 preset to 1 lpm with Oxygen Flush Assemblies:


- Faster inductions, chamber saturation takes only 1 or 2 minutes
- Flushing saturated chambers complies with OSHA waste gas recommendations



System Status Indicator tells if the RC<sup>2</sup> is on or off at a glance

Specific flow to each circuit, controlled by individual on/off switch

Oxygen flush Assembly for OSHA waste gas compliance




- Perfect for teaching procedures, without teaching specifics of inhalation anesthesia
- Ideal for multiple, simultaneous procedures
- Perfect for multiple investigators doing complex procedures with no time for adjusting oxygen flows
- Customize this system by selecting your choice of chamber sizes and nosecone sizes



## IMPAC<sup>6</sup>

- Designed specifically for high volumes of procedures
- 2 Breathing Circuits and 4 Induction Chambers allow investigator flexibility in designing custom “assembly line” protocols
- Each of the 6 circuits is equipped with a simple “ON-OFF-FLUSH” switch to reduce operator errors and speed procedures
- Separate O<sub>2</sub> and CO<sub>2</sub> flowmeters and chamber make euthanasia, O<sub>2</sub> assisted recovery and O<sub>2</sub> / CO<sub>2</sub> anesthesia accessible and precise
- Line pressure gauges create constant easy oversight of in-house gas supply
- Single waste gas port and tube for less counter top clutter



901808 IMPAC<sup>6</sup> - Integrated Multi Patient Anesthesia Center

22" wide, 18" deep, and 28" high  
Weight 50 lbs.

## Laboratory Animal Anesthesia Systems

### Safety, Quality and Economy

- Every component is chosen for durability, strength and safety
- All systems guaranteed for 7 years
- All vaporizers completely new and agent specific
- Aluminum frames with baked-on surface for scratch and stain resistance
- Clear acrylic induction chamber for clear views of the patient
- Positive pressure seals eliminate chamber leakage
- The Dual Procedure Circuit allows assembly line productivity
- Used with proper protocols, badge monitoring results in "no detectable amounts found"
- You add only: carrier gas, liquid agent, and animals

Every "Lab Animal System" comes complete with a ten foot high-pressure oxygen hose, frame, flowmeter, oxygen flush assembly, vaporizer, breathing circuit, chamber, nosecones, waste gas evacuation tubing and two VaporGuard filters.



Included with all  
Laboratory Animal Anesthesia Systems



901807 **V-10 Mobile Cart Lab Animal System**  
19" x 19" x 50.5" high, 60 lbs.  
9" x 19" shelf with 3/4" lip



901806 **V-1 Table Top Lab Animal System**  
8" x 11" x 14" high, 15 lbs.



## Laboratory Animal Anesthesia Systems





901810 V-9 Wall Mount Lab Animal System

Includes all accessories shown on page 10

901809 V-10/A Mobile Cart Lab Animal System  
With CO<sub>2</sub> Absorber

Includes all accessories shown on page 10



901811 V-9/A Wall Mount Lab Animal  
System With CO<sub>2</sub> Absorber

Includes all accessories shown on page 10

## Portable and Mobile

**All Table Top, Mobile Cart, and Wall Mount Machines have an aluminum frame with a baked on finish that resists scratching and staining, a 0 to 4 lpm oxygen flowmeter, an oxygen flush assembly, vaporizer inlet and outlet, predrilled holes to fit most vaporizer models, a 10 foot, high pressure, oxygen hose with DISS (female) fittings, all necessary tubing for component connections and a standard common outlet.**

**These units do not include vaporizers, regulators, breathing circuits or chambers.**



### 901801 Table Top (V-1)

One-handed carrying with weight-balanced handle  
Wide, flat, rubber bottomed feet eliminate tipping  
Use with "E" cylinder or "H" cylinder  
Ideal for limited floor space  
Easy storage in a locked cabinet

### 901804 Mobile Cart (V-10)

Shelf 19" by 9" with a 3/4" lip  
4 wheel, 2 leg design for more stability  
Handle doubles as a hose and breathing circuit hanger  
Cross bar holds E cylinder regulator assemblies




### 901802 Wall Mount (V-9)

Compact, saves floor and counter space  
Easily bolts to most surfaces  
Easy to reach and see, but out of the way

### 901805 Mobile Cart (V-10/A)

Same as V-10 with a 1500cc capacity CO<sub>2</sub> Absorber  
Suitable for larger animals and rebreathing circuits  
Easily adapted to ventilators



### 901803 Wall Mount (V-9/A)

Same as V-9 with a 1500cc CO<sub>2</sub> Absorber  
Suitable for larger animals and rebreathing circuits  
Easily adapted to ventilators




## Vaporizers

We sell only completely new vaporizers, specifically manufactured for the agent for which they are labeled.

We do not sell vaporizers that have been "remanufactured", "converted", "rebuilt" or "re-serviced" from old enflurane vaporizers.

"Tec3" style vaporizers all look the same on the outside, different manufacturers have different quality standards for internal components. We use only the highest quality vaporizer manufacturer.



# Vaporizers




911103 Isoflurane - Funnel-fill  
911106 Sevoflurane - Funnel-fill




911104 Isoflurane - Pin-indexed  
911107 Sevoflurane - Pin-indexed  
Requires pin index fill device

*Vaporizers for other agents available – please contact us.*



Anti-Spill Bottle Top


- 911111 Isoflurane
- 911112 Halothane
- 911113 Sevoflurane



Pin Index Bottle Fill Device

- 911117 Isoflurane
- 911118 Halothane
- 911121 Sevoflurane

## Tec3 Style Vaporizer Inlet and Outlet



911124 Inlet

911125 Outlet

## Carrier Gas Supplies

"Carrier gas," usually oxygen, is supplied from a cylinder, large or small, in the procedure room, on the anesthesia cart, or piped through the walls.

If your facility has "piped in" oxygen, you will see either wall plate outlets, gas faucets or drop lines in your procedure area. See page 15 for the appropriate hose fitting.


If you do not have piped in oxygen, you will need to have a cylinder in your procedure area. Standard sizes are "E" and "H." Most facilities do not own the cylinders, they rent them from a gas supplier. The supplier will rent the cylinders but not the regulators that are required to lower the cylinder pressures to 50 psi. Anesthesia machines need non-adjustable regulators that are preset to 50 psi.

The chart below compares E and H cylinders:

	E-Cylinder	H-Cylinder
Height	30"	54"
Diameter	3.5"	8.5"
Volume - Liters	651	7107
Initial Pressure	2200 PSI	2200 PSI
Approx. Hrs Use @ 1 LPM	11	118

## Regulators

All regulators are gas specific, preset to 50 psi and have a DISS, male, gas specific, outlet.




D & E Cylinder

- 901305 Oxygen
- 901302 Carbon Dioxide
- 901304 Nitrous Oxide
- 901309 Medical Air



H & M Cylinder

- 901306 Oxygen
- 901308 Carbon Dioxide
- 901304 Nitrous Oxide
- 901307 Medical Air




201744 "E" Cylinder Wrench with security Chain




901301 Gasket, "E" Cylinder  
Sold individually

## E Cylinder yoke/regulator assemblies for V-10 and V-10/A machines



- 932021 Oxygen, single
- 932081 CO<sub>2</sub>, single



- 932022 Oxygen, double
- 932032 CO<sub>2</sub>, double



## Quick Disconnects

Quick disconnects provide a fast and easy method for moving an anesthetic machine without having to turn the gas supply on or off. They can be installed anywhere in the high pressure hose between the machine and the gas supply. The female (coupler) end is a check valve and is connected to the supply side of the hose. The male (adapter) end is connected to the machine side of the hose.

Find the connection type you need from the photos below and match that to the part number list.



Ohmeda Coupler



Ohmeda Adapter




Oxyquip Coupler



Oxyquip Adapter



Chemtron Coupler



Chemtron Adapter



Puritan Coupler



Puritan Adapter

### Ohmeda

- 931506 Oxygen Coupler x 1/8" FPT
- 931507 Oxygen Adapter x 1/8" MPT
- 931512 CO<sub>2</sub> Coupler x 1/4" Hose Barb
- 931513 CO<sub>2</sub> Adapter x 1/8" MPT

### Oxyquip

- 931511 Oxygen Coupler x DISS Male
- 931514 Oxygen Adapter x DISS Male

### Chemtron

- 931517 Oxygen Coupler x 1/8" FPT
- 931509 Oxygen Adapter x 1/8" MPT


### Puritan

- 931516 Oxygen Coupler x 1/8" FPT
- 931515 Oxygen Coupler x 1/4" MPT
- 931510 Oxygen Adapter x 1/8" MPT

### Hose

- Color coded, nylon reinforced, 1/4" ID, conductive high pressure gas supply hose. Hose assemblies can be custom made to exact specifications.

- 931503 Oxygen, Green
- 931508 Carbon Dioxide, Grey
- 931504 Evacuation, Clear
- 931502 Nitrous Oxide, Blue
- 931501 Medical Air, Yellow



### Outlet Stations

- Recess mounted or surface mounted
- Single or multiple gas configurations
- Gas Specific
- Connection type specific

- 931505 Outlet Station for Gas Supply  
Specify gas type(s), connection type and mounting style



## Breathing Circuits


*Breathing circuits carry the agent saturated gas from the anesthesia machine to the patient and from the patient to the waste gas evacuation system. They are divided into 2 categories, rebreathing and non-rebreathing.*

## Non-rebreathing Circuits

Non-rebreathing circuits are not used with CO<sub>2</sub> absorbers or standard bellows ventilators. One tube carries the gas from the machine to the patient where a connection holds a facemask, nosecone or endotracheal tube. A second, larger diameter "exhaust down tube" also acts as a reservoir of breathable gas for the patient. Larger animals, 5 lbs and up, can be safely maintained on a non-rebreathing circuit by adding a "reservoir bag" between the down tube and the evacuation system.

Modified Bain - 1/16" ID fresh gas tube inside the 1/4" waste gas tube. Nose end flared to create a built-in nosecone.

- 921461 9mm flared end
- 921462 12.5mm flared end
- 921463 14mm flared end




- 921411 Adult Bain, 60"
- 921410 Pedi Bain, 45"
- 1/4" fresh gas tube within 22mm waste gas tube



- 921412 Mapelson E
- 1/4" fresh gas tube outside the 22mm waste gas tube



- 921425 Complete Non-Rebreathing Circuit
- For use with animals up to 15 lbs.
- With 1.0 liter reservoir bag and evacuation valve



- 921401 Bag Evacuation Valve




## Bubble Tubing Nosecones

- Create a snug seal around the animal's nose
- Cut dead space to a minimum
- Allows access to eyes, ears and top of head
- All sizes can be stored in sterile instrument solution
- Fits all standard non-rebreathing circuits

921609 **Small to Medium Mouse, 9mm**

921612 **Large Mouse, Small to Medium Rat, 12mm**

921614 **Large Rat, 14mm**



## Finger Cots

921610 **Finger Cot Nose Cone Kit**

Includes 2 19x22mm adapters  
and 12 finger cots

190501 **Finger Cot, Large,  
Black Latex**



Finger Cot installation on Standard Bain  
Non-rebreathing Circuit



Finger Cot installation on Standard  
Mapleson E Non-rebreathing Circuit

## Facemasks

- Clear cone for full visualization of muzzle
- Flexible and replaceable rubber diaphragms
- Standard 15mm connector fits all breathing circuits

921431 **Small Rodent Mask**

921430 **Medium Rodent Mask**

921429 **Large Rodent Mask**

921428 **Small Canine Mask**

921427 **Large Canine Mask**

921415 **Replacement Diaphragm**

921416 **Replacement Diaphragm**

921417 **Replacement Diaphragm**


921418 **Replacement Diaphragm**

921419 **Replacement Diaphragm**



## Stereotaxic Frame Nosecones

- Fits most stereotaxic frames
- Allows use of bite bar and ear bars
- Locks onto frame



921432 **Mouse Size**

921434 **Rat Size**

## Anesthesia Induction Chambers

- Induction Chambers that don't pollute the procedure area
- Positive pressure gaskets at lids
- Clear, durable, acrylic construction
- Large, easy to use latches
- Non-skid feet
- Standard anesthesia size inlet and outlet ports
- Both ports on one side to reduce tubing in procedure area
- Ports positioned above level where sleeping animals could create occlusions




941443	<b>1 Liter Chamber</b> (ID 3.75"W x 4.50"D x 3.75"H)	705103	<b>Replacement Gasket</b>
941444	<b>2 Liter Chamber</b> (ID 3.75"W x 9.00"D x 3.75"H)	705104	<b>Replacement Gasket</b>
941448	<b>7 Liter Chamber</b> (ID 8.25"W x 8.75"D x 5.75"H)	705105	<b>Replacement Gasket</b>
941454	<b>15 Liter Chamber</b> (ID 15.00"W x 7.50"D x 7.7"H)	705106	<b>Replacement Gasket</b>


## Vented Induction Chambers

Constant exhaust over the opening of these chambers assure investigator safety.\* Excess anesthetic vapor is completely captured during all phases of the induction process, even while removing the animal. Connect to non-recirculating exhaust, or pair with one of our evacuation units on pg. 26 for the latest in anesthesia safety. Available in the standard 2-Liter size and introducing our new 7-Liter model for larger animals or groups of animals. Convenient slide top lid and standard inlet and outlet fittings. 6ft of 30mm evacuation tubing included.

\*Independent test data available upon request @ [800-466-6463](tel:800-466-6463) or via email at [info@vetequip.com](mailto:info@vetequip.com)




942102 **2-Liter, OD 5.5"W x 13.0"D (slopes) 3.2" to 4.2" H**  
Fresh gas inlet—standard 15mm x 22mm  
Waste gas outlet—30mm OD



941447 **7-Liter Vented Induction Chamber**  
Fresh gas inlet—standard 15mm x 22mm  
Waste gas outlet—30mm OD


## Replacement Chambers

IMPAC<sup>6</sup>



941456 **IMPAC<sup>6</sup>, Divided**  
941455 **IMPAC<sup>6</sup>, Undivided**

COMPAC<sup>5</sup>




941446 **COMPAC<sup>5</sup> Chamber Complete, 2 dividers, 5 lids**  
704306 **COMPAC<sup>5</sup> Chamber Divider**  
704305 **COMPAC<sup>5</sup> Chamber Lid, large**  
704304 **COMPAC<sup>5</sup> Chamber Lid, medium**  
704303 **COMPAC<sup>5</sup> Chamber Lid, small**




## Rebreathing Circuits

Rebreathing circuits require a CO<sub>2</sub> Absorber and can be used with a ventilator. Two hoses or tubes are connected to the anesthesia machine. They meet in a “Y” at the patient’s nosecone, endotracheal tube or face mask. The gas is delivered to the patient via one hose and returned to the machine via the other. The gas then passes through the CO<sub>2</sub> absorber where CO<sub>2</sub> is removed, fresh gas is then added to the flow, and the gas goes back to the patient again.

Open ends of hoses are standard 22mm, designed to push on to the inhalation and exhalation ports of any standard anesthesia machine. The “Y” piece at the patient end connects to any elbow, endotracheal tube, face mask or nosecone.



921413 Pediatric set, 36" hoses,  
5/8" ID, 22mm fittings



921414 Adult set, 48" hoses,  
1" ID, 22mm fittings

## Reservoir Bags

Reservoir bags are used on both rebreathing and non-rebreathing circuits. With rebreathing circuits, the bag hangs from the CO<sub>2</sub> Absorber, creating an extra reservoir of “breathable” air. With non-rebreathing circuits, the bag attaches to the end of the waste gas tube, again creating an extra reservoir of “breathable” air. Since the waste gas tube itself is a reservoir of breathable air, smaller animals, such as rodents, may not need the additional reservoir of a bag.

Standard Reservoir Bags  
22mm neck with a flat oval shape

- 921402 1/2 liter
- 921403 1 liter
- 921404 2 liter
- 921405 3 liter
- 921406 4 liter
- 921407 5 liter

Reinforced Reservoir Bags  
22mm stretch neck and reinforced shoulder

- 921471 1 liter
- 921472 2 liter
- 921473 3 liter
- 921475 5 liter



Standard

Reinforced

## Rebreathing Machine Components

A “rebreathing circuit” allows the oxygen (or medical air) and anesthetic vapor mixture from the anesthesia machine to be delivered to the patient, contained and collected as it passes (or is exhaled by) the patient, has the carbon dioxide removed from it, and then is re-delivered to the patient.

The reasons for using this anesthesia delivery and re-delivery method are: the total amount of anesthetic and oxygen (or medical air) used is reduced; some heat lost with exhaled air is re-introduced to the patient, reducing the overall heat loss; some moisture lost with the exhaled air is re-introduced to the patient, lessening dehydration.

The “rebreathing” method is rarely used with small animals as total air volume passing through the patient’s lungs must be significant for the physiologic advantages (heat and moisture retention) to be effective. For those animals in excess of 35lbs to 50lbs (or 13kg to 20kg,) the advantages can be effective.

## CO<sub>2</sub> Absorber

For either the Mobile Cart or Wall Mount Machine. Holds the CO<sub>2</sub> absorbent, and with its uni-directional valves in the inhalation and exhalation ports, forces the moving gas to follow a circular pattern to and from the patient.

- Direct reading manometer displays pressure within the breathing circuit
- Anti-pollution (pop-off) valve allows breathing circuit pressure adjustments
- 1500cc absorbent container, easily removable, durable, clear



901751 CO<sub>2</sub> Absorber



401733 Anti-Pollution (Pop Off) Valve

The Anti Pollution Valve (Pop-Off) interfaces between the anesthesia machine and the waste gas evacuation system. It collects the waste gas from the CO<sub>2</sub> absorber, allows in-circuit pressure adjustments and has a 19mm port for a waste gas tubing attachment.



## 3 Stage Solution System

As more facilities opt for beginning CO<sub>2</sub> euthanasia with measured percentages of O<sub>2</sub> and CO<sub>2</sub>, economically priced systems become necessary. The Complete 3 Stage Solution System, with its individual O<sub>2</sub> and CO<sub>2</sub> flowmeter can be used for:

- Oxygen supported recovery
- Short term anesthesia
- Humane euthanasia

901703 **3 Stage Solution System, E or D Cylinders**  
(cylinders not included)



901702 **3 Stage Solution System, H or M Cylinders**  
(cylinders not included)



Both Setups Include

12lpm CO<sub>2</sub> Flowmeter/Regulator assembly, 15lpm O<sub>2</sub> Flowmeter/Regulator assembly, Recovery/CO<sub>2</sub> Mouse Housing Cage, Universal lids for rat and mouse cages, 8' tubing which connects to Flowmeter and universal lids. E or D Cylinder set up includes cylinder wrench (see pg. 14) and dual cart.



941441 **Chamber, Recovery/CO<sub>2</sub> Mouse Housing Cage**



941442 **Universal Lid, Mouse Cage**  
11.5" x 7.25"



941445 **Universal Lid, Rat Cage**  
18.75" x 9.75"



211200 **Dual "E" Cylinder Cart**

Universal lids fit on any standard mouse or rat housing cage. Excellent for either O<sub>2</sub> recovery or CO<sub>2</sub> euthanasia. No waste gas evacuation system is necessary for CO<sub>2</sub> and/or O<sub>2</sub>. These chambers do not have positive pressure seals and consequently cannot be used for anesthesia induction.



401722 **CO<sub>2</sub> flowmeter 0 to 12 lpm**  
(As Shown)

401745 **12lpm CO<sub>2</sub> Flowmeter/Regulator assembly, E or D Cylinders**

401740 **12lpm CO<sub>2</sub> Flowmeter/Regulator assembly, H or M Cylinders**



401720 **O<sub>2</sub> flowmeter 0 to 15 lpm**  
(As Shown)

401746 **15lpm O<sub>2</sub> Flowmeter/Regulator assembly, E or D Cylinders**

401748 **15lpm O<sub>2</sub> Flowmeter/Regulator assembly, H or M Cylinders**

## Anesthetic Waste Gas Overview

Anesthetic waste gases include any gas from the anesthesia machine that passes the patient without being inhaled, as well as any gases the patient exhales. Trace amounts can pollute the procedure area due to leaking equipment, spills while filling vaporizers, saturated charcoal canisters and inappropriate system use. It is assumed that prolonged and repetitive exposure to anesthetic gases can be toxic. Consequently, NIOSH has recommended Permissible Exposure Limits (PEL's) for anesthetic agents. The PEL established for halogenated anesthetics is 2 parts per million (PPM).

### Containing Anesthetic Gas

With respect to specific methods or products, OSHA makes no recommendations as to how to comply with these limits. All anesthetic gases must be contained as they travel to and from the patient. Since gas is constantly being delivered to the patient, there is always an excess that needs to be either evacuated from the area or filtered such that all the anesthetic agent is removed.

### Evacuating Gases

Evacuation of waste gas from the area without filtering means that the gas must continue to be contained until it can be released to the outside atmosphere. Channeling to the outside requires oversight to ensure that all connections are leak free. Testing must be done to assure that the number of bends in the evacuation line, as well as the overall length, does not cause resistance to the flow of gas from the anesthesia machine. Evacuation lines must be checked periodically to assure they are clear of obstructions.

### Filtering Gases

Filtering means the gas is run through an adsorbent material such as activated charcoal, which collects and holds onto the anesthetic molecules. Once the charcoal is saturated, anesthetic gases which pass through will NOT be adsorbed. Protocols should include weighing small charcoal canisters before and after each use.

### Passive vs. Active Evacuation

Passive evacuation occurs when the gas flow from the anesthesia machine pushes the waste gases to the outside atmosphere or a filtering canister. Active evacuation involves applying negative pressure (suction) to the evacuation circuit. In the passive system, a leak in the evacuation circuit could cause the waste gases to leak "outward" into the procedure area. In the active system, if negative pressure reaches the patient's nosecone, it may pull the anesthetic gases past the nosecone, depriving the patient of the delivered concentration of gas. If operating with a loose fitting mask, room air may be pulled into the nosecone area, diluting the set concentration of anesthetic.

### In-facility Vacuum, Fume Hoods & Down Draft Tables

There are two common ways to use active evacuation. The first utilizes an in-facility exhaust system and a scavenging interface device. Tubing from the in-house exhaust system and exhaust tubing from the breathing circuit are connected directly to the interface. This device constantly attenuates pressure such that negative pressure is never applied directly to the anesthesia machine or patient. It should also have a provision to avoid accidental positive pressure build-up in the breathing circuit. This is the standard method in human hospital surgery suites. The second method is to apply the negative pressure to the procedure area (e.g. fume hoods and down draft tables) in a manner that allows any waste gas to be pulled into the stream of moving air, but does not apply negative pressure to the patient or anesthesia system. This method is very popular in rodent surgery as it lends itself to the use of less obtrusive breathing circuits.

### Monitoring Anesthesia Exposure

Halogenated anesthetics, including desflurane, sevoflurane, isoflurane, halothane, and methoxyflurane are subject to NIOSH's Recommended Exposure Limits (REL's). The REL for these gases is 2 parts per million (PPM). This limit is based on a 1978 NIOSH study examining the effects of prolonged and repeated exposure to methoxyflurane and halothane. Although no similar study has been done with more current agents, we must still comply with the standard of 2 PPM.



## Testing Methods

OSHA enforces the NIOSH recommendation without outlining any specific method for monitoring these levels. OSHA inspectors use canisters to collect an air sample and send the sample back to a laboratory for analysis. Safety assessment personnel at some larger facilities are using expensive real-time monitors. Analytical chemical laboratories, which issue and analyze badges, use a time-weighted average or other method of analysis to determine exposure. AAALAC's site inspection committee has requested proof of compliance. Without a specific guideline from OSHA, some type of monitoring is always better than none. It has been our experience that OSHA has been willing to accept all of these methods.

## Real-time monitors

Real-time monitors show exactly how many PPM's are being detected as the test is being conducted. While precise, the spiked readings tend to be higher than the time-weighted average from the badge readings. The monitor sensors can be placed directly at nosecones or induction chamber lids where readings will be higher than the exposure at the investigator's breathing area, many inches away. The expense of real-time monitors may make them cost-prohibitive.

## Badge Analysis

Badges are the monitoring method of choice for most facilities. The badge is clipped to the user's shirt collar where it adsorbs samples of the air near the breathing area. The badge is then sent to a laboratory for analysis and reporting.


## Frequency

Monitoring via badges, air samples or real-time devices is not done daily or weekly. Most facilities opt to have one person in each procedure area wear a badge once every six months. You and your safety assessment committee may opt for more frequent monitoring if results exceed PEL's, procedure protocols change, if there is a change in personnel or equipment, or as needed.

## Monitoring Badge



931443 Waste Gas Monitoring Badge and Analysis Report



**"I've worked with VetEquip for years."**

**I recommend VetEquip whenever anyone is setting up a lab because they know the anesthesia end.**

**I feel VetEquip partners with me and has an interest in what I'm doing, as opposed to wanting to just sell me something."**

– Diane, RVT, LATG  
Research Coordinator  
Radiological Sciences  
Laboratory


## Waste Control Valves

The scavenging interface serves to prevent potentially dangerous increases or decreases of pressure in the anesthetic waste gas disposal system from reaching the patient's breathing circuit.

Positive pressure relief must be provided to protect the equipment and patient if occlusion of the scavenging system outlet occurs. If the scavenging system outlet becomes occluded, the positive-pressure relief vent opens to prevent transmission of high pressure to the breathing circuit.

Negative pressure relief is needed to prevent negative pressure (suction) generated by the disposal system from reaching the patient's breathing circuit.

A reservoir is necessary to allow the scavenging system to accommodate an increased volume of excess anesthetic gas, which may transiently exceed the per-minute removal capacity of the system. It may also serve as a monitor of the scavenging system if the reservoir is a distensible bag. Overdistension of the bag could indicate inadequate function of the system and the need to adjust the needle valve to allow more gas to flow through.



### 931440 Scavenging Interface Kit


Includes: 3L reservoir bag with 19mm neck, 6 each of 1/4" ID suction hose and 19mm evacuation tubing.

### 931439 Scavenging Interface Valve only

### 921405 3.0 Liter Scavenging Reservoir Bag



## VaporGuard Filter




931401 VaporGuard Filters  
Sold in cases of 8

\*Independent test data available upon request @ [800-466-6463](tel:800-466-6463) or via email at [info@vetequip.com](mailto:info@vetequip.com)



## DSx Vented, Warmed Surgery Table

- Connects to any standard warm water pump for maintenance of body heat without danger of burns
- Connects to in-house or freestanding evacuation system for collection of any waste gas leaking around facemasks or nosecones
- Stainless steel surface is an excellent heat conductor ensuring immediate and adequate surface heat
- Stainless steel construction stands up to disinfectants, hand scrubbing and cage washer heat
- Replaceable filters and removable under-table channel cover guarantee easy and thorough cleaning
- 16"W x 11.5"D x 1.25" H; Weight - 9.5 lbs



965500 DSx Vented, Warmed Surgery Table

\*Independent test data available upon request @ [800-466-6463](tel:800-466-6463) or via email at [info@vetequip.com](mailto:info@vetequip.com)

## DSx Accessories



983105 H<sub>2</sub>O Pump



933101 Extract-All Bench Top



933105 Extract All-Cabinet

## Self Contained Evacuation Systems

Self-contained recirculating air cleaning systems plug in to any 110 outlet to remove contaminated waste gas at its source.

### Filter Mate Portable Extractor



- 935101 **Filter Mate Portable Extractor**  
14" x 23" x 22½" Up to 250 CFM  
Houses a carbon based filter cell,  
surrounded by a negative pressure  
plenum to capture and re-filter any leaks

### Extract-All Cabinet



or

- 933105 **Extract-All Cabinet**  
Cabinet has double filters  
18" x 42½" x 13 not incld.  
rear filter

or

### Extract-All Bench Top



- 933101 **Extract-All Bench Top**  
12" x 10½" x 7" not incld. filter

## Accessories For Evacuation Systems

Filter Mate Portable Extractor (935101)

- 935105 **FilterMate Replacement Filter**  
935106 **FilterMate Thimble Connection**  
connects FilterMate to existing  
house exhaust

Extract-All Cabinet (935105)

- 933106 **Primary Filter**  
933107 **Secondary Filter**  
933108 **Primary Filter refill**  
933109 **Secondary Filter refill**


Extract-All Bench Top (933101)

- 933106 **Primary Filter**


## Vented Induction Chambers

Constant exhaust over the opening of these chambers assure investigator safety.\* Excess anesthetic vapor is completely captured during all phases of the induction process, even while removing the animal. Connect to non-recirculating exhaust, or pair with one of our evacuation units on pg. 26 for the latest in anesthesia safety. Available in the standard 2-Liter size and introducing our new 7-Liter model for larger animals or groups of animals. Convenient slide top lid and standard inlet and outlet fittings. 6ft of 30mm evacuation tubing included.

\*Independent test data available upon request @ [800-466-6463](tel:800-466-6463) or via email at [info@vetequip.com](mailto:info@vetequip.com)



- 942102 **2-Liter, OD 5.5"W x 13.0"D (slopes) 3.2" to 4.2" H**  
Fresh gas inlet—standard 15mm x 22mm  
Waste gas outlet—30mm OD



- 941447 **7-Liter Vented Induction Chamber**  
Fresh gas inlet—standard 15mm x 22mm  
Waste gas outlet—30mm OD



## Downdraft Anesthesia Work Station

The Downdraft Anesthesia Work Station is specifically designed to provide a small bench-mounted unit with unrestricted access for those procedures that are difficult to perform in a fume cabinet. The downflow action takes the contaminated air away from the operator and an audible alarm will alert the operator when the airflow falls to an unacceptable level.

### Features

- No ducting or installation costs
- Integral lighting illuminates the work surface
- All mechanisms in the head are on the clean side of the filter, with the switches and electrical components being isolated from any contamination
- The filter clamping mechanism enables the filter to be easily located and ensures an even seal at the filter face
- The main filter block can be chosen from fourteen different types which include speciality filters for vapors of acids, mercury and formaldehyde
- Total airflow = 145 CFM
- Work surface face velocity = 80 FPM
- Charcoal filter is easily removed for regular weight/saturation checks

966001 [Downdraft Anesthesia Work Station, 24"](#)

966002 [Replacement carbon filter for 24" DWS](#)


966003 [Downdraft Anesthesia Work Station, 36"](#)

966004 [Replacement carbon filter for 36" DWS](#)

Other filters available upon request.




*Cut away sides give extra access to the work surface.*



*Clean air recirculates into the laboratory, so there's no heat loss or costly ducting.*

### Easy to use Control Panel

- Includes visual alarm for low airflow to indicate pre-filter blockage
- Automatic "Time-in-Use" counter
- Optional electronic monitoring of both airflow and filter saturation



*Stainless steel tray is removable for easy cleaning and access to the pre-filter.*

## Anesthesia Work Station

- Basic Respiration or Volume Ventilator during anesthesia
- Animals from 150gm to 7kg
- Adjustable high pressure limit
- Airway pressure monitor
- Oxygen flush locked with expiratory phase
- Heated rebreathing chamber
- 0 to 100 ml tidal volume displacement tube
- 300 ml CO<sub>2</sub> absorber
- Control rate and volume

921255 **Anesthesia Work Station**



## HEMC 2000 Small Animal Ventilator




- Animals 1-200kg
- Tidal volume range from 20-3000cc
- Time Cycled Respiratory Rates 6-40 bpm
- Maximum working pressure adjustable from 10-60cm/H<sub>2</sub>O
- Volume constant to compensate for lung compliance changes
- Dual airway pressure alarms monitor peak inspiratory pressure

- 921761 **HEMC 2000 with 300-1600ml bellows**
- 921762 **Bellows & Housing 0-300ml**
- 921764 **Bellows 1600-3000ml**
- 921765 **Housing 1600-3000ml**
- 921766 **Bellows 300-1600ml**
- 921767 **Housing 300-1600ml**

## MicroVent 1

- Animals with tidal volumes 0-10ml
- IPPV and High Frequency Oscillatory Ventilation
- IPPV mode allows control of:
  - Oxygen flow 50-100mls per minute
  - Breaths per minute 75-2400
  - Breath volume 0-10ml
  - Mean Airway Pressure
- Oscillatory Mode:
  - Lungs are expanded and held inflated as high frequency vibrations (750-2400 per minute) cause osmotic O<sub>2</sub> - CO<sub>2</sub> exchange without lung movements

921760 **MicroVent 1 Ventilator**





## Interfacing Ventilators with Anesthesia Machines

Interface devices have air reservoirs allowing changes in piston intake volumes to occur unaffected by constant output volumes. There is no lack of air to be pulled into the ventilator and there is no back pressure within the anesthesia system.




The Universal Control Arm acts as a pressure stabilizing interface, provides an enclosed, adjustable pop off for excess volumes (easily connected to your in-house evacuation) and a manometer for constant viewing of the pressure within the interface circuit.

The Arm can also be used with any non-rebreathing circuit for the added advantages of the enclosed pop off valve and the manometer.


921451 **Universal Control Arm**

As diagramed below, a Bain non-rebreathing circuit combined with tubing, filters and ET adapters can be used to create an interface. This system can be constructed from components found in most labs or purchased as a Complete Interface System.




921452 **Complete Interface System**  
For piston-style ventilator

## Rodent Intubation Aids



984102 **Rodent WorkStand**




984102 **Rodent WorkStand**, tilted view

The HEMC Rodent WorkStand was developed in collaboration with George A. Vogler, D.V.M., DACLAM of Saint Louis University Medical School Comparative Medicine Department to facilitate endotracheal intubation and intratracheal dose instillation and other similar procedures. The WorkStand will support rats, mice and similar small rodents in a stable and comfortable position.



984104 **Mouse Intubation Pack**



984106 **Rat Intubation Pack**

Used with the Rodent WorkStand, each pack is designed to neatly hold all the tools and supplies needed to perform quick, efficient and minimally traumatic orotracheal intubation of small rodents. Each compartment is labeled for contents, reordering information and component instructions.

Both Packs include; a specially designed and molded autoclavable intubation speculum, a lidocaine applicator, and endotracheal tube guide wire, an incisor loop and a brief tutorial video. In addition, the Rat Pack contains scissors, umbilical tape, and a mirror for the verification of tracheal intubation. The Mouse Pack includes HEMC's NEW mouse endotracheal tubes in two sizes.

Otosopes in 3 quality ranges are available by special order:

**MDS, non-rechargeable NICAD batteries.** Special Order #000A3748


**Welch Allyn, rechargeable NICAD batteries.** Special Order #000A3749

**Light and easy to handle Welch Allyn, rechargeable LI Ion batteries** Special Order #000A3754



## Modification and Replacement Components

Table Top, Mobile Cart and Wall Mount anesthesia machines can be modified with additional carrier gas flowmeters and/or with a "low flow" fine calibrated oxygen flowmeter. Other flowmeters are available via special order.



**901751 CO<sub>2</sub> Absorber**  
Removable components for easy cleaning. Enclosed anti-pollution valve to capture and vent waste anesthetic gases. Direct-reading manometer to monitor pressure in the breathing circuit.



**211200 Dual "E" Cylinder Cart**


- 411723 Oxygen, 0 to 1 lpm
- 411722 Oxygen, 0 to 4 lpm
- 401726 Oxygen, 0 to 12 lpm
- 401735 Carbon Dioxide, 0 to 10 lpm
- 401730 Carbon Dioxide, 0 to 2 lpm
- 401728 Medical Air, 0 to 5 lpm
- 401721 Nitrous Oxide, 0 to 4 lpm



**401719 Flow Control**



**983701 Click-Tite Adapters**  
Short pieces of water hose with a threaded (male) connection on one end and click-tite adapters on the other.



**983701 Warm Water Pad, 15" X 22"**  
Translucent; flexible; with leak proof, click-tite water hose connections.



**941441 Chamber, Recovery/CO<sub>2</sub> Mouse Housing Cage**




**211100 Single "E" Cylinder Cart**



## Modification and Replacement Components



314525 Adapters ET, 6mm Clear




304519 Adapters ET, 6mm, Blue




304789 Y Piece, 1/4" tube adapter



304502 Adapter, Elbow only, Mapleson E Breathing circuit




921412 Mapleson Non Rebreathing E Circuit



941444 2 liter Induction Chamber



921433 Stopcocks




931435 19mm Evacuation Tubing, Blue



931436 22mm Evacuation Tubing, Clear



901301 Gasket, "E" Cylinder  
Sold individually



172104 Silicon Tubing, 1/4" ID



171434 Latex Conductive Tubing, 1/4" ID




931530 O<sub>2</sub> Hose, 10', DISS X DISS  
Included with every LAAS



931401 VaporGuard Filters



921400 Dual Procedure Circuit, Complete



921609 9mm Nosecone

921612 12mm Nosecone

921614 14mm Nosecone




304507 Adapter, 22mm Female X 19mm Male



304501 Adapter, 15mm ID X 22mm OD

# VETEQUIP gives me gas...



"I'm not as high as some thinkle peep"

Every complete Lab Animal System comes with a VetEquip Tee Shirt. Tee shirts are also available for sale!



**999997** VetEquip Tee Shirt  
Hanes® Heavyweight 100% cotton, preshrunk, coverseamed, double needle sleeve and bottom hem. Specify size— L or XL



PO Box 10785, Pleasanton, CA 94588-0785  
Telephone: 925-463-1828 • Toll free: 800-466-6463 • Fax: 925-463-1943  
Email: [info@vetequip.com](mailto:info@vetequip.com) • [www.vetequip.com](http://www.vetequip.com)



# VetEquip, Inc.

## Satisfaction Guarantee

### One Year Money Back Guarantee:

*Try any System for a full year. Not satisfied?*

*We'll make arrangements for a full refund!*

### Seven Year Repair or Replacement:

*Problems with your System?*

*We'll repair or replace it at no charge.*



**Have Questions?**

*Call 800-466-6463 for...*

***Real People with  
Honest Answers!***

**VetEquip, Inc.**

PO Box 10785, Pleasanton, CA 94588-0785

Telephone: 925-463-1828 • Toll free: 800-466-6463 • Fax: 925-463-1943

Email: [info@vetequip.com](mailto:info@vetequip.com) • [www.vetequip.com](http://www.vetequip.com)