

Vitality® hrGFP Mammalian Expression Vectors pIRES-hrGFP-1a, pIRES-hrGFP-2a, phrGFP-C, phrGFP-N1, phrGFP-1, and phrGFP

INSTRUCTION MANUAL

Catalog #240031 (pIRES-hrGFP-1a)

#240032 (pIRES-hrGFP-2a)

#240035 (phrGFP-C)

#240036 (phrGFP-N1)

#240059 (phrGFP-1)

#240062 (phrGFP)

Revision #064007g

For In Vitro Use Only



LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

United States and Canada

Stratagene

11011 North Torrey Pines Road

La Jolla, CA 92037

Telephone (858) 535-5400**Order Toll Free** (800) 424-5444**Technical Services** (800) 894-1304**Internet** tech_services@stratagene.com**World Wide Web** www.stratagene.com

Stratagene European Contacts

| Location | Telephone | Fax | Technical Services |
|-----------------------|-----------------|--------------------|--------------------|
| Austria | 0800 292 499 | 0800 292 496 | 0800 292 498 |
| Belgium | 00800 7000 7000 | 00800 7001 7001 | 00800 7400 7400 |
| | 0800 15775 | 0800 15740 | 0800 15720 |
| France | 00800 7000 7000 | 00800 7001 7001 | 00800 7400 7400 |
| | 0800 919 288 | 0800 919 287 | 0800 919 289 |
| Germany | 00800 7000 7000 | 00800 7001 7001 | 00800 7400 7400 |
| | 0800 182 8232 | 0800 182 8231 | 0800 182 8234 |
| Netherlands | 00800 7000 7000 | 00800 7001 7001 | 00800 7400 7400 |
| | 0800 023 0446 | +31 (0)20 312 5700 | 0800 023 0448 |
| Switzerland | 00800 7000 7000 | 00800 7001 7001 | 00800 7400 7400 |
| | 0800 563 080 | 0800 563 082 | 0800 563 081 |
| United Kingdom | 00800 7000 7000 | 00800 7001 7001 | 00800 7400 7400 |
| | 0800 917 3282 | 0800 917 3283 | 0800 917 3281 |

All Other Countries

Please contact your local distributor. A complete list of distributors is available at www.stratagene.com.

Vitality® hrGFP Mammalian Expression Vectors

pIRES-hrGFP-1a, pIRES-hrGFP-2a, phrGFP-C, phrGFP-N1, phrGFP-1, and phrGFP

CONTENTS

| | |
|--------------------------------------------------------------------------------|-----------|
| Materials Provided | 1 |
| Storage Conditions | 1 |
| Additional Materials Required | 1 |
| Notices to Purchaser | 2 |
| Introduction | 5 |
| Description of the Vectors | 6 |
| The pIRES-hrGFP-1a Vector | 7 |
| The pIRES-hrGFP-2a Vector | 8 |
| The phrGFP-C Vector | 9 |
| The phrGFP-N1 Vector | 10 |
| The phrGFP-1 Vector | 11 |
| The phrGFP Vector | 12 |
| Overview of Cre-mediated Site-Specific Recombination | 13 |
| Module Descriptions | 14 |
| Preparing the hrGFP Vectors | 15 |
| Ligating the Insert | 16 |
| Transformation | 17 |
| Verification of Insert Percentage, Size, and Orientation | 17 |
| PCR Amplification of DNA from Individual Colonies..... | 17 |
| Cre-Mediated Recombination Protocol | 19 |
| Confirming Recombination | 19 |
| Mammalian Cell Transfection | 20 |
| Specifications for hrGFP and EGFP Excitation and Emission Spectra | 20 |
| Troubleshooting | 21 |
| Preparation of Media and Reagents | 21 |
| References | 22 |
| Endnotes | 22 |
| MSDS Information | 22 |

Vitality® hrGFP Mammalian Expression Vectors

pIRES-hrGFP-1a, pIRES-hrGFP-2a, phrGFP-C, phrGFP-N1, phrGFP-1, and phrGFP

MATERIALS PROVIDED

| Material provided | Concentration | Quantity |
|----------------------------------|---------------|----------|
| pIRES-hrGFP-1a (Catalog #240031) | 1.0 µg/µl | 20 µg |
| pIRES-hrGFP-2a (Catalog #240032) | 1.0 µg/µl | 20 µg |
| phrGFP-C (Catalog #240035) | 1.0 µg/µl | 20 µg |
| phrGFP-N1 (Catalog #240036) | 1.0 µg/µl | 20 µg |
| phrGFP-1 (Catalog #240059) | 1.0 µg/µl | 20 µg |
| phrGFP (Catalog #240062) | 1.0 µg/µl | 20 µg |

STORAGE CONDITIONS

All Components: -20°C

ADDITIONAL MATERIALS REQUIRED

Sterile Media and Reagents

T4 DNA ligase
Ligase buffer[§]
Taq DNA polymerase[†]
Taq DNA polymerase buffer^{††}
TE buffer[§]

Equipment

Falcon® 2059 polypropylene tubes (15-ml)
Water baths (37°C and 42°C)

OPTIONAL ADDITIONAL MATERIALS

Cre recombinase [Stratagene Catalog #600270]
pExchange module EC-Hyg [Stratagene Catalog #211181]
pExchange module EC-Puro [Stratagene Catalog #211182]
pExchange module EC-Neo [Stratagene Catalog #211183]
XL1-Blue supercompetent cells [Stratagene Catalog #200236]

[§] See *Preparation of Media and Reagents*.

^{†,††} See *Endnotes*.

NOTICES TO PURCHASER

CMV Promoter

The use of the CMV Promoter is covered under U.S. Patent Nos. 5,168,062 and 5,385,839 owned by the University of Iowa Research Foundation and licensed FOR RESEARCH USE ONLY. For further information, please contact UIRF at 319-335-4546.

FLAG[®] License Agreement

The enclosed DNA expression vector and/or antibody are specifically adapted for a method of producing selected protein molecules covered by one or more of the following patents owned by Sigma-Aldrich Co.: U.S. Patent Nos. (5,011,912, 4,703,004, 4,782,137 and 4,851,341; EP Patent No. 150,126 (Austria, Belgium, Switzerland, France, United Kingdom, Italy, Netherlands and Sweden); EP Patent No. 335,899 (Belgium, Switzerland, Germany, France, United Kingdom, Italy, Luxembourg and Sweden); German Patent No. P3584260.1; Canadian Patent No. 1,307,752; and Japanese Patent Nos. 1,983,150 and 2,665,359. Your payment includes a limited license under these patents to make only the following uses of these products:

A. Vector License: You may use the enclosed vector to transform cells to produce proteins containing the amino acid sequence DYKDDDDK for research purposes provided, however, such research purposes do not include binding an unlicensed antibody to any portion of this amino acid sequence nor using such proteins for the preparation of antibodies having an affinity for any portion of this amino acid sequence.

B. Antibody License: You may only use the enclosed antibody for research purposes to perform a method of producing a protein in which the protein is expressed in a host cell and purified by use of the antibody in accordance with a claim in one of the above patents in force in a country where the use actually occurs so long as: (1) you perform such method with a DNA expression vector licensed from Sigma-Aldrich Co.; and (2) you do not bind (or allow others to bind) an unlicensed antibody to any DYKDDDDK epitope of any fusion protein that is produced by use of the method.

This license does not include any rights under any other patents. You are not licensed to use the vector and/or antibody in any manner or for any purpose not recited above. As used above, the term “unlicensed antibody” means any antibody which Sigma-Aldrich Co. has not expressly licensed pursuant to Paragraph B, above. Sigma-Aldrich Co. hereby expressly retains all rights in the above listed patents not expressly licensed hereunder.

If the terms and conditions of this License Agreement are acceptable to you, then you may open the vessel(s) containing the vector and/or antibody and, through such act of opening a vessel, will have shown your acceptance to these terms and conditions.

If the terms and conditions of this License Agreement are not acceptable to you, then please return the vessel(s) unopened to Stratagene for a complete refund of your payment.

For additional licensing information or to receive a copy of any of the above patents, please contact the Sigma-Aldrich Co. licensing department at telephone number 314-771-5765.

IRES Sequence

Use of the translation enhancer of the pIRES-hrGFP-1a and pIRES-hrGFP-2a vectors is covered by U.S. Patent No. 4,937,190 and is limited to use solely for research purposes. Any other use of the translation enhancer of the pIRES-hrGFP-1a and pIRES-hrGFP-2a vectors requires a license from WARF. WARF can be reached at P.O. Box 7365 Madison, WI 53707-7365.

Non-Commercial Research Use License For Nonprofit Entities

Stratagene agrees to sell, and Licensee agrees to purchase, Stratagene's Vitality® fluorescent protein products provided herewith (referred to as the "Products") on the following terms and conditions. (For purposes of this License, "Licensee" shall include any person or entity which ordered the Products or at any time uses the Products.) LICENSEE'S ACCEPTANCE OF DELIVERY AND/OR USE OF THE PRODUCTS SHALL CONSTITUTE LICENSEE'S BINDING AGREEMENT TO THE FOLLOWING TERMS AND CONDITIONS. IF LICENSEE IS UNWILLING TO ACCEPT SUCH TERMS AND CONDITIONS, STRATAGENE IS WILLING TO ACCEPT RETURN OF THE PRODUCTS PRIOR TO ANY USE OF THE PRODUCTS, FOR A FULL REFUND.

1. The Products, containing DNA sequences encoding for fluorescent protein or variants thereof, are proprietary or exclusively licensed to Stratagene and licensed hereunder for non-commercial research purposes only. Licensee may modify only the non-coding region outside of the nucleic acid encoding the fluorescent protein of the Products to facilitate non-commercial research. Licensee shall not have any rights to (i) modify the coding region of the nucleic acid encoding the fluorescent protein of the Products, (ii) offer the Products, or any component, derivative or modification thereof, for resale, or (iii) distribute, transfer, or otherwise provide access to, the Products, or any component, derivative or modification thereof, to any third party for any purpose or use.
2. Except as set forth above, no other rights, express or implied, are conveyed to Licensee. No rights are granted to Licensee to use the Products for (i) the provision of services to any for-profit third party (e.g., screening and profiling), (ii) diagnostic applications, (iii) methods employed in screens to evaluate compounds (e.g., high throughput screening ("HTS")), (iv) profiling chemicals for selectivity, bioavailability, drug metabolism or toxicity, (v) use *in vivo* in multicellular organisms (and methods therein) for gene therapy, (vi) quality control or quality assurance processes, including food and environmental testing, or (vii) use in manufacturing.
3. The Products shall be used solely on premises under the control of Licensee, and in compliance with all laws, regulations, rules and guidelines applicable to the Products and their use, testing, handling, or other disposition thereof, or otherwise applicable to Licensee's activities hereunder.
4. Title to the Products shall not transfer to Licensee.
5. Stratagene warrants that, at the time of shipment, the Products will conform to the specifications which accompany the Products. This warranty limits Stratagene's liability to replacement of the Products. STRATAGENE MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, WITH RESPECT TO THE PRODUCTS, INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR THAT THE PRODUCTS DO NOT INFRINGE ANY PROPRIETARY RIGHTS OF ANY THIRD PARTY. Licensee hereby waives, releases and renounces any and all warranties, guarantees, obligations, liabilities, rights and remedies, express or implied, arising by law or otherwise, with respect to the usefulness or freedom from defects of the Products, including, but not limited to, (a) any implied warranty or merchantability or fitness for a particular purpose, (b) any implied warranty arising from course of performance, course of dealing or usage in the trade, and (c) any obligation, right, liability, claim or

remedy for (1) loss of use, revenue or profit, or any other damages, (2) infringement of third party intellectual property rights, and (3) incidental or consequential damages.

6. Licensee agrees to bear all risks associated with the Products and their use, testing, handling or other disposition thereof, and all risks associated with Licensee's use of the Products purchased hereunder. Licensee hereby assumes all risks of damage or injury to Licensee's facilities, employees or agents and to any third party arising from possession or use of the Products. Stratagene shall have no liability to Licensee, its employees or agents or to any third party, regardless of the form or theory of action (whether contract, tort or otherwise, including, but not limited to, negligence and strict liability), for any direct, indirect, consequential, incidental or other damages arising out of or relating to the Products or this License.

7. Licensee shall indemnify, defend and hold Stratagene, its affiliates, distributors, suppliers, directors, officers, employees and agents, harmless from and against any and all claims, actions, demands, liabilities, damages and expenses (including attorneys' fees) relating to or arising out of any damage or injury, including, but not limited to, product liability and intellectual property infringement claims of any nature, alleged to have been caused by the Products or the use, testing, handling or other disposition thereof or Licensee's activities hereunder.

8. Licensee may at any time properly dispose of the Products in a manner which ensures their prompt destruction and is consistent with all applicable laws, regulations, rules and guidelines.

9. No modification or waiver of any terms or conditions of this License shall be effective unless in writing signed by Licensee and an authorized representative of Stratagene. For information on acquiring a license to use the Products for commercial purposes, including commercial research purposes, please contact Stratagene, Business Development, 11011 North Torrey Pines Road, La Jolla, California 92037, telephone number (858) 535-5400, facsimile number (858) 535-0071.

INTRODUCTION

The green fluorescent protein (GFP) has become an extremely versatile tool for tracking and quantifying biological entities in the fields of biochemistry, molecular and cell biology, as well as high throughput screening and gene discovery.^{1, 2} GFPs have been identified in a wide range of coelenterates, and while recently the number of cloned GFPs has expanded, to date the best characterized proteins are those from the jellyfish *Aequorea victoria*. *Aequorea* GFP forms weak homodimers at moderate to low concentrations, and is often cytotoxic when expressed at low levels. Due to this latter characteristic, researchers have often been frustrated in their attempts to produce stable GFP-expressing cell lines using the *Aequorea* protein.³

At Stratagene we have isolated a cDNA clone for GFP from a novel marine organism, and have fully humanized the gene using codons preferred in highly expressed human genes. The fluorescence spectrum for the cloned GFP protein is essentially identical to the published spectrum for the purified native protein, with the major excitation peak at 500 nm and the emission peak at 506 nm. We have expressed the protein in a wide range of human, rodent, and simian cell lines, and observed levels of fluorescence comparable to that for the red-shifted, humanized variant of *Aequorea* GFP (EGFP) in all cell-types tested. In viability experiments, we find that high level expression of functional fluorescent protein in retrovirus-transduced cells is substantially more consistent and less toxic over time and passage number for the humanized recombinant GFP (hrGFP) than for EGFP. (See http://www.stratagene.com/vectors/gfp/vitality_gfp.htm.) Thus the stable GFP-expressing cell lines are produced much more efficiently using Stratagene's hrGFP compared with EGFP.

Stratagene offers six Vitality® hrGFP vectors for mammalian expression. These vectors support a variety of expression configurations, thus providing ideal expression options for each specific application. The hrGFP allows expressed genes to be easily visualized using fluorescence microscopy or fluorescence-activated cell sorting (FACS).

Description of the Vectors

pIRES-hrGFP-1a and -2a Vectors

The pIRES-hrGFP vectors (see Figures 1 and 2) contain a dicistronic expression cassette in which the multiple cloning site (MCS) is followed by the EMCV-IRES linked to the hrGFP coding sequence. This design allows the expression of a gene of interest to be monitored at the single-cell level due to expression of hrGFP on the same transcript. The gene of interest may be fused to three contiguous copies of either the FLAG® epitope (pIRES-hrGFP-1a) or the HA epitope (pIRES-hrGFP-2a). The vectors are derived from the vector pExchange-1, and thus take advantage of *Cre*-mediated site-specific recombination to allow quick and efficient directional insertion of prefabricated modules. With the pIRES-hrGFP-based vectors, a wide variety of drug-resistant markers may be readily substituted into the core hrGFP vector bearing the gene of interest.

Figures 1 and 2 show circular maps and locations of important features for pIRES-hrGFP-1a and pIRES-hrGFP-2a, respectively.

phrGFP-C and -N1 Vectors

The vectors phrGFP-N1 and phrGFP-C allow fusion of hrGFP at either the N-terminus or the C-terminus of a protein of interest. phrGFP-C contains a copy of the hrGFP gene downstream of the MCS, allowing fusion of hrGFP to the C-terminus of the protein of interest (Figure 3). phrGFP-N1 contains a copy of the hrGFP gene that lacks a translational termination codon inserted upstream of a versatile MCS to allow fusion of hrGFP to the N-terminus of the protein of interest (Figure 4). The vectors are derived from the vector pExchange-1, and thus take advantage of *Cre*-mediated site-specific recombination to allow quick and efficient directional insertion of prefabricated modules.


phrGFP-1 Vector

The phrGFP-1 vector contains the hrGFP gene, which includes a Kozak consensus sequence and termination codons directly between two multiple cloning sites for easy transfer of the hrGFP module to new vectors (Figure 5). This vector is derived from the vector pExchange-1, and thus takes advantage of *Cre*-mediated site-specific recombination to allow quick and efficient directional insertion of prefabricated modules.

phrGFP Vector

The phrGFP vector contains the hrGFP gene and the SV40 polyadenylation signal. A *LoxP* site allows easy integration of hygromycin, neomycin, or puromycin genes via pre-existing modules from the Exchanger system in an easy 30-minute reaction using *Cre* recombinase. The phrGFP vector lacks a eukaryotic promoter. Desired promoter/enhancer elements are inserted upstream of the hrGFP gene via the extensive multiple cloning site. The hrGFP gene itself contains a Kozak consensus sequence and termination codon. See Figure 6.

The pIRES-hrGFP-1a Vector



pIRES-hrGFP-1a Multiple Cloning Site Region (sequence shown 651–727)

Sac I* Sac II* Not I* Srf I Sma I/Xma I BamH I EcoR I
 GA GCT CCA CCG CGG TGG CGG CCG CTC TAG CCC GGG CGG ATC CGA ATT C ...
 STOP*


Sph I Sal I Xho I start of FLAG tag
 ...GC ATG CGT CGA CTC GAG GAC TAC AAG GAT

* The presence of a stop codon (TAG) in-frame with the 3×FLAG tag must be considered when inserting genes into the MCS. Do not use the Sac I, Sac II, or Not I sites for cloning unless the cloning strategy removes the stop codon by double-digestion using one of these upstream sites plus a site downstream of the stop codon.

| Feature | Nucleotide Position |
|----------------------------------------------------------------------------|---------------------|
| CMV promoter | 1–602 |
| T3 primer binding site (for 5' end of insert) [5' AATTAACCCTCACTAAAGGG 3'] | 620–639 |
| multiple cloning site | 651–715 |
| stop codon (in frame with 3× FLAG tag) | 677–679 |
| 3× FLAG tag | 716–787 |
| primer binding site (for 3' end of insert) [5' GTCCTTATCATCGTCGTCTT 3'] | 747–766 |
| internal ribosome entry site (IRES) | 823–1397 |
| hrGFP ORF | 1407–2123 |
| T7 primer binding site [5' GTAATACGACTCACTATAGGGC 3'] | 2155–2176 |
| SV40 polyA signal | 2188–2571 |
| f1 origin of ss-DNA replication | 2709–3015 |
| LoxP sequence | 3178–3211 |
| ampicillin resistance (<i>bla</i>) ORF | 3256–4113 |
| pUC origin of replication | 4260–4927 |

FIGURE 1 Features of the pIRES-hrGFP-1a Mammalian Expression Vector.

The pIRES-hrGFP-2a Vector



pIRES-hrGFP-2a Multiple Cloning Site Region
(sequence shown 651-727)

Sac I* **Sac II*** **Not I***

| | |

GA GCT CCA CCG CGG TGG CGG CCG CTC TAG CCC GGG CGG ATC CGA ATT C ...

STOP*

Srf I Sma I/Xma I BamH I EcoR I

| | | |

...GC ATG CGT CGA CTC GAG TAC CCA TAT GAC


start of HA tag

* The presence of a stop codon (TAG) in-frame with the 3×HA tag must be considered when inserting genes into the MCS. Do not use the Sac I, Sac II, or Not I sites for cloning unless the cloning strategy removes the stop codon by double-digestion using one of these upstream sites plus a site downstream of the stop codon.

| Feature | Nucleotide Position |
|----------------------------------------------------------------------------|---------------------|
| CMV promoter | 1–602 |
| T3 primer binding site (for 5′ end of insert) [5′ AATTAACCCTCACTAAAGGG 3′] | 620–639 |
| multiple cloning site | 651–715 |
| Stop codon (in-frame with 3×HA tag) | 677–679 |
| 3×HA tag | 716–796 |
| Primer binding site (for 3′ end of insert) [5′ TAAGCGTAGTCAGGTACATC 3′] | 779–798 |
| internal ribosome entry site (IRES) | 832–1406 |
| hrGFP ORF | 1416–2132 |
| T7 primer binding site [5′ GTAATACGACTCACTATAGGGC 3′] | 2164–2185 |
| SV40 polyA signal | 2197–2580 |
| f1 origin of ss-DNA replication | 2718–3024 |
| LoxP sequence | 3187–3220 |
| ampicillin resistance (<i>bla</i>) ORF | 3265–4122 |
| pUC origin of replication | 4269–4936 |

FIGURE 2 Features of the pIRES-hrGFP-2a Mammalian Expression Vector.

The phrGFP-C Vector



phrGFP-C Multiple Cloning Site Region (sequence shown 646–761)

AG CTG GAG CTC CAC CGC GGT GGC GGC CGC TCT AGC CCG GGC GGA TCC CCC GGG CTG...
 ...CAG GAA TTC GAT ATC AAG CTT ATC GAT ACC GTC GAC CTC GAG ACC...
 ...ATG GTG AGC AAG CAG
 START


Restriction sites: Sac I, Sac II, Not I, BamH I, EcoR I, EcoR V, Hind III, Cla I, Sal I, Xho I.

start of hrGFP


| Feature | Nucleotide Position |
|----------------------------------------------------------------------------|---------------------|
| CMV promoter | 1–602 |
| T3 primer binding site (for 5' end of insert) [5' AATTAACCCTCACTAAAGGG 3'] | 620–639 |
| multiple cloning site | 651–743 |
| primer binding site (for 3' end of insert) [5' ACCTTGAAGCTCATGATCTC 3'] | 786–805 |
| hrGFP ORF | 747–1463 |
| T7 primer binding site [5' GTAATACGACTCACTATAGGGC 3'] | 1495–1516 |
| SV40 polyA signal | 1528–1911 |
| f1 origin of ss-DNA replication | 2049–2355 |
| LoxP sequence | 2518–2551 |
| ampicillin resistance (<i>bla</i>) ORF | 2596–3453 |
| pUC origin of replication | 3600–4267 |

FIGURE 3 Features of the phrGFP-C Mammalian Expression Vector.

The phrGFP-N1 Vector



phrGFP-N1 Multiple Cloning Site Region (sequence shown 1308–1412)




* If restriction using Xba I is desired, the phrGFP-N1 plasmid DNA must first be purified from a *dam*[−] *E. coli* strain.

| Feature | Nucleotide Position |
|--------------------------------------------------------------------------------|---------------------|
| CMV promoter | 1–597 |
| hrGFP ORF | 606–1322 |
| 5' N1 primer binding site (for 5' end of insert) [5' CAGCTGACCAGCCTGGGCAAG 3'] | 1275–1295 |
| multiple cloning site | 1323–1406 |
| T7 primer binding site (for 3' end of insert) [5' TAATACGACTCACTATAGGG 3'] | 1439–1458 |
| SV40 polyA signal | 1471–1854 |
| f1 origin of ss-DNA replication | 1992–2298 |
| LoxP sequence | 2461–2494 |
| ampicillin resistance (<i>bla</i>) ORF | 2539–3396 |
| pUC origin of replication | 3543–4210 |

FIGURE 4 Features of the phrGFP-N1 Mammalian Expression Vector.

The phrGFP-1 Vector



phrGFP-1 Multiple Cloning Site 1 Region (sequence shown 645–704)

AAG CTG GAG CTC CAC CGC GGT GGC GGC CGC TCT AGC CCG GGC GGA TCC ACC ATG GTG AGC
 Sac I Sac II Not I Srf I Sma I/Xma I BamH I start of hrGFP
KOZAK


phrGFP-1 Multiple Cloning Site 2 Region (sequence shown 1407–1478)

end of hrGFP EcoR I EcoR V Hind III Sal I Xho I
 TGG GTG TAA TAG GAA TTC GAT ATC AAG CTT ATC GAT ACC GTC GAC CTC GAG...
 STOP STOP
 ...GGG GGG CCC GGT ACC AGG TAA
 Kpn I

| Feature | Nucleotide Position |
|------------------------------------------------------------------------------------|---------------------|
| CMV promoter | 1–602 |
| T3 primer binding site (for 5' end of insert in MCS1) [5' AATTAACCCTCACTAAAGGG 3'] | 620–639 |
| multiple cloning site 1 (Sac I to BamH I) | 651–692 |
| Kozak sequence | 691–699 |
| 5' hrGFP primer (for 3' end of insert in MCS1) [5' ACCTTGAAGCTCATGATCTC 3'] | 735–754 |
| hrGFP ORF | 696–1412 |
| 3' hrGFP primer (for 5' end of insert in MCS2) [5' CAGCTGACCAGCCTGGGCAAG 3'] | 1365–1385 |
| stop codons | 1413–1418 |
| multiple cloning site 2 (EcoR I to Kpn I) | 1419–1472 |
| T7 primer binding site (for 3' end of insert in MCS2) [5' TAATACGACTCACTATAGGG 3'] | 1493–1512 |
| SV40 polyA signal | 1525–1908 |
| f1 origin of ss-DNA replication | 2046–2352 |
| LoxP sequence | 2515–2548 |
| ampicillin resistance (<i>bla</i>) ORF | 2593–3450 |
| pUC origin of replication | 3597–4264 |

FIGURE 5 Features of the phrGFP-1 Mammalian Expression Vector.

The phrGFP Vector



phrGFP Multiple Cloning Site Region (sequence shown 1–128)

Nsi I Bgl II Xba I Sac I Sac II Not I Srf I
AT GCA TAG ATC TTC TAG ATA GAG CTC CAC CGC GGT GGC GGC CGC TCT AGC CCG GGC...

BamH I EcoR I EcoR V Hind III Sal I
...GGA TCC CCC GGG CTG CAG GAA TTC GAT ATC AAG CTT ATC GAT ACC GTC GAC...

Xho I start of hrGFP
...CTC GAG ACC ATG GTG AGC AAG

| Feature | Nucleotide Position |
|----------------------------------------------------------------------------|---------------------|
| multiple cloning site | 1–113 |
| hrGFP ORF | 117–833 |
| 5' hrGFP primer binding site [5' CCAGGTTACCTTGAAGCTCAT 3'] | 162–183 |
| T7 primer binding site [5' TAATACGACTCACTATAGGG 3'] | 866–885 |
| SV40 polyA signal | 898–1281 |
| f1 origin of ss-DNA replication | 1419–1725 |
| LoxP sequence | 1888–1921 |
| ampicillin resistance (<i>bla</i>) ORF | 1966–2823 |
| pUC origin of replication | 2970–3637 |
| Primer binding site (for 5' end of insert) [5' TCACATGTTCTTCTCGGTTATCC 3'] | 3635–3659 |

FIGURE 6 Features of the phrGFP Mammalian Expression Vector.

OVERVIEW OF CRE-MEDIATED SITE-SPECIFIC RECOMBINATION

Stratagene's Vitality hrGFP reporter vectors utilize *Cre*-mediated site-specific recombination to allow the quick and efficient directional insertion of prefabricated modules. The site-specific recombination catalyzed by *Cre* recombinase is dependent on the presence and orientation of two *LoxP* sites. Each *LoxP* site has two 13-bp inverted repeats, which serve as recognition sites for the recombinase, on opposite sides of an 8-bp nonpalindromic region. The relative orientations of these *LoxP* sites determine the effects of the recombination event. Two sites oriented as tandem repeats results in the excision of the intervening sequence, while two sites in opposite orientation lead to an inversion of the DNA sequence between the two sites. In addition to these intramolecular reactions, site-specific recombination can be intermolecular, leading to the formation of a co-integrate, which is a new molecule formed by the insertion of one DNA molecule into another. A schematic diagram of the recombinant molecule formed by inserting a pExchange module into the pIRES-hrGFP-1a vector is shown in Figure 7. The same type of recombination reaction is possible using any of the three pExchange modules (see Figure 8) and any of the Vitality vectors.




FIGURE 7 Cre-mediated recombination for introduction of new resistance cassettes.

MODULE DESCRIPTIONS

The Vitality hrGFP vectors allow for the introduction of a desired eukaryotic resistance gene into the core expression vectors. Each module, shown in Figure 8, contains a *LoxP* site for insertion of either the hygromycin-, puromycin-, or neomycin-resistance module into the *LoxP* site of any of the Vitality hrGFP vectors. The hygromycin- and puromycin-resistance modules also have a chloramphenicol-resistance marker for selection. Further sequence information for the pExchange modules is available at www.stratagene.com.




FIGURE 8 The three pExchange modules, available separately from Stratagene.

TABLE I

Features of pExchange Modules

| pExchange Module EC-Hyg | | pExchange Module EC-Puro | | pExchange Module EC-Neo | |
|-------------------------|-----------|--------------------------|-----------|-------------------------|-----------|
| Feature | Location | Feature | Location | Feature | Location |
| <i>LoxP</i> sequence | 1–35 | <i>LoxP</i> sequence | 1–35 | <i>LoxP</i> sequence | 1–35 |
| HSV-TK promoter | 91–339 | SV40 polyA signal | 168–137 | SV40 promoter | 94–483 |
| hygromycin | 343–1376 | puromycin | 1031–432 | neomycin | 484–1278 |
| HSV-TK polyA signal | 1375–1937 | SV40 promoter | 1273–1138 | HSV-TK polyA signal | 1266–1724 |
| chloramphenicol | 2895–2236 | chloramphenicol | 2394–1735 | <i>LoxP</i> sequence | 1763–1775 |
| <i>LoxP</i> sequence | 3104–3116 | <i>LoxP</i> sequence | 2603–2615 | | |

PREPARING THE HRGFP VECTORS

Important When cloning into *pIRES-hrGFP-1a* and *pIRES-hrGFP-2a*, note the presence of a stop codon (TAG) in the MCS that is in-frame with the fusion tags (see the circular maps for these vectors in Description of the Vectors). Do not use the *Sac I*, *Sac II*, or *Not I* sites upstream of the stop codon for cloning unless the cloning strategy removes the stop codon by double-digestion using one of these upstream sites plus a site downstream of the stop codon.

- ♦ Ensure that the coding sequence of the insert is in the correct reading frame and that it contains an initiation codon or Kozak sequence.⁴ For gene fusions using the *phrGFP-N1* vector, ensure that the gene of interest is inserted in frame with the *hrGFP* coding sequence. If the insert lacks its own termination codon, termination codons at the 3' end of the MCS may be used (see the MCS sequence in Figure 4). For gene fusions using the *phrGFP-C* vector, ensure that the gene of interest lacks a termination codon, and reads in frame with the *hrGFP* sequence (see the MCS sequence in Figure 3).
- ♦ Stratagene suggests cloning the gene of interest into the appropriate *Vitality* vector before introducing a eukaryotic resistance marker into the vector by *Cre*-mediated recombination (see Figure 8). Introduction of the drug resistance marker to the vector may result in the loss of unique cloning sites in the vector MCS.
- ♦ Stratagene suggests dephosphorylation of the digested vector with *CIAP* prior to ligation with the insert DNA. If more than one restriction enzyme is used, the background can be reduced further by electrophoresing the DNA on an agarose gel and removing the desired plasmid band through electroelution, leaving behind the small fragment that appears between the two restriction enzyme sites.
- ♦ After purification and ethanol precipitation of the DNA, resuspend in a volume of TE buffer (see *Preparation of Media and Reagents*) that will allow the concentration of the plasmid DNA to be the same as the concentration of the insert DNA (~0.1 µg/µl).

LIGATING THE INSERT

The ideal insert-to-vector molar ratio of DNA is variable; however, a reasonable starting point is a 2:1 insert-to-vector ratio. The ratio is calculated using the following equation:

$$X \text{ } \mu\text{g of insert} = \frac{(\text{Number of base pairs of insert}) (0.1 \text{ } \mu\text{g of Vitality hrGFP vector})}{\sim Y \text{ bp of Vitality hrGFP vector}}$$

where X is the quantity of insert (in micrograms) required for a 1:1 insert-to-vector molar ratio. Multiply X by 2 to get the quantity of insert required for a 2:1 ratio. Y is the size (in base pairs) of the Vitality hrGFP vector being used; consult the circular vector maps for sizes.

1. Prepare three control and two experimental 10- μ l ligation reactions by adding the following components to separate sterile 1.5-ml microcentrifuge tubes:

| Ligation reaction components | Control | | | Experimental | |
|-----------------------------------------------------|----------------|----------------|----------------|----------------|----------------|
| | 1 ^a | 2 ^b | 3 ^c | 4 ^d | 5 ^d |
| Prepared vector (0.1 μ g/ μ l) | 1.0 μ l | 1.0 μ l | 0.0 μ l | 1.0 μ l | 1.0 μ l |
| Prepared insert (0.1 μ g/ μ l) | 0.0 μ l | 0.0 μ l | 1.0 μ l | Y μ l | Y μ l |
| rATP [10 mM (pH 7.0)] | 1.0 μ l | 1.0 μ l | 1.0 μ l | 1.0 μ l | 1.0 μ l |
| Ligase buffer (10 \times) ^e | 1.0 μ l | 1.0 μ l | 1.0 μ l | 1.0 μ l | 1.0 μ l |
| T4 DNA ligase (4 U/ μ l) | 0.5 μ l | 0.0 μ l | 0.5 μ l | 0.5 μ l | 0.5 μ l |
| Double-distilled (ddH ₂ O) to 10 μ l | 6.5 μ l | 7.0 μ l | 6.5 μ l | Z μ l | Z μ l |

^a This control tests for the effectiveness of the digestion and the CIAP treatment. Expect a low number of transformant colonies if the digestion and CIAP treatment are effective.

^b This control indicates whether the plasmid is cleaved completely or whether residual uncut plasmid remains. Expect an absence of transformant colonies if the digestion is complete.

^c This control verifies that the insert is not contaminated with the original plasmid. Expect an absence of transformant colonies if the insert is pure.

^d These experimental ligation reactions vary the insert-to-vector ratio. Expect a majority of the transformant colonies to represent recombinants.

^e See *Preparation of Media and Reagents*.

2. Incubate the reactions for 2 hours at room temperature (22°C) or overnight at 4°C. For blunt-end ligation, reduce the rATP to 5 mM and incubate the reactions overnight at 12–14°C.

TRANSFORMATION

Transform competent bacteria with 1–2 μL of the ligation reaction, and plate the transformed bacteria on LB agar plates (see *Preparation of Media and Reagents*) containing the appropriate antibiotic. Refer to references 5 and 6 for bacterial transformation protocols.

Note Competent cells with transformation efficiencies $\geq 5 \times 10^9$ cfu/ μg are available from Stratagene.

VERIFICATION OF INSERT PERCENTAGE, SIZE, AND ORIENTATION

Individual colonies can be examined to determine the percentage of vectors with inserts and the insert size and orientation by PCR directly from the colony or by restriction analysis.

PCR Amplification of DNA from Individual Colonies

The presence and size of a DNA insert in any of the Vitality hrGFP vectors may be determined by PCR amplification of DNA from individual colonies.

1. Prepare a PCR amplification reaction containing the following components:

4.0 μL of 10 \times *Taq* DNA polymerase buffer^{††}
0.4 μL of dNTP mix (25 mM each dNTP)
40.0 ng of the appropriate 5' primer (see table below for sequence)
40.0 ng of the appropriate 3' primer (see table below for sequence)
0.4 μL of 10% (v/v) Tween[®] 20
1.0 U of *Taq* DNA polymerase[†]
dH₂O to a final volume of 40 μL

Primers for hrGFP Vectors

| Vector | Primer | Nucleotide sequence (5' to 3') |
|----------------|------------------------------------|--------------------------------|
| pIRES-hrGFP-1a | Forward (T3 promoter primer) | AATTAACCCTCACTAAAGGG |
| | Reverse | GTCCTTATCATCGTCGTCCTT |
| pIRES-hrGFP-2a | Forward (T3 promoter primer) | AATTAACCCTCACTAAAGGG |
| | Reverse | TAAGCGTAGTCAGGTACATC |
| phrGFP-C | Forward (T3 promoter primer) | AATTAACCCTCACTAAAGGG |
| | Reverse (5' hrGFP Primer) | ACCTTGAAGCTCATGATCTC |
| phrGFP-N1 | Forward (3' hrGFP primer) | CAGCTGACCAGCCTGGGCAAG |
| | Reverse (T7 primer) | TAATACGACTCACTATAGGG |
| phrGFP-1 | MCS 1 Forward (T3 promoter primer) | AATTAACCCTCACTAAAGGG |
| | MCS 1 Reverse (5' hrGFP Primer) | ACCTTGAAGCTCATGATCTC |
| | MCS 2 Forward (3' hrGFP primer) | CAGCTGACCAGCCTGGGCAAG |
| | MCS 2 Reverse (T7 primer) | TAATACGACTCACTATAGGG |
| phrGFP | Forward | TCACATGTTCTTTCCTGCGTTATCC |
| | Reverse (5' hrGFP Primer) | CCAGGTTACCTTGAAGCTCAT |

^{†,††} See *Endnotes*.

2. Stab the transformed colonies with a sterile toothpick and swirl the colony into reaction tubes. Immediately following inoculation into the reaction mixture, remove the toothpick and streak onto antibiotic-containing patch plates for future reference.
3. Gently mix each reaction, overlay each reaction with 30 μ l of mineral oil and perform PCR using the following cycling parameters:

Cycling Conditions for RoboCycler® Temperature Cyclers and Single-Block Thermal Cyclers

| Number of cycles | Temperature | Length of time |
|------------------|-------------|------------------|
| 1 cycle | 94°C | 4 minutes |
| | 50°C | 2 minutes |
| | 72°C | 2 minutes/kb DNA |
| 30 cycles | 94°C | 1 minute |
| | 56°C | 2 minutes/kb DNA |
| | 72°C | 1 minute |
| 1 cycle | 72°C | 5 minutes |

4. Analyze the PCR products for the sizes of the gene inserted into the expression construct using standard 1% (w/v) agarose gel electrophoresis. Additional information can be obtained by further restriction analysis of the PCR products.
5. Following identification of colonies containing the correct insert, return to the patch plates made in step 2 above and pick a portion of each of the positive colonies into 5-ml aliquots of LB broth (see *Preparation of Media and Reagents*) containing the appropriate antibiotic. Grow overnight at 37°C with shaking.
6. The next morning, purify the plasmid DNA from the liquid cultures. You will need to prepare enough plasmid DNA so that *Cre*-mediated recombination may be performed if desired. Miniprep DNA or DNA purified using CsCl gradient centrifugation is of suitable quality.

CRE-MEDIATED RECOMBINATION PROTOCOL

1. To prepare the recombination reaction, add the following components to a 0.5-ml microcentrifuge tube:

X μ l of Vitality hrGFP vector containing the gene of interest
(500 ng)

Y μ l of pExchange module (100 ng)

1 μ l of *Cre* recombinase

1 μ l of 10 \times *Cre* recombinase reaction buffer

Z μ l of distilled water (dH₂O) to a final volume of 10 μ l

2. Incubate reaction for 30 minutes at 37°C. The reaction reaches equilibrium in ≤ 30 minutes.
3. Heat the reaction at 65°C for 20 minutes to denature the *Cre* recombinase. If not denatured, *Cre* recombinase will reduce the transformation efficiency.

Transform competent bacteria with 2.5 μ l of the *Cre*-recombination reaction, and plate the transformed bacteria on LB agar plates containing the appropriate antibiotic. Refer to references 5 and 6 for bacterial transformation protocols.

Note *Competent cells with transformation efficiencies $\geq 5 \times 10^9$ cfu/ μ g are available from Stratagene.*

Confirming Recombination

The presence of recombinant plasmids can be confirmed by analyzing the miniprep DNA of a few colonies by restriction digestion. Ninety-five percent of the colonies are expected to contain recombinant plasmids.

MAMMALIAN CELL TRANSFECTION

When the correct recombinant plasmids are confirmed, prepare enough DNA of appropriate purity for the mammalian cell transfection procedure to be carried out. Protocols for transfection of mammalian cell lines can be found in Sambrook, *et al.* (1989).⁵

The efficiency of transfection will vary depending on the host cell line used. In most cases, mammalian host cell lines transfected with plasmids should show expression of hrGFP 24–72 hours after transfection. Fluorescing cells growing in tissue culture dishes can be observed using an inverted fluorescence microscope. Fluorescence of populations of harvested cells can also be measured using FACS analysis or fluorometer assays. The table below lists excitation and emission spectra for Stratagene's hrGFP as compared to EGFP.

SPECIFICATIONS FOR HRGFP AND EGFP EXCITATION AND EMISSION SPECTRA

| GFP Form ^a | Excitation/Emission Spectra Maxima (nm) |
|-----------------------|-----------------------------------------|
| hrGFP | 500/506 |
| EGFP | 488/509 ^b |

^a Both forms of GFP compared in this table have been codon-optimized for maximum expression in human cells.

^b The emission spectrum for EGFP also shows a shoulder at 540 nm.

Note *Filter sets compatible with the detection of hrGFP and EGFP are sold by Omega Optical, Inc. (Phone: 802 254 2690, or see www.omegafilters.com):*

Exciter filter: XF1073

Emitter filter: XF3084

Beam splitter: XF2010

Microscope cube set with the exciter filter, emitter filter and beam splitter: XF100-2

TROUBLESHOOTING

| Observation | Suggestion |
|-----------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Very few colonies following transformation with Cre-recombined vector | Transformation efficiency is inhibited by Cre recombinase that is not denatured by heat treatment. Heat the recombination reaction at 65°C for 20 minutes to denature the Cre recombinase |
| | Ensure that the Cre-recombinant transformations are plated on the appropriate agar plates |
| | Transformation is inefficient. Check transformation efficiency with a control plasmid |
| | DNA is degraded. Electrophorese the DNA to check its quality |
| Low percentage of colonies with recombinant plasmids | Increase the vector:module ratio in the Cre-mediated recombination protocol |
| Western analysis does not detect fusion protein | Insert is cloned out of frame. Sequence to ensure correct reading frame. Reclone if insert is out of frame |
| | Assay is not sufficiently sensitive or is being performed incorrectly. Use a positive control |

PREPARATION OF MEDIA AND REAGENTS

| | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| LB Agar (per Liter) 10 g of NaCl 10 g of tryptone 5 g of yeast extract 20 g of agar Add deionized H ₂ O to a final volume of 1 liter Adjust pH to 7.0 with 5 N NaOH Autoclave Cool to 55°C. Add appropriate antibiotic. Pour into petri dishes (~25 ml/100-mm plate) | LB Broth (per Liter) 10 g of NaCl 10 g of tryptone 5 g of yeast extract Add deionized H ₂ O to a final volume of 1 liter Adjust pH to 7.0 with 5 N NaOH Autoclave |
| 10× Ligase Buffer 500 mM Tris-HCl (pH 7.5) 70 mM MgCl ₂ 10 mM dithiothreitol (DTT) Note <i>rATP is added separately in the ligation reaction</i> | TE Buffer 10 mM Tris-HCl (pH 7.5) 1 mM EDTA |

REFERENCES

1. Chalfie, M., Tu, Y., Euskirchen, G., Ward, W. W. and Prasher, D. C. (1994) *Science* 263(5148):802–5.
2. Tsien, R. Y. (1998) *Annu Rev Biochem* 67:509–44.
3. Hanazono, Y., Yu, J. M., Dunbar, C. E. and Emmons, R. V. (1997) *Hum Gene Ther* 8(11):1313–9.
4. Kozak, M. (1991) *J Biol Chem* 266(30):19867–70.
5. Sambrook, J., Fritsch, E. F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.
6. Hanahan, D. (1983) *J Mol Biol* 166(4):557–80.

ENDNOTES

- [†] Purchase of these products is accompanied by a license to use them in the Polymerase Chain Reaction (PCR) process in conjunction with an Authorized Thermal Cycler. Stratagene's PCR products are sold under licensing arrangements with Roche Molecular Systems, Inc., F. Hoffmann-La Roche Ltd., and The Perkin-Elmer Corporation.
- ^{††} These products are optimized for use in the Polymerase Chain Reaction (PCR) process covered by patents owned by Hoffmann-La Roche. Purchase of these PCR-related products does not convey any rights under the PCR patents owned by Hoffmann-La Roche. A license to use the PCR process accompanies the purchase of certain reagents from Stratagene when used in conjunction with an Authorized Thermal Cycler.

Vitality[®] and RoboCycler[®] are registered trademarks of Stratagene in the United States.
Falcon[®] is a registered trademark of Becton Dickinson and Company.
FLAG[®] is a registered trademark of Sigma-Aldrich Co.
Tween[®] is a registered trademark of ICI Americas, Inc.

MSDS INFORMATION

The Material Safety Data Sheet (MSDS) information for Stratagene products is provided on Stratagene's website at <http://www.stratagene.com/MSDS/>. Simply enter the catalog number to retrieve any associated MSDS's in a print-ready format. MSDS documents are not included with product shipments.