Code No. 3365

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TaKaRa

pCold™ TF DNA

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Purchaser's Agreement

Customer's order of pCold DNAs will be accepted only when the Purchaser's Agreement is signed by a customer and is attached with an order.

- PRODUCTS are for research or laboratory use only. PURCHASER understands and agrees that PRODUCTS shall not be administered to humans or animals, and shall not be used for pharmaceutical, in vitro diagnostic, or any commercial purposes other than internal research.
- PURCHASER shall not modify pCold Vector DNA sequences located between and including the CspA 3'UTR and CspA Promoter. The adjacent multicloning site (MCS) is exempt from this restriction.
- PURCHASER shall not utilize any partial sequences from PRODUCTS for the purpose of new plasmid construction using the Cold Shock Expression System.
- PURCHASER shall not transfer or sell copies of PRODUCTS, components of PRODUCTS, derivatives of PRODUCTS, and/or products obtained through the use of PRODUCTS (collectively "DERIVATIVES") to any third parties.

Notwithstanding foregoing, PURCHASER may transfer DERIVATIVES solely to a third party which has already been granted a license to use PRODUCTS for research purpose through purchase of PRODUCTS from TAKARA BIO, its subsidiaries or its local distributors, provided that PURCHASER shall enter into a prior separate agreement with TAKARA BIO for the transfer of DERIVATIVES to said third party.

For this Agreement, please contact to TAKARA BIO's subsidiaries or local distributors.

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Digested pCold TF DNA; 100 ng (~ 0.03 pmol)	1 μ l
Insert DNA fragment (0.1-0.3 pmol)	4 μ l
Ligation Mix (from DNA Ligation Kit<Mighty Mix>)	5 μ l
Total	10 μ l

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E. coli JM109 Electro-Cells (Code No.9022)

< 其它 >

IPTG (Isopropyl- β -D-thiogalactopyranoside) (Code No.9030)

Chaperone Plasmid Set (Code No.3340)

pCold™ DNA vector series (Code No.3360 - 3364)

pCold™ ProS2 DNA (Code No.3371)

NOTICE TO PURCHASER: LIMITED LICENSE

[L13a] pCold™ vectors

This product is covered by the claims of U.S. Patent No.5,981,280, 6,686,174 and their foreign counterpart patent claims, assigned to the UMDNJ. This product is covered by the claims of U.S. Patent No.6,479,260, 6,897,042 and their foreign counterpart patent claims.

[M9] pCold™ vectors

This product is covered by the claims of U.S. Patent No. 6,479,260 and its foreign counterpart patent claims.

[M10] pCold™ TF DNA

This product is the subject of the pending U.S. patent application and its foreign counterparts.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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0411-87641685 , 87641686

4006518761 , 4006518769

TAKARA BIO INC.

URL: <http://www.takara.com.cn>

v201211Da