
Replacing the NPD Bead

Instruction Guide

Agilent Technologies

Notices

© Agilent Technologies, Inc. 2004

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G1534-90318

Edition

First Edition, March 2004

Printed in USA

Agilent Technologies, Inc.
2850 Centerville Road
Wilmington, DE 19808-1610 USA

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.


Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Replacing the NPD Bead

6820 and 6890 Gas Chromatographs

This instruction guide reviews the procedure for installing a bead on a nitrogen-phosphorus detector (NPD). Before following this procedure, refer to the safety information on the inside front cover.

CAUTION

The ceramic bead is delicate. Be careful not to break or crack the bead. When you perform maintenance on the NPD, avoid touching the bead with your fingers, and prevent it from coming in contact with other surfaces.

WARNING

Be careful! The oven or detector fittings may be hot enough to cause burns.

Parts list

- 1 NPD bead assembly
- Cap for the bead
- Test chromatogram

Required tools

- T-10 Torx screwdriver

Replacing the NPD Bead

Steps

- 1 Preparing the GC
- 2 Removing the old bead and installing a new bead
- 3 Turning the bead on


Figure 1 NPD on a 6890 GC

Preparing the GC

WARNING

Hazardous voltages are present in the mainframe when the GC power cord is plugged in. Avoid a potentially dangerous shock hazard by unplugging the power cord before removing the side panels.

NOTE

ChemStation Users: The bead voltage cannot be set to zero from the ChemStation. You must use the GC Keypad. In order to do this the keypad must be unlocked and you must close the "Edit GC Parameters" screen of the ChemStation. After you have set the voltage to 0.0 save the method on the ChemStation and then shutdown the instrument session.

- 1 Set the NPD Bead Voltage to 0.0. (Do not set it to "off".)

Temp	100	24
H ₂ flow	3.0	3.0
FRONT DET (NPD)		
Air flow	60.0	60.0
Mkup (He)	5.0	5.0
Adjust offset		Off
Output (Off)		0
Bead voltage		0.0

- 3 Remove the necessary GC covers.
 - **6820:** Remove the detector and right side covers.

- **Early 6890 models:** Raise the gray plastic top cover (with the holes and ventilation slots) to the vertical position. The hinge is a metal bracket attached to the oven top. Pull the clip at its top toward you to release the hinge pin. Push the pin to the left to release the cover. Raise the cover and remove it. Remove the top electronics cover. Push the clips underneath the cover (shown in [Figure 2](#)) toward the center to disengage them.
- **Current 6890 models:** Raise the plastic top cover and remove it. Remove the top electronics cover. Push the clips underneath the cover toward the center to disengage them.


Figure 2 Removing the top electronic cover (6890)

Removing and Installing the Bead

- 1 Disconnect the NPD bead connector from the blue cable by twisting the ring and pulling the ends apart as shown in [Figure 3](#).

Push and turn the ring to move the locking pin into the slot.

Pull apart the connector.


Figure 3 Disconnecting the bead connector

- 2 Remove the three T-10 Torx screws from the NPD bead assembly.
- 3 Lift up and remove the old NPD bead as shown in [Figure 4](#). Avoid bumping the bead against the sides of the collector.


Figure 4 Removing the bead assembly

- 4 Uncap the new bead by pushing the cap off from the cable side (Figure 5). Make sure not to bump the bead on the sides of the cap.


Figure 5 Uncapping the bead

- 5 Mount the new bead assembly on the NPD collector as shown in Figure 6. Be careful not to bump the bead on the sides of the collector. Use the three new Torx T-10 screws supplied. Tighten the first screw only finger-tight; tighten the remaining two screws normally and then completely tighten the first screw.


Figure 6 Mounting the bead assembly

- 6** Carefully bend the bead assembly cable 90 degrees (Figure 7). You should support the cable as shown here.


Figure 7 Bending the bead assembly cable

- 7** Reconnect the bead assembly cable to the NPD cable as shown in Figure 8. Twist the ring to lock the connection.

Push the two parts of the cable together.

Push and twist the ring to slide the locking pin into place and lock.


Figure 8 Connecting the bead assembly cable


- 8** Reinstall the GC covers.
- 9** With all gases on, start heating the detector to 150 °C. Allow it to remain at 150 °C for about 15 minutes, then proceed to heat it to 250 °C and allow it to remain here for about 15 minutes. Finally, heat the detector to 325 °C.
- 10** After temperature equilibration is reached, the bead voltage may be adjusted (0.0 V-2.00 V should be the range to start). If the auto adjust feature is used, set the detector equilibration time to 0.00 minutes (configure detector from keyboard).
- 11** The default offset is 30 picoamps (pA). A lower offset may be used to give longer bead life.

Turning on the bead

- 1** Turn on the GC power. (ChemStation users: restart the instrument session.)
- 2** Check that the NPD bead voltage is set to 0.0.
- 3** Verify that all other GC parameters are set correctly and wait for "Ready".
- 4** Set NPD "Adjust Offset" to desired output current (pA). Normally 30 pA is a good working output current. The auto adjust process will start. Once auto adjust is completed save your method.

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