Digi-Extend® Exercise System: ENGLISH

CE

INDICATIONS

- Develops individual finger strength, flexibility and coordination.
- Also builds hand and forearm strength.

PRECAUTIONS

- If you've had an injury or a hand disorder consult with your healthcare professional prior to beginning the exercises.
- The position of the Digi-Extend® in the hand, the intensity of the exercise, level of resistance and number of repetitions should be approved by a healthcare professional.
- If unusual swelling, skin discoloration or pain occurs, discontinue use and consult a healthcare professional.

GUIDELINE FOR USE

- Exercises should be done at a slow comfortable pace.
- Use the other hand to help stabilize the device during the exercises, as needed.
- Position Digi-Extend® with the hook openings facing away from your hand during use. This will help prevent bands from slipping out.
- For some exercises it is suggested that you utilize a tabletop or rest the Digi-Extend® on your lap to perform exercises effectively and comfortably.
- Suggested exercises follow.

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


Latex-Free Color-coded bands are stretched around the circumference of the Digi-Extend® under the hooks.

FLEXION = Bending Fingers
EXTENSION = Straightening Fingers

BAND TWIST ADJUSTMENT - All the exercises can be modified to help control the amount of band length and tension by twisting the band around the hooks.

ENGLISH

1. Full Finger Extension

Place one finger into each finger guide, under the band until the mid phalanx of each finger is resting on the Digi-Extend®. Extend all five fingers. Relax to starting position and repeat.

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2. Isolated Full Finger Extension

Starting position as in #1. Flex fingers to hold Digi-Extend®. Extend one finger at a time (all three joints), and relax slowly back to the starting position. Repeat with each finger.

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3. Isolated PIP/DIP Extension

Place fingertips in the finger guides. Extend finger, keeping MP joint slightly flexed, extend only the PIP and DIP joints. Relax and let the band flex the joints and extend again. This is a small controlled motion. Repeat with each finger.

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4. Isolated MP Extension

Slide fingers all the way forward so that the proximal phalanges are resting on the device and the PIP and DIP joints are flexed. While keeping the PIP and DIP joints

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Resistive bands develop strength in each finger individually or in entire hand!

This exercise system has been created in association with:

Joanne Kassimir, OTR, CHT
Certified Hand Therapist

KASSIMIR HAND THERAPY, OT, PLLC
2375-B New York Avenue
Huntington Station, NY 11746

Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com

EC

REP

AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)

ver 0821