PRODUCT INFORMATION & MANUAL

Mouse IFN- γ Valukine™ ELISA

VAL607

For the detection of mouse IFN- γ

For research use only. Not for diagnostic or
therapeutic procedures.

Bio-Techne China Co. Ltd
P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Novus kits are guaranteed for 3 months from date of receipt

TABLE OF CONTENTS

I. BACKGROUND.....	2
II. OVERVIEW.....	2
III. ADVANTAGES.....	3
IV. EXPERIMENT	4
V. KIT COMPONENTS AND STORAGE	5
VI. PREPARATION.....	6
VII. ASSAY PROCEDURE.....	7
VIII. RERERENCES	9
IX. TROUBLESHOOTING	10

I. BACKGROUND

Interferon-gamma (IFN- γ , also known as type II interferon) is an important immunoregulatory cytokine that was originally identified through its anti-viral activity (1, 2). It plays key roles in host defense by exerting antiviral, antiproliferative, and immunoregulatory activities (3, 4). On many cell types, IFN- γ induces the production of cytokines and upregulates the expression of various membrane proteins including class I and II MHC antigens, Fc receptors, leukocyte adhesion molecules, and B7 family antigens. IFN- γ is a potent activator of macrophage effector functions. It directs the synthesis, class switching, and secretion of immunoglobulins by B cells. IFN- γ also influences T-helper cell phenotype development by inhibiting Th2 differentiation and stimulating Th1 development (3, 4). IFN- γ plays a central role in the progression of inflammatory diseases such as autoimmunity and atherosclerosis (5, 6).

Biologically active IFN- γ consists of a noncovalently linked homodimer of 20-25 kDa variably glycosylated subunits (7). Mature mouse IFN- γ shares 86% amino acid (aa) sequence identity with rat IFN- γ and 38-44% aa identity with bovine, canine, cotton rat, equine, feline, human, porcine, and rhesus IFN- γ . IFN- γ dimers bind to transmembrane IFN- γ RI (alpha subunits) which then interact with transmembrane IFN- γ RII (beta subunits) to form the functional receptor complex of two α and two β subunits (8, 9). Inclusion of IFN- γ RII in the receptor complex increases the ligand binding affinity as well as the efficiency of signal transduction (9, 10). Whereas the α -chain is expressed constitutively on many cell types, the cellular regulation of the β -chain correlates with an IFN- γ responsive state and is tightly regulated (8).

IFN- γ is produced by a number of cell types under inflammatory conditions, including dendritic epidermal/ $\gamma\delta$ T cells (11), keratinocytes (12), peripheral blood $\gamma\delta$ T cells (13), mast cells (14), neurons (15), CD8⁺ T cells (16), macrophages (17), B cells (18), neutrophils (19), NK cells (20), CD4⁺ T cells (21), and testicular spermatids (22).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. Monoclonal antibody specific for mouse IFN- γ has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse IFN- γ present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked polyclonal antibody specific for mouse IFN- γ is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of mouse IFN- γ bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Diluent 1× and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Two samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision		Inter-assay Precision		
Sample	1	2	1	2	3
n	20	20	20	20	20
Mean (pg/mL)	22.9	345	25.3	128	344
Standard Deviation	1.0	9.9	3.8	11.2	39.7
CV%	4.4	2.9	15.2	8.8	11.5

B. RECOVERY

The recovery of mouse IFN-γ spiked to different levels in cell culture media was evaluated. The recovery ranged from 89-110% with an average of 96%.

C. SENSITIVITY

The minimum detectable dose (MDD) of mouse IFN-γ is typically less than 2.1 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant mouse IFN-γ produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, six cell culture media samples were spiked with high concentrations of mouse IFN-γ and diluted with Diluent 1× to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	105	103 - 117
1:4	109	107 - 112
1:8	105	100 - 110
1:16	101	93- 107

F. SAMPLE VALUES

Cell Culture Supernates - Two spleen organ tissues from a mouse were homogenized and seeded in 100 mL of RPMI1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 g/mL streptomycin sulfate, and 10 µg/mL Con A for 2 days. The cell culture supernate was assayed for mouse IFN-γ and measured 3248 pg/mL. Mouse thymoma cells (EL-4, 2×10⁵ cells/mL) were cultured in DMEM plus 10% Horse Serum and stimulated with 10 µg/mL PHA and 10 ng/mL PMA for 4 days. The cell culture supernate was assayed for mouse IFN-γ and measured 51 pg/mL.

CTLL-2 cells (5×10⁴ cells/mL) were cultured for 3 days in RPMI plus 10% FBS, 2 mM L-glutamine, 10 ng/mL rIL-2, 50 mM 2-mercaptoethanol and stimulated with 2.5 ng/mL LPS, and 100 ng/mL rmGM-CSF. Unstimulated is 95 pg/mL and stimulated is 306 pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant mouse IFN-γ. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant mouse IFN-γ control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant mouse:


IL-10	IFN-γ R1
IFN-γ R2	IFN-kappa

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.042 0.043	0.042	—
9.4	0.077 0.077	0.077	0.035
18.8	0.113 0.112	0.112	0.070
37.5	0.183 0.186	0.184	0.142
75	0.327 0.329	0.328	0.286
150	0.580 0.594	0.587	0.545
300	1.130 1.124	1.127	1.085
600	2.020 2.042	2.031	1.989

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse IFN- γ Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with monoclonal antibody against mouse IFN- γ	1 plate
Mouse IFN- γ Conjugate	21mL/vial of a polyclonal antibody against mouse IFN- γ conjugated to horseradish peroxidase with preservatives	1 vial
Mouse IFN- γ Standard	Recombinant mouse IFN- γ in a buffered protein base with preservatives; lyophilized.	2 vial
Calibrator Diluent RD5P(5 $\times$)	21 mL/vial of a 5 $\times$ concentrated buffered protein base with preservatives. For cell culture supernate samples.	1 vial
Wash Buffer Concentrate (25 $\times$)	21 mL/vial of a 25 $\times$ concentrated solution of buffered surfactant with preservatives.	1 vial
Color Reagent A	12.5 mL/vial of stabilized hydrogen peroxide.	1 vial
Color Reagent B	12.5 mL/vial of stabilized chromogen (tetramethylbenzidine)	1 vial
Stop Solution 2	23 mL/vial of diluted hydrochloric acid solution; 0.075M Citric Acid	1 vial
Plate Covers	adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Diluted Wash Solution	May be stored for up to 1 month at 2-8 °C.*
	Stop Solution 2	
	Dilution 1 $\times$	
	Conjugate	
	Unmixed Substrate A	
	Unmixed Substrate B	
	Standard	Use a new standard for each assay.
	Microplate Wells	Return unused wells to the pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squir bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

The Stop solution 2 provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles.

B. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.


Wash Solution - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Solution Concentrate (25 $\times$) into deionized or distilled water to prepare 500 mL of Wash Buffer.

Substrate Solution - Substrates A and B should be mixed together in equal volumes within 15 minutes of use. Protect from light. 100 μL of the resultant mixture is required per well.

Diluent 1 $\times$ - Add 20 mL of Diluent Concentrate 5 $\times$ into 80 mL of deionized or distilled water to prepare 100 mL of Diluent 1 $\times$.

Mouse IFN- γ Standard - Refer to the vial label for reconstitution volume. Reconstitute the Standard with Calibrator Diluent 1 $\times$. This reconstitution produces a stock solution of 3000 pg/mL. Allow the standard to sit for a minimum of 5 minutes with gentle agitation prior to making dilutions.

Pipette 400 μL of Diluent 1 $\times$ into the 600 pg/mL tube. Pipette 200 μL of Diluent 1 $\times$ into each remaining tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 600 pg/mL standard serves as the high standard. The Diluent 1 $\times$ serves as the zero standard (0 pg/mL).

C. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Substrate Solution should remain colorless until added to the plate. Keep Substrate Solution protected from light. Substrate Solution should change from colorless to gradations of blue.
- Stop solution 2 should be added to the plate in the same order as the Substrate Solution. The color developed in the wells will turn from blue to yellow upon addition of the Stop solution 2. Wells that are green in color indicate that the Stop solution 2 has not mixed thoroughly with the Substrate Solution.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 μ L of Diluent 1 $\times$ to each well.
4. Add 50 μ L of Standard, sample, or control per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature. A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Solution (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 μ L of mouse IFN- γ conjugate to each well. Cover with a new adhesive strip. Incubate for 2 hours at room temperature.
7. Repeat the aspiration/wash as in step 5.
8. Add 100 μ L of Substrate Solution to each well. Incubate for 30 minutes at room temperature. **Protect from light.**
9. Add 100 μ L of Stop Solution 2 to each well. The color in the wells should change from blue to yellow. If the color in the wells is green or if the color change does not appear uniform, gently tap the plate to ensure thorough mixing.

10. Determine the optical density of each well within 30 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
11. **CALCULATION OF RESULTS.** Average the duplicate readings for each standard, control, and sample and subtract the average zero standard optical density.

Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the IFN- γ concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Billiau, A. and P. Matthys (2009) Cytokine Growth Factor Rev. 20:97.
2. Wheelock, E.F. (1965) Science 146:310.
3. Schoenborn, J.R & C.B. Wilson (2007) Adv. Immunol. 96:41.
4. Pestka, S. et al. (2004) Immunol. Rev. 202:8.
5. Kelchtermans, H. et al. (2008) Trends Immunol. 29:479.
6. McLaren, J.E. and D.P. Ramji (2009) Cytokine Growth Factor Rev. 20:125.
7. Gray, P.W. and D.V. Goeddel (1983) Proc. Natl. Acad. Sci. 80:5842.
8. Bach, E.A. et al. (1997) Annu. Rev. Immunol. 15:563.
9. Marsters, S.A. et al. (1995) Proc. Natl. Acad. Sci. 92:5401. 10. Krause, C.D. et al. (2000) J. Biol. Chem. 275:22995.
11. Sugaya, M. et al. (1999) J. Invest. Dermatol. 113:350.
12. Howie, S.E.M. et al. (1996) J. Invest. Dermatol. 106:1218.
13. Battistini, L. et al. (1997) J. Immunol. 159:3723. 14. Gupta, A.A. et al. (1996) J. Immunol. 157:2123.
15. Neumann, H. et al. (1997) J. Exp. Med. 186:2023.
16. Hoiden, I. & G. Moller (1996) Scand. J. Immunol. 44:501. 17. Puddu, P. et al. (1997) J. Immunol. 159:3490.
18. Yoshimoto, T. et al. (1997) Proc. Natl. Acad. Sci. 94:3948.
19. Yeaman, G.R. et al. (1998) J. Immunol. 160:5145.
20. Asea, A. et al. (1996) Clin. Exp. Immunol. 105:376.
21. Briscoe, D.M. et al. (1997) J. Immunol. 159:3247.
22. Dejuco, N. et al. (1995) Endocrinology 136:4925.

IX. TROUBLESHOOTING

Problem	Probable Cause	Suggestion
No signal	Failure to add all components	Prepare a check-list and add the components in the correct order
Low signal	Not enough supernatant per well	Check the protein concentration. Add more supernatant.
High background	Wells are not washed enough.	Wash plates thoroughly after incubation with detecting.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H

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