

AFRL-RH-WP-TR-2020-0063

COG-PACK[™] DASHBOARD USERS MANUAL

Mr. Allen W. Dukes 2Lt. Matt T. Rommel Airman Biosciences Division

Mr. Courtney J. Downs Mr. Scott J. Duberstein Mr. Ethan B. Blackford Ball Aerospace & Technologies

> May 2020 Interim Report

DISTRIBUTION A. Approved for public release:

AIR FORCE RESEARCH LABORATORY 711TH HUMAN PERFORMANCE WING, AIRMAN SYSTEMS DIRECTORATE, WRIGHT-PATTERSON AIR FORCE BASE, OH 45433 AIR FORCE MATERIEL COMMAND UNITED STATES AIR FORCE

NOTICE AND SIGNATURE PAGE

Using Government drawings, specifications, or other data included in this document for any purpose other than Government procurement does not in any way obligate the U.S. Government. The fact that the Government formulated or supplied the drawings, specifications, or other data does not license the holder or any other person or corporation; or convey any rights or permission to manufacture, use, or sell any patented invention that may relate to them.

This report was cleared for public release by 88th Air Base Wing Public Affairs Office and is available to the general public including foreign nationals. Copies may be obtained from the Defense Technical Information Center (DTIC) (http://www.dtic.mil).

AFRL-RH-WP-TR-2020-0063 HAS BEEN REVIEWED AND IS APPROVED FOR PUBLICATION.

DOUGLAS L. FRANCK, DR-III, DAF Chief, Performance Optimization Branch Airman Biosciences Division MCKINLEY.RICHA PD 202009710 -Date: 2020.09.03 10:08:15 -04'00'

R. ANDY MCKINLEY, PhD Core Research Area Lead Performance Optimization Branch Airman Biosciences Division

This report is published in the interest of scientific and technical information exchange, and its publication does not constitute the Government's approval or disapproval of its ideas or findings.

REPORT DOCUMENTATION PAGE					Form Approved OMB No. 0704-0188	
The public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Artington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number. PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS .						
1. REPORT DAT 04 08 20	E (DD-MM-YY)	2. REP	ORT TYPE erim			3. DATES COVERED (From - To) January 2019 to August 2020
4. TITLE AND SU	JBTITLE					5a. CONTRACT NUMBER
COG-PACK TM Dashboard Users Manual					FA8650-16-C-6610	
coo men	Dubiloouru	e seis manaai				5b. GRANT NUMBER
5c. PROGRAM ELEMENT NUM 62202F						5c. PROGRAM ELEMENT NUMBER 62202F
6. AUTHOR(S) 5d. PROJECT NUMBER						5d. PROJECT NUMBER
						5329
COURTNEY J. DOWNS+, ALLEN W. DUKES*, ETHAN B. BLACKFORD+,						5e. TASK NUMBER
SCOTT J. DUBERSTEIN+, MATT T. ROMMEL*						08
						5f. WORK UNIT NUMBER
						H0KG (53290819)
7. PERFORMING	ORGANIZATION	NAME(S) AND ADDR	RESS(ES)			8. PERFORMING ORGANIZATION
+Ball Aerosp	ace & Technol	ogies				REPORT NUMBER
2675 Presider	ntial Drive, Fai	rborn OH 4532	4			
9. SPONSORING	MONITORING AG	ENCY NAME(S) ANI	D ADDRESS(ES)			10. SPONSORING/MONITORING
*Air Force M	[ateriel Comm	and	Air Force M	ateriel Comma	nd**	AGENCY ACRONYM(S)
Air Force Re	search Laborat	orv	Air Force Re	esearch Laborat	orv	711HPW/RHBC
711 Human P	erformance W	ing	711 Human	Performance W	'ing	11. SPONSORING/MONITORING
Airman Syste	ens Directorate		USAF Schoo	of Aerospace	ing	AGENCY REPORT NUMBER(S)
Airman Bios	ciences Divisio	'n	Medicine	or or recospuee	, ,	
Performance	Ontimization I	Branch	Wright-Patte	erson AFR OH		AFRL-RH-WP-TR-2020-0063
Cognitive Ne	uroscience Sec	rtion	45433	13011711 D, 011		
Wright-Patter	rson AFR OH	45433	-5-55			
12. DISTRIBUTIO	wright-Patterson AFB, OH 45455 12. DISTRIBUTION/AVAILABILITY STATEMENT					
DISTRIBUT	ION A. Appro	oved for public r	elease,			
13. SUPPLEMEN		1.00.4				
88ABW-2020-2718, cleared 28 August 2020						
14. ABSTRACT The design of the COG Pack Dashboard allows a user to set up, configure, and manage multiple Sensors by						
providing a standard look and feel for each Sensor. Generally, when manufacturers develop hardware sensors, they provide a						
custom application designed to only work with their particular device. From an integration standpoint, this provides a high						
barrier of entry for researchers and users when an experiment requires multiple disparate sensors. Each independent						
approvides unque functionality for actions such as saving files, managing sensor state connectivity, and observing the data-streams in real-time. One of the objectives of the COG Pack Dashboard is to simplify these steps and						
consolidate them into a single interface. Additionally by leveraging a database back-end for data storage. COG Pack						
automatically manages data timestamps and organization. By not requiring a user to execute multiple mouse actions and						
keyboard inputs for file management. COG Pack streamlines this process and reduces the potential for 'fat-fingered' input						
errors.						
The goal of this document is to provide a user or researcher with the instructions they need to work with the COG Pack						
Dashboard. This document contains descriptions of the processes and procedures required to interact with the supported						
Sensors and Si	gnals.					
15. SUBJECT TE	ERMS COG Pack	™, Dashboard, U	User Interface, U.	II, Manual		
16. SECURITY C	LASSIFICATION O	F:	17. LIMITATION	18. NUMBER OF	19a. NA	ME OF RESPONSIBLE PERSON (Monitor)
a. REPORT	b. ABSTRACT	c. THIS PAGE	ABSTRACT:	18	1	Armando Soto
U	U		SAR	10	19b. ТЕІ N	LEPHONE NUMBER (Include Area Code) /A



COG PACKTM DASHBOARD

COURTNEY J. DOWNS, ALLEN W. DUKES, ETHAN B. BLACKFORD, SCOTT J. DUBERSTEIN, MATT T. ROMMEL

User's Manual

TABLE OF CONTENTS

1.0	OVERVIEW	1
2.0	STARTING THE COG PACK DASHBOARD	2
3.0 3.1 3.2 3.3 3.4	ASSIGNING SENSORS TO PARTICIPANTS HSI PHYSIO Configuration Profusa Lumee Configuration Intelligent Optical Systems (IOS) MASES Configuration Smart Eye Pro Configuration	2 3 4 5
4.0 4.1 4.2	ATTACHING ADDITIONAL SIGNAL PROCESSING STEPS Attaching a Threshold Calculator Attaching a "Core Temperature Estimator"	8 8 . 10
5.0	GRAPHS AND VISUALIZATIONS	. 11
6.0	LOADING AND SAVING CONFIGURATIONS	. 12
7.0	CONCLUSION	13

LIST OF FIGURES

2
2
3
4
5
7
8
9
9
0.
.1
.1
.2
.3
-

1.0 OVERVIEW

The design of the COG Pack Dashboard allows a user to set up, configure, and manage multiple Sensors by providing a standard look and feel for each Sensor. Generally, when manufacturers develop hardware sensors, they provide a custom application designed to only work with their particular device. From an integration standpoint, this provides a high barrier of entry for researchers and users when an experiment requires multiple disparate sensors. Each independent application often provides unique functionality for actions such as saving files, managing sensor state connectivity, and observing the data-streams in real-time. One of the objectives of the COG Pack Dashboard is to simplify these steps and consolidate them into a single interface. Additionally, by leveraging a database back-end for data storage, COG Pack automatically manages data timestamps and organization. By not requiring a user to execute multiple mouse actions and keyboard inputs for file management, COG Pack streamlines this process and reduces the potential for 'fat-fingered' input errors.

The goal of this document is to provide a user or researcher with the instructions they need to work with the COG Pack Dashboard. This document contains descriptions of the processes and procedures required to interact with the supported Sensors and Signals.

2.0 STARTING THE COG PACK DASHBOARD

The COG Pack system requires two main software components for execution: the back-end web service and the front-end Dashboard UI. These are the necessary components required, regardless of the configuration. In the most commonly deployed configuration, a single icon on the desktop loads both components. (Figure 1.)



Figure 1 - COG Pack Desktop Shortcut

3.0 ASSIGNING SENSORS TO PARTICIPANTS

After loading the COG Pack processes, the initialization screen is presented. (Figure 2.) This screen presents the first step in the process required for set up and configuration, allowing a user to associate particular Sensors and Signals with a participant.

Add Participant]
Study Configuration	Add Sensor(s)
Admin Settings 0000001	SmartSye PHYSIO ThresholdCalculator CoreTempEstimator MASES RegionCalculator

Figure 2 - COG Pack Initialization Screen

Initialize

To begin the association, complete the following steps:

- 1. Add a participant by clicking the "Add Participant" button on the left.
 - a. By default, the system generates a Participant ID
- 2. Click the generated Participant ID in the list.

- 3. Select the desired type of Sensor to assign to this participant by selecting the name of the Sensor from the array of available Sensors found in the "Add Sensor(s)" section.
- 4. Enter the values and settings required for each Sensor.

Proceed to the following sections for examples showing how to configure the individual Sensors.

3.1 HSI PHYSIO Configuration

Figure 3 displays the card for the configuration of an HSI PHYSIO Sensor. To assign an HSI PHYSIO Sensor to a participant, first label it with a Sensor ID by entering a value into the Sensor ID text field. Then select all of the Signals to capture. These selected Signals are then available to the Signal processing pipeline for viewing, logging, and further analysis. For this Sensor, COG Pack has provided the ability to toggle between live and simulated data. To configure this option, toggle the button labeled "Emulated" to tell the system if this device is live or simulated. Finally, the user has the option to attach additional Signal processing algorithm(s). (See also <u>Attaching Additional Signal Processing Steps</u>)

	COG Pack
Add Participant]
study Configuration	Add Sensor(s)
Admin Settings	SmartEye PHYSIO ThresholdCalculator CoreTempEstimator MASES RegionCalculator
0000001	Lumee
	Physio
	Enter Sensor ID: Emulated Assign Signal Processing:
	Mario Signal Processing
	Select Desired Signals:
	Heart Rate Battery 181 Respiration Rate SpO2
	Movement Chest Temperature Underarm Temperature Skin Temperature
	Initialize

Figure 3 - Configuring an HSI PHYSIO Sensor

3.2 Profusa Lumee Configuration

To add a Profusa Lumee Sensor to a participant, locate the 4 character alpha-numeric ID on the Sensor, and enter that value into the Sensor ID text field. Then select all of the signals to capture. Finally, attach any additional Signal processing algorithms. (See also <u>Attaching</u>

Additional Signal Processing Steps) Figure 4 demonstrates the configuration of a Profusa Lumee Sensor.

	~	
Add Participant		
Study Configuration	Add Sensor(s)	
Admin Settings	SmartEve PHYSIO Threshold	alculator CoreTemoEstimator MASES RegionCalculator
0000001		
0000002	Lumee	
0000003	Lumas	
	Lumee	
	Enter Sensor ID:	Assign Signal Processing:
	E858	Signal Processing •
	Select Desired Signals: Character Oxygen Concentration Diagnostic fits O Battery D Basel	✓ Accelerometer ✓ Motion Classification ne ✓ Intensity ✓ Gain ✓ Lib Brightness

Figure 4 - Configuring a Profusa Lumee Sensor

3.3 Intelligent Optical Systems (IOS) MASES Configuration

Adding an IOS MASES Sensor follows the same process as described for other Sensors. The IOS MASES Sensor requires only a few configuration options. First, a Sensor ID to uniquely describe this instance, then the IP Address of the interface connected to the IOS MASES Sensor. Additionally, specify the port number that the IOS MASES device uses to send data. COG Pack also supports emulation for this Sensor and can be toggled with the "Emulated" button. Finally, select all of the Signals to capture visually by clicking their respective check boxes. Figure 5 displays the configuration of an IOS MASES Sensor.

Add Participant]			
itudy Configuration	Add Sensor(s)			
dmin Settings	SmartEye PHYSIO	ThresholdCalculator CoreTempEstin	nator MASES RegionCalculator	
000001				
000002				
	Mases Enter Sensor ID:	Enter Sensor Address:	Assign Signal Processing:	
	MASES1	192.168.2.25:4023	Signal Processing	

Figure 5 - Configuring an IOS MASES Sensor

3.4 Smart Eye Pro Configuration

Currently, Smart Eye Pro is the most complicated Sensor to configure in COG Pack. To start, fill out all of the listed fields below with the settings formatted as described in Table 1.

Setting Name	Description	Data Type
Sensor ID	Unique Identifier for this Smart Eye Instance	String
IP Address	IP Address of the COG Pack system connected to Smart Eye	String
Machine Name	Name of the COG Pack system	String
Port Number	Port Number used by Smart Eye to connect to COG Pack	Integer
World Model File	Full path to *.sew file used by Smart Eye	String
Minimum Fixation Dwell Time	For AdvancedFixations calculations, define the minimum number of milliseconds required for a Fixation	Double
Vergence Angle	For AdvancedFixations calculations, define the Visual Angle, in degrees, describing the maximum dispersion of points possible for a Fixation	Double
Monitor Size in Pixels	A comma-separated list of 4 Doubles that define the pixel space of a single monitor. DEPRECATED – replaced by values in the World Model File	DEPRECATED
Monitor Size in Meters	A comma-separated list of 4 Doubles that define the physical dimensions of a single monitor. DEPRECATED – replaced by values in the World Model File	DEPRECATED
Sample Rate	The frequency at which Smart Eye is configured. Most commonly 60 or 120 Hz.	Integer
Playback File	Full path to a *.sme file recorded previously by COG Pack software to replay previously recorded binary data. (Leave empty for live data capture)	String

Table 1. Smart Eye Settings, Descriptions, and Data Types

These are the settings required for the configuration of the interface between COG Pack and Smart Eye. For additional settings of the visual elements of COG Pack, see the report titled "COG PACK Eye/Head/Gaze/Heatmap 3D Display Documentation". Figure 6 shows the configuration of a Smart Eye Sensor.

	Kan tang tang tang tang tang tang tang ta	Pack
Add Participant		
Study Configuration	Add Sensor(s)	
Admin Settings	SmartEve PHYSIO ThresholdCalculator	CoreTempEstimator MASES RegionCalculator
0000001		
0000002	Lumee	
0000003	Constant	
0000004	SmartEye	
	Binary Logging Binary Logging Di Choose File No	irectory file chosen
	Sensor ID	Minimum Eixation Dwell Time
	RIHPWPWS-4ZV3Z	100
	IP Address	Vergenance Angle
	10.0.20.102	2
	Machine Name	Monitor Size In Pixels
	RIHPWPWS4ZV3Z	0,0,3840,2160
	Port Number	Monitor Size In Meters
	35466	0,0,0 59,0 33
	World Model File	Sample Rate
	munications\F35 UDL AFRL.sew	Playback Eile
		Playback File Directory:
	Assign Signal Processing: Signal Processing	
	Select All Signals System Performance Saze Intersecti	on 🗸 Blink 🗸 Perclose
	Head Position Flation Flation Flation Gaze Direction Quality Raw Gaze Flatered Gaze Flatered Gaze Flatered All World Intersection Flatered All World Intersection Flatered All World Intersection Gaze Structure Head Head All World Intersection Left Closest World Intersection Flatered All World Intersection Flatered Kight Closest World Intersection Flatered Kight Closest World Intersection	Sye Orientation Eye Position Advanced Fixations Raw Estimated Gaze Filtered Estimated Gaze Filtered Estimated Gaze Filtered Closest World Intersection ated Closest World Intersection Closest World Intersection Closest World Intersection Usoft Intersection World Intersection World Intersection Wire Closest World Intersection Wire Closest World Intersection World Intersection Wire Closest World Intersection
	Estimated Left Closest World Intersection Estimated Right Closest World Intersection	Estimated Left All World Intersections Estimated Right All World Intersections

Figure 6 - Configuring a Smart Eye Sensor

4.0 ATTACHING ADDITIONAL SIGNAL PROCESSING STEPS

Each of the above cards for a Sensor's configuration has an option to add additional signal processing algorithm(s). To view all of the available algorithms for the selected Sensor, select the drop-down directly under the title **Assign Signal Processing (Figure 7)**.

Add Participant	
Study Configuration Admin Settings	Add Sensor(s) MASES SmartEye ThresholdCalculator RegionCalculator CoreTempEstimator PHYSIO
0000001	
	Physio Enter Sensor ID: Assign Signal Processing:
	Physio Enter Sensor ID: Emulated Mario Signal Processing
	Mario Assign Signal Processing: Select Desired Signals: ThresholdCalculator

Figure 7 – Assigning Signal Processing

4.1 Attaching a Threshold Calculator

To attach "Threshold Calculator" to a Sensor you must perform the following actions:

- 1. Select the Signal to Threshold for the current Sensor.
- 2. Select the "Threshold Calculator" option for Signal processing as shown in Figure 8.
- 3. Fill in the fields with the desired range labels and values.
- 4. Click the button labeled "Attach".

Add Participant]
udy Configuration Imin Settings 00001	Add Sensor(s) SmartSye PHYSIO ThresholdCalculator CoreTempEstimator MASES RegionCalculator Lumee
00002	
00004	Physio
00005	Enter Sensor ID: Emulated Assign Signal Processing: Mario ThresholdCalculator × • CoreTemp5timator
	Select Desired Signals: Heart Rate Battery BI Respiration Rate Sp02 Movement Chest Temperature Underarm Temperature Skin Temperature
	Thresholding
	Range Labels Range Values Attach New Label +

Figure 8 – Assigning a Threshold Calculator

For example, to set the following values for the Signal "Respiration Rate," enter the values in the same order as shown in **Figure 9.**

Low : < 60	Medium : >= 60 and < 120	High : >= 120
-------------------	---------------------------------	----------------------

Range Labels	Range Values			Attach	
Low	60				
Medium	60	120]		
High	120 +]		

Figure 9 – Assigning Threshold Ranges to a Signal

4.2 Attaching a "Core Temperature Estimator"

The process for attaching a "Core Temperature Estimator" follows a similar method:

- 1. Select the Signal (Heart Rate) used to derive the Core Temperature Estimate (Figure 10.)
- 2. Select the "Core Temperature Estimator" option for Signal processing.
- 3. Fill in the value for "Starting Core Temperature" as seen in Figure 11.
- 4. Click the button labeled "Attach".

	COG Pack
Add Participant Study Configuration Admin Settings 0000001	Add Sensor(s) PHYSIO SmartEye ThresholdCalculator CoreTempEstimator MASES
	Enter Sensor ID: Emulated Assign Signal Processing: Mario CoreTempEstimator X +
	Select Desired Signals: Heart Rate Battery IBI Respiration Rate SpO2 Movement Chest Temperature Underarm Temperature Skin Temperature
	Core Temperature Estimation Starting Core Temperature
	Initialize

Figure 10 – Assigning a Core Temperature Estimator

Core Temperature Estimation				
Starting Core Temperature	37.1		Attach	
°F °C				

Figure 11 – Configuring the Starting Temperature for Core Temperature Estimation

5.0 GRAPHS AND VISUALIZATIONS

When Sensor and Signal configuration is complete, click the button at the bottom of the screen labeled **Initialize**. This action will start the communcation process between the frontend UI and back-end **COG_Pack.Communcations.exe** processes. All of the settings and configurations will be sent to the backend, and the Sensors will begin connecting. The UI will transition to a screen that renders all of the graphs configured in the previous screen. **Figure 12** is an example of a single participant and the output of the configured Sensors and Signals. This example contains the output of various Threshold Signals derived from the Signals produced by a single HSI PHYSIO Sensor and 2 separate Profusa Lumee Sensors.



Figure 12 – Graphical Display of Signals

The different colored sections on these graphs and gauges represent the configured values for each Signal's respective "Threshold Calculator" range values. In the COG Pack Threshold Signal definition, these values can change over time, and the graphs and gauges react accordingly.

6.0 LOADING AND SAVING CONFIGURATIONS

Another very important feature of this application is its ability to save the current setup and suit of sensors that you have configured for a group of participants into a single .JSON file. This file captures the settings and setup that you have entered into the application and can represent configuration for a specific trial.

For example, if a researcher wishes to use COG Pack for a study and it always requires the same 5 participants with the same set of sensors then you can just set it up in the application one time and hit save. Then, you will be able to load that saved file into the application later and it will restore the setup without you having to go through and reenter all of that data every time.

Therefore, on the main screen of the COG Pack Dashboard is an option on the left labeled **Save/Load (Figure 13).**

COG Pack					
Add Participant Study Configuration Save/Load 0000001 0000002 0000003	Upload Sensor Configuration File Choose File No file chosen	Save Current Sensor Configurations Save Filename: Save			

Figure 13 – Dashboar	d Configuration	Save and Recall
----------------------	-----------------	-----------------

Initialize

Loading a Configuration

The process for loading an existing configuration for the Dashboard is the following:

- 1. Navigate to the Save/Load tab on the left
- 2. Click the button labeled "Choose File"

- 3. Navigate to the .JSON file that you've previously saved
- 4. Click the open button as seen in **Figure 14** below.
- 5. Click the Initialize button.

After following the above steps and clicking the initialize button at the bottom of the screen you will see the screen displaying all of the visualizations and graphs associated with the sensors that were saved to the save file.



Figure 14 – Loading a Dashboard Configuration

Saving a new Configuration

To save your existing dashboard configuration, do the following:

- 1. Configure the sensors for your participants
- 2. Navigate to the Save/Load page in the application
- 3. Enter a name for the configuration file that will be generated
- 4. Click the save button

7.0 CONCLUSION

The COG Pack dashboard is a centralized place for multiple sensor configurations and visualizations. Our unique approach of using "Reflection Based" user interface generation allows this application to be very versatile when it comes to adding even more disparate sensors to its existing library of integrated sensors. The choice that we made to use Electron along with the React.js web framework allows this application to also be able to run on many different operating systems. We are constantly working on improving the design and user experience of the Dashboard to meet the needs and suggestions of our user base.