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AUTOMATIC AUTOCLAVES FOR STEAM STERILIZATION

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Índice

1 Safety	30
1.1 Safety Icons	30
1.2 Risks to which the operator is subjected to:	30
1.3 Qualified users	30
2 General Information for Users.	31
2.2 Autoclave reception	31
2.3 Documentation	31
2.4 Other information	31
2.5 Warranty period	31
3 Technical specifications	32
4 Contents list and accessories	33
4.1 Standard contents	33
4.2 Accesories	33
5 Introduction	34
6 Equipment description	35
6.1 Design features decisive for the equipment life	35
6.3 Operating mode and range	36
6.4 Pre-established cycles	36
6.5 Test programs	36
6.6 Control panel elements	37
6.7 Main elements location	37
7 Installation	38
7.1 Location	38
7.2 Installation without water source nor wastepipe.	38
7.3 Installation with external water source inlet and wastepipe	38
7.4 Drain system selection:	39
7.5 Electrical source connection:	39
7.6 Installation qualification.	39
8 Operation	40
8.1 Preliminary operations for switch on	40
8.1.1 Opening and closing the loading and unloading door:	40
8.1.2 Cautions during usage:	40
9 Configuration, parameters and cycles selection menus	41
9.1 Cycles selection	42
9.2 Program modification	42
9.3 Language modification	43
9.4 Change °C to °F configuration	43
9.5 Printer configuration	44
9.6 External water source configuration	44
9.7 Temperature exterior sensor configuration.	45
9.8 Date and time configuration	45
9.9 Vacuum test	46
9.10 Bowie-Dick test	46
10 What to do when getting an error during the cycle operation.....	47
11 Recommendations for perfect sterilisation	48
11.1 Instruments sterilisation in bulk:	48
11.2 Instrument sterilisation in bags:	48
11.3 Tubes sterilisation:	48
11.4 Vessels sterilisation:	48
11.5 Liquids sterilisation:	48
12 Safety devices	49
12.1 Safety pressurestat: Adjusted to 2.35bar	49
12.2 Safety valve: Adjusted to 2.5 bar	49
12.3 Safety thermostat: No adjust possibility.	49
12.4 Door blockage with temperature > 98°C	49
13 Calibration	50
14 Maintenance	51
14.1 Basic maintenance operations	51
14.2 Maintenance operations on demand by U.E. legislation	51
14.2.1 Annual reviews	51
14.2.2 Periodic hydraulic tests:	51
15 Spare parts	52
16 Schematic	53

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
16.2 Esquema hidráulico

1 Safety

The equipment includes the appropriate safety features.

Risk situations that must be respected are indicated in this manual.

1.1 Safety Icons

They identify risk situations and safety measures that should be taken.

The icons make reference to the paragraph marked with a grey line.

Danger risk

Indicates danger risk.

Please, respect the indicated instructions to carry out the described operation.

Electrical risk

Indicates risk of electric shock when accessing zones shown by this icon or when taking actions indicated in this manual with this icon added.

Please, respect the indicated instructions to carry out the described operation.

Burn risks due to high temperature area's contact.

Temperature in the area indicated by this icon can exceed 60°C. Use please the thermal protective gloves to carry out the described operation.

Please, respect the indicated instructions to carry out the described operation.

Important information

- Important information to obtain good results or an optimum equipment performance.
- Important information to extend the equipment life or to avoid degradation of some of its components.

1.2 Risks to which the operator is subjected to:

- Possibility of touching hot surfaces over 60°C.
- Possibility of steams scape.
- Electrical risk.

1.3 Qualified users

This equipment can only be used by qualified staff.

This equipment can only be used by users who have read and understood these instructions or those that have been qualified for the use of this equipment.

2 General Information for Users.

2.1 Legal requirements for instalation, use and maintenance of this autoclave:

Valid reference regulation over pressure equipments:

- Directive 97/23/CE
- Pressure equipments regulation RD 1244/1979 (Spain)

Instalation requirements are established in section 7 of this manual.

Maintenance requirements are established in section 14 of this manual.

2.2 Autoclave reception

Handle the equipment in a suitable way depending on its weight and consider it as FRAGILE.

Unpack and check if the contents match up with the «packing list». If you see any damage or something missing, please notify it to your distributor as soon as possible.

Keep the original package for some days. Then, throw away in the right way, by separating its components into carton, wood and plastic.

2.3 Documentation

The following documents are supplied with the autoclave :

- User's Manual:

Do not install nor use this equipment before reading carefully this user's manual.

This instructions are an important part of this equipment and must be available to all the users.

Any doubt should be asked to the autoclave distributor.

- CE Declaration of Conformity:

This is a necessary document for authenticating the autoclave instalation.

- Operation Test:

A registry that proves the equipment correct operation is given together with the autoclave.

2.4 Other information

Autoclaves sent to J.P. SELECTA, S.A. for repairing or maintenance should be disinfected.

Any modification, elimination or lack of maintenance of any autoclave device is forbidden.

The use of this equipment in explosive environment or with steamed substances or those that can be explosive or inflammable it's forbidden.

The guaranteed protection and the functionality of the equipment can be committed if the autoclave is used in a way not specified by J.P. SELECTA, S.A..



2.5 Warranty period

Autoclaves guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA,S.A.

Any manipulation of the apparatus by personnel not authorized by J.P. SELECTA,S.A. automatically cancels the guarantee.

3 Technical specifications

	Autoclave volume			
	30	50	75	140 Liters
Network voltage:	I 230V		III 400V	
Electrical power: (W)	2800	3000	3600	7500
Max. intensity: (A)	12.5	14	16	16
Fuses: (fast. type H, dimensions: 10x38mm)	20	20	20	32 (x3)
Network connection cable:	Fase, Neutral, Earth (section 1mm ²)		Three-pin + Earth	
Weight (Net Kg):	90	108	136	250
Useful measures (cm Ø x depth)	Ø30x40	Ø30x70	Ø40x60	Ø50x70
Ext. measures: (cm Depth x Width x Height)	62x48x117	62x48x117	72x58x107	95x80x118
Volume: (Litre)	30	50	75	140
Maximum load: (Kg) (Metal)	8	12	20	35
Weight per support (Kg/cm ²) Vacia	2	2.5	3	2
Weight per support (Kg/cm ²) Load + Water	7.5	2.5	3	2
Total space + doors	Add 1 m to the height			
Max. pressure of ext. water (bar)	2	2	2	2
Min. pressure of ext. water (bar)	0.5	0.5	0.5	0.5
Water volume per cycle (Litres)	1	1	1	0
Water quality	Demineralized between 15 to 200 microS/cm; pH between 5 & 7			
Thermal energy to environment: (W/h)	525	600	650	800
Water tank volume: (Litres)	10	10	10	50
External water inlet record	3/4"	3/4"	3/4"	3/4"
Air inlet filter	Porosity 0.20 micres			
Acoustic level:	Less than 70 dBA			
Operation system:	Continuous with 20min intervals between cycles.			
Environmental conditions:	Inner use. Altitude till 2000m Temperature between 5°C & 40°C. Max. relative humidity 80% per temperatures till 31°C, reducing in a linear way till 50% of relative humidity to 40°C.			
Overpressure level:	Categorie II			
Pollution degree:	2			
Vapour contact material	Stainless steel, copper, teflon, brass			
Frame material:	Steel with cover			
Vessel material:	Stainless steel			
Vapour generator material:	Stainless steel			
Pipes material:	Copper, brass			
RS-232 outlet:	9 pin male connector. Pin 2 & 3 crossed, pin 5 common. Speed 9600, no parity, data 8 bits, stop 1 bit.			

4 Contents list and accessories

4.1 Standard contents

The standard equipment consists of the following items:

Descripción	Description	Código/Code		
Autester ST DRY PV-II	Autester ST DRY PV-II	4002514	4002515	4002425
Bandeja	Shelf		5601	5600
Racord vaciado	Empyng racord		34100	
Manual de instrucciones	Instruction manual		80206	

4.2 Accesories

- Thermal printer with continuos paper roll, also with indication of temperature, pressure, time, mode and condition. Code 4002417
- Temperature sensor for internal temperature mesurements within the chamber, with two applications (Code 4001219):
 - Measured temperature values for a register or a thermometer
 - Used in liquid programs or in packed items inside sterilization.
- Wired baskets made of AISI 316 stainless steel.

Código cestillo Basket code	Medidas (cm) Dimensions (cm)		Para autoclave cod. To autoclave code	Capacidad cestillos Number of baskets
	Ø	Altura/Height		
1000495	28	20	4002514	2
1000495	28	20	4002515	3
1000496	38	28	4002425	2

- Basket dividers made of AISI 316 stainless steel.
 - Code: 1001217 Divider of 4 compartments for basket 1000495
 - Code: 1001218 Divider of 6 compartments for basket 1000496
- Sterilisation drums

Código Code	Ø	Alto	Para autoclave cod. To autoclave code	Capacidad tambores Drums capacity
	(cm)			
1002418	28	18	4002514	2
1002418	28	18	4002515	3
1002419	35	27	4002425	2

5 Introduction

The Autester ST DRY PV II autoclaves are versatile and suitable for a wide range of applications in the field, health-care, industrial processes and quality control. These autoclaves allow solid sterilisation processes, with or without packaging, liquids, culture media and prions. They also allow disinfection and drying cycles.

5.2 Sterilisation

Sterilisation is the destruction or elimination of all forms of microorganisms including spores present in inanimate objects.

5.3 Disinfection

Disinfection is the process of destruction of infectious agents. Normally it's made at a lower temperature than sterilisation. Disinfection only achieves the elimination of some vegetative forms.

5.4 Culture media

Culture media is a nutritional material in which microorganisms can be recovered, multiplied and isolated, as well as used in susceptibility tests. They are normally shown as dry powder or granular but also can be shown hydrated and prepared. They must not be used before being sterilised.

5.5 Prions

Prions are transmissible acellular pathogen particles. They are characterized by producing diseases which affect the central nervous system (SNC). Prions are not living beings but only protein particles without nucleic acid. Its pathogen action consists in being a modified form of an existing natural protein in an organism which leads original proteins to have prion forms when coming into contact. This prion form used to be an abnormal malfunctioning form. And this happens as a chain reaction which destroys operativity of all proteins sensible to prions.

5.6 Purge (extracting the vessel air)

Purge is the cycle which extracts the air from the interior of the vessel to achieve saturated steam vapour.

5.7 Saturated vapour

Water vapour at a certain temperature corresponding to the boiling point of the original liquid.

5.8 Atmospheric purge

In the atmospheric purge, the air is expelled by vapour in a gravimetric way for a determined period of time.

5.9 Pulsed vacuum purge

Pulsed vacuum purge combines the vessel warming with the air extraction by means of a vacuum pump. During the purge cycle, the vessel can reach twice or three times a high vacuum level (depending on the configuration). Pulse vacuum purge is reliable than atmospheric purge.

5.10 Bowie Dick test

The Bowie Dick test serves as a verification of the sterilisation by water vapour applied to wrapped products and porous loads.

5.11 Vacuum test

This is a test assigned to verify the vacuum pump function and also the tightness of the vessel and the hydraulic circuit for a determined period of time.

6 Equipment description

The Autester ST DRY PV II 30, 50 and 75 litres belong to a series of saturated vapour autoclaves which are distinguished by the following features:

- External main chamber steam generator.
- Not suitable for fixing a sterilisation unit.
- Cycles for sterilisation of:
 - Metal load unpacked: Clamps, Scalpels, ...
 - Metal load packed and double packed.
 - Rubber load: Tubes, pipes, catheter, ...
 - Liquid containers load.
- Used water wastepipe:
 - Towards the tank for its reuse.
 - Towards the instalation wastepipe.
- Prepared to be used by two sensors. The instalation of the second sensor (code 4001219) must be fit in the factory.
- This function is useful for liquid sterilisation.
- Printer data registry or RS-232 external interface for temperature, pressure, time and phase readings.
- Bowie-Dick and vacuum tests, for operation and validation check.
- Operation system.
- Defective cycle indication.

6.1 Design features decisive for the equipment life

The maximum service temperature for this equipment is lower than the one for which used materials show damage for slow fluidity.

- Fatigue: This equipment has been designed to resist the maximum admissible pressure PS = 2,5 bar, as well as the fatigue in materials due to pressure cycles (from vacuum to maximum pressure). The theoretical life expected for this equipment is about 20 years, with a work system of 4 operations /day during 365 days a year. However, the equipment real life is subjected to the periodical tests every 10 years and the annual checks.
- Corrosion: The use of genuine stainless steel in the autoclave vessel guarantees resistance to corrosion of water, steam and products to be sterilised.

6.3 Operating mode and range

The autoclave allows the configurations of the sterilisation cycle in accordance with the features to be sterilised:

Table. 6.3

Program type	With external internal sensor	Without sensor
Solids atmospheric purge:	105-134°C	105-130°C
Solids two fraction purge:	115-134°C	115-130°C
Solids three fraction purge:	115-134°C	115-130°C
Liquids atmospheric purge:	105-134°C	105-130°C

Sterilisation time: from 3 to 59 minutes

Drying time: from 0 to 99 minutes

6.4 Pre-established cycles

The Autester ST DRY PV II have 10 preset programs which can be modified depending on the users requirements.

Table. 6.4

Program	Temperature °C	Sterilisation	Drying	Vacuum
		Time (minutes)	Time	purge fractions *
0 Not wrapped instruments	134	4	10	2
1 Wrapped instruments	134	12	30	3
2 Prions	134	20	20	3
3 Not wrapped delicate	121	15	15	2
4 Wrapped delicate	121	30	30	3
5 Instruments	126	12	25	3
6 Disinfection	105	25	15	1
7 Liquids	121	30	-	1
8 Solids	121	15	30	1
9 Solids	134	4	30	1

* Purge types:

1. Atmospheric purge
2. Two phase vacuum purge
3. Three phase vacuum purge

6.5 Test programs

The Autester ST DRY PV II have two test programs:

Bowie Dick test	134°C	3' 30"
Vacuum test	0.7 bar	5+ 10'

Table. 6.5


Fig. 6.1 Control panel

6.6 Control panel elements

1. Parameter indicator LCD display
2. Single button selector
3. Led indicator
4. Start / Stop button
5. Printer (optional)


Fig. 6.7a

6.7 Main elements location

Fig. 6.7a

1. Open / Close door handle
2. Main power switch

Fig. 6.7b

Fig. 6.7b

1. Water tank inlet
2. Door joint

Fig. 6.7c

Fig. 6.7c

1. RS232 connector (serial port)
2. external sensor connector (Pt100)

Fig. 6.7d

Fig. 6.7d:

1. Safety valve
2. Vapor realease valve
3. External water entry
4. Air filter
5. Steam generator safety thermostats
6. Water tank outlet
7. Condensed outlet
8. Wastepipe valve

7 Installation

7.1 Location

Locate the autoclave near an appropriate power source for the machine.

In accordance with the current legislation (see 2) the autoclave should be located in a way that the safety valve discharge outlet do NOT point to any people or could not reach them in case of steam outlet.

Fig. 7.1 Wheel block

Fig. 7.2 Condensed steam outlet

Fig. 7.3 External water source inlet

The safety valve outlet in these equipments is on the rear part.

The autoclave needs to be installed over a flat stable horizontal surface, suitable for the machine weight (see table 3), leaving a free space of 20cm around the equipment.

To easily provide maintenance operations, it is recommended to have a direct access all around the autoclave.

Fix the autoclave, adjusting it to the floor by means of the front feet.

7.2 Installation without water source nor wastepipe.

In such a case, one must provide water to the autoclave tank in a manual way.

To collect the water from the steam condensation, fix the supplied hose from the condensed outlet to a container with the following features:

- Suitable to support hot water to 80°C.
- Maximum height of the container: 30cm
- It must be situated on the floor.

7.3 Installation with external water source inlet and wastepipe

The autoclave will be filled automatically from the external source when it detects a low water level in the tank.

The connection to the external source must be done by a 3/4" connection. See photo.

The water pressure must not be higher than 2bar.

This operation mode must be activated through the control panel. See chapter 9.6.

ATTENTION

Do not use the external waste water option until a pressure hose is installed.

Fix the hose in a way that it cannot move by the effect of the pressure at the end of the cycle.

Do not hold the hose with your hands. This hose can reach temperatures higher than 60°C.

Fig. 7.4a Water recycle

Fig. 7.4b Water waste

**7.4 Drain system selection:**

With the selected valve indicated in 7.4a and b.

7.4.1 Water recycling (Fig 7.4a)

Set the valve lever towards the inside of the machine.

7.4.2 Water waste to the outside (Fig 7.4b)

At the end of each sterilisation process, set the valve lever away from the equipment.

Fix a hose to the outlet valve with the adapter 34100 (see 4.1) and also fix the other end to a container or drain, taking care **NOT TO OBSTRUCT THE TUBE**.

ATTENTION!

The steam outlet from the safety valve and the steam and water outlet from the valve adapter(Fig 7.4b) are at HIGH PRESSURE and TEMPERATURE.

Take precautions to avoid SERIOUS BURNS.

DO NOT HOLD THE HOSE IN THE WATER OUTLET (Fig 7.4b) WITH THE HANDS, AS IT MAY CAUSE SERIOUS BURNS.

7.5 Electrical source connection:

The autoclaves are supplied without the power connection plug.

It is highly recommended to connect it directly to the source by means of a connection box.

In front of it, one must install an automatic bipolar ICP overdimensioning 1.5 times the autoclave maximum power.

Please consult the type and level of the supply voltage and the intensity maximum value in the features plate of this equipment.

The automatic switch must be located in a fast and easy accessible point.

7.5.1 Protection earth wire

The installation of power supply for the autoclave **MUST** have a protection earth wire.

The autoclave protection for electrical risks needs that the network cable, marked in yellow-green colour, will be connected to the installation earth wire.

**7.6 Installation qualification.**

Autoclave commissioning must according the country regulation.

On some countries a regulated installation qualification is mandatory.

One can obtain a list of authorized installers by consulting the government industrial department of its community.

User should contact the competent organisms to consult about the installation requirements.

8 Operation

8.1 Preliminary operations for switch on

- Switch on by means of the main power switch (Fig 6.7).
- Fill the tank (see Fig 8.1) with about 8 litres of water.
Check the water height in the tank by means of the water level in the left front side.

We recommend to change the water after 50 uses or once a month.

- Check the steam outlet valve is closed (12).
- Put the tray at the bottom of the vessel.
- Place the items to be sterilised.
- Close the door.

Note:

When the autoclave is cold, do a prewarming-up for 10 - 15 minutes before switching on.

To do the prewarming-up:

Plug the machine, switch the main power switch on and let the machine on for 10-15 minutes.

8.1.1 Opening and closing the loading and unloading door:

There exist a couple of safety devices which block the opening door if it's not safe to do so.

- Door opening and closing must be done by previously informed staff.
- Before opening the door, the user will make sure that the autoclave cycle is already finished and that there no exists pressure inside the vessel.

8.1.2 Cautions during usage:

The LCD display indications will be periodically checked. Pressure must not be higher than acceptable maximum pressure (PS = 2,5 bar)


ATTENTION!

Danger of serious burns

When the cycle finishes, the vessel load could be at a high temperature.

Please use the thermal protective gloves to remove the material from the vessel.

Fig 9.1: The single button selector has two movements:


9 Configuration, parameters and cycles selection menus

Once pushed the general power switch (see Fig. 6.7), the display shows a message indicating the controller software version for some seconds.

And then, it will appear a message corresponding to the last program used by default.

All functions will be selected by means of the single button selector: See Fig 9.1

9.1 Cycles selection

P0 134 C DRY 10
SOLID 2 PRESS START

Push the single button selector, and it will appear «PROGRAMS» in the display.

PROG/
PROGRAMS

Push again to enter the cycles selection.

PROG/P0
134 C 04
SOLID2
DRY 10

Turn the single button till you find the required cycle.

PROG/P1
134 C 12
SOLID3
DRY 38

Once found, push the single button selector.

SELECT

«SELECT» will appear, push again to validate.
The confirmation message «SELECTED» will appear for a few seconds.

And finally, the display with the selected program will be shown.

P1 134 C DRY 38
SOLID 3 PRESS START

Push the Start / Stop button, the LED will light up and the equipment will start the cycle.

9.2 Program modification

P0 134 C DRY 10
SOLID 2 PRESS START

This menu allows the modification of the cycle parameters. All cycles can be modified. (See 6.3 to know the limits of every parameter in a cycle)

PROG/
PROGRAMS

Through the «PROGRAMS» menu, we push the single button to access different programs.

PROG/P0
134 C 04
SOLID2
DRY 10

By turning the button we can select the program number which better fits our needs and also to modify it.

PROG/P1
134 C 12
SOLID3
DRY 38

Push the single button selector and the message «SELECT» will appear.

SELECT

Turn the single button till the message «MODIFY» appears, push again the button and the program we want to modify will appear in the display.

MODIFY

To modify a parameter: select the parameter by turning the single button selector, confirm it by pushing, modify it by turning the button and again confirm it by pushing.

P1 134 C DRY 38
SOLID 3 PRESS START

Once the program is modified, push the button and the message «SELECT» will appear. Push again to validate, and the message «SELECTED» will appear for a few seconds.

Finally, the display will show the selected program.

Push the Start / Stop button, the LED will light up and the equipment will start the cycle.

9.3 Language modification

**CONF/
CONFIGURATION**

**CONF/LNG/
LANGUAGE**

**CONF/UNIT/ENG>
ENGLISH**

EXIT

Turn the single button selector till the message «CONFIGURATION» appears in the menu display; Push again the single button and the message «LANGUAGE» will appear.

Push the single button and the display will show the last language memorized. Turn the button till you see the required language and push the button to validate it; the display will show «SELECTED» for a few seconds.

Turn the button again till the message «EXIT» appears.

Push the single button selector, the message «CONFIGURATION» will appear in the display.

Turn the single button selector till the message «EXIT» appears in the display.

9.4 Change °C to °F configuration

**CONF/
CONFIGURATION**

**CONF/LNG/
LANGUAGE**

**CONF/UNIT/ T (C)
UNITS**

EXIT

Turn the single button selector till the message «CONFIGURATION» appears in the display.

Push the single button. The message «LANGUAGE» will appear. Turn the button till the message «UNITS» appears.

Push the single button selector and it will change from Celsius to Fahrenheit degrees.


Turn again the button till the message «EXIT» appears.

9.5 Printer configuration

- CONF/ CONFIGURATION** 
- To activate and configure the printer, begin from the «CONFIGURACIÓN» menu and push the single button selector.
- CONF/LNG/ LANGUAGE** 
- The message «LANGUAGE» will appear in the display. Turn the button till the message «PRINTER».
- CONF/PRN/ PRINTER** 
- Pulsh the single button selector, and the display will show «PRINTER STA-TUS».
- CONF/PRN/SET OFF PRINTER STATUS** 
- The display will show the printer status (ON or OFF). By pushing the single button, one sets the printer On / Off.
- CONF/PRN/INT 01 PRINT INTERVAL** 
- Turn the single button selector and the message «PRINT INTERVAL» will appear in the display.
- Push the single button and turn it to change the time in minutes. Accept by pushing again.
- CONF/PRN/FED/FEED PRESS BUTTON** 
- Turn the single button. The display will show «PRESS BUTTON».
- This function advances the printer paper everytime one pushes the single button.
- This is used when changing the paper roll.
- EXIT** 
- Turn the single button selector till the message «EXIT» appears in the display.
- Push again to exit the «CONFIGURATION» menu.

9.6 External water source configuration

- CONF/ CONFIGURATION** 
- If you have installed the autoclave the way the tank is filled automatically from an external source (see 7.4) you must configure this function.
- CONF/LNG/ LANGUAGE** 
- Beginning from the «CONFIGURATION» menu, push the single button selector to enter.
- CONF/WTR/ OFF EXTERNAL WATER** 
- The message «LANGUAGE» will appear. Turn it till «EXTERNAL WATER» is shown.
- The display will show the external water inlet valve status: Open (ON) or closed (OFF).
- By pushing the single button, the printer will be Open / Closed.
- EXIT** 
- Turn the single button selector till «EXIT» and then push. Repeat it to exit to the main display.

9.7 Temperature exterior sensor configuration.

You must configure this function if you have installed a second sensor, for liquid sterilisation cycles (see 7).


Beginning from the «CONFIGURATION» menu, push the single button selector to enter.

The message «LANGUAGE» appears. Turn the button till «EXTERNAL PROBE».

This display shows the external sensor status: On or Off.

By pushing the single button, the external sensor is set On /Off.

Turn the single button till «EXIT» and then push.
Repeat this action to exit till the main display.

9.8 Date and time configuration

Use this menu to update time and date in the autoclave clock.

Beginning from the «CONFIGURATION» menu, push the single button selector to enter.

The message «LANGUAGE» appears.

Turn the button till «CLOCK»

Push the single button selector, the display will show the message «CHANGE DATE».

The date and the time will appear in the display.

To modify the day, month and year, select it by turning the single button, then confirm by pushing it, modify it by turning again and then confirm by pushing the button once again.

**CONF/
CONFIGURATION**

9.9 Vacuum test

To complete the vacuum test, please follow the following instructions:

**TEST/
TEST**

Push the single button selector till the «TEST» menu appears in the display.

**TEST/VAC> OFF
VACUUM TEST**

Push the single button and the display will show «VACUUM TEST».

**TEST/VAC> ON
VACUUM TEST**

Push it again to select vacuum test «ON».

**VACUUM TEST
PRESS START**

Push the button and the message «PRESS START» will appear in the display. The machine will be ready to complete the vacuum test by pushing the start/stop button (4), and lighting up the LED indicator (3).

**PROG/
PROGRAMS**

9.10 Bowie-Dick test

To complete the vacuum test, please follow the following instructions:

**CONF/
CONFIGURATION**

Push the single button till the «PROGRAMS» menu appears in the display.

**TEST/
TEST**

Turn the single button selector till the message «TEST» appears in the display.

**TEST/VAC> OFF
VACUUM TEST**

Push the single button selector and the display will show the message «VACUUM TEST», turn then the single button till «BOWIE-DICK» appears.

**TEST/DCK> OFF
BOWIE DICK TEST**

Push the BOWIE-DICK test «ON».

EXIT

Turn the single button till the message «EXIT» appears.

**BOWIE DICK TEST
PRESS START**

Push the button and the display will show «PRESS START». The machine is ready to complete the Bowie-Dick test by pushing the start/stop button(4), and then lighting up the LED indicator (3).

10 What to do when getting an error during the cycle operation

The autoclave control system can identify an incorrect cycle operation.

In case the cycle interrupts, for instance, due to power supply, lack of water, etc... one must consider that the load has not been sterilised and sterilisation must be repeated again.

As general rule: switch off the autoclave, using the main switch, wait for some seconds and the switch on again.

If the door is blocked, wait about 20 minutes for the vessel temperature to fall.

You can accelerate the process by opening the depressurization manual valve. (Please, close the valve before switching it on again).

Display	Error	Description
ERROR 01! EEPROM FAIL	01	Defective EEPROM error.
ERROR 02! PRESSURE SENSOR FAIL	02	Pressure sensor error. Defective sensor or cut / disconnected cable.
WATER TANK EMPTY	10	Water tank empty.
DOOR OPEN	12	Door open.
ERROR 13! OVERPRESSURE	13	Overpressure. The maximum admissible pressure of +2.5bar has been exceeded.
ERROR 14! VACUUM FAIL	14	Vacuum error. There's a vacuum leak or the pump is damaged.
ERROR 15! PRESSURE FAIL	15	Pressure error during fractions. There's no water or the resistance is defective.
ERROR 16! INTERNAL OVERTEMP.	16	Tank overtemperature. The internal temperature has exceeded the 150°C.
ERROR 17! EXTERNAL OVERTEMP.	17	External overtemperature. The external temperature has exceeded the 150°C.
ERROR 18! STERILIZATION FAIL	18	Sterilisation error. The margin error has been exceeded or there's still air in the inside.
ERROR 19! SETPOINT UNREACHABLE	19	Value not reached. There's no water omsoode or the resistance is defective.
ERROR 20! EMPTY TOO LONG	20	It takes too long to empty the pressure and the water in the tank. Valves or pipes blocked.
ERROR 21! PRINTER PAPER END	21	There's no paper in the printer.
ERROR 22! PRINTER ERROR	22	Printer error. Possible paper jam.
ERROR 23! POWER FAIL	23	Power failure error (detected when switching on).
ERROR 26! TEST FAILED	26	Vacuum test error.

11 Recommendations for perfect sterilisation

The material to be sterilised must be perfectly clean, free from any encrusted or residual material. We recommend to wash it with a good detergent and distilled water and after that, rinse with plenty of water.

Do not overload the baskets, the trays or the grids, trying to leave always a space of 1 or 2 cm between them, allowing a free way for vapour and easy drying.

11.1 Instruments sterilisation in bulk:

Place the items over the water resisting paper leaving a space in between them.

Do not place different kind of material instruments in the same tray.

11.2 Instrument sterilisation in bags:

There must not be more than one instrument per bag.

Place the bags over the water resisting paper, in a vertical position, in a way they do not touch each other neither the autoclave walls.

11.3 Tubes sterilisation:

Tubes must be perfectly clean, rinsed and drained. Place the tubes, if possible open on both ends.

Try the tubes not to touch the autoclave walls nor block the inner diameter.

11.4 Vessels sterilisation:

Never place the vessels hermetically sealed.

Place them face downwards to avoid any residual water left inside.

11.5 Liquids sterilisation:

Place the liquids to be sterilised over the trays in prepared containers to withstand sterilisation temperature, so one can pick any possible spillages.

The containers will be fulfilled about 2/3 of its capacity. DO NOT SEAL THEM HERMETICALLY. Just close them with a cotton or without anything to let the containers air pressure release.

Try to avoid using narrow neck containers.

When finishing the sterilisation cycle, let the autoclave cool down slowly (till pressure is 0 kg/cm²).

12 Safety devices

The autoclave has the following safety devices:

12.1 Safety pressurestat: Adjusted to 2.35bar

It cuts the electrical power of vapour generator (or electrical resistance) to stop the vapour production when an unusual overpressure is produced in the sterilised vessel.

12.2 Safety valve: Adjusted to 2.5 bar


In case the vessel pressure increases over 2.5bar, this valve is automatically open, releasing the vapour from the inside of the vessel and, consequently, depressurizing the vessel.

In case the safety valve is triggered, stop the autoclave and call the technical service department.

12.3 Safety thermostat: No adjust possibility.

In case of unusual increase of temperature or a considerable decrease of the water level (which produces an overtemperature of the vapour generator), the safety thermostat activates (14) and stops the heater.

This thermostat cannot be reclosed. Stop the autoclave and call the technical service department.


12.4 Door blockage with temperature > 98°C

In case of door blockage due to the safety system break, please follow the protocol below:

- 1 Stop the machine.
- 2 Unplug the machine from the power supply.
- 3 Put the tube on the air valve (11).
- 4 Open carefully the valve to release the vapour pressure. Take care to prevent burns.
- 5 Wait till all the vapour is released.
- 6 Proceed the following way, using the hook provided with the machine (see photos):

Look through the lid grid with the aid of a torch and find the closing mechanism. You could see the opening lever.

Once located, insert carefully the hook through the grid and pull the lever outwards. Hold the lever this position while simultaneously opening the lid with the door handle.

13 Calibration

This section describes the parameters that can be changed through the calibration menu. These parameters should only be changed when the initial selected values do not satisfy the applications requirements.

How to access the calibration menu:

Push the single button selector and select the «CONFIGURATION» menu inside «PROGRAMS».

Push again and the message «LANGUAGE» will appear; then turn the button till the display shows «PASSWORD».

Push the button to select «ENTER PASSWORD». The defective code is «0000» (the user can change the password by selecting «CHANGE PASSWORD»); push again and the message «ENTER PASSWORD OK» will appear, then turn the single button selector till «EXIT».

Push the single button till «CONFIGURATION», turn it till «CALIBRATION», and push the button to enter the following submenus:

«INTERNAL TEMP»: Inserts a (constant) correction in the temperature measure.

«EXTERNAL TEMP»: Inserts a (constant) correction in the second sensor of the temperature measure.

«PRESSURE»: Inserts a (constant) correction in the pressure measure.

«PREHEAT»: Inserts a (constant) correction in the boiler heater temperature (belt heating element).

«COUNTER»: This is the counter for the autoclave executed cycles.

These submenus allow the modification of the zero of pressure and temperature sensors, change the fixed temperature of each program and see the sterilisations quantity completed by the machine.

Turn the single button selector till «EXIT», then push it till you see the message «CONFIGURATION», and turn again till the display shows «EXIT».

Then push the button and we'll return to the main display.

14 Maintenance

Before removing the equipment's lid, please disconnect it from the power supply.

Any repair or maintenance operation should only be made by authorized personnel of J.P. SELECTA S.A.

14.1 Basic maintenance operations

14.1.1 Cleaning:

We recommend the following products to clean the different parts of the equipment:

To clean the vessel and the stainless steel outside: Alcohol.

To clean the covers and plastic areas: Alcohol with a cotton or with a non abrasive cloth.

It's recommended to change the water after 50 cycles or once a month.

14.2 Maintenance operations on demand by U.E. legislation

Concerning the current legislation of pressure equipments:

- Directive 97/23/CE (Europe Union)
- Pressure Equipments Regulations RD 1244/1979 (Spain)

14.2.1 Annual reviews

Once a year, the autoclave must have its safety elements revised. This check must be done by competent staff. Please contact J.P. SELECTA, S.A. to get a list of authorized technical services.

Check once a year the safety valve operation, the pressure switch, the levels control and the temperature and pressure sensors.

Please be advised that the autoclave is protected by the enclosure and so, the user cannot have access, as there could be very hot parts and casual pressure hot leaks could happen. It is recommended to take the following cautions when accessing these parts:

1. Disconnect and cool the equipment.
2. Use suitable protections
3. Do not operate devices you don't know, especially manual valves, taps or atmospheric connections, where a leak could happen.

14.2.2 Periodic hydraulic tests:

Check country regulations if periodic test are mandatory.

On some countries each 10 years from the equipment installation, and once made an annual review (see 14.2.1.1), the pressure vessel must have a hydraulic pressure test according to regulated procedure.

The test will be made by a personnel authorized by the administration.

J.P. SELECTA can provide the data of some accredited companies in Spain.

15 Spare parts

15520 Air filter


- 21199 Autoclave joint 30 & 50 litres
21200 Autoclave joint 75 litres
21206 Autoclave joint 140 litres
21193 Autoclave joint 25 litres



- 29309 Autoclave electronic main board 25 L
29343 Autoclave electronic main board 30 L
29390 Autoclave electronic main board 50 L
29344 Autoclave electronic main board 75 L
29345 Autoclave electronic main board 140 L

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