***RECOVER YOUR ENERGY &
ENHANCE YOUR VITALITY***

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THANK YOU FOR CHOOSING OUR PRODUCT. WE REMAIN AT YOUR ENTIRE DISPOSAL FOR ANY ASSISTENCE OR ADVICE YOU MAY NEED.

El Electroestimulador Prostim 2 es fabricado para:
The electrostimulator ProStim 2 is manufactured for:

REHAB EUROPA SL
Berenguer de Palou, 86-90, Local 1
08027. Barcelona, España
Tlf.: (+34) 933527650
Fax: (+34) 933527595
E-mail: info@rehabmedic.com

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This product has been manufactured according to the technical regulations in force and is certified according to Directive 93/42/EEC updated by 2007/47 directive for medical devices, by Kiwa Cermet Italia S.p.A. No. 0476, in order to ensure the product's safety.

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tensión/flicker Voltage fluctuation/flicker emissions IEC 61000-3-3	are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
---	--

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***GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC
IMMUNITY – FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING***

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The Prostim 2 electrostimulator are intended for use in the electromagnetic environment specified below. The customer or the user of the Prostim 2 electrostimulator should assure that it is used in such an environment.

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			<i>no closet to any part of the Prostim 2 , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</i> Distancia de separación recomendada Recommended separation distance
RF conducida <i>Conducted RF</i> IEC 61000-4-6	3 Veff_Vrms de 150 kHz a 80 MHz 150 kHz to 80 MHz 6 Veff_Vrms de 150 kHz a 80 MHz banda ISM 150 kHz to 80 MHz for ISM band	3 Veff_Vrms ([V₁] V)	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P} =$ $d = \left[\frac{12}{V_1} \right] \sqrt{P} \text{ for ISM band}$
RF radiada <i>Radiated RF</i> IEC 61000-4-3	10 V/m de 80 MHz a 2,7 GHz 80mhz to 2,7 GHz	10V/m [E₁] V/m	$d = \left[\frac{12}{E_1} \right] \sqrt{P} \text{ de 80 MHz a 800 MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \text{ de 800 MHz a 2,7 GHz}$

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Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Prostim 2 electrostimulator is used exceeds the applicable RF compliance level above, the Prostim 2 electrostimulator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the A3R electrostimulator.

- b La intensidad de campo en el intervalo de frecuencia de 150 kHz a 80 MHz debería ser menor que [V₁] V/m Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V₁] V/m.

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**RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF
COMMUNICATIONS EQUIPMENT FOR PROSTIM 2 THAT ARE NOT LIFE-
SUPPORTING**

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The Prostim 2 series electrostimulator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Prostim 2 series electrostimulator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Prostim 2 series electrostimulator as recommended below, according to the maximum output

<i>power of the communications equipment.</i>					
Potencia de salida máxima del transmisor especificada <i>Rated maximum output power of transmitter</i> W	Distancia de separación a la frecuencia del transmisor (m) <i>Separation distance according to frequency of transmitter (m)</i>				
	150 kHz to 80 MHz	150 kHz to 80 MHz (ISM band)	De 80 MHz to 800 MHz	80 MHz to 2,7 GHz	De 2.7 GHz to 6 GHz (to RF wireless radio communication equipment)
0,01	0,12	0,2	0,12	0,23	–
0,1	0,38	0,63	0,38	0,73	–
0,2	–	–	–	–	0,9
1	1,20	2,0	1,20	2,30	–
1,8	–	–	–	–	2,7
2	–	–	–	–	2,8
10	3,80	6,3	3,80	7,30	–
100	12,00	20	12,00	23,00	–

PROstim 2

USER'S MANUAL

ENGLISH

39-67

TECHNICAL FEATURES	42
Device	42
Use conditions	42
Technical features of the currents	42
Charger	42
EQUIPMENT	42
INTENDED USE	43
CONNECTIONS	44
Accessories	44
Electrode application	44
LABELLING AND SYMBOLS	44
PANEL/KEYBOARD	46
DISPLAY	46
WARNINGS AND CONTRAINDICATIONS	46
Mandatory behavior	46
Warnings before use	47
Warnings during the use	47
Preventative measures to take using the device for incontinence treatments	47
SIDE EFFECTS AND CONTRAINDICATIONS	48
Contraindications	48
MAINTENANCE AND CLEANING	48
Dispositivo	48
Accessorios	49
Battery: how to charge the batteries	49
Disposal of the device	49
INSTRUCTIONS FOR USE	50
Programming	51
ELECTROSTIMULATION	52
TENS	54
INCONTINENCE PROGRAMS	55
PROGRAM LIST	56
ELECTRODE PLACEMENT	61
FREQUENTLY ASKED QUESTIONS	65
WARRANTY	66
EMC TABLES	28

TECHNICAL FEATURES

DEVICE

Size:	13 x 8 x 2 cm
Weight:	220 g.
Case:	in Food Grade ABS
Protection level:	IP 22
Storage and transportation temperature:	from -10°C to 45°C
Max. relative humidity:	30% - 75%
Made in PRC	

USE CONDITIONS

Temperature:	from 0°C to 35°C
Max. relative humidity:	from 15% to 93%
Atmospheric pressure:	from 700 hPa to 1060 hPa

The values indicate the limits allowed if the product or its accessories are not in the original package.

TECHNICAL FEATURES OF THE CURRENTS

EMS and TENS:

Channels available:	Channel 1-2
Constant current:	Yes
Intensity:	0-100 mA p.p. per channel
Wave form:	Rectangular, biphasic, symmetric, compensated
Working frequency:	1-150 Hz*
Recovery frequency:	1-150 Hz*
Pulse amplitude:	50-450 µs*
Working time:	from 1 to 30 seconds
Recovery time:	from 0 to 1 minute
Frequency modulation range:	continuous variation from 1 to 150 Hz
Min. modulation time:	3 seconds
Amplitude modulation range:	continuous variation from 50 to 400 µs


CHARGER

Model:	DK7-065-0200-EU
INPUT:	100-240Vac - 50-60Hz - 0,06 A
OUTPUT:	5,8Vdc - 200 mA
Polarity:	
Battery pack:	Ni-MH AAA 4,8V 800 mAh

EQUIPMENT

The electrostimulator is supplied complete of cables and electrodes to use: therefore, opening the package, it is necessary to check that the basic equipment is complete. If some elements should be missing, contact immediately the authorized retailer where you purchased the product. Control carefully the integrity of the device and its electrodes.

* Depending on the type of stimulation

- a. Device.
- b. 2 colored electrode connection cables (for EMS and TENS treatments).
- c. 1 bag containing 4 reusable self-adhesive electrodes (50 x 90 mm) (It is suggested to use these electrodes for big areas such as thighs, abdomen and glutei...).
- d. 1 bag containing 4 reusable self-adhesive electrodes (50 x 50 mm) (It is suggested to use these electrodes for small areas such as upper limbs, calves, cervical areas...).
- e. Charger (see technical features).
- f. Carrying bag.
- User manual and warranty

All the supplied information can be modified without previous notice.

The device can be used with some optional accessories. If you are interested in buying these accessories, please contact the retailer.

ACCESSORIES NOT INCLUDED (available on charge)

- Find Motor point pens
- Kit of 8 elastic bands for legs and thighs
- Kit of 4 elastic bands for thighs
- Face electrodes
- Kit Y cables

INTENDED USE

The after sales service is guaranteed for 5 years. We suggest having a check of the device every 2 years for the maintenance and to ensure the safety. The numbers of treatments depends on the battery charge.

The electrostimulators are designed to be used at home environment to make aesthetic, sporting and pain treatments.

The use of this device is permitted to the patient itself (properly informed about the use conditions) and physicians.

CONNECTIONS

ACCESSORIES

Cable connections: It is possible to use one or two cables according as you have chosen to use one or two channels of the unit.

To connect the cable connector to the device plug it into the appropriate inlet in the upper part of the unit.

The inlets are situated exactly under the correspondent channel (left or right). Insert the cable at the channel you want to use.

ELECTRODE APPLICATION

Take the electrodes from the original package; all new electrodes have a seal on the package. Be sure that the device is off. To start, connect the two cable plugs to the electrodes, then disconnect the electrodes from their position and apply them on the skin. To place the electrodes correctly, see the pictures included in this manual.

After the use, place the electrodes in their original position again.

ATTENTION: Do not unplug the electrodes if the unit is working.

LABELLING AND SYMBOLS

It refers to the manufacturer

Warning

This symbol on your device indicates that it complies with the directives on medical devices (93/42/CEE 47/2007CEE). The number of the notified unit is 0476.

It indicates that this is a II class device

It indicates that this device has type BF parts

WEEE symbol (Waste of Electrical and Electronic Equipment). Recycling symbol. The WEEE symbol used for this product indicates that the device may not be disposed of as a household product. Properly dispose of the product to help protect the environment. For more information on recycling this product, contact the local competent department, the household waste management company or the store in which the product was purchased

It indicates the optimal temperatures for the storage and transportation of the product

It indicates that the product has been produced in compliance with the directive 2011/65/EEC.

It informs the operator that she must read the manual before using the device.

IP22

It indicates the water protection degree

It informs the operator of a compulsory conduct

It refers to the pressure of the environment where the device is used

It refers to the humidity of the environment where the device and its accessories are used

It refers to the expiry date

It refers to the production lot


It refers to the manufacturing date

Indoor use

CHARGER

INPUT	INPUT: It indicates the value of the mains voltage for the power supply unit.
OUTPUT	OUTPUT: - Indicates the power supply unit outbound voltage - Indicates the maximum power value of the magnetic field emitted by the device - Indicates the range of frequencies of the magnetic field emitted by the device
MODEL	Indicates the device type; power supply model unit of the device; battery pack inside the device

ELECTRODES


The electrodes supplied are designed for use with a single patient adhesive, pregelados with female cable. They are CE marked according to directive 93/42 / EEC on medical devices. All information provided may be subject to change without notice.

	It indicates the dimensions of the electrode
	It indicates the number of electrodes contained in the package
REF	It indicates the product code
CE	It refers to product certification and indicates that it complies with directive 2001/95/EC updated as 2014/357/EU
	It indicates the storage temperature of the electrodes


- Clean and degrease the skin.
- Do not apply the electrode on wounds or injured skin.
- Connect the cable connector to the electrode connector.
- Remove the electrode.
- Apply on the skin.
- Start the program.
- At the end, turn off and put the electrode back in the package.
- Electrodes are for personal use.
- Do not remove the electrode by grabbing the connector.
- Electrodes should not touch each other.
- Do not apply the electrodes on the temples, the neck and in a transthoracic way.
- Do not leave the electrodes in the car.

PANEL/KEYBOARD

DISPLAY

WARNINGS AND CONTRAINDICATIONS

MANDATORY BEHAVIOR

For safety reasons, the device must be used in the prescribed manner and within the limits of use explained in this manual.

The manufacturer declines any responsibility due to a different use from what it is indicated in this manual.

No part of this manual (texts and photos) may be reproduced by electronic or mechanical manner without the prior written authorization of the manufacturer.

It is suggested not to carry out treatments in presence of skin wounds.

If the package, the cable or the connector of the power supply unit show signs of wear or damage, replace it instantly.

The unit should be connected to the mains by its power supply unit. Before starting the treatment, make sure that the power system specifications comply with the directives in force within your country. Make sure that the power supply unit will be in a comfortable position and that it will be easy to be removed.

WARNINGS BEFORE USE

We suggest avoiding the use of the device together with other electronic devices, especially with those used to maintain vital functions. For a correct use of the Electromedical device, make reference to the tables attached. If it is necessary to use the device near or together with other devices, pay attention to its working.

- It is recommended to read carefully the entire operating manual before using the device; keep carefully this operating manual.
- The device can emit current value above 10mA.
- Before each use always check the integrity of the unit. This is a fundamental requirement for carrying out therapies, do not use the unit if either the buttons or cables are defective or malfunctioning.
- It must be used only by people over 18 who are able to understand and take action.
- It should not be used for purposes other than neuromuscular stimulation as described in this manual.
- It must be used following the indications and under the physician or physiotherapist's control;
- It must be used with the supplied electrodes intended for transcutaneous neuromuscular stimulation.
- It must be kept out of the reach of children.
- ECG monitoring devices may not operate properly when electrostimulation is working.
- Do not use the device in trans-thoracic modality because it may cause cardiac arrhythmia, imposing its frequency over the heart one. Do not stimulate the pectoral and dorsal muscles simultaneously.
- If there is any health problem, it must be used only after seeing a doctor.
- A simultaneous connection of the patient to a high frequency electrosurgery device can cause burns where the electrodes are placed and the electrostimulator can result damaged.
- After turning on the device, make sure that in the display the software version and the model of the device appear, it means that the device is working and it is ready to be used.

If not, or in the display all the segments do not appear, turn it off and on again. If the problem persists, contact the service center and do not use the device.

An unexpected switching-off of the device means that the battery has run-down. Charge it according to the instructions in the paragraph: how to charge the batteries.

WARNINGS DURING THE USE

While using the electrostimulator some warnings should be followed:

- In case of damaged cables, they must be replaced with original parts and no longer used;
- Use only the electrodes supplied by the company;
- The device must be kept out of the reach of any pet that could damage it and contaminate it with parasites;
- Pay particular attention when the current density for every electrode is above 2mA/cm² (effective value).
- The cables of the electrostimulation should not be wrapped around people's neck to avoid any risk of strangulation and suffocation.
- The mobile and fixed radio communications devices could influence the functioning of the electromedical device: see the tables attached to this manual.

PREVENTATIVE MEASURES TO TAKE USING THE DEVICE FOR INCONTINENCE TREATMENTS

- Patients with extra-urethral incontinence should not be treated with the electrostimulator.
- Patients suffering excessive incontinence due to evacuation problems should not be treated with the electrostimulator.
- Patients with severe urinary retention to the upper urinary passages should not be treated with the stimulator.
- Patients with total peripheral denervation of the pelvic floor should not be treated with the stimulator.

- Patients suffering of a total/subtotal prolapse of the uterus/vagina should be stimulated with extreme care.
- Patients with infections to the urinary passages should be treated for these symptoms before starting the stimulation treatment.
- Before removing or touching the probe it is necessary to turn off the stimulator or to regulate the intensity of both the channels on 0,0 mA.
- The treatment is a personalized medical prescription: do not borrow the stimulator to other persons.

SIDE EFFECTS AND CONTRAINDICATIONS

Isolated incidents of skin irritation may occur in subjects with high epidermal sensitiveness. In case of an allergic reaction to the electrode gel, suspend the treatment and contact a specialist. If during the treatment signs of tachycardia and extrasystole appear, suspend the treatment and contact your physician.

CONTRAINDICATIONS

Do not use the device in the following cases:

- Stimulation of the front part of the neck (carotid sinus).
- Pacemaker weavers.
- Patients with tumor diseases (see your oncologist).
- Stimulation of the brain region.
- Pains whose etiology is unknown.
- Sores and dermatological diseases.
- Severe traumas.
- Stimulation on recent scars.
- Pregnancy.
- It is strictly forbidden to use the electrostimulator in the ocular area.
- Near body areas with metallic implants or infratissue metals (prostheses, osteosynthetic devices, coils, screws, orthopedic plates), when using monophasic current, interferential, or continuous current, ionophoresis.

It is recommended to use the device with caution in people with capillary fragility, as an excessive stimulation could cause capillary ruptures.

MAINTENANCE AND CLEANING

DEVICE

Maintenance and cleaning of the device

- In case of real or alleged malfunction, do not tamper with the device or try to repair it yourself. Do not intervene on the device and do not open it. Only specialized and authorized centers can repair it.
- Avoid violent impacts that may cause damage and malfunctions to the device even if undetectable immediately.
- Use this device in a dry environment and in an open space (not wrapped in any materials).
- Clean the device and accessories only with disinfectant with sodium hypochlorite or quaternary ammonium salt diluted with distilled water equal to 0,2-0,3%. After cleaning/disinfecting it, dry the device and its accessories with a clean cloth.
- Always use the device and its accessories with clean hands.
- It is recommended to use the device in a clean room, to avoid the contamination of the device with dust and dirt.

- It is recommended to use the device in a ventilated space, with regular air change.
- It is recommended to clean/disinfect the parts after every use.

ACCESSORIES

Use and storage of the electrodes and the cables.

After using the multi-purpose single patient and/or single-use electrodes, they must be stored using their plastic film and placed in the plastic bag.

Avoid that the electrodes touch each other or that they overlie one over the other.

Once the package has been opened, the electrodes can be used for 25-30 applications.

The electrode must be always used with clean hands and they must be replaced if they are not perfectly in contact with the skin.

If using non self-adhesive electrodes it is suggested to clean the surface with proper cleansers that respect the requirements described in the manual.

The electrodes must be stored in their bag and in an environment that respects the requirements described in the manual.

After the end of a treatment, unplug the cables from the connectors and clean them carefully with proper cleansers that respect the requirements described in the manual.

After cleaning and drying them, they must be folded up and placed in the plastic bags supplied with the cables.

BATTERY: HOW TO CHARGE THE BATTERIES

The device is supplied with a set of rechargeable nickel-metal hydrate batteries (4,8V 800mAh).

When in the display the battery to be charged icon () appears, we suggest recharging the device, as after the appearance of this symbol the electrostimulator autonomy is limited and could not allow the execution of a complete program.

To charge the batteries, turn off the electrostimulator and disconnect the electrodes, then connect the electrostimulator to the charger provided by inserting the plug in the appropriate inlet.

We suggest charging the device for 8-12 hours to obtain the greater duration of the battery pack.

A possible increase in the temperature in the battery area during the charge is normal.

Do not use a different charger from the one provided with the device. To replace the electrostimulator batteries, contact an authorized service center.


DISPOSAL OF THE DEVICE

Do not throw the device or part of it in the fire, but dispose of the product in the specialized centers and respecting the directives in force within your country.

When the product has to be disposed of, the user can give it back to the retailer when purchasing a new unit.

A correct separate waste collection or following what above mentioned contribute to avoid possible negative effects on environment and health and promote the reuse and/or recycle of materials of which the device is composed. The illegal disposal of the product by the user involves the application of the administrative fines according to the current regulations.

INSTRUCTIONS FOR USE

PROGRAMMING

Press Fn for 2 seconds in the start screen to program the device with 6 modes (PED1 - PED6).

Press OK to pass from a parameter to the next one.

It is possible to change the duration of all programs from 5' to 100' with steps of 5'.

Default duration 20'.

PED1:

It is possible to set pulse frequency and duration.

From 1 to 150 Hz with steps of 1 Hz.

From 50 to 250 μ s with steps of 10 μ s.

Default values: 100 Hz / 200 μ s.

PED2:

It will be possible to set pulse duration with a fixed frequency of 100 Hz.

From 50 to 250 μ s with steps of 10 μ s.

Default values: 200 μ s.

PED3:

It is possible to set frequency values.

From 1 to 150 Hz with steps of 1 Hz.

Default values: 100 Hz

Modulated pulse duration from 100 μ s to 250 μ s.

PED4:


It is possible to set pulse frequency and duration.

From 1 to 150 Hz with steps of 1 Hz.

From 50 to 320 μ s with steps of 10 μ s.

It is possible to change the following values: t_{up} = 1 s; t_{hold} = 5 s; t_{down} = 1 s; t_{off} = 5 s.

Default values: 30 Hz / 250 μ s.

PED5:

It is possible to set frequency values.

From 1 to 150 Hz with steps of 1 Hz.

Default values: 30 Hz

Modulated pulse duration from 100 μ s to 320 μ s.

PED6:

It is possible to set tHOLD and tOFF values.

From 1 to 30 s with steps of 1 s.


It is possible to set the frequency value.

From 1 to 150 Hz with steps of 1 Hz.

Pulse duration fixed at 300 μ s.

Default values: tHOLD = 2 s; tOFF = 10 s; 30 Hz.

tUP and tDOWN cannot be changed and have a duration of 1 s.

ELECTROSTIMULATION

ELECTROSTIMULATION ADVANTAGES

The electro stimulation does not want to replace physical activity but it must be considered as an integrating treatment.

According to the aims it will be useful:

- To the sportsmen/sportswomen, in order to integrate the normal training and to increase performances.
- To who wants to intervene on some imperfections.
- To who suffers from some pathologies.
- To who has suffered from a trauma or in the rehabilitation.

The electro stimulation is a technique that provokes a muscular contraction completely similar to the volunteer one, using electrical TENS impulse. There are two different ways of use:

- Muscular stimulation (ideal for the development of the force and for aesthetic treatments).
- The stimulation to nervous ends (ideal for the treatments against the pain).

TYPES OF MUSCLE

The muscle can be subdivided in three different types: streaked or voluntary muscle; cardiac muscle and smooth or involuntary muscle.

The voluntary muscle includes the skeletal muscle that gives movement to the skeleton and the cutaneous muscle that acts on facial expressions.

The cardiac muscle and the smooth muscle are not voluntarily controlled. The most part of the muscles of the human body belongs to the category of streaked or voluntary muscles, with approximately 200 muscles for every side of the body (approximately 400 in total).

Skeletal muscles are the target of the EMS.

MECHANISM OF THE MUSCULAR CONTRACTION

The skeletal muscle exercises its functions through the mechanism of the contraction.

When the muscular contraction happens, the movement of the articulations is produced and consequently the movement of the skeleton, too.

The muscle contracted in the following way: when a person decides to make a movement, an alteration of the electric potential in the motor center of the brain and an electric impulse is transmitted to the muscle that has to move.

The electrical impulse is transmitted as a change of potential (voltage); an electrical impulse crosses the motor nerve to the muscle that must be stimulated. After the reception of the impulse,

Ca ions are freed by the reservoirs causing the approach of actin and myosin. The result is that the distance between the bands Z diminishes causing therefore the muscular contraction.

The energy demanded for the contraction is supplied by the supply of sugars and fats stored in the human body. In other words, the electrical stimulation is not a directed energy resource but it works as an instrument that triggers the muscular contraction. The same type of mechanism is activated when the muscular contraction is produced by the EMS.

They assume in other words the same role of a natural impulse transmitted by the motor nervous system. Normally the muscle relaxes and returns to its original state them at the end of the contraction.

ISOTONIC AND ISOMETRIC CONTRACTIONS

The isotonic contraction reveals when, in a movement of our body, the interested muscles produce a state of constant tension. The isotonic contraction reveals when, in a movement of our body, the interested muscles produce a state of constant tension. Instead, when the joint heads of a physical segment are blocked and the musculature produces a tension, we talk about an isometric contraction.

In the event of the electro stimulation an isometric contraction is preferred because it fosters to obtain a more effective control.

THE DISTRIBUTION OF THE DIFFERENT TYPES OF FIBERS IN THE MUSCLE

The relationship between the two main categories (type I and the type II) can vary in sensitive way. There are muscular groups that are typically constituted by fibers of type I, like the soleus, and there are muscles that have only fibers of type II like the orbicular muscle. But in the most part of the cases we have a simultaneous presence of various types of fibers.

The studies lead on the distribution of fibers in the muscle have put in evidence the strait relationship that exists between the motoneuron (tonic or phasic) and the functional features of fibers from it innervate and have demonstrated as to a specific motor activity (and sports in particular) can determine a functional adaptation of the fibers and a modification of the metabolic features of the same ones.

Motor unit	Concentration type	Fiber type	stimulation frequencies
Tonic	ST slow contraction	I a	10 - 50 Hz
Phasic	FTa fast contraction	II a	50 - 70 Hz
Phasic	FTb fast contraction	II b	80 - 120 Hz

In order to make a tissue pass from the phase of rest to the one of excitation, by means of an induced electrical stimulus (impulse of the electrostimulation), some conditions are necessary:

- the current amount.
- the duration of the stimulus that must be adapted to the corporeal district that is wanted to be stimulated.

From this consideration, for the excitation of a tissue, the relationship between the duration and the intensity of the stimulus is particularly important and not only the value of the peak of the intensity. This relationship varies according to the muscular district.

STIMULATION INTENSITY

The current intensity necessary to obtain muscular contraction is personal and depends on the position of the electrodes, the underlying adipose tissue layer, sweating, the presence of hair on the area to be treated, etc. For these reasons, the same current intensity can generate different feelings from person to person, from day to day, and from the right side to the left side of the body. During different sessions, it will be necessary to regulate the intensity in order to obtain the same level of contraction because of the accommodation phenomena. The current intensities recommended in the different phases are proposed as indicative values, and each person should

modify these levels according to his/her personal needs.

- Moderate intensity. The muscle does not tire, not even during prolonged treatments. The contraction induced is tolerable and pleasant. This is the first level on the graphic representation of intensity.
- Intermediate intensity. The muscle is visibly contracted but the stimulation does not cause the joints to move. This is the second level on the graphic representation of intensity.
- Elevated intensity. The muscle is contracted substantially. The muscular contraction will cause the extension or bending of the limb if this is not blocked. This is the third level on the graphic representation of intensity.
- Maximum intensity. The muscle is contracted maximally. This is an intense treatment that should be performed only after having executed different applications at lower.

In the descriptions of the treatments, the best levels of intensities are recommended.

NOTE: The recommended levels of current are only indicative.

NOTE: For Microcurrent programs, it is not necessary to set an intensity value (in mA) because this is preset and automatically activated for all phases.

TENS

Transcutaneous Electrical Nerve Stimulation (TENS) is a selective stimulation of the large fibers of the peripheral nerves favoring the closing of the gate entrance for the pain pulses and increasing the release of endorphin substances, reducing in this way the pain intensity. Therefore with TENS we want to treat the severe and chronic pain due to the main musculoskeletal pains.

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The pain decrease following to the TENS currents application is due to these factors:

- Gate control theory.
- Endorphin secretion.
- Different sedative effects in relation with the frequency.

GATE THEORY

If the electric signals that lead to the brain information about pain are stopped, also the pain perception is eliminated. If, for instance, we hit our head into an object the first thing we do is massaging the area affected by the trauma. In this way we stimulate the receptors related to touch and pressure. TENS in continuous mode and in frequency modulation can be used to generate signals similar to the ones of touch and pressure. If their intensity is enough, their priority is so high that it prevails on the pain signals. Once the priority is gained, the gate related to the sensory signals is opened and the pain one is closed, impeding in this way the passage of these signals to the brain.

ENDORPHIN SECRETION

When a nervous signal proceeds from the pain area to the brain, it spreads through a chain of connections joined together called synapse. The synapse can be seen as the space between the end of a nerve and the start of another. When an electric signal reaches the end of a nerve, it produces

some substances called neurotransmitters that pass through the synapse and activate the start of the next nerve. This process repeats for all the length needed to the signal to reach the brain. The opioids involved in the pain reduction have the task to insinuate in the synapse space and impede the neurotransmitter propagation. In this way a chemical block of the pain signals occurs. The endorphins are opioids naturally produced by the body to tackle the pain and they can act both on the marrow and on the brain, in this way they are effective analgesics. Tens can increase the natural production of endorphins and, thereby, they act decreasing the pain perception.

INCONTINENCE PROGRAMS

The urology programs require the use of specific endovaginal and endorectal electrode probes, certified according to the Directive for Medical Devices 93/42/EEC. These are bipolar probes with a 2-mm female adapter which attaches to 2-mm male cables.

WARNINGS

Urological electrostimulation is a medical application, which must be carried out under medical supervision.

USE

To correctly use the probe electrode, follow the instructions provided by the manufacturer and given by the physician.

MAINTENANCE

For cleaning, sterilization and disinfection, refer to the manufacturer's instructions.

SUGGESTIONS

In case of deterioration of the probe electrode, replace it immediately.

GENERAL NOTES ON ELECTRODE POSITIONING

The correct electrode positioning and the correct choice of their size are critical aspects for the electrostimulation to be effective.

To choose the size of the electrodes and as for their positioning it is necessary to refer to the images at the end of this manual.

NOTE For all the programs that cause an important muscle contraction (such as, for example, strength, hypertrophy, toning and firming programs ...) it is important to position the electrode on the muscle motor point that is the most sensitive point to stimulation.

If the electrode is not positioned exactly on the motor point, the contraction could be small and/or annoying. In this case it is necessary for the positive electrode to be moved a few millimeters up to feel an effective and comfortable contraction.

PROGRAM LIST

SPORT

DEMO	P1
Motor point pen	P2
Warm-up lower limbs	P3
Pre-competition warm-up lower limbs	P4
Maximum strength lower limbs	P5
Endurance strength lower limbs	P6
Explosive strength lower limbs	P7
Reactivity lower limbs	P8
Capillarization lower limbs	P9
Active recovery lower limbs	P10
Muscle cool down lower limbs	P11
Endurance lower limbs	P12
Decontracting lower limbs	P13
Abdominal maximum strength	P14
Abdominal endurance strength	P15
Abdominal active recovery	P16
Abdominal cool down	P17
Abdominal capillaries	P18
Abdominal decontracting	P19
Warm-up upper limbs	P20
Pre-competition warm-up upper limbs	P21
Maximum strength upper limbs	P22
Endurance strength upper limbs	P23
Explosive strength upper limbs	P24
Reactivity upper limbs	P25
Capillarization upper limbs	P26
Active recovery upper limbs	P27
Muscle cool down upper limbs	P28
Endurance upper limbs	P29
Decontracting upper limbs	P30
TOTAL	30

FITNESS & BEAUTY

Leg anti-cellulite lipolysis	P31
Hip-gluteus anti-cellulite lipolysis	P32
Abdominal area anti-cellulite lipolysis	P33
Upper limb anti-cellulite lipolysis	P34
Leg drainage	P35
Hip-gluteus drainage	P36
Abdominal area drainage	P37
Upper limb drainage	P38
Leg firming	P39
Hip-gluteus firming	P40
Abdominal area firming	P41
Upper limb firming	P42
Leg toning	P43
Dorsal toning	P44
Hip-gluteus toning	P45
Abdominal area toning	P46
Upper limb toning	P47
Leg sculpting	P48
Hip-gluteus sculpting	P49
Abdominal area sculpting	P50
Pectoral sculpting	P51
Leg vascularization	P52
Upper limb vascularization	P53
Woman abdominal area toning	P63
Man abdominal area definition	P64
Man pectoral definition	P65
Woman leg mass building	P66
Woman upper limb mass building	P67
Swollen upper limbs	P68
Swollen legs	P69
Leg toning massage	P70
Face beauty 1	P71
Face beauty 2	P72
Breast toning 1	P73
Breast toning 2	P74
TOTAL	35

SPECIAL SPORTS

Downhill skiing	P54
Soccer strength endurance	P55
Cross-country skiing	P56
Running: lower-limb endurance	P57
Swimming: upper-limb endurance	P58
Cycling: aerobic workout	P59
Swimming: trunk endurance	P60
Triathlon: maintenance of physical condition	P61
Triathlon: recovery of physical condition (after a period of inactivity)	P62
TOTAL	9

TENS

Conventional antalgic Tens*	P75
Endorphinic Tens*	P76
Rotator cuff tendinitis*	P77
Muscle pain*	P78
Knee osteoarthritis*	P79
Post-surgery pain*	P80
Muscle injury*	P81
Cervical pain*	P82
Trapezius pain*	P83
Sciatica*	P84
Lumbago*	P85
Epicondylitis*	P86
Scapulo-humeral periarthritis (shoulder)*	P87
Carpal tunnel*	P88
Osteoarthritis*	P89
Menstrual pain*	P90
TOTAL	16

PREVENTION

Ankle prevention	P91
Knee prevention	P92
Quadriceps muscle prevention	P93
Lower limb reathletization	P94
Shoulder reathletization	P95
Elbow prevention	P96
TOTAL	6

REHAB

Quadriceps atrophy with prosthesis*	P97
Shoulder subluxation prevention*	P98
Multiple sclerosis - Ankle flexor extensor mm*	P99
TOTAL	3

INCONTINENCE

Mixed Incontinence*	P100
TOTAL	1

*Because of the presence of clinical programs, this product is a medical device.

The devices are certified by the Kiwa Cermet Italia S.p.A. Body n. 0476 in compliance with the 93/42/EEC EU directive for medical devices. The certification refers only to medical, rehabilitation and functional recovery applications.

THE BODY POSITION DURING THE STIMULATION

The body position during the electrostimulation session depends on the body part involved and on the program type that is being carried out. During the treatment execution with high intensities, we suggest blocking the limbs in order to work in isometry. For instance, if you want to treat the quadriceps with a strength program, we suggest carrying out the treatment while sitting with the foot blocked, in order to avoid an involuntary leg extension during the contraction phase.

For all the programs that do not imply high intensity (massages, decontracting, drainage programs) the body position is not important, as long as it is comfortable.

ELECTRODE POSITIONING FOR TENS PROGRAMS

In the following pages of this manual you can find the images with the correct electrode positioning for tens treatments. If the localization of your pain is not included in the images represented, you can position the electrodes by forming a "square" on the painful area. Here's an example.

TREATMENT DESCRIPTIONS

Following is a short description of the programs provided by the electrostimulator.

SPORT

WARM-UP: Increases muscular temperature for better training preparation. Increases the blood flow to and metabolism of the muscle. Recommended before training sessions or competitions.

Intensity: low or intermediate.

MAXIMUM STRENGTH: Improves maximal force capability and increases muscular mass. Allows training with maximum load and reduced trauma risk as compared to classical training. Recommended for activities requiring great physical strength. Also useful during the early preparation period as it trains the muscles to receive significant stress.

Intensity: maximum, at pain threshold limit. During active Rest, adjust the intensity to obtain comfortable contractions.

ENDURANCE STRENGTH: Improves the ability to display a high level of force for an extended period of time and to resist toxin accumulation, thus delaying muscular fatigue. The program involves a long series of lengthy contractions, alternated with brief active rest. Recommended for sports requiring intense muscular work for an extended period of time.

Intensity: intermediate or maximum. Active Rest should induce muscular vibrations.


EXPLOSIVE STRENGTH: The Explosive Force program increases the ability to rapidly achieve maximum levels of force. The objective is to recruit the maximum number of muscle fibers in the shortest time possible. It is recommended to begin with a Maximal Force program to prepare the muscle for the Explosive Force program. Recommended in all sports requiring maximal force under natural load conditions, such as sprints, ball sports, and alpine skiing. Intensity: maximum muscular contraction at pain threshold limit. Active Rest should induce muscular vibrations.

DECONTRACTING: Enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programs can be used whenever muscular relaxation is needed. Intensity: low to intermediate.

CAPILLARIZATION: Increases the blood flow to the muscle and improves muscle resistance and recovery. An increased arterial capacity keeps the primary and secondary capillary network active, thus improving the tissue oxygenation system. This delays fatigue during very intense work. Capillarization programs are particularly recommended for sports requiring resistance force and endurance.

Intensity: intermediate.

POST-COMPETITION/TRAINING ACTIVE RECOVERY: Improves recovery from muscular fatigue after training or competitions and enables easier training resumption. The stimulation varies from muscular vibrations to light contractions. This allows greater circulation and reduced muscular

hardening. It also alleviates pain thanks to the release of endorphins. This program is valid for all sports as it improves and speeds up the recovery of maximum muscular capacity. This program is recommended in the hours immediately after training or competitions.

Intensity: low to intermediate.

BEAUTY

DRAINAGE: The coordinated stimulation obtained by using all channels of the electrostimulator produces an in-depth muscular massage effect. The massage activates venous and lymphatic systems by means of a mechanical effect. The purpose is to improve blood circulation and enable the toxin elimination. Intensity: low to intermediate.

LIPOLYSIS: The low frequency of this program positively affects the metabolism and increases the elimination of toxins and excess fluid. This is a method for reshaping critical areas, modifying the metabolism of lipocytes. Improves blood circulation. Helps to eliminate the orange peel skin. Intensity: low to intermediate.

VASCULARIZATION: Enables effective muscle relaxation thanks to the pleasant contractions that increase blood circulation. These programs can be used whenever muscular relaxation is needed. Intensity: intermediate.

FIRMING: This program is suitable for sedentary individuals with flaccid muscles that require reactivation of muscular tone. It prepares the muscles for more intense stimulations. It can be performed every day on large muscle groups such as thighs, abdominal and dorsal muscles. Intensity: intermediate to high.

TONING: The toning program allows to develop and maintain muscular tone. This type of treatment may be combined with aerobic activity and is suitable for users of any level. The program is especially directed at all those who, mindful of their own physical condition, wish to integrate their own regular training. The program can be used up to three times a week. Intensity: intermediate to high.

SCULPTING: The sculpting program is a merely aesthetic treatment and is ideal for those who wish to better define their already toned muscle system. The work performed is very intense, so slight muscle soreness may be experienced the day after the treatment. However, this effect may be avoided by performing stretching exercises before and after the session or a winding-down program. The program is suitable for those who wish to achieve muscular definition in those areas where traditional muscular development exercises are insufficient. Intensity: intermediate to high.

MASS BUILDING: A specific program for those who wish to increase muscle mass and volume. Intensity: maximum muscular contraction at pain threshold limit.

A muscle warm-up program should be carried out before performing any of the following programs: Firming, Toning, Sculpting and Mass Development.

SPECIAL SPORTS

RUNNING: it improves the ability to maintain the effort for an extended period of time.

The training programs for running are intended to stimulate slow fibers with very long contraction periods.

SWIMMING: it permits to train the swimmer's musculature with the aim to improve the effort endurance. Therefore, the athlete's musculature could become accustomed to the training sessions with extended stimulations, obtaining a lot of benefits during the trainings in water.

CYCLING: it improves the ability to maintain the effort for an extended period of time. The training programs for the biker are intended to stimulate slow fibers with very long contraction periods.

TRIATHLON: these programs permit to train the musculature to bear at best the prolonged and intense efforts, by improving the ability of the muscles to consume oxygen.

SKIING: These programs improve the muscle aerobic endurance and therefore the ability to bear an extended physical effort.

SOCCER: This program is especially intended to improve the specific conditional abilities of soccer and especially strength and endurance. The program effects are to increase the strength and the endurance in carrying out sprints and goal kicks.

TENS

Uses mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. This analgesic effect occurs during program execution and can last for a few hours. However, its duration may vary. The TENS programs also improves circulation. Place the electrodes on the painful area. These programs can be used in any situation requiring pain relief. Caution! TENS programs alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice.

TENS (Antalgic-Endorphinic): Allows pain relief and effective muscle relaxation. TENS programs use mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. Muscular relaxation improves blood circulation. It can be used in any situation requiring both pain relief and muscle relaxation. Caution! TENS programs alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice.

Intensity: ANALGICTENS: tingling but not painful sensation. Should not induce muscular contractions. ENDORPHINIC TENS: visible muscular vibrations.

TENS (Cervical pains, Lumbago, Periarthritis, Epicondylitis, Sciatica,...): Uses mechanisms of the nervous system to alleviate both acute and chronic pain. Impulses travel along the nerves to block the pain propagation paths. This program combines stimulation with light contractions that stimulate the release of endorphins. Caution! TENS programs alleviate the pain but do not eliminate the cause. If the pain persists, seek medical advice. Intensity: tingling but not painful sensation.

REHAB

Resumption of physical activity after a period of interruption. After a trauma or other interruptions of physical activity, muscles may lose their tone and mass. In these cases, normal muscular capacity must be regained as soon as possible. The program enables training of specific muscles concerned.

It can be used in any situation requiring activity resumption. Intensity: low to intermediate.

PREVENTION PROGRAMS

The stimulation is aimed at neuromuscular facilitation, which especially improves the mobilization of the injured articulation. Prepares the limb to resume its normal functionality. It is recommended to use this program during the first phases of the recovery of muscular trophism. Intensity: intermediate.


INCONTINENCE

Programs recommended for the stimulation of perianal muscles.


These programs are recommended in case of urological pain or incontinence.

Programs to be used with probes in compliance with the law.


ELECTRODE PLACEMENT - TENS


BRACHIAL BICEPS


BRACHIAL BICEPS


FLEXOR CARPI


EXTENSOR CARPI

DELTOID


PECTORAL


PECTORAL


TRAPEZIUS


TRAPEZIUS


DORSAL MUSCLES


DORSAL MUSCLES


SUBSPINATUS


LUMBAR MUSCLES


LUMBAR/DORSAL MUSCLES


ABDOMINALS

ABDOMINALS


ABDOMINALS


RECTUS ABDOMINI


GLUTEUS


GLUTEUS


FEMORAL BICEPS


ADDUCTORS


RECTUS FEMORIS

QUADRICEPS

GASTROCNEMIUS


GASTROCNEMIUS

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

ELECTRODE PLACEMENT - TENS


SHOULDER PAIN


ELBOW PAIN


CERVICAL PAIN


LUMBAR PAIN


LUMBAR PAIN


SCIATICA

KNEE PAIN

KNEE PAIN

KNEE PAIN

FREQUENTLY ASKED QUESTIONS

Why, when I increase the current intensity above 1 mA, it returns immediately to 0 mA?

Before starting any program, check that the cables are connected and that the electrodes are placed on the area to be treated. In addition, be sure that the electrodes are not worn, because this could inhibit their ability to conduct current.

What kind of electrodes should be used for electrostimulation?

Use self-adhesive electrodes, which are practical and improve the quality of stimulation. If used with care, they will last for 25-30 applications. The electrodes should be replaced when they have no longer a good connection with the skin.

Where should the electrodes be placed?

In the back of this manual, there is a comprehensive electrode placement guide (it is not necessary to respect the polarities indicated). You can follow these instructions. To verify the correct placement of the electrodes, use the special Find Motor Point Pen program or follow this empirical method: place the electrodes as indicated in the pictures in the back of this manual; start the stimulation; with your hand, move the electrode by sliding it along the muscle without removing it from the skin. You will notice a change in contraction according to the different positions. Once you locate the point where the stimulation is greatest, decrease the channel intensity to zero (0,0 mA), replace the electrode and increase gradually the intensity.

Use of Y cables. This permits to use more electrodes on the same channel.

This permits, for instance, to stimulate vastus medialis and vastus lateralis of the quadriceps with one single channel. They are not recommended for medical applications.

Does the power decrease using Y cables?

The power intensity for each channel does not vary. However, when Y cables are used to split one single channel in two, the current is distributed on a wider muscle area, therefore contraction will be less pronounced. Increase the intensity to obtain the same contraction level.

Can electrostimulation hurt me?

It is very unlikely that electrostimulation damages muscles. In this case it is important to increase intensity gradually while observing the muscle behavior and avoiding to keep the limb completely outstretched. In case of doubt, please contact a specialist.

Is it possible to use the electrostimulator during menstruation cycle?

Electrostimulation may interfere in some way with menstruation, causing anticipation, delay, accentuation or reduction of the cycle; however, these effects are subjective and highly variable. It is recommended to avoid treatments in the abdominal zone during menstruation cycle and immediately before or after it.

Is it possible to use the electrostimulator during lactation?

Until now, no collateral effects regarding lactation have been observed. Nonetheless, during lactation, it is recommended not to stimulate the thoracic region.

Are dermatological diseases (e.g. psoriasis, urticaria) contraindications for electrostimulation?

Yes, do not treat areas affected by dermatological diseases.

When are the first results visible?

Aesthetic results of electrostimulation are always subjective. For Toning program, with a regular program of 3-4 sessions per week, a noticeable result may be observed after 15 days. For Lipolysis and Drainage programs, 40 days of treatment are necessary. Results are obtained more quickly if treatments are combined with good physical activity and a correct life style.

How many sessions can I weekly perform?

For fitness and aesthetics programs, the number of sessions depends on the type of treatment: 3-4 sessions per week on alternate days are suggested for toning; daily treatments are permitted for Lipolysis and Drainage programs.

WARRANTY

The device is guaranteed to the first user for a period of twenty-four (24) months from the date of purchase against defects in materials or of the manufacturing (if the user uses the device for professional purposes, is twelve, 12, months), provided that it is used properly and maintained under normal operating conditions.

Warranty coverage is limited in the following cases:

- Six (6) months for supplied accessories such as batteries, chargers, power supplies, cables, etc.
- No warranty for extendable accessories and materials such as, for example, electrodes, etc.

The warranty is valid and enforceable in the country where the product was purchased. In the event that the product is purchased in any country of the European Community, the warranty is still valid in all its countries.

In order to take advantage of the warranty service, the user must comply with the following warranty clauses:

1. The products, and all accessories, to be repaired must be sent by and at the expenses of the customer in their original packages.
2. The product's warranty is subject to the production of a fiscal document (fiscal receipt, receipted bill or sale invoice), attesting the product's purchase date.
3. The repair work shall have no effect on the original warranty expiry date and shall neither renew nor extend it.
4. If, during the repair work, no defects are found, the costs related to inspection times shall in any case be charged.
5. The warranty becomes void where the fault has been caused by: impacts, falls, erroneous or improper use of the product, use of non-original power supply unit or charger, accidental events, alteration, replacement/detachment of the warranty seals and/or tampering with the product. The warranty does not cover damages caused during transportation when unsuitable packages are used.
6. The warranty does not cover the inability to use the product, other incidental or consequent costs or other expenses incurred by the purchaser. The features and dimensions reported in this manual are not binding. The manufacturer reserves the right to make changes without previous notice.

PROstim2

if you like it, share it...

[@rehab_medic](https://twitter.com/rehab_medic)

www.facebook.com/rehabmedic

[rehabmedic](https://www.pinterest.com/rehabmedic)

www.rehabmedic.com