# Primado2

Control Unit and Foot Controls

## OPERATION MANUAL

This product is a surgical system used for bone operations (removal, cutting, perforation and formation), or cutting implants. We recommend that prior to use, you carefully read this Operation Manual regarding "Precaution for Handling and Operation", "Connection method", "Operation method", "Check before Operation" and "Periodical Maintenance Check" so that you can continue using the Attachments in the future. In addition, keep this Operation Manual in a place where a user can refer to it at any given time.

\* In this operation manual Bur, Drill, Blade and Disk are collectively called Cutting Accessory(ies).

\* For reprocessing of the Control Unit, Foot Control, Motor and Attachment or operation method of the Attachment, refer to each manual (see Table-1).

Table-1

| Operation Manual Title                | Operation Manual No. |
|---------------------------------------|----------------------|
| Primado2 REPROCESSING MANUAL          | OM-SZ0911MA          |
| Primado2 OPERATION MANUAL Attachments | OM-SH0913MA          |

## CONTENTS

|                                              |    |
|----------------------------------------------|----|
| Precautions for Handling and Operation ..... | 1  |
| 1. Component Names .....                     | 4  |
| 2. Part Names .....                          | 5  |
| 3. How to Connect Each Component .....       | 12 |
| 4. Check before Operation .....              | 16 |
| 5. Operation .....                           | 16 |
| 6. Option Set Up .....                       | 20 |
| 7. Fuse Replacement .....                    | 25 |
| 8. Periodical Maintenance Checks .....       | 25 |
| 9. Reprocessing .....                        | 25 |
| 10. Error Codes .....                        | 26 |
| 11. Troubleshooting .....                    | 29 |
| 12. Specifications .....                     | 31 |
| 13. Disposing of Product .....               | 32 |
| 14. Warranty .....                           | 32 |
| 15. Symbols .....                            | 32 |
| 16. EMC Information .....                    | 33 |

## Precautions for Handling and Operation

- Read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

| Class                                                                                              | Degree of Risk                                                                                                                      |
|----------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
|  <b>DANGER</b>  | Hazard that could result in personal death or serious injury if the safety instructions are not correctly followed.                 |
|  <b>WARNING</b> | Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.           |
|  <b>CAUTION</b> | Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed. |
| <b>NOTICE</b>                                                                                      | General product specification information highlighted to avoid product malfunction and performance reduction.                       |

Before use, read carefully the Operation Manual of the Motor and the Attachments.

## User and Intended Purpose

Primado2 is intended for the following applications:

Cutting, drilling, removal, and shaping of bones and implants in the fields of Neuro, Spine, Orthopedic, ENT, Oral maxillofacial, and Plastic and reconstructive surgery.

### DANGER

- Do not use this system in a room where there is danger of explosion or in the vicinity of flammable substances. It is especially important not to use this system on a patient to whom a flammable anesthetic agent (dinitrogen monoxide) has been administered.
- Should a short circuit occur inside this product, an electric shock or fire may result if the product is not connected to properly grounded electric supply.
- After use, turn off the Power Switch. If the unit is not being used for a prolonged period, disconnect the Power Cord Plug from the outlet. Otherwise dust or deteriorated insulation might cause electric shock or fire.
- Never disassemble or modify the product. Personal injury, electric shock or fire may result.

### WARNING

- This is a medical system; do not allow anybody to use it other than surgeons who are familiar with the procedures. Use this system in accordance with its intended purpose and method of use described in this Operation Manual.
- Do not subject this unit to strong impacts, such as dropping it on the floor. Personal injury, abnormal heating and fire may result caused by the internal damage.
- This product may be affected when used in rooms with electromagnetic interference and portable or mobile RF communication equipment. Do not use this product around any equipment that generates electromagnetic waves such as an ultrasonic generator and electrical scalpel. When using equipment generating electromagnetic waves, turn off the power to the Control Unit.
- Do not directly spray detergent/liquid on the LCD Touch Panel of the Control Unit when cleaning the Control Unit. Otherwise, it may cause failure.
- Except for the Foot Control, this system is not designed to be water-proof. Handle the Control Unit and Motor with care so that they will not be exposed to water or sterile water (saline solution). An electric shock may result due to a short circuit.
- The Irrigation Tube and Cutting Accessory are for single use only. Do not reuse or reprocess it. If you reuse or reprocess it, the following risks may occur:
  - Infection.
  - Personal injury due to damaged Cutting Accessory.
  - System defect due to damaged/deformed Cutting Accessory.
- Irrigate appropriately while cutting in order to prevent necrosis of bones and tissues.
- After purchase or repair of the Motor (optional), be sure to clean and sterilize the Motor before use (refer to OM-SZ0911MA Primado2 REPROCESSING MANUAL).
- After purchase or repair of the Attachment (optional), be sure to clean, lubricate and sterilize the Attachment before use (refer to OM-SZ0911MA Primado2 REPROCESSING MANUAL).
- Do not lubricate the Motor. It may result in an abnormal heating and failure.

### CAUTION

- When operating the product, always consider the safety of the patient.
- Users shall be responsible for any judgement that relates to the application of this product to a patient.
- Users are responsible for the operational control, maintenance and continual inspection of this System.
- Be sure to read this Operation Manual and the Operation Manuals for "Primado2 Motor" and each Attachment in order to become completely familiar with the function of each part before use.
- Connect the product to an electrical outlet for medical use or a properly grounded outlet.
- Be sure to perform inspection by operating the product with the Motor, Attachment and Cutting Accessory mounted before use (see "4 Check before Operation"). If any abnormal condition is detected, stop using the product and contact your Authorized NSK Dealer.
- Use these products after performing "4 Check before Operation". When you feel an abnormality such as vibration, overheating or noise while in use, immediately stop using them and perform inspection (see "11 Troubleshooting"). If the issue does not improve after troubleshooting, contact your Authorized NSK Dealer.
- Do not bend or kink the Irrigation Tube while the Irrigation Pump is operating. Damage or disconnection of the Irrigation Tube may result.
- Wear Personal Protective Equipment (PPE) during operation.
- It is recommended that a spare product be kept on hand in case of a breakdown during surgery.
- When the Control Unit and Foot Control are stained with blood, turn off the power and wipe the stains off using a lint-free cloth moistened with neutral detergent for medical use. Then, remove the detergent using a lint-free cloth moistened with distilled, deionized or demineralized water.


**⚠ CAUTION**

- If chemicals, solutions or sterilizing liquid is left on the surface of the Control Unit and/or Foot Control, wipe it away immediately. If it is left, it may lead to discoloration, deformation, rust and/or corrosion.
- Do not use the following fluids to wipe or clean the product: strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, and solvents such as benzene or thinner.
- Do not disconnect the Motor Cord from the Motor. It may cause failure.
- Dispose of the used Irrigation Tube and Cutting Accessories as medical waste (see "13 Disposing of Product").
- Connect or disconnect the cords or Attachments only after turning off the Power Switch. If the Power Switch is ON, an unintentional movement may occur, which could result in injury to users.
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information.
- The use of accessories, Motors and Attachments other than those specified, with the exception of accessories, Motors and Attachments sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- This product should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, observe this product and other equipment in order to verify a normal operation prior to use.
- Regarding the operation time of the Motor, follow specification table in step "12. Specifications". Continuous use of the Motor for a prolonged period might cause overheating, which can cause an accident.
- This device is for indoor use only.
- Perform regular function and maintenance checks for the devices and parts (see "8 Periodical Maintenance Checks").
- If the product is not used for a long period of time, check if it functions correctly and safely before use.
- The temperature at the tip of the Attachment may reach up to 60°C depending on the environment.
- Check the correct connection of the following before use:
  - Connection of the Motor Cord Plug and Foot Control Plug to the Control Unit.
  - Connection of the Power Cord Plug to the Control Unit.
  - Connection of the Power Cord Plug to the outlet.
- Only use the Irrigation Tube (Order Code: Y1003915) recommended by NSK. Otherwise, it may cause failure.
- Disconnect the Power Cord Plug, Motor Cord Plug and Foot Control Plug from the Control Unit when the product is not used. If excessive impact is applied, such as hitting the connections without disconnecting them, failure of the Control Unit, Power Cord Plug, Motor Cord Plug or Foot Control Plug may occur.
- Please report any serious incident that may occur in connection with this device and its accessories to the manufacturer as well as to your national authority.

**NOTICE**

- After use, turn OFF the Power Switch. If the unit is not being used for a prolonged period, disconnect the Power Plug from the outlet. If the unit is left plugged in and unattended for too long dust or deteriorated insulation might cause electric shock or fire.
- The user shall be responsible for maintenance and operation.
- Only use irrigation tubes recommended by NSK.

# 1 Component Names

| No. | Model       | Name of Component                 |
|-----|-------------|-----------------------------------|
| 1   | P200-CU-120 | Primado2 Control Unit (120V/230V) |
| 2   | P200-CU-230 | Power Cord                        |
| 3   |             | Irrigation Pole                   |
| 4   | FC-74       | Primado2 Multi Foot Control       |
| 5   | FC-73       | Primado2 Single Foot Control      |
| 6   | P200-HMH    | Primado2 High Torque Motor        |
| 7   | P200-HMH-HS | Primado2 High Torque Motor/HS     |
| 8   | P200-SMH-S  | Primado2 Slim Motor Short         |
| 9   | P200-SMH    | Primado2 Slim Motor               |
| 10  | P200-SMH-HS | Primado2 Slim Motor/HS            |
| 11  | P200-BMH    | Primado2 Micro Bone Saw Motor     |
| 12  | P200-BMH-HS | Primado2 Micro Bone Saw Motor/HS  |
| 13  | P200-WPD    | Primado2 Wire Pin Driver          |

These components are sold separately, according to the client's specific needs.

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





## 2-3 Motor

P200-HMH, P200-SMH-S, P200-SMH, P200-BMH

P200-HMH-HS, P200-SMH-HS, P200-BMH-HS (Motor with Hand Control)

1) Hand Switch\*


2) Adjustable Lever\*

3) Safety Lock\*

4) Motor Cord Plug

\*Motor with Hand Control only.

P200-WPD

1) Trigger (Upper) (A)

2) Safety Lock


3) Trigger (Lower) (B)

4) Mounting/Removal Button

5) Motor Cord Plug

## 2-4 LCD Touch Panel: Operation Key





HOME screen


**NOTICE** • Whenever the icon  is touched, the HOME screen appears, regardless of the current screen.

## 2-5 Buttons &amp; Icons


|             |            | Icon | Name                              | Descriptions                                                                                                                                |
|-------------|------------|------|-----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| HOME Screen | Main Panel |      | Motor Select Button A             | Selects the Motor which is connected to Motor Connector A.<br>If the Motor is not connected, "Disconnection" is displayed.                  |
|             |            |      | Motor Select Button B             | Selects the Motor which is connected to Motor Connector B.<br>If the Motor is not connected, "Disconnection" is displayed.                  |
|             |            |      | Speed Adjustment Area             | Rotational speed can be adjusted by the Track Bar and UP/DOWN Button.                                                                       |
|             |            |      | Irrigation Volume Adjustment Area | The irrigation volume can be adjusted by the Track Bar and UP/DOWN Button.                                                                  |
|             |            |      | Irrigation ON/OFF Button          | Irrigation is OFF by default. Irrigation is turned ON by touching the button.                                                               |
|             |            |      | Track Bar                         | The rotational speed or irrigation volume can be adjusted by moving the Track Bar to the left or right.                                     |
|             |            |      | UP Button                         | Increases the rotational speed or irrigation volume.                                                                                        |
|             |            |      | DOWN Button                       | Decreases the rotational speed or irrigation volume.                                                                                        |
|             |            |      | Rotation Mode Select Button       | Changes the Rotation Mode of the Motor connected to Motor A or B between Forward, Reverse and Oscillate (P200-WPD only).                    |
|             |            |      | "Rotating" display                | Shows that the Motor is rotating.                                                                                                           |
|             |            |      | Priming Button                    | Whether the Irrigation Mode at ON or OFF, as long as the Priming Button is kept pushed, the irrigation volume is at 100% capacity.          |
|             | Sub Panel  |      | Foot Control Panel                | Displays the function set up by the individual buttons of the Foot Control Panel. If the Foot Control is not connected, only  is displayed. |
|             |            |      | Foot Control ON/OFF Button        | Default setting: OFF. The Foot Control is turned ON by pressing the button (see Figure above).                                              |
|             | Mini Panel |      | Option Set Up Button              | Displays the "Option Set Up" screen.                                                                                                        |

|                                     |            | Icon                                                                                | Name                                                                       | Descriptions                                                                                     |
|-------------------------------------|------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| Option Set Up Screen                | Main Panel |    | Control Unit Set Up Button                                                 | Displays the “Control Unit Set Up” screen.                                                       |
|                                     |            |    | Motor Set Up Button                                                        | Displays the “Motor Set Up” screen.                                                              |
|                                     |            |    | Attachment Set Up Button                                                   | Displays the “Attachment Set Up” screen.                                                         |
|                                     |            |    | Foot Control Set Up Button                                                 | Displays “Foot Control Set Up” screen.                                                           |
|                                     | Sub Panel  |    | Back Button                                                                | Displays previous screen.                                                                        |
|                                     | Mini Panel |    | HOME Button                                                                | Displays the HOME screen regardless of the current screen.                                       |
| Control Unit Set Up Screen          | Main Panel |    | Background Brightness Adjustment Area                                      | By using the Track Bar and the UP/DOWN Button, the brightness of the back light can be adjusted. |
|                                     |            |    | Volume Adjustment Area                                                     | By using the Track Bar and the UP/DOWN Button, the volume can be adjusted.                       |
|                                     |            |   | Language Set Up Button                                                     | Displays the “Language Set Up” screen.                                                           |
|                                     |            |  | Information Button                                                         | Shows the version of the hardware and software of the Control Unit.                              |
| Language Set Up Screen              | Main Panel |  | Language Select Button (English)                                           | The language in which message alerts are to be displayed.                                        |
|                                     |            |  | Language Select Button (German)                                            |                                                                                                  |
|                                     |            |  | Language Select Button (French)                                            |                                                                                                  |
|                                     |            |  | Language Select Button (Spanish)                                           |                                                                                                  |
|                                     |            |  | Language Select Button (Italian)                                           |                                                                                                  |
|                                     |            |  | Language Select Button (Japanese)                                          |                                                                                                  |
| Language Set Up Confirmation Screen | Main Panel |  | “Yes” Button<br>English / German<br>French / Spanish / Italian<br>Japanese | The language selected in the Language Set Up screen is set as the main language.                 |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  | “No” Button<br>English / Italian / Spanish<br>German / French<br>Japanese  | Cancels the language selected in the Language Set Up screen.                                     |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  |                                                                            |                                                                                                  |
|                                     |            |  |                                                                            |                                                                                                  |

|                     |            | Icon | Name                                    | Descriptions                                                                                                                |
|---------------------|------------|------|-----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Motor Set Up Screen | Main Panel |      | Acceleration Area                       | The acceleration time for the Motor to reach its maximum rotational speed can be adjusted by the UP/DOWN Button (10%~100%). |
|                     |            |      | Brake Area                              | The braking time for the Motor to stop can be adjusted by the UP/DOWN Button (10%~100%).                                    |
|                     |            |      | Torque Area                             | The torque of the Motor can be adjusted by the UP/DOWN Button (30%~100%).                                                   |
|                     |            |      | Oscillating Turns Area                  | The Oscillating Turns can be adjusted by the UP/DOWN Button (1~5) (P200-WPD only).                                          |
|                     |            |      | Factory Default Button                  | Returns the settings back to the original factory setting made before delivery.                                             |
|                     |            |      | Motor Select Button                     | Changes the connection of the Motor between Motor Connector A and B.                                                        |
|                     | Main Panel |      | Attachment Select Button (Direct Drive) | Indicates the rotational speed suited for the Constant Speed Handpiece.                                                     |
|                     |            |      | Attachment Select Button (Perforator)   | Indicates the rotational speed suited for the Perforator Handpiece.                                                         |
|                     |            |      | Attachment Select Button (High Speed)   | Indicates the rotational speed suited for the High Speed Handpiece.                                                         |
|                     |            |      | Attachment Select Button (Super Slim)   | Indicates the rotational speed suited for the Super Slim Handpiece.                                                         |
|                     |            |      | Attachment Select Button (Contra 4:1)   | Indicates the rotational speed suited for the 4:1 Speed-Reduction Contra Angle Handpiece.                                   |
|                     |            |      | Attachment Select Button (Contra 16:1)  | Indicates the rotational speed suited for the 16:1 Speed-Reduction Contra Angle Handpiece.                                  |
|                     |            |      | Attachment Select Button (Contra 64:1)  | Indicates the rotational speed suited for the 64:1 Speed-Reduction Contra Angle Handpiece.                                  |
|                     |            |      | Attachment Select Button (Drill)        | Indicates the rotational speed suited for the Constant Speed Attachment.                                                    |
|                     |            |      | Attachment Select Button (Bone Saw)     | Indicates the rotational speed suited for the Bone Saw Attachment.                                                          |
|                     |            |      | Attachment Select Button (Reamer)       | Indicates the rotational speed suited for the Reamer Attachment..                                                           |

|                            |            | Icon                                                                                | Name                                                          | Descriptions                                                                                                                                                                  |
|----------------------------|------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Foot Control Set Up Screen | Main Panel |    | Foot Control Button (Blue) (Only for the Multi Foot Control)  | Displays the “Foot Control Set Up” screen.                                                                                                                                    |
|                            |            |    | Foot Control Button (Silver)                                  |                                                                                                                                                                               |
|                            |            |    | Foot Control Button (Green) (Only for the Multi Foot Control) |                                                                                                                                                                               |
|                            |            |    | Foot Control Pedal                                            |                                                                                                                                                                               |
|                            |            |    | ON/OFF Select                                                 | Turns ON/OFF the operation of the Foot Control by the Foot Control. When turned OFF, the operation of the Foot Control becomes unavailable.                                   |
|                            |            |    | Motor Select                                                  | Sets up the Motor connected to Motor Connector A or B so that it can be switched over to another Connector by the Foot Control.                                               |
|                            |            |    | Speed Select                                                  | Sets up the rotational speed of the Motor connected to Motor Connector A or B so that it can be switched over to another speed (30%→Default, Value→100%) by the Foot Control. |
|                            |            |  | Rotation Mode Select                                          | Sets up the Rotation Mode of the Motor connected to Motor Connector A or B so that it can be switched over to another by the Foot Control.                                    |
|                            |            |  | Irrigation ON/OFF Select                                      | Sets up the Irrigation so that it can be turned ON/OFF by the Foot Control.                                                                                                   |
|                            |            |  | Flush                                                         | Makes 100% irrigation available as long as this button is pressed.                                                                                                            |
|                            |            |  | Disable                                                       | Sets the status as no function.                                                                                                                                               |
|                            |            |  | “Foot Control Button” Display (Only for the Multiple Model)   | Highlights the Function Button currently used.                                                                                                                                |
| Error Code Display Screen  | Sub Panel  |  | Error Code Button                                             | When an error has occurred simultaneously, a blinking light announces the error. Touch the Error Code Button to display the error code for that error.                        |
|                            |            |  | Motor                                                         | When an error has occurred in the Motor, this indicates which Connector is connected to the Motor.                                                                            |
|                            | Mini Panel |  | Close Button                                                  | Closes the Error display.                                                                                                                                                     |

## 2-6 Flow Chart of the LCD Touch Panel

### 3 How to Connect Each Component


#### CAUTION

- Do not remove the Power Cord or Motor Cord by holding the cord portion or failure might occur. Hold the plug.
- Connect or disconnect the Power Cord or Motor Cord only after turning OFF the Power Switch. Otherwise, personal injury may result caused by an unexpected rotation.

#### 3-1 Connection of the Power Cord

- 1) Connect the Power Cord into the end into the Inlet Box located at the back of the Control Unit.

- 2) Plug the Power Cord to the medical-use outlet.

Power Supply plug type to be used in North America area.

Plug Type NEMA 5-15P  
(Hospital Grade Type)

#### CAUTION


- Be sure to connect the Power Cord into the medical-use outlet.
- Be sure to confirm that the grounding is established. Should a short circuit occur inside this product, an electric shock or fire may result.
- Do not connect/disconnect the Power Cord with wet hands: an electric shock might occur.
- If the Power Cord cannot be plugged into the outlet easily, do not use force. A foreign substance might be obstructing the outlet or a deformation of the connection portion may be present. Contact your dealer immediately.
- If you have any suspicion that the connection portion might be deformed, the connector might be cracked or otherwise damaged, or the wire in the Power Cord might be broken, do not plug the Power Cord into the wall outlet. Contact your dealer immediately.

#### 3-2 Connection of the Foot Control

By matching the mark ● on the Foot Control Plug and the mark ▼ on the Foot Control Connector at the front of the Control Unit, plug the Foot Control into the Connector all the way until you hear a clicking sound.

To disconnect the Foot Control, hold the Foot Control Plug and pull it out with a straight movement.

#### CAUTION

- When connecting/disconnecting the Foot Control, do not turn the Foot Control Plug. Otherwise, the inside pins might become damaged, leading to failure.

### 3-3 Connection of the Motor

By matching the mark ● on the Motor Cord Plug and the mark ● on the Motor Connector (A or B) at the front of the Control Unit, plug the Motor into the Connector all the way.

#### CAUTION

- There are two Motor Connectors. By connecting two Motors to the two connectors, you can select the one you want to use by the Motor Select Button (however, both Motors cannot be operated simultaneously). The selected Motor can be confirmed in the Main Panel.
- Do not pull the Motor Cord to disconnect the Motor.


Fig. 4

To disconnect the Motor, hold the Motor Cord Plug and pull it out with a straight movement.

#### CAUTION

- Whenever you use two Motor Connectors alternately, be sure to confirm which Motor is under operation.
- Do not connect anything else other than the Motor exclusively used for this system. Otherwise, the Motor Connector might become deformed or failure could occur.
- If the Motor Cord cannot be connected to the Control Unit easily, do not insist. Deformation or cracks might occur. A foreign substance might be obstructing the connection or the connection portion might be deformed. Contact your dealer.
- Whenever connecting/disconnecting the Motor with Hand Control and Wire Pin Driver, be sure to switch OFF the Safety Lock.
- When connecting/disconnecting the Motor, do not turn the Motor Cord Plug. Otherwise, the inside pins might become damaged, leading to failure.
- Prevent the Motor Cord Plug from falling onto the floor or receiving a strong impact; otherwise, deformation might occur and it might become impossible to establish a proper connection.
- Do not connect/disconnect the Motor Cord with wet hands; otherwise, electric shock might occur.
- Be sure to hold the Motor Cord Plug when disconnecting the Motor from the Control Unit to prevent failure caused by breaking wires.

### 3-4 Mounting of the Attachment

#### CAUTION

- Whenever mounting/removing the Attachment, be sure to switch OFF the Safety Lock of the Wire Pin Driver / Motor with Hand Control.
- Prior to mounting/removing the Attachment, remove the Cutting Accessories. Otherwise, injury might occur.
- After mounting the Attachment, confirm that the Attachment and Motor are connected properly.

#### ◆ P200-HMH and P200-HMH-HS Handpieces

##### <Mounting>

- 1) Push the Attachment into the Motor straight.
- 2) Turn the Attachment rightward or leftward so that the Alignment Pin of the Motor matches the hole of the Attachment with a clicking sound to connect firmly.

##### <Removing>

Pull out the Attachment from the Motor straight.


Fig. 5

#### ◆ P200-SMH-S, P200-SMH and P200-SMH-HS Handpieces

##### <Mounting>

- 1) Push the Attachment into the Motor straight.
- 2) Align the "▼" Symbol of the Attachment and "●" Symbol of the Motor, then push the Attachment in until it clicks to connect firmly.

##### <Removing>

Pull out the Attachment from the Motor straight.


Fig. 6

### ◆ P200-BMH and P200-BMH-HS Attachments

#### <Mounting>

- 1) Push the Attachment into the Motor straight.
- 2) Turn the Attachment rightward or leftward so that the Alignment Pin of the Motor matches the hole of the Attachment with a clicking sound to connect firmly.
- 3) Tighten the Removal Ring of the Attachment manually to fix it firmly.

#### <Removing>

Loosen the Removal Ring and pull out the Attachment from the Motor Straight.


Fig. 7

### ◆ P200-WPD Attachment

#### <Mounting>

Insert the Attachment straight into the Wire Pin Driver until it clicks.

- \* Wire Collet or Pin Collet should be mounted with the Attachment Lever facing down (Motor Cord side).

#### <Removing>

Pull the Attachment straight out while pressing the Mounting/Removal Button of the Wire Pin Driver.


Fig. 8

### 3-5 Mounting of the Irrigation Pole

Match the position of the convex part on the Irrigation Pole with the position of the groove in the Irrigation Pole Stand located at the back of the Control Unit, and insert the Irrigation Pole into the hole.


Fig. 9

### 3-6 Mounting of the Irrigation Tube

- 1) By rotating the Knob towards OPEN, open the Pump Cover.
- 2) Insert the Irrigation Tube (part with no blue line) into the Irrigation Pump in such a manner that the needle side of the Irrigation Tube is positioned on the right-hand side of the Pump, and securely place the stopper of the Irrigation Tube into the Guide Parts located at both ends.

Fig.10


Fig.11

- 3) Rotate the Knob towards CLOSE to close the Pump Cover.

Fig. 12

**CAUTION**

- After having confirmed that the Irrigation Tube is properly set on the roller, CLOSE the Knob of the Irrigation Pump. If the Knob is turned to CLOSE with the Irrigation Tube displaced from the proper position, the Irrigation Tube might become damaged (Fig. 13).
- Do not mount the Irrigation Tube in the opposite direction. The saline solution bottle might burst if the Control Unit keeps operating with the Irrigation Tube mounted opposite.
- Prior to inserting the Irrigation Tube into the Irrigation Pump, confirm that there is no foreign substance in the roller portion. When closing the Pump Cover, do not force the Knob excessively and be careful not to get your fingers caught.
- Only use the Irrigation Tube (Order Code: Y1003915) recommended by NSK. Otherwise, it may cause failure.


Fig.13

### 3-7 Mounting the Saline Solution

- 1) Close the Tube Clamp located between the Irrigation Tube Needle and the Irrigation Pump (Fig. 14).
- 2) Hang a saline solution on the Irrigation Pole, remove the cover of the Irrigation Tube Needle and insert the Needle into the bottle outlet. (Fig. 15.)
- 3) Mount an Irrigation Nozzle for the Attachment onto the other end of the Irrigation Tube (Fig. 16).
- 4) Open the Tube Clamp before operating the Control Unit.


Fig.14


Fig.15


Fig.16

**CAUTION**

- If the Tube Clamp is opened while the Pump Cover is open, the saline solution might flow out from the end of the Irrigation Tube. Therefore, prior to opening the Tube Clamp, confirm that the Knob of the Irrigation Pump is in the CLOSE position.
- If using an Irrigation Tube which is bent or kinked or if operating the Irrigation Pump in a state where the saline solution cannot flow out of the bottle properly, the Irrigation Tube might become damaged or the Irrigation Tube Needle might become disconnected from the bottle.

## 4 Check before Operation

Follow the check list below, before use. If any abnormalities are found, stop using the product immediately and contact your dealer.

### 4-1 Control Unit

- 1) Check for any abnormalities or damage on the outer surface.
- 2) Check that the Power Cord can be securely connected to the Control Unit and the medical-use outlet.
- 3) Turn the power ON and check for any abnormalities on the LCD display.
- 4) Touch the LCD display to check each setup can be customized.
- 5) Check that the Irrigation Pole can be securely inserted into the Pole Stand of the Control Unit.
- 6) Adjust the Irrigation Volume. Open the Pump Cover and check that the Roller operates normally.
- 7) Check that the Irrigation Tube can be securely mounted to the Irrigation Pump.

### 4-2 Foot Control

- 1) Check for any abnormalities or damage on the outer surface and for any damage or deformation on the Pins in the Foot Control Plug.
- 2) Connect the Foot Control to the Control Unit to check its secure connection.
- 3) Depress the Pedal and Buttons of the Foot Control to check that they operate normally.

### 4-3 Motor

- 1) Check for any abnormalities or damage on the outer surface and for any damage or deformation on the Pins in the Motor Cord Plug.
- 2) Check that the Motor can be securely connected to the Control Unit.
- 3) Depress the Pedal of the Foot Control to check for any abnormal vibration, noise or heating in the Motor. Check in the Main Panel of the Control Unit that the rotational speed can reach the maximum speed you have set.
- 4) Check that Forward, Reverse, and Oscillate (P200-WPD only) Modes operate normally.
- 5) Mount the Attachment to the Motor to check its secure connection.

\* If any abnormalities are found, stop using the product immediately and contact your dealer.

\* Operate the product if none of the above abnormalities are present.

\* When checking Attachments before treatment, refer to the Operation Manual of each Attachment.

## 5 Operation

### CAUTION

- The LCD Display of this product is a touch panel. Refrain from doing the following, which might result in failure or malfunction of the touch panel.
  - Operating the touch panel using a sharp object, such as a ballpoint pen.
  - Operating the touch panel with wet hands.
- Do not push the touch panel with force, otherwise failure might occur.
- Refer to the Operation Manuals of the Motor and Attachment to be mounted.

### 5-1 How to turn on the power

Turn ON the Power Switch. The HOME screen is displayed after the startup screen appears and the Control Unit finishes an automatic self-check of the system.

### CAUTION

- Do not turn OFF the Power Switch when the startup screen appears, because the Control Unit is conducting an automatic self-check of the system.
- If the Power Switch is turned ON while the Hand Switch/Trigger or Foot Control is depressed, the Motor does not start for safety reasons, and a message alert is displayed in the Main Panel. In order to deactivate the message, release the Hand Switch/Trigger or Foot Control.

## 5-2 Operation of the Main Panel

### 5-2-1 How to set the rotational speed and irrigation volume

Push the UP/DOWN Button or move the Track Bar to set the speed or volume.

(The Track Bar moves to the position which you have touched, within the adjustment area. Since it is rather difficult to accomplish a fine setting, use the UP/DOWN button for fine tuning.)


Fig.17

### 5-2-2 How to set the rotation mode

Touch the Rotation Mode Select Button to change the rotation mode sequentially.

Forward (clockwise)

Reverse (counterclockwise)

Oscillate (Clockwise and counterclockwise rotations are alternately repeated) Available only when P-200 WPD is mounted.

### 5-2-3 Switching between Motor A and B

Touch Motor Select Button A or B to switch to the desired Motor.


Fig.18

### 5-2-4 Priming Button

While holding the Motor, keep the Priming Button pushed, and confirm that the saline solution is being fed from the Irrigation Nozzle of the Attachment.


Fig.19

### 5-3 How to operate the Motor

This can be done in three ways: by the Motor with Hand Control, the Primado2 Wire Pin Driver or the Foot Control.

#### CAUTION

- The Motor with Hand Control and the Primado2 Wire Pin Driver can also be operated by the Foot Control. Connect it with the Foot Control and depress the Pedal. The Motor will work in the normal mode of the Foot Control, regardless of the Hand Control.
- If operating the Motor only with the Motor with Hand Control and the Primado2 Wire Pin Driver, it is recommended to disconnect the Foot Control from the Control Unit. Alternatively, set the Foot Control to OFF by the Foot Control ON/OFF Button (to disable the Foot Control).
- The Hand Switch of the Motor with Hand Control and the Pedal of the Foot Control are variable. The rotational speed increases/decreases within the preset rotation range, depending on the extent/duration of depressing.

#### ◆ When using the Motor with Hand Control (P200-HMH-HS, P200-SMH-HS, P200-BMH-HS)

- 1) Confirm that the Safety Lock of the Hand Control is OFF.
- 2) Slide the Adjustable Lever to the desired length.
- 3) Hold the Hand Switch and slide the Safety Lock to ON. (Fig. 20)
- 4) Depress the Hand Switch to start the Motor. When the Hand Switch is depressed or released, there is a short alarm sound.

When the Reverse Mode has been selected, the Bur rotates in reverse. During the reverse rotation in this Mode, the alarm continues to sound.


Fig.20

#### CAUTION

- When operating the Safety Lock, be sure to hold the lever with your finger. If operating it without holding the lever, the lever might inadvertently be depressed when the Safety Lock is being turned ON/OFF, thereby triggering a rotation of the Motor. This might cause an accident.

#### ◆ When using the Primado2 Wire Pin Driver (P200-WPD)

- 1) Confirm that the Safety Lock is working (  /  ).
- 2) Release the Safety Lock of the Wire Pin Driver from LOCK (  /  ) and switch it to Reverse OFF (FWD /  ) or Both ON (FWD/REV).
  - In the Reverse OFF (FWD /  ) Mode, when depressing the lower Trigger (B) of the Wire Pin Driver, the Attachment rotates forward. When depressing the Trigger or releasing it, a short alarm sounds. In this Mode, the upper Trigger (A) is locked, and therefore the Trigger cannot be depressed (i.e. it cannot be put into operation).
  - In the Both ON (FWD/REV) Mode, when depressing the lower Trigger (B), the Attachment rotates forward. When depressing the Trigger or releasing it, a short alarm sounds. When depressing the upper Trigger (A), the Attachment rotates in reverse. During the reverse rotation, the alarm continues to sound.
  - In the Both ON (FWD/REV) Mode, when depressing both the upper(A) and lower(B) Triggers simultaneously, the forward and reverse rotations of the Attachment are repeated alternately (i.e., the rotation direction oscillates). The rotation direction changes at every 360 degrees of rotation. If one of the Triggers has been depressed earlier, it will take priority in the function and rotation direction. When the other Trigger is depressed with the first Trigger already functioning, the forward and reverse rotations will start to alternate: namely, Oscillate Mode will start.


Fig.21

| Trigger                                        | Rotation Direction                                                   |
|------------------------------------------------|----------------------------------------------------------------------|
| Upper (A)                                      | Reverse Rotation                                                     |
| Lower (B)                                      | Forward Rotation                                                     |
| Upper and Lower depressed simultaneously (A B) | Oscillation (Forward and reverse rotations are repeated alternately) |

#### CAUTION

- The Wire Pin Driver operates in the rotation direction set not by the Control Unit but by the Wire Pin Driver.

### ◆ When Using the Foot Control (FC-73, FC-74)

- 1) When depressing the Pedal of the Foot Control, the Motor starts to operate. When depressing the Pedal or releasing it, a short alarm sounds.
- 2) In the Reverse Mode, the Bur starts rotating in reverse. During this reverse rotation, the alarm continues to sound.
- 3) Push the Foot Control Button (Multi Foot Control: blue, silver and green buttons; Single Foot Control: silver button) to set the Control Unit to the pre-determined function. For any further information, refer to Step “6-4 Setting the Foot Control”


Fig.22


Fig.23

#### CAUTION

- Before using the Foot Control Button, be sure to confirm what function is assigned to that Button in the Sub Panel.

#### How to confirm the Sub Panel

##### <Multiple Model>


Fig.24

##### <Single Model>

Fig.25

The function setting done by each Button is displayed. For more information on these functions, refer to Step “6-4 Setting the Foot Control”.

##### <Foot Control not connected>

Fig.26

##### <Foot Control OFF>


Fig.27

This shows that the Motor is switched OFF by the ON/OFF switching function operated from the Foot Control. Operation by the Foot Control is not available in this case. In order to operate it by the Foot Control, touch the Sub Panel  or turn the Foot Control ON by the Foot Control Button.

## 6 Option Set Up

It is possible to change the setting of the Control Unit and Motor by setting up the options.

The components which can be set up are as follows:

- Control Unit
- Motor
- Attachment
- Foot Control

- 1) Touch the Option Set Up Button in the Mini Panel in the HOME screen to display the Option Set Up screen.
- 2) Touch the Set Up Button for each component to display the Set Up screen.
- 3) Touch the HOME Button or BACK Button to return to the HOME screen.


Fig.28

Option Set Up Screen

Fig.29

### CAUTION

- The individual settings made for the Control Unit, Motor and Foot Control in the Option Set Up screen are stored in the internal memory. Therefore, even after turning OFF the power, the settings will remain. Exception: the Attachment settings are reset whenever the power is turned OFF or the Motor is disconnected from the Control Unit. Therefore, prior to use, be sure to confirm the settings of the Attachment.
- If the Motor or the Foot Control is not connected and the Option Set Up is attempted, a message alert will be displayed announcing that connection is not made.

### 6-1 Setting the Control Unit

Touch the Control Unit Set Up Button in the Option Set Up screen to display the Control Unit Set Up screen.

Control Unit Set Up Screen


Fig.30

### 6-1-1 Adjusting the background brightness and volume

The background brightness and volume level can be adjusted by either touching the UP/DOWN Button, or by moving the Track Bar to the desired level (Touch a position within the Adjustment Area to move the Track Bar to that position. For fine tuning, use the UP/DOWN Button).

Fig.31

### 6-1-2 Selecting the language

The language used for message alerts can be selected from amongst several choices.

- 1) Touch the Language Set Up Button in the Control Unit Set Up screen to display the Language Set Up screen.
- 2) Touch the button labeled with the desired language.
- 3) The Language Set Up Confirmation screen is displayed.
  - Touch the YES Button in the Language Set Up Confirmation screen to set the selected language as the main language; the Language Set Up screen will be displayed.
  - Touch the NO Button to cancel the selected language; the Language Set Up screen will be displayed.

Language Unit Set Up Screen

Fig.32

Language Set Up Confirmation Screen


Fig.33

### 6-1-3 Information display

Information on the version of the Control Unit is displayed.

## 6-2 Adjusting the Motor

Touch the Motor Set Up Button in the Option Set Up screen to display the Motor Set Up screen.

### 6-2-1 Setting methods for Acceleration, Brake, Torque and Oscillating Turns

Push the UP/DOWN Button to set the desired value.

When two Motors are connected to Motor Connector A and B, the Motor to be set up can be selected by the Motor Select Button in the Sub Panel.

- **Acceleration adjustment:**  
The acceleration time for the Motor to reach its maximum rotational speed can be set at any level between 10% and 100%.
- **Brake adjustment:**  
The time it takes for the Motor to stop can be set at any level between 10% and 100%.
- **Torque adjustment:**  
The torque level of the Motor can be set between 30% and 100%.
- **Oscillating turns adjustment:**  
The oscillating turns can be set between 1 and 5 (This is available only when the P200-WPD Wire Pin Driver is connected).


- While the Motor Set Up screen is displayed, the Motor can be started. Prior to using the Motor, be sure to return to the HOME screen.

### 6-2-2 How to return to the Factory Default settings

Touch the Factory Default Button and then the YES Button to return to the Factory Default settings provided before delivery.

Touch the NO Button to keep the current settings.

## 6-3 Setting Up the Attachment

By setting up the connected Attachment, the optimal rotational speed for the particular Attachment becomes available. After setting up the Attachment in the Attachment Set Up screen, return to the HOME screen to adjust the rotational speed and the Irrigation volume.

When two Motors are connected to Motor Connector A and B, the Attachment to be set up can be selected by the Motor Select Button in the Sub Panel.

Touch the Attachment Set Up Button in the Option Set Up screen to display the Attachment Set Up screen.

- The settings of the Attachment are reset whenever the power is turned off or the Motor is disconnected from the Control Unit. Prior to using the Attachment, be sure to confirm the current settings.
- If the Motor is not connected, a message alert announcing that the connection is not made will be displayed in the Main Panel.

### Attachment Set Up Screen

Since the Attachments available differ depending on the type of the Motor used, the Attachment Set Up screen matching the connected Motor is displayed.


Fig.36

### Compatible Attachments

| Motor    | Attachment   | Compatible Attachments                                 |
|----------|--------------|--------------------------------------------------------|
| P200-HMH | Direct Drive | PD-1A135, PD-1A155, PD-1A175                           |
|          |              | PD-2AES, PD-2AS, PD-2AM, PD-2AL                        |
|          |              | PD-2SES, PD-2SS, PD-2SM                                |
|          |              | PD-2SD, PD-2SDM                                        |
|          |              | PD-2AD, PD-2ADM                                        |
|          |              | PD-CRA                                                 |
|          | Perforator   | PD-PER                                                 |
|          | High Speed   | P200-3GAS, P200-3GAM, P200-3GAL                        |
|          | Super Slim   | PD-RA330, PD-RA330-L                                   |
|          | Contra 16:1  | PD-CA-16                                               |
|          | Contra 64:1  | PD-CA-64                                               |
| P200-SMH | Direct Drive | P200-1A135, P200-1A155, P200-1A155-C                   |
|          |              | P200-1A175, P200-1A175-C                               |
|          |              | P200-2AMS, P200-2AES, P200-2AS, P200-2AM, P200-2AL     |
|          |              | P200-2SMS, P200-2SES, P200-2SS, P200-2SM, P200-2SL     |
|          |              | P200-2SD                                               |
|          |              | P200-CRA                                               |
|          |              | P200-3MC                                               |
|          |              | P300-2S10, P300-2S20, P300-2S40, P300-2S60, P300-2S80  |
|          |              | P300-2A10, P300-2A20, P300-2A40, P300-2A60, P300-2A80  |
|          |              | P300-1AHS                                              |
|          |              | P300-1AHA                                              |
|          | Perforator   | P200-PER                                               |
|          | Super Slim   | P200-RA330, P200-RA330-L                               |
|          | Contra 4:1   | P200-CA-4                                              |
| P200-BMH | Bone Saw     | P200-REC, P200-SAG, P200-OSC                           |
| P200-WPD | Drill        | PD-DA-J4, PD-DA-J6, PD-DA-A, PD-DA-K, PD-DA-W, PD-DA-P |
|          | Bone Saw*    | PD-SA-S, PD-SA-R                                       |
|          | Reamer       | PD-RA-M, PD-RA-H, PD-RA-T, PD-RA-A                     |

\* P200-WPD is compatible with the Oscillating Mode, but the Mode does not start when this type of Attachment is mounted.

## 6-4 Setting Up the Foot Control

1) Touch the Foot Control Set Up Button in the Option Set Up screen to display the Foot Control Set Up screen.


Fig.37

2) Touch the Foot Control Button (Multi Foot Control: blue, silver and green buttons; Single Foot Control: silver button) which you want to set up to display the Foot Control Set Up screen. Select the desired functional Button; then, that Button will be assigned to the selected function, and the Foot Control Set Up screen will be displayed.


Fig.38

- **ON/OFF Select:**  
The ON/OFF selection of operations by the Foot Control can be set up. Touch the OFF Button to cancel operation by the Foot Control; the ON/OFF selection itself remains.
- **Motor Select:**  
Switching between the two Motors connected to Motor Connector A and B can be done.
- **Speed Select:**  
The rotational speed of the Motor can be switched from 30% to the default value, then to 100%.

**CAUTION** • Beware that previous rotational speed is not stored.

- **Rotation Mode Select:**  
The rotation mode of the Motor can be switched.
- **Irrigation ON/OFF Select:**  
Irrigation can be switched ON/OFF.
- **Flush:**  
Push the Foot Control Button continuously to keep the irrigation at 100% capacity.
- **Disable:**  
Disables function to the Foot Control Button.
- **“Foot Control Button” Display (only for the Multi Foot Control):**  
The Button which is being set up is highlighted.

## 7 Fuse Replacement

- 1) Remove the Fuse Box located in the Inlet Box at the back of the Control Unit, while keeping the nails at both ends pushed inward.
- 2) Replace the fuses (2 pcs) with a new one and push the Fuse Box into the Inlet Box until you hear a clicking sound.

### CAUTION

- Do not use any type of fuse other than ones NSK recommends. Otherwise, failure might occur.

| Electrical fuses |              | Order Code  |
|------------------|--------------|-------------|
| 100V/120V        | T3.15AL/250V | D1200131510 |
| 230V             | T1.6AL/250V  | D1200116010 |

- If the reason for the fuse burning is unknown, contact your dealer.


Fig.39

## 8 Periodical Maintenance Checks

Every three months, perform periodical maintenance checks on this product according to the following list. If any abnormalities are found, contact your Authorized NSK Dealer.

| Points to Check      | Details                                                                                                                                                                                                                                                                                                                                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Outer Surface        | Check for any scratches, damage or abnormalities on the outer surface.                                                                                                                                                                                                                                                                                                                                                         |
| Each Connection Part | Check for any scratches, damage or deformation of each connection part.                                                                                                                                                                                                                                                                                                                                                        |
| Control Unit         | Check for any damage, breakage or deformation on each cord plug.                                                                                                                                                                                                                                                                                                                                                               |
|                      | Make sure that the Motor Cord Plug, Foot Control Plug and Power Cord Plug can be securely connected to the Control Unit.                                                                                                                                                                                                                                                                                                       |
|                      | Connect the Power Cord to the medical-use outlet and turn on the power switch.                                                                                                                                                                                                                                                                                                                                                 |
|                      | Make sure that the HOME screen appears on the LCD Touch Panel.                                                                                                                                                                                                                                                                                                                                                                 |
|                      | Touch the LCD Touch Panel to check each setup can be customized.                                                                                                                                                                                                                                                                                                                                                               |
| Foot Control         | Depress the Pedal and Buttons of the Foot Control to check that they operate normally.                                                                                                                                                                                                                                                                                                                                         |
| Motor                | Depress the Pedal of the Foot Control or Hand Switch of the Motor to check for any abnormal vibration, noise or heating in the Motor.                                                                                                                                                                                                                                                                                          |
|                      | Check that Forward, Reverse, and Oscillate (P200-WPD only) Modes operate normally.                                                                                                                                                                                                                                                                                                                                             |
|                      | Mount the Attachment to the Motor to check its secure connection.                                                                                                                                                                                                                                                                                                                                                              |
|                      | After connecting each part, rotate the Motor for approximately 1 minute by pressing the Foot Control Pedal or Hand Switch of the Motor. Check for any abnormal vibration, noise or heat generation. Check on the Main Panel of the LCD Touch Panel that the rotational speed can reach the maximum speed you have set. Check that the Motor stops normally after releasing the Foot Control Pedal or Hand Switch of the Motor. |

### Additional safety information:

Periodical safety testing shall be performed every two years following IEC 62353 (EN 62353), in accordance with test procedure OM-SE0904EN. As periodical safety testing is limited exclusively to professional who have been trained and qualified under the control of NSK, contact your authorized NSK dealer.

## 9 Reprocessing

Perform reprocessing appropriately according to OM-SZ0911MA Primado2 REPROCESSING MANUAL included in the package of this product.

## 10 Error Codes

When the Motor has stopped operating due to failure, overloading, breaking of a wire, misuse etc., the Control Unit displays Error Codes. Then, an error code is displayed in the Mini Panel. If an error code is displayed, follow the check and measure steps shown in the display to resolve the error. Touch the Close Button in the Sub Panel or operate the Motor or the Foot Control after resolving the error to close the error display. In case the error is unresolved, contact your dealer for repair.

Fig.40

### NOTICE

- If two or more error codes are involved, the error code displayed in the Mini Panel will blink. Touch the Mini Panel to display the other error codes in sequence.

| Error codes | Error contents                | Causes                                                                         | Checks and measures                                                                                                                                                                                                                                                                          |
|-------------|-------------------------------|--------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| E00         | Error of the Control Unit.    | Abnormality or failure of the internal memory of the Control Unit.             | Turn ON the power source again. If the Control Unit operates normally, there is no problem.                                                                                                                                                                                                  |
| E01         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                        | Fasten the Twist Collet of the Attachment to lock the Bur. If the Motor Handpiece works normally, there is no problem.<br>Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem.                                                    |
|             |                               | Use under high load for a long period (excessive electric current - hardware). | Re-start. If it works normally, the error was caused by temporary overloading.                                                                                                                                                                                                               |
|             |                               | Failure of the Motor Handpiece.                                                | Replace with a spare Motor Handpiece to check whether or not it works normally. If the spare Motor Handpiece works normally, this means that the original one is damaged. Contact your dealer for repair.                                                                                    |
|             |                               | Failure of the internal circuits in the Control Unit.                          | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                                                                               |
| E02         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                        | Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem.                                                                                                                                                                              |
|             |                               | Failure of the Motor Handpiece (excessive voltage).                            | Replace with a spare Motor Handpiece to check whether or not it works normally. If the spare Motor Handpiece works normally, this means that the original one is damaged. Contact your dealer for repair.                                                                                    |
|             |                               | Failure of the internal circuits in the Control Unit.                          | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                                                                               |
| E03         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                        | Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem.                                                                                                                                                                              |
|             |                               | Failure of the Motor Handpiece (sensor signals).                               | Replace with a spare Motor Handpiece to check whether or not it works normally. If the spare Motor Handpiece works normally, this means that the original one is damaged. Contact your dealer for repair.                                                                                    |
|             |                               | Failure of the internal circuits in the Control Unit.                          | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                                                                               |
| E04         | Error of the Control Unit.    | Temperature rise in the Control Unit under high load for a long period.        | After letting the Control Unit rest in a cooler place for approximately 10 minutes, confirm the operation. If the Control Unit works normally, there is no problem.<br>Check whether or not the environment where the Control Unit is used or the installation place has a high temperature. |
|             |                               | Failure of the internal circuits in the Control Unit.                          | Turn ON the power source again. If the Control Unit operates normally, there is no problem.                                                                                                                                                                                                  |

| Error codes | Error contents                | Causes                                                                                           | Checks and measures                                                                                                                                                                                                                       |
|-------------|-------------------------------|--------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| E05         | Error of the Control Unit.    | Occurrence of abnormal voltage in or failure of the startup / stop circuits of the Control Unit. | Turn ON the power source again and start-stop several times.<br>If the Control Unit works normally, there is no problem.                                                                                                                  |
| E06         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                                          | Fasten the Twist Collet of the Attachment to lock the Bur. If the Motor Handpiece works normally, there is no problem.<br>Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem. |
|             |                               | Use under high load for a long period (rotor lock).                                              | If the re-starting proceeds normally, the error was caused by temporary overloading. There is no problem.                                                                                                                                 |
|             |                               | Failure of the Motor Handpiece.                                                                  | Replace with a spare Motor Handpiece to check whether or not it works normally. If the spare Motor Handpiece works normally, this means that the original one is damaged. Contact your dealer for repair.                                 |
|             |                               | Failure of the internal circuits of the Control Unit.                                            | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                            |
|             |                               |                                                                                                  |                                                                                                                                                                                                                                           |
| E07         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                                          | Fasten the Twist Collet of the Attachment to lock the Bur. If the Motor Handpiece works normally, there is no problem.<br>Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem. |
|             |                               | Use under high load for a long period (excessive electric current - hardware).                   | Re-start. If it works normally, the error was caused by temporary overloading. There is no problem.                                                                                                                                       |
|             |                               | Failure of the Motor Handpiece.                                                                  | Replace with a spare Motor Handpiece to check whether or not it works normally. If the spare Motor Handpiece works normally, this means that the original one is damaged. Contact your dealer for repair.                                 |
|             |                               | Failure of the internal circuits of the Control Unit.                                            | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                            |
|             |                               |                                                                                                  |                                                                                                                                                                                                                                           |
| E08         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                                          | Fasten the Twist Collet of the Attachment to lock the Bur. If the Motor Handpiece works normally, there is no problem.<br>Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem. |
|             |                               | Use under high load for a long period (excessive electric current – software).                   | Re-start. If it works normally, the error was caused by temporary overloading. There is no problem.                                                                                                                                       |
|             |                               | Failure of the Motor Handpiece.                                                                  | Replace with a spare Motor Handpiece to check whether or not it works normally. If the spare Motor Handpiece works normally, this means that the original one is damaged. Contact your dealer for repair.                                 |
|             |                               | Failure of the internal circuits of the Control Unit.                                            | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                            |
|             |                               |                                                                                                  |                                                                                                                                                                                                                                           |
| E09         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                                                          | Re-connect the Motor Handpiece to the Control Unit to check whether or not the Hand Switch/Trigger works normally. If it works normally, there is no problem.                                                                             |
|             |                               | Failure of the Motor Handpiece (Hand Switch/Trigger).                                            | Replace with a spare Motor Handpiece to check whether or not the Hand Switch/Trigger works normally. If the Hand Switch/Trigger works normally, this means that the original Motor Handpiece is damaged. Contact your dealer for repair.  |
|             |                               | Failure of the internal circuits of the Control Unit.                                            | Turn ON the power source again. If the Motor Handpiece operates normally, there is no problem.                                                                                                                                            |
| E10         | Error of the Foot Control.    | Poor connection of the Foot Control.                                                             | Re-connect the Foot Control to the Control Unit to check whether or not the Foot Control works normally. If it works normally, there is no problem.                                                                                       |
|             |                               | Failure of the Foot Control.                                                                     | If you have a spare Foot Control, replace the Foot Control with the spare one to check whether or not the spare one works. If it works, the original Foot Control is damaged. Contact your dealer for repair.                             |
|             |                               | Failure of the internal circuits of the Control Unit.                                            | Turn ON the power source again. If the Foot Control operates normally, there is no problem.                                                                                                                                               |

| Error codes | Error contents                | Causes                                                               | Checks and measures                                                                                                                                                                                           |
|-------------|-------------------------------|----------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| E11         | Error of the Irrigation Pump. | Foreign substance caught in the Irrigation Pump.                     | Check whether or not any foreign substance is caught in the Irrigation Pump. If any, remove it. Re-mount the Irrigation Tube. If the Irrigation Pump works normally, there is no problem.                     |
|             |                               | Failure of the Irrigation Pump (excessive electric current).         | Turn ON the power source again. If the same error code is displayed, the Irrigation Pump might be damaged. Contact your dealer for repair of the Control Unit.                                                |
|             |                               | Failure of the internal circuits of the Control Unit.                | Turn ON the power source again. If the Irrigation Pump works normally, there is no problem.                                                                                                                   |
| E12         | Error of the Irrigation Pump. | Failure of the Irrigation Pump (excessive voltage).                  | Turn ON the power source again. If the same error code is displayed, the Irrigation Pump might be damaged. Contact your dealer for repair of the Control Unit.                                                |
|             |                               | Failure of the internal circuits of the Control Unit.                | Turn ON the power source again. If the Irrigation Pump works normally, there is no problem.                                                                                                                   |
| E13         | Error of the Control Unit.    | Signal communication abnormality or failure of the Display.          | Turn ON the power source again. If the Control Unit works normally, there is no problem.                                                                                                                      |
|             |                               | Failure of the internal circuits of the Control Unit.                | Turn ON the power source again. If the Control Unit works normally, there is no problem.                                                                                                                      |
| E14         | Error of the Foot Control.    | Poor connection of the Foot Control.                                 | Re-connect the Foot Control to the Control Unit. If the Foot Control works normally, there is no problem.                                                                                                     |
|             |                               | Failure of the Foot Control.                                         | If you have a spare Foot Control, replace the Foot Control with the spare one to check whether or not the spare one works. If it works, the original Foot Control is damaged. Contact your dealer for repair. |
|             |                               | Failure of the internal circuits of the Control Unit.                | Turn ON the power source again. If the Foot Control works normally, there is no problem.                                                                                                                      |
| E15         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                              | Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece is recognized and displayed normally, there is no problem.                                                                         |
|             |                               | Failure of the Motor Handpiece (recognition of the Motor Handpiece). | Replace the Motor Handpiece with a spare one to check whether or not the spare one works normally. If it works normally, the original Motor Handpiece is damaged. Contact your dealer for repair.             |
|             |                               | Failure of the internal circuits in the Control Unit.                | Turn ON the power source to check whether or not the Motor Handpiece works normally. If it does, there is no problem.                                                                                         |
| E16         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                              | Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem.                                                                                               |
|             |                               | Failure of the Motor Handpiece (speed of the Motor Handpiece).       | Replace the Motor Handpiece with a spare one to check whether or not the spare one works normally. If it works normally, the original Motor Handpiece is damaged. Contact your dealer for repair.             |
|             |                               | Failure of the internal circuits in the Control Unit.                | Turn ON the power source to check whether or not the Motor Handpiece works normally. If it does, there is no problem.                                                                                         |
| E17         | Error of the Control Unit.    | Abnormality or failure of the internal memory of the Control Unit.   | Turn ON the power source to check whether or not the Control Unit works normally. If it does, there is no problem.                                                                                            |
| E18         | Error of the Motor Handpiece. | Poor connection of the Motor Handpiece.                              | Re-connect the Motor Handpiece to the Control Unit. If the Motor Handpiece works normally, there is no problem.                                                                                               |
|             |                               | Failure of the Motor Handpiece (Sensor signals).                     | Replace the Motor Handpiece with a spare one to check whether or not the spare one works normally. If it works normally, the original Motor Handpiece is damaged. Contact your dealer for repair.             |
|             |                               | Failure of the internal circuits in the Control Unit.                | Turn ON the power source to check whether or not the Control Unit works normally. If it does, there is no problem.                                                                                            |
| E19         | Error of the Control Unit.    | Signal communication abnormality or failure of the Display.          | Turn ON the power source to check whether or not the Control Unit works normally. If it does, there is no problem.                                                                                            |
|             |                               | Failure of the internal circuits in the Control Unit.                |                                                                                                                                                                                                               |

## 11 Troubleshooting

If you are suspecting any faults, check the following once again prior to asking for repair work. If none of these cases is applicable or if the situation does not improve despite your efforts or measures, the damage could be serious; call your dealer to fix the problem.

### <Control Unit>

| Symptoms                                                                   | Causes/Points of checking                                                     | Countermeasures                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------|-------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The HOME screen does not appear.                                           | Is the power plug disconnected?                                               | Connect the power plug to the wall outlet properly (Refer to 3-1).<br>* If the screen still does not appear, contact your dealer.                                                                                                                                                                                            |
|                                                                            | Is the fuse blown out?                                                        | Replace it with the fuse specified by us (Refer to 7).<br>* Be sure to replace it with a fuse specified by us.                                                                                                                                                                                                               |
|                                                                            | Failure of the Power Switch                                                   | Contact your dealer.                                                                                                                                                                                                                                                                                                         |
| An error code appears.                                                     | Refer to 9. Error Codes.                                                      | Refer to the column of Checks and measures.                                                                                                                                                                                                                                                                                  |
| The Motor does not rotate.                                                 | Is the Foot Control Cord disconnected?                                        | Plug the Cord into the Connector of the Control Unit properly (Refer to 3-2, 3-3).                                                                                                                                                                                                                                           |
|                                                                            | Is the Motor Cord disconnected?                                               | * If it still does not function, contact your dealer.                                                                                                                                                                                                                                                                        |
|                                                                            | Is the Safety Lock of the Motor with Hand Control/Wire Pin Driver set to OFF? | If the Safety Lock of the Motor with Hand Control/Wire Pin Driver is set to OFF, release the Safety Lock (Refer to 5-3).                                                                                                                                                                                                     |
|                                                                            | Is the Foot Control set to OFF?                                               | If the Foot Control is set to OFF, cancel that setting (Refer to 5-3).<br>* If the Motor still does not function, contact your dealer.                                                                                                                                                                                       |
|                                                                            | Is the load too high?                                                         | If excessively loaded by the Attachment or the Cutting Accessory, the Control Unit activates the safety mechanism to stop the Motor from rotating. Remove the root cause, and then depress the Pedal of the Foot Control or re-push the Hand Switch/Trigger.<br>* If the Motor still does not function, contact your dealer. |
| Cannot increase the rotational speed.                                      | Is the rotational speed restricted?                                           | Adjust the Speed Adjustment Area in the Main Panel to set the desired rotational speed (Refer to 5-2).                                                                                                                                                                                                                       |
| Cannot reverse the rotation.                                               | Is the setting for Forward?                                                   | Select the Rotation Mode Select Button in the Main Panel or the Foot Control Button to switch (Refer to 5-2, 6-4).<br>* If the rotation direction does not change, contact your dealer.                                                                                                                                      |
| Cannot irrigate.                                                           | Are there kinks, signs of plugging or cracks in the Irrigation Tube?          | Fix the kinks in the Irrigation Tube first and then use it. If there are signs of plugging or cracks, replace the Tube with a new one (Refer to 3-6, 3-7).                                                                                                                                                                   |
|                                                                            | Is the irrigation setting at 0% or OFF?                                       | Set the desired irrigation level at the Irrigation Volume Adjustment Area in the Main Panel of the Control Unit (Refer to 5-2).<br>If the irrigation is set to OFF, push the ON Button displayed in the Main Panel of the Control Unit or push the Foot Control Button to cancel that setting.                               |
|                                                                            | Is the Irrigation Tube mounted in the wrong direction?                        | Re-check the mounting orientation of the Irrigation Tube installed in the Irrigation Pump. If the Tube is mounted in the wrong direction, remove it and re-install (Refer to 3-6).                                                                                                                                           |
| Weak torque; slow responsiveness in rotation; takes time to stop rotating. | Is the Motor set at a low value?                                              | Re-check the setting of the Motor in the Option Set Up screen. If the set value is low, adjust it.<br>Alternatively, return the setting of the Motor to the factory default (Refer to 6-2).                                                                                                                                  |
| Dark screen.                                                               | Is the Control Unit set to a low value?                                       | Adjust the brightness of the background by setting up the Control Unit in the Option Set Up screen (Refer to 6-1-1).                                                                                                                                                                                                         |
| No sound.                                                                  | Is the Volume set at 0?                                                       | Adjust the sound volume by setting the Control Unit in the Option Set Up screen (Refer to 6-1-1).                                                                                                                                                                                                                            |

<Motor (without Hand Control)>

| Symptoms                                                    | Causes/Points of checking                                                                                                                                        | Countermeasures                                                                                                                                                                       |
|-------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No rotation after a Cutting Accessories has been attached.  | The Twist Collet on the Attachment is OPEN.                                                                                                                      | After inserting the Cutting Accessory, fasten the Twist Collet (Turn it in the direction of '  '). |
|                                                             | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> | Contact your dealer.                                                                                                                                                                  |
| Abnormal heating during rotation.                           | Continuous use over a long period.                                                                                                                               | Stop operating the Attachment and Motor immediately. Wait until the Attachment and Motor cool down before continuing to operate.                                                      |
|                                                             | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> | Contact your dealer.                                                                                                                                                                  |
| Abnormal vibration and noise.<br>Cutting Accessory run-out. | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> |                                                                                                                                                                                       |

<Motor (with Hand Control)>

| Symptoms                                                                     | Causes/Points of checking                                                                                                                                        | Countermeasures                                                                                                                                                                       |
|------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No rotation when the Hand Switch is depressed.                               | Twist Collet of the Attachment is OPEN.                                                                                                                          | After inserting the Cutting Accessory, fasten the Twist Collet (Turn it in the direction of '  '). |
|                                                                              | Safety Lock is "ON".                                                                                                                                             | Release the safety Lock.                                                                                                                                                              |
|                                                                              | Safety Lock is "OFF".                                                                                                                                            | Contact your dealer.                                                                                                                                                                  |
| Rotation when the Hand Switch is not depressed.                              | The Hand Switch cannot get back to the specified position due to foreign substances.                                                                             |                                                                                                                                                                                       |
| No rotation even after the Twist Collet of the Attachment has been fastened. | The bearing has seized; Foreign substances are obstructing the bearing.                                                                                          |                                                                                                                                                                                       |
| Abnormal heating during rotation.                                            | Continuous use over a long period.                                                                                                                               | Stop operating the Attachment and Motor immediately. Wait until the Attachment and Motor cool down before continuing to operate.                                                      |
|                                                                              | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> | Contact your dealer.                                                                                                                                                                  |
| Abnormal vibration and noise.<br>Cutting Accessory run-out.                  | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> |                                                                                                                                                                                       |

<Wire Pin Driver>

| Symptoms                                                    | Causes/Points of checking                                                                                                                                        | Countermeasures                                                                                                                  |
|-------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| No rotation when the Trigger is depressed.                  | The Safety Lock is "ON".                                                                                                                                         | Release the Safety Lock.                                                                                                         |
|                                                             | The Safety Lock is "OFF".                                                                                                                                        | Contact your dealer.                                                                                                             |
| Rotation when the Trigger is not depressed.                 | The Trigger cannot get back to the specified position due to foreign substances.                                                                                 |                                                                                                                                  |
| The Attachment cannot be mounted.                           | The insertion is deformed.                                                                                                                                       |                                                                                                                                  |
| Abnormal heating during rotation.                           | Continuous use over a long period.                                                                                                                               | Stop operating the Attachment and Motor immediately. Wait until the Attachment and Motor cool down before continuing to operate. |
|                                                             | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> | Contact your dealer.                                                                                                             |
| Abnormal vibration and noise.<br>Cutting Accessory run-out. | <ul style="list-style-type: none"> <li>• The bearing has seized.</li> <li>• Foreign substance inside the bearing.</li> <li>• The bearing is worn out.</li> </ul> |                                                                                                                                  |

## 12 Specifications

### <Control Unit>

|                                      |                                      |
|--------------------------------------|--------------------------------------|
| Model                                | P200-CU-120<br>P200-CU-230           |
| Power source                         | AC 120V 50/60 Hz<br>AC 230V 50/60 Hz |
| Electricity consumption              | 65 VA                                |
| Maximum feed rate of Irrigation Pump | 75 ml/min                            |
| Dimensions                           | W278×D268×H160 mm                    |
| Weight                               | 6.9Kg                                |
| Length of Power Cord                 | 3.6 m                                |

### <Motor>

| Model                | P200-HMH                | P200-HMH-HS | P200-SMH-S              | P200-SMH     | P200-SMH-HS    |
|----------------------|-------------------------|-------------|-------------------------|--------------|----------------|
| Rotational speed     | 60,000min <sup>-1</sup> |             | 80,000min <sup>-1</sup> |              |                |
| Rated input          | DC 39.5V                |             | DC 23.7V                |              |                |
| Dimensions           | Ø23×L98.9mm             | Ø23×L98.9mm | Ø18×L84.2mm             | Ø18×L102.7mm | Ø19.4×L120.5mm |
| Weight               | 361g                    | 390g        | 289g                    | 305g         | 340g           |
| Operation time       | 3 min. ON / 10 min. OFF |             |                         |              |                |
| Length of Motor Cord | 3.5m                    |             |                         |              |                |

| Model                | P200-BMH                | P200-BMH-HS     | P200-WPD                |
|----------------------|-------------------------|-----------------|-------------------------|
| Rotational speed     | 48,000min <sup>-1</sup> |                 | 1,200min <sup>-1</sup>  |
| Rated input          | DC 21.8V                |                 | DC 39.1V                |
| Dimensions           | Ø18 x L100.8mm          | Ø20.4x L102.3mm | W74.9 x D28.0 x H143.1  |
| Weight               | 302g                    | 359g            | 790g                    |
| Operation time       | 1 min. ON / 15 min. OFF |                 | 1 min. ON / 10 min. OFF |
| Length of Motor Cord | 3.5m                    |                 |                         |

### <Foot Control>

| Model                       | FC-73              | FC-74              |
|-----------------------------|--------------------|--------------------|
| Dimensions                  | W153×D223×H175.5mm | W220×D223×H175.5mm |
| Weight                      | 1.4Kg              | 2.0Kg              |
| Length of Foot Control Cord | 3.5 m              |                    |

|                                      | Temperature              | Humidity  | Atmospheric pressure |
|--------------------------------------|--------------------------|-----------|----------------------|
| Use Environment                      | 0-40°C (no condensation) | 30-75% RH | 700-1,060hPa         |
| Store and Transportation Environment | -10-50°C                 | 10-85% RH | 500-1,060hPa         |

### Classifications of equipment

- Type of protection against electric shock:
  - Class I equipment
- Degree of protection against electric shock:
  - Type BF applied part  : Motor, Attachment
- Method of sterilization or disinfection recommended by the manufacturer:
  - Refer to OM-SZ0911MA Primado2 REPROCESSING MANUAL.
- Degree of protection against ingress of water as detailed in the current edition of IEC 60529:
  - Foot Control: IPX8 (Protected against immersion in water)
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:
  - Foot Control: AP Equipment 
- Mode of operation:
  - Intermittent operation

### Operation Principal

This system is composed of the Control Unit, the Foot Control, the Motor and the Attachment.

Power is supplied from the Control Unit to the Motor by operating the Foot Control or the Hand Switch of the Motor. Due to the motor rotary, the Cutting Accessory attached to the Attachment rotates.

## 13 Disposing of Product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose of the product for you.

## 14 Warranty

Our products are subjected to rigorous quality control and inspections. In the unlikely event of a product proving to be faulty within the guarantee period when it has been used normally, it will be repaired free of charge subject to terms and conditions below. For such repair requests, please contact your dealer.

Guarantee period ··· One year

- In the following cases, the repair has to be paid for even if it is within the guarantee period.
  - Fault and damage caused by misuse, or inappropriate repair or modification.
  - Fault and damage caused during transit or dropping the product after it was purchased.
  - Fault and damage attributable to fire, earthquake, flooding, power surge, environmental pollution or other act of providence.
  - Fault and damage caused by the use of consumables, or methods other than those specified by us.
  - When the product has been repaired using components other than genuine components supplied by us.
- For repairs after the guarantee period, too, please contact your dealer.

## 15 Symbols

|                                                                                     |                                                                                                                                                         |                                                                                     |                                                                                                 |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
|    | Manufacturer                                                                                                                                            |    | Medical Device                                                                                  |
|  | Conforms to CE European Directive of "Medical device directive 93/42/EEC"                                                                               |  | EOG sterilization                                                                               |
|  | This product is designed not to become an ignition source in flammable air and anesthetic gas environment                                               |  | Single use only. Do not reuse.                                                                  |
|  | Protected against the effects of continuous immersion in dust and water                                                                                 |  | Use by                                                                                          |
|  | Type BF applied part                                                                                                                                    |  | Batch number                                                                                    |
|  | Authorized representative in the European community                                                                                                     |  | Catalog number                                                                                  |
|  | Date of manufacture                                                                                                                                     |  | Serial No.                                                                                      |
|  | See operation manual                                                                                                                                    |  | "ON" state of Main Power Switch                                                                 |
|  | Caution, Consult operation instructions.                                                                                                                |  | "OFF" state of Main Power Switch                                                                |
|  | Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment. |  | Keep dry                                                                                        |
|  | Dispose of this device and its accessories via methods approved for electronic device and in compliance with the Directive 2012/19/EU                   |  | This is the correct upright position of the distribution packages for transport and/or storage. |
|  | This product can be sterilized in a steam sterilizer at 135°C.                                                                                          |  | Fragile / handle with care                                                                      |
|  | This product can be cleaned using a washer-disinfector                                                                                                  |  | Temperature limitation                                                                          |
|                                                                                     |                                                                                                                                                         |  | Humidity limitation                                                                             |
|                                                                                     |                                                                                                                                                         |  | Atmospheric pressure limitation                                                                 |

## 16 EMC Information

| Guidance and manufacturer's declaration - electromagnetic emissions                                                                                                                 |                              |                                                                                                                                                                                   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The Primado2 is intended for use in the electromagnetic environment specified below. The customer or the user of the Primado2 should assure that it is used in such an environment. |                              |                                                                                                                                                                                   |
| Emissions test                                                                                                                                                                      | Compliance                   | Electromagnetic environment - guidance                                                                                                                                            |
| RF emissions CISPR11/ EN55011                                                                                                                                                       | Group 1<br>Class B           | The Primado2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR11/ EN55011                                                                                                                                                       | Group 1<br>Class B           |                                                                                                                                                                                   |
| Harmonic emissions<br>EN/IEC61000-3-2                                                                                                                                               | Class A<br>(exception 120V)  |                                                                                                                                                                                   |
| Voltage fluctuations/flicker emissions<br>EN/IEC61000-3-3                                                                                                                           | Complies<br>(exception 120V) |                                                                                                                                                                                   |

| Guidance and manufacturer's declaration - electromagnetic immunity                                                                                                                  |                                                                                                                                                        |                                                                                                                                                        |                                                                                                                                                                                                                                                                                  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The Primado2 is intended for use in the electromagnetic environment specified below. The customer or the user of the Primado2 should assure that it is used in such an environment. |                                                                                                                                                        |                                                                                                                                                        |                                                                                                                                                                                                                                                                                  |
| Immunity test                                                                                                                                                                       | EN/IEC60601 test level                                                                                                                                 | Compliance level                                                                                                                                       | Electromagnetic environment - guidance                                                                                                                                                                                                                                           |
| Electrostatic discharge (ESD)<br>EN/IEC61000-4-2                                                                                                                                    | ±6kV contact<br>±8kV air                                                                                                                               | ±6kV contact<br>±8kV air                                                                                                                               | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                    |
| Electrical fast transient/burst<br>EN/IEC61000-4-4                                                                                                                                  | ± 2kV for power supply lines<br>± 1kV for input/output lines                                                                                           | ± 2kV for power supply lines<br>± 1kV for input/output lines                                                                                           | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                              |
| Surge<br>EN/IEC61000-4-5                                                                                                                                                            | ±1kV line(s) to line(s)<br>±2kV line(s) to earth                                                                                                       | ±1kV line(s) to line(s)<br>±2kV line(s) to earth                                                                                                       | Mains power quality should be that of a typical commercial or hospital environment                                                                                                                                                                                               |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>EN/IEC61000-4-11                                                                            | <5% Ut (>95% dip in Ut)for 0.5 cycle<br>40% Ut (60% dip in Ut)for 5 cycles<br>70% Ut (30% dip in Ut)for 25 cycles<br><5% Ut (>95% dip in Ut)for 5 secs | <5% Ut (>95% dip in Ut)for 0.5 cycle<br>40% Ut (60% dip in Ut)for 5 cycles<br>70% Ut (30% dip in Ut)for 25 cycles<br><5% Ut (>95% dip in Ut)for 5 secs | Mains power quality should be that of a typical commercial or hospital environment. If the user of the Primado2 requires continued operation during power mains interruptions, it is recommended that the Primado2 be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>EN/IEC61000-4-8                                                                                                                         | 3A/m                                                                                                                                                   | 3A/m                                                                                                                                                   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                        |
| NOTE 'Ut' is the AC mains voltage prior to application of the test level                                                                                                            |                                                                                                                                                        |                                                                                                                                                        |                                                                                                                                                                                                                                                                                  |

| Guidance and manufacturer's declaration - electromagnetic immunity.                                                                                                                 |                           |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The Primado2 is intended for use in the electromagnetic environment specified below. The customer or the user of the Primado2 should assure that it is used in such an environment. |                           |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Immunity test                                                                                                                                                                       | EN/IEC60601 test level    | Compliance level | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Conducted RF<br>EN/IEC61000-4-6                                                                                                                                                     | 3Vrms<br>150 kHz to 80MHz | 3Vrms            | <p>Portable and mobile RF communications equipment should be used no closer to any part of the Primado2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance<br/> <math>d = 1.2\sqrt{P}</math><br/> <math>d = 1.2\sqrt{P}</math> 80MHz to 800MHz<br/> <math>d = 2.3\sqrt{P}</math> 800MHz to 2.5GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters as determined by an electromagnetic site survey<sup>(a)</sup>, should be less than the compliance level in each frequency range<sup>(b)</sup>. Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| Radiated RF<br>EN/IEC61000-4-3                                                                                                                                                      | 3V/m<br>80MHz to 2.5 GHz  | 3V/m             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Primado2 is used exceeds the applicable RF compliance level above, the Primado2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Primado2.

b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.

| Cables and accessories | Maximum length    | Complies with                                                                                                                                                                                                                        |
|------------------------|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Motor Cord             | 3.5m (Unshielded) | RF emissions, CISPR11: Class B/Group 1<br>Harmonic emissions (exception 120V): EN/IEC61000-3-2<br>Voltage fluctuations/ flicker emission (exception 120V): EN/IEC61000-3-3<br>Electrostatic discharge (ESD): EN/IEC61000-4-2         |
| Foot Control Cord      | 3.5m (Unshielded) | Electric fast transient / burst: EN/IEC61000-4-4<br>Surge: EN/IEC61000-4-5                                                                                                                                                           |
| Power Cord             | 3.6m (Unshielded) | Voltage dips, short interruptions and voltage variations on power supply input lines: EN/IEC61000-4-11<br>Power frequency (50/60Hz) magnetic field: EN/IEC61000-4-8<br>Conducted RF: EN/IEC61000-4-6<br>Radiated RF: EN/IEC61000-4-3 |

| Recommended separation distances between portable and mobile RF communications equipment and the Primado2.                                                                                                                                                                                                                                                                                                                             |                                                                |                                    |                                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|------------------------------------|-------------------------------------|
| The Primado 2 is intended for use in an electromagnetic environment in which radiated RF disturbances are control LCD. The customer or the user of the Primado 2 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Primado 2 as recommended below, according to the maximum output power of the communications equipment. |                                                                |                                    |                                     |
| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                                                                                                                         | Separation distance according to frequency of transmitter<br>m |                                    |                                     |
|                                                                                                                                                                                                                                                                                                                                                                                                                                        | 150kHz to 80MHz<br>$d=1.2\sqrt{P}$                             | 80MHz to 800MHz<br>$d=1.2\sqrt{P}$ | 800MHz to 2.5GHz<br>$d=2.3\sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                   | 0.12                                                           | 0.12                               | 0.23                                |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                    | 0.38                                                           | 0.38                               | 0.73                                |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                      | 1.2                                                            | 1.2                                | 2.3                                 |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                     | 3.8                                                            | 3.8                                | 7.3                                 |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                    | 12                                                             | 12                                 | 23                                  |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.                                                                                                   |                                                                |                                    |                                     |
| <b>NOTE 1</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.                                                                                                                                                                                                                                                                                                                                   |                                                                |                                    |                                     |
| <b>NOTE 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                              |                                                                |                                    |                                     |

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**NSK Europe GmbH**