Store this Addendum together with your ventilator Operator's Manual.

This Addendum summarizes changes and corrections to the HAMILTON-C2 Operator's Manual. Information is organized according to the affected section of the manual.

Preface: Training

Training

Hamilton Medical AG provides online training for users via the Hamilton Medical College. See <https://www.hamilton-medical.com/E-Learning-and-Education/College.html>.

Safety message updates

The following safety messages apply to HAMILTON-C2 use:

WARNING

- Do *not* place an HMEF between the flow sensor and the patient as doing so limits the ventilator's ability to identify disconnection at the patient, including displacement of a mask or nasal interface.
- Do *not* connect nitric oxide to the oxygen inlet; it is *not* permitted to use the ventilator with nitric oxide or mixtures of nitric oxide.
- You must remove the low-pressure oxygen adapter before using high-pressure oxygen.
- The use of a pneumatic nebulizer adds gas to the ventilator breathing system, which can affect the accuracy of volume or flow measurements.

- When connecting the ventilator to primary power, it is *not* permitted to use any type of extension cable or multi-socket adapter.

CAUTION


- *If you place an additional component, such as an HMEF, between the flow sensor and the patient, the additional resistance limits the ventilator's ability to identify disconnection at the patient.*
*To correctly identify a patient disconnection, be sure to appropriately set the lower limit of the **Pressure** alarm, as well as the **Volume** alarm limits, and carefully monitor the patient's SpO₂ and PetCO₂ values, if available.*
- *Do NOT cover the ventilator or position it in such a way that the operation or performance of the ventilator is adversely affected.*
- *Do NOT place the CO₂ sensor on the patient. It can burn the skin as the sensor may reach a temperature of 46°C (115°F).*
- *If you connect a distributed alarm system (DAS) to the ventilator, the DAS cannot be relied upon for the receipt of ventilator alarms.*

NOTICE

- It is recommended that additional independent monitoring devices, including pulse oximeters measuring SpO2, be used during mechanical ventilation. The operator of the ventilator must still maintain full responsibility for proper ventilation and patient safety in all situations.
- It is not possible for additional gases to be present at the O2 sampling site, which would thereby affect gas measurements.
- Alarm conditions, including technical faults/events, that are *not* directly related to a physiological sensor (CO2) do *not* affect the function of any attached physiological sensor, including the values of any associated CO2 measurements. Real-time waveforms on the ventilator provide a method for assessing the displayed numeric values.
- The factory default alarm limit settings are set in line with the selected patient group, allowing for unattended monitoring. These settings, however, can *never* replace individual review of the patient and adjustment of alarm limits based on their condition.
- Not all features or products are available in all markets.
- Product description and order number may differ depending on region.

Section 1.3.1 Breathing circuits and accessories


Figure 1. HAMILTON-C2 with accessories

- | | |
|---------------------------------|---------------------|
| 1 Support arm and infusion arm | 4 Breathing circuit |
| 2 Display and controls | 5 Humidifier |
| 3 Breathing circuit connections | 6 Trolley |

Section 2.11 Working with the trolley

The trolley warning labels described in Table 2-3 have been updated as follows:

Label	Description
	Make sure the wheel brakes are unlocked when moving the trolley.
	Do not lean on the trolley.
	Weight The maximum safe working load applies to a stationary properly load-balanced trolley.

Section 2.11.2 Preparing the trolley for intrahospital transport

WARNING

- Only the components listed in this section are approved for intrahospital transport.
- Use of additional items, such as a support arm, can result in the trolley tipping over.
- The ventilator must be attached to the trolley. Ensure the device is securely attached before use.

If using a HAMILTON-C2 trolley, the ventilator and its components, as well as the trolley, **must be** configured and positioned as follows during transport within the hospital:

- The oxygen cylinders must be securely attached to the trolley.
- Only the following components are allowed to be connected during transport:
 - Breathing circuit
 - Flow sensor
 - CO2 sensor (mainstream or sidestream)
 - O2 cylinder
 - Aerogen nebulizer
 - Humidifier
 - Water bottle
 - Infusion arm (water bottle holder)

Section 8.5 Alarm troubleshooting table

Note that the alarm system is always active when the ventilator is turned on. For additional details about alarms, see Chapter 8 of your ventilator *Operator's Manual*.

Table 1. Alarms and other messages

Alarm	Definition	Action needed
Battery low	The low battery alarm has different levels of priority, depending on how much charge is left, and which power supply is in use. High priority. The ventilator is running on battery power, and the battery charge is critically low. You have a minimum of 5 minutes operating time left. If the high-priority Battery low alarm occurs when starting up the ventilator, you may have less than 5 minutes of operating time remaining. Medium priority. The ventilator is running on battery power and the battery charge is low. Low priority. The ventilator is running on primary power and the battery charge is low.	• Connect the ventilator to a primary power source. • Install charged battery. • If necessary, be prepared to provide alternative ventilation.
High pressure	<i>High priority, Low after Audio pause is activated.</i> The measured inspiratory pressure exceeds the set high Pressure alarm limit. The ventilator immediately closes the inspiratory valve to stop gas flow to the patient and opens the expiratory valve to reduce pressure to the PEEP/CPAP level. If the pressure reaches 15 cmH2O above the high Pressure alarm limit for longer than 5 seconds, the ventilator opens the release valve. If the pressure reaches 15 cmH2O above the high Pressure alarm limit for longer than 7 seconds, the ventilator enters the Ambient state.	• Check patient condition. • Adjust the Pressure alarm limit. • Check the artificial airway of the patient for kinks and occlusions. • Check the breathing circuit limbs and flow sensor tubes for kinks and occlusions. • Provide alternative ventilation once the ventilator enters the Ambient state.

Section 9.5 Suctioning maneuver

The following section provides information about performing closed-suctioning maneuvers.

NOTICE

When performing a closed suctioning maneuver, follow your institution's protocols.

Verify alarm limit settings and consider whether O2 enrichment should be used prior to performing a closed-suctioning maneuver.

If the Suctioning tool is enabled on your device, ensure O2 enrichment is not active when performing the closed-suctioning maneuver.

When performing a closed-suctioning maneuver, ventilation continues and the current settings do not need to be adjusted.


You can perform a closed-suctioning maneuver with the following pressure-controlled ventilation modes: APVcmv, APVsimv, PCV+, PSIMV+, DuoPAP, APRV, SPONT, or ASV.

Section A.1 Physical characteristics

Table 2. HAMILTON-C2 physical characteristics

Dimension	Specifications
Weight	Ventilation unit: 9.5 kg (21 lb) Ventilation unit with standard trolley: 37 kg (81.6 lb) The trolley can accommodate a maximum safe working load of 80 kg (176 lb). ¹
Dimensions	See Figure 2.

Figure 2. HAMILTON-C2 dimensions

¹ The maximum safe working load applies to a stationary, properly load-balanced trolley.

Section A.4 Electrical specifications

Table 3. Electrical specifications

Element	Specifications
Battery	Hamilton Medical provides a high-capacity ² battery. An optional second battery is available.
	Electrical specifications: 14.4 V DC, 6.8 Ah, 98 Wh ²
	Type: Lithium-ion, supplied by Hamilton Medical only
	Recharge time: While ventilator is connected to primary power, approximately 2.45 h to fully recharge one battery, approximately 5 h to fully recharge two batteries. At temperatures over 43°C, the charge time is doubled (a minimum of 5 hours to charge one battery, 10 hours to charge two).
	Storage: -20°C to 50°C, ≤ 95% relative humidity. The storage location should be free from vibration, dust, direct sunlight, moisture, and corrosive gases, and with a recommended temperature range < 21°C. Extended exposure to temperatures above 45°C can degrade battery performance and life.

² Battery revision 3 and later.

Section A.6 Monitored parameters

The monitored parameters displayed on the ventilator are rounded to the nearest whole number, when required.

Waveforms displayed on the ventilator are not filtered and represent the actual monitored values.

The following table lists corrections and updates to the monitoring parameter specifications.

Table 4. Monitored parameters, ranges, and accuracy

Parameter (units)	Range: Adult/Ped	Range: Neonatal	Accuracy ³
Pinsp ⁴ (cmH ₂ O)	0 to 50	--	± (2 cmH ₂ O + 4% of actual reading)
PetCO ₂ ⁵ (mmHg)	0 to 150	0 to 150	CO ₂ (BTPS): 0 to 40 mmHg: ±2 mmHg 41 to 70 mmHg: ±5% of reading 71 to 100 mmHg: ±8% of reading 101 to 150 mmHg: ±10% of reading For sidestream CO ₂ sensor above 80 b/min: ±12% of reading
Vt/IBW (ml/kg)	2 to 20	--	--
Vt/Weight (ml/kg)	--	2 to 20	--

³ The stated accuracy includes the tolerance interval for each measurement.

⁴ Inspiratory pressure displayed in the Vent Status panel.

⁵ If the PetCO₂ value falls below 10 mmHg, the display shows dashes (---).

Sections B.3 and B.5 APVcmv and APVsimv modes

NOTICE

The minimum inspiratory pressure ($P_{peak} - PEEP$) in APVcmv and APVsimv modes is 5 cmH₂O. Be aware that a small set tidal volume with high lung compliance may lead to higher-than-expected tidal volumes.

Section B.5.2.2 IntelliSync/PSync

The name IntelliSync has been changed to PSync.

Chapter G Parts and Accessories

Corrected part numbers for the HAMILTON-C2 trolley and its accessories

Table 5. Ventilator trolley parts and accessories

Description	PN
Trolley	160170
Humidifier mount	160091
Oxygen cylinder holder	160090

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