NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to www.agilent.com/chem.

Application Note 02110

Agilent Technologies

# Direct Aqueous Injection and Automated On-Line Solid Phase Extraction (SPE) of Phenoxy Herbicides in Drinking Water Using the Varian 320-MS Triple Quadrupole LC/MS/MS

Richard Robinson and Tiffany Payne Varian, Inc.

#### Introduction

Phenoxy herbicides are members of a class of chemicals related to the growth hormone indoleacetic acid. Contact with broadleaf plants induces rapid, uncontrolled growth, killing the affected plants. Phenoxy herbicides are widely used to selectively destroy broad-leaved weeds amongst wheat or corn crops.

The best-known phenoxy herbicides are 2,4-dichlorophenoxyacetic acid (2,4-D) and 2,4,5-trichlorophenoxyacetic acid (2,4,5-T, shown in Figure 1). 2,4,5-T itself is only moderately toxic, with oral LD<sub>50</sub> of 389 mg/kg. However, the manufacturing process for 2,4,5-T contaminates this chemical with trace amounts of 2,3,7,8-tetrachlorodibenzo-p-dioxin (TCDD). TCDD is reported to be extremely toxic to humans. 2,4-D has not been reported as being contaminated with TCCD, however it is one of the most widely used herbicides in the world and is a regulated drinking water contaminant in the United States and other countries.



Figure 1. Molecular structure of the phenoxy herbicide 2,4,5-trichloropenoxyacetic acid (2,4,5-T).

The trace-level quantitation of herbicides in drinking water requires analytical speed, selectivity, and sensitivity. Fully online automated SPE analysis provides sample preparation and detection in a single step, saving time and labor for the analyst. In addition, the "closed" nature of an on-line system prevents contamination of the sample and/or losses due to evaporation. Further, in on-line SPE analysis, the complete sample — rather than an aliquot of the sample — is analyzed, which improves the ability for detection of the analyte. As a result, water sample sizes are reduced, saving sample shipping costs and storage space. A method is detailed below using the Varian 320-MS LC/MS/MS to provide automated direct aqueous injection and on-line SPE for routine surveillance of 17 acid herbicides in drinking water. A limit of detection (LOD) of less than 20 ng/L for all compounds studied is reported with a five point linear calibration from 20 to 500 ng/L.

Reproducibility over 10 injections at less than 11% RSD is obtained at the routine surveillance level of 100 ng/L. These typical results demonstrate the performance of 320-MS LC/MS/MS method for monitoring drinking water.

Instrumentation

- Varian 320-MS LC/MS/MS with ESI Source (Varian Part Number 0293795212)
- Varian 212-LC Binary Solvent Delivery Modules (Varian Part Number LCMSHPLC02)
- Varian ProStar<sup>™</sup> 210 Loading Pump (Varian Part Number 0393550001)
- 5 mL Ti Pump Head for 210 Pump (Varian Part Number 0393594091)
- 8700 PSI Ti Pressure Module (Varian Part Number 0393552501)
- Varian ProStar 410 AutoSampler (Varian Part Number 0393529001)
- Varian ProStar 510 Column Oven (Varian Part Number 0393572102)

#### **Sample Preparation**

Samples were prepared in a North Oxfordshire drinking water matrix. Sodium thiosulfate (18 mg) was added to 1 L of tap water. Samples were made up in this water and the water acidified with 0.1% with formic acid.

#### **HPLC Conditions**

Trapping Column:	Pursuit <sup>®</sup> XRs C18, 5 $\mu$ m, 50 x 2.0 mm (Varian Part Number A6000050X020)
Analytical Column:	Polaris $^{\rm \tiny TM}$ C18, 3 $\mu m$ , 100 x 2.0 mm (Varian Part Number A2001100X020)

#### Varian 212-LC Pumps:

Solvent A: 0.02% formic acid in water Solvent B: 0.02% formic acid in acetonitrile

### Varian ProStar 210 Pump:

Solvent C: 0.1% formic acid in water

#### Programs:

	Time			Flow
Analytical LC:	(min:sec)	%A	%B	(mL/min)
	00:00	100	0	0.250
	07:00	100	0	0.250
	17:00	0	100	0.250
	21:00	0	100	0.250
	21:06	100	0	0.250
	26:00	100	0	0.250
	Time		Flow	
Loading Pump LC:	(min:sec)	%C	(mL/min)	
	00:00	100	1.0	
	10:00	100	1.0	
	10:10	100	0.0	
	25:00	100	1.0	
	26:00	100	1.0	
	Time			
Divert Valve Program:	(min:sec)	Load/I	Load/Inject	
	00:00	Loa	d	
	06:00	Inje	ct	
	21:00	Loa	d	

Injection Volume: 10 mL

#### Flow Diagram



Table 2. 320-MS segment parameters.					
Compound	MS/MS Transitions	Capillary Voltage (V)	Collision Energy (V)		
2,3,6-TBA	$223 \rightarrow 179$ $225 \rightarrow 181$	-30 -30	7 7		
Dicamba-d₃	224 <del>→</del> 180	-30	6.5		
Dicamba	219 → 175	-30	6.5		
Benazolin	242 <b>→</b> 170	-30	12.5		
	242 <b>→</b> 198	-30	5.5		
Fluroxypyr	$253 \rightarrow 196$ $253 \rightarrow 233$	-30 -30	12 6		
Bentazone-d <sub>6</sub>	245 → 132	-30	25		
Dantaran	239 → 132	-30	25		
Bentazone	239 <b>→</b> 197	-30	19		
<sup>13</sup> C <sub>6</sub> 2,4-D	225 <b>→</b> 167	-30	12.5		
2.4-D	219 <b>→</b> 161	-30	12.5		
2,7 0	219 <b>→</b> 125	-30	25		
MCPA-d <sub>6</sub>	205 <b>→</b> 147	-30	15		
Bromoxynil	276 <b>→</b> 276	-30	3		
MCPA	$199 \rightarrow 141$ $199 \rightarrow 155$	-30 -30	15 8.5		
2,4-DP-d <sub>6</sub>	239 <b>→</b> 164	-30	15		
2,4-DP	233 <b>→</b> 161	-30	15		
loxynil	$\begin{array}{c} 370 \rightarrow 126 \\ 370 \rightarrow 370 \end{array}$	-30 -30	22 4.5		
MCDD	213 → 141	-30	15		
IVICFF	215 <b>→</b> 143	-30	15		
2,4,5-T	253 <b>→</b> 195	-30	10		
2,4-DB	247 <b>→</b> 161	-20	7.5		
MCPB	227 <b>→</b> 141	-20	10		
2,4,5-TP	267 <b>→</b> 195	-20	10		
Dichlorophen	267 <b>→</b> 127	-30	19.5		
	$\begin{array}{c} 267 \rightarrow 267 \\ 269 \rightarrow 127 \end{array}$	-30 -30	7 19.5		
<sup>13</sup> C <sub>6</sub> PCP	$\begin{array}{c} 271 \rightarrow 271 \\ 273 \rightarrow 273 \end{array}$	-30	7 7		
РСР	$263 \rightarrow 263$ $265 \rightarrow 265$ $267 \rightarrow 267$	-30	7 7 7		

### **MS** Parameters

Ionization Mode: ESI Negative API Drying Gas: 18 psi at 200 °C API Nebulizing Gas: 65 psi Detector: 1550 V Needle: -4500 V Shield: -600 V CID Gas: 1.5 mTorr

## Results and Discussion

This LC/MS/MS method separates and quantifies 17 phenoxy herbicides at trace levels using a 26-minute run time.

Figures 2 and 3 show overlaid chromatograms of 2,4-D and bentazone at concentrations ranging from 20 to 500 ng/L.



Figure 2. Overlaid chromatogram showing injections of 2,4–D from 20 to 500 ng/L with the automated on-line SPE technique.



Figure 3. Overlaid chromatogram showing injections of bentazone from 20 to 500 ng/L.

Calibration curves for all analytes were found to be linear from 20 to 500 ng/L. Figure 4 shows the calibration curve for bentazone, with an r<sup>2</sup> value equal to 0.9997 from 20 to 500 ng/L. Table 3 shows the r<sup>2</sup> values for all calibration curves as well as the RSD values for 10 replicate injections at 100 ng/L.



Figure 4. Standard calibration curve for bentazone on the 320-MS,  $r^2 = 0.9997$ .

Table 3. Calibration curve r<sup>2</sup> values and %RSD values for 10 replicate injections.

	10 Deviliante	Calibratian
	TO Replicate	Calibration
Compound	%RSD	Curve (r <sup>2</sup> )
2,3,6-TBA	7.54	0.9978
Dicamba	10.42	0.9993
Benazolin	6.72	0.9992
Fluroxypyr	10.59	0.9989
Bentazone	5.44	0.9997
2,4-D	6.15	0.9994
Bromoxynil	8.91	0.9993
MCPA	6.17	0.9993
2,4-DP	5.15	0.9999
loxynil	7.99	0.9977
MCPP	5.83	0.9995
2,4,5-T	7.97	0.9992
2,4-DB	5.44	0.9998
MCPB	10.49	0.9984
2,4,5-TP	6.31	0.9989
Dichlorophen	5.62	0.9991
РСР	4.51	0.9987

#### Conclusion

This method, including on-line SPE and MS/MS provides an accurate and sensitive analysis for phenoxy herbicides. The calibration curves were found to be linear from 20 to 500 ng/L with r<sup>2</sup> values ranging from 0.9977 to 0.9999, and 10 replicate injections at 100 ng/L produced %RSD values ranging from 4.51% to 10.59%. The on-line technique reduces labor, sample handling, and provides very precise and accurate results.

ProStar, Pursuit, Polaris, Varian and the Varian logo are trademarks or registered trademarks of Varian, Inc. in the U.S. and other countries. © 2009 Varian, Inc.

These data represent typical results. For further information, contact your local Varian Sales Office.

> Varian, Inc. www.varianinc.com North America: 800.926.3000, 925.939.2400 Europe The Netherlands: 31.118.67.1000 Asia Pacific Australia: 613.9560.7133 Latin America Brazil: 55.11.3238.0400

NOTICE: This document contains references to Varian. Please note that Varian, Inc. is now part of Agilent Technologies. For more information, go to www.agilent.com/chem.

Agilent Technologies