

Effective Polymer Analysis

Agilent PLgel 5 μ m MIXED-D Columns

Technical Overview

Introduction

The PLgel 5 μ m MIXED-D column is specifically designed for the analysis of polymers, paints and resin systems where material above 400,000 MW is unlikely to be present. High pore volume, concentrated in this operating range, combined with the 5 μ m efficiency, provides excellent resolution for low MW polymers and oligomers. Two, or even three, PLgel 5 μ m MIXED-D columns are the perfect replacement for the popular 10⁴/500 Å or 10⁴/10³/500 Å/100 Å column combinations.



Agilent Technologies

Polymers produced by condensation polymerization rarely contain material above 400,000 molecular weight and are therefore ideally analyzed using PLgel 5 μ m MIXED-D columns.

Column: 2 x PLgel 5 μ m MIXED-D, 300 x 7.5 mm
(part number PL1110-6504)
Eluent: THF
Flow Rate: 1.0 mL/min
Detection: RI




Figure 1. Analysis of condensation polymers using a PLgel 5 μ m MIXED-D column set

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 21, 2010

SI-01829



Agilent Technologies