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CONTINUATION

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Description in the manual is based on the newest products. Owing to improvement or other changes, contents of manual may differ from practical situation. Our company will reserve the right of recension at any moment. Please kindly forgive not to notify the revised place one by one.

The copyright of this manual belongs to Wenzhou Sundoo Instruments Co., Ltd. Any companies or individuals have no right to copy or plagiarize part or whole contents of the manual without permission of our company.

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MANUFACTURER: WENZHOU SUNDOO INSTRUMENTS CO.,LTD

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Add: Building E, No. 10 High-tech Industrial Zone WENZHOU CHINA

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SVM-210 metallurgical microscope measuring system includes video digital microscope and professional image measuring software. It is integrated traditional optic and digital video technology. The microscope is inspected with not only optical inspection but also connect with video device, that provide the image for inspection of many persons in same time. The software can be used in catching the image, saving the image, measuring sample, printing data, and so on. The microscope is applied with metallurgic inspection in mining industry, research, and teaching especially, It is used in metallurgic inspection of tiny granule, powder and silica wafer.

装箱清单
Package List

部件名称 Name		数量 Quantity
落射光照器 Inside illumination		1
镜座 Lens base		1
平场消色差物镜 Flat- field Achromatic object lens	4x	1
	10x	1
	40x	1
平场目镜 Flat-field ocular	16x	2
载物板 Carry object platform		1
滤色片、淡蓝、淡黄 Achromatic lens ,Light blue, Light yellow		1
装箱单 Package list		1
CCD成像系统 CCD Photographic System		1
1A保险管(备件) Fuse 1A		1
检验合格证 Certificate		1
使用说明书 User' s Manual		1

Maintains

1. When the microscope is not used
Please fix the protector of the lens and cover the dustproof cover. You are advised that stock the lens and base of lens into damp proof container, and place drier and prevent them to be mould.
2. When the microscope is cleaned
Don't touch the optical component with hands. Use soft cotton or brush to clean the lens. You can clean the obstinate dirt as fingerprint, grease and so on with the cotton or lens tissue, which tipped alcohol.
3. You can blow the dust on the lens with blow bloom.
The design of products is being improved, specification and subject to change without notice.

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c. Achromatic lens

There are light blue, light yellow and light green achromatic lens applied with the microscope. The achromatic lens can off yellow, green light, which is sensitivity of eyes.

7. Changing the light bulb.

Please put off the power before changing light bulb. After the lamp be cold, remove the slant light cover, disassembled the halogen bulb and change it. Take off the window on the baseboard, you can change the bulb (halogen and fluoresce light) through the window.

8. Adjusting the center of bulb

The center of bulb has been adjusted, before the microscope is sold. That can keep illumination be nice. You must adjust it after change the bulb. Select the max of field of view raster and diaphragm raster, and select the lower times of ocular, adjust the focus. Loose the screw on the back of the lamp room, and move it up or down to move the bulb up and down. Turn in or turn out the screw on the right of back of lamp room. Then fix the bulb on correct location for nice illumination.

Optional appurtenance

1. AV/VGA, XGA-transfer card + CRT monitor

The XGA-transfer card can provide hi-definition image, display on CRT, without flash and with bright and flamboyant picture.

Connecting Step:

- Attach the power, Connect the CRT's signal plus with the transfer card's connector.
- Inert MMI signal provide with the instrument into the reciprocal connector on the transfer card, attach Ac cable With Photographic lens' output connector, insert another AV terminal into the connector on the yellow MMI cable.

2. Monitor, Computer, Video Card, Measuring Software

The Measuring Software can provide measuring function, that can measure width, radii, diameter, distance between two lines and angle.

Connecting Step:

Please insert the Video Card into the PCI connector, install Video Card's drive program and measuring program. Insert the AV plug into the CCD output connector, connect another AV terminal with input connector of VAG.

In the course of inspection, select a piece of metal for the sample thickness less than 35mm, and put it on the carried object platform. The surface of the metal will be treated with rubdown, polishing and corrode. And move the sample with the motive platform, locate the inspected area under the object lens correctly.

3. Use the ocular base.

a. Pull the board on the ocular base to adjust the distance of pupil for the inspector. The image in the field of view must be overlapping, and record the value for future inspection.

b. Turn the angle of view adjust loop on the right ocular, make the line correct with your pupil-gap value. The length of cylinder is 160mm.

c. Inspect on left ocular with left eye, and adjust the focus adjust knob still the image be clear.

d. Inspect on right ocular with right eye, and adjust the nice focus adjust knob still the image be clear.

4. The binocular base are adjusted by the way, that can keep the same focus in different magnification. Find the inspected object with 10X object lens. Then select the magnification adjusting knob4 to get the higher magnification for inspecting the details of object. The image is little clear. When the magnification is changed from lower to higher, turning the fine adjusting focus knob slightly, you can get clear image.

5. The strain ness of the focus adjusting knob is adjusted by the strain knobs. The strainness has been adjusted in the factory.

6. adjusting brightness

a. diaphragm raster

Fiddle the role of diaphragm raster to change the magnitude of diaphragm raster. To get good definition and contrast of object lens, you will change the diaphragm raster to smaller. For good brightness, you should increase the brightness of illumination.

b. Field of view raster

Fiddle the role of field of view raster to change the magnitude of diaphragm raster. Select the small field of view raster to eliminate the interfere beam. The small field of view raster is for the contrast. You can adjust the screws on illumination base to change the center of the field of view raster. In general condition, you will adjust the field of view raster when you select new magnification.

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bright light will affect the image of samples.

d.Prevent operate the microscope in the dusty environment.The microscope should be covered with dustproof cover,when it is not used.

e.The microscope should be place on the area without shock.

6.Operate in carefulness

a.Microscope is precision optical instrument. The violent shock and obstinate operation will damage the microscope. especially, do not turn the left and right adjust knobs in against direction. That will make the focus adjusting in malfunction.

b.Do not touch optical components with hands, which will cause the image to be blurred.

Installation(Refer to instructure diagram)

1. Please take binocular ocular base and main body out of the package.

2. Take the inside illumination out of the foam. Take off the dustproof cover, lean the V connector of inside illumination little, put it into to the connector of the main body, which is on the top of the main body. And fasten the screw on the left side of the connector to fix the illumination. Then connect the power wire with the power output connector on the base.

3. Fix the ocular lens base on the illumination base, and fasten the screw on the right of the illumination base to tight the lensbase.

4. Pull the 10X or 16X ocular into the lens base.

5. Fix three object lens into the turntable, from the lower to higher magnification. Keep the object lens tight with the screw of turntable.

6. Place the carrier object board onto the carrier object platform and fix it with clamps.

7. Please insert the Achromatic lens into the slots on the illumination base.

Operation

1. Connect power

Make sure the voltage of power supply is same as the rate input voltage. Then connect the microscope with the power supply. Turn the brightness controller on the base of main body. The lamp will be bright. Adjust the controller to get the suitable brightness.

2. Locate the sample

Cautions

1. Do not disassemble the microscope.
To reduce the risk of electric shock, do not remove the parts, which are not described to be removed in the manual. If the microscope is in trouble, refer servicing of quality service personnel.
2. Check the Power
 - a. The rate input voltage is signed on the pedestal of Microscope. The input voltage is according with the rate input voltage. If the voltage of power supplied is not accord with the rate voltage, please ask the experience personnel to open the cover on the pedestal, turn the switch, which have a symbol 110V~220V. Select the voltage same as the power supply, then fix the cover on the pedestal. Or refer service with the service personnel. That will damage the microscope; the microscope is operated with the rate voltage.
 - b. The socket must be with earth.
3. Prevent the Fire or Scald
 - a. During the bulb provides illumination, it is in high temperature. Pay more attention to do not touch the bulb to prevent scald.
 - b. When you change the bulb, please wait for the bulb is cold.
 - c. To prevent fire, do not take it in close to the bulb.
4. Change the Light Bulb
 - a. Please cut off the power before change bulb, prevent electric shock or damage the microscope.
 - b. The bromine tungsten lamp is only applied with the microscope, specification 6V20W
 - c. When you change the bulb, please wear gloves or pack the bulb with cover. Do not touch the glass part of the bulb with hand. After change the bulb, please sweep the finger mark or dirty with clean soft cloth dip in pure alcohol, on the bulb. Otherwise, the bulb should be easy to crack or reduce the brightness.
5. Delivery and Stock.
 - a. The power must be cutoff, before remove the microscope.
 - b. Operation Environment
Temperature: 0°C-40 °C; Relative Humidity: 85%.
 - c. The microscope is prevented the direct sunlight and bright light. The

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
Picture 2


Parameters

- 1.Parameters View p17 form.
- 2.Conjuquate distance of object and image;
- 3.Object lens(view P2 form 1) ;
- 4.Ocular (view P2 form 2);
- 5.Length of focus adjust: 40mm graduation of nice focus adjusting is 0.002mm;
- 6.The area of mechanic carry platform is 130×122mm²Distance of horizontal motion is 70mm distance of vertical motion is 30 mm;
- 7.Illumination : Input voltage :220V or 110V;
Output : 6V 20W halogen bulb
inside brightness adjusting power supply.

Summary

SVM-210 Metallurgical Microscope Measuring System includes Video Digital Microscope and professional image measuring software. It is integrated Traditional Optic and Digital Video Technology. The brightness of Inside Same Axis Illumination is continuous adjusted. The Electric Magnifications is 400x-4000x and the Optic Magnifications is 40X-400X. The accuracy, which is 300 Nanometers, make the imagination be more clear. The professional measuring software can capture the image quickly and exactly. It has the picture measurement, management, treatment and so on.

Principle

Picture 1

Diagram 1 is the optical principle of the microscope, the emitted light after gathered by converging lens and reflect by transparent reflect lens, shine to sample by affixture object lens, the light emit from sample pass object lens, affixture object lens, pass transparent reflect lens, through affixture object lens to transfer image prism, the image be in ocular and be magnified again, inspector will inspect the magnified image. When the image transfer to prism, the beam will be divided into 2 beams, one beam is refracted to vertical direction for CCD camera, another beam is refracted to the 30° angle direction for person's eye observation.

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