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The screenshot shows the 'Light Curve' dialog box in the Dual PAM v1.19 software. The main window displays a graph of ETR (umol/m²/s) versus PAR (μE). The 'Light Curve' dialog is open on the right, showing a table of 20 steps with PAR, Intensity, and Time values. The 'Model' section is set to 'EP', and the 'Fit' button is highlighted. The 'Edit' button is also highlighted. The 'Start' button is visible at the bottom of the dialog.

Step	PAR	Intens.	Time/10s
1	5	1	3
2	13	2	3
3	31	3	3
4	89	5	3
5	167	7	3
6	209	8	3
7	325	10	3
8	496	12	3
9	754	14	3
10	1173	16	3
11	1450	17	0
12	1804	18	3
13	496	12	0
14	325	10	3
15	209	8	3
16	167	7	3
17	89	5	3
18	31	3	3
19	13	2	3
20	5	1	3

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1. Nishikawa Y, Yamamoto H, Okegawa Y, et al. PGR5-dependent cyclic electron transport around PSI contributes to the redox homeostasis in chloroplasts rather than CO₂ fixation and biomass production in rice[J]. *Plant and Cell Physiology*, 2012, 53(12): 2117-2126.
2. Alric J, Lavergne J, Rappaport F. Redox and ATP control of photosynthetic cyclic electron flow in *Chlamydomonas reinhardtii* (l) aerobic conditions[J]. *Biochimica et Biophysica Acta (BBA)-Bioenergetics*, 2010, 1797(1): 44-51.
3. Suzuki K, Ohmori Y, Ratel E. High root temperature blocks both linear and cyclic electron transport in the dark during chilling of the leaves of rice seedlings[J]. *Plant and Cell Physiology*, 2011, 52(9): 1697-1707.
4. Wang C, Yamamoto H, Shikanai T. Role of cyclic electron transport around photosystem I in regulating proton motive force[J]. *Biochimica et Biophysica Acta (BBA)-Bioenergetics*, 2015, 1847(9): 931-938.

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1. Benson S L, Maheswaran P, Ware M A, et al. An intact light harvesting complex I antenna system is required for complete state transitions in Arabidopsis[J]. *Nature Plants*, 2015, 1: 15176.
2. Pinnola A, Cazzaniga S, Alboresi A, et al. Light-Harvesting Complex Stress-Related Proteins Catalyze Excess Energy Dissipation in Both Photosystems of *Physcomitrella patens*[J]. *The Plant Cell*, 2015, 27(11): 3213-3227.
3. Pinnola A, Dall'Osto L, Gerotto C, et al. Zeaxanthin binds to light-harvesting complex stress-related protein to enhance nonphotochemical quenching in *Physcomitrella patens*[J]. *The Plant Cell*, 2013, 25(9): 3519-3534.
4. Iwai M, Yokono M, Kono M, et al. Light-harvesting complex Lhcb9 confers a green alga-type photosystem I supercomplex to the moss *Physcomitrella patens*[J]. *Nature Plants*, 2015, 1(2).
5. Allahverdiyeva Y, Suorsa M, Rossi F, et al. Arabidopsis plants lacking PsbQ and PsbR subunits of the oxygen-evolving complex show altered PSII super-complex organization and short-term adaptive mechanisms[J]. *The Plant Journal*, 2013, 75(4): 671-684.
6. Pietrzykowska M, Suorsa M, Semchonok D A, et al. The light-harvesting chlorophyll a/b binding proteins Lhcb1 and Lhcb2 play complementary roles during state transitions in Arabidopsis[J]. *The Plant Cell*, 2014, 26(9): 3646-3660.
7. Belgio E, Kapitonova E, Chmeliov J, et al. Economic photoprotection in photosystem II that retains a complete light-harvesting system with slow energy traps[J]. *Nature communications*, 2014, 5.
8. Sukhov V, Surova L, Sherstneva O, et al. Influence of variation potential on resistance of the photosynthetic machinery to heating in pea[J]. *Physiologia plantarum*, 2014, 152(4): 773-783.
9. Xu, D., Y. Wang, et al. (2014). Long-Term Experiment on Physiological Responses to Synergetic Effects of Ocean Acidification and Photoperiod in the Antarctic Sea Ice Algae *Chlamydomonas* sp. ICE-L. *Environmental Science & Technology* 48(14): 7738-7746. (黄海水产研究所)
10. Xu, L., Q. Wang, et al. (2014). Improvement of hydrogen yield of lba-transgenic *Chlamydomonas reinhardtii* caused by increasing respiration and impairing photosynthesis. *International Journal of Hydrogen Energy*. (上海师范大学)
11. Yuan, X., L. Zhang, et al. (2014). Photosynthetic Physiological Response of *Radix Isatidis* (*Isatis indigotica* Fort.) Seedlings to Nicosulfuron. *PLoS ONE* 9(8): e105310. (山西农业大学)
12. Zhang, G., Y. Liu, et al. (2014). Exogenous Calcium Alleviates Low Night Temperature Stress on the Photosynthetic Apparatus of Tomato Leaves. *PLoS ONE* 9(5): e97322. (沈阳农业大学)
13. Zhang, J., F. Gao, et al. (2014). NdhP Is an Exclusive Subunit of Large Complex of NADPH Dehydrogenase Essential to Stabilize the Complex in *Synechocystis* sp. Strain PCC 6803. *Journal of Biological Chemistry* 289(27): 18770-18781. (上海师范大学)
14. Zhang, X., S. Mou, et al. (2014). Roles of the transthylakoid proton gradient and xanthophyll cycle in the non-photochemical quenching of the green alga *Ulva linza*. *Estuarine, Coastal and Shelf Science*. (黄海水产研究所)
15. Zhao, J., F. Gao, et al. (2014). Deletion of an electron donor-binding subunit of the NDH-1 complex, NdhS, results in a heat-sensitive growth phenotype in *Synechocystis* sp. PCC 6803. *Chinese Science Bulletin*: 1-7. (上海师范大学)
16. Zhao, J., F. Gao, et al. (2014). NdhO, a Subunit of NADPH Dehydrogenase, Destabilizes Medium Size Complex of the Enzyme in *Synechocystis* sp. Strain PCC 6803. *Journal of Biological Chemistry* 289(39): 26669-26676. (上海师范大学)
17. Zhao, P., W. Gu, et al. (2014). Silicon enhances the growth of *Phaeodactylum tricornutum* Bohlin under green light and low temperature. *Scientific Reports* 4. (中科院海洋所)
18. Zhou, W., L. He, et al. (2014). *Pyropia yezoensis* can utilize CO₂ in the air during moderate dehydration. *Chinese Journal of Oceanology and Limnology* 32: 358-364. (中科院海洋所)

19. An, Y., X. Han, et al. (2014). Poplar GATA transcription factor PdGNC is capable of regulating chloroplast ultrastructure, photosynthesis, and vegetative growth in *Arabidopsis* under varying nitrogen levels. *Plant Cell, Tissue and Organ Culture (PCTOC)*: 1-15. (北京林业大学)
20. An, Y., S. Wang, et al. (2014). O₂ evolution and cyclic electron flow around photosystem I in long-term ground batch culture of *Euglena gracilis*. *Advances in Space Research*. (中科院武汉水生所)
21. Deng, C.-N., D.-Y. Zhang, et al. (2014). Toxic effects of erythromycin on photosystem I and II in *Microcystis aeruginosa*. *Photosynthetica*: 1-7. (云南师范大学)
22. Deng, C., X. Pan, et al. (2014). Cu²⁺ Inhibits Photosystem II Activities but Enhances Photosystem I Quantum Yield of *Microcystis aeruginosa*. *Biological Trace Element Research* 160(2): 268-275. (云南师范大学)
23. Deng, C., X. Pan, et al. (2014). Influence of ofloxacin on photosystems I and II activities of *Microcystis aeruginosa* and the potential role of cyclic electron flow. *Journal of bioscience and bioengineering*. (云南师范大学)
24. Dong, Y., C. Wang, et al. (2014). A novel bHLH transcription factor *PebHLH35* from *Populus euphratica* confers drought tolerance through regulating stomatal development, photosynthesis and growth in *Arabidopsis*. *Biochemical and Biophysical Research Communications*. (北京林业大学)
25. Fan, J., Y. Cui, et al. (2014). Lipid accumulation and biosynthesis genes response of the oleaginous *Chlorella pyrenoidosa* under three nutrition stressors. *Biotechnology for biofuels* 7(1): 17. (华东理工大学)
26. Fan, X., H. Qiao, et al. (2014). Short-term retention of kleptoplasty from a green alga (*Bryopsis*) in the sea slug *Placida* sp. YS001. *Biologia* 69(5): 635-643. (黄海水产研究所)
27. Fan, X., D. Xu, et al. (2014). The effect of nutrient concentrations, nutrient ratios and temperature on photosynthesis and nutrient uptake by *Ulva prolifera*: implications for the explosion in green tides. *Journal of Applied Phycology* 26(1): 537-544. (黄海水产研究所)
28. Gao, S., W. Gu, et al. (2014). Desiccation enhances phosphorylation of PSII and affects the distribution of protein complexes in the thylakoid membrane. *Physiologia plantarum*. (中科院海洋所)
29. Gao, S., Z. Zheng, et al. (2014). Photosystem I shows a higher tolerance to sorbitol-induced osmotic stress than photosystem II in the intertidal macro-algae *Ulva prolifera* (Chlorophyta). *Physiologia plantarum*. (中科院海洋所)
30. Gao, Z., D. Xu, et al. (2014). The green tide-forming macroalga *Ulva linza* outcompetes the red macroalga *Gracilaria lemaneiformis* via allelopathy and fast nutrients uptake. *Aquatic Ecology* 48(1): 53-62. (黄海水产研究所)
31. Huan, L., S. Gao, et al. (2014). Specific photosynthetic and morphological characteristics allow macroalgae *Gloiopeltis furcata* (Rhodophyta) to survive in unfavorable conditions. *Photosynthetica*: 1-7. (中科院海洋所)
32. Huan, L., X. Xie, et al. (2014). Positive Correlation Between PSI Response and Oxidative Pentose Phosphate Pathway Activity During Salt Stress in the Intertidal Macroalgae. *Plant and Cell Physiology*: pcu063. (中科院海洋所)
33. Jia, X.-H., P.-P. Zhang, et al. (2014). Regulation of *pepc* gene expression in *Anabaena* sp. PCC 7120 and its effects on cyclic electron flow around photosystem I and tolerances to environmental stresses. *Journal of Integrative Plant Biology*: n/a-n/a. (上海海洋大学)
34. Lei, Y.-b., Y.-l. Zheng, et al. (2014). Different responses of photosystem I and photosystem II in three tropical oilseed crops exposed to chilling stress and subsequent recovery. *Trees* 28(3): 923-933. (中科院西双版纳植物园)
35. Ou, H.-S., C.-H. Wei, et al. (2013). Principal component analysis to assess the efficiency and mechanism for enhanced coagulation of natural algae-laden water using a novel dual coagulant system. *Environmental Science and Pollution Research* 21(3): 2122-2131. (华南理工大学)
36. Shi, Y., R. Wang, et al. (2014). Molecular Cloning and Functional Characterization of the Lycopene ϵ -Cyclase Gene via Virus-Induced Gene Silencing and Its Expression Pattern in *Nicotiana tabacum*. *International journal of molecular sciences* 15(8): 14766-14785. (郑州大学)
37. Song, Y., Q. Diao, et al. (2014). Polyamine metabolism and biosynthetic genes expression in tomato (*Lycopersicon esculentum* Mill.) seedlings during cold acclimation. *Plant Growth Regulation*: 1-12. (沈阳农业大学)

38. Teng, K., J. Li, et al. (2014). Exogenous ABA induces drought tolerance in upland rice: the role of chloroplast and ABA biosynthesis-related gene expression on photosystem II during PEG stress. *Acta Physiologiae Plantarum*: 1-9. (河南农业大学)
39. Wang, B., Z. Zhang, et al. (2014). Cellular Capacities for High-Light Acclimation and Changing Lipid Profiles across Life Cycle Stages of the Green Alga *Haematococcus pluvialis*. *PLoS ONE* 9(9): e106679. (厦门大学)
40. Wang, S., X. Pan, et al. (2014). PSI Showed Higher Tolerance to Sb (V) than PSII Due to Stimulation of Cyclic Electron Flow Around PSI. *Current microbiology*: 1-8. (中国科学院新疆生地所)
41. Cao, S., X. Zhang, et al. (2013). A transthylakoid proton gradient and inhibitors induce a non-photochemical fluorescence quenching in unicellular algae *Nannochloropsis* sp. *FEBS Letters* 587(9): 1310-1315. (青岛农业大学)
42. Dai, H., L. Zhang, et al. (2013). Identification of a cyanobacterial CRR6 protein, Slr1097, required for efficient assembly of NDH-1 complexes in *Synechocystis* sp. PCC 6803. *The Plant Journal* 75(5): 858-866. (上海师范大学)
43. Deng, C., D. Zhangb, et al. (2013). Toxic effects of mercury on PSI and PSII activities, membrane potential and transthylakoid proton gradient in *Microsorium pteropus*. *Journal of Photochemistry and Photobiology B: Biology* 127: 1-7. (云南师范大学)
44. Gao, S., J. Niu, et al. (2013). The physiological links of the increased photosystem II activity in moderately desiccated *Porphyra haitanensis* (Bangiales, Rhodophyta) to the cyclic electron flow during desiccation and re-hydration. *Photosynthesis Research* 116(1): 45-54. (中科院海洋所)
45. Gao, Z., D. Li, et al. (2013). Survival and proliferation characteristics of the microalga *Chlamydomonas* sp. ICE-L after hypergravitational stress pretreatment. *Icarus* 226(1): 971-979. (山东理工大学)
46. Gu, W., X. Xie, et al. (2013). Comparison of different cells of *Haematococcus pluvialis* reveals an extensive acclimation mechanism during its aging process: From a perspective of photosynthesis. *PLoS ONE* 8(7): e67028. (中科院海洋所)
47. He, Z. H., H. W. Li, et al. (2013). Comparative analysis of the chloroplast proteomes of a wheat (*Triticum aestivum* L.) single seed descent line and its parents. *J Plant Physiol* 170(13): 1139-1147. (中科院植生所)
48. Huang, W., P.-L. Fu, et al. (2013). Differences in the responses of photosystem I and photosystem II of three tree species *Cleistanthus sumatranus*, *Celtis philippensis* and *Pistacia weinmannifolia* exposed to a prolonged drought in a tropical limestone forest. *Tree physiology* 33(2): 211-220. (扬州大学)
49. Jiang, Y. P., L. F. Huang, et al. (2013). Brassinosteroids accelerate recovery of photosynthetic apparatus from cold stress by balancing the electron partitioning, carboxylation and redox homeostasis in cucumber. *Physiologia plantarum* 148(1): 133-145. (浙江大学)
50. Hu, Y.-Y., D.-Y. Fan, et al. (2013). Whole-tissue determination of the rate coefficients of photoinactivation and repair of photosystem II in cotton leaf discs based on flash-induced P700 redox kinetics. *Photosynthesis Research* 117(1-3): 517-528. (石河子大学)
51. Lv, Z., L. Wei, et al. (2013). Moderate Heat Pretreatment Alleviates the Inhibition of Photosystem II Activity Caused by the Response of *Cyanobacterial* Cells to High Red Light. *Photosynthesis Research for Food, Fuel and the Future*, Springer: 502-504. (上海师范大学)
52. Mou, S., X. Zhang, et al. (2013). Photoprotection in the green tidal alga *Ulva prolifera*: role of LHCSR and *PsbS* proteins in response to high light stress. *Plant Biology* 15(6): 1033-1039. (黄海水产研究所)
53. Ou, H.-S., C.-H. Wei, et al. (2013). Principal component analysis to assess the composition and fate of impurities in a large river-embedded reservoir: Qingcaosha Reservoir. *Environmental Science: Processes & Impacts* 15(8): 1613-1621. (华南理工大学)
54. Pfündel, E. E., C. Klughammer, et al. (2013). Deriving fluorometer-specific values of relative PSI fluorescence intensity from quenching of F₀ fluorescence in leaves of *Arabidopsis thaliana* and *Zea mays*. *Photosynthesis Research* 114(3): 189-206. (德国 WALZ)
55. Qi, M., Y. Liu, et al. (2013). Nano-TiO₂ Improve the Photosynthesis of Tomato Leaves under Mild Heat Stress. *Biological*

- Trace Element Research 156(1-3): 323-328. (沈阳农业大学)
56. Qu, C., C. Liu, et al. (2013). Improvement of Cerium on Photosynthesis of Maize Seedlings Under a Combination of Potassium Deficiency and Salt Stress. *Biological Trace Element Research* 155(1): 104-113. (苏州大学)
 57. Shao, N., G. Y. Duan, et al. (2013). A mediator of singlet oxygen responses in *Chlamydomonas reinhardtii* and *Arabidopsis* identified by a luciferase-based genetic screen in algal cells. *The Plant Cell Online* 25(10): 4209-4226. (华南植物园)
 58. Li, D.-D., X.-L. Xia, et al. (2013). Two poplar calcineurin B-like proteins confer enhanced tolerance to abiotic stresses in transgenic *Arabidopsis thaliana*. *Biologia Plantarum* 57(1): 70-78. (北京林业大学)
 59. Lin, A., G. C. Wang, et al. (2013). Simultaneous measurements of H⁺ and O₂ fluxes in *Zostera marina* and its physiological implications. *Physiologia plantarum* 148(4): 582-589. (中科院海洋所)
 60. Liu, J., P. Wang, et al. (2013). A deficiency in chloroplastic ferredoxin 2 facilitates effective photosynthetic capacity during long-term high light acclimation in *Arabidopsis thaliana*. *The Plant Journal* 76(5): 861-874. (中山大学)
 61. Song, R., C.-Y. Zhao, et al. (2013). Effect of sulphate nutrition on arsenic translocation and photosynthesis of rice seedlings. *Acta Physiologiae Plantarum* 35(11): 3237-3243. (河南农业大学)
 62. Wang, J., M. R. Sommerfeld, et al. (2013). Combined effect of initial biomass density and nitrogen concentration on growth and astaxanthin production of *Haematococcus pluvialis* (Chlorophyta) in outdoor cultivation. *Algae* 28(2): 193-202. (中科院植物所)
 63. Wang, J. H., S. C. Li, et al. (2013). Differences in the stimulation of cyclic electron flow in two tropical ferns under water stress are related to leaf anatomy. *Physiologia plantarum* 147(3): 283-295. (云南省农业科学院)
 64. Wang, L. Y., D. Li, et al. (2013). Antisense-mediated depletion of tomato GDP-L-galactose phosphorylase increases susceptibility to chilling stress. *J Plant Physiol* 170(3): 303-314. (山东农业大学)
 65. Wang, P., J. Liu, et al. (2013). Evidence for a Role of Chloroplastic m-Type Thioredoxins in the Biogenesis of Photosystem II in *Arabidopsis*. *Plant Physiology* 163(4): 1710-1728. (中山大学)
 66. Wang, S., F. Chen, et al. (2013). Simultaneous analysis of photosystem responses of *Microcystis aeruginosa* under chromium stress. *Ecotoxicology and Environmental Safety* 88: 163-168. (中科院生地所)
 67. Wang, S., D. Zhang, et al. (2013). Effects of cadmium on the activities of photosystems of *Chlorella pyrenoidosa* and the protective role of cyclic electron flow. *Chemosphere* 93(2): 230-237. (中科院生地所)
 68. Wang, Y., D. Xu, et al. (2013). Variation of photosynthetic performance, nutrient uptake, and elemental composition of different generations and different thallus parts of *Saccharina japonica*. *Journal of Applied Phycology* 25(2): 631-637. (青岛农业大学)
 69. Wei, L., X. Li, et al. (2013). A simple approach for the efficient production of hydrogen from Taihu Lake *Microcystis* spp. blooms. *Bioresource Technology* 139: 136-140. (上海师范大学)
 70. Wang, J., D. Han, et al. (2013). Effect of initial biomass density on growth and astaxanthin production of *Haematococcus pluvialis* in an outdoor photobioreactor. *Journal of Applied Phycology* 25(1): 253-260. (中科院植物所)
 71. Xie, X., S. Gao, et al. (2013). Desiccation Induces Accumulations of Antheraxanthin and Zeaxanthin in Intertidal Macro-Alga *Ulva pertusa* (Chlorophyta). *PLoS ONE* 8(9): e72929. (中科院海洋所)
 72. Xu, D., F. Li, et al. (2013). Facilitative interactions between the green-tide macroalga *Monostroma arctium* and the red macroalga *Porphyra yezoensis*. *Journal of Experimental Marine Biology and Ecology* 444: 8-15. (黄海水产研究所)
 73. Xu, D., N. Ye, et al. (2013). Variation in morphology and PSII photosynthetic characteristics of *Macrocystis pyrifera* during development from gametophyte to juvenile sporophyte. *Aquaculture Research*. (黄海水产研究所)
 74. Xu, J., X. Zhang, et al. (2013). Activities of principal photosynthetic enzymes in green macroalga *Ulva linza*: functional implication of C₄ pathway in CO₂ assimilation. *Science China Life Sciences* 56(6): 571-580. (中科院海洋所)
 75. Yuan, X., P. Guo, et al. (2013). Safety of herbicide Sigma Broad on *Radix Isatidis* (*Isatis indigotica* Fort.) seedlings and their photosynthetic physiological responses. *Pesticide Biochemistry and Physiology* 106(1): 45-50. (山西农业大学)
 76. Zhang, Q., C. Zhou, et al. (2013). Stereoselective Phytotoxicity of HCH Mediated by Photosynthetic and Antioxidant Defense

Systems in *Arabidopsis thaliana*. PLoS ONE 8(1): e51043. (浙江大学)

77. Zhang, X., N. Ye, et al. (2013). Occurrence of the *Psb S* and *Lhc SR* products in the green alga *Ulva linza* and their correlation with excitation pressure. *Plant Physiology and Biochemistry* 70: 336-341. (黄海水产研究所)
78. Zhang, Y.-J., Q.-Y. Yang, et al. (2013). Extended leaf senescence promotes carbon gain and nutrient resorption: importance of maintaining winter photosynthesis in subtropical forests. *Oecologia* 173(3): 721-730. (中国科学院西双版纳热带植物园)
79. Zhang, Y.-M., H. Chen, et al. (2013). Nitrogen starvation induced oxidative stress in an oil-producing green alga *Chlorella sorokiniana* C3. PLoS ONE 8(7): e69225. (中科院武汉水生所)
80. Zhang, Y., K. Cao, et al. (2013). Winter Photosynthesis of Evergreen Broadleaf Trees from a Montane Cloud Forest in Subtropical China. *Photosynthesis Research for Food, Fuel and the Future*, Springer: 812-817. (中国科学院西双版纳热带植物园)

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Bilger W and Schreiber U (1990) Chlorophyll luminescence as indicator of stress-induced damage to the photo-synthetic apparatus. Effects of heat-stress in isolated chloroplasts. *Photosynth Res* 25: 161-171

Bilger W, Björkman O, Thayer SS (1989) Light-induced spectral absorbance changes in relation to photosynthesis and the epoxidation state of xanthophyll cycle components in cotton leaves. *Plant Physiol* 91: 542-551

Bilger W, Heber U and Schreiber U (1988) Kinetic relationship between energy-dependent fluorescence quenching, light scattering, chlorophyll luminescence and proton pumping in intact leaves. *Z Naturforsch* 43c: 877-887

Cruz JA, Avenson TJ, Kanazawa A, Takizawa K, Edwards GE and Kramer DA (2004) Plasticity in light reactions of photosynthesis for energy production and photo-protection. *J Exp Botany* 56: 395-406

Deamer DW, Crofts AR and Packer L (1966) Mechanisms of light-induced structural changes in chloroplasts. I. Light-scattering increments and ultrastructural changes mediated by proton transport. *Biochim Biophys Acta* 131: 81-96

Heber U (1969) Conformational changes of chloroplasts induced by illumination of leaves in vivo. *Biochim Biophys Acta* 180: 302-319

Junge W and Witt HT (1968) On the ion transport system in photosynthesis: investigations on a molecular level. *Z Naturforsch* 23b: 244-254

Kanazawa A and Kramer DA (2002) In vivo modulation of non-photochemical quenching (NPQ) by regulation of the chloroplast ATP synthase. *Proc Natl Acad Sci USA* 99: 12798-12794 Zealquest Scientific Technology Co., Ltd.

Klughammer C, Kolbowski J and Schreiber U (1990) LED array spectrophotometer for measurement of time resolved difference spectra in the 530–600 nm wave-length region. *Photosynth Res* 25: 317-327

Kramer DA and Sacksteder CA (1998) a diffused-optics flash kinetic spectrophotometer (DOFS) for measurements of absorbance changes in intact plants in the steady-state. *Photosynth Res* 56: 103-112

Ruban AV, Pascal AA, Robert B, Horton P (2002) Activation of zeaxanthin is an obligatory event in the regulation of photosynthetic light harvesting. *J Biol Chem* 277: 7785–7789

Sacksteder CA and Kramer DA (2000) Dark interval relaxation kinetics of absorbance changes as a quantitative probe of steady-state electron transfer. *Photosynth Res* 66: 145-158

Velthuys BR (1978) A third site of proton translocation in green plant photosynthetic electron transport. *Proc Natl Acad Sci USA* 76: 2765-2769

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