

	Made in USA
Catalog Number	230325
Product Name	SoloPack Gold Competent Cells
Materials Provided	SoloPack Gold competent cells, 15 single-tube transformations (25 µl each) pUC18 control plasmid (0.1 ng/µl in TE buffer), 10 µl
Certified By	****
Quality Controlled By	****
Shipping Conditions	Shipped on dry ice.
Storage Conditions	Ultrapotent cells must be placed immediately at the bottom of a -80°C freezer directly from the dry ice shipping container. Do not store the cells in liquid nitrogen. Ultrapotent cells are sensitive to even small variations in temperature. Transferring tubes from one freezer to another may result in a loss of efficiency.
Guaranteed Efficiency	$\geq 1 \times 10^8$ cfu/µg pUC18 DNA
Test Conditions	Transformations are performed both with and without plasmid DNA using 25-µl aliquots of cells and 1 µg of pUC18 control plasmid (10 pg/µl) (the pUC18 control plasmid is diluted 1:10 in sterile water) following the protocol outlined below. Following transformation, 5-µl samples of the culture are plated in duplicate on LB agar plates with 100 µg/ml ampicillin. The plates are incubated at 37°C overnight and the efficiency is calculated based on the average number of colonies per plate.
Antibiotic Resistance	SoloPack Gold competent cells are tetracycline and chloramphenicol resistant.
Genotype and Background	Tet ^r Δ(<i>mcrA</i>)183 Δ(<i>mcrCB</i> - <i>hsdSMR</i> - <i>mrr</i>)173 <i>endA1 supE44 thi-1 recA1 gyrA96 relA1 lac</i> Hte [F' <i>proAB lacI</i> ^q Δ <i>M15</i> Tn10 (Tet ^r) Amy Cam ^r]. (Genes listed signify mutant alleles. Genes on the F' episome, however, are wild-type unless indicated otherwise.) SoloPack Gold competent cells are packaged in single-reaction aliquots, providing high performance and convenience. The single-tube transformation format eliminates the need to thaw, aliquot, and refreeze cells and allows for fewer pipetting steps because the entire transformation reaction occurs in the tube supplied. SoloPack Gold cells are deficient in all known <i>E. coli</i> restriction systems (McrA-, McrCB- McrF-, Mrr-, HsdR-), preventing the cleavage of methylated DNA. The cells are also endonuclease deficient (<i>endA</i>), improving the quality of miniprep DNA, and are recombination deficient (<i>recA</i>), improving insert stability. The Hte phenotype increases the transformation efficiency of ligated and large supercoiled DNA and causes the production of large colonies. The <i>lacI</i> ^q Δ <i>M15</i> gene is present on the F' episome, allowing blue-white screening for recombinant plasmids.
Transformation Protocol	<ol style="list-style-type: none"> 1. Thaw two aliquots of SoloPack Gold cells on ice (one tube for each of the experimental and control transformations). Preheat SOC broth to 42°C. 2. When the cells have thawed, swirl the tubes to gently mix the cells. 3. Add 0.1-50 ng of the experimental DNA to one aliquot of cells. Dilute the pUC18 control DNA 1:10 with sterile dH₂O, then add 1 µl of the diluted pUC18 DNA to the other aliquot of cells. 4. Swirl the tubes gently, then incubate the tubes on ice for 30 minutes. 5. Heat-pulse the tubes in a 42°C water bath for 60 seconds. The duration of the heat pulse is critical. 6. Incubate the tubes on ice for 2 minutes. 7. Add 175 µl of preheated (42°C) SOC broth and incubate the tubes at 37°C for 1 hour with shaking at 225-250 rpm. 8. Plate ≤200 µl of the transformation mixture on LB agar plates containing the appropriate antibiotic (and containing IPTG and X-gal if color screening is desired). For the pUC18 control transformation, plate 5 µl of the transformation on LB-ampicillin agar plates. 9. Incubate the plates at 37°C overnight. If performing blue-white color screening, incubate the plates at 37°C for at least 17 hours to allow color development (color can be enhanced by subsequent incubation of the plates for 2 hours at 4°C). 10. For the pUC18 control, expect 25 colonies ($\geq 1 \times 10^8$ cfu/µg pUC18 DNA). For the experimental DNA, the number of colonies will vary according to the size and form of the transforming DNA, with larger and non-supercoiled DNA producing fewer colonies.
Blue-White Color Screening	Blue-white color screening for recombinant plasmids is available when transforming this host strain (containing the <i>lacI</i> ^q Δ <i>M15</i> gene on the F' episome) with a plasmid that provides α-complementation (e.g. the Stratagene pBluescript II vector). When lacZ expression is induced by IPTG in the presence of the chromogenic substrate X-gal, colonies containing plasmids with inserts will be white, while colonies containing plasmids without inserts will be blue. If an insert is suspected to be toxic, plate the cells on media without X-gal and IPTG. Color screening will be eliminated, but lower levels of the potentially toxic protein will be expressed in the absence of IPTG.

Critical Success Factors and Troubleshooting

Quantity and Volume of DNA: The greatest efficiency is obtained from the transformation of 1 µl of 0.01 ng/µl supercoiled pUC18 DNA per aliquot of cells. A greater number of colonies may be obtained by transforming up to 50 ng DNA, although the resulting efficiency (cfu/µg) may be lower. The volume of the DNA solution added to the reaction may be increased to up to 2.5 µl, but the transformation efficiency may be reduced.

Heat Pulse Duration and Temperature: Optimal transformation efficiency is observed when cells are heat-pulsed at 42°C for 60 seconds in the tubes provided. Efficiency decreases sharply when either the temperature or the duration of the heat pulse is altered.

Plating the Transformation Mixture: If plating <100 µl of cells, pipet the cells into a 200 µl pool of medium and then spread the mixture with a sterile spreader. If plating ≥100 µl, the cells can be spread on the plates directly. Tilt and tap the spreader to remove the last drop of cells. If desired, cells may be concentrated prior to plating by centrifugation at 1000 rpm for 10 minutes followed by resuspension in 200 µl of SOC medium.

Preparation of Media and Reagents

SOB Medium (per Liter)

20.0 g of tryptone
5.0 g of yeast extract
0.5 g of NaCl

Add dH₂O to a final volume of 1 liter and then autoclave

Add 10 ml of filter-sterilized 1 M MgCl₂ and 10 ml of filter-sterilized 1 M MgSO₄ prior to use

SOC Medium (per 100 ml)

Prepare immediately before use

2 ml of filter-sterilized 20% (w/v) glucose or 1 ml of filter-sterilized 2 M glucose

SOB medium (autoclaved) to a final volume of 100 ml

LB Agar (per Liter)

10 g of NaCl
10 g of tryptone
5 g of yeast extract
20 g of agar

Add deionized H₂O to a final volume of 1 liter

Adjust pH to 7.0 with 5 N NaOH and then autoclave

Pour into petri dishes (~25 ml/100-mm plate)

LB-Ampicillin Agar (per Liter)

1 liter of LB agar, autoclaved and cooled to 55°C

Add 10 ml of 10 mg/ml filter-sterilized ampicillin

Pour into petri dishes (~25 ml/100-mm plate)

Plates for Blue-White Color Screening

Prepare the LB agar and when adding the antibiotic, also add 5-bromo-4-chloro-3-indolyl-β-D-galactopyranoside (X-gal) to a final concentration of 80 µg/ml [prepared in dimethylformamide (DMF)] and isopropyl-1-thio-β-D-galactopyranoside (IPTG) to a final concentration of 20 mM (prepared in sterile water). Alternatively, 100 µl of 10 mM IPTG and 100 µl of 2% X-gal may be spread on solidified LB agar plates 30 minutes prior to plating the transformations. (For consistent color development across the plate, pipet the X-gal and the IPTG into a 100-µl pool of SOC medium and then spread the mixture across the plate. Do not mix the IPTG and the X-gal before pipetting them into the pool of SOC medium because these chemicals may precipitate.)

Limited Product Warranty

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

Endnotes

*U.S. Patent Nos. 6,706,525, 5,512,468 and 5,707,841 and equivalent foreign patents.

For in vitro use only. This certificate is a declaration of analysis at the time of manufacture.