




Use of Short, Small-Particle Columns for Fast and Efficient Separations


Application
Agrichemical
Robert Ricker

Herbicide Mixture

ZORBAX SB-C8
(4.6 x 75 mm) 3.5 μ m
(P/N: 866953-906)
Flow Cell 5mm, 2.5 μ l
20-60% B in 10 min
2 mL/min.



ZORBAX SB-C8
(4.6 x 150 mm) 5 μ m
(P/N: 883975-906)
Flow Cell 10mm, 8 μ l
20-60% B in 30 min.
1 mL/min.



1. Tebuthiuron
2. Prometon
3. Prometryne
4. Atrazine
5. Bentazon
6. Propanil
7. Propazine
8. Metolachlor

Conditions:
Column: ZORBAX SB-C8
Mobile Phase:
A H2O to pH 2.0 with TFA
B ACN

Highlights

- High-speed, high-resolution separation of a herbicide mixture in less than 8 minutes.
- Smaller particle size (3.5 μ m) leads to separations of similar efficiency when compared to 5 μ m particles in columns of twice the length.
- Increased flow rate on the 3.5 μ m, 75 mm column further decreases the run-time. Efficiencies are nearly equivalent to those of the 5 μ m 150 mm column run at a slower flow rate.



Agilent Technologies

*Robert Ricker is an application chemist
based at Agilent Technologies, Wilmington,
Delaware.*

For more information on our products and
services, visit our website at:
www.agilent.com/chem

Copyright© 2002 Agilent Technologies, Inc.
All Rights Reserved. Reproduction,
adaptation or translation without prior
written permission is prohibited, except as
allowed under the copyright laws.

Agilent shall not be liable for errors
contained herein or for incidental or
consequential damages in connection with
the furnishing, performance, or use of this
material.

Information, descriptions, and specifications
in this publication are subject to change
without notice.

Printed in the USA
April 25, 2002
5988-6288EN



Agilent Technologies