

# <sup>1</sup>H-<sup>19</sup>F Analyses on the Agilent 400-MR

Data Sheet

#### Introduction

Agilent's 400-MR, when equipped with an Automated Triple Broadband (ATB) or AutoSwitchable (ASW) probe, is ideal for analyses of compounds containing <sup>19</sup>F, such as drug precursors. The following experiments are available: 1-dimensional (1-D) fluorine decoupled proton spectra, 1-D proton decoupled fluorine spectra, and 2-D protonfluorine through-bond and through-space correlation.

#### **Key Benefits**

- **Simplified spectra**–Reduced complexity <sup>19</sup>F and <sup>1</sup>H NMR data analyses give you faster fluorine moiety confirmation.
- **Easy to configure**—These applications are included on a standard 400-MR when equipped with a standard ATB or ASW probe. No special hardware or cabling required.
- **Easy to use**–All calibrations and experimental parameters are included in the standard Agilent MR Workstation software available with the 400-MR.
- **Robust and reliable**—The ATB and ASW probes have dual tuned H/F coils and require minimal mechanical operation.





Agilent Technologies

#### Studying Compounds Containing <sup>19</sup>F

The use of <sup>19</sup>F in place of <sup>1</sup>H nuclei can radically alter the ADME properties of a drug and is widely used in the drug discovery process [1]. NMR is a valuable tool to ensure that a fluorine moiety remains intact and in the correct position through a synthetic procedure. The standard Agilent 400-MR is capable of analyzing compounds containing <sup>19</sup>F constituents, when equipped with an ATB or ASW probe.

### **Simplified Spectra**

In NMR data acquisition, <sup>1</sup>H and <sup>19</sup>F couple to each other, often creating complex splitting patterns. Using the 400-MR, spectra can be simplified by decoupling techniques resulting in clearer assignments, more visible impurites and more reliable confirmations.

### **Information Rich Data**

The 400-MR is also capable of performing experiments that define the position of a <sup>19</sup>F within a molecule by correlating through-space and through-bond H-F information.

#### **Applications Include**

- Drug discovery, organic synthesis confirmation
- · Drug development, purity analyses
- Metabolite identification and quantitation

The Agilent 400-MR can be used to provide useful structural information on compounds containing <sup>19</sup>F moieties. For the molecule shown, the vertical lines indicate the spectral elements attributed to each nucleus. The bottom spectrum in each group of four shows the coupled data; the top spectrum shows the simplification obtained through decoupled data.

#### Reference

1. For example, see Kitteringham et al, Annu. Rev. Pharmacol. Toxicol., 2001, 41, 443



<sup>1</sup>H observe

#### www.agilent.com

Product specifications and descriptions in this document are subject to change without notice.

© Agilent Technologies, Inc., 2011 Published in USA, February 1, 2011 Publication Number SI-0984



## **Agilent Technologies**