Notices

© Agilent Technologies, Inc. 2006

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G3100-90000

Edition

10/06

Printed in Germany

Agilent Technologies Hewlett-Packard-Strasse 8 76337 Waldbronn

Microsoft [®] is a U.S. registered trademark of Microsoft Corporation.

Software Revision

This guide is valid for A.01.xx revisions of the Agilent 3100 OFFGEL Fractionator software, where xx refers to minor revisions of the software that do not affect the technical accuracy of this guide.

Warranty

The material contained in this document is provided "as is," and is subiect to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

In This Guide...

This book describes how to operate the Agilent 3100 OFFGEL Fractionator including the control module.

Chapter Overview

1 Product overview

This chapter gives an overview of the Agilent 3100 OFFGEL Fractionator. It describes the intended use and the fractionation principles of the instrument, highlights its features and outlines the different instrument components.

2 Site requirements and specifications

This chapter describes information about site requirements that have to be met so that the Agilent 3100 OFFGEL Fractionator performs within specifications.

3 Installation

This chapter describes how to install the Agilent 3100 OFFGEL Fractionator.

4 **Operation**

This chapter details how to operate the Agilent 3100 OFFGEL Fractionator. It familiarizes the user with all control module functions and describes the necessary steps to perform a fractionation.

5 Maintenance

This chapter describes cleaning procedures that should be used to do maintenance on your instrument at regular intervals or before every use. In addition, it provides instructions on how update the firmware.

6 Troubleshooting

This chapter provides support for troubleshooting the instrument, the control module and the application.

7 Safety

This chapter contains safety information relevant to the Agilent 3100 OFFGEL Fractionator.

8 Appendix

This chapter contains details on non-condensing condition as well as ordering information for consumables and spare parts.

Related Documents

The Agilent 3100 OFFGEL Fractionator can be used with different consumable kits which each are supplied with a *Kit Manual*:

• 5969-1582 Kit Manual

1 Product overview 11

Intended Use 12

Fractionation principles of the Agilent 3100 OFFGEL Fractionator 12

Features 14

Instrument and Accessories 15

2 Site requirements and specifications 19

Site Requirements for the Agilent 3100 OFFGEL Fractionator 20 **Power Considerations** 20 **Bench Space** 20 Environment 20 The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) 21 Do not dispose of in domestic household waste 21 Sound Emission 21 22 Instruction for Use Physical Specifications of the Agilent 3100 OFFGEL Fractionator 23

3 Installation 25

Unpacking the Agilent 3100 OFFGEL Fractionator26Setting up the Agilent 3100 OFFGEL Fractionator28

4 Operation 31

Control Module Keys 32 The Esc key 33 The Enter key 33 Navigation keys 33

Action keys 33 Selection keys 33 **Direction keys** 33 Numeric/Alphanumeric keys 34 The Control Module User Interface 35 Start Screen 35 **Tray Status Screen** 36 The Plot Screen 38 The Method Screen 39 **Timetable Screen** 40 Save As Screen 40 Load Screen 41 Logbook Screen 41 Last Result Screen 42 **Configure Screens** 43 43 LAN Screen and Service **Date & Time Screen** 44 **Display Screen** 45 About Screen 45 47 **General Functions Setting Cooling Platform Temperature** 47 Setting Date and Time 47 Setting Display Contrast 48 Working with Methods 49 Loading a Method 49 Modifying a Method 49 Saving a Method 51 52 **Deleting a Method** Transferring a Method 52 **Timetable Programming** 53

55 **Displaying Data Graphically** Selecting Signals 55 **Rescaling the Plot Screen** 56 X (time) axis 56 Y (signal unit) axis 56 57 **Toggling Signals Displaying Exact Signal Values** 57 **Exporting Results** 59 How to find the exported files on the USB stick 60 How to display exported results 60 Performing a Fractionation 65 Preparing a tray for the OFFGEL mode 66 Preparing a tray for the standard IPG IEF mode 66 Starting a run 67 **Maintenance** 69 Cleaning the 3100 OFFGEL Fractionator 70 Cleaning the electrodes 70 70 **Regular cleaning procedure** Cleaning procedure for heavily contaminated electrodes 70 **Cleaning the Trays** 71 Cleaning the cooling platform 72 Cleaning the anode connector 72 Cleaning the top cover and the housing of the instrument 75 **Firmware Update** 76 Installation of the Firmware Update Tool 76 How to un-install the Firmware Update Tool 77 77 Installation of the Firmware on your PC Connect your PC to the Agilent 3100 OFFGEL Fractionator 77 Setup of the Firmware Update Tool 78

5

Configure the Firmware Update Tool80Doing the Firmware Update83The Log-file86

6 Troubleshooting 89

Troubleshooting the Instrument 90 **Status Indicators** 90 Status LED 90 Status bars 90 **Error Messages** 91 Hardware Diagnosis 93 Checking the LED and the cover sensor 93 Checking temperature control of the cooling platform 94 Checking high voltage, electrodes and electrode connectors 94 **Troubleshooting the Control Module** 97 Troubleshooting the Application 98 Looking at Run Parameters 101 Example 1. Typical results of the protein test sample 101 Example 2. Current drop during fractionation 103 104 Example 3. Reduced electrical charge in one lane

7 Safety 107

Safety Information 108 Power Cords 108 Operation 109 Safety Symbols 110 Chemical and Biological Safety 110

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Product overview

1

Intended Use 12 Features 14 Instrument and Accessories 15

This chapter gives an overview of the Agilent 3100 OFFGEL Fractionator. It describes the intended use and the fractionation principles of the instrument, highlights its features and outlines the different instrument components.

1 Product overview Intended Use

Intended Use

The Agilent 3100 OFFGEL Fractionator system fractionates proteins or peptides based on their isoelectric points (pI) in liquid phase (OFFGEL mode). In addition, the system performs standard in-gel isoelectric focusing (IEF) using immobilized pH gradient gel (IPG) strips (in-gel mode). The system consists of a standalone instrument and a handheld control module. The control module provides complete local control and monitoring of the Agilent 3100 OFFGEL Fractionator. The system is capable of running a total of 16 samples in parallel (2×8).

OFFGEL pre-fractionation of complex protein and peptide sample provides additional sensitivity of downstream analysis tools such as LC-MS. In addition, the pI is a useful parameter that can be used for validation of peptide MS hits.

Fractionation principles of the Agilent 3100 OFFGEL Fractionator

The Agilent 3100 OFFGEL Fractionator performs isoelectric focusing of proteins or peptides in immobilized pH gradient (IPG) gel strips. The 3100's OFFGEL electrophoresis differs from conventional gel electrophoresis in that the sample components do not remain in the gel. Instead, they are recovered from a buffer solution, making the recovery much easier than with conventional gels.

The three-part diagram in Figure 1 shows starting and ending states, as well as the OFFGEL fractionation process itself. After rehydration of gel strips, the immobilized pH-gradient (IPG) gel seals tightly against the compartment frame. Then, the diluted sample is distributed across all wells in the strip (Figure 1A). When a voltage is applied to the ends of the gel strip, the proteins or peptides move through the gel until the molecules reach a well covering a portion of the gel where the pH equals the pI of the molecule (Figure 1B). After fractionation, the liquid fractions containing pI-based separated proteins or peptides can easily be removed with a (multichannel-) pipette and processed for downstream experiments (Figure 1C).

Figure 1 Fractionation principle of the Agilent 3100 OFFGEL Fractionator.

1 **Product overview** Features

Features

- pI-based OFFGEL fractionation with in-solution recovery for easy transfer to downstream methods such as LC/MS.
- Reproducible pI-based OFFGEL fractionation down to 0.1 pH resolution depending on the consumable kit used.
- Preconfigured methods for OFFGEL and traditional in-gel IEF protein and peptide fractionation.
- Compatible with other upstream or downstream methods (immunodepletion, LC-MS, gel-based analysis)
- pI obtained acts as an additional validation parameter.
- No MS-incompatible additives used in fractionation buffer
- The instrument can also be operated in the traditional in-gel mode (IEF).
- Fractionation of up to 16 samples in parallel on 2 trays, accommodating 8 samples each.
- Broad loading capacity range from 50 μg sample load up to 5 mg sample load.
- Two independent power supplies allow two different methods to be run at the same time.
- Easy method setup and monitoring of important fractionation parameters and settings such as voltages and currents during a run.
- Online current control for each individual sample or gel-strip during fractionation.
- Viewing of up to three run parameters in an online plot.
- Easy export of important run parameters of a fractionation to a USB stick for further analysis
- Automated storage of up to 1200 data files, including important run parameters, on the compact flash card
- Monitoring of all operations and error events using the self-updating logbooks.

Instrument and Accessories

Figure 2 gives an overview of the complete Agilent 3100 OFFGEL Fractionator system:

The instrument has two individual power supplies which provide voltages to the two trays. On the left side of the cooling platform is the anode connector; across the cooling platform are two cathode connectors (see Figure 3).

1 Product overview

Instrument and Accessories

Figure 3 View inside the instrument.

Attached to each tray are a fixed electrode and a movable electrode that establish the connection to the anode and the cathode connectors (see Figure 4). Each tray accommodate up to eight samples. Each OFFGEL fractionation requires a frame to build up the well compartments as depicted schematically in Figure 1. There are two different frame lengths available, depending on the consumable kit and application chosen. The 12-well frame is typically used for fractionations requiring lower resolution, the 24-well frame is used for fractionations requiring higher resolution. Trays and electrodes can be ordered as spare parts if needed (for ordering details refer to Table 7 on page 150).

Figure 4 Electrodes, trays and frames

The frames are included in the consumable kits. The available consumable kits are listed in Table 1.

Product	Order Number
3100 OFFGEL, High Res Kit, pH 3-10	5188-6424
3100 OFFGEL, Low Res Kit, pH 3-10	5188-6425
3100 OFFGEL, High Res Kit, pH 4-7	5188-6426
3100 OFFGEL, Low Res Kit, pH 4-7	5188-6427
3100 OFFGEL, Starter Kit	5188-6444
OFFGEL, Protein Test Sample	5188-6428

1 **Product overview**

Instrument and Accessories

Agilent 3100 OFFGEL Fractionator User Manual

2 Site requirements and specifications

Power Considerations 20 Bench Space 20 Environment 20 The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) 21 Sound Emission 21 Instruction for Use 22 Physical Specifications of the Agilent 3100 OFFGEL Fractionator 23

This chapter describes information about site requirements that have to be met so that the Agilent 3100 OFFGEL Fractionator performs within specifications.

Site Requirements for the Agilent 3100 OFFGEL Fractionator

Power Considerations

The Agilent 3100 OFFGEL Fractionator power supply has wide ranging capabilities and accepts any line voltage in the range 100–240 V with a line frequency of 50–60 Hz (see also on page 23 for more details). Consequently, there is no voltage selector in the rear of the instrument.

CAUTION

2

Shock hazard or damage to your Agilent 3100 OFFGEL Fractionator may result if the Agilent 3100 OFFGEL Fractionator is connected to a line voltage higher than specified.

Bench Space

The Agilent 3100 OFFGEL Fractionator requires sufficient desk or laboratory bench space; dimensions and weights are shown in Table on page 23. Approximately 15 cm (6 inches) of clear space at the rear of the instrument for electric connections is needed.

Environment

The Agilent 3100 OFFGEL Fractionator will work within the specifications of ambient temperatures and relative humidity described in Table on page 23.

CAUTION

Do not store, ship or use the Agilent 3100 OFFGEL Fractionator under conditions where temperature fluctuations could cause condensation within the Agilent 3100 OFFGEL Fractionator. Condensation will damage the electronics. If the Agilent 3100 OFFGEL Fractionator was shipped in cold weather, leave it in its box and allow it to warm up slowly to room temperature to avoid condensation.

NOTE Put the Agilent 3100 OFFGEL Fractionator only on a flat and horizontal work top. Inclinations of more than 0.5% will cause the instrument to malfunction.

NOTE

The Agilent 3100 OFFGEL Fractionator is designed to operate in a controlled electromagnetic environment (EN61326/A1) where RF transmitters, such as mobile phones, should not be used in close proximity.

The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC)

The Waste Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC), adopted by EU Commission on 13 February 2003, is introduced producer responsibility for all Electric and Electronic appliances from 13 August 2005.

This product complies with the WEEE Directive (2002/96/EC) marking requirements. The affixed label indicates that you must not discard this electrical/electronic product in domestic household waste.

Product Category:

With reference to the equipment types in the WEEE Directive Annex I, this product is classed as a "Monitoring and Control instrumentation" product.

Do not dispose of in domestic household waste

To return unwanted products, contact your local Agilent office, or see www.agilent.com for more information.

Sound Emission

This statement is provided to comply with the requirements of the German Sound Emission Directive of January 18th 1991.

2 Site requirements and specifications

Site Requirements for the Agilent 3100 OFFGEL Fractionator

This product has a sound pressure emission (at the operator position) < 70 dB (A).

- Sound Pressure Lp < 70 dB (A)
- At Operator Position
- Normal Operation
- According to ISO 7779:1988/EN 27779/1991 (Type Test).

Instruction for Use

The Agilent 3100 OFFGEL Fractionator is designed and tested as general purpose electrical equipment for laboratory use. It is designed and tested for research use only and not for use in diagnostic procedures.

CAUTION

If the equipment is used in a manner not specified in this manual, the protection provided by the equipment may be impaired.

Physical Specifications of the Agilent 3100 OFFGEL Fractionator

Туре	Specification	Comment	
Weight	13 kg (29 lbs)		
Dimensions	157 x 355 x 427 mm	Height × width × depth	
Line voltage	100–240 VAC	Wide-ranging capability	
Line frequency	50–60 Hz		
Power consumption	140 VA		
Ambient operating temperature	5–40 °C (41–104 °F)	see also humidity specification	
Ambient non-operating temperature	-40-70 °C (-40-158 °F)	Humidity 0 - 90%	
Humidity	15 - 80 %, at 5–40 °C (41–104 °F) (non-condensing; e.g. max humidity is 55 % if the cooling platform is used 10 °C below ambient temperature)	see details in Appendix 8, "Non-condensing conditions"	
Cooling platform temperature	10 - 35 °C (50-95 °F) maximal 10 °C below ambient temperature		
High voltage power supply	500 V - 10 kV, maximal 150 μA per gel strip, maximal 1 W per gel strip		
Operating altitude	2000 m (6666 ft)		
Non-operating altitude	4500 m (15000 ft)		
Safety Standards: IEC, CSA, UL	Installation Category II Safety Class I, Pollution Degree 2		

Table 2Physical Specifications

2 Site requirements and specifications

Site Requirements for the Agilent 3100 OFFGEL Fractionator

Agilent 3100 OFFGEL Fractionator User Manual

Installation

3

Unpacking the Agilent 3100 OFFGEL Fractionator 26 Setting up the Agilent 3100 OFFGEL Fractionator 28

This chapter describes how to install the Agilent 3100 OFFGEL Fractionator.

3 Installation

Unpacking the Agilent 3100 OFFGEL Fractionator

Unpacking the Agilent 3100 OFFGEL Fractionator

Before you begin the installation, ensure that the place you have chosen meets the site requirements as specified in "Site Requirements for the Agilent 3100 OFFGEL Fractionator" on page 20. The shipment consists of three boxes, one containing the handheld controller, one containing the instrument and one containing the power cable.

- **1** Unpack the handheld controller. The box includes the control module and a *Declaration of Conformity*.
- **2** Open the box in which the Agilent 3100 OFFGEL Fractionator was shipped. On the top you will find a box containing the accessories.

Figure 5 Opened instrument shipment box.

3 Take out the accessory box and lift the instrument with its foam packing out of the outer box.

Installation 3

Unpacking the Agilent 3100 OFFGEL Fractionator

Figure 6 Instrument inner shipment box.

- **4** Remove the foam packing and the inner box and place the instrument in a suitable place.
- **5** Check that you have the complete set of accessories available using Table 3.

De	escription	Part Number	
1	User Manual	G3100-90000	
2	Hex key set 1 – 5 mm	8710-0641	
3	Compact Flash Card 512 MB	01100-13001	
4	RS232 cable	RS232-61601	
5	USB-Memory Stick Kit	G4208-68700	

Table 3Accessories

Setting up the Agilent 3100 OFFGEL Fractionator

Setting up the Agilent 3100 OFFGEL Fractionator

After you have unpacked the Agilent 3100 OFFGEL Fractionator and positioned it in a suitable place on your bench, you need to install the flash card and connect the power cord and the control module.

The Agilent 3100 OFFGEL Fractionator has the following connectors at the rear:

Figure 7 Connections at the rear of the instrument.

- **1** Insert the compact flash card into the dedicated position.
- **2** Connect the control module to one of the CAN connectors.

Be sure to not to confuse the CAN connector with the LAN connector. The control module will not work when connecting it to the LAN connector.

- **3** Connect the power cable to the instrument and to the line power
- **4** Switch on the instrument.

NOTE

3

On the control module, you should see the following welcome screen:

Figure 8 Control module welcome screen.

After the welcome screen, the control module displays the following start screen:

<u>Tray I</u> Idle		Tray II Idle					
Method	def <i>a</i> ult	Method	def <i>a</i> ult				
Volt-Hours	0.000 kVh	Volt-Hours	0.000 kVh				
Voltage 0 V		Voltage	0 V				
Temperature 20.00 °C							
Ready							
Tray I 📲	Config 👔 Met	hods 👔 Abou	t 📔 🛛 Tray II				

Figure 9 Control module start screen.

NOTE

If you have problems with the control module please refer to "Troubleshooting the Control Module" on page 97 in this manual.

The instrument LED is usually yellow, which indicates that the cooling platform is cooling down to the specified temperature (typically 20°C). Once the temperature has been reached, the LED is switched off.

NOTE

To test the instrument for proper functionality, we recommend that you perform the tests as specified in chapter "Troubleshooting the Instrument" on page 90.

3 Installation

Setting up the Agilent 3100 OFFGEL Fractionator

Agilent 3100 OFFGEL Fractionator User Manual

Operation

Control Module Keys 32 The Control Module User Interface 35 General Functions 47 Working with Methods 49 Displaying Data Graphically 55 Exporting Results 59 Performing a Fractionation 65

This chapter details how to operate the Agilent 3100 OFFGEL Fractionator. It familiarizes the user with all control module functions and describes the necessary steps to perform a fractionation.

Control Module Keys

The display will show you a variety of menu buttons [**F1-F5**] (in the lower section) or function buttons [**F6-F8**] (in the right hand section) that can be accessed with the corresponding Navigation (for menus) and Action (for functions) keys.

accepts current input entry or action

Figure 10 3100 Control module keys.

NOTE

The i-key and the m-key have no function in the 3100 control module.

The Esc key

Esc (Escape) allows you to exit the current window or screen and leads you back to the last window or screen you were working with.

In an edit field, the previous value can be restored by pressing ESC.

The Enter key

With **Enter**, you accept a current entry or action. When entering a parameter into a certain field, **Enter** leads you on to the next accessible entry field. In this case, it has the same function as the right Direction key.

Navigation keys

These five keys [F1 - F5] allow you to switch between the menus. Within these menus, the relevant parameters can be set, and certain functions can be accessed. The Navigation keys always correspond to a button displayed above them on the screen. The menus accessed via the buttons vary according to the screen you are working with.

Action keys

The three Action keys [**F6-F8**] trigger a variety of functions. The available functions depend on the screen you are working with.

Selection keys

The Selection keys (**arrow up/down**) allow you to select settings in various list boxes. You can also change values in certain parameter entry fields, or enter alphanumeric characters.

Direction keys

The Direction keys (**arrow left/right**) allow you to move back and forth between the entry fields.

Numeric/Alphanumeric keys

These keys allow you to enter numeric values in parameter entry fields. In certain fields where alphabetical characters may be entered, you can use the **Numeric/Alphanumeric** keys to do so. Pressing them several times in sequence changes the current value according to Figure 11.

Figure 11 Keypad of the Control Module

4

The Control Module User Interface

The **Start Screen** of the control module represents the main screen of the application. All other screens can be accessed from the start screen by using the navigation keys [**F1 - F5**]. An outline of the application structure is depicted in Figure 12.

Figure 12 Application Structure: Screens of the Control Module.

Start Screen

The Start screen (Figure 13) gives an overview about the actual state of the instrument.

4 Operation

The Control Module User Interface

Tray status							
	Tray I Idle			Tray II Idle			
Actual cooling platform temperature	Method Volt-Hours Voltage	default 0.000 kVh 0 V		Met Volt-He Volt	thod ours tage	def <i>a</i> ult 0.000 0	kVh V
p		Temperature	20.00 Ready	0 °C /			
Instrument status	Tray I	Config 👔	Method	s	About		Tray II

A tray status is displayed for each tray individually. Since space on the display is limited, you see only the most important parameters on the start screen; details for each tray can be monitored on the tray status screens Tray 1 [**F1**] and Tray 2 [**F5**]. The following four parameters are shown:

- status **Idle** if no fractionation is running, status **Focusing** if a fractionation is running, status **"cover open"** when the top cover is opened.
- Method: the name of the method that was loaded.
- Volt-hours: the volt hours that the fractionation is running.
- Voltage: the actual voltage that is currently applied to the tray.

Below the tray status are two status bars; one displays the actual temperature of the cooling platform and the other the status of the instrument itself.

From the status screen, you can move to all other screens using the navigation keys [**F1-F5**] (see Figure 12).

Tray Status Screen

The Tray status screen (Figure 14) displays all actual and set run parameters for each tray individually during focusing (Tray 1 [F1] and Tray 2 [F5]). Current is even measured for each IPG-strip individually.
Operation 4 The Control Module User Interface

Figure 14 Tray status screen information.

The status bar on the top of the screen displays the screen name, the elapsed volt hours (the integral of voltage over time), the name of the loaded method, the instrument status and the time of day. This status bar is also visible on subsequent screens when navigating through them. An asterisk (*) occurs next to the method name denotes that the method was modified during the run. The icon in the top right corner indicates whether an USB stick is connected to the instrument or not. Black denotes no USB stick inserted and white denotes a USB stick inserted.

The main panel of the status screen contains information about the current run parameters. The current method step (e.g. **Focusing** or **Hold**) is displayed on the top left corner of the main panel; next to it is displayed the elapsed time since run start in hours and minutes (h:m).

The run parameter table contains the specified limits for volt hours, voltage, power and current of the loaded method and the actual measured values. These parameters apply to each strip individually.

The bar graph on the right side visualizes the actual currents for each strip online; in addition, the numeric value is displayed within each bar.

The screen provides access to:

- Plot [F1] displays important parameters graphically,
- Method [F2] gives access to method management functions (load, save, delete),

The Control Module User Interface

- Logbook [F3] shows a log of the events for this specific tray,
- Other Tray [F4] switches to tray status screen of the other tray,
- Results [F5] shows a summary of the last performed fractionation.
- Start/Stop [F8] starts or stops a fractionation.

The Plot Screen

The **Plot screen** (Figure 15) can display online important parameters graphically. It is possible to display up to three signals at the same time.

You can view the following signals:

- current for each strip individually,
- voltage,
- volt hours,
- limits of current and voltage.

You can use a cursor to pinpoint information from the x-axes directly from the Plot screen; the value is then displayed in the upper right corner. You can use the Direction and Selection keys of the control module to change the scale of the plot window. You can use the Rescale button [**F7**] you to optimize the screen according to the selected signal and parameter ranges. You can modify the range of the signals using the Select button [**F6**].

4

The Method Screen

The **Method screen** (Figure 16) is available either from the start screen or from the tray status screen (see Figure 12 on page 35). Navigating to the method screen via the tray status screen allows you to work in the online mode, whereas navigating directly from the start screen allows you to work in the offline mode. The main difference between online and offline mode is that changes that are made in the method and timetable online are directly applied to a currently running fractionation without storing the changes in the method file. If you want to permanently change a method, you can either work in the offline mode or save the method from the online mode under a different file name. Detailed instructions on how to work with methods can be found under "Working with Methods" on page 49.

Method	0.4	98 kVh	default	Ready	14	:00:05 📖
	Volt.Hour	Voltage	Current	Power [mW]	Time [h:m]]
Focusing	100.0	8000	50	200	1:00	X
Hold		8000	50	200		Cancer
Timetable <empty></empty>						
Timetable				Load	1 8	ave

You use **Cancel** [**F7**] to exit the screen, or **Done** [**F6**] either to apply the changes to the running fractionation (online mode) or to save the changes to the file (offline mode).

The screen provides access to:

- Timetable [F1] editing, inserting, and deleting of focusing steps,
- Save [F4] saving a method on the internal disk or on the USB stick,
- Load [F5] loading a method for execution or editing,

See "Timetable Screen", "Save As Screen" and "Load Screen" for further details.

The Control Module User Interface

Timetable Screen

You use the **Timetable screen** (Figure 17) to implement time-based steps in the focusing phase of your method. The values from the focusing step (method screen) become part of the timetable and are displayed as startup values at the bottom of the timetable screen.

Figure 17 Timetable screen

You use **Insert** [**F7**] to call the timetable edit dialog and insert a new row in the time table. You can use **Edit** [**F8**] to change existing timetable entries, or **Delete** [**F6**] to delete existing timetable entries. You use **Exit** [**F5**] to confirm the changes and exit the screen.

Save As Screen

The **Save As** screen (Figure 18) allows you to create a name for your method and save it [**F6**] either to the internal disc or to an USB stick. You can save the method using the same name when working in the offline mode but you have to use a different name when working in the online mode.

The Control Module User Interface

Save As	3	default	Ready	17:54:00 📖
Drive Internal Disk Internal Disk Internal Disk Internal Disk	Method Name 111 123 222 333	Ŷ Nam Driv	e <mark>default</mark> e Internal Disk	
Internal Disk Internal Disk	50 51			Save

Figure 18 Save as screen

Load Screen

You use the **Load screen** (Figure 19) to load, delete or transfer a method which was saved on the internal disk or a USB stick.

Load		default	Ready	17:53:13 📖
Drive	Method Name	Date	Time	
Internal Disk	333		±	Delete
Internal Disk	50			
Internal Disk	51			
Internal Disk	default			Taister
USB Stick	0			
USB Stick	1		+	
	ak.	Ani		

Figure 19 Load screen

You use **Load** [**F6**] to call a method for editing (offline mode) or execution (online mode). You use **Transfer** [**F7**] to transfer a method from the internal disc to a USB stick or vice versa, or you can use **Delete** [**F8**] to delete existing methods.

Logbook Screen

The **Logbook screen** (Figure 20) shows a log of all events for the selected tray, including the events that have no specific tray context. The table displays the message text, the event Id, date and time.

The Control Module User Interface

Tray I	0.000 kVh	OG12PR00	Ready	10:1	8:31 🗆
Message			ld Date	Time	
Focusing end			5TRTE 09/26	10:18:11	
Profile start			5TRTE 09/26	10:17:06	
Focusing start			5TRTE 09/26	10:17:06	
Parameter changed			NFD 09/26	10:16:58	
Ready			5TRTE 09/26	10:16:46	4
*** Control module n	eady ***		INFD 09/26	10:18:18	
USB not mounted	-		5TRTE 09/26	10:16:37	Done
Ack. Error	<u> </u>				

You use **Ack.Error** [**F1**] to acknowledge permanent errors, such as a leak, after corrective actions have been performed (see Table on page 92).

The logbook can store up to 2000 events. Events are deleted when the instrument is power cycled.

Last Result Screen

The **Last Result** screen (Figure 21) shows a summary of the last performed run. If a fractionation is running on the selected tray, the screen shows the current status of that fractionation. In contrast to the tray status screen, the bar graph shows ampere hour for each strip separately and thus allows quick review of the focusing performed.

Last Result I 1.012 kVh	default	* Ready	10:2	21:41 🗆
ti	Fectric Char	ge		
Started 07/28/06 09:58		A-19.0 R-0.0		
Ended 07/28/06 10:21		C • 18.9		
Method default		E - 18.8		Export
Volt Hours 1.00 kVb		G-18.9		
	0.0	H • 0.0 [mAh]	19.0	✓ Exit
		1		

Figure 21 Last result screen

You use **Export** [**F7**] to export all acquired run parameters to an USB stick, and **Exit** [**F6**] to exit the screen.

4

Configure Screens

On the **Configure main screen** (Figure 22) you can define the temperature of the cooling platform. The temperature applies to both trays. The default temperature of the cooling platform is 20°C. The status bar shows the actual cooling platform temperature, the instrument status and the time.

Figure 22 Configure main screen

The screen provides access to:

- LAN [F1] information about the LAN settings,
- Date & Time [F2] configuring date and time,
- **Display** [F5] adjusting display contrast of the control module.

You use **Done** [F6] to confirm the changes and exit the screen, and **Cancel** [F7] to exit the screen without changing the temperature.

LAN Screen and Service

The LAN screen (Figure 23) and the subsequent service screen provide the possibility to connect the instrument to the LAN. This feature is currently not usable but is reserved for future applications.

The Control Module User Interface

Figure 23 LAN and service screen

Date & Time Screen

The **Date and Time screen** (Figure 24) allows you to adjust date and time to your local environment. This also corrects the time of day on the status bar as well as the date used for the logbook entries.

Figure 24 Date & Time screen.

You confirm the changes and exit the screen with Done [F6].

4

NOTE

It is important to set the correct time and date because the results files are stored by time when they were created and not by name.

Display Screen

You use the **Display screen** (Figure 25) to adjust the contrast of the control module screen.

You use **Done** [**F6**] to confirm the changes and exit the screen, and **Default** [**F7**] to set the display contrast back to the default value of 11.

About Screen

The **About screen** (Figure 26) shows information about the currently loaded firmware and controller module software, as well as the serial numbers of the instrument and controller. G3100A specifies the instrument and firmware; in the example in Figure 26 the serial number is DE00000037 with firmware revision U.00.12 build 0002. G3101A specifies the control module; in the example, the serial number is CN23605354 with software revision T.01.05. The build is reflected by the date and time.

The Control Module User Interface

[Agilent	Technologie © 2006 Agilent Tec	es OFFGEL chnologies. All righ	Fractionator ts reserved.
	G3100A	Serial Number	DE00000037	
		Revision	U.00.12	
		Build	[0002]	
000000	G3101A	Serial Number	CN23605354	
		Revision	T.01.05	
		Build	07/13/2006*1	5:00*2MB
L	L		1	

Figure 26 About screen

General Functions

Setting Cooling Platform Temperature

- 1 From the Start screen, press **Configure** [F2].
- 2 Enter the temperature for the cooling platform with the Numeric Keys.
- **3** Press **Done** [**F6**] to apply the changes.

NOTE

The temperature applies to both trays.

Setting Date and Time

- 1 From the Start screen, press **Configure** [F2].
- 2 In the Configure screen select Date & Time [F2].
- **3** Use the **Selection Keys** to specify a month; press **Enter** to confirm the month.
- 4 Use the Numeric Keys to specify the day; press Enter to confirm the day.
- **5** Use the **Numeric Keys** to specify the year; press **Enter** to confirm the year.
- **6** Proceed the same way with the time.
- 7 Press **Done** [F6] to confirm your changes.

When date and time are changed in the control module, the settings are automatically stored in the instrument. The date and time are backed up by a battery in the instrument. The control module does not have a clock of its own.

NOTE

Upon startup, the control module synchronizes the clock with the instrument.

Setting Display Contrast

- 1 From the start screen, press **Configure** [F2].
- 2 In the Configure screen, select **Display** [F3].
- **3** Use the **Selection Keys** to enter values between 0 and 31 to adjust the display contrast according to your personal and location requirements.
- 4 Press Done [F6] to confirm the settings or press Default [F7] to restore the default settings.

Working with Methods

A method consists of a focusing step and a hold step which contain all run parameters except temperature for a fractionation. The focusing step can include a timetable if desired. You can work with methods either in the offline or in the online mode. You work in the offline mode if you navigate to the method screen directly from the start screen. You work in the online mode if you navigate to the method screen via the tray status screens (see Figure 12). The following section describes the offline mode and makes separate notes for differences to the online mode.

NOTE

You will find several default methods on the control module which are ready to use. These can be modified to your needs if desired.

Loading a Method

- 1 From the Start screen press Method [F3].
- 2 Press Load [F4]; the Load screen is displayed.
- **3** Use **Selection Keys** to select the method from the list
- **4** Press **Load** [**F6**] to load the method; the run parameters of the loaded method are shown, and the method name is displayed on the status bar of the method screen.

NOTE

If you press the **Done** button [**F6**] in the online mode, the method gets loaded to the selected tray for execution, in the offline mode you exit the method screen.

Modifying a Method

A Method can be modified by changing its run parameters on the Method screen.

Working with Methods

- 1 Load the method you want to modify (see "Loading a Method").
- **2** Use the **Direction Keys** to highlight the parameter you want to change.
- **3** Use the **Numeric Keys** to enter the new values.
- **4** Proceed with step 2 and step 3 until you have modified all parameters.

NOTE

NOTE

Voltage, current and power are mandatory values required for each method; they define the upper limits during a run. Whatever parameter is reached first limits the others. Voltage hour and time define the upper limit of the duration of the focusing step; the run will stop once one of these values has been reached.

5 Press **Done** [**F6**] to save the method changes under the current method name.

If the method you are modifying is open in the online mode at the same time, a pop-up screen will asks you:if you want to continue (Figure 27).

Figure 27 Pop-up window when changing a method that is loaded on tray I.

Only if you confirm with **OK** by pressing **Enter** will the method changes be saved under the current method name. Otherwise, you need to save the modified method under a different name (step 6).

It is not possible to change a method offline while it is executing.

6 Press Save As [F5] to save the method under a different name (see "Saving a Method".

NOTE If you change a method online and press **Done** [**F6**], an asterisk (*) appears next to the method name in the status bar to indicate that the current method has been modified. This can be done either before or even during a run. The changes apply to the fractionation but

The time-programmable settings can be modified in the Timetable screen (see "Timetable Programming".

Saving a Method

- 1 From the Start screen, press Method [F3].
- 2 Select the Save As button [F5].

are not saved in the method.

- **3** If you want to save the method under the same name, continue with step 8. If you want to specify a new method name, press the **Selection Key up** to enter the letter input mode.
- **4** Select the required letter or number in the method name entry box using the **Selection Keys**. Numbers and the minus sign (-) can be entered directly.
- **NOTE** Alphanumeric characters can be entered also by pressing the keys 1 to 9 several times (e.g. 1 A B C, 2 D E F, 3 G H I,). See also Figure 10 on page 32.

Default methods on the internal disk are part of the system and cannot be overwritten.

- **5** Move the cursor to the next entry position using the **right Direction Key** and repeat step 5.
- 6 Repeat step 5 and step 6 until you have specified your method name.

NOTE The left direction key can be used as backspace, the right direction key is used to move to the next character position

- **7** Press **Enter** to exit the letter input mode.
- **8** Use the **right Direction Key** to move to the drive field.

Working with Methods

- **9** Press the **Selection Key down** and select the drive (Internal Disk or USB Stick).
- **10** Confirm your selection by pressing **Enter**.
- 11 Select Save [F6] to accept and save the changes.

Deleting a Method

- 1 From the Start screen, press Method [F3].
- 2 Press Load [F4]; the load screen is displayed.
- **3** Use the **Selection Keys** to select the method from the list
- 4 Select **Delete** [**F8**] to delete the method. You are asked if you really would like to delete the method.
- **5** Select **Yes** to confirm the deletion by pressing the **Enter** key.

NOTE

A method that is currently loaded to one of the trays or is being executed cannot be deleted from the method list.

Default methods on the internal disk are part of the system and cannot be deleted.

Transferring a Method

You can transfer methods from the internal disk of the instrument to the USB stick and vice versa.

- 1 From the Start screen, press Method [F3].
- 2 Press Load [F4]; the load screen is displayed.
- **3** Use the **Selection Keys** to select the method from the list
- **4** Select **Transfer** [**F7**] and a copy of the selected method is transferred to the other drive.

NOTE Default methods on the internal disk are part of the system and cannot be transferred to a USB stick.

Timetable Programming

The focus step of a method can be sub-divided into different time-based steps if desired. You can create a time-based program using the timetable screen.

- 1 From the Start screen, press Method [F3].
- 2 Press Timetable [F1]; the timetable screen is displayed.

NOTE

The values from the focusing step (Method screen) become part of the timetable and are displayed as startup values at the bottom of the Timetable screen.

- **3** Select **Insert** [**F7**] to insert a new row to the timetable.
- **4** Use the **Numeric Keys** to enter a time (hours and minutes) that specifies the time-span between startup value and the desired parameter change. Move from hours to minutes by using the **right Direction Key**
- **5** Use the **right Direction Key** to move to the Type field.
- **6** Press the **Selection Key down** and select the parameter type (Voltage, Power or Current).
- 7 Confirm your selection by pressing **Enter**.
- 8 Use the Numeric Keys to enter a value for the selected parameter type. This value specifies the upper limit of the selected parameter.
- 9 Confirm your selection by pressing Enter.
- **10** Use the **right Direction Key** to move to the Mode field.
- **11** Press the **Selection Key down** and select the Mode (Linear or Step)

NOTE

If you select the **Step** mode, the instrument approaches the set value from the current value stepwise. If you select the **Linear** mode the instrument approaches the set value from the current value in a linear manner, e.g. a gradual increase of voltage.

- **12** Confirm your selection by pressing **Enter**.
- **13** Press **Done** [**F6**] to confirm all entries; this brings you back to the Timetable screen.

Working with Methods

The timetable becomes part of the current method only when the method is saved.

You can edit an existing timetable line by pressing Edit [F8]. Use Delete [F6] to delete the selected line from the timetable.

Displaying Data Graphically

In the Plot screen, you have many opportunities to display up to three signals on a graphic display while the analysis is performed.

Additionally, all run parameters of a fractionation are stored in a .csv file on the compact flash card and can be visualized in Microsoft Excel after exporting (see "How to display exported results" on page 60

Selecting Signals

- 1 From the Start screen, press Tray Status [F1 or F5].
- 2 Select Plot [F1]; the plot screen is displayed.
- **3** From the Plot screen press the **Select** button [**F6**] to show the **Plot selection** menu.

- **4** Use the **Direction** and **Selection** keys to navigate within and between the Available Signals and Selected Signals list boxes.
- 5 Exchange signals between the list boxes by pressing the **Move** button [F8] or the **Enter** key.

NOTE

At the right of the Selected Signals list box you is the legend to the signals.

Displaying Data Graphically

- **6** You can also change the time range (X axis) for the plot in this screen by using the **Direction** and **Numeric** keys.
- **7** The different signals can be set up by pressing the **Setup** button [**F7**]. Depending on which signal is highlighted, you can enter an individual Y-Range setting here.

Figure 29 Signal Setup menu.

8 When the signals and their X (time) and Y (signal unit) ranges have been specified, press the **Done** button [F6] to switch to the graphic view.

Rescaling the Plot Screen

X (time) axis

There are several possibilities to rescale the X (time) axis:

- Enter a time in the Plot Selection windows (available from the Plot screen via the **Select** button [**F6**])
- Perform a rescale directly in the Plot screen by pressing the Direction left/right keys. Each press of the right key halves the time range, each press of the left key doubles it. Press the keys several times to set up the appropriate time frame.

The time range is indicated at the bottom of the **Plot** screen. This setting is independent of the active signal.

Y (signal unit) axis

There are several possibilities to rescale the Y (signal unit) axis:

4

- From the Plot Selection windows (available from the Plot screen via the **Select** button [**F6**]) choose a signal from either list box and press the **Setup** button [**F7**]. You can specify a Y range separately for each signal. This setting can also be made for signals that are no part of the Selected Signals list box. Rescaling directly from the Plot screen overwrites these settings.
- Use the **Rescale** button [**F7**] in the **Plot** screen to adjust the Y axis according to the minimum and the maximum signal value within the set time range. Using this function provides the optimum signal display. It refers only to the active signal indicated at the top of the screen.
- Use the selection keys to change the scaling of the Y axis. Each press of the right key halves the time range, each press of the left key doubles it.

Toggling Signals

You can monitor up to three different signals in real-time (updates every second) on the **Plot** screen. Although all three signals are shown in the display, the **Rescale** button [**F7**] only refers to the active signal. The active signal is shown in the **Plot** screen title and can be toggled by pressing the 1, 2 and 3 keys on the numeric keypad.

Displaying Exact Signal Values

Selecting the **Cursor** button [**F8**] in the **Plot** screen displays the X and Y value of the current cursor position for the active signal. In this mode, you can rescale the Y axis using the Selection keys. You can use the Direction left/right keys to move the cursor along the graph in an X direction in order to find the position you want.

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Figure 30 Signal plot with cursor.

NOTE

Pressing the keys 1, 2, 3 switches the displayed parameter in the top right corner of the graph and provides a way to figure out which run parameter a graph is showing.

Exporting Results

You can export all run parameters of a fractionation as .csv file to a USB stick for further analysis.

NOTE

Since USB memory sticks may vary from vendor to vendor or from type to type, incompatibilities can occur. In general, USB-memory sticks from Sandisk and Kingston should work. The USB-memory sticks must be FAT-16 formatted and without encryption, maximum size is 2 GB. See "USB-Memory Stick Kit" on page 150.

A compatible USB stick is included with the shipment.

- 1 From the Start screen, press Tray Status [F1 or F5].
- 2 Select **Result** [**F5**]; the result screen is displayed.
- **3** Select the **Export** button [**F7**] to navigate to the Export menu.

Figure 31 Export menu.

- **4** Use the right **Direction** keys and the **Selection** keys to specify the result set you would like to export. You can either export the data for the complete year or specify a single fractionation that you would like to export.
- 5 Select the Export button [F7] to export the results to a USB stick.

NOTE

The compact flash card can store up to about 3 years' worth of results. The storage capacity is dependent on the file size and is thus method dependent. When the compact flash card is full, the oldest data file is overwritten automatically.

How to find the exported files on the USB stick

When exporting results or transferring methods to an USB stick, the following directories are generated:

Figure 32 Folder directories on the USB stick

The METH folder contains all exported methods. The RES folder contains all exported results. The subdirectory of the RES folder starts with an instrument folder in the name of the instrument serial number. The following directories follow the structure of the export menu: year, month, day, hour, minute (see Figure 31). The .csv file in the minute folder finally contains all run parameters of a fractionation.

How to display exported results

The results of all run parameters can be exported as .csv file and thus can be visualized with Microsoft Excel. Agilent Technologies has designed the *Agilent 3100 OFFGEL fractionator plotter* which is a Excel add-on and facilitates the analysis of the run parameters. The tool is provided via the Agilent

Technologies web page http://www.agilent.com (Life Sciences and Chemical Analysis Solutions - Technical Support - Downloads & Utilities) or direct via http://www.chem.agilent.com/scripts/DownloadsUtilities.asp and can be downloaded free of charge.

NOTE

Note that the tool was designed for Microsoft Excel 2003 and US-American system locale settings only and has not been tested with other versions of Microsoft Excel. However, it is expected that the tool also runs under Excel 2003 international versions and international system locale settings.

- 1 In Microsoft Excel 2003, open the downloaded file **Agilent 3100 OFFGEL Fractionator.xla**. You may also want to copy the file in the XLSTART folder (for office 2003 it can be C:\Program Files\Microsoft Office\OFFICE11\ XLSTART) so that the program is loaded automatically when Excel is started.
- 2 Click **Enable Macros** when the security warning window pops up.

NOTE If you cannot run the program, check that Excel security settings under tools>macro>security allow macros to be run.

The following button is generated in your tool bar:

Figure 33 Agilent 3100 OFFGEL Fractionator toolbar button.

3 Click on the Agilent 3100 OFFGEL Fractionator toolbar button to activate the following user interface:

NOTE

Exporting Results

Agilent	3100 OFFGEL	ractionator	
	Open File Ma	ke Graphic	_
~	v		
Γ	\checkmark		
	~	Γ	
	V		
	~		
v	V	Г	
Unch	eck All	Γ	

Figure 34 Agilent 3100 OFFGEL Fractionator user interface.

You can close the Agilent 3100 OFFGEL Fractionator user interface by clicking the close window button. You can reopen it by clicking the Agilent 3100 OFFGEL Fractionator toolbar button.

4	Click Open File,	select an exporte	d file with the	browser and	open the file.
---	------------------	-------------------	-----------------	-------------	----------------

NOTE The file must be an unchanged Agilent 3100 OFFGEL Fractionator .csv file.

- **5** The file name of the loaded file appears in the Agilent 3100 OFFGEL Fractionator list box (Figure 35)
- **NOTE** Multiple files can be loaded simultaneously. You can toggle between them by clicking on the file name.

Open File	Make	Graphic
O TII 40.	csv	
1		
✓ kVHours	CurrA	🗖 uAhA
VoltTT	CurrB	🗐 uAhB
🔽 Volt	CurrC	🗐 uAhC
PowerTT	CurrD	🗐 uAhD
Power	CurrE	🗐 uAhE
Temp	CurrF	IT uAhF
Top Cover	CurrG	🖵 uAhG
		1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.

- **6** Select all run parameters which should be plotted.
- 7 Click Make Graphic. The plot is displayed on a separate worksheet

Figure 36 Plotted graph of run parameters (y-axes) versus time [h].

NOTE

If a worksheet is manually deleted or the data file is reloaded, the corresponding graph is deleted. If a plot to a certain worksheet is repeatedly generated, the previous plot version is deleted.

Table gives details of the different run parameters:

Run Parameter	Description
VoltTT, PowerTT	Fixed settings of Volt and Power of the method that was used.
kVHours, Volt, Power	kVHours, Volt, Power measured during the fractionation for the complete tray.
Temp	Temperature of the cooling platform during the fractionation.
Top Cover	The events or time intervals when the top cover was opened during a fractionation are marked with the red diamond symbol on the X-axis. If there are no Top Cover events registered, the Top Cover check box is automatically unchecked.
CurrA, CurrB, CurrC,	Currents measured for individual tray lanes (A-H).
µAhA, µAhB, µAhC,	Measured µA-hours acquired for individual tray lanes (A-H).

Table 4Run parameters

Performing a Fractionation

Fractionations using the Agilent 3100 OFFGEL Fractionator can be performed in the OFFGEL mode or in the standard IPG IEF mode (in-gel mode).

The table lists the default methods delivered with the system and ready to use with a cooling platform temperature of $20\,^{\circ}\mathrm{C}$

Method name	recommended use
OG24PR00	OFFGEL default method for proteins for use with 3100 OFFGEL High Res Kit, pH 3-10, 24 well frames OR for use with 3100 OFFGEL High Res Kit, pH 4-7, 24 well frame
OG24PE00	OFFGEL default method for peptides for use with 3100 OFFGEL High Res Kit, pH 3-10, 24 well frames OR for use with 3100 OFFGEL High Res Kit, pH 4-7, 24 well frame
OG12PR00	OFFGEL default method for proteins for use with 3100 OFFGEL Low Res Kit, pH 3-10, 12 well frames OR for use with 3100 OFFGEL Low Res Kit, pH 4-7, 12 well frame
OG12PE00	OFFGEL default method for peptides for use with 3100 OFFGEL Low Res Kit, pH 3-10, 12 well frames OR for use with 3100 OFFGEL Low Res Kit, pH 4-7, 12 well frame
IG_24_00	IPG-IEF in-gel default method for protein and peptides for use with 24 cm IPG-strips pH 3-10 OR pH 4-7.
IG_13_00	IPG-IEF in-gel default method for protein and peptides for use with 13 cm IPG-strips pH 3-10 OR pH 4-7.
DIAGNOSE	use for hardware diagnosis as described in "Checking high voltage, electrodes and electrode connectors" on page 94.

Table 5Default methods

Preparing a tray for the OFFGEL mode

To prepare a tray for the OFFGEL mode, refer to the *Kit Manual* (5969-1582) for detailed instructions. The manual is included in every consumable kit (see Table 1 on page 17).

Preparing a tray for the standard IPG IEF mode

You can use any IPG-strip for fractionations in the standard IPG IEF mode (in-gel mode). Follow the recommendations of the IPG-strip supplier about rehydration loading of the sample to the IPG-strip. The tray of the 3100 OFFGEL fractionator system can be used as a rehydration tray.

The following workflow uses an IPG-strip already loaded with sample and describes an example of the tray preparation for one IPG strip.

- 1 Place the IPG-strip in a new tray lane with the gel side up. The low pH side (Anode), marked with a "+", must be at the left, where the fixed electrode will be attached.
- **2** Take two electrode pads with tweezers, wet them with rehydration buffer, and place one at each end of the IPG-strip so that it connects to the swelled gel.
- **3** Attach the electrodes to the tray exactly as described in the OFFGEL *Kit Manual* (5969-1582). If necessary, reposition the IPG-strip with the electrode pads so that the fixed electrode connects to the left electrode pad.
- **4** Place the tray into the instrument.
- **5** Cover the complete tray lane with mineral oil using 1.5 ml mineral oil per IPG-strip.

NOTE Distribute the mineral oil evenly across the whole tray lane by pipetting droplets across the lane. Make sure that no mineral oil is spilled out of the tray.

6 Start the run as described in "Starting a run" on page 67.

4

Starting a run

- **1** To start the fractionation, make sure you have configured the cooling platform according to your needs. For details, refer to "Setting Cooling Platform Temperature" on page 47.
- 2 From the Start screen, press Tray Status [F1 or F5].
- 3 Select Method [F2].
- 4 Press Load [F4]; the load screen is displayed.
- 5 Use the Selection Keys to select the method from the list
- **6** Press the **Load** button [**F6**] to load the method; the run parameters of the loaded method are shown and the method name is displayed on the status bar of the method screen.
- 7 Press the Done button [F6] to load the method to the selected tray
- 8 Press Start [F7] so start the run.

NOTE

A blue light below the anode connector of the respective tray indicates that the instrument is running. A blinking blue light indicates that the fractionation is finished.

Both trays can be used independently. If the top cover is opened, the fractionation pauses until the top cover is closed again.

Performing a Fractionation

Agilent 3100 OFFGEL Fractionator User Manual

Maintenance

5

Cleaning the 3100 OFFGEL Fractionator 70 Firmware Update 76

This chapter describes cleaning procedures that should be used to do maintenance on your instrument at regular intervals or before every use. In addition, it provides instructions on how update the firmware. Only perform maintenance on the Agilent 3100 OFFGEL Fractionator as specified in this chapter. Other maintenance and repairs must be done by Agilent Technologies trained personnel. Unauthorized maintenance can be dangerous and damage is not covered by warranty.

5 Maintenance

Cleaning the 3100 OFFGEL Fractionator

Cleaning the 3100 OFFGEL Fractionator

Materials used in the instrument that may come in contact with liquid reagents are resistant to chemicals typically used for isoelectric focusing (IEF).

WARNING

Always turn off the instrument and unplug the power cord before doing any cleaning procedures.

Always wear powder-free gloves when handling parts that are in contact with sample or mineral oil.

CAUTION Do not use ketones, such as acetone, or strong acids to clean the instrument

Cleaning the electrodes

Regular cleaning procedure

Both the fixed and the movable electrodes have to be cleaned before every use.

- **1** Wipe the electrode pins with water using a lint-free tissue.
- **2** Wipe the electrode pins with isopropanol using a lint-free tissue to remove the mineral oil.
- **3** Thoroughly air dry the electrodes or dry well with a lint-free tissue prior to use.

Cleaning procedure for heavily contaminated electrodes

If for any reason not only the electrode pins but the complete electrodes are contaminated use the following cleaning procedure:

5

- 1 Rinse the complete electrodes with hand-warm tap water. Do not exceed a water temperature of 50°C.
- **2** Use a mild detergent or a 10% SDS solution and a brush to clean the electrodes.
- **3** Carefully rinse the electrodes with hand-warm tap water to remove the detergent.
- **4** Carefully rinse the electrodes with distilled water to remove remaining salts.
- **5** Thoroughly air dry the electrodes prior to use to avoid any electrical short circuit.

CAUTION Do not use heat such as a hair dryer to dry the electrodes since the plastic parts of the electrodes may be deformed.

NOTE

Additional electrodes are also available as spare parts see Table 7 on page 150

Cleaning the Trays

Trays have to be cleaned before every use to avoid cross-contamination of protein or peptides from one sample to the other. This is especially necessary when the instrument is used in the IPG IEF mode (in-gel mode).

- 1 Rinse the complete tray with hand-warm tap water. Do not exceed a water temperature of 50°C.
- 2 Use a mild detergent or 10% SDS solution and a brush to clean the tray.
- 3 Carefully rinse the tray with hand-warm tap water to remove the detergent.
- **4** Carefully rinse the tray with distilled water to remove remaining salts.
- **5** Thoroughly air dry the tray or dry well with a lint-free tissue prior to use.

5 Maintenance

Cleaning the 3100 OFFGEL Fractionator

CAUTION

Do not use heat such as a hair dryer to dry the trays since the material may be deformed.

NOTE

Additional trays are also available as spare parts see Table 7 on page 150.

Cleaning the cooling platform

If for any reason the cooling platform is contaminated, use the following cleaning procedure:

- 1 Wipe the cooling platform with water using a lab tissue or paper towel.
- **2** Wipe the cooling platform with isopropanol using a lab tissue or paper towel to remove the mineral oil.
- **3** Thoroughly air dry the cooling platform or dry well with a lab tissue or paper towel prior to the next fractionation.

NOTE

A leak sensor is installed on the left side of the rear of the cooling platform. The sensor and its surroundings has to be completely dry before it is possible to start the next fractionation.

Cleaning the anode connector

If you see heavy contamination of the fixed electrode, such as salts or mineral oil, it is strongly recommended to also detach and clean the anode connector.
5

1 Loosen the 4 screws on top of the anode connector with a hex key.

Figure 37 Disassembling the anode connector.

2 Carefully lift the anode connector from the front (Figure 38). Take care to not damage the cables at the rear of the anode connector.

5 Maintenance

Cleaning the 3100 OFFGEL Fractionator

Figure 38 Detached anode connector.

- **3** Wipe the anode connector with water using a lab tissue or paper towel.
- **4** Wipe the anode connector with isopropanol using a lab tissue or paper towel to remove the mineral oil.
- **5** Thoroughly air dry the anode connector or dry well with a lab tissue or paper towel prior to the next fractionation.
- **6** Return the anode connector to its original position and fasten the screws with a hex key.

Cleaning the top cover and the housing of the instrument

The top cover and the housing can be cleaned with a damp cloth either using water or a mild detergent.

Do not use alcohols or ketones for the top cover, since the material is not resistant to these solvents.

5 Maintenance Firmware Update

Firmware Update

At regular intervals, Agilent Technologies releases new firmware revisions of its products to provide new features or bug fixes. These revisions and the Firmware Update Tools are provided via the Agilent Technologies web page http://www.agilent.com (Life Science and Chemical Analysis Solutions -Technical Support - Downloads & Utilities) or direct via http://www.chem.agilent.com/scripts/DownloadsUtilities.asp

This section describes how you can update the firmware of the Agilent 3100 OFFGEL Fractionator instrument and the control module. It requires

- installation of the Firmware Update Tool,
- installation of the firmware onto a PC,
- connection of the PC to the instrument using a RS232 cable.

NOTE

Use a 9-to-9 pin (female) connector RS-232 cable (Agilent part number G1530-60600 or RS232-61601). Such a cable is included in the accessory kit (see Table 3 on page 27) that is delivered with the system

This is a so-called "Null Modem Cable" with full handshaking, where the connection are made between pins 1-1, 2-3, 3-2, 4-6, 5-5, 6-4, 7-8, 8-7, 9-9.

Installation of the Firmware Update Tool

1 Download the Firmware Update Tool from the Agilent web site (http://www.chem.agilent.com/scripts/DownloadsUtilities.asp) onto your PC.

NOTE

Do not use older versions of the Firmware Update Tools as provided via the Agilent web. Versions below 2.3 [001] (introduced 07/2006) should no longer be used.

2 Double-click the exe-file to install the program files. Follow the instructions. The program is installed into a folder (default is *C:\Program Files\Agilent\LC Firmware Update Tool*). You may change this according your needs. A Firmware Update Tool icon is installed onto your desktop, and a program folder is installed under Start - Programs (see Figure 39).

How to un-install the Firmware Update Tool

The Firmware Update Tool has an un-install routine. If required, remove the Firmware Update Tool via the "**Uninstall Firmware Update Tool**" function within the program folder, see Figure 39.

Installation of the Firmware on your PC

- 1 Download the Agilent 3100 OFFGEL Fractionator firmware from the Agilent web site (http://www.chem.agilent.com/scripts/DownloadsUtilities.asp) onto your PC.
- **2** Double-click the exe-file to extract the content into a folder (default is *c:\Fwupgrad\lc*).

NOTE

The Control Module firmware should be in the same folder as the Agilent 3100 OFFGEL Fractionator instrument firmware This will ensure that the firmware is found for complete system updates.

Connect your PC to the Agilent 3100 OFFGEL Fractionator

Your PC has to be connected by RS-232 cable to the Agilent 3100 OFFGEL Fractionator. For a description of connections at the rear of the instrument see Figure 7 on page 28.

Firmware Update

The FW Update Tool uses the following port settings when performing the update:

- Baudrate 19200 (overrides the PC settings)
- No Parity
- 8 Data Bits
- 1 Stop Bit
- COM1, COM2 etc.

The settings are set automatically for the time of the firmware update.

NOTE

Note that the upgrade process can take up to 2 hours. No fractionations can be performed during this time.

Setup of the Firmware Update Tool

1 Start the Firmware Update Tool via the shortcut (see Figure 39). The following initial screen is displayed:

When you move with the cursor over a information is given about its function. Alternatively, you may use the menu items.

File	Open Protocol File - Exit
Action	Connect - Disconnect - Select for update - Identify module - Start/Cancel update
Options	Configure directories - Add/Delete/Modify Instrument - Support G1323B (active/not active)
Help	Help - About (revision of program)

Table 6 Functions of Button

Button	Description
₽	Connect to module/system
₹	Disconnect from module/system
✓	Select for update of all modules of a system
	Start the update for all selected modules
.	Add instrument to the list
ā	Delete instrument from the list
*-	Modify instrument (name or the connection can be changed)
C	Configure directories of the firmware files and report
	Open protocol file (logbook with details of the actions made during the update)

Configure the Firmware Update Tool

1 Select the button **Configure Directories** and browse for the location of the firmware files (e.g. E:\Fwupgrad\LC).

🔆 Firmware Update Tool	
<u>File Action Options H</u> elp	
	Ĩ

Figure 41 Select the button Configure Directories

Configure Directories	×
Set the directory that contains the firmware files and specify the location of the protocol file.	
Firmware directory: E:\Fwupgrad\LC	
Protocol file: C:\Program Files\Agilent\LC Firmware Update Tool\FWReport.txt	
<u>Q</u> K <u>C</u> ancel	

Figure 42 Configure Directories

NOTE

If you want to update both the instrument and control module firmware simultaneously, place the firmware files for the Control Module als in this firmware directory.

The Protocol file is a text file that gives the complete status of the firmware update process. This file is useful as a compliance document to prove that the module firmware files were updated successfully or as diagnostics information in the case of problems, see "The Log-file" on page 86. Therefore, name the protocol file appropriately.

2 Press OK.

3 Select the **Add Instrument** button.

Figure 43 Select the Add Instrument button.

- 4 Enter a name for the system (e.g. System 1)
- **5** Enter the COM port connection information in the instrument address field, e.g. COM1 .

Add Instrument	×
Select an instrument name/address from the lists below or edit the text fields directly.	
Instrument Name:	_
System I	-
Instrument Address:	
СОМ1	•
<u>Q</u> k <u>C</u> ancel	

Figure 44 Add instrument screen

Firmware Update

6 Press OK.

ጅ Firmware Update Tool				
File Action Options Help				
₽ ₽ ✓ ► × 3 5 5 4 9 1				
C C All Instruments				
System1	All Instruments:			
	Instrument	Address		
	System1	COM1		

- 7 Double-click on **All Instruments** to expand the view.
- 8 If you want to update the Control Module, select **Options** and check **Support Control Module**.

Firmware Up tote Tool				
File Action Options Help				
로 🛒 🖓 Configure Directories				
🖃 🚐 All Ins				
📖 🖻 👝 Remove Instrument	All Instruments:			
📛 Modify Instrument				
	Instrument	Address		
Support Control Module	System1	COM1		

9 Select **System 1** and click on the **Connect to Instrument** button.

🖗 Firmware Update Tool		
File Action Options Help		
2// / / × 1555 4		
E - E All Instruments		
System1		
	Instrument	System1
	Address	COM1

Figure 47 Select system 1

All modules within System 1 are listed with serial number, firmware revision and build number. G3100A refers to the Agilent 3100 OFFGEL Fractionator instrument and G3101A refers to the control module.

🔆 Firmware Update Tool			
Program Action Options Help			
₽₽ √ ► + - ½			
🖻 🚎 System 1			
G3100A:DE00000038	Instrument System 1		
G3101A:CN53314946	COM1		
	Address COMI		
	Module	Revision	Build
	G3100A:DE0000038	U.00.20	[0001]
	G3101A:CN53314946	T.01.08	

Figure 48 Modules of System 1 are listed

Doing the Firmware Update

You have the possibility to update the firmware for the control module and instrument separately (e. g. if both do not need updating) or at the same time. The following describes the general workflow for an individual update and highlights the differences to a complete system update.

1 Select a module.

🔆 Firmware Update Tool		
Program Action Options Help		
코 로 √ ▶ + - ½ 🛱 🗎		
E-œ⊐ All Instruments È-œ⊒ System 1	Current Firmware Dodate Information	n Revision History
G3100A:DE00000038 G3101A:CN53314946	Module: G3100A:DE0000	0038 Q Identify Module
	Main System:	Resident System:
	Revision Build #	Revision Build #
	U.00.20 [0001]	U.01.04 [0151]

Figure 49 Specific module is selected - Current Firmware Information screen

- **2** Above screen shows the current firmware information
- **3** Select the **Revision History** tab. This window shows details of the used main board and the update history in the module, see Figure 50.

Figure 50 Specific module is selected - Revision History Information screen

4 Select the **Update Information** tab. In this window, you can select the firmware for the main and resident part in the module.

If no firmware is listed when you click on the file selection button, correct the firmware folder locations, see Figure 42 on page 80.

5 Activate the check boxes for the firmware. This firmware part is added to the update list.

🔆 Firmware Update Tool		
Program Action Options Help		
	Current Firmws e Update Information Dievis Main System:	sion History
	Update Filename I 3100A_U020_001.dlb	Revision Build #
	Resident System:	Revision Build #
	Res_U104_151.db	U.01.04 151

Figure 51 Update Information screen

- **NOTE** You may now select the G3101control module in addition for multiple module update. The firmware is added to the update list.
 - **6** Select the **Start Update** button. All selected modules are updated (resident first, then main firmware). The update status is shown at the bottom of the screen.

When the main firmware is being updated, the Status LED of the instrument is red and blinking.

If the Control Module has also been selected for update, the Control Module update starts first.

Module	System	Instrument	Progress [%]	A
G3100A:DE00000038	Main	DE_PP00038		Start Update
G3100A:DE00000038	Resident	DE_PP00038		
G3101A:CN53314946	G1323B	DE_PP00038		CancelUpdate
				<u></u>
				-
Download in progress	nlease wait			

If you cancel an update, the module remains in the current mode, resident or main, depending when the process was stopped. It cannot be used until the update process has been restarted.

At the end of the update, a report is shown. The update status is either PASS (successful) or FAILED (not successful), see Figure 53. If the update failed, the Protocol file contains more information, see "The Log-file" on page 86.

C:\data\Tanja\OGE_Prome	theus\Userinfo\sw relaese\F	WReport.txt	
Module	System	Result	^
G3100A:DE00000037	Main	PASSED	
G3100A:DE00000037	Resident	PASSED	

- 7 If desired, open and print the Protocol file.
- 8 Press **OK** to continue.
- **9** Perform a power cycle of the instrument.

The Log-file

This file is the verification that the firmware update was completed successfully or, in the case of problems, the log-file may give some specific information.

Press the button **Open Protocol File** to open the report. Information can be copied/pasted into an e-mail for further help from Agilent support.

<mark>Ж</mark> ғі	irmware	e Update	Tool			
Eile	<u>A</u> ction	Options	Help			
2	杲、	()	×	i	ä	

Figure 54 Select the button Open Protocol File

📕 FWReport.txt - Notepad 📃 🗖 🛛
File Edit Format View Help
[8/8/2006 13:47:28] Download for G1323B-System of module G3101A:CN23605354 (REV="T.01.05", Build="") with file C:\data\Tanja\OGE_Prometheus\Userinfo\sw relaese\ogt106en.bin (REV="1.01.06" Build="") started
[8/8/200614:20:41] Download for G1323B-System of module G3101A:CN23605354 succeeded [8/8/2006144:21:52] Download for Decident System of module
G3100A:DE00000037 (REV="U.01.04", Build="[0151]") with file C:\data\Tanja\OGE_Prometheus\Userinfo\sw relaese\Res_U104_151.dlb (REV="U.01.04" = Build="[151") started
[8/8/2006 14:26:54] Download for Resident-System of module G3100A:DE00000037 succeeded [8/8/2006 14:27:41] Download for Main-System of module
[5]002-000 +:27.04-] REV="1.00.12", Build="[0002]") with file (S100A-R:DE00000037 (REV="1.00.12", Build="[0002]") with file (C:\data\Tanja\OGE_Prometheus\Userinfo\sw relaese\Prometheus_U12_02.dlb
[8/8/2006 14:29:04] Download for Main-System of module [331004-R:DE00000037 succeeded

Figure 55 Example of log-file with successful update of both modules

If you encounter any problems with the firmware update, feel free to send an e-mail with your findings to our technical support specialists

(http://www.chem.agilent.com/scripts/cag_techsupport.asp) or contact your local Agilent Technologies representative. Attach your log-file and details of your system, revision and errors.

5 Maintenance

Firmware Update

Agilent 3100 OFFGEL Fractionator User Manual

Troubleshooting

Troubleshooting the Instrument90Troubleshooting the Control Module97Troubleshooting the Application98

This chapter provides support for troubleshooting the instrument, the control module and the application.

Troubleshooting the Instrument

Status Indicators

Status LED

The instrument has a status LED on the front top right position (see Figure 2 on page 15). With the power switched on, the LED can have the following states:

LED Status	Description	Comment
Green	The instrument is performing a fractionation.	
Yellow	 The instrument is currently not ready. Possible reasons: the set temperature of the cooling platform has not been reached top cover is open 	A running fractionation will be paused if the top cover is opened.
Red	The instrument has detected an error.	A running fractionation will be aborted. see Table on page 92
Off	The instrument is not performing a fractionation. The set temperature of the cooling platform has been reached.	

Table 7	LED	status	overview

Status bars

The control module displays several status bars (see "Start Screen" on page 35)

Status bar title	Description	Comment
Cover open	The top cover of the instrument is open.	A running fractionation will be paused if the top cover is opened
Temp	The set temperature of the cooling platform has not been reached.	You can start a fractionation with the control module; however, the status bar will switch to Wait Temp I/II.
Wait Temp I or Wait Temp II	The set temperature of the cooling platform has not been reached; however, a fractionation was already started on one of the trays.	Once the cooling platform temperature is reached the fractionation will be started automatically.

Table 8 Status bar title overview - LED is yellow

Error Messages

Error message	Solution
USB drive full. Exchange USB drive or delete files.	Exchange the USB drive or delete files on the USB drive.
File export failed. The selected files cannot get copied to USB drive. Check if drive is full.	Exchange the USB drive or delete files on the USB drive.
Result Storage Empty. No result data to copy. CF-Card is empty or not present	Insert a compact flash card (for details refer to Figure 7 on page 28)
No USB drive. USB drive not present or incompatible drive format.	Insert a USB drive. Make sure the format is compatible (for details refer to page 59).
Cannot overwrite default method. A default method cannot be overwritten. Change method name.	Change the method name to save modification of default methods.
The method file <name> on <drive> cannot be overwritten. It is currently running on <tray>.</tray></drive></name>	Change the method name to save modification to this method or save the method once the fractionation is finished.
Delete Denied. The method <name> from <drive> is currently in use on <tray>.</tray></drive></name>	Wait until the fractionation is finished to delete this method.

Table 9 Error messages due to wrong operation

Troubleshooting the Instrument

Error message	Comment	Action
Cooling platform temperature error (EE 20302).	The cooling platform temperature has not been reached.	 Make sure the ambient temperature and the set cooling platform temperature are within instrument specifications. A Peltier element might be defective. Call Agilent Technologies.
A leak was detected in the module. Use the logbook to determine date and time the fault occurred (EE 00064)	A leak has been detected by the leak sensor.	 Check if there is liquid spillage or condensation on the cooling platform and near the leak sensor and wipe it away. For details refer to "Cleaning the cooling platform" on page 72. Acknowledge the error in the logbook to reset the red LED status.
Cover sensor failed (EE 20300).	Opening and closing the top cover is no longer recognized.	 Check that the top cover closes tightly. The cover sensor might be defective. Call Agilent Technologies.
Temperature sensor failed (EE 20301).		 The temperature sensor is defective. Call Agilent Technologies.
Fan failed (EE 00068).		 The fan might be defective. Call Agilent Technologies.
Leak (EE 00080, EE 00081, EE 00082 or EE 00083).		• The leak sensor is defective. Call Agilent Technologies.

Table 10 Error messages hardware related - LED is red

Error message	Comment	Action
High voltage (EE 20305).	The voltage did not reach the set point.	 Check for condensation or dirt on anode connector, electrodes or tray. Clean the the parts as described in "Cleaning the 3100 OFFGEL Fractionator" on page 70. The high voltage module might be defective. Call Agilent Technologies.
Peltier Error (EE 20304).	The voltage or currents of the Peltier elements are out of limits.	 A Peltier element might be defective. Call Agilent Technologies.

 Table 10
 Error messages hardware related - LED is red (continued)

Hardware Diagnosis

With the following few steps, you can verify if your instrument is performing to its specifications. Note that a meaningful hardware diagnosis requires a fully functional control module (see "Troubleshooting the Control Module" on page 97).

Checking the LED and the cover sensor

1 Switch on or restart the instrument and watch the Status LED during startup.

The Status LED **must** show all three LED colors for approximately one second during startup in the following order: yellow-green-red.

If this is not the case, a LED is broken and needs repair. Call Agilent Technologies for further help. Otherwise, proceed with step 2.

- **2** Set the cooling platform temperature to 20 °C and wait until the cooling platform temperature is reached and the LED is turned off.
- **3** Open the top cover.

Troubleshooting the Instrument

The LED **must** turn yellow. The control module must show **Cover Open** in the status bar. This confirms that the cover sensor is working properly. A failed cover sensor would be recognized by the instrument and result in an error, see Table on page 92.

Checking temperature control of the cooling platform

NOTE

For this test, the ambient temperature must be within specification (maximum 30°C if the cooling platform is set to 20 °C) and conditions should be non-condensing.

- **1** Set the cooling platform temperature to 20 °C and wait until the cooling platform temperature is reached and the LED is turned off.
- **2** Set the cooling platform temperature to 30 °C and measure the time until the LED is turned off.

The temperature ramp up should not take longer than 5 minutes.

3 Set the cooling platform temperature back to 20 °C and measure the time until the LED is turned off.

The temperature ramp down should not take longer than 15 minutes.

Checking high voltage, electrodes and electrode connectors

- 1 Make sure your two electrode sets and trays are completely dry and clean.
- **2** Assemble both trays with the electrodes (without frames) and insert them to the instrument. Attach the movable electrodes at the most right position.
- **3** Load the DIAGNOSE method to both trays and press **Start** [**F8**] on both trays to start the runs.
- 4 Check that the actual voltage values reach 7500 8500 V after one minute and that the currents of all lanes are < 1 μ A.

If this is not the case, check once again that the tray and the electrodes are completely dry and clean. Inspect the cathode connector for contamination and clean it if necessary. Make sure that no sparking occurs. If you cannot solve this problem call Agilent Technologies for further help. Otherwise, proceed with step 5.

- **5** Press **Stop** [**F8**] to stop both runs.
- **6** Open the top cover and pipette 4 ml deionized water (18 MOhm/m) into lanes A, C, E and G on each tray.

7 Close the top cover.

NOTE

Make sure the water is evenly distributed and the electrode pins are in contact with the water. However, avoid spillage of water across the electrodes into neighboring lanes and the electrode connectors. This might result in sparking.

- 8 Press Start [F8] on both trays to start the DIAGNOSE method again.
- **9** Look at both tray status screens. You should see a regular current pattern as shown in Figure 56.

l Tra	ay I	0.077 kVh	DIAGNOSE	Ready	10:	15:38 🗆
Focusi	ng	0:03h:m	Current	A - 45.7		
	Limit	Actual	1	B - 0.0		
Volt.Hours	20	0.076 kVh		D • 0.0		
Voltage	0008	1333 v		E • 50.0 F • 0.0		
Power	200	66.7 mW		G • 43.7 H • 0.0		
Current	50.0	50.0 μA	ö .o	[µĂ]	50.0	
Plot		Method	Logbook	Other Tray	l Re	sults

Figure 56 Typical current pattern on tray I when every second lane is filled with water.

No currents should flow in the empty tray lanes (<1.0 μ A). If this is the case proceed with step 10. Otherwise stop here because you may have contaminated electrodes, electrode connectors or trays. Clean and carefully dry the electrodes, trays and the electrode connectors as described in "Cleaning the 3100 OFFGEL Fractionator" on page 70 and start this test again at step 6.

Typical values for currents are between 20 μ A and 50 μ A and one lane has to be at maximum current of 50 μ A. A typical voltage value for this test is between 1000 volts and 3000 volts. The voltage has to be greater than 500 volts.

NOTE

The absolute value for the voltage is dependent on the conductivity of the deionized water used and is not a criterion for the proper function of the instrument. The more salt ions the water contains, the lower is the voltage value.

Troubleshooting the Instrument

- **10** Open the top cover and pipette 4 ml deionized water (18 MOhm/m) into lanes B, D, F and H on each tray.
- **11** Close the top cover.
- **12** You should now see the a typical current pattern for all eight lanes (Figure 57).

Tra	ay I	0.132 kVh	DIAGNOSE	Ready	10:2 <i>′</i>	1:39 🗆
Focusi	ng	0:05h:m	Current	A • 40.1		Stop
1	Limit	Actual		B • 48.0		
Volt.Hours	20	0.131 kVh		C - 45.6 D - 50.0		
Voltage	8000	1306 v		E - 45.5 F - 45.4		
Power	200	65.3 mW		G - 35.7 H - 46 1		ĺ
Current	50.0	50.0 μΑ	0.0	[µA]	50.0	
Plot	Ì	Method	Logbook	Other Tray	Resu	ılts

13 Press **Stop** [**F8**] on both trays to stop the runs.

NOTE

If you require further assistance, call Agilent Technologies and send the corresponding results file of this diagnose test to the customer call center.

Troubleshooting the Control Module

	Internal diagnostics continuously monitor the module's condition and record any unusual events in an electronic logbook. The logbook is self-updating, where the newest entry replaces the oldest entry.
NOTE	Only entries that take place during a fractionation are permanently saved with the fractionation log file, all other entries are deleted with every power cycle of the instrument.
	If your control module does not work correctly, disconnect the module CAN connector from the rear of the Agilent 3100 OFFGEL Fractionator and reconnect it.
NOTE	Make sure you are using the CAN connector and not the LAN connector.
	If the problem still remains, then
	• power down the Agilent 3100 OFFGEL Fractionator and wait 1 minute and then restart.
	If the problem still remains, then
	• try to connect the control module to an Agilent 1200 Series module if available. A visible start screen on the Agilent 3100 OFFGEL Fractionator indicates the general function of the control module. This result would point to a defective CAN connector or a firmware/software problem.
	Call Agilent Technologies if the problem cannot be solved.

Troubleshooting the Application

Troubleshooting the Application

Impaired focusing results can have many different root causes. In addition to hardware errors, the impaired focusing results can be sample and/or handling related. Therefore, always make sure that you set up your focusing experiment exactly as described in the kit manual. The sample concentration should be between 50 μ g/ml and 5 mg/ml and the salt content of your sample buffer should be at most 10 mM. Run only samples similar in protein content and buffer composition on the same tray, otherwise it is highly recommended that you use separate trays.

With every focusing, a results file containing all electrical parameters during a run are stored to the compact flash card and can be exported. The traces also provide useful information for troubleshooting (see "How to display exported results" on page 60), and some examples are discussed in this chapter (see "Looking at Run Parameters" on page 101)

The following table gives an overview of symptoms you might experience when performing OFFGEL experiments, their possible causes and actions required for remedy.

Symptom	Possible Cause	Action
No current or very low current (less than 20% of the current limit) during the first 5 minutes after start of a run in one lane compared to other lanes.	1 No electrical contact or impaired contact.	 check if there is any damage to the electrodes . If so, exchange the electrodes. check if there is any damage to the anode connector. If so, call Agilent Technologies. Make sure that the electrodes are in place correctly and the electrode pins are connecting to the electrode pads. Make sure that the electrode pads are damp and provide a good contact between the electrode pins and the gel ends of the IPG-strip.
	2 IPG strip is improperly rehydrated.	 Make sure that the IPG-strip was rehydrated completely and that all wells in the frame are filled with rehydration buffer during rehydration.
One lane shows very high current signal compared to other lanes on the same tray. This might impair the focusing result of all other lanes on that tray.	1 the corresponding lane is leaking.	 inspect the corresponding frame for any damage. If damage is present, exchange the frame. Make sure that the IPG-strip is inserted with the gel side up. Make sure that the corresponding frame is snapped into the tray lane as described in the kit manual.
	2 the corresponding IPG-strip was inserted in the opposite direction.	 Make sure that the "+" on the IPG-strip is in the direction of the fixed electrode/anode connector.

Table 11	Troubleshooting the OFFGEL application
	indubication and the off off off officiation

Troubleshooting the Application

Symptom	Possible Cause	Action
	3 the corresponding sample has higher salt concentration or higher protein/peptide content than the other samples.	 Run only samples similar in protein content and buffer composition on the same tray.
After fractionation has finished, the external wells of the frame have reduced liquid or are empty.	During IEF, a water transport is observed due t migrating ions that drag th water in their hydration shell to the electrodes. Furthermore, water is dragged from wells at the basic strip end (cathodic end of strips pH 3-10) to wells in the neutral pH range due to the ionization of DTT. Therefore, after OFFGEL electrophoresis th wells close to the electrodes are sometimes filled with less liquid than the inner wells or are even empty.	 If low recovery in the extreme wells is a problem, refill the corresponding wells with diluted OFFGEL Stock Solution without ampholytes (4 parts OFFGEL Stock Solution without ampholytes plus one part water) one hour prior to the end of the focusing phase.
Sparking at the electrodes.	1 short circuits at the electrodes.	 make sure that the electrodes and the electrode connectors are completely dry.
After fractionation has finished, the cooling platform is contaminated with mineral oil.	 Mineral oil tends to creep across the electrodes and tray wells once they are moistened. However, if the electrodes do not dry out during fractionation this does not impair the focusin results. 	 carefully pipette the mineral oil and avoid spillages. when inserting the tray into the instrument, make sure that no oil is spilled out of tray. Insert the tray before adding the mineral oil to the instrument.

Table 11 Troubleshooting the OFFGEL application (continued)

Looking at Run Parameters

The results of all run parameters can be exported as a .csv file and thus can be visualized with Microsoft Excel. Agilent Technologies has designed the *Agilent 3100 OFFGEL fractionator plotter* which is a Excel add-on and facilitates the analysis of the run parameters (see "How to display exported results" on page 60). The following examples show typical results for good runs as well as some impaired results where the root cause of the problem is known.

Example 1. Typical results of the protein test sample

The fractionation of the protein test sample is a way to check that the Agilent 3100 OFFGEL Fractionator in combination with the consumables kit is working in your hands, and functions as a system checkout test. You find detailed information on how to run the protein test sample in the kit manual.

Figure 58 shows a typical result for the electrical run parameters of the protein test sample using the 3100 OFFGEL Low Res Kit, pH 3-10. The method used for this fractionation was the default method OG12PR00, which has a current limit of 50 μ A, a voltage limit of 8000 V, a power limit of 200 mW and a volt hour limit of 20 kVh. With this setup the fractionation is running in current limiting mode. Therefore, the current signal in Figure 58 is stable and at 50 μ A.

Troubleshooting the Application

Figure 58 Example 1- typical results for the Protein test sample using the 3100 OFFGEL Low Res Kit, pH 3-10.

The run stopped when the volt hour limit of 20 kVh was reached. Both volt signal and power signal are increasing throughout the run, since the sample is depleted of ions and thus the resistance is increasing. After fractionation, the protein test samples should be focused in fraction 3 (see Figure 59).

Troubleshooting the Application

Figure 59 Beta-lactoglobulin Protein Test Sample in well 3 of a 12 cm pH 3-10 IPG strip. (Gel-like image from an Agilent 2100 bioanalyzer analysis using the Protein 80 kit of wells 1-9.)

Example 2 . Current drop during fractionation

Figure 60 shows an example where the current (see arrow) dropped close to zero after approximately 30 minutes of the run. In addition, the volt and power signal is highly fluctuating.

Troubleshooting the Application

Figure 60 Example 2- current drop during fractionation.

The root cause of this focusing failure was that the IPG strip was dried out at one position so that the conductivity was lost.

Example 3. Reduced electrical charge in one lane

Figure 61 shows an example of a fractionation of eight samples in parallel and just highlights the different electrical charges (in μ Ah) that were measured throughout the focusing for each lane. The signal of electrical charge for lane F (see arrow) is clearly lower than for the other lanes. Indeed, it could be verified in a follow-up experiment that the focusing of lane F was impaired.

Troubleshooting the Application

Figure 61 Example 3- one lane shows reduced electrical charge.

The root cause of this failure was a conductivity problem in lane F.

Troubleshooting the Application

This chapter contains safety information relevant to the Agilent 3100 OFFGEL Fractionator.

7 Safety Safety Information

Safety Information

The following general safety precautions must be observed during all phases of operation, service, and repair of the Agilent 3100 OFFGEL Fractionator.

All safety instructions should be read and understood before installation, operation and maintenance of the instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture and intended use of the instrument.

Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

WARNING A warning alerts you to a situation that could cause physical injury or damage to the equipment. Do not proceed beyond a warning until you have fully understood and met the indicated conditions.

CAUTION

A caution alerts you to situations that could impair Agilent 3100 OFFGEL Fractionator functionality or performance and lead to loss of data. Do not proceed beyond a caution until you have fully addressed the indicated conditions.

NOTE

Instrument Safety

This is a Safety Class I instrument (provided with a terminal for protective grounding) and has been manufactured and tested according to international safety standards.

Power Cords

Different power cords are offered as options with the Agilent 3100 OFFGEL Fractionator. The female end of all power cords is identical. It plugs into the power-input socket at the rear of the instrument. The male end of each power cord is different and designed to match the wall socket of a particular country or region.
WARNING

Never operate your instrument from a power outlet that has no ground connection. Never use a power cord other than the Agilent Technologies power cord designed for your region.

WARNING

Never use cables other than those supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.

Operation

Before applying power, follow the installation instructions. Additionally the following must be observed.

Do not remove instrument covers when operating. Before the instrument is switched on, all protective ground terminals, extension cords, auto-transformers, and devices connected to it must be connected to a protective ground socket. Any interruption of the protective grounding will cause a potential shock hazard that could result in serious personal injury. Whenever it is likely that the protection has been impaired, the instrument must be made inoperative and be secured against any intended operation.

Any maintenance of the instrument under voltage should be avoided.

Do not replace components of the instruments with power turned on.

Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.

Do not install substitute parts or make any unauthorized modification to the instrument.

Safety Symbols

Table 12	Safetv S	vmbols us	ed on Inst	ruments and	in Manuals
	ourou, o	y	Ju Jii 1110 0	unionito unu	in manaalo

Symbol	Description
Â	The apparatus is marked with this symbol when the user should refer to the <i>User Manual</i> in order to prevent risk of harm to the operator and to protect the apparatus against damage.

Chemical and Biological Safety

WARNING

When working with solvents please observe appropriate safety procedures (for example goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

WARNING

If pathogenic, toxic, or radioactive samples are intended to be used in this instrument, it is the responsibility of the user to ensure that all necessary safety regulations, guidelines, precautions and practices are adhered to accordingly. This includes also the handling of genetically modified organisms. Ask your laboratory safety officer to advise you about the level of containment required for your application and about proper decontamination or sterilization procedures to follow if fluids escape from containers.

Observe all cautionary information printed on the original solution containers prior to their use.

Because leaks, spills, or loss of sample may generate aerosols, observe proper safety precautions.

Agilent 3100 OFFGEL Fractionator covers have not been designed as bioseals for aerosol or liquid containment.

Handle body fluids with care because they can transmit disease. No known test offers complete assurance that they are free of micro-organisms. Some of the most virulent – Hepatitis (B and C) and HIV (I-V) viruses, atypical mycobacteria, and certain systemic fungi – further emphasize the need for aerosol protection.

Always follow local state and federal biohazard handling regulation when disposing of biohazardous waste material e.g. contaminated compartment frames.

Handle all infectious samples according to good laboratory procedures and methods to prevent spread of disease.

Dispose of all waste solutions and products according to appropriate environmental health and safety guidelines.

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Seguridad

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This chapter contains details on non-condensing conditions as well as ordering information for consumables and spare parts.

Non-condensing conditions

Non-condensing conditions

The Agilent 3100 OFFGEL Fractionator is specified to be operated only under non-condensing conditions and with a cooling platform temperature of maximum 10 $^{\circ}$ C below ambient temperature.

For example, assuming a typically experimental setup of 20 °C on the cooling platform and an ambient temperature of 25 °C, the humidity must be below 73 %. A plot of humidity versus ambient temperature for the cooling platform at 20 °C, 25 °C and 30 °C is depicted in Figure 1. For non-condensing conditions the humidity has to be below the plotted curves.

Figure 1 Plot of humidity versus ambient temperature assuming the cooling platform to be set at 20 °C, 25 °C or 30 °C. For non-condensing condition the humidity has to be below the plotted curves.

Ordering Information

List of Consumables

Product	Order Number	Product Description
3100 OFFGEL, 3-10, Low Res Kit	5188-6425	includes frames (12 well), IPG-strips pH 3-10 (13 cm), cover seals and fractionation reagents for 12 samples.
3100 OFFGEL, 3-10, High Res Kit	5188-6424	includes frames (24 well), IPG-strips pH 3-10 (24 cm), cover seals and fractionation reagents for 12 samples.
3100 OFFGEL, 4-7, Low Res Kit	5188-6427	includes frames (12 well), IPG-strips pH 4-7 (13 cm), cover seals and fractionation reagents for 12 samples.
3100 OFFGEL, 4-7, High Res Kit	5188-6426	includes frames (24 well), IPG-strips pH 4-7 (24 cm), cover seals and fractionation reagents for 12 samples.
3100 OFFGEL, Starter Kit	5188-6444	includes frames (6x 12 well, 6x 24 well), IPG-strips pH 3-10 (3x 13 cm, 3x 24cm), IPG-strips pH 4-7 (3x 13 cm, 3x 24cm), cover seals and fractionation reagents for 12 samples
OFFGEL, Protein test sample	5188-6428	includes beta-lactoglobulin as a protein test sample, enough for more than 100 tests.

 Table 6
 Agilent 3100 OFFGEL Fractionator consumable kits

8 Appendix Ordering Information

List of Spare Parts

Product	Order Number	Product Description
Tray Set	G3100-60007	includes 4 replacement trays
Electrode Assy	G3100-60002	includes 1 tray, 1 fixed electrode and 1 movable electrode
USB-Memory Stick Kit	G4208-68700	includes 256 MB USB flash memory

Table 7 Agilent 3100 OFFGEL Fractionator spare parts

A

about screen, 46 Absatzmarken, 124 accessories, 27 achtung, 124 action key, 33 add instrument, 81 alphabetic characters, 34 alphanumeric characters, 33, 51 alphanumeric keypad, 34 analysis time programming, 53 anode, 17 anode connector, 16 cleaning, 72 aviso, 142

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baudrate, 78 bench space, 20

C

É^ÉO íióé. 114 cables de alimentación, 142 íióéÉ^ÉO. 114 CAN connection. 28 cathode, 17 cathode connector. 16 caution, 108, 114, 120, 130, 136 cleaning, 70 anode connector, 72 cooling platform, 72 electrodes. 70 top cover, 75 tray, 71 COM ports, 78 compact flash card, 28, 60 compatibilidad electromagnética, 143 condensation, 20

configure directories. 80 firmware update tool, 80 configure main screen, 43 connection CAN, 28 LAN. 28 USB. 28 consumable kits, 17, 149 contrast. 48 control module, 15 about screen, 46 configure main screen, 43 contrast. 48 date & time screen, 44 display, 48 display screen, 45 error messages, 91 export menu, 59 keypad, 34 keys, 32 LAN and service screen, 44 last result screen, 42 load screen, 41 logbook screen, 42 method screen, 39 plot screen, 38, 55 save as screen, 41 serial number, 46 software, 35 start screen, 29, 36 support, 82 tray status screen, 37 user interface, 35 welcome screen, 29 controlmodule timetable screen, 40 cooling platform, 16 cleaning, 72 temperature, 47, 94

cover seals, 17

D

Dämpfe, 125 date & time screen, 44 date and time settings, 47 default methods, 52 deleting a method, 52 dimensions, 23 direction keys, 33 display contrast, 45, 48 settings, 48 display screen, 45

Ε

electrode assy, 150 electromagnetic environment, 21 EMC regulations, 109, 115, 121, 131, 137 enter key, 33 entzündlicher Gase, 125 environment, 20 error messages, 91, 92 esc key, 33 export menu, 59 exporting results, 59

F

firmeware revision, 46 firmware doing the update, 83 log file, 87 update tool, 76 fixed electrode, 17 flammable gases, 115, 121, 131, 137 focusing step, 49 fractionation principle, 13 frame, 17 fumes, 115, 121, 131, 137

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hardware diagnosis, 93 hold step, 49 humidity, 23, 148

I

in-gel mode, 12, 14, 66, 71 installation instrument, 28 RS-232 firmware update tool, 76 instrucciones de seguridad, 142 Instrument top cover. 15 instrument accessories. 27 anode connector, 16 cathode connector. 16 cooling platform, 16 errors, 92 hardware diagnosis, 93 installation. 28 serial number, 15, 46 setup, 28 spare parts, 150 status LED, 15, 90 troubleshooting, 90 IPG IEF mode, 66 IPG strips, 12 IPG-strip, 66 isoelectric points, 12

Κ

keypad, 34 keys action, 33 alphanumeric, 34 direction, 33 enter, 33 esc, 33 keypad, 34 navigation, 33 numeric, 34 selection, 33

L

LAN and service screen, 44 LAN connection, 28 last result screen, 42 line frequency, 20, 23 line voltage, 20, 23 load screen, 41 loading capacity, 14, 98 logbook screen, 42 log-file, 86

Μ

Marken Absätze, 124 method deleting, 52 DIAGNOSE, 94 focusing step, 49 hold step, 49 loading, 41, 49 modifying, 49 offline mode, 39 online mode. 39 save as. 41 saving, 51 time-based steps, 53 timetable, 49 transfer. 52 method screen, 39 modifying a method, 49 movable electrode, 17

Ν

navigation key, 33 Netzanschluss, 125 Netzkabel, 125 non-condensing conditions, 23, 94, 148 non-operating altitude, 23 non-operating temperature, 23 numeric keypad, 34 numeric keys, 34

0

OFFGEL fractionator plotter, 60 OFFGEL mode, 12, 66 operating altitude, 23 operating temperature, 23 Operation, 115, 121, 125, 131, 137, 143

Ρ

paragrafi tag, 136 paragraph tags, 120, 130 pH resolution, 14 physical specifications, 23 pl. 12 plot screen, 38, 55 rescale, 56 signal unit axis, 56 time axis. 56 toggling signals, 57 x axis, 56 v axis. 56 plotter, 60 graph. 63 toolbar button, 61 user interface. 62 power considerations, 20 power consumption. 23 power cords, 108, 115, 120, 131, 137 power outlet, 109, 115, 121, 131, 137 power supply, 20 precaución, 142 precauciones de seguridad, 142 protein test sample, 17, 101, 103, 149 protocol file, 87

R

report screen, 86 rescale plot screen, 56 RF transmitters, 21 riesgo para la seguridad, 143 RS-232 baudrate, 78 connection, 28 firmware update tool installation, 76 firmware update tool setup, 78 run parameters, 101

S

safety hazard, 109, 115, 121, 131, 137 safety instructions, 108, 114, 120, 130, 136 safety precautions, 108, 114, 120, 130, 136 safety standards IEC. CSA. UL. 23 Safety Symbols, 116, 121, 126, 132, 138, 144 safety symbols, 110 sample concentration, 14, 98 save as screen, 41 screens system, 35 timetable, 51 select signals, 55 selection keys, 33 serial number, 15 setup firmware update tool, 78 instrument, 28 Sicherheitsanweisungen, 124 Sicherheitsmaßnahmen, 124 Sicherheitsrisiko, 125 signal unit axis, 56 software revision, 46 sound emission, 21 spare parts, 150

start a run, 67 start screen, 29, 36 status bar titles, 91 status LED, 15 status overview, 90 support control module, 82

T

tag di paragrafo, 136 tags paragraphs, 120, 130 temperature fluctuations, 20 time axis. 56 time programming, 53 timetable, 40, 49, 53 linear mode, 53 screen, 40 step mode, 53 toggling signals, 57 toma de alimentación, 143 top cover, 15, 75 tray, 17, 66, 71 cleaning, 71 preparing, 66 tray set, 150 tray status screen, 37 troubleshooting application, 98 control module, 97 instrument, 90

U

un-install firmware update tool, 77 unpacking, 26 update info screen, 85 USB connection, 28 USB stick, 59, 60, 150 folder directories, 60 user interface, 35

V

vapores inflamables, 143 voltage selector, 20

W

warning, 108, 114, 120, 130, 136 warnung, 124 weight, 23 welcome screen, 29

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x axis, 56

Υ

y axis, 56

www.agilent.com

In This Book

© Agilent Technologies 2006

Printed in Germany 10/06