Agilent 3200EA
Electrode Holder
电极支架

Operating Guide
用户手册

Agilent Technologies

Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

5973-1782

Edition

First edition, March 2012

Printed in China

Agilent Technologies, Inc.
412 Ying Lun Road
Shanghai 200131 P.R. China

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Introduction

The 3200EA Electrode Holder

3200EA Electrode Holder accessory can support up to five electrodes or probes during analysis and measurement tasks. The accessory consists of a base and a flexible arm with two double-rockers. This construction provides the ability to keep the electrode holder horizontal at any height. The Electrode Holder accessory includes five (5) 16.5 mm diameter holes for holding electrodes or probes.

3200EA Electrode Holder is suitable for holding electrodes in laboratories of factories and academic organizations. The holder is especially useful to support measurement electrodes (pH, DO, temperature, and so forth) used with Agilent 3200 series meters and a stirrer (for example the Agilent 3200SA Stirrer).

Agilent customer portal

Agilent also provides customized information for the products you own through a customer portal. This web service provides many customizable services as well as information related directly to your Agilent products and orders. Log onto the portal at <http://www.agilent.com/chem>.

What's in the Kit

Agilent provides any necessary special tools needed for installation. The electrode holder shipping case includes:

- 3200EA Electrode holder base
- 3200EA Flexible arm

Installing the 3200EA Electrode Holder

- 1 Open the 3200EA Electrode Holder shipping case.
- 2 Remove the electrode holder base, flexible arm, and other parts.

- 3 Place the electrode holder base onto a stable, clean, dry, and non-slip table surface.
- 4 Insert the flexible arm into the base as shown in [Figure 1](#).


Figure 1 Installing the flexible arm into the base

- 5 Secure the arm to the base by tightening the setscrew at the bottom of the arm. See [Figure 2](#).


Figure 2 Securing the arm to the base

Positioning the Flexible Arm

First adjust the overall height by moving the lower arm, then level the upper arm. Refer to [Figure 3](#).


Figure 3 Flexible arm adjustment controls

To adjust the overall height, loosen the large and small thumbscrews on the lower arm. Move the lower arm to the desired position, then tighten the thumbscrews.

To level the upper arm, loosen the large and small thumbscrews on the upper arm. Move the upper arm to the desired position, then tighten the thumbscrews.

Rotating the flexible arm: The flexible arm normally rotates 360 degrees in the base. If needed, loosen the setscrew. See [Figure 2](#) on page 4.

Installing electrodes/probes: Slide the electrode into the electrode holder. Route the electrode/probe cable under the clips in the upper arm. If needed for measurement, loosen the

thumbscrews in the upper arm and adjust the electrode holder slightly to put the electrode/probe(s) into the sample solution. Tighten the thumbscrews.

CAUTION

Do not completely remove thumbscrews when adjusting the flexible arm's position.

CAUTION

Do not move the flexible arm if the thumbscrews are tight. Always loosen the thumbscrews before moving the flexible arm. Otherwise, the movement can damage the screw or arm threads.

CAUTION

When using the electrode holder with a stirrer, keep the electrode/probe far enough from the stirrer to avoid damage.

For more information and for related supplies, contact your local Agilent distributor or log onto the Agilent website at: <http://www.agilent.com>.

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