**Agilent 3200I
Ion Meter
离子计**

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Agilent Technologies

Notices

© Agilent Technologies, Inc. 2012

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

5973-1778

Edition

First edition, March 2012

Printed in China

Agilent Technologies, Inc.
412 Ying Lun Road
Waigaoqiao Free Trade Zone
Shanghai 200131 P.R. China

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Safety Notices

CAUTION

A CAUTION notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a CAUTION notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

Contents

1 Installation	7
Tools and Components Needed for Installation	8
Installation of the 3200I Ion Meter	9
Installation of electrochemical probes and electrodes	9
Installing the power adapter	12
Installing the Ground Line	13
Installing the Optional Software	14
2 Operation	15
Introduction	16
Where to find information	16
Definition of terms	16
Proper use	17
Features of the 3200I Ion Meter	18
Specifications	19
Physical overview	20
Instrument Control	23
Turn on the Meter	25
Setup	26
Measuring Mode	26
Set Manual Temp	30
Set Condition	30
System Setup	31
Language Select	35
Auto Power Off	35
Set Default	36
Set Ion Mode	37
To adjust mV zero	38

To select a pH buffer group	39
Prepare Electrodes	43
To prepare a pH electrode	43
Calibrate	45
To calibrate a pH electrode	45
To calibrate an ISE electrode	47
To calibrate a blank concentration	48
Measure	50
To measure pH value	50
To measure mV	50
To measure temperature	50
To measure pX	51
To measure ion concentration	52
Managing Data	59
To save data	59
To delete data	59
To output data	59
To view data	59
To Turn Off the 3200I Ion Meter	69
Software Operation	70
To use the optional G4390A Electrochemical Data Collecting software	70
To use the Data Printing software	70
3 Troubleshooting and Maintenance	71
Troubleshooting	72
General Troubleshooting Procedure	73
Check the temperature measurement	73
Restore factory default settings	74
Troubleshoot ion measurement	74

Meter Self-Diagnosis and Messages	76
Severe errors	76
Error warning messages	76
Maintenance	79
Electrode cleaning	79
Electrode storage	79
Consumables and Accessories for the 3200I Ion Meter	80

1 **Installation**

Tools and Components Needed for Installation 8

Installation of the 3200I Ion Meter 9

Installing the Ground Line 13

Installing the Optional Software 14

Use the tools and accessories shipped with 3200I Ion Meter during installation.

Tools and Components Needed for Installation

Agilent provides the tools needed for installation. The following accessories are included in the shipping case.

- Electrode holder (G4389A)
- Electrode
- Power adaptor (5185-8389)

Installation of the 3200I Ion Meter

Open the 3200I Ion Meter shipping case. Remove the meter, electrode holder, and other accessories.

Installation of electrochemical probes and electrodes

- 1 Place the electrode holder near the meter and move the arm into position. Use the thumbscrews shown in [Figure 1](#) to secure the arm in place.


Figure 1 The electrode holder

- 2 Clip the pH/pX electrode and ATC probe into the electrode fixture shown in [Figure 1](#).

1 Installation

3 Route the electrode cable as shown **Figure 2**.


Figure 2 Installed meter and electrodes/probes

- 4 Locate the pH/pX and ATC probe sockets on the back of the meter. Plug the ion selective electrode (ISE) or pH combination electrode connector and ATC probe connector into the sockets on the rear of the meter. See [Figure 3](#).


Figure 3 The back view of the 3200I Ion Meter

Installing the power adapter

A universal power adapter is included with the 3200I. The power adapter operates at 100–240 VAC, 1 A, 50/60 Hz. The output from the power adapter is 9 VDC 1 A.

CAUTION

Do not use any other power adapter with the 3200I.

- 1 The power adapter includes several outlet adapters. Choose the appropriate outlet adapter for the power outlets in your country.
- 2 Slide the outlet adapter on to the two metal prongs on the power adapter. A click will be heard when the plug is properly engaged. See [Figure 4](#).

Figure 4 Assembling the power adapter

- 3 Connect the power adapter cord to the power socket on the back of the meter. See [Figure 3](#) on page 11.


Installing the Ground Line

Agilent provides a ground line with the meter, but normally you do not need to install it for use. If other devices, for example a constant temperature bath, causes electrical interference and unstable readings, use the ground wire. Connect the ground wire terminal to the meter and connect the other end to the interference source, for example the cover of the constant temperature bath.

Installing the Optional Software

If purchased, install the optional G4390A Electrochemical Data Collecting Software now. The software provides communications between the meter and a computer. Connect the computer to the meter with a USB cable. See the G4390A software documentation for more information.

If not using the G4390A software, Agilent provides downloadable data printing software on the Agilent Customer Portal (see “[Agilent customer portal](#)” on page 16). To use this software, download it from the portal and install it. Then connect the meter to the PC using the USB cable. Refer to the data printing software documentation for more information.

2 Operation

- Introduction 16
- Instrument Control 23
- Turn on the Meter 25
- Setup 26
- Prepare Electrodes 43
- Calibrate 45
- Measure 50
- Managing Data 59
- To Turn Off the 3200I Ion Meter 69
- Software Operation 70

Introduction

Where to find information

This document describes how to install, operate, maintain, and troubleshoot the 3200I Ion Meter.

Before operating your meter, be sure to read the meter installation and operation information.

Agilent customer portal

Agilent also provides customized information for the products you own through a customer portal. This web service provides many customizable services as well as information related directly to your Agilent products and orders. Log onto the portal at <http://www.agilent.com/chem>.

Definition of terms

pH Slope

Refers to the mV change when the pH has been changed each unit and is represented by mV/pH or percent (%) of the theoretical slope.

pH E₀

Also called *Zero Potential* and refers to the mV value at pH7.

One-point calibration

Calibration with one pH buffer.

Multi-point pH calibration

Calibration with two or more pH buffers.

Proper use

To avoid common safety issues:

- Ground the meter using the connection on the back.
- The meter can be used continuously for a long time. After each measurement, soak the electrode in distilled water. If an electrode is out of use for more than 6 hours, rinse and store it in a protective container.
- Improper positioning of the electrode can result in abnormal measurements. When immersing the electrode in a sample solution, place the electrode in a location where the solution can flow freely around the electrode.

CAUTION

The ion meter is an analytical instrument with high accuracy. To protect high resistance components from damage, place the short circuit plug (G4383-40000) shipped with the meter on the electrode connector when the meter is not connected with an electrode.

CAUTION

Do not expose to corrosive gas. Keep the sockets on the back of the meter clean and dry. Do not allow contact with acid, alkaline, or salt solutions.

CAUTION

Only use the power adapter included with the meter.

CAUTION

When the short-circuit plug is not being used, it must be kept dry and clean. The plug could be damaged by high humidity. The meter could be damaged or the measurement affected if a damaged plug is used.

CAUTION

Do not use where nearby vibrations will affect the performance. Do not use if corrosive gas is in the air. Do not use near strong electromagnetic fields.

Features of the 3200I Ion Meter

The 3200I Ion Meter is a state-of-the-art and customer-friendly benchtop analytical instrument, designed for the measurement of pH/pX, ion concentration, potential, and temperature in aqueous solutions with high accuracy.

Features include:

- Measurement of pH/pX, ion concentration, potential, and temperature.
- Automatic pH buffer recognition for buffers prepared according to NIST, DIN, and GB.
- Multi-point calibration up to five calibration points.
- Ion measurement modes for H^+ , Ag^+ , Na^+ , K^+ , NH_4^+ , Cl^- , F^- , NO_3^- , BF_4^- , CN^- , Cu^{2+} , Pb^{2+} and Ca^{2+} .
- Custom defined ion measurement modes for other ions
- A clear dot-matrix LCD display
- GLP norm. Stores 200 pH data sets, 100 pX data sets, and 100 ion concentration data sets.
- The ability to view, output, and delete stored data
- Three measuring modes to meet the needs of various users:
 - Continuous Mode
 - Timed Reading (Timed Reading Mode)
 - Auto-Lock Mode (Auto-Lock Reading Mode)
- USB PC connectivity with available communication software
- Power-off protection. When the meter is manually or automatically shut off, the stored measuring data, calibration data, and setting parameters will not be lost.
- A back-lit design that can be used in a dark environment
- A durable key pad
- An Auto Power Off function
- A self-diagnosis mode for troubleshooting

Specifications

Measuring range	
pH	-2.000–20.000 pH
pX	0.000–14.000 pX
mV	-1999.9–1999.9 mV
Unit	pX, mol/L, ppm, %, mg/L, µg/L
Temperature	-5.0–110.0 °C
Resolution	
pH/pX	0.1/0.01/0.001 pH/pX
mV	0.1 mV
Temperature	0.1 °C
Accuracy	
pH/pX	±0.002 pH/pX
mV	±0.03% FS
Concentration	±0.3%
Temperature	±0.1 °C
Normal working conditions	
Environmental temperature	0–40 degrees
Relative humidity	≤ 85%
Power supply	Power adaptor, 5185-8389
Power input	100–240 VAC, 1 A
Power output	9 VDC, 1 A
Size	(length × width × height, mm): 190 × 190 × 105
Weight	about 1 kg

Physical overview

The Agilent 3200I Ion Meter is composed of a meter and an electrode system. The electrode system contains a pH combination electrode, an ion selective electrode, and an ATC probe held in place by the electrode holder. See [Figure 5](#), [Figure 6](#), and [Figure 7](#).


Figure 5 The meter, electrode holder, and probes.

Figure 6 The front view of the 3200I Ion Meter



Figure 7 The back view of the Agilent 3200I Ion Meter

The Display

The display shows the working condition and current setting of the 3200I Ion Meter. The initial state is shown in [Figure 8](#). The left of the screen shows current system time and the right shows current measuring mode, parameters, and the most recent calibration data.

[Figure 8](#) The Initial State display

The Keypad

The keypad has 15 keys. All keys have double functions except for [**Enter**] and [**Cancel**]. Functionality automatically changes based on context. Key functions are listed in the [Table 1](#) on page 23:

[Figure 9](#) Keypad

Instrument Control

The 3200I Ion Meter is usually controlled by the keypad. Most keys have two functions. Functionality automatically changes based on context. Key functions are listed in the [Table 1](#):

Table 1 Keypad functions






Button	Function 1	Function 2
	Type the number 1 .	Output data when viewing or calibrating stored data.
	Type the number 2 .	Move the selection cursor upward when selecting.
	Type the number 3 .	Store measurement data.
	Type the number 4 .	Move the selection cursor left when selecting.
	Type the number 5 .	Opens the Setup menu from the initial screen and also is used as a general “select” button.
	Type the number 6 .	Move the selection cursor right when selecting.
	Type the number 7 .	View stored or calibrated data.
	Type the number 8 .	Move the selection cursor down when selecting.
	Type the number 9 .	Switch the display window or parameter when measuring.

Table 1 Keypad functions (continued)

Button	Function 1	Function 2
	Type the number 0 .	Begin measurement.
	Type a decimal point.	Calibrate electrode.
	Type a negative number.	Delete the data being viewed.
	Meter power switch.	

Alternately, the 3200I Ion Meter can be controlled through a computer.

Uses for the 3200I Ion Meter

The Ion Meter is suitable for measuring pH, pX and ion concentration of aqueous solutions.

Turn on the Meter

- 1 Press [On/Off] to turn on the meter. When the meter is powered on, the display shows the Agilent name, the meter model, and other information.
- 2 After the self-check is completed, the display shows the initial state. The display shows the measuring mode, pH parameter, and last calibration data. See [Figure 10](#).


Figure 10 The initial state

Setup

Measuring Mode

This meter supports three measuring modes including **Continuous Mode**, **Timed Reading**, and **Auto-Lock Mode**.

To set the measuring mode:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#).


Figure 11 The setup screen

- 2 Use the directional arrow keys to navigate to **Measuring Mode**.
- 3 Press [**Enter**] to select **Measuring Mode**. If you need to go back after selecting a menu item, press [**Cancel**] to exit.

- 4** The **Measuring Mode** screen displays. The left column displays the available measurement parameters pH, pX (ion conductivity), and CONC (ion concentration).

The right column displays the measuring modes for the selected (highlighted) parameter. For most parameters, these are: Continuous Mode, Timed Reading, and Auto-Lock Mode. For ion concentration, these are: Direct Reading, STD Addition, Sample Addition, and GRAN Method.

See [Figure 12](#). Check marks indicate the current selections in use.


Figure 12 Select a measuring mode

- 5** Use the arrow keys to highlight the desired parameter and press [**Setup**] to select it. A check appears next to the selected parameter. During measurement, you can still view, save and print the other two parameters. See “[To view other parameters during measurement](#)” on page 68.
- 6** Use the arrow keys to highlight the desired measurement mode for the selected parameter.
- 7** Press [**Setup**] to enable the mode. See “[Descriptions of the modes](#)” on page 28 for more information.
- 8** Repeat [step 5](#) through [step 7](#) to select or clear all desired parameters and modes.
- 9** After setting all desired parameters and modes: Press [**Enter**] to save the new setup, exit from setup mode and return to the initial state, or Press [**Cancel**] to exit setup mode without saving changes and return to the initial state.

The new parameters and measuring mode are displayed as shown in [Figure 13](#).


Figure 13 Example measuring mode displayed

Descriptions of the modes

Continuous Mode The meter will simply continue measuring, calculating, and displaying the data. It is up to the user to manually save data or end the measurement.

This is the most frequently used measuring mode. When measurement begins, the meter continuously measures, calculates, and displays the results. You can view the calibrated parameter, calibrate the electrode, save, or print results at any time during measurement. To end a measurement run, press **[Cancel]** followed by **[Enter]**.

Timed Reading The meter will automatically save measurement data periodically while measuring.

When you select Timed Reading, specify a time interval from 1 to 99 minutes. The default time interval is 10 minutes.

When measurement begins, the meter calculates and displays results as it normally does. When the set interval time elapses, the meter automatically stores the data and continues measuring. Data is constantly recorded at the specified time interval.

If the USB interface is connected to a PC, the meter will automatically print the data.

Read about Auto Delete before using this mode. See [“Descriptions of the System Setup settings”](#) on page 33.

To end a measurement run, press **[Cancel]** followed by **[Enter]**.

Auto-Lock Mode The meter will stop measuring when recorded values stabilize enough to fall within a set precision and time.

Before performing a measurement in Auto-Lock reading mode, set the Auto-Lock conditions (see “[Set Condition](#)” on page 30). When the measurement begins, the meter automatically measures, calculates, and displays the results. Once the measurements meet the preset Auto-Lock condition, the measurement finishes and the screen displays the final readings.

You can view the calibrated parameter, calibrate the electrode, save, or print results at any time during measurement. After measurement, you can save and print the results. Press [**Cancel**] to exit measuring mode or press [**Measure**] to begin the next measurement.

Direct Reading Mode (Direct Reading) This is the most commonly used ion concentration measuring method. When measurement begins, the meter reads, calculates and displays the potential (mV) value. When the reading becomes stable, press [**Enter**] and the meter calculates the current concentration value. See “[Direct Reading mode](#)” on page 52 for details.

Standard Addition Mode (STD Addition) In this mode, the user adds standard solution to the sample. Measure the potential before and after the addition of standard solution. The meter then calculates sample concentration. See “[To clear the blank concentration calibration](#)” on page 49 for details.

Sample Addition Mode (Sample Addition) This method is similar to the standard addition method, except that the user adds sample solution to the standard solution. Measure potential variation before and after adding sample solution. See “[Sample Addition Mode](#)” on page 57 for details.

GRAN Method The GRAN Method uses multiple additions of standard solution to the sample solution. The user adds a fixed amount of standard solution to the sample repeatedly. Measure the potential after each addition. See “[GRAN method](#)” on page 58 for details.

Set Manual Temp

If the meter is connected to an ATC probe, the probe will supply a temperature value. If an ATC probe is not connected to the meter, the user must enter the sample solution temperature manually.

To manually set the sample solution temperature:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Set Manual Temp**.
- 3 Press [**Enter**] to select **Set Manual Temp** or press [**Cancel**] to exit.
- 4 The Set Manual Temp screen appears on the display. See [Figure 14](#).


Figure 14 Set Manual Temp

- 5 With the Manual Temp value highlighted, press [**Setup**] to edit the value.
- 6 Use the numeric keys to enter a temperature.
- 7 Press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

Set Condition

Use **Set Condition** to specify Auto-Lock parameters when measuring in the Auto-Lock measuring mode. An Auto-Lock measurement ends when all measured parameters meet the conditions set on this screen. For example, if **Measure pH** is set to

0.010pH, when the measured pH holds stable within a 0.010pH tolerance for the set Auto-Lock time interval, the measurement ends and the screen displays the final reading.

To set the conditions:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Set Condition**.
- 3 Press [**Enter**] to select **Set Condition** or press [**Cancel**] to exit.
- 4 The **Set Condition** screen appears on the display. See [Figure 15](#).


Figure 15 The Set Condition screen

- 5 Use the arrow keys to select a condition to edit. Press [**Setup**] to edit the value.
- 6 Enter the desired value. The Auto-Lock time has a range of 1–200 seconds.
- 7 Press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

System Setup

The System Setup screen is where the user sets the time on the system clock, sets the calibration interval, and enters the information for the Good Laboratory Practices (GLP) standards.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **System Setup**.
- 3 Press [**Enter**] to select **System Setup** or press [**Cancel**] to exit.

- 4 The **System Setup** screen appears on the display. See [Figure 16](#).


Figure 16 The System Setup screen

- 5 Use the arrow keys to select a value to edit. Press [**Setup**] to edit the value.
- 6 Enter the desired value. See “[Descriptions of the System Setup settings](#)” on page 33 for more information.

- 7** If needed, edit the system date and time. To modify the date and time:
- Use the arrow keys to highlight the Time field shown in [Figure 16](#) on page 32 and press [**Setup**].
 - The **Set Date & Time** screen opens. This screen displays the current year, month, day, hour, minute, and second (see [Figure 17](#)). Press the arrow keys to highlight the proper cell, then press [**Setup**] to select it.


Figure 17 The Set Date & Time screen

- Enter the correct value using the numeric key pad and press [**Enter**].
- Once all changes are made, press [**Enter**] to save the settings and return to the **System Setup** screen.
- Once all of the system settings are made, press [**Enter**] to save the new setup, exit, and return to the initial state. Press [**Cancel**] to exit without saving changes.

Descriptions of the System Setup settings

Calib Interval is the number of hours between the calibration messages that the Meter periodically displays as a reminder to recalibrate the electrode. The meter begins counting the interval time from when the previous calibration ended. When the calibration interval has elapsed, the meter displays a popup window to remind the user to recalibrate the electrode. A **Calib Interval** value of **0** disables the reminder.

Operator No is a three-digit number from 000–200 used to identify the person operating the meter. The **Operator No** is recorded every time the meter saves data.

When stored data reaches the meter's memory limit **Auto Delete** allows you to automatically overwrite old data with new data. To automatically overwrite old data, select the **Auto Del** checkbox. If **Auto Delete** is not enabled, the meter does not save new data when its memory is full. For example, the meter permits up to 200 sets of data. When you want to store the 201st set of data, if **Auto Delete** is on, meter will delete the first set of data and store the 201st data. If **Auto Delete** is disabled, the new data will not be stored.

CAUTION

If **Auto Delete** is disabled, data can be lost.

Language Select

This meter supports Chinese and English languages.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Language Select**.
- 3 Press [**Enter**] to select **Language Select** or press [**Cancel**] to exit.
- 4 The **Language Select** screen appears on the display. Use the arrow keys to choose either **Chinese** or **English** and press [**Enter**].

Auto Power Off

The meter can power off automatically after a set time period.

To set the **Auto Power Off** time:

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Auto Power Off**.
- 3 Press [**Enter**] to select **Auto Power Off** or press [**Cancel**] to exit.
- 4 The **Auto Power Off** screen appears. Enter a power off time from 10–480 minutes. Enter a time of 0 to disable **Auto Power Off**.
- 5 Press [**Enter**].

When the **Auto Power Off** time elapses after the start of a measurement run, the meter shuts off. If the meter is connected to the optional G4390A Electrochemical Data Collecting Software, **Auto Power Off** is disabled.

Set Default

Use **Set Default** to reset the meter to its original factory settings.

- 1 From the initial state, press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 2 Use the directional arrow keys to navigate to **Set Default**.
- 3 Press [**Enter**] to select **Set Default** or press [**Cancel**] to exit.
- 4 The **Set Default** screen appears. Press [**Enter**] to restore parameters to their original settings.

Table 2 Default values for meter parameters

Parameter	Default value
mV Zero	Clears the value. Re-zero before use.
Buffer Group calibrations (NIST standards, DIN standards, GB standards)	Clears all values. The default electrode slope is set to 100% and the E_0 is set to 0.0 mV. Select the desired buffer group before taking measurements.
Other calibration data	Resets to default values
Ion concentration calibration units	mol/L
Ion concentration measuring units	mol/L
current Ion Measuring Method	Na ⁺
Manual temperature	25.0 °C
Measuring mode	Continuous
Timed Reading Mode interval	10 minutes
Calib Interval	message disabled (0 hours)
Operator No	000
Auto Delete	Enabled
Auto-Lock Set Condition	0.1 pH
pH resolution	0.001 pH
Auto-Lock time	5 s
Auto Power Off	Disabled (0 minutes)

Set Ion Mode

The meter provides 10 preset ion modes, used to measure different ion concentrations. To set the meter to measure concentration for a listed ion:

- 1 If needed, install the appropriate electrode for the desired ion.
- 2 Press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 3 Select **Set Ion Mode**.

Press [**Enter**] and the meter will display as shown in [Figure 18](#).


Figure 18 Set ion mode

- 4 Use the scroll keys to highlight the desired ion. The currently-selected ion name and molecular weight appear at the bottom of the screen.

Always select the correct ion mode before measurement. An improper choice leads to inaccurate measurements. For example, to measure sodium ion concentration, select **Na⁺**.

- 5 Press [**Enter**] to select the highlighted ion, or press [**Cancel**] to close this screen and return to the initial state.

Regular ions

These ion modes are supplied with the meter: H⁺, Ag⁺, Na⁺, K⁺, NH₄⁺, Cl⁻, F⁻, NO₃⁻, BF₄⁻, CN⁻, Cu₂⁺, Pb₂⁺ and Ca₂⁺. Since H⁺ is always allowed, it is not displayed in the ion mode.

Custom ions

You can create up to five custom ions (custom ion modes), named **Cus00** to **Cus04**. A custom ion mode requires an ISE electrode capable of measuring the ion.

To create a new custom ion mode:

- 1 Press [**Setup**] (at the Set Ion Mode screen) to establish a new one. (The meter supplies the name for the mode.)
- 2 Press [**Enter**] to create the new Custom ion mode.
- 3 Press [**Setup**], then select **Amend Ion Mode**.
- 4 Set the ion valence and molecular weight for this custom ion.
- 5 Press [**Enter**] to confirm.
- 6 To delete a custom ion mode, navigate to the Set ion mode screen, then press [**Setup**]. Scroll to Delete Ion Mode then press [**Enter**].

CAUTION

When deleting custom ion modes, all corresponding stored data will be deleted at the same time.

To adjust mV zero

Although the meter compensates for many factors which may affect measurement (such as temperature), it still cannot guarantee zero drift. For accurate measurement, perform **Adjust mV Zero** before measurement.

To check for drift in pH measurement:

- 1 Turn off the meter.
- 2 Install the short circuit plug (G4383-40000) into the pH/pX socket.
- 3 Turn on the meter and let stabilize for 30 seconds.
- 4 Press [**Setup**]. The setup screen appears. See [Figure 11](#) on page 26.
- 5 Select **Adjust mV Zero**, and press [**Enter**]. The meter prompts **Adjust mV zero?**.

- 6 Press [Enter] to calibrate a new mV zero. If the potential (mV) value deviates far from zero, double check the short circuit plug connection. Ensure the meter is in zero mV state to avoid an incorrect adjustment.
- 7 To abort the adjustment, press [Cancel]. The meter will return to the measuring state.

To select a pH buffer group

- 1 View the last pH calibrated data. See “To view calibration data ” on page 47.
- 2 Press [Setup] to enter Group Setup as shown in **Figure 19**. **Figure 19** indicates that the current Buffer Group Setup meets the DIN standard (DIN is checked). The upper area displays the 3 buffer groups and the bottom area displays the buffers included in the selected buffer group.

Figure 19 Group setup

- 3 Use arrow keys to highlight the desired buffer group.
- 4 Press [Setup] to select the highlighted buffer group. (Only one group can be active at a time.)

To select a buffer in a group

After selecting the desired buffer group, select the individual buffers that best correspond to the sample.

To prevent selecting buffers with overlapping pH ranges, select only buffers necessary for measurement and clear selection of other buffers.

- 1 Verify the current buffer group.
- 2 Use arrow keys to highlight the target buffer group and press [**Setup**]. The meter shows the buffers available in the group. See [Figure 20](#).


Figure 20 Selecting a buffer for use or calibration (example NIST buffer group shown)

- Icons with \checkmark indicate that a selected buffer.
 - Icons without \checkmark indicate an available buffer.
 - The bottom part of the display shows the pH value and the range of the highlighted buffer.
- 3 Use the arrow keys to select (highlight) a buffer.
 - 4 Press [**Setup**] to select or deselect this buffer.

For example, to select the 4.008 pH buffer, use arrow keys to highlight the 4.008 pH icon. Press [**Setup**] to select the icon.

Calibrate the meter to all selected buffers.

About pH buffer groups

The meter provides auto-recognition of several buffer groups prepared according to either the NIST standard, the deionizedN (DIN) standard, or the GB standard. Each group contains several buffers. Using a buffer group simplifies ion calibration. Before measurement, select the buffer solutions most applicable to the sample solution, then calibrate the electrode to those buffers. The meter uses these calibrations to provide accurate results.

Table 3 Buffer groups

Buffer pH	NIST group	DIN group	GB group
1.677	●		
1.680		●	●
3.557		●	
3.559			●
3.775		●	
4.003			●
4.008	●	●	
6.864	●		●
6.865		●	
7.000	●	●	
7.409			●
7.416	●	●	
9.182			●
9.184		●	
10.014	●	●	
12.454		●	
12.460			●
12.469	●		

2 Operation

For each Buffer Group setup, you may select up to five buffers. Since the pH range of a buffer may overlap with the pH range of another buffer in the same buffer group, the meter prevents selecting two or more buffers with overlapping pH ranges.

Prepare Electrodes

To prepare a pH electrode

A new electrode should be stored wet in the storage solution and should be in working condition when received from the factory.

- 1 Visually check for any mechanical damage. If the storage bottle is dry, soak the electrode in reference filling solution for at least 2 hours before use.
- 2 Take the storage bottle off and keep it for future use.
- 3 If there is electrolyte solidified on the measuring tip, rinse it off with distilled or deionized water.
- 4 Unplug the filling hole.
- 5 Siphon away the reference solution and add new reference solution. The level of the reference solution must be at least 20 mm higher than the level of sample. The liquid junction must be immersed completely in solution.
- 6 To guarantee mobility of the reference solution, keep the filling hole open during measurement.
- 7 Hold the electrode measuring tip downwards and swing it several times to remove air bubbles near the sensitive glass bulb.
- 8 Connect the electrode to the meter. Put the electrode in the first buffer solution with the measuring tip downwards.

To prepare an ISE electrode

- 1 Different ISE electrodes need different preparation. Ensure the unused electrode was stored in accordance with the requirements of “[Electrode storage](#)” on page 79. If not, soak the electrode in storage solution for at least 2 hours.
- 2 Take the electrode tip out of its protective sleeve or storage bottle. Preserve them for future use.
- 3 If there is electrolyte solidified on the outside of the electrode tip, rinse it off with distilled water.
- 4 Hold the electrode tip downwards and swing it several times to remove any air bubbles from the electrode.
- 5 If the electrode is refillable, open the filling hole and add filling solution. Fill to at least 20 mm higher than the sample solution level.
- 6 To guarantee filling solution mobility during measurement, leave the filling hole open.

Calibrate

To calibrate a pH electrode

Before each measurement, we recommend recalibrating the electrode to the series of buffers chosen for the sample. After calibration, any previous calibration data is overwritten.

The meter can perform calibration manually (**Manual config**) or using its auto-recognition feature (**Auto-config**).

- In **Auto-config**, the meter automatically compares the current pH measurement against the selected buffers and assigns the calibration to the correct buffer. If the meter cannot recognize a buffer, it returns a **Calibration error**.
 - In **Manual config**, the meter prompts you for the correct buffer pH at the selected temperature.
- 1 Prepare 1 to 5 buffers. They can be purchased or can be prepared by you. Place them in a constant temperature environment to equilibrate.
 - 2 View the current calibration. See [Figure 21](#). The upper area shows the current pH (in this example the slope is set to 100.00%), potential and temperature value. The bottom area shows the current calibration result.


Figure 21 Calibrating the pH electrode

- 3 Confirm the Buffer Group setup. See “To select a pH buffer group” on page 39 and “To select a buffer in a group” on page 40.
 - If you try to calibrate to a buffer that is not selected, it can cause a calibration error or a bad calibration.

- 4 Press [Setup] and scroll to **Set Config Type** to select the recognition method.


Figure 22 Calibration setup

5 Select Auto-config or Manual config.

- For regular buffers, use the Auto-Recognition method. Select **Auto-config**.
 - If using a custom pH buffer, select **Manual config**.
 - If using both regular buffers and custom buffers, use both recognition methods. You can calibrate to some buffers, change the buffer selections, then calibrate to another buffer. See “[Calibrating custom buffers and buffers with overlapping pH values](#)” on page 47.
- 6 Clean the pH electrode, reference electrode, and ATC probe. Place them in the calibration buffer. See “[To prepare a pH electrode](#)” on page 43.
- 7 When the reading becomes stable, press [**Enter**]. The meter will display **Storing.....** and save the calibration data.
- 8 After a few seconds, the meter will display a **Continue?** prompt. To calibrate other buffers, select the next buffer to continue. Repeat [step 6](#) and [step 7](#).

At any time during calibration, press [**Cancel**] to end calibration.

The meter supports up to five calibration points. When the fifth buffer is calibrated, the meter prompts to end calibration.

Calibrating custom buffers and buffers with overlapping pH values

For buffers such as 6.864 pH and 7.000 pH which have overlapping pH ranges, try the following:

- When you calibrate the 6.864 pH buffer, set the buffer group to contain only 6.864 pH, then calibrate it. Reset the buffer group to contain only 7.000 pH, then calibrate it.
- Use **Manual config** and manually input the pH of each selected buffer at the selected temperature.

To calibrate an ISE electrode

- 1 Prepare the electrode for use as described in “To prepare an ISE electrode” on page 44.
- 2 Install the electrode. See “Installation of electrochemical probes and electrodes” on page 9.
- 3 Set **pX** as the measurement parameter and set **Continuous Mode**. See “Measuring Mode” on page 26.
- 4 Press [**Measure**] to start pX measurement.
- 5 Press [**Calibrate**], then scroll to and select **Calib pX EC**.
- 6 Press [**Enter**] to begin the electrode slope calibration. See Figure 23. The upper area shows the current pX (the slope is set at 100.00 %), potential and temperature value. The bottom area shows the current calibration result.


Figure 23 pX calibration

- 7 Prepare one to five standard solutions.
- 8 Equilibrate the solutions at the required temperature.

- 9 Clean the appropriate ISE electrode or reference electrode and ATC probe, and place them in the standard solutions to be calibrated. See “To prepare an ISE electrode” on page 44.
- 10 Press [Setup] and select **Set STD Value**. Manually input the STD pX value of the selected standard solution at the selected temperature.
- 11 When the reading becomes stable, press [Enter]. The meter displays **Storing ...** and saves the calibration data.
- 12 After a few seconds, the meter displays a **Continue?** prompt.

Figure 24 Storing calibration

- 13 To calibrate another standard solution, select your new standard solution and press [Enter] to repeat the procedure. Continue until all standard solutions are calibrated.

At any time during calibration, press [Cancel] to end the process.

To calibrate a blank concentration

Both the direct reading mode and the standard addition mode have a blank calibration. Performing a blank calibration is the same in either mode.

- 1 Create a blank solution with a similar chemical composition to the sample solution (except for the concentration of the ion measured).
- 2 Set either **Direct Reading** or **STD Addition**. See “Measuring Mode” on page 26.
- 3 Press [Setup] and Scroll to **CONC Unit**. Select the correct units.
- 4 Scroll to **Blank CONC** and press [Setup].

- 5 Put the ISE and the ATC probe in the blank solution. See Figure 25.


Figure 25 Blank Concentration Calibration in Direct Reading mode

- 6 When the reading becomes stable, press [**Enter**] and the meter calculates and saves the blank concentration value.
- 7 Press [**Cancel**] to return to measuring ion concentration in direct reading mode.

To clear the blank concentration calibration

To clear the last blank concentration value:

- 1 Set either **Direct Reading** or **STD Addition**. See “[Measuring Mode](#)” on page 26.
- 2 Press [**Setup**] and select **Clear Blank CONC**.

Measure

To measure pH value

- 1 From the initial state display, press [**Setup**].
- 2 Select the pH parameter and a measuring mode. See “[Measuring Mode](#)” on page 26.
- 3 If the electrode has not been recently calibrated, calibrate it. See “[Calibrate](#)” on page 45.
- 4 Press [**Measure**] to begin measurement. [Figure 26](#). The upper area of the screen shows the current measuring mode and system time. The central area of the screen shows present pH, potential, and temperature.


Figure 26 Enter Continuous Mode

- 5 When measurement begins, adjust mV zero. See “[To adjust mV zero](#)” on page 38.
- 6 During measurement, you can recalibrate the electrode, set parameters and select the pH display resolution.
- 7 After measurement, press [**Save**] to save data, [**Output**] to output data, or [**Cancel**] to end measurement.

To measure mV

The current pH and mV values are always shown when measuring.

To measure temperature

The current temperature is always shown when measuring.

To measure pX

With a recently-calibrated ISE installed, you can begin measurements.

- 1 From the initial display, ensure the measured parameter is pX. If not, select pX measurement as described in “[Measuring Mode](#)” on page 26.
- 2 Check the current ion mode. See “[Set Ion Mode](#)” on page 37.
- 3 Adjust mV zero. See “[To adjust mV zero](#)” on page 38.
- 4 If the ISE is not calibrated, prepare and calibrate the ISE before use. See “[To calibrate an ISE electrode](#)” on page 47.
- 5 Press **[Measure]**. See [Figure 27](#). During measurement, you may recalibrate the electrode and set parameters.


Figure 27 pX measurement

- 6 After measurement, you can press **[Save]** to save measurement data, **[Output]** to print data, or **[Cancel]** to end measurement.

To measure ion concentration

There are four concentration measurement modes:

- Direct Reading of Concentration Mode (**Direct Reading**)
- Standard Addition Mode (**STD Addition**)
- Sample Addition Mode (**Sample Addition**)
- GRAN Method (**GRAN Mode**)

Direct Reading mode

The Direct Reading mode uses the following Nernst formula to calculate the concentration:

$$E_x = E_0 + S[\log(C_x + C_b)]$$

In which:

E_x = Equilibrium potential of sample

E_0 = Zero potential, potential for ion activity = 1

S = Electrode slope

C_x = Concentration value of sample

C_b = Blank concentration value

Obtain slope and zero potential through slope calibration. Calibrate the blank concentration also. Then measure concentration of the sample directly.

To measure concentration in Direct Reading mode:

- 1 Select **CONC** as the parameter and select **Direct Reading** as the mode. See “[Measuring Mode](#)” on page 26.
- 2 Select the ion mode that matches the ISE (ion selective electrode) type. For instance, to measure Ag^+ concentration, select the **Ag^+** mode. See “[Set Ion Mode](#)” on page 37.
- 3 From the initial screen, press [**Measure**]. See [Figure 28](#). During measurement, you can view calibrated data, calibrate the electrode, adjust mV zero, calibrate blank concentration, and reset blank concentration.


Figure 28 Direct Reading mode

- The upper left area shows current system time, potential, temperature, and pX value. The bottom area shows current measuring result and blank concentration value.
- 4** If necessary, perform a blank concentration correction. See “To calibrate a blank concentration” on page 48.
- 5** Rinse the ISE thoroughly and put it in sample solution. The meter displays current measurement data.
- 6** When the reading becomes stable, press [**Enter**] and the meter calculates measurement result as shown in Figure 29.


Figure 29 Calculated result

- 7** After measurement, you can press [**Save**] to save measurement data, [**Measure**] to perform another ion concentration measurement, [**Setup**] to select the proper concentration unit, [**Output**] to print data, or [**Cancel**] to end measurement.

Standard Addition Mode (STD Addition)

Standard Addition Mode is also called “known addition mode”. First, measure the potential of the sample solution after the reading becomes stable. Then add standard solution of known concentration to the sample. Measure the potential again after the reading becomes stable. Calculate the sample concentration through the change in potential using the formula:

$$C_x = \frac{\rho \times C_s}{(1 + \rho) \times 10^{(E_2 - E_1)/S} - 1} + \frac{\rho \times C_b}{(1 + \rho) \times 10^{(E_{b2} - E_{b1})/S} - 1}$$

In which:

C_x = Concentration value of sample to be measured

C_s = Concentration value of STD solution (add solution)

S = Slope of electrode

C_b = Blank STD concentration value

E_1 = Potential value of the system before adding solution

E_2 = Potential value of the system after adding solution

ρ = Ratio of added volume of the standard solution (V_s) to volume of the sample to be measured (V_x)

E_{b1} = Potential value of the system before adding solution when calibrating blank

E_{b2} = Potential value of the system after adding solution when calibrating blank

To measure concentration in STD Addition mode:

- 1 Select **CONC** as the parameter and select **STD Addition** as the mode. See “[Measuring Mode](#)” on page 26.
- 2 Select the ion mode. See “[Set Ion Mode](#)” on page 37.
- 3 Press **[Measure]**. The meter enters standard addition mode as shown in [Figure 30](#).

Figure 30 Standard Addition mode

- 4 Input the required information. Use the arrow keys to select each parameter. Press [**Setup**] to modify the parameter or to calibrate the blank concentration (if **Blank CONC** is selected).

Adding Volume Addition volume. The standard solution volume to be added.

Vol before Add The sample volume before adding any standard.

STD CONC Value Standard solution concentration

STD CONC Unit Standard solution concentration unit

Blank CONC Blank concentration value. See “[To calibrate a blank concentration](#)” on page 48.

- 5 After setup, rinse the ISE electrode thoroughly and put it in the sample solution.
- 6 Press [**Enter**] to perform the standard addition measurement. The meter displays current potential, temperature and pX value. See [Figure 31](#).


Figure 31 Before adding

- 7 When the reading becomes stable, press [Enter] to save this potential. The meter displays **STD Addition**.
- 8 Add the correct volume of standard solution and mix.
- 9 When the reading becomes stable, press [Enter]. The meter prompts to end the process and shows the calculated sample concentration. See **Figure 32**.


Figure 32 Final screen

- 10 After measurement, you can press [**Save**] to save measurement data, [**Output**] to print data, or [**Cancel**] to end measurement.

Sample Addition Mode

This mode is similar to the standard addition mode, except that the solution is added to the standard. The meter calculates the concentration from the formula:

$$C_x = C_s \times [(1 + \rho) \times 10^{(E_2 - E_1)/S} - \rho]$$

In which:

C_x = Concentration value of sample to be measured (add solution)

C_s = Concentration value of STD solution

ρ = Volume of STD solution (V_s)/volume of sample to be measured (V_x)

E_1 = Potential value of the system before adding measured sample

E_2 = Potential value of the system after adding measured sample

S = Slope of electrode

To measure concentration in Sample Addition mode:

- 1 Select **CONC** as the parameter and select **Sample Addition** as the mode. See “[Measuring Mode](#)” on page 26.
- 2 Select the ion mode. See “[Set Ion Mode](#)” on page 37.
- 3 Press [**Measure**]. The meter enters sample addition mode as shown in [Figure 33](#).

Figure 33 Sample Addition

- 4 Input the required information. Use the arrow keys to select each parameter. Press [**Setup**] to modify the parameter or to the calibrate blank concentration (if **Blank CONC** is selected).
- 5 After setup, rinse the ISE electrode thoroughly and put it in the sample solution.
- 6 Press [**Enter**] to perform the sample addition measurement. The meter displays current potential, temperature and pX value.
- 7 When the reading becomes stable, press [**Enter**] to save this potential. The meter displays **Sample Addition**.
- 8 Add the correct volume of sample solution and mix.
- 9 When the reading becomes stable, press [**Enter**]. The meter prompts to end the process and shows the calculated sample concentration.

GRAN method

The meter supports the GRAN method for measuring samples with low concentration. Solve for C_x , using this equation:

$$(V_s + V_x) \times 10^{(E)/S} = 10^{(E_0)/S} \times (C_x \cdot V_x) + 10^{(E_0)/S} \times (C_s \cdot V_s)$$

Operation is similar to standard addition mode.

- 1 Select **CONC** as the parameter and select **GRAN Mode** as the mode. See “[Measuring Mode](#)” on page 26.
- 2 Select the ion mode. See “[Set Ion Mode](#)” on page 37.
- 3 Press [**Measure**]. The meter enters GRAN mode.
- 4 Input the concentration (C_s) and volume (V_s) of the STD solution and the volume (V_x) of the sample.
- 5 Measure the electrode potential value of the sample after the standard solution is added.
- 6 Repeat the measurement 3 to 8 times. The meter calculates the concentration of the sample.

Managing Data

To save data

The Meter can save 200 sets of pH data. The procedure to save data varies with different measuring modes.

- In continuous mode and Auto-Lock reading mode, press [**Save**] to save data when readings become stable.
- In timed reading mode, you can still save data manually, but the meter automatically saves measuring data at a fixed periodic time interval. See “[Measuring Mode](#)” on page 26.

To delete data

You can delete individual data entries or delete all data at once.

- 1 In either the initial screen or while measuring, press [**View**].
- 2 Select the data and press [**Delete**].

To output data

There are two ways to output current measurements, calibration data, and saved data.

If using the optional G4390A the data collecting software with the meter, the software can collect the data, print it, and so forth.

If not using the optional G4390A software, install the optional data printing software. See “[Installing the Optional Software](#)” on page 14.

To view data

The Agilent 3200I Ion Meter allows you to view parameters such as the last calibration data and current parameter. It also allows you to set the parameter directly and view saved data. The

meter stores data according to the parameter measured. All saved data meets GLP standards. The meter can store 200 sets of data.

- 1 In the initial state, press [View].
- 2 Navigate to the type of data to view and press [Enter]. See Figure 34.


Figure 34 Viewing data

To view saved pH data (View Saved pH)

- 1 In the initial state, press [View].
- 2 Navigate to **View Saved pH** as shown in Figure 35. The upper area displays the current viewing mode and actual stored amount.


Figure 35 Viewing saved pH data

Each page can display 10 saved data at most. The data includes saved time and operator No. Use arrow keys to view each saved data. Press [Output] to print data or press [Delete] to delete data.

The output format follows:

=====

MODEL

3200I
ION METER
SOFT VERSION
VER 1.00
PRINT TIME
10:25:42
2010/06/20
OPERATOR NO
000

STORED NUM: 003

NO:001
OPERATOR NO: 000
STORED TIME: 10:19:00
2010/06/20
SLOPE: 100.00%
E0: -0.0mV

POTENTIAL: 0.0mV
pH: 7.000pH
TEMP: 25.0c
TC.TYPE: MTC

NO:002
OPERATOR NO: 000
STORED TIME: 10:19:00
2010/06/20
SLOPE: 100.00%
E0: -0.0mV

POTENTIAL: 0.0mV
pH: 7.000pH
TEMP: 25.0c
TC.TYPE: MTC

NO:003
OPERATOR NO: 000
STORED TIME: 10:19:00
2010/06/20
SLOPE: 100.00%
E0: -0.0mV

POTENTIAL: 0.0mV
pH: 7.000pH

```
TEMP:           25.0C
TC.TYPE:        MTC
=====
=====
```

To view saved pX data (View Saved pX)

The meter saves data according to the ion mode. Each ion mode can save 100 sets of pX data and 100 sets of ion concentration data. Since the meter can save data from six different ion modes, there are 600 sets of pX data and 600 sets of ion concentration data.

- 1 In the initial state, press [View].
- 2 Navigate to **View Saved pX** as shown in **Figure 36**. The upper area displays the current viewing mode and actual stored amount.


Figure 36 Viewing saved pX data

- 3 Press [**Output**] to print data or press [**Delete**] to delete data.

The meter supports six ion modes at once. To view saved data in other ion modes, press [**Setup**] (or [**Mode**] or [**View**]) to change to other ion modes. The meter will ask **View Other Ion?** Press [**Enter**] to view data of other ion modes. **Figure 37** shows **View Na⁺** ion mode is changed to **View Ag⁺** ion mode in **View saved pX**.


Figure 37 Changing ion

Press [**Output**] to print data or press [**Delete**] to delete data.

The output format follows:

```
=====
MODEL
      3200I
      ION METER
SOFT VERSION
      VER 1.00
PRINT TIME
      10:26:23
      2010/06/20
OPERATOR NO
      000
*****
STORED NUM:      003
*****
NO:001
OPERATOR NO:      000
STORED TIME:     10:34:00
                  2010/06/20
SLOPE:          100.00%
E0:             118.3mV

POTENTIAL:       59.2mV
pX:              1.000pX
TEMP:            25.0c
ION NAME:        Na
*****
NO:002
OPERATOR NO:      000
STORED TIME:     10:34:00
                  2010/06/20
```

2 Operation

```
SLOPE:          100.00%
E0:            118.3mV

POTENTIAL:      59.2mV
pX:             1.000pX
TEMP:           25.0c
ION NAME:       Na
*****
NO:003
OPERATOR NO:    000
STORED TIME:    10:34:00
                2010/06/20
SLOPE:          100.00%
E0:            118.3mV

POTENTIAL:      59.2mV
pX:             1.000pX
TEMP:           25.0c
ION NAME:       Na
=====
```

To view saved CONC data (View Saved Ion Concentration)

This is very similar to **View Saved pX**.

Figure 38 View saved ion concentration

View pH calibration data (View pH Last Calib)

- 1 In the initial state, press [View] and select **pH Last Calib**.
- 2 Press [**Enter**] to view last calibrated pH data as shown in **Figure 39**.


Figure 39 Viewing pH calibration data

The upper area shows the last calibrated data. The bottom area shows the current buffer groups. Use the arrow keys to view detailed calibrated data. To print current parameter data, connect to a PC with USB and press [**Output**].

The output format follows:

```
=====
MODEL
      3200I
      ION METER
SOFT VERSION
      VER 1.00
PRINT TIME
      10:11:27
      2010/06/20
OPERATOR NO
      000
*****
PH CALIB DATA
CALIB TIME:    08:12:00
```

2 Operation

2008/06/01
OPERATOR NO: 000

POINT 1
pH: 4.003pH
POTENTIAL: 177.3mV
TEMP: 25.0c

POINT 2
pH: 6.864pH
POTENTIAL: 8.0mV
TEMP: 25.0c

POINT 3
pH: 9.182pH
POTENTIAL: -129.1mV
TEMP: 25.0c

CALIB RESULT
SLOPE 1: 100.00%
E0 1: -0.0mV
SLOPE 2: 100.00%
E0 2: -0.0mV
=====

View pX calibration data (View pX Last Calib)

- 1 In the initial state, press [View] and select **pX Last Calib**.
- 2 Press [**Enter**] to view calibrated data in current ion mode. See **Figure 40**.


Figure 40 pX Last Calib

- 3 Press arrow keys to view detailed calibrated data, or press [**Setup**] to set new ion mode and print calibrated data.

Press [Output] to print data or press [Delete] to delete data.

The output format follows:

```
=====
MODEL
      3200I
      ION METER
SOFT VERSION
      VER 1.00
PRINT TIME
      10:23:44
      2010/06/20
OPERATOR NO
      000
*****
PX CALIB DATA
CALIB TIME:      08:12:00
                2008/06/01
OPERATOR NO:      000
ION NAME:        Na

      POINT 1
pX:          1.000pX
POTENTIAL:    59.2mV
TEMP:         25.0c
*****
      POINT 2
pX:          5.000pX
POTENTIAL:   -177.5mV
TEMP:         25.0c
*****
      POINT 3
pX:          7.000pX
POTENTIAL:   -295.8mV
TEMP:         25.0c
*****
BLANK CONC:
CONC:        0.000e+00mol/L
*****
CALIB RESULT
SLOPE 1:      100.00%
E0 1:         118.3mV
SLOPE 2:      100.00%
E0 2:         118.3mV
=====
```

To view other parameters during measurement

You can view other parameters during a measurement, even if you did not select them during setup. To view other parameter values during a measurement:

- 1 Press [**Mode**] and the display will highlight the measuring window. Press [4/◀] or [6/▶] repeatedly to switch parameters. When you stop pressing [4/◀] or [6/▶] for more than a few seconds, the meter automatically returns to the normal measuring mode display.
- 2 Press [**Output**] to print or press [**Save**] to save displayed values.

To Turn Off the 3200I Ion Meter

After using the meter, press [**Save**] to save the data, then press [**On/Off**] to turn off the meter. When not in use, all electrodes should be soaked in distilled water. If the meter is out of use for a long time, pay attention to:

- 1 Disconnect power adaptor to avoid damage to the analyzer and the power supply.
- 2 Keep the socket of the meter clean and dry. Do not allow the meter to contact acid, alkaline, or salt solutions.
- 3 If the meter is out of use, insert the short circuit plug in the back of the meter to prevent damage to the open circuit.
- 4 After measurement, store the electrode in reference filling solution. For refillable electrode, plug the filling hole. If the electrode will be out of use for a long time, put it back in the box and store it at ambient temperature.

Software Operation

See also “[Installing the Optional Software](#)” on page 14.

To use the optional G4390A Electrochemical Data Collecting software

- 1 Connect the meter to the power adapter.
- 2 Press [**On/Off**] to turn on the meter.
- 3 Connect the meter to the computer using the USB cable.
- 4 Run the Electrochemical Data Collecting Software in the computer and the software will recognize the meter model and type. The software provides the functions corresponding to the meter type. Refer to the software documentation for details.

To use the Data Printing software

- 1 Connect the meter to the power adapter.
- 2 Press [**On/Off**] to turn on the meter.
- 3 Connect the meter to the computer using the USB cable.
- 4 Launch the printing software.
- 5 With the meter taking measurements or when viewing calibration data or stored data, press [**Output**]. The software receives the data. Once received, print the data as desired. For details, refer to the Data Printing software documentation.

3

Troubleshooting and Maintenance

Troubleshooting 72

General Troubleshooting Procedure 73

Meter Self-Diagnosis and Messages 76

Maintenance 79

Consumables and Accessories for the 3200I Ion Meter 80

This chapter describes how to verify whether the meter is working properly and how to maintain the meter.

Troubleshooting

If the meter is not working properly during operation, perform trouble-shooting according to the following table.

Table 4 Troubleshooting

Number	Failure mode	Failure cause	Solutions
1	No display after meter starts up LCD display is not lighting.	1 Incorrect installment of power adaptor. 2 The power supply does not meet with the requirements. 3 Power adaptor is damaged. 4 The power socket has a poor contact. 5 LCD is damaged (After start up for a few minutes, user can hear buzzing sound after pressing [On\Off] but there is no display on LCD.)	1 Connect it again according to the instructions. 2 Use an approved power supply. 3 Replace the power adaptor. 4 Check all connections for good contact.
2	No buzzing when pressing key.	1 You may have pressed invalid keys under current setup. 2 The buzzer has been damaged.	The buzzer only sounds when you press a valid key.
3	No response when pressing key.	1 You may have pressed invalid keys under current setup. 2 The key has been damaged.	Press the valid key for operation.
4	Meter does not display 25 °C when no ATC electrode is connected with meter.	When the meter is not connected with an ATC electrode, the meter will automatically switch to manual temperature setting.	Manually input temperature.
5	The reading is not stable after a long time.	1 The electrode has been damaged or aged. 2 There is strong electrical signal interference source nearby (electrical leaking, strong electromagnetic field, and so forth).	1 Replace the electrode. 2 Remove the electrical signal interference source. Lift or move the beaker away from the source of interference. Shield the meter and beaker from the electromagnetic field. Use a wire to connect the meter with a ground line. For the wire, one terminal is connected with the meter and the other terminal is connected with interference source.

General Troubleshooting Procedure

During the use of the meter, there are many factors that may affect the measurement, including the electrode, the standard solutions used to calibrate electrode, the sample solution, the temperature during measurement, incorrect operation, or use, and so forth. When the measured results are significantly different from what was expected, first determine whether the meter itself or factors other than the meter caused the error. Follow the suggestions described below to diagnose the problem. Since there are multiple measuring methods, you will need to diagnose the meter based on your application conditions.

Diagnose whether the errors are caused by meter hardware, electrode calibration data error, parameter errors, or other errors, including the sample solution, solution preparation, buffer solutions used for calibration, and so forth.

For factors other than the meter, compare measurement results between known solutions. Put the electrode in different standard solutions. Based on the comparison, judge whether the deviation is resulted from the electrode, the solution, or something else.

Check the temperature measurement

Because temperature measurements are used for all measurements and calculations, always ensure temperature measurement is functioning well.

Connect the meter to the ATC temperature diagnostic tool (5185-8390) shipped with the meter. Turn on the meter to enter into measurement state. With the ATC temperature diagnostic tool connected, the meter should display a temperature reading between 49.0 to 51.0 °C. If yes, the meter is correctly measuring temperature. If the meter displays a temperature reading significantly different from 50 °C, there is a problem with the meter. Contact Agilent service.

Restore factory default settings

If other troubleshooting does not solve the problem, you can restore the meter to the factory default settings. This clears all user calibrations and settings so that you can verify that an erroneous calibration or setting is not the cause of the problem. After restoring the factory default settings, retest the meter.

See “[Set Default](#)” on page 36.

Troubleshoot ion measurement

Check the potential (mV) value

- 1 Disconnect the ISE electrode.
- 2 Connect the meter to the short circuit plug (G4383-40000) shipped with the meter.
- 3 Turn on the meter.
- 4 Start measuring pH.

At this moment, the meter should display a potential (mV) reading between -0.5 to 0.5 mV (also called mV Zero). If yes, this indicates that the meter is functioning properly. If the potential (mV) reading significantly deviates from zero, then use the Adjust mV Zero function to adjust mV Zero. (see “[To adjust mV zero](#)” on page 38) After adjustment, the potential (mV) reading should be near to zero. If the potential value (mV) reading cannot be adjusted near to zero by the mV Zero adjustment, this indicates a problem the meter. Contact Agilent service.

Check the pH value

Disconnect the ATC probe. Set the meter for manual temperature input. Manually set the temperature to 25 °C. When the potential (mV) is zero, the meter should display a pH reading between 6.99 to 7.01 pH. If yes, this indicates the electrode calibration data of the meter is acceptable. If not, there may be a problem with the electrode calibration data. Check the slope of electrode. Press [**View**], select **Last pH Calib**, and press [**Enter**]. Now, view the calibration data and slope from

last calibration. If the calibration data is severely out of the measurement requirement (the meter displays error message—refer to “[Error warning messages](#)” on page 76), this indicates a problem in the last electrode calibration which may lead to a large measurement error. Recalibrate the electrode. Pay attention to the selection of the Buffer Group during calibration.

NOTE

In the measurement state, the meter may display an error message simultaneously with measurement data.

Check the pX value

If user observes a large deviation in pX measurement, the deviation in pX measurement may be related to the stirring rate and reagents used to adjust ionic strength. This procedure can only diagnose the meter. To diagnose other factors, refer to professional application notes.

To diagnose the meter, first check the temperature measurement (see “[Check the temperature measurement](#)” on page 73). Next check the potential measurement (see “[Check the potential \(mV\) value](#)” on page 74), and then the pH value measurement (see “[Check the pH value](#)” on page 74).

For factors other than the meter, compare standard versus sample solutions to determine the cause of the problem. Put the electrode in different standard solutions and sample solutions. Check the corresponding potential (mV) value, pH or pX value. Based on the comparison, judge whether the deviation results from the electrode, the solution, or something else.


Meter Self-Diagnosis and Messages

The meter supports self-diagnosis, which can find some common errors caused by the meter, electrode/probe, solution or operation. It will deal with these errors differently according to their severity. There are two levels of severity: **Severe system errors** hinder further use of the meter or make the meter unable to finish necessary tasks. In this case, replacement or maintenance of the meter is required.

Minor errors, which are prompts or warnings, are caused by various reasons and can be corrected in several ways. However, you must pay attention to these errors to ensure measurement integrity and reliability. If you ignore these error messages, operation can continue.

Severe errors

When severe errors occur, the meter will shut down and display a prompt window. In this case, you cannot do any further operation except turn off the meter. The meter will display an error message as shown in [Figure 41](#).

[Figure 41](#) Severe error message

Error warning messages

User should pay attention to messages or warnings resulting from an improper electrode/probe, solution, or operation. Read these error messages carefully. Maintain and operate the meter and electrode according to their instructions to minimize the occurrence of such errors.

Figure 42 shows a warning that the temperature is out of range.


Figure 42 Temperature warning

is the icon for an error warning. The icon includes the error code. See [Table 5](#).

Table 5 Error warning codes

Code	Description	Solution	Memo
00	Conductivity measuring module error	Contact Agilent for service.	Serious error
01	D0 measuring module error	Contact Agilent for service.	Serious error
02	Temperature measuring module error	Contact Agilent for service.	Serious error
03	Data storage error	Contact Agilent for service.	Serious error
20	Potential is out of range	Replace electrode.	-1999.9 to 1999.9 mV
21	pH/pX is out of range	Replace electrode.	-3.000 to 21.000 pH
22	Conductivity is out of range	Replace probe and solution.	0 to 2000 mS/cm
23	Resistivity is out of range	Replace probe and solution.	0 to 100 MΩ•cm
24	Temperature is out of range	Replace electrode and decrease solution temperature.	-6.0 to 120.0 °C
25	D0 electric current is out of range	Replace probe.	0 to 4000 nA
26	pH electrode slope is out of range	Replace electrode and recalibrate it.	80 to 120%
27	Failed to recognize pH buffers	Replace electrode, setup proper buffer group and replace buffers.	

3 Troubleshooting and Maintenance

Table 5 Error warning codes (continued)

Code	Description	Solution	Memo
28	Temperature of pH buffer is out of range	Cool or heat buffer.	
29	Calibrate the same buffer repeatedly	Replace the buffer.	Due to incorrect operation
30	Number of pH buffers exceeds maximum.	Remove one or more buffers from the group.	
31	pH buffers conflict with each other	Remove one or more buffers with overlapping pH values.	
32	The data storage is full.	Overwrite the previous data and store new data.	Do not delete all data
33	The internal clock has a low battery.	Set time manually.	
34	Number of customer-defined ions exceeds maximum	Delete one or more unnecessary customer-defined ions.	
35	The maximum ion mode number stored	Delete all storage data of certain ion mode.	

Maintenance

Electrode cleaning

General electrode cleaning includes inorganic cleaning, organic cleaning, grease cleaning, protein precipitation cleaning and glass sensitive membrane regeneration. The type of cleaning needed depends on the contaminants and the electrode. After one or more cleaning procedures, rinse the outside of the electrode with distilled water. Siphon the filling solution away and add fresh solution. Repeat 2 or 3 times. Store the pH electrode in the storage solution.

Electrode storage

pH measuring electrode

When the electrode is not in use, store it in reference filling solution. For refillable electrodes, plug the filling hole. When the electrode is not used for a long time, put the electrode back in the box and store it a dry place at ambient temperature.

Consumables and Accessories for the 3200I Ion Meter


Table 6 Consumables and accessories for 3200I Ion Meter

Order No.	Model and Name	Description
G4388A	3200SA Stirrer	Electrode holder and magnetic stirrer: combined to stir solution with stable and precise speed with a large adjustable range.
G4389A	3200EA Electrode Holder	It is used to fix different electrodes.
5185-8389	AC Adaptor	100/240 VAC, 1 A, 50/60 Hz
G4383-40000	Short Circuit Plug	Zero potential calibration and diagnostication.
G4388-27000	Stirring Bar	It is used with stirrer.
5190-4002	I9111 Fluoride Combination ISE	Combination Electrode, plastic body and refillable.
5190-4004	I9121 Ammonia Combination ISE	Combination Electrode, plastic body and refillable.
5190-4005	I9131 Sodium Combination ISE	Combination Electrode, glass body and refillable.
5190-3988	P3211 pH Combination Electrode	Combination Electrode, glass body and refillable.
5190-3989	P3212 pH Combination Electrode	Combination Electrode, plastic body and sealed
5190-3990	P3311 pH Triode Combination Electrode	Triode Combination Electrode with built-in ATC probe, glass body and refillable
5190-3991	P3111 pH Electrode	Glass body and single electrode
5190-3992	P3213 pH Combination Electrode	Combination Electrode, plastic body, flat sensitive glass membrane and refillable
5190-3993	P3214 pH Combination Electrode	Combination Electrode, plastic body, spear-tip sensitive glass and sealed
5190-4003	R8111 Reference Electrode	Glass body, Ceramic junction, Ag/AgCl and single-junction
5190-3999	ORP8211 ORP Electrode	Combination Electrode, glass body and refillable
5190-3998	T7111 ATC Probe	Stainless Steel body Measuring Range: 0 to 100
5190-0533	pH Buffer Pouches	Bottle, 3 X 250 mL, 4.01, 7.00, 10.01 pH
5190-0534	pH Buffer Pouches	Bottle, 3 X 250 mL, 4.00, 6.86, 9.18 pH

Table 6 Consumables and accessories for 3200I Ion Meter (continued)

Order No.	Model and Name	Description
5190-0541	pH 1.68 Buffer	Bottle, 3 X 250 mL
5190-0538	pH 4.00 Buffer	Bottle, 3 X 250 mL
5190-0535	pH 4.01 Buffer	Bottle, 3 X 250 mL
5190-0539	pH 6.86 Buffer	Bottle, 3 X 250 mL
5190-0537	pH 7.00 Buffer	Bottle, 3 X 250 mL
5190-0540	pH 9.18 Buffer	Bottle, 3 X 250 mL
5190-0536	pH 10.01 Buffer	Bottle, 3 X 250 mL
5190-0542	pH 12.46 Buffer	Bottle, 3 X 250 mL
5190-0546	Reference solution	Bottle, 3 X 30 mL
5190-0543	Sodium ISE Reference solution	Bottle, 3 X 30 mL
5190-0544	Ammonia electrode membrane kit	5 pcs

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