

A Comparison of the Relative Cost and Productivity of Traditional Metals Analysis Techniques Versus ICP-MS in High Throughput Commercial Laboratories

Application

Author

Steve Wilbur
Agilent Technologies, Inc.
3380 146th Pl SE Suite 300
Bellevue, WA 98007 USA

Abstract

A financial model was developed to help the metals laboratory using graphite furnace atomic absorption and inductively coupled plasma optical emission spectroscopy calculate the potential savings by switching to inductively coupled plasma mass spectrometry. Results based on several typical laboratory examples are presented.

Introduction

The past 5 years have seen significant growth in the use of inductively coupled plasma mass spectrometry (ICP-MS) for the analysis of trace metals in many applications in the environmental, semiconductor, geological, and health sciences industries. This growth is driven by three factors. First is the need for increasingly lower limits of detection for many metals in many applications. Second is the significantly improved performance, reliability, and ease of use of modern ICP-MS instruments. And third is economics.

Traditionally, most elemental analysis has been performed by either atomic absorption (AA) or optical emission spectroscopy (OES). Generally, the ultratrace (sub-ppb) elements were measured by graphite furnace atomic absorption (GFAA), a highly sensitive single-element technique. The trace and minor (ppb to ppm) elements were measured

by inductively coupled plasma optical emission spectroscopy (ICP-OES), which is less sensitive but capable of simultaneous multi-element analysis.

As the need for sub-ppb detection limits extends to more elements in more samples, ICP-OES becomes less useful and the reliance on GFAA increases. However, GFAA, while sensitive, is slow, expensive to operate, and has limited dynamic range. Because GFAA is much slower than ICP-OES, many routine labs have a dedicated GFAA instrument for each analyte that is required to be measured by GFAA - multiple GFAAs working with one ICP-OES. Furthermore, the analysis of mercury will add the need for a third technique, either cold vapor AA or atomic fluorescence. However, in the interest of simplicity, a separate mercury analyzer was not considered in the examples used. Each of these techniques may require separate sample handling and preparation, as well as separate analysis, data processing and archival, significantly increasing the cost per sample.

The subject of this application note is to evaluate the productivity and cost effectiveness of ICP-MS as a routine, highly sensitive, multi-element technique where a single ICP-MS instrument has the potential to replace an ICP-OES, multiple GFAAs, and a mercury analyzer for most routine elemental analyses. The analytical applicability of ICP-MS to many types of samples is already well established. More recently, the introduction of the Octopole Reaction System on the 7500 Series ICP-MS instruments from Agilent has removed the final performance barriers that have prevented ICP-MS being proposed as a complete replacement for GFAA and ICP-OES.



Agilent Technologies

Methods

To facilitate this study, a spreadsheet-based sample cost comparison model was developed in Excel. This tool allows the user to provide detailed parameters related to numbers and types of samples, as well as associated costs of sample preparation, instrumentation, and analysis. Output is simply cost of analysis per sample. Also reported are the total time required for sample analysis per month, the number of analysts required, and the number of instruments. The model compares the results for GFAA, ICP-OES, and ICP-MS. While it will allow almost any values to be entered for most parameters, the results presented here are based on values obtained from several commercial laboratories doing these analyses. No model can exactly predict the results for all situations and still be simple enough to be useful. Therefore, in the interest of simplicity, a number of assumptions were made in the design of the model and in the example data entered. We feel that the assumptions are realistic and do not impart significant bias on the results. The tool is easy to use and can allow a laboratory to quickly and simply evaluate the cost effectiveness of the three techniques based on laboratory-specific information.

Assumptions

- GFAA system costs US\$30K
- ICP-OES system costs US\$100K
- ICP-MS system costs US\$180K
- Cost of funds (finance) is 6%
- General facilities costs, such as laboratory space, utilities etc., are ignored since they are difficult to estimate and do not significantly affect the results in most cases.

- An instrument operator can keep a modern, automated GFAA, ICP-OES, or ICP-MS running for two shifts (16 hours) per day. When analysis times exceed 16 hours per day for any technique, additional instrumentation and operators will be required. Instruments are added in increments of one; operators are added in fractions since it is assumed that they can be shared with other tasks in the laboratory and cost calculations are based only on the portion of time the operator spends on the specific analysis.
- GFAA is a single element technique. Instruments with multiple lamps still perform a single analysis at a time. Typical analysis time is 90 seconds per element and each element requires two replicate analyses (burns).
- ICP-OES and ICP-MS are multi-element techniques and the number of elements does not significantly effect the analysis time. This is not strictly true, but the assumption is reasonable for the sake of simplicity.
- GFAA will use pressurized argon and the consumption is 40 hours of use per cylinder (\$100).
- GFAA graphite tubes and platforms cost \$50 per set and last for 100 burns.
- ICP-MS and ICP-OES will use liquid argon and the typical consumption is 3 weeks of use per dewar (\$250).
- ICP-MS detectors last typically for 3 years and the cost per year is amortized based on 3-year lifetime.

Results

Several typical laboratory scenarios were evaluated by varying the current instrument complement of the laboratory, and by varying the current and anticipated number of samples to be analyzed per month. Also examined was the effect of the number of elements that must be analyzed by GFAA (in the case of laboratories without ICP-MS) to meet required DLs.

Scenario 1

Laboratory currently has one GFAA plus one ICP-OES, which are paid for. ICP-MS must be purchased and amortized over 3 years. See Table 1.

Table 1. Scenario 1

Samples/ month	GFAA elements	# GFAA required	Cost/sample GFAA + ICP-OES	# ICP-MS required	Cost/ sample ICP-MS	Savings/ month
400	8	1	\$41	1	\$30	\$4,536
1000	8	2	\$33	1	\$15	\$18,196
5000	8	9	\$31	2	\$9	\$112,968

Scenario 2

Laboratory currently has two GFAA plus one ICP-OES, which are paid for. ICP-MS must be purchased and amortized over 3 years. See Table 2.

Table 2. Scenario 2

Samples/ month	GFAA elements	# GFAA required	Cost/sample GFAA + ICP-OES	# ICP-MS required	Cost/ sample ICP-MS	Savings/ month
400	8	1	\$41	1	\$30	\$4,536
1000	8	2	\$32	1	\$15	\$17,283
5000	8	9	\$31	2	\$9	\$112,055

Scenario 3

Laboratory currently has no instrumentation and must decide on purchasing GFAA plus ICP-OES versus ICP-MS. See Table 3.

Table 3. Scenario 3

Samples/ month	GFAA elements	# GFAA required	Cost/sample GFAA + ICP-OES	# ICP-MS required	Cost/ sample ICP-MS	Savings/ month
400	8	1	\$51	1	\$30	\$8,491
1000	8	2	\$37	1	\$15	\$22,151
5000	8	9	\$32	2	\$9	\$116,923

Scenario 4

Comparison of costs per sample as a function of number of GFAA elements. (All instruments must be purchased.) See Table 4.

Table 4. Scenario 4

Samples/ month	GFAA elements	# GFAA required	Cost/sample GFAA + ICP-OES	# ICP-MS required	Cost/ sample ICP-MS	Savings/ month
1000	2	1	\$24	1	\$14	\$9,601
1000	4	1	\$28	1	\$14	\$12,751
1000	8	2	\$38	1	\$14	\$22,151
1000	10	3	\$42	1	\$14	\$27,490

Discussion

In all cases, even when the laboratory already owns two graphite furnaces and one ICP-OES (a common configuration) and must purchase the ICP-MS, the cost per sample is lower for ICP-MS. This is mainly due to the high cost of consumables for GFAA plus the fact that GFAA and ICP-OES requires two separate sample prep steps. Additionally, as the number of samples increases from a conservative number of 400 per month to 1000 and 5000 per month, the differential becomes much greater. This is caused by rapidly increasing labor costs for GFAA, as well as the much higher sample capacity of ICP-MS, lower consumables costs, and requirements for only a single sample prep.

Return on Investment for ICP-MS

A simple return on investment (ROI) can be calculated from the above tables. In this case, the cost per month of the new ICP-MS system is approximately US \$5500.00 (assuming purchase price of US\$180K financed for 3 years at 6%). Figure 1 shows the payback times for a laboratory that already owns two GFAAs and one ICP-OES as a function of the sample load. The y-axis represents the accumulated monthly savings of using ICP-MS versus GFAA + ICP-OES for three different sample loads compared to the unpaid balance on the ICP-MS instrument. As can be seen, the accumulated savings of ICP-MS is equal to the payoff amount after just 4 months when analyzing 2000 samples per month. Even when analyzing as few as 400 samples per month, the accumulated savings is sufficient to pay off the ICP-MS instrument in around 20 months. In this case, eight furnace elements are assumed. Other assumptions are as above.




Figure 1. Cumulative return on investment of ICP-MS purchase for three sample levels plotted against the monthly unpaid balance on the ICP-MS. In this case, it is assumed that the accumulated revenue will be used to pay off the loan when the balance equals the residual loan amount. At that point, the net monthly revenue is increased by the loan amount. In this example, laboratories running 2000 samples per month will be able to pay off the ICP-MS in about 4 months, 1000 sample laboratories in about 8 months, and 400 sample laboratories in about 20 months. At the end of 36 months (the original loan period), net revenue exceeds \$200K for the 400 sample lab, \$750K for the 1000 sample lab, and \$1.7 million for the 2000 sample lab.

Conclusions

For almost any metals laboratory, analyzing at least 100 samples per week (400 per month) and using a combination of GFAA and ICP-OES for the analysis, converting to ICP-MS will save money. Depending on the number of samples, the payback for the ICP-MS can be as short as a few months. The cost advantages are not reduced significantly, even if the laboratory already owns its GFAA and ICP-OES instruments. They are also not significantly affected by the number of GFAA elements. As Scenario 4 shows, for the laboratory analyzing at least 1000 samples per month with only two elements by GFAA, the cost savings of switching to ICP-MS is approximately \$10,000 per month. Add to this the increased confidence in results obtained by ICP-MS, the ability to analyze all analyte elements at GFAA (or better) DLs, and the robustness and simplicity of operation of modern ICP-MS instruments, and the choice becomes simple. The productivity of ICP-MS in a high-volume laboratory can quickly pay off the purchase price and increase laboratory profitability significantly.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2005

Printed in the USA
January 17, 2005
5989-1585EN



Agilent Technologies