




A Full Family of Polar PEG Columns From Agilent

Product Brief



As the frequency of polar laboratory expeditions increases and the range of exploration expands, so has our line of polyethylene glycol (PEG) capillary columns.

Added to the already powerful lineup of HP-Wax, HP-INNOWax, HP-FFAP, and HP-20M capillaries is the new family of HP-Basic Wax columns. Each of these five column types is available in a variety of diameters, lengths and film thickness including HP-INNOWax, HP-Wax, HP-Basic Wax, HP-FFAP and HP-20M.

Even though each phase is based on the polyethylene glycol polymer, strict control of the cross-linking and deactivation processes results in a variety of unique phase characteristics to meet the varying analysis needs of your laboratory.

Columns and Supplies


Versatile HP-INNOWax

Since its introduction, the HP-INNOWax column has become the workhorse of the PEG columns. Because of a compatibility with a

wide range of analytes and solvents, long lifetime, and low-bleed characteristics, the HP-INNOWax columns are used in a wide range of applications.

Industrial Solvents

1. Pentane
2. Methyl formate
3. Acetone
4. Ethyl acetate
5. Methyl ethyl ketone
6. Methanol
7. 2-Methyl-2-propanol
8. Methylene chloride
9. Benzene
10. Ethanol
11. 2-Butanol
12. Toluene
13. n-Propanol
14. Ethyl benzene
15. p-Xylene
16. m-Xylene
17. 1-Butanol
18. o-Xylene




Agilent Technologies
 Innovating the HP Way

The HP-INNOWax column has the highest upper temperature limit of any of the wax phases. With a maximum upper temperature limit for programmed operation of 270°C for standard dimensions and 250°C for 0.53 mm columns, the HP-INNOWax column is ideal for the analysis of high boiling point analytes that may be difficult or impossible to do on other "Wax" columns. For example, free fatty acids to C24 can be analyzed with this column. The versatility of the HP-INNOWax column is second to none. The superior inertness produces sharp, reproducible peaks for both basic and acidic compounds without requiring regeneration of the column with injections of acids or bases.

Aldehydes and Acids


- | | |
|---------------------|----------------------|
| 1. Butanal | 10. iso-Butyric acid |
| 2. 2-Methyl butanal | 11. Butyric acid |
| 3. Pentanal | 12. iso-Valeric acid |
| 4. Hexanal | 13. Valeric acid |
| 5. Heptanal | 14. Hexanoic acid |
| 6. Octanal | 15. Heptanoic acid |
| 7. Acetic acid | 16. Octanoic acid |
| 8. Decanal | 17. Decanoic acid |
| 9. Propanoic acid | |



Column: HP-INNOWax
30 m x 0.32 mm x 0.5 µm
(Part number 19091N-213)
Carrier: Helium, 40 cm/sec, 11.7 psi (60°C),
2.5 ml/min constant flow
Oven: 60°C (1 min) to 250°C at 10°C/min
Injection: Split (40:1), 0.5 µl, inlet at 250°C
Detector: FID 275°C

PEG Columns for Petrochemical Applications

Reformate Gasoline Analysis




A
Column: HP-Wax, 30 m x 0.32 mm x 0.5 µm
(Part number 19091X-213)
Carrier: Helium, 6.4 psi, constant pressure
Oven: 60°C (4 min) 10°C/min to 140°C
(0 min) 140–200°C (4 min) at 15°C/min
Injection: 250°C, 0.5 µl, Split 200/1
Detector: 275°C, FID

B
Column: HP-Wax, 10 m x 0.1 mm x 0.2 µm
(Part number 19091X-241)
Carrier: Helium, 37.3 psi, constant pressure
Oven: 60°C (1 min) 36.7°C/min to 140°C
(0 min) 140–200°C (2 min) at 55.1°C/min
Injection: 250°C, 0.1 µl, Split 800/1
Detector: 275°C, FID

Impurities in p-Xylene: ASTM D3798

Column: HP-INNOWax, 60 m x 0.32 mm x 0.5 µm
(Part number 19091N-216)
Carrier: He, 32 cm/sec, 19.9 psi (60°C), 2.5 ml/min constant flow
Oven: 60°C (1 min) to 92°C (4.5 min) at 4°C/min,
to 220°C (5 min) at 20°C/min
Injection: Split (split ratio 100:1), 0.5 µl Inlet 220°C
GC System: Agilent GC/FID (270°C)
Sample: Neat, 99%+

1. Non-aromatic hydrocarbon
2. Toluene
3. Ethylbenzene
4. p-Xylene
5. m-Xylene
6. o-Xylene



HP-INNOWax temperature range:

40°C to 260°C;

0.53 mm id columns:

40°C to 240°C

HP-INNOWax

Similar Phases

Stabilwax, Supelcowax-10,
CP WAX 52 CB, DB-WAXetr

Applications

Alcohols, Free Acids,
Aromatics, Essential Oils,
Solvents


The HP-INNOWax columns have the lowest bleed of any of the Wax phases. This results in lower detector noise so you can detect and quantify lower levels of compounds. This phase is also a good choice for use in mass spectrometry where column bleed can reduce the accuracy of peak identification.

The combination of a highly deactivated glass surface with the cross-linked and surface bonding results in a column that is compatible with a wide range of polar compounds. This makes the HP-INNOWax column the ideal choice when developing new methods. The HP-INNOWax family of columns have been used successfully in a variety of different industries and applications including petrochemical, solvent analysis, food and flavors, and beverage analysis.

PEG Columns for Petrochemical Applications


Impurities in Styrene: ASTM D5135

Column: HP-INNOWax, 30 m x 0.32 mm x
0.5 µm (Part number 19091N-213)
Oven: 80°C (5 min), 10°C/min to 150°C
(2 min)
Carrier: Helium 8.6 psi (28 cm/sec),
constant flow, 1.5 ml/min
Injection: 0.5 µl split 40:1, 230°C
GC System: Agilent GC/FID 240°C



Impurities in Ethylbenzene: ASTM D5060

1. Hydrocarbon
2. Benzene
3. Toluene
4. Ethylbenzene
5. p-Xylene
6. m-Xylene
7. Cumene
8. o-Xylene
9. Propylbenzene
10. Styrene



Column: HP-INNOWax, 60 m x 0.32 mm x 0.5 µm
(Part number 19091N-216)
Carrier: He, 32 cm/sec, 19.9 psi (60°C),
2.5 ml/min constant flow
Oven: 60°C (1 min) to 92°C (4.5 min) at
4°C/min, to 220°C (5 min) at 20°C/min
Injection: Split (split ratio 100:1), 0.5 µl Inlet 220°C
GC System: Agilent GC/FID 270°C
Sample: Neat, 99%+


Low Operating Temperature HP-WAX

The HP-Wax column has the lowest operating temperature of any of the Wax phases. The proprietary cross-linking lowers the usable temperature of the stationary phase to 20°C compared with 40°C for other Wax phases making this phase ideal for low boiling point compounds.

If your application is purge and trap or headspace, the HP-Wax column can be used with cryogenic cooling for cold trapping analytes. The HP-Wax column is designed to match the polarity of competitive Wax columns making substitution easy. The HP-Wax phase is the closest in polarity to the HP-20M so you can transfer your methods to a cross-linked phase without lengthy method development.

PEG Columns for Petrochemical/Chemical Applications

Impurities in Styrene




Solvent Analysis

HP-INNOWax, HP-Wax Comparison

- Pentane
- Methyl formate
- Acetone
- Ethyl acetate
- Methyl ethyl ketone
- Methanol
- 2-Methyl-2-propanol
- Methylene chloride
- Benzene
- Ethanol
- 2-Butanol
- Toluene
- n-Propanol
- Ethyl benzene
- p-Xylene
- m-Xylene
- 1-Butanol
- o-Xylene


Carrier: Hydrogen, 40 cm/sec, constant pressure
 Oven: 35°C (2.5 min), 10°C/min to 85°C
 Injection: 250°C, Split (100:1), 0.5 µl
 GC System: Agilent GC/FID, 300°C



PEG Columns for Petrochemical/Chemical Applications


Substituted Aromatics HP-Wax Column

1. Hexane	10. tert-Butylbenzene	Column:	HP-Wax, 60 m x 0.32 mm x 0.5 µm (Part number 19091X-216)
2. Benzene	11. Isobutylbenzene	GC System:	Agilent GC/FID, 300°C
3. Toluene	12. Mesitylene	Injection:	Split, 250:1, 0.5 µl injection, 250°C
4. Ethylbenzene	13. sec-Butylbenzene	Concentration	
5. 1,4-Xylene	14. Styrene	of analytes:	500 µg/ml in hexane
6. 1,3-Xylene	15. 1,2,4-Trimethylbenzene	Carrier:	Hydrogen at 40 cm/sec constant pressure
7. Cumene	16. Butylbenzene	Oven:	50°C for 1 min; 50–75 at 4°C/min, 75–100°C at 10°C/min, (hold 4 min) 100–150°C at 10°C/min
8. 1,2-Xylene	17. 1,4-Diisopropylbenzene		
9. Propylbenzene	18. trans-β-Methylstyrene		



Solvents in Paint Thinner

HP-Wax,
10 m x 100 µm x 0.2 µm
(Part number 19091X-241)



Temperature range:

20°C to 250°C;

Series 530-µm and thick-film
columns: 20°C to 230°C

HP-Wax

Similar Phases

DB-Wax, Rtx-wax,
Carbowax 20M

Applications

Alcohols, Aromatics, Essential
Oils, Solvents

The HP-Wax columns are available in 100 µm diameters for fast analysis. The 100 µm columns have the same phase chemistry as the standard dimension HP-Wax columns so translating methods to the smaller dimensions is simple. If you are using the HP-INNOWax columns, you can still translate your methods to the 100 µm HP-Wax columns. Most methods translate directly without the need for time consuming method development. Faster analysis means more productivity for your lab.

The HP-Wax columns are cross-linked and bonded to ensure long life-time. The columns can be solvent rinsed to extend the column lifetime and reduce costs. Although the HP-Wax columns are used for the same application areas as the HP-INNOWax columns, the HP-Wax has slightly different selectivity for some compounds like alcohols (see comparative chromatogram on page 4). If some analytes are difficult to separate on the HP-INNOWax columns and you need the polarity of a Wax column, try the HP-Wax column.

HP-FFAP Columns for Free Fatty Acid Analyses


The HP-FFAP columns are designed primarily for the analysis of organic acids, free fatty acids or samples that require quantitation of acidic impurities. The stationary phase is modified with acid to provide a very inert column that can accommodate the demanding analysis of acids dissolved in water. Free fatty acids up to C24 can be analyzed without costly and time consuming derivatization. If you want to analyze the derivatized free acids, the FFAP column is a good choice to resolve the derivatized saturated and unsaturated acids to C24.

The HP-FFAP is cross-linked and bonded to resist the damage that can occur when injecting water based samples and can be used at operating temperatures of 60°C to 240°C. These columns do not require several injections of acids (priming) before accurate quantitation is possible so accurate results are obtained quickly. The columns are also solvent rinsable for long column lifetimes.

HP PEG Columns for Foods and Flavors Applications


Free fatty acids (C2–C7)

Injection: 0.5 µl, split 20/1, 260°C
 Column: HP-FFAP, 25 m 0.32 mm 0.5 µm (Part number 19091F-112)
 Carrier: Helium @ 10.8 psi (42 cm/sec)
 Oven: 80°C (1 min), to 120°C (20°C/min), then to 205°C (6.13°C/min), 205°C (2 min)
 GC System: Agilent GC/FID, 260°C



Perfume Calibration Mix



1. α-Pinene 10. Anisaldehyde Column: HP-20M (Carbowax 20 M)
 2. β-Pinene 11. Methyl cinnamate 50 m x 0.32 mm x 0.3 µm
 3. δ-Limonene 12. Ethyl cinnamate (Part number 19091W-015)
 4. Citronellal 13. Eugenol Carrier: Hydrogen, 64 cm/sec
 5. Menthol 14. Jasmonal Oven: Temperature program listed
 6. Benzyl acetate 15. Diethyl phthalate Injection: Split, 210°C
 7. Citronellyl acetate 16. Coumarin Detector: FID, 230°C
 8. Geraniol 17. Vanillin
 9. Hydroxy citronellal 18. Benzyl benzoate



PEG Column	Features	Advantages	Benefits
HP-INNOWax	Highest Upper Temperature Limit Chemical Compatibility Low Bleed Bonded and Cross-linked Highly Inert	Analyze High-Boiling Point Compounds General Purpose Column Low Detector Noise Solvent Rinsable General Purpose Column	Wide Usage Wide Usage Best Choice for MS Use Low Lifetime Broad Analyte Compatibility
HP-Wax Columns	Lowest Operating Temperature Limit Similar in Polarity to HP-20M Available in 100 µm id Bonded and Cross-linked Highly Inert	Analyze Low-Boiling Point Analytes Transfer Older Methods to Cross-linked Phase Fast Analysis Solvent Rinsable Excellent for Slightly Basic Compounds	Wide Usage Easier Method Transfer Higher Sample Throughput Long Lifetime Accurate Quantitation
HP-FFAP Columns	Designed for Acids Bonded and Cross-linked	Can Inject Acids Directly Solvent Rinsable	No Need for Derivatization Long Column Lifetime
HP-Basic Wax Columns	Selective for Amines Upper Temperature Limit to 260°C Bonded and Cross-linked	Sharp Peaks Analyze High-Boiling Point Compounds Solvent Rinsable	No Need to Precondition Expanded Range of Analytes Long Column Lifetime

PEG Columns for Foods and Flavors Applications



Derivatized Olive Oil (FAMES) on Two Different Columns



Common Parameters

GC:	6890	Liner:	Deactivated with glass wool plug
Sample:	Greek Olive Oil (Derivatized) in Methylene Chloride	Detector Temp.:	260°C
Split Ratio:	20/1	Hydrogen Flow:	40 ml/min
Injector Temp.:	260°C	Air Flow:	450 ml/min
Injection Size:	1.5 µl	Makeup Flow:	45 ml/min (Nitrogen)
		Constant Flow Mode	

Jasmine Extract



HP-FFAP

Similar Phases

DB-FFAP, Stabilwax-DA, OV-351, CP WAX 58CB, Nukol

Applications


Acids, Alcohols, Aldehydes, Acrylates, Ketones, Nitriles

Temperature Range

60°C to 240°C

Ethoxyethanol HP-FFAP Column

- | | |
|--------------------------|---------------------------------|
| 1. Ethylene oxide | 7. Hydroxy acetate |
| 2. Ethyl formate | 8. Acetic acid |
| 3. Ethyl alcohol | 9. Formic acid |
| 4. Water | 10. Ethylene glycol/monoformate |
| 5. 2-Ethoxyethanol | 11. Ethylene glycol/monoacetate |
| 6. 2-Ethoxyethyl acetate | |



Column: HP-FFAP, (Cross-linked polyethylene glycol) 30 m x 0.53 mm x 1.0 µm (Part no. 19095F-123)
Carrier: Helium at 10 ml/min
Oven: Temperature program listed below
Injection: Split (10:1), 220°C
Detector: MSD (280°C)

HP-Basic Wax

Similar Phases

CAM, Carbowax Amine, Stabilwax-DB, CP-51 WAX for amines, CP-WAX for amines and diamines

Applications

Amines and other basic compounds

Temperature Range isothermal/ temperature programmed


40°C to 240/260°C

40°C to 220/240°C for 0.53 mm id

HP-Basic Wax for Amines and Other Basic Compounds

The most recent addition to the Agilent line of PEG capillary columns is the HP-Basic Wax column. This phase has been specifically developed for the analysis of amines and other basic compounds. Peak shapes for amines show excellent symmetry from the very first injection so there is no need to “prime” the column with sample before quantitative results can be obtained.

Free Fatty Acids



Column: HP-FFAP 30 m x 0.53 mm x 1.0 μ m
(Part number 19095F-123)
Carrier: Helium, 49 cm/sec, 6 ml/min
constant flow
Oven: 120°C (1 min) to 240°C at 10°C/min
Injection: Split (15:1), 1 μ l, inlet = 240°C
Detector: FID (265°C)
Sample: 1 to 2% acids in acetone

The HP-Basic Wax column can be used from 60°C to 260°C; one of the highest in the industry so you can extend the range of compounds you can analyze. Best of all, accurate quantitation is possible from the first

injection and the performance of the column does not deteriorate even after many injections. The column can be solvent rinsed (water and methane are not recommended) to extend the lifetime of the column.


PEG Columns for Foods and Flavors Applications

Analysis of Amines

Column: HP-Basic Wax, 30 m x 0.32 mm x 0.25 μ m
(Part number 19091N-613)
Carrier: Hydrogen, 6.0 psi
Oven: 120°C
Injector: 250°C, Split ratio: 100:1
Detector: 250°C, FID


Test Mixture
(0.5 mg/ml in Hexane)

- | | |
|----------------------|-------------------------|
| 1. n-Octylamine | 6. n-Heptadecane |
| 2. n-Nonylamine | 7. Trihexylamine |
| 3. n-Decylamine | 8. n-octyldecane |
| 4. Benzylamine | 9. 2,6-Dimethylaniline |
| 5. Dicyclohexylamine | 10. 2,4-Dimethylaniline |



Flavors/Fragrance Standard Mix


- | | | | |
|--------------------------|-------------------------|---------------------------|----------------------------------|
| 1. 1-Limonene | 11. Benzyl aldehyde | 20. Phenyl propyl acetate | 29. Heliotropine |
| 2. cis-3-Hexenal | 12. Linalool | 21. Maltol | 30. Helional |
| 3. Prenyl acetate | 13. Octanol | 22. Cyclamen aldehyde | 31. Idole |
| 4. cis-3-Hexenyl acetate | 14. Menthol | 23. Cinnamic aldehyde | 32. Methyl naphthyl ketone |
| 5. Rose oxide | 15. Citronellyl acetate | 24. Isopropyl myristate | 33. Musk xylol |
| 6. Hexenol | 16. alpha-Terpineol | 25. Unknown | 34. Vanillin |
| 7. Nonanal | 17. Benzyl acetate | 26. Dimethyl anthranilate | 35. Ethylene brassylate (Musk T) |
| 8. p-Methyl cresol | 18. Citronellol | 27. p-Cresol | |
| 9. Decanal | 19. Geraniol | 28. Eugenol | |
| 10. Camphor | | | |



Column: HP-INNOWax, 30 m x 0.25 mm x 0.25 μ m (Part number 19091N-133)
Carrier: Helium, 30 cm/sec, 0.9 ml/min, Constant Flow
Oven: 80°C (1 min), 5°C/min to 250°C (2 min)
Injection: 0.2 μ l, On-column
Detector: MSD, 280°C
Sample: 500 ppm in ethanol


PEG Columns for Foods and Flavors Applications

Vanillin Mix




Brandy

HP 5890 Series II with EPC and FID



Peppermint Oil

- | | |
|----------------|--------------------|
| 1. Limonene | 5. d-sio-Menthone |
| 2. Cineole | 6. Menthyl acetate |
| 3. Menthone | 7. Menthol |
| 4. Menthofuran | 8. Germacrene |




Column: HP-INNOWax (cross-linked PEG)
30 m x 0.32 mm x 0.5 µm (Part number 19091N-213)
Carrier: Helium, 40 cm/sec, 11.7 psi (60°C), 2.5 ml/min constant flow
Oven: 60°C (1 min) to 250°C at 10°C/min
Injection: Split (60:1), 0.5 µl, inlet at 220°C
Detector: FID (275°C)
Sample: Neat

Wine, Cabernet

HP 5890 Series II with EPC and FID

GC System: Agilent GC/FID
Automated Sampler: 7673 Automated Sampler
Column: HP-INNOWax, 30 m x 0.25 mm x 0.25 µm (Part number 19091N-133)
Carrier: He, 33 cm/sec, 15.5 psi (35°C), 1.5 ml/min constant flow
Oven: 35°C (5 min) to 150°C at 5°C/min, to 250°C (2 min) at 20°C/min
Injection: Split (split ratio 25:1), 1 µl, inlet temp at 220°C
Detector: FID (280°C)
Sample: Neat



HP-20M

Similar Phases

Carbowax 20M, 007-CW, BP-20

Applications

Alcohols, Free Acids, Ethers, Glycols

Nonbonded HP-20M Columns


The HP-20M is the nonbonded wax phase that was designed to support established methods that require the polarity be identical to Carbowax 20M. The columns have a temperature range of 60°C to 220°C.

PEG Columns for Paramaceuticals Applications

Alcohols, Aldehydes, and Ketones

Calibration standards

1. Acetaldehyde
2. Acetone
3. Methanol
4. Isopropanol
5. Ethanol
6. n-Propanol IS



Column: HP-INNOWax, 15 m x 0.25 mm x 0.5 µm (Part number 19091N-231)


Detector: FID, 275°C

Carrier: Helium, 1.2 ml/min, Constant Flow

Oven: 40°C (4 min), 10°C/min to 60°C; 25°C/min to 110°C

Injection: 0.5 µl, 250°C, Split 250:1

Analysis of Ethylene Oxide in a Pharmaceutical Formulation



Inlet: Split/10:1, 1 µl, 250°C
Column: HP-WAX, 60 m x 0.53 mm i.d. x 1 µm (Part number 19095X-126)

Carrier Gas: Hydrogen


Head Pressure: 25 kPa at 50°C—constant flow mode

Flow: 5.2 ml/min

Oven Temperature: 50°C, 5°C/min to 75°C, 25°C/min to 200°C

GC System: Agilent GC/FID, 250°C

Volatile Hydrocarbons and Ketones in Blood



Headspace

sampler: Agilent 7694

Column: HP-INNOWax, 30 m x 0.53 mm x 1.0 µm (Part number 19095N-123)

Inlet: Split/2:1 (180°C)

GC System: Agilent GC/FID 220°C

GC Conditions

Oven temperature: 40°C (10 min), 5°C/min to 110°C

Carrier gas

flow rate: 3.4 ml/min (25 cm/sec), constant flow

Headspace Conditions

Oven temperature: 90°C

Valve and

loop temp.: 92°C

Transfer line

temp.: 110°C

Incubation time: 30 min

Shaking: High


Injection timing: 0.5 min pres; 0.02 min vent; 1.0 min inject

Vial pressure: 10.8 psi

Loop size: 1 ml

Sample: 1.0 ml blood and 100 µl internal standard (0.01% n-butanol in deionized water) to each 22 ml headspace vial.

Organic Volatile Impurities in Pharmaceuticals



1. Methylene chloride
2. Chloroform
3. Benzene
4. Trichloroethylene
5. 1,4-Dioxane

Equipment

Column: HP-INNOWax, 30 m x 0.53 mm x 1.0 µm (Part number 19095N-123)

Detector: FID, 300°C

Sample introduction: Agilent 7694 Headspace Sampler

Experimental

Carrier gas: He, 40 cm/sec, constant flow

Oven: 50°C (hold 3.5 min) to 60°C at 10°C/min, to 160°C at 40°C/min

Injection: Splitless, 1 µl sample volume, 175°C

Sample: 100 mg pharmaceutical with 50 ppm CHCl_3 ; 100 ppm remaining components

Headspace Sampler Conditions

Setpoint:

Oven—80°C Transfer line—110°C Loop—90°C

Time events (min):

Vent time—0.50 Loop equilibrium time—0.05

Loop fill time—0.15

GC cycle time—10 Pressurize time—0.20

Vial equilibrium time—10 to 60 min with vigorous agitation

Agilent PEG Ordering Information**HP-Wax, HP-INNOWax, and HP-Basic Wax Columns**

Description*	HP-Wax Part No.	HP-INNOWax Part No.**	HP-Basic Wax Part No.
10 m x 0.1 mm x 0.1 µm	19091X-141		
10 m x 0.1 mm x 0.2 µm	19091X-241		
20 m x 0.1 mm x 0.1 µm	19091X-142		
20 m x 0.1 mm x 0.2 µm	19091X-242		
40 m x 0.1 mm x 0.2 µm	19091X-244		
25 m x 0.2 mm x 0.2 µm	19091X-102	19091N-102	
25 m x 0.2 mm x 0.4 µm	19091X-202	19091N-202	
50 m x 0.2 mm x 0.2 µm	19091X-105	19091N-105	
50 m x 0.2 mm x 0.4 µm	19091X-205	19091N-205	
15 m x 0.25 mm x 0.15 µm	19091X-031	19091N-031	
15 m x 0.25 mm x 0.25 µm	19091X-131	19091N-131	19091N-631
15 m x 0.25 mm x 0.5 µm	19091X-231	19091N-231	19091N-731
30 m x 0.25 mm x 0.15 µm	19091X-033	19091N-033	
30 m x 0.25 mm x 0.25 µm	19091X-133	19091N-133	19091N-633
30 m x 0.25 mm x 0.5 µm	19091X-233	19091N-233	19091N-733
60 m x 0.25 mm x 0.15 µm	19091X-036	19091N-036	
60 m x 0.25 mm x 0.25 µm	19091X-136	19091N-136	
60 m x 0.25 mm x 0.5 µm	19091X-236	19091N-236	
15 m x 0.32 mm x 0.15 µm	19091X-011	19091N-011	
15 m x 0.32 mm x 0.25 µm	19091X-111	19091N-111	19091N-611
15 m x 0.32 mm x 0.5 µm	19091X-211	19091N-211	19091N-711
30 m x 0.32 mm x 0.15 µm	19091X-013	19091N-013	
30 m x 0.32 mm x 0.25 µm	19091X-113	19091N-133	19091N-613
30 m x 0.32 mm x 0.5 µm	19091X-213	19091N-213	19091N-713
60 m x 0.32 mm x 0.15 µm	19091X-016	19091N-016	
60 m x 0.32 mm x 0.25 µm	19091X-116	19091N-116	
60 m x 0.32 mm x 0.5 µm	19091X-216	19091N-216	
15 m x 0.53 mm x 1.0 µm	19095X-121	19095N-121	19095N-621
30 m x 0.53 mm x 1.0 µm	19095X-123	19095N-123	19095N-623
60 m x 0.53 mm x 1.0 µm	19095X-126	19095N-126	

* Length x Internal Diameter x Film Thickness

** INNOWax™ bondable

HP-FFAP Columns

Description*	HP-FFAP Part No.
25 m x 0.2 mm x 0.3 µm	19091F-102
50 m x 0.2 mm x 0.3 µm	19091F-105
15 m x 0.25 mm x 0.25 µm	19091F-431
30 m x 0.25 mm x 0.25 µm	19091F-433
15 m x 0.32 mm x 0.25 µm	19091F-411
25 m x 0.32 mm x 0.50 µm	19091F-112
30 m x 0.32 mm x 0.25 µm	19091F-413
50 m x 0.32 mm x 0.50 µm	19091F-115
10 m x 0.53 mm x 1.0 µm	19095F-121
15 m x 0.53 mm x 1.0 µm	19095F-120
30 m x 0.53 mm x 1.0 µm	19095F-123

* Length x Internal Diameter x Film Thickness

HP-101 and HP-20M Capillary Columns

Description*	HP-101 Part No.	HP-20M Part No.
25 m x 0.2 mm x 0.1 µm		19091W-102
50 m x 0.2 mm x 0.1 µm		19091W-105
25 m x 0.32 mm x 0.3 µm	19091Y-012	19091W-012
50 m x 0.32 mm x 0.3 µm	19091Y-015	19091W-015
10 m x 0.53 mm x 1.33 µm		19095W-121
30 m x 0.53 mm x 1.33 µm		19095W-123

* Length x Internal Diameter x Film Thickness

So when you search the polar horizon, squinting to distinguish where the snow line begins and the sky ends, be sure you have Agilent PEG columns to help you find the answers you need quickly and reliably.

Column polarity, ruggedness, and lifetime second to none sets Agilent PEG columns apart from the competition.

Copyright © 1998, 2000
Agilent TechnologiesPrinted in USA 4/00
5968-0322E**Agilent Technologies**
Innovating the HP Way