

Analysis of a Brenzcatechol Additive in Styrene using HPLC

Angelika Gratzfeld-Huesgen

Polymer/ chemical industry

Abstract

Brenzcatechol (TBC) or Benzene-1,2-diol is used as an antioxidant for polymers. The concentration range of Brenzcatechol used varies between 10 and 100 ppm.

Method Performance

Figure 1 shows the HPLC chromatogram of the analyzed styrene sample and the standard chromatogram of Brenzcatechol. For additional identification purposes, spectra can be taken and a comparison with the standard spectrum can be made. For this application, 2.1 mm columns were used in order to improve sensitivity so that the detection of 1 ng with signal to noise of 2, could be achieved.



Conditions

Column 200 x 2.1 mm Hypersil ODS, 5 μm **Mobile Phase** A = Water, B = Acetonitrile **Gradient** at start 50 % B, at 10 min 99.9 %B, at 20 min 50 %B **Post Time** 6 min **Flow Rate** 0.5 ml/min **Oven Temp** 40 °C **Injection Vol** 1 μl **Diode array detector** 280/30 nm; Reference 500/50 nm

Sample preparation 1 ml styrene sample was diluted with 1 ml Tetrahydrofurane (THF)

Figure 1 Analysis of brenzcatechine (TBC) additive in styrene



Agilent Technologies



Figure 2 Brenzcatechin (Benzene-1,2-diol)

Method performance

LOD: 1 ng or 1 ppm with signal/noise = 2 red RT <0.2 % rsd area <2 %

Equipment

Agilent 1100 Series

- degasser
- binary pump
- autosampler
- thermostatted column compartment
- diode array detector Agilent ChemStation + software



Angelika Gratzfeld-Huesgen is application chemist at Agilent Technologies, Waldbronn, Germany.

For more information on our products and services, visit our worldwide website at http://www.agilent.com/chem

© Copyright 1997 Agilent Technologies Released 06/97 Publication Number 5965-9056E



Agilent Technologies Innovating the HP Way