

# Analysis of Acidulants in White Wine using HPLC

Angelika Gratzfeld-Huesgen

Food

#### Abstract

Sorbic acid and citric acids are commonly used as acidulants<sup>1</sup> and/or as preservatives. Acetic, propionic, succinic, adipic, lactic, fumaric, malic, tartaric, and phosphoric acids can serve as acidulants as well. Acidulants are used for various purposes in modern food processing. For example, citric acid adds a fresh, acidic flavor, whereas succinic acid gives food a more salty, bitter taste. In addition to rendering foods more palatable and stimulating, acidulants act as

- flavoring agents to intensify certain tastes and mask undesirable aftertastes
- buffering agents to control the pH during food processing and of the finished products
- preservatives to prevent growth of microorganisms
- synergists to antioxidants to prevent rancidity and browning
- · viscosity modifiers in baked goods
- · melting modifiers in cheese spreads and hard candy
- · meat curing agents to enhance color and flavor



Figure 1 Analysis of acidulants in white wine

# Conditions

#### Column

300 <sup>~</sup> 7.8 mm BioRad HPX 87-H, 9 μm **Mobile phase** 0.0035M H<sub>2</sub>SO<sub>4</sub> isocratic **Flow rate** 0.6 ml/min **Column compartment** 65 °C

Injection vol 10 µl Detector

UV-VWD detection wavelength 192 nm or 210 nm

Sample preparation Filtration



Agilent Technologies Innovating the HP Way

#### Sample preparation

Sample preparation depends strongly on the matrix to be analyzed, but in general steam distillation and solid-phase extraction techniques can be used.

### Chromatographic conditions

High-performance liquid chromatography (HPLC) with UV-visible diode-array detection (UV-DAD) has been applied in the analysis of citric acid in wine and in a vodka mixed drink. Retention time and spectral data were used as identification tools.

#### HPLC method performance

Limit of detection

100ng injected amount, S/N = 2 equivalent to 2 ppm with 50  $\mu$ l injected volume

#### Repeatability of

RT over 10 runs <0.1% areas over 10 runs <3 %



Figure 2 Analysis of citric acid in vodka

# References

1.

Official Methods of Analysis, Food Compositions; Additives, Natural Contaminants, 15th ed; AOAC: Arlington, VA, 1990, Vol. 2.; Official Method AOAC 986.13: quinic, malic, citric acid in cranberry juice cocktail and apple juice.



Agilent Technologies

Innovating the HP Way

# Conditions

#### Sample preparation filtration Column 300 ~ 7.8 mm BioRad HPX 87-H, 9 μm Mobile phase 0.007M H<sub>2</sub>SO<sub>4</sub> isocratic Flow rate 0.6 ml/min Column compartment 65 °C Injection vol 10 μl Detector UV-DAD

# Equipment

# Agilent 1100 Series

- degasser
- isocratic pump
- autosampler
- thermostatted column compartment
- diode array detector, variable wavelength detector

or refractive index Agilent ChemStation + software

Angelika Gratzfeld-Huesgen is application chemist at Agilent Technologies, Waldbronn, Germany.

For more information on our products and services, visit our worldwide website at http://www.agilent.com/chem

© Copyright 1997 Agilent Technologies Released 09/97 Publication Number 5966-0627E