

Analysis of Amitriptyline, Fluoxetine, Quetiapine, and Sertraline Using Dried Blood Spots

Application Note

Pharmaceutical

Authors

William Hudson, Ben Yong, and
Paul Boguszewski
Agilent Technologies, Inc.

Introduction

Dried blood spotting (DBS) has gained in popularity with pharmaceutical laboratories for pre-clinical and clinical trials for new chemical entities in the pharmaceutical industry. Fifteen microliters of blood is spotted onto a blood spotting card and allowed to dry for at least 2 hours. A 2–4 mm core is then punched from the spot. Analytes are desorbed from the punched spot using solvent, then analyzed by LC-MS/MS. Cellulose materials have typically been used for these applications, but a new non-cellulose material is now available. This Dried Matrix Spotting (DMS) media can also be used for other biological fluids such as plasma, urine, and spinal fluid.



Agilent Technologies

Experimental

Four analytes were chosen: quetiapine, amitriptyline, sertraline, and fluoxetine. Fluoxetine-D6 was readily available as an internal standard, so it was used for all compounds. Precision and accuracy recoveries were calculated based on the linear regression of the calibration standards. The accuracy was calculated at a low level (0.5 ng/mL), a mid-level (5.0 ng/mL), and high level (500 ng/mL).





LC/MS conditions

Column	Agilent Poroshell 120 EC-C18 3.0 mm x 50 mm 2.7 μ m
Mobile phase	A: 0.1% Aqueous formic acid B: ACN
Pump program	Flow rate 300 μ L/ min
	t_0 A: 50%, B: 50%
	$t_{1.5-2.0}$ A: 10%, B: 90%
	$t_{2.01-300}$ A: 50%, B: 50%
Run time	3:00 minutes
Gas temp	350 $^{\circ}$ C
Gas flow	10 L/min
Nebulizing	20 psi
Pol	Pos

Compound	Q1 ion	Product ion	CE
Quetiapine	384.0	253.1	18 V
Amitriptyline	278.2	105.1	22 V
Sertraline	306.1	159.01	26 V
Fluoxetine	310.1	163.0	26 V
Fluoxetine-D6	316.1	154.1	2 V

Results and Discussion

1 mL of human blood was spiked with 10 μ L of each working standard to create a calibration curve of 0.1, 0.2, 0.5, 1.0, 2.0, 5.0, 20, 50, 200, 500 and 2,000 ng/mL. Fifteen μ L of each of these standards was spotted onto an Agilent DMS card.



A 3 mm disk was punched from each dried spot and placed into a 96-well collection plate.

300 μ L of 0.1% formic acid in 80% methanol (with 0.066 ng/mL of deuterated internal standard mix) was added to each well and vortexed.

The samples were evaporated to dryness and reconstituted in 100 μ L of mobile phase.

Analyte Recoveries (n=6)

	0.5 ng/mL		5.0 ng/mL		500 ng/mL	
	% Rec	RSD	% Rec	RSD	% Rec	RSD
Quetiapine	91%	12%	87%	7%	99%	8%
Amitriptyline	90%	7%	95%	13%	92%	2%
Fluoxetine	105%	10%	104%	16%	101%	2%
Sertraline	90%	11%	100%	5%	97%	2%





Conclusions

Four compounds in blood were successfully desorbed and analyzed by LC-MS/MS. Good detection levels were achieved using an Agilent 1290 LC system and an Agilent 6460 mass spectrometer. Linearity was demonstrated using a 1st order regression and correlation coefficients were better than 0.995. Relative recoveries were within 10% of the true value with RSD values all less than 10%. DMS proved to be an effective means of storing plasma samples and protein removal for a simple approach to sample preparation prior to LC-MS/MS.

For More Information

Bond Elut DMS cards are intended for use in DMPK/ADME research applications only. They should not be used in diagnostic procedures.

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.



www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011
Printed in the USA
June 30, 2011
5990-8168EN



Agilent Technologies