



# Analysis of catechins in tea by HPLC with electrochemical detector

Kuniko Koizumi and  
Hiroki Kumagai

Food

## Abstract

Catechin (C), epicatechin (EC), epigallocatechin (EGC) and epigallocatechingallate (EGCG)—collectively called catechins, have many phenolic hydroxyl groups in the molecule. It is known that catechins have various physiological benefits such as anti-oxidation, anti-bacterium, and anti-tumor effects, as well as the suppression of the increase of cholesterol concentration in blood. Catechins also effect the flavor of foods such as tea, wine, and beer.

This brief demonstrates the analysis of catechins in tea using the Agilent 1100 Series modules and systems for LC with an electrochemical detector.

## Analyzed Compounds

(+)-catechin,  
(-)-epicatechin  
(-)-epigallocatechin  
(-)-epigallocatechingallate

## Sample

Green tea, tea and wulong tea




Figure 1  
Chromatogram of standard solution, 2 mg/L each

## Conditions

### Column

250 x 4.0 mm Hypersil BDS  
(Agilent Parts No: 79926 OB-584)

### Mobile phase

Methanol: (2g NaNO<sub>3</sub> + 0.05 g  
H<sub>2</sub>SO<sub>4</sub> in 1l Water) = 22:78

**Column temp** 25 °C

**Injection vol** 20 µl

### Detector

Agilent 1049 Electrochemical detector

### Mode

Amperometry Working electrode;  
glassy carbon. Applied Potential: 0.850 V

## Sample preparation

Green tea, tea and wulong tea were  
diluted with water and filtrated.



**Agilent Technologies**

Innovating the HP Way

## Method Performance

Limit of detection  
0.26–2.00 ng (S/N = 3)

RSD of peak area  
0.73–1.14 %

RSD of retention time  
0.09–0.14 %





Figure 2  
Chromatogram of catechins in various teas

## Equipment

### Agilent 1100 Series

- binary pump with
  - vacuum degasser
  - autosampler
  - thermostatted column compartment
  - programmable electrochemical detector
  - diode array detector
- Agilent ChemStation + software



Kuniko Koizumi and Hiroki Kumagai are application engineers based at Yokogawa Analytical Systems Inc., Tokyo, Japan.

For more information on our products and services, visit our worldwide website at <http://www.agilent.com/chem>

© Copyright 1997 Agilent Technologies  
Released 10/97  
Publication Number 5965-9802E



**Agilent Technologies**

Innovating the HP Way