



Analysis of Polybutadiene using Gel Permeation Chromatography

Angelika
Gratzfeld-Huesgen

Polymer
/chemical industry

Abstract

Polybutadiene is used as homo- and co-polymers for the production of tyres (70 %), moving belts and soles of shoes. In 1989 1.1 Mio.t. were used worldwide.

For example, the performance of tyres depends strongly on the molecular weight (MW) of the polybutadiene used and its additives. To ensure highest quality and consequently highest safety MW data need to be evaluated for each batch of produced polymer. Gel Permeation Chromatography (GPC) is an analytical tool used to characterize polymers which are soluble in organic solvents. In general an isocratic pump is sufficient for GPC analysis, however for ease of solvent change and rinsing, a pump with two or more channels would be advantageous. The pump should be able to pump the selected flow rate with a precision of typically < 0.15 %. Solvent degassing is recommended either offline or even better online with vacuum degassing. For high sample throughput the use of an autosampler would be beneficial. The temperature of the column oven needs to be very stable to avoid retention time shift and therefore MW errors—a peltier controlled thermostat is ideal for highest temperature stability especially at and below ambient temperatures

For detection a UV detector and/or a refractive index detector can be used. In this example we used both detectors to demonstrate, that results can be quite different, depending on the

detection system used. The use of a diode array system enables the taking of spectra as an additional identification tool. This can help to identify for example remaining monomers (figure 1).




Figure 1
Spectra of oligomers and monomers

Method Performance

In figure 2 the different signal traces of UV-DAD and refractive index detection are shown. It can be seen that the calculated MW data differ significantly. In addition to MW data like Mw and Mn, GPC evaluation software also calculates molecular weight distribution curves (MWD),



Agilent Technologies
Innovating the HP Way

which offer information about the relation between for example height percentage and log molecular weight or cumulative height percentage versus log molecular weight. (figure 3).

Method performance

Precision of weight: average molecular weight (rsd of Mw) = < 1 %

Precision of number weight average: molecular weight (rsd of Mn) = < 2 %




Figure 2
Analysis of polybutadiene using UV and refractive index




Figure 3
Molecular weight data based on refractive index detection

| MW data | Refractive index | UV Detection |
|----------------|------------------|--------------|
| Mp | 10283 | 10000 |
| Mn | 10543 | 6567 |
| Mw | 12054 | 13565 |
| Polydispersity | 1.143 | 2.066 |
| Mz | 14804 | 22037 |
| Mz + 1 | 21860 | 35784 |
| Mv | 11780 | 12579 |

Table 1
Molecular weight data refractive index versus UV detection

Conditions

Column 2 × PLgel mixed-D,
7.5 × 300 mm, 5 µm

Mobile phase

Tetrahydrofurane (THF)

Flow rate 1 ml/min

Oven Temp 20 °C

Injection vol 20 µl

UV detector DAD 254/100 nm,
reference 360/100 nm

Refractive index detector

Sample preparation

33 mg sample dissolved in 1 ml
THF; Polystyrene standards from
PSS were used for narrow
standard calibration

Equipment

Agilent1100 Series:

- isocratic pump
 - degasser (recommended)
 - autosampler
 - thermostatted column compartment
 - diode array detector and/or HP 1047A refractive index detector
- Agilent ChemStation +
software + polymer
labs GPC software

Angelika Gratzfeld-Huesgen is
application chemist at Agilent
Technologies, Waldbronn,
Germany.

For more information on our
products and services, visit
our worldwide website at
<http://www.agilent.com/chem>

© Copyright 1997 Agilent Technologies
Released 06/97
Publication Number 5965-9046E



Agilent Technologies

Innovating the HP Way