



Analysis of Quercetin and Kaempferol in Gingko Extract and Tablets (*Gingko Biloba*) by HPLC

Udo Huber

Pharmaceutical

The extract of *Gingko biloba* contains bilobaloides, ginkgoloides and other glucosides of flavones, for example, quercetin, isoquercetin, kaempferol, luteolin and sitosterole glucosides. The extract increases blood circulation and energy metabolism of the brain. It also increases the blood sugar and shows some vasodilatory effects. The ginko extract or tablets are used in treating blood circulation problems of the brain, arms and legs.

Figure 1 shows the separation of quercetin and kaempferol in the extract and tablets of *Gingko biloba* using gradient analysis on a reversed phase column and UV detection. The autosampler temperature was set to 4 °C to avoid decomposition of the samples.




Figure 1
Analysis of quercetin and kaempferol

Conditions

Column

4 x 125 mm Hypersil ODS, 5 µm

Mobile phase

A = 0.5 % H₃PO₄ in water,
B = methanol

Flow rate

2.0 ml/min

Gradient

at 0 min 38 % B
at 12 min 48 % B

Column wash

at 17 min 100 % B
at 20 min 38 % B

UV detector

variable wavelength detector
370 nm, standard cell

Column compartment temperature

25 °C

Stop time

20 min

Post time

5 min

Injection volume

10 µl



Agilent Technologies

Innovating the HP Way

Extraction

4 g of *Ginkgo biloba* extract (from Caesar & Loretz GmbH, Germany) were refluxed for 30 min in 70 ml methanol and 10 ml 25 % HCl. After cooling to room temperature the mixture was filtered and the filter washed with approximately 100 ml methanol. The solvent was partly removed *i. vac.* and diluted with methanol to 100 ml in a volumetric flask. 5 ml of this solution were filtered through a C18 disposable cartridge. The cartridge was washed with 4 ml methanol and the filtrate diluted to 10 ml in a volumetric flask. The same procedure was used to extract 10 ginkgo tablets ('Promod Ginkgo biloba L.', 40 mg *Ginkgo biloba* extract per tablet, Sine Laboratories', Shanghai, China)

Method and extraction from:

A. Hasler, O. Stichler, *J. Chromatogr.* 508 (1990), 236-240




Figure 2
Comparison of sample and standard spectra

The HPLC method presented shows an easy but reliable and precise analysis of quercetin and kaempferol in the extract and tablets of *Ginkgo biloba*.

Equipment

Agilent 1100 Series

- Quaternary pump (includes vacuum degasser)
- Thermostatted autosampler
- Thermostatted column compartment
- Variable wavelength detector, standard flow cell 10-mm path length, 13- μ l cell volume

Alternative:

- Binary pump
- Vacuum degasser
- Diode array detector standard flow cell 10-mm path length, 13- μ l cell volume
- Agilent ChemStation + 3D software

Columns

- Hypersil ODS, 5 μ m, 4 x 125 mm (Agilent part number 7992618-564)
- *Recommended:* Guard cartridges Hypersil ODS, 5 μ m, 4 x 4 mm (Agilent part number 7992618-504, 10/pk)

Note:

Since the method was specifically developed on the Agilent 1100 Series system you might not be able to reproduce this analysis on an older system or even on a new system with lower performance. To avoid sample decomposition it is necessary to use a cooled autosampler, for example, the Agilent 1100 Series thermostatted autosampler.

© Copyright 1998 Agilent Technologies
Released 12/98
Publication Number 5968-2973E



Agilent Technologies

Innovating the HP Way