



Analysis of Testosterone Metabolites using HPLC

Noriko Shimoi
Hiroki Kumagai

Biopharmaceutical

Abstract

Cytochrome P450 (CYP) is a family of haemoproteins and oxygenase. CYP plays an important role in the detoxification of drugs, mutagen, and carcinogen. They exist in multiple forms (or isoenzymes) which show broad substrate specificity and differing catalytic activity. The CYP isoenzyme shows the stereospecific hydroxylation of steroid nucleus, therefore, the metabolism of testosterone was used as a probe for CYP isoenzyme activities.

In this Application Brief we demonstrate the highly precise analysis of testosterone and its metabolites.




Figure 1
Chromatogram of testosterone and its metabolites
(standard solution - each peak contained 20 ng steroid)

Conditions

Column:

150 x 4.6 mm, 3 μ m
Cacellpak C₁₈ UG120 (Shiseido)

Mobile phase:

A = water

B = CH₃OH, C, CH₃CN

Column compartment:

40 °C

Injection vol.:

25 μ l

Detector:

diode-array detector

wavelength 245 nm

reference 450 nm



Agilent Technologies

Innovating the HP Way

Analyzed compounds

The following compounds were analyzed in this application:
testosterone (TS), androstenedione (AD), methyltestosterone (MTS),
2 α -6 β -7 α -, 16 β - hydroxytestosterone (2 α -, 6 β -, 7 α -, 16 α -,16 β - OHT)




Figure 2
Overlay of 5 consecutive runs

Equipment

Agilent 1100 Series

- vacuum degasser
- quaternary pump
- autosampler
- thermostatted column compartment
- diode array detector, Agilent ChemStation + software

HPLC method performance

Limit of detection:
0.18 - 1.33 ng with S/N = 3
RSD of peak area:
0.11 - 0.78 %
RSD of retention time:
0.04-0.11 %

The authors are application chemists at Yokogawa Analytical Inc., Tokyo, Japan.

For more information on our products and services, visit our worldwide website at <http://www.agilent.com/chem>

© Copyright 1998 Agilent Technologies
Released 04/98
Publication Number 5967-5596E



Agilent Technologies

Innovating the HP Way