



# Analysis of Triglycerides in Olive Oil and Rap Oil using HPLC

Angelika  
Gratzfeld-Huesgen

Food

## Abstract

Unsaturated triglycerides in olive oil have very characteristic patterns. Other fats and oils are also complex mixtures of triglycerides but with different patterns.

## Sample preparation

Triglycerides can be extracted from homogenized samples with petrol ether. Fats and oils can be dissolved in tetrahydrofurane.

## Chromatographic conditions

The presented HPLC method was used to analyze the unsaturated triglycerides, LnLnLn, LLL, and OOO.<sup>1</sup>




Figure 1  
Analysis of the triglyceride pattern of olive and rape oil

## Conditions

### Column

200 ~ 2.1 mm Hypersil MOS, 5 µm

### Mobile phase

acetone/ACN (30:70)

### Flow rate

0.5 ml/min

### Column compartment

30 °C

### Injection vol

2 µl

### Detector

refractive index

### Sample preparation

Samples were dissolved in tetrahydrofurane.



**Agilent Technologies**

Innovating the HP Way

## HPLC method performance

Limit of detection for ECD  
50 µg/l with S/N = 2

Repeatability of  
RT over 10 runs <0.3 %  
areas over 10 runs 5 %


## References

1.  
"Determination of triglycerides in vegetable oils", EC Regulation  
No. L248, 28ff.

## Equipment

### Agilent 1100 Series

- degasser
  - isocratic pump
  - autosampler
  - thermostatted column compartment
  - refractive index detector
- Agilent ChemStation + software



Angelika Gratzfeld-Huesgen is application chemist at Agilent Technologies, Waldbronn, Germany.

For more information on our products and services, visit our worldwide website at <http://www.agilent.com/chem>

© Copyright 1997 Agilent Technologies  
Released 09/97  
Publication Number 5966-0634E



**Agilent Technologies**

Innovating the HP Way