

Analyze Injection-Molding Polymers on Agilent PLgel 5 µm MIXED-C by GPC

Application Note

Materials Testing and Research

Introduction

Injection-molding thermoplastics are used in a wide variety of applications, including panels for plastic doors, casings for scientific instruments, refrigerator linings, and telephone cases. Ideally, the materials should display excellent weather and impact resistance, good color retention and consistency, resistance to environmental stress, cracking, and high heat resistance. A number of copolymers have been developed that display these properties, such as copolymers of acrylonitrile/styrene/acrylate (ASA) and acrylonitrile/butadiene/styrene (ABS). Typically, these materials are either layered or blended to meet the requirements of a specific application.

This application note describes the analysis of a sample of an ASA/ABS-blended material from a car door panel by gel permeation chromatography (GPC), using Agilent PLgel 5 µm MIXED-C columns.

Injection-Molding Polymer Analysis

Figure 1 shows the blended material that eluted as a single monomodal peak, indicating that the two components of the blend were very similar in molecular size, with an Mw of 160,000 g/mol and a polydispersity of 2.6.



Agilent Technologies

Authors

Greg Saunders and Ben MacCreath Agilent Technologies (UK) Ltd Essex Rd **Church Stretton** SY6 6AX UK

Conditions

Samples	Blend of ASA and ABS
Columns	2 × Agilent PLgel 5 μm MIXED-C, 7.5 × 300 mm (p/n PL1110-6500)
Eluent	THF (stabilized)
Flow rate	1.0 mL/min
lnj vol	100 μL
Detector	RI
System	Agilent PL-GPC 50

Conclusions

Gel permeation chromatography with Agilent PLgel 5 μ m MIXED-C columns permits the molecular weight distribution of polymers to be assessed, and the size of the molecular components in a polymer blend to be investigated.

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.



Figure 1. Individual components of an injection-molding polymer elute as a single peak, showing their similarity in molecular size.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011 Printed in the USA July 8, 2011 5990-8601EN

