

ASAPHMQC – Quick **Multidimensional Spectra**

Application Note

Author

Jarrett Farias, Ph.D. NMR Applications Scientist **Research Products Group** Santa Clara, CA USA

Abstract

VnmrJ 3 software provides easy-to-use, interactive tools for setting up advanced experiments. Simplified set-up procedures allow even novice users to get critical information about their research samples using the most advanced NMR experiments available. This application note is one of a series designed to provide step-by-step guidance for setting up sophisticated experiments to collect exactly the data you need for your analyses.

Introduction

The Acceleration by Sharing Adjacent Polarization Heteronuclear Multiple Quantum Coherence (ASAPHMQC)¹ experiment, similar to the commonly used Heteronuclear Single Quantum Coherence (HSQC)² experiment, provides information regarding one-bond (¹H-¹³C) connectivity and heteronuclear chemical shifts for the protonated carbons. ASAPHMQC, however, requires far less time than most two-dimensional (2D) experiments, and, as such, provides an elegant tool for the collection of inspection or survey data.

In this application note, comparative ASAPHMQC and gradient heteronuclear single quantum coherence (adiabatic version) (gHSQCAD) data were acquired with default settings on 5 mg of clindamycin HCI. All 16 ¹H-¹³C correlations are clearly identifiable in the ASAPHMQC experiment in 56 seconds, whereas the gHSQCAD required 8.25 minutes to complete. Additionally, with a few minor parameter changes, a 10 mg sample of clindamycin HCl gives a complete ASAPHMQC spectrum in 19 seconds. These examples provide an illustration of the power of the ASAPHMQC experiment.



The efficiency gain for ASAPHMQC comes from a much shorter relaxation delay as compared to an HSQC-type experiment (0.06 seconds versus 1 second for the default). The loss of longitudinal magnetization from incomplete relaxation is compensated for by a magnetization transfer from nearby ¹²C-attached protons (which are kept along z during the pulse sequence) to the observed ¹³C-attached protons. The default parameters give a total experiment time of 56 seconds, not much longer than a simple proton spectrum.

One way to exploit the efficiency of ASAPHMOC is to use the experiment, along with a one-dimensional (1D) ¹H spectrum, as an efficient and elegant strategy for the investigation of a molecular structure at the initial stage. Typically, if the ¹H spectrum does not provide sufficient information, the sample is resubmitted with strategic homonuclear or heteronuclear correlation experiments. Adding heteronuclear information in the inspection stage can often solve the structure and, if not, more targeted data can be acquired, for example, selective 1D or band-selective 2D experiments. This approach to structure elucidation will often be much faster and minimize the impact to spectrometer usage.

Another use of the efficiency of ASAPHMQC is to obtain high-resolution (in F1) data when necessary. For 2D ¹H-¹³C correlations. the resolution in F2, determined by the acquisition time, is often limited by the carbon-decoupling power requirements. Resolution in F1, however, requires additional t1 increments (ni), which can be expensive in terms of time. Because of the rapid recycle time of the ASAPHMOC experiment, a high-resolution 2D proton carbon correlation (for example, ni = 2048) can be obtained in approximately 13 minutes for a sufficiently concentrated sample. This can be extremely useful for compounds with spectral overlap issues.

An ASAPHMQC example

A sample of clindamycin (Figure 1) was used to demonstrate the speed and data quality obtained using the ASAPHMQC experiment. All data was collected on an Agilent 400-MR DD2 instrument with a OneNMR probe.



Figure 1. The chemical structure of clindamycin.

The spectrum in Figure 2 is of a ¹H-¹³C ASAPHMQC experiment. Figure 3 shows a ¹H-¹³C gHSQCAD for comparison. Using the unmodified default parameters to set up both experiments, ASAPHMQC provided signals for all 16 protonated carbons with a significantly reduced experiment time, 56 seconds as compared to 8 minutes 18 seconds for the gHSQCAD. The intensities of correlations 12 and 20 are smaller in the ASAPHMOC experiment than other correlations due to the limited number of ¹H-¹³C pairs for magnetization transfer, and can be increased with a longer relaxation delay than the default of 60 milliseconds.



Figure 2. $^{1}\text{H-}^{13}\text{C}$ ASAPHMQC on 5 mg clindamycin dissolved in 600 μL DMSO-d_6. The total experiment time was 56 seconds.



Figure 3. 1 H- 13 C gHSQCAD on 5 mg clindamycin dissolved in 600 μ L DMSO-d₆. The total experiment time was 8 minutes 18 seconds.

Depending on the sensitivity of the probe and the amount of sample present, the default number of scans per increment could be decreased to ensure the rapid acquisition of an ASAPHMOC experiment with minimal impact to spectrometer time. In Figure 4, the concentration of clindamycin was doubled, so fewer repetitions are required to achieve an adequate signal-to-noise ratio. The relaxation delay was increased to 70 milliseconds to demonstrate the impact to correlations 12 and 20, and the number of repetitions per increment was decreased to one. All 16 protonated carbons are still clearly identifiable, with correlations 12 and 20 more intense than in Figure 2. The total time for this ASAPHMQC experiment was only 19 seconds, which is similar to the time needed to acquire a short 1D ¹H experiment.



Figure 4. ¹H-¹³C ASAPHMQC on 10 mg clindamycin dissolved in 600 μ L DMSO-d₆. The total experiment time was 19 s, using 56 increments, a spectral with of 100 ppm in F1, and a delay of 70 ms.

le Edit View Experiments .		cess Iools Help		
			۳ 🕓 🖌 🐂	TOCSY
Sample Changer set VTWAI				
port ProcessPlot ArrayedS	ipectra	p:5 Sea: Submita Index: 1		
Totocols QuickSubmit	Frame			
riment Selector	JL × +			
ethylindanone Common	Service			
Std 1D 🕶 0+	HHJHomo2D 🔻			
J1(CH)corr - HSQCAD	*p(CH)corr ▼			
(H)Sel1D - gHSQCAD HSQC	Sel2D 🖛			
(HC)Crisis2 - gHSQC	C)HerToxys 🖛			
HF_Expts gHMQC	(CC)corr 🔻			
HETCOR				
Spectra V gHETCOR	adamard 🔻			
Spectra GHETCOR Dosy 2D ASAPHMQC	adamard ▼ Dosy 3D ▼ System			
Spectra gHETCOR Dosy 2D ASAPHMQC	adamand ♥ hasy 30 ♥ System			
Spectra V gHETCOR ASAPHMQC Queue Submit Queue	adamard • hasy 30 • System • a × one			
Spectra UHETCOR ASAPHMQC Submit Queue Submit Queue	adamard •			
Specra V Dosy 2D V ASAPHMQC Submit Queue Sample mpleInfo	adamard • have 20 • System) a × one			
Spectra V gHETCOR ASAPHMQC Queue Submit Queue Sample mpletnfo	adamad • hss.20 • System a × one Stat	Acquire Process		
Spectra V gHETCOR Dosy 2D V ASAPHINQC Queue Submit Queue Sample mpleInfo	adamad • han: 3) • System one Stat	Acquire Process CleaSampletine		
Specra V gHETCOR ASAPHMQC Queue Submit Queue Di Sample mpldinfo	adamad • Dase30 • yytem ore Star	Acquire Process CeaCampterre	Sample Preparation AutoBure	
Specira 2 aHETCOR Dasy 2D 4 Submit Queue D Sample mplifinfo	alamara ♥ hacv.30 ♥ yytem 0 × x one Sta	Acquire Process Cleadompleins lard Operator : vmr/1 Logout Sample Name Use My Prifs	Sample Preparation Equilibrium for 10.0. pec	
Secora V Dosy 2D V Submit Queue D Sample mpdrifo	alamar	Acquire Process CleaSamplaire lard Operator : vnmr1 Logout Sample Name Job My Profs Semptime (20119902.01)	Sample Preparation Equilibrium for <u>Do.0</u> sec MonSemptie at <u>Do.0</u> (Current Status Update	
Sectra V Desy 2D V ASAPHNQC Submit Queue D Sample	alamara •	Acquire Process CleaSamplerits Ard Operator : vmm1 Logeut Sampto Hane Sampto I Collision Collisi	Sample Preparation Bouldbard for [D-0.] see Man Sample and Sample Control of	
Specia V gHETCOR Desy 2D V ASAPHAQC Submit Queue D Sample repaires	Strem	Acquire Process CleaSamplains lard Operator : vmm1 Logeout Sample Hues (2019902.03) Notebook Concertation P Eader p	Sampla Preparation Emiliarua for 10.0 mc Mun Sonne Status Update More SP Available Location: 1 Inter Status Location: 1 Inter Status Jonathan	
Specia V altETCOM Dasy D ASPHHQC Submit Queue D Sample mpidrfo	o x system) o x o x o x star st	Acquire Process CleaSamplerits lard Sample Name (20130902).01 Hotelook Page A Concentration P Bador Page Solvert 0M50 CC13 D20 Other A	Sample Preparation Realitizator [Dio.] acc AutoBure ans.2010817.01 Current Status Update Augusto Location: Augusto Location: AutoBure United accession acces	
Specira V	Alamad Stan	Acquire Process CteaSempleirfs arc Operator : vmm1 Loggott Semple Inee Conservation Part Contentiation P Eader p Solvent DMSO CDCI3 D20 Other Comment	Sample Preparation Equilatera for [0-0] ecc RonService Annoted Some United Annoted Some United Annoted Some United Annoted Barles Location I. Loc. Provide Uds For M Barles Location I. Loc. Provide Uds For M DeQ 10456 FM Soft Ods For M	
Sectra V	o x yytem) ○ x ord □ · ord □ · ord · vity sample · · NagnQ	Acquire Process Ceasamplerrs Iard Operator : vmm1 Logeut Sampbor. Nesteok: Pape: A Concertation P Eador Selvent DMSO	Sample Preparation Balibicard for [0.0.3] sec Annosamble all Constants Update ther DP Annopeet Hard Queue Lenessage Before SL DP Queues Lenesson Annored V (1.016) and V (1.016) for Her Stable V (1.016) for Before SL DP Queues (1.016) for Her Stable V (1.016) for Her Stab	
Specira V aHETCOR Dosy 2D V ASAPHIQC Ouser Submit Queue Da Sample mplefino Add Next Selection ts: • DagQ New Study Submit to	Linging Linging Linging Linging	Acquire Process CleaSampletts lard Operator : vmrv1 Logeut Sample Name Use My Prds A Intotecker, D Esser A Schwirth DMSO Sador A Schwirth DMSO COCI3 D2O Other	Sampla Preparation Equitorara for [0:0] ec RunService an [0:0] Current Status Update der SVP Janeoton partier United (Content of the Status) partier United (Content of the Status) Partier Status Update Avantete Autoreton I Partier Status Update Partier Status Update	

Figure 5. Adding the ASAPHMQC protocol to the Study Queue. The ASAPHMQC protocol is found in the **J1(CH)corr** drop down menu and will automatically add a PROTON protocol to the queue if no PROTON is already present.



Figure 6. Customizing the ASAPHMQC experiment. The **Defaults** page contains the fields most commonly modified, such as **Scans per t1 Increment**, **t1 Increments**, and **C13 Spectral Width (ppm)**.

Experimental Method

The ASAPHMQC experiment can be run within or outside the Study Queue mode without any modification. To acquire an ASAPHMQC in the Study Queue:

- 1. Select **New Study** at the bottom of the Study Queue.
- 2. In the Experiment Selector, select the Liquids tab, then select J1(CH)corr and choose ASAPHMOC. (Figure 5). A PROTON will automatically add to the Study Queue before the ASAPHMQC. The default settings are 170 ppm (160 to -10) for the indirect dimension spectral width, a relaxation delay of 60 milliseconds, a mixing time of 25 milliseconds, 96 increments, two repetitions per increment, an acquisition time of 64 milliseconds, and a ¹H-¹³C coupling of 146 Hz.
- 3. To customize the ASAPHMQC experiment, right-click on the **ASAPHMQC** protocol in the Study Queue and select **Open Experiment**. The experiment is retrieved, the pulse sequence displayed, and the Defaults panel is displayed (Figure 6). The typical fields modified are **Scans per t1 Increment, t1 Increments**, and **C13 Spectral Width (ppm)**.

- To lengthen the relaxation delay, select the Acquisition page and enter a longer time in the Relaxation delay field (Figure 7). Lone ¹H-¹³C pairs without neighboring protons may need a longer relaxation delay than the default of 60 milliseconds.
- The experiment is now ready to acquire ASAPHMQC data. Use the Submit button in the Study Queue to initiate data collection.



Figure 7. Customizing the ASAPHMQC experiment. The **Acquisition** page contains additional parameters, such as the **Relaxation delay**, and provides information such as resolution in the indirect dimension.

Conclusions

The ASAPHMQC experiment is a superb choice for obtaining quick protonated carbon (or other heteronuclear) information. The efficiency gained from the crosspolarization period and shortening of the relaxation delay, means information can be obtained in the ASAPHMQC experiment eight times faster than in a gHSQCAD experiment. As such, the experiment is an excellent choice complimentary experiment to be collected along with a 1D ¹H spectrum to quickly provide additional structural information with no special setup requirements.

References

1. Kupce, E., and Freeman, R. Fast multidimensional NMR by polarization sharing. *Magn. Reson. Chem.*, **2007**, 45:2-4.

2. Bodenhausen, G., and Ruben, D. J. Natural abundance nitrogen-15 NMR by enhanced heteronuclear spectroscopy. *Chem. Phys. Lett.*, **1980**, 69:185-189.

www.agilent.com/chem/nmr

Information, descriptions, and specifications in this document are subject to change without notice.

© Agilent Technologies, Inc., 2011 Published in the USA, October 27, 2011 5990-9237EN



Agilent Technologies