

# **Avoid Column Mismatch when Analyzing Polymers**

## **Agilent PLgel Columns for Polymer Analysis**

### **Technical Overview**

#### **Introduction**

When low MW or molecular weight individual pore size columns are combined with mixed gel columns, artefacts can occur when polydisperse samples are analyzed. Careful column selection can avoid such artifacts.




**Agilent Technologies**

The Polymethylurethacrylate shown in the figure has a MW of approximately 73,000 (A). When a PLgel 5  $\mu\text{m}$   $10^3$  Å column is used with a PLgel 5  $\mu\text{m}$  MIXED-C column, resolution is enhanced below the exclusion limit of the  $10^3$  Å column (60,000 MW).

This results in a shoulder on the distribution, which may be wrongly interpreted (B). All of the sample components in (C) have eluted above the exclusion limit of the PLgel 5  $\mu\text{m}$  50 Å column (2,000 MW) and no artefacts are observed, but the pore volume, over which the polymer elutes, is significantly reduced.

#### Conditions

Eluent: THF  
Flow Rate: 1.0 mL/min  
Detection: RI



**Figure 1. Polymethylurethacrylate on different column combinations**

[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 23, 2010

SI-01751



**Agilent Technologies**