Generic method recommendations



Note: This user guide is a convenient starting point for any SPE method development. Further optimization may be required to adjust the method to your application needs.

Europe info_agilent@agilent.com

Asia Pacific inquiry lsca@agilent.com

India india isca marketing@agilent.com

© Agilent Technologies, Inc. 2011 Published in the USA, August 4, 2011 Part Number 7704034 Publication Number 04034-0311

Learn more: www.agilent.com/chem/samplepreparation

Buy online: www.agilent.com/chem/store

Find an Agilent office or authorized distributor: www.agilent.com/chem/contactus

> U.S. and Canada 1-800-227-9770, agilent_inquiries@agilent.com

Bond Elut Plexa SPE method guide



Accuracy Begins Here

The Bond Elut Plexa Family is a new generation of polymeric SPE products, designed for simplicity, improved analytical performance and ease-of-use.These advanced SPE sorbents offer excellent flow characteristics due to their monodisperse particle size distribution, affording superior ease-of-use, with minimal clogging of the packed bed.

Optimized surface chemistries and extraction protocols deliver ultra clean extracts with minimized ion suppression.

The Measure of Confidence



Method development and troubleshooting for plasma samples

Bond Elut Plexa PAX

Bond Elut Plexa PAX contains a strong anion exchange functionality. Simple generic methodology and excellent batch to batch reproducibility offer robust anion exchange SPE workflow.

Strong Anion Exchange SPE for Acidic Analytes			
SPE device Pre-conditioning	1. 500 μL MeOH 2. 500 μL H ₂ O		
Sample	100 µL Plasma		
Pretreatment	Dilute 1:3 with: $2\% \text{ NH}_4 \text{OH in H}_2 \text{O}$		
Washes	1. 500 μL H ₂ 0 2. 500 μL MeOH		
Elution	$2x250\mu\text{L}5\%\text{HCO}_2\text{H}$ in MeOH		

Volumes stated for all methods are for a 30 mg, 1 mL SPE format device.

pH adjustment – To improve ion exchange interactions on Plexa PAX, ionize analytes prior to loading. For acidic analytes the pH should be at least 2 pH units above the pK_a .

Bond Elut Plexa

Deley CDE for neur

Bond Elut Plexa is a non-polar divinylbenzene-based neutral polymeric sorbent. This sorbent is the best choice for non-ionic extraction of a wide range of acidic, neutral and basic analytes from different matrices.

analytes				
SPE device Pre-conditioning	1. 500 μL MeOH 2. 500 μL H ₂ O			
Sample	100 µL Plasma			
Pretreatment	Dilute 1:3 with: 2% NH₄OH <i>(neutrals and bases)</i> 1% HCO₂H in H₂O <i>(acids)</i>			
Washes	500 µL 5 % MeOH in H ₂ 0			
Elution	2 x 250 µL MeOH			

pH adjustment – To improve hydrophobic interaction on Plexa, neutralize analytes prior to loading. Basic analytes should be at least 2 pH units above the pK_a . Acidic analytes should be 2 pH units below the pK_a .

Bond Elut Plexa PCX

Bond Elut Plexa PCX is a cation exchanger with mixed mode sorbent characteristics and is therefore suitable for the extraction and clean-up of polar and non-polar bases from biofluids.

Strong Cation Exchange SPE for Basic Analytes				
SPE device Pre-conditioning	1. 500 μL MeOH 2. 500 μL H ₂ O			
Sample	100 µL Plasma			
Pretreatment	Dilute 1:3 with: $2\% H_3PO_4$ in H_2O			
Washes	1. 500 μL 2% HCO ₂ H in H ₂ O 2. 500 μL MeOH:ACN (1:1, v/v)			
Elution	2 x 250 μL 5% NH ₃ (28-30%) in MeOH: ACN (1:1, v/v)			

 $\label{eq:phi} \begin{array}{l} \textbf{pH adjustment} - \text{To improve ion exchange interactions on} \\ \text{Plexa PCX, ionize analytes prior to loading. Basic analytes} \\ \text{should be at least 2 pH units below the pK}_{a}. Acidification is \\ \text{also necessary to disrupt analyte-protein interaction.} \end{array}$

Troubleshooting	Bond Elut Plexa	Bond Elut Plexa PCX	Plexa PAX	
	Reduce volume of washing stepReduce concentration of organics in the wash step			
step(s)	 Rinse with either 2% NH₃ for basic analytes or 1% formic acid for acids to ensure hydrophobic interactions Increase sorbent bed mass 	Increase sorbent bed mass for increased ion exchange capacity		
Inadequate Elution (Eluent does not contain >90% of the	 Decrease flow rate, (1 mL/min is recommended) Check solubility of analyte in the eluent Increase strength of elution solvent Increase the eluent volume or use multiple aliquots of eluent 			
analyte.j	 Add modifier (depending on analyte type) to the elution solvent, thereby promoting ionization 	 Use up to 10% ammonia (28-30%) in solvents such as MeOH and ACN 	 Use up to 10% formic acid in MeOH for anion exchange elution 	