




Agilent Cerity Networked Data System for Pharmaceutical QA/QC



Getting Started Guide



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2003

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G4000-90011

Edition

5/2003

Printed in Germany

Agilent Technologies Deutschland GmbH
Hewlett-Packard-Strasse 8
76337 Waldbronn, Germany

Microsoft® is a U.S. registered trademark of Microsoft Corporation.

Software Revision

This guide is valid for A.02.xx revisions of the Agilent Cerity Networked Data System for Pharmaceutical QA/QC software, where xx refers to minor revisions of the software that do not affect the technical accuracy of this guide.

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as

defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

Before you start 5

Running Routine Samples 9

**Basic Exercise #1
Equilibrate the instrument 11**

**Basic Exercise #2a
Run a single sample to produce an example chromatogram 17**

**Basic Exercise #2b
Run a group of single samples to identify compounds 23**

**Basic Exercise #3a
Run a sequence to quantify compounds with single-level
calibration 29**

**Basic Exercise #3b
Reintegrate and reprocess the results 39**

**Advanced Exercise #4a
Run a sequence to quantify compounds with multi-level calibration 45**

**Advanced Exercise #4b
Change sample variables in the method and reprocess 53**

**Advanced Exercise #5a
Run a sequence to quantify impurities 61**

**Advanced Exercise #5b
Use a different method to reprocess 67**

Setting Up Methods 71

Basic Exercise #1

Set up an equilibration method 73

Basic Exercise #2

Set up a method for single samples to identify compounds 81

Basic Exercise #3


Set up a single-level calibrated method for a sequence 91

Advanced Exercise #4

Set up a multi-level calibrated method for a sequence 105

Advanced Exercise #5

Set up a method for a sequence to quantify impurities 117



Before you start

The Getting Started Exercises provide a quick way to learn the Certy Pharmaceutical QA/QC application. Use the *Certy Concepts Guide* to help you do the tasks in these exercises.

Setting Up Methods

If you develop methods for your laboratory, you should go through these exercises. You can use these methods to run samples and sequences with the Running Routine Samples exercises.

Running Routine Samples

If you run samples but do not develop methods, you can do these exercises with the default methods that come with the Certy Networked Data System, or you can use the methods set up with the Setting Up Methods exercises.

Before you start

Make sure that you or your administrator transfer the default methods and example chromatogram from the Certy CD-ROM to the database. For details to transfer the methods and make them usable for your system, turn to the next page.



Step 1. Restore the default methods

The default methods for the Basic and Advanced exercises are located on Cerity software CD in **\GettingStarted\DefaultMethods**.

- 1 Restore the default methods.

The default methods for the Basic and Advanced exercises are located on Cerity software CD in **\GettingStarted\DefaultMethods**.

- 2 Select **Start > Programs > Agilent Cerity > Administration and Maintenance > Archive and Restore**.
- 3 Enter logon information and click **OK**.
- 4 Select **Restore**, and click **Next**.
- 5 Click the ... button.
- 6 Select **\GettingStarted\DefaultMethods\Basic** (or **\Advanced**) on the CD-Drive.
- 7 Click **OK**, click **Next**, and click **Yes** to the messages.
- 8 Click the >> button to move the default methods to the **Restore Objects** list.
- 9 Click **Next**, click **Start**, and click **OK** for each message that appears.

The following message appears: "These tables contain duplicates".

Step 2. Resolve database duplicates

- 1 Click **Next**.
- 2 Make sure that the **Select instruments to enable** check box is clear.
- 3 Click **Next** and select the second Administrator role.
- 4 Click **Rename**, enter the new role name Admin and click **OK**.
- 5 Click **Next**, click **Start**, and click **OK**.
- 6 Click **OK** and any **Close** buttons.

Step 3. Restore the example chromatogram

The example chromatogram is on Cerity-CD-1 in **\GettingStarted\DefaultResults**. Make sure that the default example chromatogram has been restored.

- 1 Repeat step 1 through step 4 in “[Step 1. Restore the default methods](#)” on page 6.
- 2 Select **\GettingStarted\DefaultResults** on the CD-ROM drive, click **OK**, and click **Next**.
- 3 Select **defexchrom2a**, click **>**, and click **Next**.
- 4 Click **Start**, click **OK** to the messages that appear, and click **Close**.
- 5 Select **Start > Programs > Agilent Cerity > Cerity Pharmaceutical QA/QC**.
- 6 Enter logon information and click **OK**.
- 7 Select **Result** from the Current View list.
- 8 Select **AllResultsRestored** from the Query list.

Step 4. Copy the default method to use with your instrument

Refer to “Basic Exercise #2 Set up a method for single samples to identify compounds” on page 81 if you need to.

- 1 Select **Method** from the **Current View** list.
- 2 Select **AllMethodsRestored** from the **Query** list.
- 3 For each default method:
 - a Select **File > New > Method**.
 - b Click **Browse**, select **defaultmethodN** for Basic exercises, or **AdvdefaultmethodN** for Advanced exercises, and click **OK**.
 - c Name the new method **defexerN**, and click **Next**.
 - d Select the instrument where the method will be used, and click **Next**.
 - e Click **Next** until you reach the New Method Review panel.
 - f Click **Finish**, and click **Save** when the Save to the database message appears.
- 4 Select **AllMasterMethods** from the **Query** list.
- 5 Expand **defexerN**.
- 6 Expand **Instrument Setup**, and adapt the settings.
- 7 Adapt the instrument settings for the non-matching LC modules.

NOTE

The first time that you copy and rename **Advdefaultmethod4**, name it **defexer4a**. The first user will alter this method in Exercise 4b. You must then copy **Advdefaultmethod4** and rename it **defexer4b** for the second user to use the method.

You can use the default methods ONLY on instruments with an Agilent VWD detector. Your other LC modules do NOT have to match the modules on which the default methods were set up (autosampler, quaternary pump, thermostatted column compartment).

If you have no instrument available with a VWD detector to use with these exercises, then the administrator or advanced user should set up the methods using the Setting Up Methods sections in this guide.



Running Routine Samples

These exercises help you learn how to run routine samples. You can use the default methods for the “a” exercises or set up methods in the Setting Up Methods exercises. You must have results from the “a” exercises to do the “b” exercises. The set of basic and advanced exercises includes the topics below:

Basic **Exercise 1 – Equilibrate the instrument** Learn how to equilibrate the instrument with the instrument panel or with a method.

Exercise 2a – Run a single sample to produce an example chromatogram Learn how to produce an example chromatogram that you can use to set up integration and identification in a method.

Exercise 2b – Run a group of single samples to identify compounds Learn how to enter and run a group of single samples with a method to identify the compounds in the sample.

Exercise 3a – Run a sequence to quantify compounds with single-level calibration Learn how to run a sequence with single-level, single-update calibration, ESTD quantitation, and fixed amounts.

Exercise 3b – Reintegrate and reprocess the results Learn how to manually reintegrate the sequence results, and reprocess the results with the original method revision. For more information on running routine samples, see the *Concepts Guide*, “Sample Analysis”.



Advanced **Exercise 4a – Run a sequence to quantify compounds with multi-level calibration** Learn how to run a sequence set up for multi-level, overall calibration, variable compound amounts and sample variables.

Exercise 4b – Change sample variables in the method and reprocess Learn how to reprocess the results with the most current version of the method and a version with new sample variables.

Exercise 5a – Run a sequence to quantify impurities Learn how to create and run a sequence set up for ISTD quantitation, custom calculations, limits, bracketed calibration and system suitability.

Exercise 5b – Use a different method to reprocess Learn how to reprocess with a new method.

Before you start Read “**Before you start**” on page 5.

If you plan to use default methods in these exercises, make sure that these methods are in your database. From the Query list, select AllMethodsRestored to view defexer1-5 or AllResultsRestored to view defexchrom2a.


Your system administrator must have configured an Agilent 1100 LC series liquid chromatograph for your system.

If you choose to do the Running Routine Samples exercises with the default methods, you must use an instrument with a VWD detector. If you use the methods created in the Setting Up Methods exercises, you need only an autosampler, pump (quaternary or binary) and UV-Vis detector (VWD, MWD, DAD).

Solvent A is water. Solvent B is methanol or acetonitrile.

Use Agilent Technologies column Eclipse XDB-C8 (or C-18), 4.6MM X 15 CM (5µM).

Prepare the following three vials of the isocratic standard, Agilent Part # 01080-68704: undiluted, diluted by factor 2, and diluted by factor 4.



Basic Exercise #1

Equilibrate the instrument

This exercise contains a series of tasks to help you learn how to:

- Equilibrate the instrument with the instrument panel in the Cerity Pharmaceutical QA/QC application
- Enter and run an equilibration sample (blank run) with a method created to equilibrate the instrument

You can use a copy of the default method that comes with the system to equilibrate the instrument, or you can use the method created in the [“Basic Exercise #1 Set up an equilibration method”](#) on page 73.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Make sure that the pump is on standby and the VWD lamp is off.


Make sure that the methods for this exercise have been set up or restored.



Basic Exercise #1 Equilibrate the instrument

Task 1. Purge the pump from the Instrument Panel

Task 1. Purge the pump from the Instrument Panel

Steps	Detailed Instructions		
<div>1 Disengage pump and purge line B.</div> <div>Flow rate: 5ml/min</div> <div>%B = 100%</div>	<div><div>a Turn the black valve on the pump counterclockwise two full turns.</div><div>b Select Instrument from the Current View list.</div><div>c Select the instrument that you intend to equilibrate.</div></div> <div>The Instrument Panel appears, along with the Online Plot.</div> <div></div> <div><div>d Click the pump module on the Instrument Panel.</div><div>A menu appears.</div></div> <div><div>e Select Set Pump.</div><div>f Enter a Flow of 5ml/min and %B=100, and click OK.</div></div> <tr><td><div>2 Purge line A and engage pump.</div><div>%A = 100</div></td><td><div><div>a When there are no more bubbles in the line, repeat steps d and e from step 1.</div><div>b Set %B = 0, and click OK.</div><div>c When there are no more bubbles in the line, click the pump module, and select Standby.</div><div>d Tighten the black valve.</div></div></td></tr>	<div>2 Purge line A and engage pump.</div> <div>%A = 100</div>	<div><div>a When there are no more bubbles in the line, repeat steps d and e from step 1.</div><div>b Set %B = 0, and click OK.</div><div>c When there are no more bubbles in the line, click the pump module, and select Standby.</div><div>d Tighten the black valve.</div></div>
<div>2 Purge line A and engage pump.</div> <div>%A = 100</div>	<div><div>a When there are no more bubbles in the line, repeat steps d and e from step 1.</div><div>b Set %B = 0, and click OK.</div><div>c When there are no more bubbles in the line, click the pump module, and select Standby.</div><div>d Tighten the black valve.</div></div>		

Task 2. Equilibrate the instrument from the Instrument Panel

Steps

Detailed Instructions

1 Enter the pump parameters

Methanol as Solvent B:

- Flow rate: 2ml/min.
- Solvent composition: 80%MeOH/20%H₂O

Acetonitrile as Solvent B:

- Flow rate: 1.5ml/min
- Solvent composition: 65%ACN/35%H₂O

a Click the pump module on the Instrument Panel.

b Select **Set Pump**.

The Set Pump dialog box appears.

c Enter the pump parameters as shown in the left column, and click **OK**.

The screenshot shows the 'Set pump LC_18' dialog box. The 'Flow' field is set to 2 ml/min. The 'Solvents' section shows Solvent A at 20% and Solvent B at 80% (checked). The 'Act. Fill (liters)' and 'Max. Fill (liters)' fields are also visible.

Solvents	Act. Fill (liters)	Max. Fill (liters)
A: 20 %	0.097	3.5
B: <input checked="" type="checkbox"/> 80 %	0.597	3.3
C: <input type="checkbox"/> Off	0	5
D: <input type="checkbox"/> Off	0	5

d Click the pump module, and select **On**.

2 Turn the detector lamp on


a Click the detector module on the Instrument Panel.

b Select **Lamp On**.

Wait until baseline has stabilized.



Basic Exercise #1 Equilibrate the instrument
Task 2. Equilibrate the instrument from the Instrument Panel

Steps	Detailed Instructions
3 Monitor the baseline until it appears stable. After this step, you are ready to do the remaining exercises, or you can move on to the next task to learn to equilibrate the instrument with a method.	<ul style="list-style-type: none">a Click Change at the bottom of the Online Plot. The Edit Signal Plot dialog box appears.b Select the detector signal you need from the Available Signals list, and click the Add button to put the signal in the Selected Signals list. (You can also select the pump pressure).c Set the Predictable Range (Y-axis) as -10 to +10.d Set the X-Axis range as 10 min.e Click OK.



- f** Click the detector module after the lamp has been on for a few minutes.
- g** Select **Balance**.
When the baseline stays at zero for a few minutes after the balance, the baseline is considered stable.



Task 3. Equilibrate the instrument with a method—Enter an equilibration sample

Steps	Detailed Instructions
1 Enter the sample information Sample Name: equilsamp <i>iii</i> , where <i>iii</i> are your initials Method: defexer1 or equilmeth <i>iii</i> See “Before you start” on page 5 for instructions on how to restore and copy the default methods.	<ol style="list-style-type: none"> Select Instrument from the Current View list. Expand the Sample Entry folder for the instrument that you need to equilibrate. Select Single Samples. Enter the Sample Name as equilsamp<i>iii</i>. Select the Method as equilmeth<i>iii</i> or defexer1. Select the Sample Type as Blank Run. Click Apply. <p>You can also enter the sample in the Sample View when you need to enter samples and sequences during a run.</p>
2 Enter the tasks that the system will do during the analysis.	<ol style="list-style-type: none"> Clear the Quantify and Report check boxes. Click Apply.
	
3 Save the sample to the database	<ol style="list-style-type: none"> On the Standard toolbar, click . Review the list of changes Under Reason for changes, enter a reason or select a reason from the list. Enter your electronic signature if required. Click the Save button.

Basic Exercise #1 Equilibrate the instrument

Task 4. Equilibrate the instrument with a method—Run the equilibration sample

Task 4. Equilibrate the instrument with a method—Run the equilibration sample

Steps	Detailed Instructions
1 Run <i>equilsamp111</i>	<p>a Select the sample, <i>equilsamp111</i>, in the Sample Table. The Run button is now active.</p> <p>b Click the Run button  on the Actions toolbar.</p>
2 Monitor the baseline until it is stable.	<p>a Select the instrument that you want to equilibrate. The Instrument Panel appears, along with the Online Plot.</p> <p>b Click Change at the bottom of the Online Plot. The Edit Signal Plot dialog box appears. (See the figure on page 14.)</p> <p>c Select the detector signal you need from the Available Signals list, and click the Add button to put the signal in the Selected Signals list.</p> <p>d Set the Predictable Range as -10 to +10.</p> <p>e Set the X-Axis range as 10 min.</p> <p>f Click OK.</p> 



Basic Exercise #2a

Run a single sample to produce an example chromatogram

This exercise contains a series of tasks to learn how to:

- Enter a sample to produce an example chromatogram
- Run the sample
- Review the results

An example chromatogram can be any chromatogram that you produce. Use the example chromatogram to test new integration parameters and identify peaks as compounds.

You can use one of the following methods for this exercise:

- A copy of the default method provided with the Cerity Networked Data System.
- The method saved in [“Task 3. Save and audit method changes”](#) on page 86 in the Setting Up Methods section.
- An equilibration method that you created in [“Basic Exercise #1 Set up an equilibration method”](#) on page 73.

For the tasks on the next pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Running Routine Samples”](#) on page 9 for running routing samples.



Equilibrate the instrument. See [“Basic Exercise #1 Equilibrate the instrument”](#) on page 11. Make sure that the methods for this exercise have been set up or restored.



Basic Exercise #2a Run a single sample to produce an example chromatogram


Task 1. Enter a single sample

Task 1. Enter a single sample

Steps	Detailed Instructions
1 Start the Instrument View to find the sample table for single samples.	<p>a Select Instrument from the Current View list.</p> <p>b Expand the folder for the instrument that will produce the example chromatogram.</p> <p>c Select Single Samples.</p> <p>The sample table and sample entry panel appear in the workspace.</p>
2 Enter a sample with the following information:	<p>a Enter <i>exchromiii</i> in the Sample Name box.</p> <p>b Select a method from the Method list.</p> <p>The instrument associated with the method appears in the Instrument box.</p> <p>c Select Sample from the Sample Type list.</p> <p>d Enter the vial number for the sample in the Vial Number box.</p> <p>e Click Apply to put the sample information in the sample table.</p> <p>Use the default values for all other parameters</p>
3 Enter the tasks to perform during the run.	<p>a Clear the Quantify and Report check boxes.</p>
	
4 Save the sample.	<p>a On the Standard toolbar, click .</p> <p>The Save Changes To The Database dialog box appears.</p> <p>b Review the List of changes.</p> <p>c Under Reason for changes, enter a reason or select a reason from the list.</p> <p>d Enter your electronic signature if required.</p> <p>e Click the Save button.</p>

Task 2. Run the sample

Steps	Detailed Instructions
1 Check that the instrument is ready for use.	<p>a On the selection tree, select your instrument.</p> <p>b Click the Online Plot tab.</p> <p>c Click the Change button.</p> <p>The Edit Signal Plot dialog box appears.</p> <p>d Select the detector signal you need from the Available Signals list.</p> <p>e Click the Add button to put the signal in the Selected Signals list.</p> <p>f Select the Predictable Range option and set the predictable range from -20mAU to 300mAU.</p> <p>g Under Window Properties, enter 5 min in the X-Axis range box.</p> <p>h Click the OK button.</p>



Basic Exercise #2a Run a single sample to produce an example chromatogram

Task 2. Run the sample

Steps

2 Run the sample.

Detailed Instructions

a

On the selection tree, expand your instrument folder.


b

Select **Single Samples**.


c

Select the sample, *exchromiii*.

The Run button



becomes available on the Tools toolbar.




	NAME	METHOD	INSTRUMENT
1	equilsampeme	equilmetheme [Older version]	ee-lc3
2	ssexchromeme	ssexchromeme	ee-lc3
3	ssexchromdec	ssexchromeme	ee-lc3


d

Click the **Run** button.

You can also run the sample from the Sample View.

3 Monitor the signal, and track the status of the sample.	<p>a On the selection tree, select your instrument.</p> <p>b Click the Online Plot tab to view the signal.</p> <p>Change the axes if necessary.</p> 
---	--

c Click the **Worklist** tab to track the status of the sample.



The screenshot shows the Agilent Cerity NDS software interface. The title bar reads "Agilent Cerity NDS for Pharmaceutical QA/QC - AGILENT\MuskService - Administrator - Cerity for Pharma QA/QC". The menu bar includes "File", "Edit", "View", "Go", "Tools", "Actions", and "Help". Below the menu bar is a toolbar with various icons. To the right of the toolbar is a dropdown menu labeled "Instrument" with a downward arrow. Below the toolbar is a row of icons for file operations, including "Fill Down". On the left side, there is a tree view under "AllInstruments" showing a folder structure with "e-vdt" and "e-lc3". On the right side, the "Instrument Panel" is visible, with the "Worklist" tab selected. The Worklist table has the following data:

	Sample Name	Status	Sample Type	Method	Pri
1	ssexchromeme2	Completed	Sample	metsscpcdsme	500

After you click the **Worklist** tab, the **Abort**, **Pause** and **Resume** buttons become available.


Task 3. Review the chromatogram

Steps	Detailed Instructions
1 Review the sample result and make sure all four peaks are integrated.	<div>a Select Result from the Current View list.</div> <div>b Select MySamplesRunLast24h from the Query list.</div> <div>c Expand the Samples folder.</div> <div>d Expand the exchromiii folder.</div> <div>e Select the exchromiii #1 injection.</div> <div>f View the chromatogram and results.</div>

RT	Compound Name	Amount	RF (Resp/Amt)	Peak Area	Peak Height	P
0.56		N/A	N/A	0.5678	0.1215	
0.76		N/A	N/A	0.7701	0.3293	
0.94		N/A	N/A	419.6385	153.4289	
1.11		N/A	N/A	374.5102	126.7572	
1.44		N/A	N/A	2.6036	0.7431	
1.75		N/A	N/A	0.2067	0.0663	
1.89		N/A	N/A	357.0249	98.3153	
3.09		N/A	N/A	623.8801	90.8962	

Basic Exercise #2a Run a single sample to produce an example chromatogram

Task 3. Review the chromatogram



Basic Exercise #2b

Run a group of single samples to identify compounds

This exercise contains a series of tasks to learn how to:

- Enter a sample
- Run and track groups of single samples
- Review the results to check compound identification

You can use one of the following methods for this exercise:

- A copy of the default method provided with the Cerity Networked Data System (NDS).
- The method completed in [“Basic Exercise #2 Set up a method for single samples to identify compounds”](#) on page 81.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Running Routine Samples”](#) on page 9.

Equilibrate the instrument. See [“Basic Exercise #1 Equilibrate the instrument”](#) on page 11.

Make sure that the methods for this exercise have been set up or restored.



Basic Exercise #2b Run a group of single samples to identify compounds

Task 1. Enter three single samples

Task 1. Enter three single samples

Steps

Detailed Instructions

- 1 Start the Instrument View and find the sample table for single samples.**

- Select **Instrument** from the **Current View** list.
- Expand your instrument folder
- Select **Single Samples**.

The sample table and **Sample Entry** tab sheet appear in the workspace.

- 2 Enter a sample with the following information:**


- Name the sample *exer2biii1*, where *iii* are your initials.
- Select the method for the sample: *defexer2* or *exer2iii*
- Select the Vial # that contains the the full-strength isocratic standard.

- Enter *exer2biii1* in the **Sample Name** box.
- Select the *exer2* method from the **Method** list (or copy of *defexer2b*).
The instrument associated with the method appears in **Instrument** box.
- Select **Sample** from the **Sample Type** list.
- Enter the **Vial Number** that contains the standard.
- Click **Apply** to put the sample information into the sample table.

- 3 Enter the tasks that you want the system to do during the run**

- Mark the **Quantify** check box, and clear the **Report** check box.
You must mark the **Quantify** check box to identify the compounds, even though Calibration and Quantitation are not set up in the method.
- Click **Apply**.

- 4 Save the sample**

- On the Standard toolbar, click .
The **Save Changes To The Database** dialog box appears.
- Review the **List of changes**.
- Under **Reason for changes**, enter a reason or select a reason from the list.
- Click the **Save** button.

Basic Exercise #2b Run a group of single samples to identify compounds

Task 1. Enter three single samples

Steps

Detailed Instructions

5 Repeat Steps 2 through 4 for the next two samples.

Name these samples, exer2biii2 and exer2biii3.

- Select the empty row.
- Start with Step 2a and finish with Step 4d for exer2biii2.
- Repeat steps a and b for exer2biii3.

1

	INSTRUMENT NAME	METHOD NAME	SAMPLE NAME	NUM OF INJECTIONS
1	EMELC3	exer2dec	exer2bdec3	1
2	EMELC3	exer2dec	exer2bdec2	1
3	EMELC3	exer2dec	exer2bdec1	1
4				

Sample Entry | Sample Logbook

Sample Name:
exer2bdec3

Method:
exer2dec

Sample Type:
Sample

Instrument:
EMELC3

Vial Number Injections Volume [µl]
1 1 as method

Apply

Run | Amounts | Identification | Description | Report Destination

Run with:
Priority: Medium Schedule: Ready for Analysis

Task(s) to perform:
☒ Acquire ☒ Quantity
☒ Integrate ☐ Report


Analyst
SCHEIDERER,ROBIN

Basic Exercise #2b Run a group of single samples to identify compounds

Task 2. Run the samples


Task 2. Run the samples

Steps	Detailed Instructions
1 Check that the instrument is ready.	<p>a Select Instrument from the Current View list.</p> <p>b Click the Online Plot tab.</p> <p>c Click the Change button.</p> <p>The Edit Signal Plot dialog box appears.</p> <p>d Select the detector signal you need from the Available Signals list.</p> <p>e Click the Add button to put the signal in the Selected Signals list.</p> <p>f Select the Predictable Range option and set the range from -20mAU to 300mAU.</p> <p>g Under Window Properties, enter 15 min in the X-Axis range box.</p> <p>h Click the OK button.</p>



Basic Exercise #2b Run a group of single samples to identify compounds

Task 2. Run the samples

Steps	Detailed Instructions
2 Run the samples.	<div><div><div>a Expand your instrument folder.</div><div>b Select Single Samples.</div><div>c Select the sample, exer2biii1.</div><div>d Click the Run button .</div><div>e Select the sample, exer2biii2.</div><div>f Click the Run button.</div><div>g Select the sample, exer2biii3.</div><div>h Click the Run button.</div></div><p>The samples run in the order started, unless exer2biii3 is of a higher priority than exer2biii2. Then, exer2biii3 runs before exer2biii2. The first sample started will always run first even if the sample is a lower priority than the other samples.</p></div>

3 Monitor the signal, and track the status of the samples.	<div><div>a Click the Online Plot tab to view the signal.</div><div>Change the axes if necessary.</div><div>b Click the Worklist tab, and track the status of the three samples</div></div>
--	---

1

Instrument Panel		Worklist						
	Name	Status	Type	Method	Priority #	Vial #	Injections #	Description
1	exer2bdec1	Running(1)	Sample	exer2dec	500	1	1	
2	exer2bdec2	Queued	Sample	exer2dec	500	1	1	
3	exer2bdec3	Queued	Sample	exer2dec	500	1	1	

Basic Exercise #2b Run a group of single samples to identify compounds

Task 3. Review the chromatogram

Task 3. Review the chromatogram

Steps	Detailed Instructions
1 Review the sample results and make sure all the compounds are identified in each sample.	<p>a Select Result from the Current View list.</p> <p>b Expand the Calibration - exer2biii folder or defexer2 folder. Even though calibration was not set up in the method, the result appears in a Calibration folder.</p> <p>c Expand the Samples folder.</p> <p>d Expand the exer2biii1 folder.</p> <p>e Select the exer2biii1 #1 injection.</p> <p>f View the result.</p> <p>g Repeat steps d through f for the following samples:</p> <ul style="list-style-type: none">• exer2biii2• exer2biii3.

The screenshot displays the Agilent Certify NDS software interface. On the left, a tree view shows the project structure: 'All Samples Not Approved Run Last 7 Days' > 'Samples' > 'Calibration - exer2dec Calb Rev 2' > 'Samples' > 'exer2bdec3 [Rev 2]' > 'exer2bdec2 [Rev 2]' > 'exer2bdec1 [Rev 2]' > 'exer2bdec3 #1 [Rev 2]'. The main window shows two chromatograms. The top one is titled 'Example Chromatogram' and the bottom one is 'Signal'. Both show peaks at retention times 0.93, 1.10, 1.88, and 3.08 minutes. Below the chromatograms is a 'Results' table.

RT	Compound Name	Amount	RF (Rsp/Amt)	Peak Area	Peak Height
0.93	dimethylphthalate	N/A	N/A	402.2817	153.2890
1.10	diethylphthalate	N/A	N/A	377.3364	126.5203
1.88	biphenyl	N/A	N/A	358.5160	98.0923
3.08	o-tolphenyl	N/A	N/A	506.3579	97.2367



Basic Exercise #3a

Run a sequence to quantify compounds with single-level calibration

This exercise contains a series of tasks to help you learn to:

- Create a sequence with a method set up for single-level, single-update calibration, ESTD quantification and fixed compound amounts
- Select report types and set up a directory for reports
- Run and track the sequence
- Review the results to make sure the compounds have been identified and quantified correctly
- Review the reports

You can choose between two methods to use with this exercise:

- a copy of the default method provided with the system.
- method that you created in [“Basic Exercise #3 Set up a single-level calibrated method for a sequence”](#) on page 91.

For the Basic exercises, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before You Start

Read [“Running Routine Samples”](#) on page 9.

Equilibrate the instrument. See [“Basic Exercise #1 Equilibrate the instrument”](#) on page 11.



Place all the vials of prepared samples into the ALS tray. Make sure that the methods for the exercise have been set up or restored.



Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration

Task 1. Create a new sequence

Task 1. Create a new sequence

Steps	Detailed Instructions
<p>Create a new sequence.</p> <p>Name the sequence <code>exer3seqiii</code>, where <i>iii</i> are your initials.</p> <p>Use one of the two methods:</p> <ul style="list-style-type: none">• <code>defexer3</code>• <code>exer3iii</code> (created with Exercise 3 of Setting Up Methods)	<p>a Click the New button,  , in the Standard toolbar, and select Sequence. The Create New Sequence dialog box appears.</p> <p>b Enter the Sequence Name as <code>exer3seqiii</code>.</p> <p>c Select the Instrument that will run the sequence.</p> <p>d Select the Method for the sequence.</p> <p>e Click OK.</p> <div><p>The image shows the 'Create New Sequence' dialog box. It has a title bar with a question mark and a close button. Inside, there are three text input fields: 'Sequence Name' with the value 'exer3seqeme', 'Instrument' with the value 'GetStartLC', and 'Method' with the value 'exer3singlevel'. To the right of the 'Instrument' and 'Method' fields are 'Browse...' buttons. At the bottom are 'OK' and 'Cancel' buttons.</p></div> <p>f If the Save Changes to the Database dialog box appears, select the Reason for changes, if present, and click Save.</p>

Task 2. Enter sample and sequence information

Steps

Detailed Instructions

1 Review the Sequence Table

Note how the sequence table matches the sequence template setup in the method.

- Select Instrument from the Current View list.
- Expand the instrument you are using, and select the sequence you just created.
- Review the table.

Sequence Table		Sequence Options						
	Sample Name	Sample Type	Cal. Level	Summary Group	Vial #	Injections #	Injection Volume [μl]	Sample Amount [mg/ml]
1	Cal1	Calibration	1		2	1	as method	0
2	sample 1_2	Sample			5	1	as method	0
3	sample 1_4	Sample			9	1	as method	0
4	Cal1	Calibration	1		2	1	as method	0
5	sample 1_2	Sample			5	1	as method	0
6	sample 1_4	Sample			9	1	as method	0
7	Cal1	Calibration	1		2	1	as method	0
8	sample 1_2	Sample			5	1	as method	0
9	sample 1_4	Sample			9	1	as method	0
10								

2 Enter the tasks to be performed during the run:

Quantify, Report.




Acquire and Integrate are always marked.

- Click the **Sequence Options** tab.
- Make sure that the **Quantify** and **Report** check boxes are marked for the Task(s) to perform.

Sequence	Identification	Description	Report Destination
Run with			
Priority:	Schedule:		
Medium	Ready for Analysis		
Calibration Mode:			
Single Update Calibration			
Sequence Created by			
Task(s) to perform			
<input checked="" type="checkbox"/> Acquire		<input checked="" type="checkbox"/> Quantify	
<input checked="" type="checkbox"/> Integrate		<input checked="" type="checkbox"/> Report	
<input type="checkbox"/> Allow Online Editing			

Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration

Task 2. Enter sample and sequence information

Steps	Detailed Instructions
<p>3 Enter the destination path for, but do not print, the reports:</p> <p>Enter Exercise3iii, where "iii" are your initials.</p>	<p>a Click the Report Destination tab.</p> <p>b Clear the Printer check box, if necessary.</p> <p>c Mark the Path check box, and enter the directory, Exercise3iii.</p> <p>The system automatically creates this directory if it does not exist and places the generated reports into the directory Agilent\Cerity\Reports\Pharmaqc\Reports</p>
	
<p>4 Select the following reports to be generated:</p> <p>Single Injection</p> <p>Standard Injection</p> <p>Sequence</p>	<p>a Mark the Print check box to the left of the Report Types noted on the left margin.</p> <p>b Clear all the Print check boxes that are not those noted on the left margin.</p>
	
<p>5 Save the sequence</p>	<p>a Click , and enter reasons for changes and your password, if necessary.</p>

Task 3. Run and track the sequence


Steps

1 Make sure that the instrument is ready.


- Use the same conditions as set in the method.
- Online Plot settings:
Y-Axis range: -20 to 300
X-Axis range: 15 minutes

Detailed Instructions

- Select the instrument for the sequence from the selection tree.
- Make sure the instrument and column are equilibrated, and the conditions are the same as those set in the method for the sequence.





- Click **Change** at the bottom of the Online Plot.
The Edit Signal Plot dialog box appears.
- Select the detector signal you need from the Available Signals list, and click **Add** to place this signal on the right.
- Set the **Predictable Range** as -20 to 300.
- Set the **X-Axis range** as 15 min.
- Click **OK**.




Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration

Task 3. Run and track the sequence

Steps	Detailed Instructions
2 Run the sequence.	<p>a Expand the instrument folder.</p> <p>b Select the sequence that you just set up.</p> <p>The Run button, , appears.</p> <p>c Click the Run button.</p>



3 Monitor the signal, and track the status of the sequence.	<p>a Select the instrument.</p> <p>b Observe the signal in the Online Plot tab, and change the axes if you need to.</p>
---	--



c Click the Worklist tab, and observe the status of the sequence.
--

Instrument Panel Worklist								
	Name	Status	Type	Method	Priority #	Vial #	Injections #	Descripti
1	defseq3	Running(3.1)	Sequence	defaultmethod3	500	N/A	N/A	

Note that the Abort, Pause and Resume buttons appear when you enter the Worklist.


Task 4. Review the results and reports

Steps

Detailed Instructions

- 1 Review the calibration table and curve for each revision of the calibration.

- a Select **Result** from the Current View.
 - b Select **AllSeqNotApprovedRunLast7Days** from the Query list.
 - c Expand the **exer3seqiii** folder.
 - d Select the **Calibration - exer3seqiii Calib Rev 2** folder.
- The calibration table and curve appear in the workspace.



- e Select the **Calibration - exer3seqiii Calib Rev 3** folder.
- f Select the **Calibration - exer3seqiii Calib Rev 4** folder.

Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration

Task 4. Review the results and reports

Steps

Detailed Instructions

2 Review the results for each calibration standard in each revision.

Note the different response factors used to quantify the samples.

- Expand the **Calibration - exer3seqiii Calib Rev 2** folder.
- Expand the **Calibrations** folder.
- Expand the Cal1 folder.
- Select Cal1 #1.
- Observe the response factor in the workspace.

The screenshot displays the Agilent NDS software interface for Pharmaceutical QA/QC. The left pane shows a tree view with the following structure:

- AllSeqNotApprovedRunLast7Days
 - defseq3 [Rev 3]
 - Calibrations
 - Calibration - defaultmethod3 Calib Re
 - cal1 [Rev 1]
 - cal1 #1 [Rev 1] (Selected)
 - Samples
 - Calibration - defaultmethod3 Calib Re
 - Calibration - defaultmethod3 Calib Re
 - Old Revisions

The main window shows a chromatogram titled 'Signal' with 'VWD: Absorbance' on the y-axis (0 to 120 mAU) and time on the x-axis (0 to 4 minutes). Three peaks are labeled: '0.94 - dimethylphthalate', '1.11 - diethylphthalate', and '3.05 - biphenyl'.

Below the chromatogram is a 'Results' table:

RT	Compound Name	Amount	RF (Rsp/Amt)	Peak Ar
0.94	dimethylphthalate	10.0000	43.3062	433.062
1.11	diethylphthalate	7.1673	43.3062	387.986
1.87	biphenyl	15.0000	24.4655	366.982
3.05		19.8666	27.1839	540.052

- Expand the **Calibration - exer3seqiii Calib Rev 3** folder.
- Repeat steps b-c.
- Select the second Cal1 standard.
- Observe the response factor.
- Expand the **Calibration - exer3seqiii Calib Rev 4** folder.
- Repeat steps b-c.
- Select the third Cal1 standard.
- Observe the response factor.

Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration

Task 4. Review the results and reports

Steps

Detailed Instructions

3 Review the sample results for each revision.

Note the response factor used for the quantitation.

- a Expand the **Calibration - exer3seqiii Calib Rev 2** folder.
- b Expand the **Samples** folder.
- c Expand the **Sample1_2** folder.
- d Select **Sample1_2 #1**.
- e Observe the response factor in the workspace.
- f Repeat steps c-e for **Sample1_4**.

The chromatogram displays four distinct peaks. The first peak at 0.94 minutes is the tallest, followed by peaks at 1.11, 1.87, and 3.05 minutes. The baseline is stable throughout the run.

RT	Compound Name	Amount	RF (Rsp/Amt)	Peak Area	Peak Height
0.94	dimethylphthalate	2.4563	43.5878	107.0642	38.4523
1.11	diethylphthalate	1.7347	43.5878	94.5155	31.6319
1.87	biphenyl	3.4917	24.6553	86.0869	23.8970
3.05	biphenyl	4.6752	27.3947	128.0765	22.6674

- g Expand the **Calibration - exer3seqiii Calib Rev 3** folder.
- h Repeat steps b-f.
- i Expand the **Calibration - exer3seqiii Calib Rev 4** folder.
- j Repeat steps b-f.


4 Review the reports.

Hint: Use the Report Viewer to open the reports.

- a Select **Start > Programs > Agilent Cerity > Report Viewer**.
- b Select **File > Open**.
- c Open **Cerity > Agilent > Reports > PharmaQC > Reports > Exercise3iii**.
- d Open and view each report.

Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration

Task 4. Review the results and reports



Basic Exercise #3b

Reintegrate and reprocess the results

This exercise contains a series of tasks to help you learn to:

- Manually reintegrate the calibration standard results
- Change sample variable values
- Reprocess the sequence with the original method revision

You use the data produced in Exercise #3a.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Running Routine Samples”](#) on page 9.




Task 1. Make changes to the results and sample information

Steps

Detailed Instructions

- 1 Find the single injection result for the third quantitation of sample1_4 in the sequence exer3seqiii.

- a Select **Result** from the Current View.
- b From the Query List, select **MySeqNotApprovedRunLast7days**.
- c Expand the **exer3seqiii** folder.
- d Expand the **Calibration - exer3iii Calib Rev 4** folder.
- e Expand the **Samples** folder.
- f Expand the **sample 1_4** folder.
- g Select **sample 1_4#1**.




The chromatogram displays absorbance (mAU) on the y-axis (0 to 20) against time (min) on the x-axis (0 to 4). Four peaks are identified and labeled with their retention times: 0.94 (dimethylthalate), 1.11 (diethylthalate), 1.87 (biphenyl), and 3.05. The peak at 3.05 is the most prominent, reaching an absorbance of approximately 20 mAU. The peaks at 0.94 and 1.11 are clustered together and highlighted with a blue vertical bar. The peak at 1.87 is also labeled. The x-axis is labeled 'min' and the y-axis is labeled 'mAU'. The signal is labeled 'Signal' and 'VWD: Absorbance'.

RT	Compound Name	Amount	RF (Rsp/Amt)	Peak Area	Peak Height
0.94	dimethylthalate	2.4390	43.4855	106.0625	38.1547
1.11	diethylthalate	1.7108	43.4855	92.9914	31.3272
1.87	biphenyl	3.4804	24.5864	85.5719	23.7081
3.05		4.6565	27.3183	127.2068	22.5348

Basic Exercise #3b Reintegrate and reprocess the results

Task 1. Make changes to the results and sample information

Steps	Detailed Instructions
<p>2 Manually reintegrate the dimethylphthalate peak.</p> <p>Draw the baseline from the bottom left corner of the peak to the inflection point on the bottom right of the peak.</p> <p>Note that the Amount and RF values disappear.</p>	<p>a Click Manual Integration, and select Draw Peak Baseline.</p> <p>A mouse pointer in the shape of a bell curve appears on the chromatogram.</p> <p>b Place the pointer at the bottom left of the peak at the intersection between the baseline and peak, and click once.</p> <p>c Hold the mouse button down, and move the pointer to the inflection point at the bottom right of the peak.</p> <p>d Release the mouse button.</p> <p>The new baseline appears, but the bell curve pointer remains.</p> <p>e Click the Select Objects button to change the pointer from a bell curve pointer to a normal pointer.</p>



Basic Exercise #3b Reintegrate and reprocess the results

Task 1. Make changes to the results and sample information

Steps	Detailed Instructions
3 Change sample variable values. <ul style="list-style-type: none">• Dilution = 5• Purity = .9	<p>a Select the sequence, <i>exer3seq.iii</i>.</p> <p>The sequence table and Sample Entry panel appear in the workspace.</p> <p>b Select the first sample 1_4 in the sequence.</p> <p>c Click the Amounts tab, and enter a default value for the Dilution factor of 5.</p> <p>d Enter a default value for the Purity of .9, and click Apply.</p> <p>e Repeat steps c and d for every sample 1_4 in the sequence.</p>

Sequence TableSequence Options

	Sample Name	Sample Type	Cal. Level	Custom Sample Group	Vial #	Inject #
1	cal1	Calibration	1		2	1
2	sample 1_2	Sample			5	1
3	sample 1_4	Sample			9	1
4	cal1	Calibration	1		2	1
5	sample 1_2	Sample			5	1
6	sample 1_4	Sample			9	1
7	cal1	Calibration	1		2	1
8	sample 1_2	Sample			5	1
9	sample 1_4	Sample			9	1

Sample EntrySequence Logbook

Sample Name:
sample 1_4

Sample Type:
Sample

Custom Sample Group:

Apply

RunAmountsIdentificationDescription

Sample variables

Sample Amount: 0

Sample Amount U mg/ml


Multiplier: 1

Dilution: 5

Purity: .9

Task 2. Reprocess the sequence results

Steps	Detailed Instructions
1 Open the Reprocess window. See Chapter 3, "Sample Analysis", in the <i>Concepts Guide</i> for a chart that helps you select the correct reprocessing option.	a Select the sequence, <i>exer3seqiii</i> . The Save Reasons for Changes dialog box appears. b Enter any information requested, and click Save . c Select Actions > Reprocess in the top menu bar.
2 Select the reprocessing option that uses all other method settings of the original method, except for the integration settings and default sample variable values. In the Cerity system all sample, sequence, method and instrument information is attached to the result.	a Select Use the method revision now attached to the result . b Click OK . The Cerity system uses the settings of the method originally used to run the sequence, the new manual integration setting and the new sample variable values to process the sequence.



Basic Exercise #3b Reintegrate and reprocess the results

Task 2. Reprocess the sequence results

Steps	Detailed Instructions
-------	-----------------------

- 3 Track reprocessing until complete.
- a Select the sequence, exer3seqiii.

b Click the **Sequence Options** tab.

Sequence Table

Sequence Options

Sequence Name:
defseq3 - Reprocessed

Instrument:
EMELC3

Sequence Template

Apply

Sequence

Identification

Description

Report Destination

Run with

Priority:
Medium

Schedule:
Running Reprocessing

Calibration Mode:
Single Update Calibration

Sequence Created by

When the system has completed reprocessing, the message “Completed Reprocessing” appears on the Sequence Options panel

Sequence Table

Sequence Options

Sequence Name:
defseq3 - Reprocessed

Instrument:
EMELC3

Sequence Template

Apply

Sequence

Identification

Description

Report Destination

Run with

Priority:
Medium

Schedule:
Completed Reprocessing

Calibration Mode:
Single Update Calibration

Sequence Created by



Advanced Exercise #4a

Run a sequence to quantify compounds with multi-level calibration

This exercise contains a series of tasks to help you learn how to:

- Create a sequence with a method set up for multi-level, overall calibration, ESTD quantitation and variable compound amounts
- Enter new information for an individual sample or standard
- Edit a sequence during a run
- Review the results to view the multi-level, overall calibration process.
- View the early quantitation single injection reports and the sequence report

You can choose between two methods to use with this exercise:

- A copy of *defexer4iii*, the instrument method copied from the default method provided with the system.
- *Exer4iii*, the method that you created in “[Advanced Exercise #4 Set up a multi-level calibrated method for a sequence](#)” on page 105.

For tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read “[Running Routine Samples](#)” on page 9.

Equilibrate the instrument. See “[Basic Exercise #1 Equilibrate the instrument](#)” on page 11.



Advanced Exercise #4a Run a sequence to quantify compounds with multi-level calibration

Task 1. Create a new sequence, and enter sample and sequence information

Task 1. Create a new sequence, and enter sample and sequence information

Steps	Detailed Instructions
<p>1 Create a new sequence.</p> <p>Name the sequence <code>exer4seqiii</code>, where <i>iii</i> are your initials.</p> <p>Use one of the two methods:</p> <ul style="list-style-type: none">• <code>defexer4iii</code>• <code>exer4iii</code> (created with Exercise 4 of Setting Up Methods)	<ul style="list-style-type: none">• For detailed instructions, see “Task 1. Create a new sequence” on page 30. <p>After you create a new sequence, the revision number is set to 1.</p>
<p>2 Enter values for sample amounts and variables.</p> <p>For the first sample <code>1_2</code>, enter:</p> <ul style="list-style-type: none">• Sample Amount - 2.5 mg• Dilution Factor - 2• Purity - .93	<ul style="list-style-type: none">a Select Instrument from the Current View list.b Expand the instrument folder.c Select <code>exer4seqiii</code>.d Select the first sample <code>1_2</code> in the Sequence Table.e Click the Amounts tab.f Enter 2.5 for the Sample Amount.g Change the Dilution Factor value to 2.h Change the Purity value to .93.
<p>3 Enter compound amounts.</p> <p>To quantify a compound in a sample, you must select to use the compound amount for the standard.</p> <p>For the second set of calibration standards for dimethyl phthalate, enter compound amounts:</p> <ul style="list-style-type: none">• Cal1 - 10.17 µg• Cal2 - 37.62 µg	<ul style="list-style-type: none">a Click the Sequence Table tab, and select Cal1 from the second set of standards.b Enter 10.17 for the Compound amount.c Select Cal2 from the second set of standards.d Enter 37.62 for the Compound amount.

Sequence Table | Sequence Options

	Sample Name	Sample Type	Cal. Level	Immediate Quantitation	Custom Sample Group	Vial #	Injections #	Injection Volume [µl]
4	sample 1_4	Sample		YES		9	1	as method
5	cal1	Calibration	1	NO		9	1	as method
6	cal2	Calibration	2	NO		2	1	as method
7	sample 1_2	Sample		NO		5	1	as method
8	sample 1_4	Sample		NO		9	1	as method
9	cal1	Calibration	1	NO		9	1	as method
10	cal2	Calibration	2	NO		2	1	as method

Sample Entry | Sequence Logbook

Sample Name:
cal2

Sample Type:
Calibration Standard

Custom Sample Group:

Run | Amounts | Identification | Description

Sample variables



Sample Amount:	0
Sample Amount U	mg/ml
Multiplier:	1
Dilution Factor:	5

Compound amounts

Use	Name	Amount
<input checked="" type="checkbox"/>	dimethylphthal	39.75
<input checked="" type="checkbox"/>	biphenyl	60
<input type="checkbox"/>	diethylphthalat	0

Advanced Exercise #4a Run a sequence to quantify compounds with multi-level calibration

Task 1. Create a new sequence, and enter sample and sequence information

Steps	Detailed Instructions
4 Enter the tasks to be performed during the run: Quantify, Report, Allow Online Editing	<ul style="list-style-type: none"> a Select the sequence that you just created. b Click the Sequence Options tab. c Make sure that the Quantify and Report check boxes are marked for the Task(s) to perform. d Mark the Allow Online Editing check box. 
5 Enter the destination path for, but do not print, the reports: Enter Exercise4 <i>iii</i> , where <i>iii</i> are your initials.	<ul style="list-style-type: none"> • For the detailed instructions, see step 3 on page 32.
6 Save the sequence	<ul style="list-style-type: none"> • On the Standard toolbar, click , and enter reasons for changes and your password, if necessary. <p>After you save the sequence, the revision increments by one. Here, the revision number is set to 2.</p>

Advanced Exercise #4a Run a sequence to quantify compounds with multi-level calibration

Task 2. Edit the sequence during the run

Task 2. Edit the sequence during the run

Steps

Detailed Instructions

1

Run the sequence after the instrument is ready.

2

Edit the sequence during the run:

For detailed instructions, see “Task 3. Run and track the sequence” on page 33, Steps 1 and 2.

Note that the sequence disappears from beneath the instrument folder.

After you run the sequence, the revision number increments by one. Here, the revision number is set to 3.

After the last peak comes off during the run of the first standard, select to immediately quantify the first sample 1_4 in the sequence.

a

On the selection tree, select the instrument.


b

On the instrument workspace, click the **Worklist** tab.

c

Select the sequence.

d

After the last peak comes off during the run of the first standard, click  , in the toolbar.

The sequence in the worklist now says “Preparing to edit”. When the sample run is complete, the sequence is stopped and the status says “Editable”.

Instrument Panel

Worklist

	Name	Status	Type	Method	Priority #	Vial #	Injections #	Description
1	exer4seqdec	Preparing to edit...(1:1)	Sequence	exer4dec	500	N/A	N/A	

Instrument Panel

Worklist

	Name	Status	Type	Method	Priority #	Vial #	Injections #	Description
1	exer4seqdec	Editable	Sequence	exer4dec	500	N/A	N/A	

e

Expand the instrument folder. (Note that the sequence has reappeared.) If you do not see the sequence, click the **Redo Query** button or F5.

f

Select the sequence, and select the first sample 1_4 in the Sequence Table.

g

Double-click the **Immediate Quantitation** cell.

h

Double-click **Yes**.

i

Save and run the sequence.

The revision number increments to 4 after you save the sequence.

The revision number increments to 5 after you run the sequence.

j

Select the instrument and click the **Worklist** tab. (The sequence starts with the second standard.)

Instrument Panel

Worklist

	Name	Status	Type	Method	Priority #	Vial #	Injections #	Description
1	exer4seqdec	Running(2:1)	Sequence	exer4dec	500	N/A	N/A	

Task 3. Review the calibration results

Steps

1 Review the calibration table and curve.


If you ran the sample more than 7 days ago, you must modify the query to retrieve older results from the database. See the online *How To* help, "Define a query."

Note that when you first view the sequence result in the Result View, the revision number equals the number of saves that you made plus the number of run executions. In this exercise, the revision number for the sequence result is 5.

See Chapter 5, "Sample Analysis", in the *Concepts Guide* for information on sequence and calibration revisioning.

Detailed Instructions

- Select **Result** from the Current View.
- Select **AllSeqNotApprovedRunLast7Days** from the Query list.
- Expand the **exer4seqiii** folder.
One folder appears that contains the calibration and single injection results.
- Select any one of the **Calibration - exer4seqiii Calib Rev 5** folder.
The calibration table and curve appear in the workspace.



- View how the system uses the standards in overall calibration to quantify the samples compared to single-level calibration in Exercise 3a.

Advanced Exercise #4a Run a sequence to quantify compounds with multi-level calibration

Task 3. Review the calibration results

Steps


Detailed Instructions

2 Review the single injection results for both sample 1_2 injections.

Note that the Amount is different for the first sample 1_2. Why?

The Amount is the compound amount in the sample. The value for this exercise represents the compound amount in the injection times the values of the dilution factor and purity. When you entered this sample, you changed these values.


- Expand any one of the Calibration folders.
- Expand the **Samples** folder.
- Expand the first **sample 1_2** folder.
- Select the single injection.
- Note the value in the Amount column.



Chromatogram showing peaks for dimethylphthalate, diethylphthalate, and biphenyl. The x-axis is labeled 'min' and the y-axis is labeled 'mAU'. The peaks are labeled 1, 2, and 3.

RT	Compound Name	Amount	RF (Rsp/Amt)	Peak Area	Peak Height
0.93	dimethylphthalate	15.6528	10.9313	195.7383	69.9187
1.10	diethylphthalate	11.6001	10.9230	170.3073	57.6104
1.85	biphenyl	23.2052	6.1513	153.4873	42.9007
2.98		31.9479	6.8819	236.4102	42.5603

- Expand the second **sample 1_2** folder.
- Select the single injection.
- Compare the Amounts from the first and second sample1_2's.



Chromatogram showing peaks for dimethylphthalate, diethylphthalate, and biphenyl. The x-axis is labeled 'min' and the y-axis is labeled 'mAU'. The peaks are labeled 1, 2, and 3.

RT	Compound Name	Amount	RF (Rsp/Amt)	Peak Area	Peak Height
0.93	dimethylphthalate	44.7873	10.9132	195.5095	69.9118
1.10	diethylphthalate	31.1959	10.8913	169.8817	57.5508
1.86	biphenyl	62.4277	6.0983	152.2814	42.5782
2.91		96.0701	6.0980	226.6047	42.0420

Task 4. Review the reports


Steps	Detailed Instructions
1 Review the two single injection reports for the first sample 1_2 and sample 1_4. Note that there is only one folder for each of the second set of samples because they were not marked for Immediate Quantitation.	a Select Start > Programs > Agilent Cerity > Report Viewer . b Select File > Open , or click the Open button. c Expand the Exercise4iii folder. d Expand the 003 Multi-Injection Summary Group folder, and expand the 01 Sample single injection folder. e Double-click default.htm . Note the compound amounts. f Expand the 003 Multi-Injection Summary Group 0001 folder, and expand the 01 Sample single injection folder. g Double-click default.htm . Note the compound amounts. h Repeat steps d-g for the 004 Multi-Injection Summary Group folders.
2 View the sample amount for the first sample 1_2 in the Sequence Report.	a Click the Open button, and expand the Exercise4iii folder. b Expand the Sequence folder, and double-click default.htm.

Sequence samples

	Name	Position	Modified inj. volume	Amount	Unit	Cal. level
1	cal1	9	(As Method)	0.0000	mg/ml	1
2	cal2	2	(As Method)	0.0000	mg/ml	2
3	sample 1_2	5	(As Method)	2.5000	mg/ml	1
4	sample 1_4	9	(As Method)	0.0000	mg/ml	1
5	cal1	9	(As Method)	0.0000	mg/ml	1
6	cal2	2	(As Method)	0.0000	mg/ml	2
7	sample 1_2	5	(As Method)	0.0000	mg/ml	1
8	sample 1_4	9	(As Method)	0.0000	mg/ml	1
9	cal1	9	(As Method)	0.0000	mg/ml	1
10	cal2	2	(As Method)	0.0000	mg/ml	2

Advanced Exercise #4a Run a sequence to quantify compounds with multi-level calibration

Task 4. Review the reports



Advanced Exercise #4b

Change sample variables in the method and reprocess

This exercise contains a series of tasks to help you learn to:

- Change an integration setting in the method.
- Remove a calibration point.
- Change the sequence so that no sample is immediately quantified after processing
- Reprocess the sequence with the most current method revision.
- Add a new sample variable to the method.
- Reprocess the sequence after you add the new variable
- Regenerate the reports

You use the data produced in Exercise #4a.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Running Routine Samples”](#) on page 9.



Task 1. Update the method and result

Steps	Detailed Instructions
<p>1 Change the integration setting in the method.</p> <p>Set the Height Reject to 0.</p> <p>If you are using a copy of the defexer4iii method, make sure that no one else has modified it. Check the old revisions. If it has been modified, make another copy of the default method. See "Before you start" on page 5.</p>	<p>a Select Method from the Current View.</p> <p>b Expand the exer4iii folder.</p> <p>c Expand the Data Analysis folder.</p> <p>d Select Integration.</p> <p>e Click the Height Reject cell, and enter 0.</p> <p>f Save the method.</p>
<p>2 Remove the second Cal2 calibration point for dimethylphthalate.</p>	<p>a Select Result from the Current View.</p> <p>b Expand the exer4seqiii folder.</p> <p>c Select the Calibration - Exer4iii folder.</p> <p>d Click the Calibration cell for the second Cal2 calibration.</p> <p>e Click the .. button, and double-click the cell to change Yes to No.</p>

Compound Summary

dimethylphthalate

Peak Area

0 100

0 20

Weighed Amount

Peak Area

0 20

0 20

Weighed Amount

Sample Name	Valid Calibration	Weighed Amount
cal1 #1	Yes	10.0000
cal1 #1	Yes	10.1700
cal1 #1	Yes	10.0000
cal2 #1	Yes	40.0000
cal2 #1	No	37.6200
cal2 #1	Yes	40.0000

Advanced Exercise #4b Change sample variables in the method and reprocess

Task 1. Update the method and result

Steps

3

Change the sequence so that no sample is immediately quantified during processing

Detailed Instructions

a

Select the exer4seq*iii* sequence.

b

In the sequence table, double-click the **Immediate Quantitation** cell for the first Sample1_2.

c

Double-click **No**.

d

Repeat steps b and c for the first Sample1_4.

e

Save the changed result.


Note that the revision is incremented by 1.

Sequence Table

Sequence Options

	Sample Name	Sample Type	Cal. Level	Immediate Quantitation	Custom Sample Group	Vial #	Injections #
1	cal1	Calibration	1	NO		2	1
2	cal2	Calibration	2	NO		3	1
3	sample 1_2	Sample		NO		5	1
4	sample 1_4	Sample		NO		9	1
5	cal1	Calibration	1	NO		2	1

Task 2. Reprocess and review the result

Steps	Detailed Instructions
<p>1 Reprocess the sequence with the most current revision of the method.</p> <ul style="list-style-type: none">• Use the integration settings in the method.• Set up to print (regenerate) reports	<p>a Select the exer4seqjws2 sequence.</p> <p>b Select Actions > Reprocess.</p> <p>c Select Use the most current revision of the method that is attached to the result.</p> <p>d Mark the Use integration settings in the method check box.</p> <p>e Mark the Print Reports check box.</p> <p>f Click OK.</p> <p>g To follow reprocessing, click the Sequence Options tab.</p>
	
<p>2 Make sure the integration change appears in the reprocessed result.</p> <p>If you cannot see the calibration standard chromatogram because of the example chromatogram, click the Layout button and clear the Display Example Chromatogram check box.</p>	<p>a Expand the second Calibration - Exer4jws2 folder.</p> <p>b Expand the Calibrations folder and the Cal1 folder.</p> <p>c Select Cal1 #1.</p> <p>Note that one or more peaks are now integrated and appear in the Results Table.</p>

Advanced Exercise #4b Change sample variables in the method and reprocess

Task 2. Reprocess and review the result

Steps

Detailed Instructions

3 Review the calibration summary.

Select the second **Calibration - Exer4iii** folder.

Note that the calibration point that you removed before reprocessing is gone.

Compound Summary

dimethylphthalate

Peak Area

Weighed Amount

Area

Weighed Amount

Sample Name	Valid Calibration	Weighed Amount
cal1 #1	Yes	10.0000
cal1 #1	Yes	10.1700
cal1 #1	Yes	10.0000
cal2 #1	Yes	40.0000
cal2 #1	Yes	40.0000

4 Review the reports for the first set of samples to make sure that they were quantified with all the calibration standards.


a Select **Start > Programs > Agilent Cerity > Report Viewer**.

b Select **File > Open**.

c Expand **Exercise4iii-0001**.

Note that only one report exists for all the samples. After initial processing two reports existed for the first Sample1_2 and Sample1_4.

Task 3. Add a new sample variable to the method and reprocess

Steps	Detailed Instructions
1 Add a new variable to the method. Add a divisor called "attenuation factor" with a default value of 3.	<ul style="list-style-type: none">a Select Method from the Current View list.b Expand the current revision of exer4iii.c Select Sample Variables.d Type "attenuation factor" into a Divider cell of the System Sample Variables table.e Enter a Default Value of 3.f Save the method.
2 Reprocess the sequence with the revised method. <ul style="list-style-type: none">Enter a new value for the Attenuation Factor of 7 for the first Sample 1_2.Set up to print (regenerate) reports.	<ul style="list-style-type: none">a Select Result from the Current View list.b Select exer4seqiii.c Select Actions > Set up reprocessing for new sample entry fields. <div></div> <ul style="list-style-type: none">d Click OK. The new Sample Entry panel appears.e Click the Amounts tab, and enter 7 for the "Attenuation Factor".f Select Actions > Reprocess.g Select Use the method revision now attached to the result.h Mark the Print Reports check box.i Click OK.

Advanced Exercise #4b Change sample variables in the method and reprocess

Task 3. Add a new sample variable to the method and reprocess

Steps

Detailed Instructions

3 Find the report for the first Sample1_2.

Note that the quantitation value is different after reprocessing. The software used the “Attenuation Factor” in the calculation.

a Select **Start > Programs > Agilent Cerity > Report Viewer.**

b Select **File > Open.**

c Expand **Exercise4iii-0002.**

d Expand the 003Multi-Injection Summary folder.

e Expand the 01Sample Single Injection folder.

f Double-click Default.htm.

The report appears with the new amount for Sample1_2.

The chromatogram displays absorbance (mAU) on the y-axis (0 to 30) against time (min) on the x-axis (0 to 4). Five peaks are identified with their retention times: 0.76, 0.94, 1.11, 1.89, and 3.11. The peak at 1.89 min is labeled 'biphenyl'. The plot title is 'VWD: Absorbance'.

Sample single injection compounds

RT	Compound	Peak area	Amount	Unit	Resp. f.	Tailing f.
0.76	N/A	0.8370	0.4500	N/A	0.2471	N/A
0.94	dimethylphthalate	124.1833	2.4779	ug	6.6582	N/A
1.11	diethylphthalate	109.6416	1.7791	N/A	6.5501	N/A
1.89	biphenyl	106.8904	3.7001	ug	3.8380	N/A
3.11	N/A	153.0533	4.5837	N/A	4.4362	N/A

Advanced Exercise #4b Change sample variables in the method and reprocess

Task 3. Add a new sample variable to the method and reprocess



Advanced Exercise #5a

Run a sequence to quantify impurities

This exercise contains a series of tasks to help you learn to review results and reports of a sequence run with a method set up for multi-level, bracketed calibration, ISTD quantitation and variable compound amounts. You learn how to:

- Recognize the results of an overall calibration
- Find the system suitability calculations that were selected for the review layout in the method
- Find the custom calculations that were set up in the method
- Review the reports for the calculations that were set up in the report template

You can choose between two methods to use with this exercise:

- instrument method copied from the default method provided with the system, defexer5.
- method that was created in [“Advanced Exercise #5 Set up a method for a sequence to quantify impurities”](#) on page 117.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Running Routine Samples”](#) on page 9.

Equilibrate the instrument. See [“Basic Exercise #1 Equilibrate the instrument”](#) on page 11.



Advanced Exercise #5a Run a sequence to quantify impurities

Task 1. Set up and run the sequence

Task 1. Set up and run the sequence

Steps	Detailed Instructions
<p>1 Create a new sequence.</p> <p>Name the sequence <code>exer5seq<i>iii</i></code>, where <i>iii</i> are your initials.</p> <p>Use one of the two methods:</p> <ul style="list-style-type: none">• <code>defexer5</code>• <code>exer5<i>iii</i></code> (created with Exercise 5 of Setting Up Methods)	<ul style="list-style-type: none">• For detailed instructions, see “Task 1. Create a new sequence” on page 30.
<p>2 Make sure that Quantify and Report are selected.</p>	<ul style="list-style-type: none">• For detailed instructions, see “Task 2. Enter sample and sequence information” on page 31, Step 2.
<p>3 Enter the destination path for, but do not print, the reports, and save the sequence.</p> <p>Enter <code>Exercise5<i>iii</i></code>, where “<i>iii</i>” are your initials.</p>	<ul style="list-style-type: none">• For detailed instructions, see “Task 2. Enter sample and sequence information” on page 31, Step 3.
<p>4 Run and track the sequence.</p>	<ul style="list-style-type: none">• For detailed instructions, see “Task 3. Run and track the sequence” on page 33.

Task 2. Review the results and reports

Steps

- 1 Compare the response factors for dimethylphthalate for the first set of bracketed samples with the second set.

Hint: If you can't see the RFs, click the bottom of the Compound Summary panel to make the scroll bar appear.

Note that the RFs for the second Cal1 and Cal2 for the first set of bracketed samples are the same as for the first Cal1 and Cal2 for the second set of bracketed samples.

Detailed Instructions

- a Select **Result** from the Current View.
- b Select **AllSeqNotApprovedRunLast7Days** from the Query list.
- c Expand the **exer3seq** folder.
- d Select the second **Calibration - exer3seq** folder.
The first calibration folder contains the blank run.
- e Scroll to see the RFs if not visible.
- f Select the third **Calibration - exer5seq** folder.
- g Scroll to see the RFs if not visible.
- h Compare the RFs.

Sample Name	Weighed Amount	RF (Rsp/Amt)
cal1 #1	10.0000	1.7832
cal1 #1	10.0000	1.7784
cal2 #1	40.0000	1.7247
cal2 #1	40.0000	1.7271

Weighed Amount	RF (Rsp/Amt)
10.0000	1.7784
10.0000	1.7727
40.0000	1.7271
40.0000	1.7248

- 2 Review the system suitability calculations for Cal1 #1 under the second calibration folder.

Note the values for the Average Percent Specified Impurity and the Average Percent Unspecified Impurity calculations that were set up as custom calculations in the method.

- a Expand the second **Calibration - exer3seq** Calib folder.
- b Expand the **Calibrations** folder.
- c Expand the Cal1 folder.
- d Select Cal1 #1.
- e Review the Results Table for the system suitability calculations.

You may have to click the bottom of the Results table to see the scroll bar.

Results					
RT	Compound Name	Peak Width	TailingFactor	SignalToNoise	Peak resolution USP
0.94	dimethylphthalate	0.0424	1.144	97.300	N/A
1.11	diethylphthalate	0.0443	1.050	79.413	2.303
1.89	biphenyl	0.0560	0.887	1041.299	9.108
3.10		0.0905	0.666	607.791	9.690

Summary Results	
Percent Specified Impurity:	13.42
Percent Unspecified Impurity:	37.91

Advanced Exercise #5a Run a sequence to quantify impurities

Task 2. Review the results and reports

Steps	Detailed Instructions
3 Review the percent impurity results for the first Sample1_2 and for the sample group. Note that the percent impurity values exceeded their limits.	a Expand the second Calibration - exer3seq<i>iii</i> folder. b Expand the Samples folder. c Select the Sample1_2 folder. Note that the average percent specified and unspecified impurities for both injections appear here.

Results Table	
Compound Name	Injection#
dimethylphthalate	1
	2
diethylphthalate	1
	2
biphenyl	1
	2
Not Identified Peaks	1
	2

Summary Results	
Avg Percent Specified :	13.65
Avg Percent Unspecified :	37.80

- d** Expand the **Group Results** folder.
- e** Select **Samples**.
- The results for the average of the percent impurities over all samples appear here, as do the results of the limit checks for these impurities.

Summary Results	
Avg % S All Samples :	13.73
Avg % S All Samples Limit Check :	Not Passed
Avg % U All Samples :	37.72
Avg % U All Samples Limit Check :	Not Passed

Advanced Exercise #5a Run a sequence to quantify impurities

Task 2. Review the results and reports

Steps

Detailed Instructions

4

Review the sample single injection report for the first Sample1_2 and the report for the sample group.

a

Select **Start > Programs > Agilent Cerity > Report Viewer**.

b

Select **File > Open**.

c

Expand **Exercise5iii**.

d

Expand 003Multi-InjectionSummary.

e

Expand 01Sample Single Injection, and double-click default.htm.

Note the system suitability calculation values in the table that was set up in the method.

Retention Time	Compound Name	Amount	Response Factor	Tailing Factor	Peak resolution USP	SignalToNoise
0.93	dimethylphthalate	24.8892	0.1169	1.178	N/A	237.192
1.10	diethylphthalate	17.5561	0.1169	1.135	2.308	194.383
1.89	biphenyl	37.5000	0.0667	1.090	9.129	2554.088
3.11	N/A	48.6177	0.0741	1.043	9.713	1489.322

f

Expand **Exercise5iii**.

g

Expand **Sample Group**, and click default.htm.

Note the percent impurity calculations and limits that were set up in custom calculations and the report template in the method.

Avg % S All Samples:	13.73
Avg % U All Samples:	37.72

Sample group limit results


#	Sample name	Compound	Limit (Compound)	Limit (Sample)
1	sample 1_2	dimethylphthalate	XXXXXXXXXX	XXXXXXXXXX
2	sample 1_4	dimethylphthalate	XXXXXXXXXX	XXXXXXXXXX
3	sample 1_2	dimethylphthalate	XXXXXXXXXX	XXXXXXXXXX
4	sample 1_4	dimethylphthalate	XXXXXXXXXX	XXXXXXXXXX
1	sample 1_2	diethylphthalate	XXXXXXXXXX	XXXXXXXXXX
2	sample 1_4	diethylphthalate	XXXXXXXXXX	XXXXXXXXXX
3	sample 1_2	diethylphthalate	XXXXXXXXXX	XXXXXXXXXX
4	sample 1_4	diethylphthalate	XXXXXXXXXX	XXXXXXXXXX
1	sample 1_2	biphenyl	XXXXXXXXXX	XXXXXXXXXX
2	sample 1_4	biphenyl	XXXXXXXXXX	XXXXXXXXXX
3	sample 1_2	biphenyl	XXXXXXXXXX	XXXXXXXXXX
4	sample 1_4	biphenyl	XXXXXXXXXX	XXXXXXXXXX

Avg % S All Samples Limit Check: Not Passed

Avg % U All Samples Limit Check: Not Passed

Advanced Exercise #5a Run a sequence to quantify impurities

Task 2. Review the results and reports



Advanced Exercise #5b

Use a different method to reprocess

This exercise contains a series of tasks to help you learn to:

- set up a different method with a new calibrated compound
- set up reprocessing for a different method
- reprocess the sequence with the different method

You use the data produced in Exercise #5a.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before You Start

Read [“Running Routine Samples”](#) on page 9.



Advanced Exercise #5b Use a different method to reprocess

Task 1. Set up a different method

Task 1. Set up a different method

Steps	Detailed Instructions
1 Copy exer5iii and rename it exer5iii2. <ul style="list-style-type: none">• Or, copy defexer5.• Or, use defexer5iii2 to reprocess.	<ul style="list-style-type: none">a Select File > New > Method.b Click the Browse button in the Method Wizard.c Select exer5iii.d Enter a New Method Name of exer5iii2, and click Next.e Click Next until you reach the New Method Review panel.f Click Finish and click Save.
2 Add diethylphthalate as a calibrated compound. <ul style="list-style-type: none">• Cal Level 1 - 8 µg• Cal Level 2 - 32 µg• Set biphenyl as the ISTD for this compound.	<ul style="list-style-type: none">a Expand the exer5iii2 folder.b Expand the Data Analysis folder.c Select Calibration.d Right-click the calibration table, and select Insert Compound.e Select diethylphthalate, click > and click OK.f In the calibration table, select diethylphthalate.g Click on the Level 1 Use Default Amount cell, and click the .. button.h Select the + sign, and enter 8 µg into the Weighed Amount and Unit cells.i Repeat steps g and h for Level 2 and 32 µg.j Select Quantitation.k Select diethylphthalate.l Mark the Use ISTD Compound check box, and select biphenyl.m Save the method.

Task 2. Reprocess the sequence result

Steps


Detailed Instructions

1 Set up reprocessing for a different method.

Select **exer5iii2** or **defexer5iii2**.


See Chapter 3, "Sample Analysis", in the *Concepts Guide* for a chart that helps you select the correct reprocessing options.

- Select **Result** from the Current View.
- From the Query List, select **MySeqNotApprovedRunLast7days**.
- Select the **exer5seqiii** folder.
- Select **Actions > Set up reprocessing for a different method**.



- Click **Browse**, select **exer5iii2** and click **OK**.
- Click **OK**, and click **Save**.

A copy of the sequence appears in the selection tree, ready for reprocessing. This copy is now attached to the new method but has no folders underneath until it is reprocessed.



2 Enter amounts for the new calibrated compound, diethylphthalate, for every calibration standard.

Level 1 - 8

Level 2 - 32

- Select this copy (note the date and time after it).
- Click the **Amount** tab on the Sample Entry panel in the sequence workspace.
- For each Level 1 standard, mark the **Use** checkbox for diethylphthalate, and enter 8.
- For each Level 2 standard, mark the **Use** checkbox for diethylphthalate, and enter 32.
- Save the result.

Advanced Exercise #5b Use a different method to reprocess


Task 2. Reprocess the sequence result

Steps	Detailed Instructions
3 Reprocess the copy.	<div>a Select Actions > Reprocess</div> <div>b Make sure that Use the method revision now attached to the result is marked.</div> <div>c Click OK.</div> <div>d Monitor the reprocessing from the Sequence Options panel.</div> <div>e Click the Redo Query button.</div> <div>f Expand the copy.</div> <div>g Select a calibration folder.</div> <div>h Make sure that diethylphthalate is now included as a calibrated compound.</div>

The screenshot displays the software interface with a file tree on the left and a main panel on the right. The file tree shows a sequence named 'AllSeqNotApprovedRunLast7Days' with several sub-items, including 'Copy Of exer5seqv2 - Reprocessed (4/12)' and 'Calibrations'. The main panel shows a 'Calibration Table' with columns for 'Compound Name', 'Weighed Amount', 'Comment', 'Signal Short Description', and 'RF (Rsp/Amt)'. The table lists compounds: dimethylphthalate, diethylphthalate, and biphenyl. Below the table is a 'Compound Summary' section for 'diethylphthalate' showing a graph of 'Area Ratio (STD)' vs 'Weighed Amount' and a table of sample data.

Compound Name	Weighed Amount	Comment	Signal Short Description	RF (Rsp/Amt)
dimethylphthalate	10.0000		Vw/D1 A	1.7832
	40.0000			1.7271
diethylphthalate	8.0000		Vw/D1 A	1.9196
	32.0000			1.9107
biphenyl	15.0000		Vw/D1 A	1.0000
	60.0000			1.0000

Sample Name	Valid Calibration	Weighed #
cal1 #1	Yes	8.000
cal1 #1	Yes	8.000
cal2 #1	Yes	32.000
cal2 #1	Yes	32.000



Setting Up Methods

These exercises help you learn how to set up methods for your laboratory. See Chapter 4, “Method Setup”, in the *Concepts Guide* for background information that can help you use these exercises. The set of basic and advanced exercises includes the following topics:

Basic

Exercise 1 – Set up an equilibration method Learn how to set up a method template and enter instrument parameters to equilibrate the instrument.

Exercise 2 – Set up a method for single samples to identify compounds Learn how to use an example chromatogram to set up integration and compound identification for single samples.

Exercise 3 – Set up a single-level calibrated method for a sequence Learn how to set up single-level, single-update calibration, ESTD quantitation, and fixed compound amounts.

Advanced

Exercise 4 – Set up a multi-level calibrated method for a sequence Learn how to set up multi-level, overall calibration, ESTD quantitation, variable compound amounts, and sample variables.

Exercise 5 – Set up a method for a sequence to quantify impurities Learn how to setup ISTD quantitation, custom calculations, limits, bracketed calibration, and system suitability.

After you set up the methods in Exercises 1-5, you can use them to run the samples and sequences in Exercises 1-5 of the section—[“Running Routine Samples”](#).



Agilent Technologies

Before you start Read **“Before you start”** on page 5!

Your system administrator must have configured an Agilent 1100 LC series liquid chromatograph for your system.

If you plan to copy a default method to create a new method as in Exercises 3 and 5, make sure that the default methods are in your database. From the Query list, select AllMethodsRestored to view defexer1-5. If they do not appear, see the instructions in **“Before you start”** to transfer these methods from the CD-ROM to your database.



Basic Exercise #1

Set up an equilibration method

This exercise provides a series of tasks to learn how to:

- Create a method template to set up instrument parameters
- Set up instrument parameters
- Save and audit method changes
- View the history of method changes

A *method template* is a framework to let you enter only those conditions and parameters that you need to acquire and process data. A *method* is a method template that contains entered parameter values.

Use this method to equilibrate the instrument as outlined in the chapter “[Basic Exercise #1 Equilibrate the instrument](#)” on page 11.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start


Read “[Setting Up Methods](#)” on page 71 for setting up methods.



Basic Exercise #1 Set up an equilibration method

Task 1. Create a method template to enter instrument parameters

Task 1. Create a method template to enter instrument parameters


Steps	Detailed Instructions
<p>1 Create a new method template for a single sample.</p> <ul style="list-style-type: none">Name the method template, equilmeth<i>iii</i>, where <i>iii</i> are your initials.	<p>a Select File > New > Method or click  and select Method. The Method Wizard appears.</p> <p>b On the New Method panel, enter the Method Name, equilmeth<i>iii</i>.</p> <p>c Select Single Sample.</p> <div data-bbox="534 519 1195 1036"></div> <p>d Click Next to scroll to the Instrument panel.</p>

Basic Exercise #1 Set up an equilibration method

Task 1. Create a method template to enter instrument parameters


Steps	Detailed Instructions
-------	-----------------------

- 2 Select the instrument to equilibrate.**
- a** On the Instrument panel, select the instrument you need to equilibrate. The instruments that appear in the **Available Instruments** list depends on your configuration of the Cerity Networked Data System.



- b** Click **Next** to scroll to the Data Analysis panel.


- 3 Clear all data analysis selections.**
- a** On the Data Analysis panel, clear the **Compound Identification** check box.



- b** Click **Next** to scroll to the New Method Review panel.

Basic Exercise #1 Set up an equilibration method

Task 1. Create a method template to enter instrument parameters

Steps	Detailed Instructions
4 Review and save the method template.	<div><div>a On the New Method Review panel, review the setting in the Method Wizard Settings section.</div><div>b Add the words "Test Comment" in the Comment section.</div><div>c Click Finish.</div></div> <div></div> <div>d Click Save if the Save Changes to the Database dialog box appears.</div>


- 5 View the Method Wizard settings in the method.

After you save the method template, the Method View appears.

a Select the method you just created - *equilmethiii*.

b View the **Method Description** in the workspace.

Notice that the Method Description corresponds to the Comment section of the New Method Review panel in the Method Wizard.



Task 2. Enter the instrument conditions for the equilibration

Steps	Detailed Instructions
1 Set the pump parameters: Methanol as Solvent B: <ul style="list-style-type: none"> Flow rate: 2ml/min. Solvent composition: 80%MeOH/20%H₂O Stoptime: 10 min. Acetonitrile as Solvent B: <ul style="list-style-type: none"> Flow rate: 1.5ml/min Solvent composition: 65%ACN/35%H₂O Stoptime: 10 min. 	<ul style="list-style-type: none"> a On the selection tree, expand the <i>equilmethiii</i> method folder. b Expand the Instrument Setup folder and select Quaternary Pump or Binary Pump. c Enter the Flow as 2 ml/min. d Under Solvents, mark the B check box and enter 80 in the % box. The percent of solvent A is automatically set to 20 %. e Under Stoptime, select the min option and enter 10. f Under Posttime and Pressure Limits, accept the default values.
2 Set the autosampler (ALS) injection volume to zero	<ul style="list-style-type: none"> a Select the ALS folder. b Click the Setup tab. c Under Injection, select Standard Injection. d Set the Injection Volume to zero.

Setup | Timetable | Auxiliary & Data Curves

Flow: 2 ml/min

Solvents:

A: 20 %

B: ☒ 80 %

C: ☐ Off

D: ☐ Off

Stoptime: ☐ no Limit ☒ 10 min

Posttime: ☒ Off ☐ 0 min

Pressure Limits: Min: 0 bar Max: 400 bar

Setup | Auxiliary & Time

Injection:

☒ Standard Injection Injection Volume: 0 µl

☐ Injection with Needle Wash Wash Vol: 1

☐ Use Injector Program

Basic Exercise #1 Set up an equilibration method

Task 2. Enter the instrument conditions for the equilibration

Steps	Detailed Instructions
3 Set the same stoptime for all modules. Stoptime: 10 min.	<ul style="list-style-type: none">a Select the ALS folder,b Click the Auxiliary & Time tab.c Under Stoptime, select the as Pump option.d Select the DAD, MWD, or VWD folder that appears in your detector configuration.e Under Stoptime, select the as Pump/Injector option.f Select the TCC folder.g Under Stoptime, select the as Pump/Injector option.h Accept default values for all other module parameters


Task 3. Save and audit method changes

Steps


Detailed Instructions

1 Save the method.

The Cerity administrator must set up auditing for the **Save Changes To The Database** dialog box to appear. The Cerity administrator can provide a list of reasons and require you to enter your electronic signature to complete the dialog box. These requirements can only appear when a Cerity GMP license is installed and Auditing is set up by the Cerity administrator.

- a** On the Standard toolbar, click .

The **Save Changes To The Database** dialog box appears.



The dialog box titled "Save Changes To The Database" contains a list of changes, a reason for changes field, and buttons for Save, Discard, and Cancel.

List of changes

- Change the 'Stoptime' from 'no Limit' to 'as Pump/Injector' for the TCC Setpoint.
- Change the 'Stoptime' from 'no Limit' to 'as Pump/Injector' for the VWD Setpoint.
- Change the 'Flow' from '0' to '2' for the Quaternary Pump Setpoint.
- Change the 'Stoptime' from 'no Limit' to '5' for the Quaternary Pump Setpoint.
- Change the 'Solvent D Ratio' from '0' to 'off' for the Quaternary Pump Setpoint.
- Change the 'Solvent C Ratio' from '0' to 'off' for the Quaternary Pump Setpoint.
- Change the 'Solvent B Ratio' from '0' to '80' for the Quaternary Pump Setpoint.
- Change the 'Solvent A Ratio' from '100' to '20' for the Quaternary Pump Setpoint.
- Change the 'Injection Volume' from '5' to '0' for the ALS Setpoint.
- Change the 'Stoptime' from '10' to 'as Pump' for the ALS Setpoint.

Reason for changes


Save Discard Cancel

- b** Review the **List of changes**.
- c** Under **Reason for changes**, enter a reason or select a reason from the list.
- d** Click the **Save** button.

2 View the history of changes to the method.

If you need to use this method before you set up other methods, use the method with Running Routine Samples, Basic Exercise #1, Equilibrate the instrument.

- a** On the selection tree, select the method, *equilmethiii*.
- b** View the list of changes to the method.



Description	Item	Comment	E-Sig	Timestamp
Change the 'Stoptime' from 'no Limit' to 'as Pump/Injector' for the TCC Setpoint.	TCC Setpoint	Initial configuration	None	03/17/2002, 16:31:51
Change the 'Stoptime' from 'no Limit' to 'as Pump/Injector' for the VWD Setpoint.	VWD Setpoint	Initial configuration	None	03/17/2002, 16:31:51
Change the 'Flow' from '0' to '2' for the Quaternary Pump Setpoint.	Quaternary Pump Setpoint	Initial configuration	None	03/17/2002, 16:31:51
Change the 'Stoptime' from 'no Limit' to '5' for the Quaternary Pump Setpoint.	Quaternary Pump Setpoint	Initial configuration	None	03/17/2002, 16:31:51

Individual setpoint changes can only appear in the history of changes when a Cerity GMP license is installed and Auditing is set up by the Cerity administrator.

Basic Exercise #1 Set up an equilibration method

Task 3. Save and audit method changes



Basic Exercise #2

Set up a method for single samples to identify compounds

This exercise contains a series of tasks to learn how to:

- Create a method template for single samples to include only compound identification in the method
- Set up and save the method to produce an example chromatogram
- Use an example chromatogram to set up integration
- Set up compound identification

A *method template* is a framework to let you enter only conditions and parameters that you need to acquire and process data.

Use the method created in the first part of this exercise to enter and run a single sample to produce an example chromatogram. You can use the completed method to enter and run a group of samples to identify compounds. See [“Basic Exercise #2a Run a single sample to produce an example chromatogram”](#) on page 17 and [“Basic Exercise #3b Reintegrate and reprocess the results”](#) on page 39.

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Setting Up Methods”](#) on page 71 for setting up methods.



Basic Exercise #2 Set up a method for single samples to identify compounds

Task 1. Create a method template to identify compounds only

Task 1. Create a method template to identify compounds only

Steps

Detailed Instructions

1 Create a new method template for a single sample.


- Name the method template, *exer2iii*, where *iii* are your initials.

- c** Select **File > New > Method** or click  and select **Method**.

The Method Wizard appears.

- d** Enter *exer2iii* in the **Method Name** box.


- e** Select **Single Sample**.



- f** Click **Next** to scroll to the Method Wizard Instrument panel.

2 Select an instrument for the method.


- a** On the Instrument panel, select the instrument that will run the sample.



- b** Click **Next** to scroll to the Data Analysis panel.

Basic Exercise #2 Set up a method for single samples to identify compounds



Task 1. Create a method template to identify compounds only

Steps	Detailed Instructions
3 Mark only Compound Identification.	<p>a On the Data Analysis panel, Clear the Calibration and Quantification and Include System Suitability Calculations check boxes.</p>  <p>b Click Next to scroll to the Identification panel.</p>
4 Complete the setup of the method template. Do not mark any check box on the Method Wizard Identification panel.	<p>a Click Next, and click the Finish button.</p> <p>b Click Save if the Save Changes to the Database dialog box appears.</p>

Basic Exercise #2 Set up a method for single samples to identify compounds

Task 2. Enter the instrument conditions for the equilibrium

Task 2. Enter the instrument conditions for the equilibrium

Steps	Detailed Instructions
<p>1 Enter the pump parameters:</p> <p>Methanol as Solvent B:</p> <ul style="list-style-type: none">Flow rate: 2ml/min.Solvent composition: 80%MeOH/20%H₂OStoptime: 5 min. <p>Acetonitrile as Solvent B:</p> <ul style="list-style-type: none">Flow rate: 1.5ml/minSolvent composition: 65%ACN/35%H₂OStoptime: 6 min.	<p>a On the selection tree, expand the exer2iii method folder.</p> <p>b Expand the Instrument Setup folder and select the Quaternary Pump or Binary Pump.</p> <p>c Enter the Flow as 2 ml/min.</p> <p>d Under Solvents, mark the B check box and enter 80 in the % box. The percent of solvent A is automatically set to 20 %.</p> <p>e Under Stoptime, select the min option and enter 5.</p>
	
<p>2 Enter the injection volume and stop time for the autosampler.</p> <p>Injection Volume: 1µl</p> <p>Stop Time: same as pump</p>	<p>a On the selection tree, select the ALS folder.</p> <p>b Click the Auxiliary & Time tab.</p> <p>c Under Stoptime, select the as Pump option.</p> <p>d Click the Setup tab and select Standard Injection.</p> <p>e Enter 1µl for the Injection Volume.</p>
	

Basic Exercise #2 Set up a method for single samples to identify compounds

Task 2. Enter the instrument conditions for the equilibrium

Steps

Detailed Instructions

3 Make sure the stop time is the same for all instrument modules.

Stop Time: same as pump

- a On the selection tree, select the **VWD** folder.
- b Under **Stoptime** select **as Pump/Injector**.
- c On the selection tree, select the **TCC** folder.
- d Under **Stoptime**, select **as Pump/Injector**.

Signal & Time | Timetable | Options | Special Setpoints

Signal

Wavelength: 254 nm

Peakwidth (Responsetime): >0.10 min (2.0 s)

Stoptime:

☒ as Pump / Injector

☐ no Limit

☐ 0 min

Posttime:

☒ Off

☐ 0 min

Basic Exercise #2 Set up a method for single samples to identify compounds

Task 3. Save and audit method changes

Task 3. Save and audit method changes

Steps


1 Save the method.

After you save the method here, you can use the method to produce an example chromatogram.


See “Basic Exercise #2a Run a single sample to produce an example chromatogram” on page 17.

Continue with Task 4 after you produce an example chromatogram.

Detailed Instructions

- a On the Standard toolbar, click .

The **Save Changes To The Database** dialog box appears.



- b Review the **List of changes**.
c Under **Reason for changes**, enter a reason or select a reason from the list.
d Click the **Save** button.

The Cerity administrator must set up auditing for the **Save Changes to the Database** dialog box to appear. The Cerity administrator can provide a list of reasons and require you to enter your electronic signature to complete the dialog box.

Task 4. Select an example chromatogram and set up integration


Steps


1 Select an example chromatogram.

If no chromatogram of the isocratic sample exists, you must run a sample to produce the example chromatogram. See “Basic Exercise #2a Run a single sample to produce an example chromatogram” on page 17.

You do not need the example chromatogram to set up integration and identification, but it is recommended.


Detailed Instructions

- On the selection tree, expand the *exer2iii* method folder, if necessary.
- Expand the **Data Analysis** folder.
- Select **Example Chromatogram**.
- On the **Tools** toolbar, click .




- Expand the **Samples** folder.
- Expand the *exchromiii* or *defexchrom2a* folder.
- Select the sample name with the injection number.
- Click the **Select** button.

The example chromatogram appears in the workspace.




Basic Exercise #2 Set up a method for single samples to identify compounds

Task 4. Select an example chromatogram and set up integration

Steps	Detailed Instructions
2 Change the initial event values so that there are only four integrated peaks.	<div>a On the selection tree, select Integration under Data Analysis. The example chromatogram appears with the integration events tables.</div> <div>b Change the Height Reject event value to 1 (or the lowest value that will still integrate the four main peaks).</div> <div>c Click  on the Actions toolbar</div>

Example Chromatogram

VWD: Absorbance



Initial Events Timed Events

VWD Select...

Initial Event Name	Initial Event Value
Area Reject	0.0000
Slope Sensitivity	1.0000
Peak Width	0.0400
Shoulder Detection Mode	Disabled
Height Reject	1.0000


For All Signals

Tail Peak Skim Height Ratio	0.0000
Front Peak Skim Height Ratio	0.0000
Skim Valley Ratio	20.0000
Baseline Correction	Classical
Tangent Skim Mode	Standard
Peak to Valley Ratio	500.0000

Results


RT	Signal Short Description	Peak Area
0.9349	VWD1 A	419.5843
1.1044	VWD1 A	374.2865
1.8794	VWD1 A	356.2544
3.0739	VWD1 A	523.9493

Task 5. Set up compound identification

Steps	Detailed Instructions
<div><div>1</div><div>Set up the compound table for the following compounds: RT=.9 to 1.1, dimethylphthalate RT=1.1 to 1.2, diethylphthalate RT=1.8 to 2.1, biphenyl RT=3 to 3.2, o-terphenyl</div></div>	<div><div>a</div><div>On the selection tree, select the Identification item for Data Analysis.</div><div>b</div><div>On the Tools toolbar, click  . The peaks appear with the names New CompoundN in the compound table, where N = 1 - 4.</div><div>c</div><div>Under Compound Name, select the first cell and enter dimethylphthalate. After you select the cell, enter the name. The previous entry is overwritten.</div><div>d</div><div>Under Compound Name, select the second cell and enter diethylphthalate.</div><div>e</div><div>Under Compound Name, select the third cell and enter biphenyl.</div><div>f</div><div>Under Compound Name, select the fourth cell and enter o-terphenyl.</div></div>

CompoundOptions

VWD: Absorbance



Compound Name	Expected Time	Peak Signal	Time Reference Peak	Use Default Time Window	Low Time Limit	High Time Limit
dimethyl phthalate	0.9908	VwD1 A	-	+	0.9661	1.01
diethyl phthalate	1.1668	VwD1 A	-	+	1.1376	1.19
biphenyl	1.9700	VwD1 A	-	+	1.9207	2.01
o-terphenyl	3.1861	VwD1 A	-	+	3.1065	3.26

Basic Exercise #2 Set up a method for single samples to identify compounds


Task 5. Set up compound identification

Steps


Detailed Instructions

2 Save the method.

If you need to run this method to identify compounds before you set up the other methods in these exercises, use the method with “[Basic Exercise #2b Run a group of single samples to identify compounds](#)” on page 23.


- a On the Standard toolbar, click .

The **Save Changes To The Database** dialog box appears.



- b Review the **List of changes**.
c Under **Reason for changes**, enter a reason or select a reason from the list.
d Click the **Save** button.

The Cerity administrator must set up auditing for the **Save Changes to the Database** dialog box to appear. The Cerity administrator can provide a list of reasons and require you to enter your electronic signature to complete the dialog box.



Basic Exercise #3

Set up a single-level calibrated method for a sequence

This exercise contains a series of tasks to learn how to:

- Create a method template for a sequence to include single-level, single-update calibration and ESTD quantitation
- Set up calibration and quantitation with fixed compound amounts
- Set up a sequence template

A *sequence template* is a sequence table containing the order of calibration standards and samples that you need to run with this method. A sequence template is useful if the order, sample names and characteristics are similar every time you run a sequence with this method.

You can use this method with [“Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration”](#) and [“Basic Exercise #3b Reintegrate and reprocess the results”](#).

For the tasks on the following pages, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Setting Up Methods”](#) on page 71 for setting up methods.



Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 1. Copy a method to create a method template for a sequence

Task 1. Copy a method to create a method template for a sequence


Steps


1 Create a new method template from an existing method.

- Name the method template, *exer3iii*, where *iii* are your initials.
- Use *exer2iii* or *defexer2* as the template for the new method template.
- Make sure that only Compound Identification and Calibration and Quantitation are marked.


You copy a method when you want to keep the instrument and data analysis settings from the old method. You do not have to take the time to enter the settings into the new method.

Detailed Instructions

- Select **File > New > Method** or click  and select **Method**. The **Method Wizard New Method** panel appears.
- Click the **Browse** button and select *exer2iii* or *defexer2* from the **Method Template Selection** dialog box.
- Enter *exer3iii* in the **New Method Name** box.
- Select **Sequence**.





- Click **Next** until you reach the **Data Analysis** panel.
- Mark the **Calibration and Quantitation** check box.




Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 1. Copy a method to create a method template for a sequence

Steps	Detailed Instructions
2 Make selections to keep the compound table and set up new calibration. <ul style="list-style-type: none">Make selections to set up:<ul style="list-style-type: none">single-level calibrationfixed compound amountssingle-update calibrationsequence-specific calibration	<p>a Click Next to scroll to the Compound Table panel.</p> <p>b Select the Keep Compound Calibration from Method template option. This option lets you keep the compound table from the previous method (even though no calibration was set up).</p>  <p>c Click Next to scroll to the Identification panel.</p> <p>d Do not mark any check boxes on the Identification panel.</p> <p>e Click Next to scroll to the Calibration panel.</p> <p>f Select Fixed Amount and use the default options.</p> 

Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 1. Copy a method to create a method template for a sequence

Steps	Detailed Instructions
3 Set up Quantitation and review your new method.	<div><div>a Click Next to scroll to the Quantitation panel.</div><div>b Make sure that the Limit checks check box is clear and the ESTD option is selected.</div></div> <div></div> <div><div>c Click Next to scroll to the New Method Review panel.</div><div>d Review the Method Wizard Settings.</div><div>e Click the Finish button to save your new method.</div></div>

Task 2. Select an example chromatogram


Steps


1 Select an example chromatogram.

- Use the example chromatogram you produced with Basic Exercise 2a or 2b of the “Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration” and “Basic Exercise #3b Reintegrate and reprocess the results”.
- Or, use defexchrom2a.

You do not need to select an example chromatogram. However, it is easier to modify compound identification if you do.


Detailed Instructions

- On the selection tree, expand the exer3iii folder.
- Expand the **Data Analysis** folder.
- Select the **Example Chromatogram** item.
- On the **Tools** toolbar, click .



- Select the sample name with the injection number to produce the example chromatogram.
- Click the **Select** button.

The example chromatogram appears in the workspace.



After you have selected the example chromatogram, you can see the integration and identification settings that belong to the original method.

Task 3. Modify compound identification

Steps	Detailed Instructions
<div>1 Remove a compound from the compound table. Remove the o-terphenyl compound.</div>	<div>a On the selection tree, select Identification under the Data Analysis folder. b Select the o-terphenyl cell. c Right-click the o-terphenyl cell and select Remove Compound.</div>

CompoundOptions

VWD: Absorbance

200
100
0

01234min

Compound Name	Expected Time	Peak Signal	Time Reference Peak	Use Default Time Window
dimethylphthalate	1.0160	VwD1 A	-	+
diethylphthalate	1.1910	VwD1 A	-	+
biphenyl	1.9944	VwD1 A	-	+

Task 4. Set up calibration

Steps

1 Set up calibration for dimethylphthalate.

dimethylphthalate - 10µg

Detailed Instructions

- a** On the selection tree, select **Calibration** under the Data Analysis folder.
- b** On the calibration table, select dimethylphthalate.
- c** On the **Options** tab, enter 10 in the **Weighed Amount** box and µg in the **Amount Unit** box.

Compound Name	Expected Time	Weighed Amount	Amount Unit	Quantitation Based On	RF (Rsp/Amt)
dimethylphthalate	0.9349	10.0000	µg	area	0.0000
diethylphthalate	1.1044	0.0000		area	N/A
biphenyl	1.8794	15.0000	µg	area	0.0000

Options	
Compound Name :	<input type="text" value="dimethylphthalate"/>
Weighed Amount :	<input type="text" value="10"/>
Amount Unit :	<input type="text" value="µg"/>
Comment :	<input type="text"/>

Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 4. Set up calibration

Steps

2 Set up calibration for biphenyl.

Biphenyl - 15µg

Detailed Instructions

a On the calibration table, select biphenyl.

b On the **Options** tab, enter 15 in the **Weighed Amount** box and µg in the **Amount Unit** box.

Compounds

Compound Name	Expected Time	Weighed Amount	Amount Unit	Quantitation Based On	RF (Rsp/Amt)
dimethylphthalate	0.9349	10.0000	µg	area	0.0000
diethylphthalate	1.1044	0.0000		area	N/A
biphenyl	1.8794	15.0000	µg	area	0.0000

Options

Compound Name :

biphenyl

Weighed Amount :

15

Amount Unit :


µg

Comment :

Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 4. Set up calibration

Steps	Detailed Instructions
3 Remove diethylphthalate from the calibration table.	<p>a On the calibration table, right-click anywhere and select Remove Compound from the shortcut menu.</p> <p>The Select Compound(s) dialog box appears.</p> <p>b In the Calibration Table list, select diethylphthalate.</p> <p>c Click the < button to put diethylphthalate in the Available Compounds list.</p> <p>d Click the OK button.</p>



Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 5. Set up quantitation for all four peaks

Task 5. Set up quantitation for all four peaks

Steps	Detailed Instructions
<p>1 Base the quantitation of diethylphthalate on dimethylphthalate.</p> <p>Use an amount multiplier of .8.</p>	<p>a On the selection tree, select Quantitation Setup under the Data Analysis folder.</p> <p>b Click the Uncalibrated Compounds tab.</p> <p>c Under Compound Calibration Type, select the Use Compound option.</p> <p>d Select dimethylphthalate from the Use Compound list.</p> <p>e Enter .8 in the Amount Multiplier (Compound) box.</p>

Calibrated CompoundsUncalibrated CompoundsUnidentified Peaks

Compound Name	Expected Time	Compound Calibration Type	Amount Multiplier (Compound)	RF (Rsp/Amt)	Compound Group
diethylphthalate	1.1044	dimethylphthalat ..	1.0000	N/A	

Compound Name

diethylphthalate

Compound Calibration Type

☒ Use Compound

dimethylphthalate

☐ Manual Response Factor

N/A

☐ No Quantification

Compound Group

None

New...

Compound Info

Amount Multiplier (Compound)

.8

Task 5. Set up quantitation for all four peaks

Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 6. Set up the sequence template

Task 6. Set up the sequence template

Steps

1

Enter the following calibration standards and samples into the sequence template:

Cal1- full-strength isocratic standard

Sample 1_2 - isocratic standard diluted 1/2 with methanol

Sample 1_4 - isocratic standard diluted 1/4 with methanol

NOTE

You cannot set up a sequence template with calibration standards until you have set up calibration in Data Analysis.

Detailed Instructions

a

On the selection tree, select **Sequence Template** for the method.

b

On the sample table, enter the calibration standard for row one.

Enter Cal1 in the **Sample Name** box.

Select **Calibration Standard** from the **Sample Type** list.

Enter the **Vial#** where this standard is located in the ALS.

Click the **Apply** button to put the sample information into the sample table.

c

Enter sample 1_2 for row two.

Select Row 2 in the sample table.

Enter sample 1_2 in the **Sample Name** box.

Select **Sample** from the **Sample Type** list.

Enter the **Vial#** where this sample is located in the ALS.

Click the **Apply** button to put the sample information into the sample table.

d

Enter sample 1_4 for row three.

Select Row 3 in the sample table.

Enter sample 1_4 in the **Sample Name** box.

Select **Sample** from the **Sample Type** list.

Enter the **Vial#** where this sample is located in the ALS.

Click the **Apply** button to put the sample information into the sample table.

	Sample Name	Sample Type	Cal. Level	Summary Group	Vial #	Injections #	Injection Volume [ul]	Samp Amou [mg/n]
1	Cal1	Calibration	1		2	1	as method	0
2	sample 1_2	Sample			5	1	as method	0
3	sample 1_4	Sample			9	1	as method	0

Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 6. Set up the sequence template

Steps

- 2 Enter two more sets of Cal1, sample1_2 and sample 1_4 into the template.**


Hint: Use the Fill Down Wizard and Copy command.

The standards and samples in the final template appear in the following order:

- calibration standard
- two samples,
- calibration standard
- two samples,
- calibration standard
- two samples

Detailed Instructions

- a** Click **Fill Down** on the Edit toolbar, and select **Fill Down Wizard**.
The Fill Down Wizard appears.
- b** Under **Range**, select **Append**, enter 6, and click **Next**.




- c** On the **Sample Names** panel, enter cal1 in the **Name** box, and click **Next**.
- d** On the **Vial Numbers** panel, clear the **Define Vial numbers?** check box, and click **Finish**.
- e** When the **Apply Sample Changes** dialog box appears, click **Yes**.
See that the six new rows display copies of the first row of the template.
- f** Select the two samples on lines 2 and 3, and click the **Copy** button on the Edit toolbar.
- g** Select rows 5 and 6, and click the **Paste** button on the Edit toolbar.
- h** Select rows 8 and 9, and click the **Paste** button on the Edit toolbar.

	Sample Name	Sample Type	Cal. Level	Summary Group	Vial #	Injections #	Injection Volume [μl]	Samp Amou [mg/n]
1	Cal1	Calibration	1		2	1	as method	0
2	sample 1_2	Sample			5	1	as method	0
3	sample 1_4	Sample			9	1	as method	0
4	Cal1	Calibration	1		2	1	as method	0
5	sample 1_2	Sample			5	1	as method	0
6	sample 1_4	Sample			9	1	as method	0
7	Cal1	Calibration	1		2	1	as method	0
8	sample 1_2	Sample			5	1	as method	0
9	sample 1_4	Sample			9	1	as method	0
10								

Basic Exercise #3 Set up a single-level calibrated method for a sequence

Task 6. Set up the sequence template

Steps	Detailed Instructions																																																																																																														
<p>3 Set how the calibration will be updated:</p> <p>First Cal1- Replace (both RF and RT)</p> <p>Second Cal1 - Average for RF and Floating average for RT (Weighted 60% after RT)</p> <p>Third Cal1 - Average for RF and Floating average for RT (Weighted 75% after RT)</p>	<p>a On the sequence table, select the first Cal1.</p> <p>b Click the Run tab.</p> <p>c Under Calibration, select Replace from the Response Factor Update list and select Replace from the Retention Time Update list.</p> <p>d Select the second Cal1 in the sequence table.</p> <p>e Select Average from the Response Factor Update list and Floating average from the Retention Time Update list.</p> <p>f Select 60%.</p> <p>g Repeat steps d and e for the third Cal1.</p>																																																																																																														
	<table> <tr> <th></th><th>Sample Name</th><th>Sample Type</th><th>Cal. Level</th><th>Custom Sample Group</th><th>Vial #</th><th>Injections #</th><th>Injection Volume [µl]</th><th>Sample Amount</th><th>Multipl</th></tr> <tr><td>1</td><td>cal1</td><td>Calibration</td><td>1</td><td></td><td>2</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>2</td><td>sample 1_2</td><td>Sample</td><td></td><td></td><td>5</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>3</td><td>sample 1_4</td><td>Sample</td><td></td><td></td><td>9</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>4</td><td>cal1</td><td>Calibration</td><td>1</td><td></td><td>2</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>5</td><td>sample 1_2</td><td>Sample</td><td></td><td></td><td>5</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>6</td><td>sample 1_4</td><td>Sample</td><td></td><td></td><td>9</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>7</td><td>cal1</td><td>Calibration</td><td>1</td><td></td><td>2</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>8</td><td>sample 1_2</td><td>Sample</td><td></td><td></td><td>5</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>9</td><td>sample 1_4</td><td>Sample</td><td></td><td></td><td>9</td><td>1</td><td>as method</td><td>0</td><td>1</td></tr> <tr><td>10</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> </table> <div> <div> <div>Sample Name: cal1</div> <div>Sample Type: Calibration Standard</div> <div>Custom Sample Group: <div></div> <div>New</div></div> <div> <div>Vial Number</div>2 <div>Injections</div>1 <div>Volume [µl]</div>as method </div> </div> <div> <div>Run Amounts Identification Description</div> <div> <div>Calibration</div> <div>Calibration Mode: Single Update</div> <div>Calibration Level: 1</div> <div>Response Update: Average</div> <div>Retention Time Update: Floating Average 60 %</div> </div> </div> </div>		Sample Name	Sample Type	Cal. Level	Custom Sample Group	Vial #	Injections #	Injection Volume [µl]	Sample Amount	Multipl	1	cal1	Calibration	1		2	1	as method	0	1	2	sample 1_2	Sample			5	1	as method	0	1	3	sample 1_4	Sample			9	1	as method	0	1	4	cal1	Calibration	1		2	1	as method	0	1	5	sample 1_2	Sample			5	1	as method	0	1	6	sample 1_4	Sample			9	1	as method	0	1	7	cal1	Calibration	1		2	1	as method	0	1	8	sample 1_2	Sample			5	1	as method	0	1	9	sample 1_4	Sample			9	1	as method	0	1	10									
	Sample Name	Sample Type	Cal. Level	Custom Sample Group	Vial #	Injections #	Injection Volume [µl]	Sample Amount	Multipl																																																																																																						
1	cal1	Calibration	1		2	1	as method	0	1																																																																																																						
2	sample 1_2	Sample			5	1	as method	0	1																																																																																																						
3	sample 1_4	Sample			9	1	as method	0	1																																																																																																						
4	cal1	Calibration	1		2	1	as method	0	1																																																																																																						
5	sample 1_2	Sample			5	1	as method	0	1																																																																																																						
6	sample 1_4	Sample			9	1	as method	0	1																																																																																																						
7	cal1	Calibration	1		2	1	as method	0	1																																																																																																						
8	sample 1_2	Sample			5	1	as method	0	1																																																																																																						
9	sample 1_4	Sample			9	1	as method	0	1																																																																																																						
10																																																																																																															
<p>4 Save the method.</p> <p>After you complete this method, you may use it to run a sequence. See “Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration” on page 29 and “Basic Exercise #3b Reintegrate and reprocess the results” on page 39.</p>	<p>a Click , and enter your reasons for changes and electronic signature, if necessary.</p>																																																																																																														



Advanced Exercise #4

Set up a multi-level calibrated method for a sequence

This exercise contains a series of tasks to learn how to:

- Use an existing method to create a new method template for a sequence
- Include multi-level, overall calibration and ESTD quantitation in the method
- Set up calibration and quantitation with variable compound amounts for a calibration table with two levels
- Set up system sample variables
- Set up a sequence template for overall calibration
- Select a new report template for a single standard injection report

See [“Basic Exercise #3 Set up a single-level calibrated method for a sequence”](#) on page 91 to learn what a sequence template is.

You can use this method with [“Advanced Exercise #4a Run a sequence to quantify compounds with multi-level calibration”](#) on page 45 and [“Advanced Exercise #4b Change sample variables in the method and reprocess”](#) on page 53.

For the tasks in this exercise, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Setting Up Methods”](#) on page 71 to set up methods.



Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 1. Copy a method to create a new method template for a sequence

Task 1. Copy a method to create a new method template for a sequence


Steps


1 Copy the method to create a new template.

- Copy either *exer3iii* or *defexer3*.
- Name the method template, *exer4iii*, where *iii* are your initials.
- Change nothing until you reach the Compound Table panel.

Note that the Method Wizard panels contain the method selections from Exercise 3.

Detailed Instructions

- Select **File > New > Method** or click  and select **Method**.
The Method Wizard appears.
- On the New Method panel, click the **Browse** button, and select *exer3iii* or *defexer3*.
- Enter *exer4iii* in the **New Method Name** box.




- Click **Next** until you reach the Compound Table panel.

2 Set up the Compound Table panel.

Because you will set up a multi-level calibration, you set up a new calibration table.


- On the Compound Table panel, select **Set up a new Compound Calibration**.



- Click **Next** until you reach the **Calibration** panel.

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 1. Copy a method to create a new method template for a sequence

Steps	Detailed Instructions
3 Set up the calibration panel. Choose to set up: <ul style="list-style-type: none">• multi-level calibration (2 levels)• variable compound amounts• overall calibration• sequence-specific calibration	<p>a Select Variable Amount.</p> <p>b Mark the Multi Level check box, and enter 2 levels.</p> <p>c Select Overall Calibration.</p>  <p>d Click Next until you reach the New Method Review panel.</p>
4 Review your new method template.	<p>a On the New Method Review panel, review the Method Wizard Settings.</p> <p>b Click the Finish button to save your new method.</p> <p>c Save all changes to the database, with a reason if necessary.</p>

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 2. Set up example chromatogram and compound identification

Task 2. Set up example chromatogram and compound identification

Steps

1 Select an example chromatogram.


Use the example chromatogram that you produced with “Basic Exercise #3a Run a sequence to quantify compounds with single-level calibration” and “Basic Exercise #3b Reintegrate and reprocess the results”.

Or, use defexchr2a. (To use this chromatogram, use an instrument with a VWD detector.)


If you cannot see the sample whose chromatogram that you want to select, select another query.

Hint: The result, defexchr2a, is a restored result.

Detailed Instructions


- a On the selection tree, expand the new method template, **exer4iii**.
- b Expand the **Data Analysis** folder and select **Example Chromatogram**.
- c On the **Tools** toolbar, click .

The **Select example chromatogram** dialog box appears.



- d Select the injection from the analysis that contains the example chromatogram for the new method. If you do not see the defexchr2a under the Samples folder, select the query, **AllResultsRestored**.
- e Click the **Select** button.

The example chromatogram appears in the workspace.



The integration parameters are retained from the Exercise 3 method. You do not have to set up integration.

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 2. Set up example chromatogram and compound identification

Steps

2 Set up the compound table for these compounds:

RT=0.9-1.1 min, dimethylphthalate


RT=1.1-1.3 min, diethylphthalate

RT=1.8-2.0 min, biphenyl

Do not identify the fourth peak. In another exercise, you will set up the fourth peak as an unspecified impurity that is not identified based on retention time.

Detailed Instructions

a On the selection tree, select **Identification** under the Data Analysis folder.

b On the Tools toolbar, click .

The peaks appear with the names New Compound one through four in the compound table.

c Under **Compound Name**, select the first cell and enter dimethylphthalate.

After you select the cell, enter the name. The previous entry is overwritten.


d Under **Compound Name**, select the second cell and enter diethylphthalate.

e Under **Compound Name**, select the third cell and enter biphenyl.

f Under **Compound Name**, right-click the fourth cell.

g Select **Remove Compound**.

On the Identification workspace, view the three identified peaks and one unidentified peak.



Task 3. Set up calibration and quantitation

Steps	Detailed Instructions
<div><p>1 Set up calibration for dimethylphthalate and biphenyl.</p><p>Default amounts for dimethylphthalate:</p><ul style="list-style-type: none">• Level 1 - 10µg• Level 2 - 40µg<p>Default amounts for biphenyl:</p><ul style="list-style-type: none">• Level 1 - 15µg• Level 2 - 60µg<p>When you set up a method with variable compound amounts, the application lets you enter the actual weight (concentration) of the standard compounds in sample entry.</p></div>	<div><p>a On the selection tree, select Calibration under the Data Analysis folder.</p><p>b On the Compounds table, select dimethylphthalate.</p><p>c On the Options sheet, click the Use Default Amount cell and select +.</p><p>When you make this selection, the amount that you enter in the Weighed Amount cell for each level appears in the Amounts sheet in Sample Entry.</p><p>d For level 1, enter 10 in the Weighed Amount box and µg in the Amount Unit box.</p><p>e For level 2, enter 40 in the Weighed Amount box.</p><p>f Repeat steps c-e for biphenyl.</p></div>

Compounds

Default Calibration Curve

Compound Name	Level Id	Weighed Amount	Use Default Amount	Amount Unit	Quantitation Based On
dimethylphthalate	1	10.0000	+	ug	area
	2	40.0000			
diethylphthalate	1	0.0000	-		area
	2	0.0000			
biphenyl	1	15.0000	+	ug	area
	2	60.0000			

Options


Calibration Curve

Compound Namebiphenyl

Level Id	Weighed Amount	Use Default Amount	Amount Unit	Low Amount Limit	Use Low L
1	15.0000	+	ug	14.2500	-
2	60.0000			57.0000	

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 3. Set up calibration and quantitation

Steps	Detailed Instructions
<p>2 Remove diethylphthalate from the calibration table.</p> <p>The system has automatically added all compounds from the compound identification table to the calibration table.</p> <p>In this step, remove diethylphthalate to use it as an uncalibrated compound that is quantified based on the response factors of a different compound.</p>	<p>a On the calibration table, right-click anywhere and select Remove Compound from the shortcut menu.</p> <p>The Select Compounds dialog box appears.</p> <p>b In the Calibration Table list, select diethylphthalate.</p> <p>c Click the < button to put diethylphthalate in the Available Compounds list.</p> <p>d Click the OK button.</p> 
<p>3 Set up quantitation as you did in Exercise 3.</p>	<p>See “Task 5. Set up quantitation for all four peaks” on page 100.</p>

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 4. Set up system sample variables

Task 4. Set up system sample variables

Steps

1

Set up a multiplier called “dilution factor”.

Use a default value of 5.

2

Set up a divisor called “correction factor”.

Use a default value of 2.

Detailed Instructions

a

On the selection tree, select **Sample Variables**.

b

Double-click the Dilution cell, and add the word Factor.

c

Enter a default value of 5.

a

Click the Divisor cell once, and enter the name, Correction Factor.

b

Enter a default value of 2.

System Defined Sample Variables (Set by the user in Sample Entry and used in quantification)

	Variable ID	Display Name	Default Value
1	Multiplier_1	Multiplier	1
2	Multiplier_2	Dilution Factor	5
3	Multiplier_3	Purity	1
4	Multiplier_4		1
5	Multiplier_5		1
6	Divider_1	Correction Factor	2
7	Divider_2		1
8	Divider_3		1
9	Divider_4		1

Task 5. Edit the sequence template

Steps

Detailed Instructions

1 Edit the template to look like this:

- two calibration standards (Lev1,2)
- two samples,
- two calibration standards
- two samples,
- two calibration standards

NOTE

You cannot set up or edit a sequence template with calibration standards until you set up calibration in Data Analysis.

2 Set up to quantify the first sample, Sample 1_2, immediately.

When you make this selection, Sample 1_2 will be quantified using the first set of calibration standards. Sample 1_2, along with the other samples, will also be quantified at a later time using the average of all the calibration standards.

Note that the sequence template still contains the information for the method from Exercise 3 but no longer identifies calibration standards.

- On the selection tree, select **Sequence Template**.
- On the sample table, select the calibration standard for row one.
- Select **Calibration Standard** from the **Sample Type** list.
- Move to another row or click the **Apply** button.
- Repeat steps b-d for the next two standards.
- Select the standard in the first row.
- Click the **Insert** button in the toolbar.
- Change the **Sample Name** of the second standard to Cal2.
- Set the **Vial#** to 3 and the **Calibration Level** to 2.
- Click **Apply**.
- Repeat steps g-j for the next two standards.
- Select the last two sample rows, and click the **Delete** button.

a Double-click the cell for Sample 1_2 under the heading, **Immediate Quantitation**.

b Double-click the **Yes** that appears.

	Sample Name	Sample Type	Cal. Level	Immediate Quantitation	Custom Sample Group	Vial #	Injections #
1	cal1	Calibration	1	NO		2	1
2	cal2	Calibration	2	NO		3	1
3	sample 1_2	Sample		YES		5	1
4	sample 1_4	Sample		NO		9	1
5	cal1	Calibration	1	NO		2	1
6	cal2	Calibration	2	NO		3	1
7	sample 1_2	Sample		NO		5	1
8	sample 1_4	Sample		NO		9	1
9	cal1	Calibration	1	NO		2	1
10	cal2	Calibration	2	NO		3	1
11							

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 5. Edit the sequence template

Steps

Detailed Instructions

3 Use the default compound amounts for all standards.

a Click the **Amounts** tab on the Sample Entry panel

b For each calibration standard:

- Select the standard in the sequence table.
- Under Compound amounts, mark the **Use** check boxes for dimethylphthalate and biphenyl.

	Sample Name	Sample Type	Cal. Level	Immediate Quantitation	Custom Sample Group	Vial #	Injections #	Injection Volume [µl]	Sample Amount
1	cal1	Calibration	1	NO		2	1	as method	0
2	cal2	Calibration	2	NO		3	1	as method	0
3	sample 1_2	Sample		YES		5	1	as method	0
4	sample 1_4	Sample		NO		9	1	as method	0
5	cal1	Calibration	1	NO		2	1	as method	0
6	cal2	Calibration	2	NO		3	1	as method	0
7	sample 1_2	Sample		NO		5	1	as method	0
8	sample 1_4	Sample		NO		9	1	as method	0
9	cal1	Calibration	1	NO		2	1	as method	0
10	cal2	Calibration	2	NO		3	1	as method	0
11									

Sample Name:
cal2

Sample Type:
Calibration Standard

Custom Sample Group:

New

Vial Number
3

Injections
1

Volume [µl]
as method

Run

Amounts

Identification

Description

Sample variables

Sample Amount: 0

Sample Amount Unit: mg/ml

Multiplier: 1

Dilution Factor: 5

Purity: 1

Correction Factor: 2

Compound amounts

Use	Name	Amount
<input checked="" type="checkbox"/>	dimethylphthalate [u]	40
<input type="checkbox"/>	diethylphthalate:	0
<input checked="" type="checkbox"/>	biphenyl [ug]:	60

114

Cerity NDS Getting Started Guide


Task 6. Select a new report template for a report

Steps

Detailed Instructions

- 1 Select a report template for a single standard injection report

- a On the selection tree, select **Reporting**.
- b On the Reporting table, select the Standard single injection report type.
- c Click the **Select Template...** button.
The **Select Report Template** dialog box appears.
- d On the **Select Report Template** dialog box, select the template for the Standard Single Injection Detailed report.
- e Click **OK**.



- 2 Select these report types to print:

- Sample single injection
- Standard single injection
- Sequence

- a Double-click the **Print** cell for the Multi-Injection Summary Group report to change **Yes** to **No**.
- b Repeat step a for the Calibration Standards Group report to change **Yes** to **No**.


Print	Report Type	Report Template
Yes	Sample single injection	exer5injdec.html
Yes	Standard single injection	sin_d.html
Yes	Multi-Injection Summary Group	Smp_short.htm
No	Calibration Standards Group	Cal_short.htm
No	QC Sample Group	QC_short.htm
Yes	Sample Group	exer5sqdec.html
No	Custom Sample Groups	Sum_short.htm
Yes	Sequence	Seq_short.htm
No	Customer Report 1	Composite_1.xml
No	Customer Report 2	Composite_2.xml
No	Customer Report 3	Composite_3.xml

- 3 Save the method.

- a On the Standard toolbar, click , and enter your reasons for changes and electronic signature, if required.

Advanced Exercise #4 Set up a multi-level calibrated method for a sequence

Task 6. Select a new report template for a report



Advanced Exercise #5

Set up a method for a sequence to quantify impurities

This exercise contains a series of tasks to learn how to:

- Include custom, noise and system suitability calculations in the method for a sequence
- Include bracketed calibration and ISTD quantitation in the method
- Set up a custom calculation to average the percent impurities of all the samples in the sequence over multiple injections
- Set limits for custom and system suitability calculations
- Set up a sequence template for bracketing, multiple injections and a blank run for a S/N calculation
- Set up the Result View layout to see the custom and system suitability calculations.
- Edit a report template for a sample group report to include the custom and system suitability calculations

You can use this method with [“Advanced Exercise #5a Run a sequence to quantify impurities”](#) on page 61 and [“Advanced Exercise #5b Use a different method to reprocess”](#) on page 67.

In the tasks in this exercise, try the steps on the left without the detailed instructions. If you need more help, follow the detailed instructions on the right.

Before you start

Read [“Setting Up Methods”](#) on page 71 to set up methods.



Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 1. Copy a method to create a new method template for a sequence

Task 1. Copy a method to create a new method template for a sequence


Steps


1 Copy the method to create a new template.

- Copy either *exer4iii* or *defexer4iii*. You can use the original method from Exercise 4 or the modified method from Exercise 4b.
- Name the method template, *exer5iii*, where *iii* are your initials.

Note that the Method Wizard panels contain the method selections in Exercise 4.

Detailed Instructions


- Select **File>New>Method** or click  and select **Method**. The Method Wizard appears.
- Click the **Browse** button, and select *exer4iii* or *defexer4iii*.
- Enter *exer5iii* in the **New Method Name** box.



- Click **Next** until you reach the Data Analysis panel.

2 Include the capability to set up custom calculations and system suitability calculations



- On the Data Analysis panel, mark the **Custom Calculations** check box.
- Mark the **Include Noise Calculations** check box. Note that when you mark the **Include Noise Calculations** check box, the **Include System Suitability** check box appears marked and dimmed.



- Click **Next** to scroll to the Compound Table panel.

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 1. Copy a method to create a new method template for a sequence

Steps	Detailed Instructions
3 Select a Compound Table option. Even though you are changing the mode of calibration to Bracketing, you can keep the calibration setup from Exercise 4.	<p>a On the Compound Table panel, Select Keep Compound Calibration from Method template.</p>  <p>b Click Next until you reach the Calibration panel.</p>
4 Select Calibration options. Select Bracketing and keep all other options the same.	<p>a On the Calibration panel, select Bracketing.</p>  <p>b Click Next to scroll to the Quantitation panel.</p>


Task 1. Copy a method to create a new method template for a sequence

Task 1. Copy a method to create a new method template for a sequence

Method Wizard

?

Quantitation



Do you want to include limit checks on the calculated results ?

☒ Limit checks

Which Calibration Mode do you want to use in your Method ?

☐ ESTD
☒ ISTD

Task 2. Edit quantitation for an internal standard

Steps

1 Set up the ISTD quantitation.

Set biphenyl as the internal standard and use it for the quantitation of dimethylphthalate.

Detailed Instructions

- a Expand the method that you just created, and expand the Data Analysis folder.
- b On the selection tree, select **Quantitation Setup**.
- c Click the Calibrated Compounds tab.
- d On the calibration table, select biphenyl.
- e Under Internal Standard, mark **Set this Compound as the ISTD**.
- f Select dimethylphthalate.
- g Under Internal Standard, mark **Use ISTD compound**.
- h Click the down arrow, and select biphenyl from the list.

Calibrated Compounds		Uncalibrated Compounds		Unidentified Peaks	
Compound Name	Expected Time	Compound Group	ISTD	ISTD Name	Com
dimethylphthalate	0.9349			biphenyl	
biphenyl	1.8902		ISTD		

Compound Name

Internal Standard

☐ Set this Compound as the ISTD

☒ Use ISTD compound

Compound Group

None

New..

Compound Info

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 3. Set up a custom calculation to average the percent impurities of all the samples in a sequence

Task 3. Set up a custom calculation to average the percent impurities of all the samples in a sequence

Steps	Detailed Instructions
1 Set up to calculate the percent of impurity in each single injection. The isocratic standard is a well-defined sample with known compounds. To help you learn how to set up a custom calculation, pretend that the composition of the isocratic standard is the following: Main compound - dimethylphthalate Specified impurity - diethylphthalate ISTD - biphenyl Unspecified impurity - unknown peak You can also point and drag the cell reference to specify the cells in the calculation.	<ol style="list-style-type: none"> In the selection tree, select Custom Calculations under Data Analysis. Click the Single Injection tab, if necessary. Add a column that contains the Amount variable for all compounds/peaks. <ul style="list-style-type: none"> Right-click the table, and select Add Column. In the Existing Column sheet, expand Compounds and select Amount. Click Apply. Add a column for the percent specified impurity calculation. <ul style="list-style-type: none"> Click the Add a New Custom Calculation Column tab. Enter the Variable ID for the specified impurity as anything you want, e.g. PercentSpecifiedImpurity (no spaces). Enter the Display Name, e.g., Percent Specified Impurity. Select the Level as Single Inj. Variables, then click Apply. Add a column for the percent unspecified impurity calculation. <ul style="list-style-type: none"> Enter the Variable ID, the Display Name, and select the Level as Single Inj. Variables, and click OK. Enter the formula for the percent specified impurity calculation into the Single Inj. Variables cell. <ul style="list-style-type: none"> Enter the syntax $=D8 / \text{SUM}(D7 : D13) * 100$, which represents the amount of diethylphthalate divided by the sum of the amounts of all the peaks x 100. You can use the f_x button to find the SUM function, or you can type SUM. Enter the formula for the percent unspecified impurity calculation into the Single Inj. Variables cell. (Use same syntax as for the specified impurity.)

1			New	New
2		Amount	Percent Specified Impurity	Percent Unspecified Impurity
3	-			
4	Single Injection			
5	Single Inj. Variables		9.48	19.07
6	- Identified Compounds			
7	dimethylphthalate	0.9993		
8	diethylphthalate	1.9968		
9	biphenyl	3.0126		
10	- Not Identified Peaks			
11	Unknown 1	4.0158		
12	..	4.9725		
13	Unknown n	6.0583		

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 3. Set up a custom calculation to average the percent impurities of all the samples in a sequence

Steps	Detailed Instructions
2 Set up to calculate the average percent of impurity for all injections of a sample. Do this for both the specified and unspecified impurity.	<p>a On the Custom Calculations workspace, click the Multi-injection tab.</p> <p>b Add a column for percent specified impurity.</p> <ul style="list-style-type: none"> Right-click the table, and select Add Column. On the Existing Column sheet, expand User Defined, and select Percent Specified Impurity. Click Apply. <p>c Add a column for the percent unspecified impurity.</p> <ul style="list-style-type: none"> Select Percent Unspecified Impurity. Click Apply. <p>d Add a column for the average of the percent specified impurity for all injections.</p> <ul style="list-style-type: none"> Click the Add a New Custom Calculation Column tab. Enter the Variable ID as anything you want, e.g., AvgPercentSpecified. Enter the Display Name as a variant of the ID, e.g., Avg Percent Specified. Enter the Level as Multiple Inj. Variables, and click Apply. <p>e Add a column for the average of the percent unspecified impurity for all injections of a sample.</p> <ul style="list-style-type: none"> Enter the Variable ID, Display Name and Level as Multiple Inj. Variables. Click OK. <p>f Enter the formula for the average of the percent specified impurity into the Multiple Inj. Variable cell.</p> <ul style="list-style-type: none"> Enter the syntax <code>=AVERAGE(D6:D8)</code>, which represents the average of the percent impurity calculation for each sample or all injections. You can use the f_x button to access the AVERAGE function, or you can type AVERAGE. <p>g Enter the formula for the average of the percent unspecified impurity.</p>

	A	B	C	D	E	F	G
1						New	New
2				Percent Specified Impurity	Percent Unspecified Impurity	Avg Percent Specified	Avg Percent Unspecified
3	-						
4	Multi-Injection Summary						
5	-	Multiple Inj. Variable				2.00	2.00
6		Single Inj. #1		1.00	0.99		
7		..		2.00	2.02		
8		Single Inj. #n		3.01	2.98		
9	-	dimethylphthalate					
10		Single Inj. #1					
11		..					
12		Single Inj. #n					

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 3. Set up a custom calculation to average the percent impurities of all the samples in a sequence

Steps	Detailed Instructions
3 Set up to calculate the average percent of impurity for all samples. Do this for both the specified and unspecified impurity.	<p>a Click the Sample Group tab in the Custom Calculations workspace.</p> <p>b Add a column for average percent specified impurity.</p> <ul style="list-style-type: none"> Right-click the table, and select Add Column. Expand User Defined, and select Avg Percent Specified. Click Apply. <p>c Add a column for the average percent unspecified impurity.</p> <ul style="list-style-type: none"> On the Existing Column sheet, expand User Defined, and select Avg Percent Unspecified. Click Apply. <p>d Add a column for the average of the percent specified impurity for all samples.</p> <ul style="list-style-type: none"> Click the Add a New Custom Calculation Column tab. Enter the Variable ID as anything you want, e.g., AvgPercentSAIISamples. Enter the Display Name as a variant of the ID, e.g., Avg % S All Samples. Enter the Level as Sample Group Variables, and click Apply. <p>e Add a column for the average of the percent unspecified impurity for all samples, e.g., AvgPercentUAIISamples.</p> <ul style="list-style-type: none"> Enter the Variable ID, Display Name and Level as Sample Group Variables. Click OK. <p>f Enter the formula for the average of the percent specified impurity.</p> <ul style="list-style-type: none"> Enter the syntax =AVERAGE(F6:F8), which represents the average of the percent impurity calculation for all samples. You can use the f_x button to access the AVERAGE function, or you can type AVERAGE. <p>g Enter the formula for the average of the percent unspecified impurity for all samples.</p>

	A	B	C	D	E	F	G
1						New	New
2				Avg Percent Specified	Avg Percent Unspecified	Avg % S All Samples	Avg % U All Samples
3	-						
4			Samples				
5	-		Sample Group Variable:			1.99	=AVERAGE (E6:E8)
6			Sample #1	0.99	1.01		
7			..	2.01	1.98		
8			Sample #n	2.97	3.01		
9	-		dimethylphthalate				
10			Sample #1				
11			..				

Task 4. Set up limits for the custom and system suitability calculations

Steps

Detailed Instructions

1 Set up limits for system suitability calculations

- If tailing factor > than 1.7, then say Not Passed - all samples and only dimethylphthalate
- If USP resolution < than 1.5, then say Not Passed - all samples and all compounds
- If signal to noise is less than 5, then say Not Passed.

- Select **Limits** under Data Analysis.
- Make sure the Single Injection sheet appears.
- Right-click the Limits table, and select **Insert New Limit**.
- Expand the **Peak** folder, and select TailingFactor.
- From the **Condition** list, select >, and for **Value**, enter 1.7.
- From the **Apply to** list, select dimethylphthalate, and click **OK**.
- Repeat steps c and d for Peak resolution USP.
- From the **Condition** list, select <, and for **Value**, enter 1.5.
- Click **OK**.
- Repeat steps c and d for SignalToNoise.
- From the **Condition** list, select <, and for **Value**, enter 5.
- Click **OK**.

Limit Options for:

Single Injection				
Variable ID	Header	Units	Condition	Value
SignalToNoise	SignalToNoise		<	5
TailingFactor	TailingFactor		>	1.7
USP_Resolution	Peak resolution USP		<	1.5

2 Set up limits for both the average of the specified impurity and the average of the unspecified impurity for all samples.

- If specified impurity > 10%, not passed
- If unspecified impurity > 5%, not passed

Hint: The tab Summary Groups let you set limits for all the variables and calculations associated with sample-type groups, such as sample group, calibration standard group, custom sample group and QC group.

- Click the **Summary Groups** tab.
- Right-click the table, and select **Insert New Limit**.
- In the Insert New Limit dialog box, expand the **Single Values** folder and select Avg % S All Samples.
- From the **Data Set** list, select Sample.
- From the **Condition** list, select >.
- Enter a value of 10, and click **OK**.
- Repeat steps b-f for the Avg % U All Samples and a value of 5.

Limit Options for:

Summary Groups				
Variable ID	Header	Units	Data Set	Apply To
AvgPercentKAllSamples	Avg % K All Samples		All	Selected Variable ID
AvgPercentUAllSamples	Avg % U All Samples		All	Selected Variable ID

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 5. Edit the sequence template for bracketing and multiple injections

Task 5. Edit the sequence template for bracketing and multiple injections

Steps

Detailed Instructions

1 Set up the brackets

- Quantify the first set of samples with the average RFs of the first and second sets of standards.
- Quantify the second set of samples with the average RFs of the second and third set of standards.

a Select **Sequence Template** in the selection tree.

b Double-click the **Bracketing** cell for Cal1 in row 1, and double-click **Open**.

c Double-click the **Bracketing** cell for Cal1 in row 5, and double-click **Open**.

d Double-click the **Bracketing** cell for Cal2 in row 6, and double-click **Close**.

e Double-click the **Bracketing** cell for Cal2 in row 10, and double-click **Close**.

2 Enter a blank sample in the first row and enter two injections for each sample.

a Select row 1, and click the **Insert** button. (Use tooltip.)




b Enter NoiseBlank for the **Sample Name**, and select Blank Run for the **Sample Type**.

c Enter a different Vial#, and click **Apply**.

d Enter 2 for the Injections # for each sample in the sequence.


	Sample Name	Sample Type	Cal. Level	Bracketing	Custom Sample Group	Vial #	Injections #	Injection Volume [µl]	Sample Amount
1	NoiseBlank	Blank Run				4	1	as method	0
2	cal1	Calibration	1	Open		2	1	as method	0
3	cal2	Calibration	2	None		3	1	as method	0
4	sample 1_2	Sample				5	2	as method	0
5	sample 1_4	Sample				9	2	as method	0
6	cal1	Calibration	1	Open		2	1	as method	0
7	cal2	Calibration	2	Close		3	1	as method	0
8	sample 1_2	Sample				5	2	as method	0
9	sample 1_4	Sample				9	2	as method	0
10	cal1	Calibration	1	None		2	1	as method	0
11	cal2	Calibration	2	Close		3	1	as method	0

Task 6. Set up the Result View layout to see custom and system suitability calculations

Steps	Detailed Instructions
1 Set up to view the percent specified impurity and the percent unspecified impurity.	<p>a On the selection tree, expand the Data Review Layout folder.</p> <p>b Select Single Injection in the selection tree.</p> <p>c Select Summary Table in the workspace.</p> <p>d Select Percent Specified Impurity from the Available Items list, and click > to move it to the Display Items list.</p> <p>e Repeat step d for Percent Unspecified Impurity, and click Apply.</p>
	
2 Set up to view the tailing factor, USP resolution and the S/N for each compound.	<p>a Select the Results Table.</p> <p>b Select Tailing Factor from the Available Items list, and click > to move it to the Display Items list.</p> <p>c Repeat step b for Peak resolution USP and SignalToNoise, and click Apply.</p>
	
3 Set up to view the average of the specified impurity and the average of the unspecified impurity for each sample.	<p>a In the selection tree, select Multiple Injection.</p> <p>b Select the Summary Table in the workspace.</p> <p>c Select Avg Percent Specified from the Available Items list, and click > to move it to the Display Items list.</p> <p>d Repeat step b for Avg Percent Unspecified, and click Apply.</p>
	


Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 6. Set up the Result View layout to see custom and system suitability calculations

Steps	Detailed Instructions
4 Set up to view the average of the percent specified and unspecified impurities in all the samples and their limit checks.	<div><div>a Select Samples in the selection tree.</div><div>b Select Summary Table in the workspace.</div><div>c Select Avg % S All Samples from the Available Items list, and click > to move it to the Display Items list.</div><div>d Repeat step c for Avg % U All Samples, Avg % S All Samples Limit Check and Avg % U All Samples Limit Check.</div><div>e Click Apply.</div></div> <div></div>

Task 7. Edit a report template for the sample group

Steps	Detailed Instructions
<p>1 Edit a report template for a sample single injection report.</p> <ul style="list-style-type: none">Edit the inj.html report.Add a column for USP resolution and Signal to Noise to the existing compounds table under the chromatogram.Save the template as exer5inj<i>iii</i>, where <i>iii</i> is your initials.	<p>a On the selection tree, select Reporting.</p> <p>b Select the Sample single injection report type, and click Edit Template....</p> <p>c Double-click Individual Report Templates, and double-click inj.html.</p> <p>d Place the cursor in the last column of the compounds table located beneath the chromatogram.</p> <p>e Right-click the table, and select Table Properties. The Compound Table Properties dialog box appears.</p> <p>f In the Select Column Fields list, mark the Peak resolution USP and SignalToNoise check boxes, and click OK.</p>



The compound table in the resulting template looks like this:

Retention Time	Compound Name	Amount	Response Factor	Tailing Factor	Peak resolution USP	SignalToNoise
#####.##	×	####.##	×.DDDD	#####.###	##.###	##.###

g Select **File > Save As**, enter exer5inj*iii*, and click **OK**.

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 7. Edit a report template for the sample group

Steps	Detailed Instructions
2 Edit the sample group detailed report template (sus_d.html). <ul style="list-style-type: none">Insert an html table under the Sample group variables table.Enter the text for the Avg. % S Impurity All Samples and Avg% U Impurity All Samples.Enter the placeholder for the values for the % impurities.Under the Sample Group Limits table, enter the Limit check information for the sample group.Save the template as exer5sg<i>iii</i>, where <i>iii</i> are your initials.	<ul style="list-style-type: none">a Exit the Report Template Editor.b Select the Sample Group report type, and click Edit Template....c Double-click Individual Report Templates, and double-click sus_d.html.d Insert a line below the Sample group variables table, and click the Insert HTML table button.e In the Insert Table dialog box, select the Classic Table Style and click OK.f Click the Fields tab and expand the Sample Group folder.g Expand the Sample Group Variables Results folder.h Place the cursor into the first cell of the HMTL table, press the Alt key and double-click Avg % S All Samples.i Place the cursor into the second cell in the first row and double-click Avg % S All Samples.j Repeat steps h and i for Avg % U All Samples, using the second row.k Place the cursor below the Sample group limit results table.l Press the Ctrl key and double-click Avg % S All Samples Limit Check.m Do the same for Avg % U All Samples Limit Check.n Select File > Save As, enter exer5sg<i>iii</i> and click Save.

After you finish, the template displays as the Sample group template



Sample group (detailed)

Sequence name:	XXXXXXXXXXXXXXXXXXXX
Sequence Start:	sys_Date, sys_Time
Sequence End:	sys_Date, sys_Time
Method (rev):	XXXXXXXXXXXXXXXXXXXX (###)

Number of unidentified peaks: ##

Sample group variables				
#	Sample name	Amount	Position	Inj. vol.
##	XXXXXXXXXXXXXXXXXXXX	## DDDD	XXXXXXXX	### DD

Avg % S All Samples:	## DD
Avg % U All Samples:	## DD

Sample group limit results				
#	Sample name	Compound	Limit (Compound)	Limit (Sample)
##	XXXXXXXXXXXXXXXXXXXX	XXXXXXXXXXXXXXXXXXXX	XXXXXXXX	XXXXXXXX

Avg % S All Samples Limit Check: XXXXXXXX

Avg % U All Samples Limit Check: XXXXXXXX

Task 8. Select report templates and report types

Steps

Detailed Instructions

1

Select report templates for report types.

• Use exer5inj*iii* for the Sample single injection report.

• Use exer5sg*iii* for the Sample group report.

a

Exit the Cerity Report Template Editor.

b

Select the Sample single injection report type and click **Select Template....**

c

Select exer5inj*iii* and click **OK**.

d

Select the Sample group report type and click **Select Template....**

e

Select exer5sg*iii* and click **OK**.

2

Select these report types to print.

• Sample single injection

• Standard single injection

• Multi-injection summary

• Sample group

• Sequence

a

Double-click the **Print** cell for the Multi-Injection Summary Group report to change **No** to **Yes**.

b

Repeat instruction (a) for the Sample Group report to change **Yes** to **No**.


Print	Report Type	Report Template
Yes	Sample single injection	exer5injdec.html
Yes	Standard single injection	sin_d.html
Yes	Multi-Injection Summary Group	Smp_short.htm
No	Calibration Standards Group	Cal_short.htm
No	QC Sample Group	QC_short.htm
Yes	Sample Group	exer5sgdec.html
No	Custom Sample Groups	Sum_short.htm
Yes	Sequence	Seq_short.htm
No	Customer Report 1	Composite_1.xml
No	Customer Report 2	Composite_2.xml
No	Customer Report 3	Composite_3.xml

Select Template...

Edit Template...

3

Save the method

On the Standard toolbar, click  , and enter your reasons for changes and electronic signature, if necessary

Advanced Exercise #5 Set up a method for a sequence to quantify impurities

Task 8. Select report templates and report types

In This Book

This Getting Started Guide is a collection of basic and advanced exercises that provide a quick way to learn the Cerity Pharmaceutical QA/QC application.

The exercises are grouped into two groups:

Running Routine Samples
exercises help lab technicians learn how to run routine samples.

Setting up methods
exercises help Chemists learn how to set up methods for the laboratory.

© Agilent Technologies Deutschland GmbH 2003

Printed in Germany
5/2003



G4000-90011



Agilent Technologies