

# INFORMATICS GENOMICS PROTEOMICS METABOLOMICS


0 1 1 0 1 0 0 0 1 0 1 1 1 0 0 0 1 1 0 0 0 0 1 1 0 1 1 0 1 0 0 0 1  
0 0 1 0 1 1 1 0 0 1 1 0 1 0 0 0 0 1 0 1 1 0 1 0 1 1 1 0 1 1 0 0 1 1

## Agilent ChIP Analytics 1.2 Software

ChIP Analytics 1.2 software is an easy-to-use application for analyzing data from Agilent ChIP-on-chip microarrays. It accepts data output from Agilent Feature Extraction software or AXON GenePix® software and analyzes the significance of protein binding events without extensive processing and manipulation. The software takes advantage of the high-speed statistical analysis functions that include user-configurable heuristics for binding event identification based on p-values and adjacent probes, as well as inter- and intra-array intensity normalization and error modeling.

### Features and Benefits at a Glance

- New visualization capabilities with reference genome, QC plots, and binding event data
- Peak detection algorithm for convenient visualization of binding sites
- Simultaneous analysis of sets of microarrays and experiments
- Output reports containing probe, sequence, and gene detail
- UCSC track report option allows visualization of binding event data in the UCSC Genome Browser
- Support for multiple biological and technical replicates




**Figure 1.** The peak detection visualization feature allows inspection of the regions of interest and a display of nearby genomic features.



Agilent Technologies

## Simultaneous Analysis of Sets of Microarrays with Annotations

Overview of ChIP Analytics workflow and inputs and outputs. Software supports use of multiple biological and technical replicates.




## Output Reports Containing Probe, Sequence, and Gene Detail

Output allows identification of probes, genes, and genomic loci that exhibit significant binding. ChIP Analytics reports its findings in file formats that can be easily loaded into other analysis software.

Microsoft Excel 050310_E2F4_HES_HP01a_251314910018_1_013149_20050507.Probe.tsv																			
File Edit View Insert Format Tools Data Window Help																			
Type a question for help																			
A1 Name																			
Name	ID	Block	Row	Column	Chromosome	Start	Local Control	Type	Foreign IP	Background	WCE Fore	WCE Back	IP Background	WCE Norm	IP Norm	WCE Norm	Log2	P-value	
chr1:14866.A_76_P02	1	267	78	chr1	1.47E+08	data	19145	103	1081	73	16	11	20591.85	951.6588	4.435485	2.75919	9.22E-14		
chr1:14866.A_76_P01	1	286	1	chr1	1.47E+08	data	12431	91	1120	66	31	15	13339.04	997.3867	3.741473	2.668752	5.10E-13		
chr1:44874.A_76_P01	1	245	33	chr1	4.874E12	data	4827	111	596	80	40	40	5088.443	453.4251	3.456816	2.616833	3.38E-12		
chr1:21921.A_76_P02	1	352	48	chr1	2.19E+08	data	2101	100	345	73	15	11	2150.305	221.2929	3.280512	2.580355	2.73E-12		
chr1:15335.A_76_P02	1	234	16	chr1	1.53E+08	data	9465	96	1240	70	22	11	10123.86	1112.419	3.189988	2.559241	1.05E-12		
chr1:14785.A_76_P02	1	388	92	chr1	1.48E+08	data	4635	91	701	64	13	9	4902.307	583.4888	3.070659	2.53029	6.91E-12		
LACC-Hwi	LACC-Hwi	1	324	78	NONE	NONE	common	8178	98	1298	70	17	10	8728.92	1169.975	2.899326	2.482591	1.65E-11	
chr1:54371.A_76_P02	1	303	15	chr1	5.437791E	data	3942	91	725	67	14	8	4152.351	604.3381	2.7805	2.445628	3.19E-11		
chr1:19904.A_76_P01	1	207	69	chr1	1.99E+08	data	3850	98	723	73	31	38	4045.215	596.3994	2.761866	2.438167	3.64E-11		
LACC-Hwi	LACC-Hwi	1	3	34	NONE	NONE	common	7033	92	1274	66	13	10	7496.309	1150.120	2.704386	2.420385	4.98E-11	
chr1:14662.A_76_P02	1	59	87	chr1	1.47E+08	data	23485	87	4048	64	14	9	25305.86	3904.878	2.696122	2.417592	5.23E-11		
chr1:91671.A_76_P01	1	103	31	chr1	9.167852E	data	6891	96	1281	70	21	12	7338.31	1153.105	2.669924	2.408357	6.15E-11		
chr1:22405.A_76_P02	1	282	46	chr1	2.24E+08	data	2879	100	621	73	16	15	2992.247	495.1802	2.595204	2.380916	9.92E-11		
chr1:11731.A_76_P02	1	427	11	chr1	1.17E+08	data	5289	94	1054	68	58	34	5506.811	929.8273	2.592145	2.378779	1.03E-10		
chr1:16002.A_76_P01	1	111	69	chr1	1.6E+08	data	5475	92	1129	70	21	26	5810.262	1002.269	2.535334	2.358415	1.46E-10		
chr1:93255.A_76_P02	1	310	24	chr1	9.325679E	data	1877	103	471	78	29	24	1904.649	341.3667	2.480131	2.332846	2.26E-10		
chr1:21565.A_76_P02	1	349	22	chr1	2.156501E	data	25679	101	5329	71	67	57	27665.03	5169.125	2.420071	2.311871	3.22E-10		
chr1:42551.A_76_P02	1	182	10	chr1	4.255782E	data	2091	98	527	67	14	9	2152.47	407.8538	2.398869	2.302858	3.75E-10		
chr1:19384.A_76_P02	1	261	45	chr1	1.94E+08	data	7710	107	1711	77	32	45	8212.717	1572.867	2.384463	2.2964	4.18E-10		
chr1:91671.A_76_P01	1	279	55	chr1	9.167827E	data	3946	100	961	73	15	11	4146.94	832.5775	2.316391	2.266384	6.89E-10		
chr1:35994.A_76_P02	1	209	34	chr1	3.599477E	data	10013	95	2292	68	20	18	10717.58	2158.35	2.312032	2.264535	7.10E-10		
chr1:14932.A_76_P01	1	249	90	chr1	1.493203E	data	1878	87	513	64	14	9	1923.646	386.538	2.276408	2.247351	9.42E-10		
chr1:40666.A_76_P02	1	322	38	chr1	4.066637E	data	2297	98	619	74	29	35	2364.578	492.2031	2.264257	2.238867	1.08E-09		
chr1:20855.A_76_P02	1	88	37	chr1	2.08E+08	data	2439	91	664	66	14	11	2525.824	544.7974	2.212963	2.217171	1.54E-09		
chr1:24055.A_76_P01	1	51	27	chr1	2.405224E	data	1748	93	520	71	19	22	1775.869	396.530	2.161539	2.189442	2.41E-09		
LACC-Hwi	LACC-Hwi	1	393	66	NONE	NONE	common	1991	94	590	67	13	9	2037.758	470.3715	2.11511	2.17529	3.42E-09	
chr1:88861.A_76_P02	1	348	91	chr1	8.886165E	data	1356	93	431	67	12	9	1351.652	312.5887	2.112386	2.16526	3.54E-09		
chr1:20822.A_76_P02	1	202	80	chr1	2.08E+08	data	5116	92	1445	72	13	9	5421.756	1313.865	2.049444	2.128866	6.16E-09		
chr1:84175.A_76_P02	1	13	19	chr1	8.417582E	data	1935	86	601	67	14	10	1985.813	481.2873	2.04476	2.228144	6.23E-09		
chr1:32855.A_76_P01	1	373	60	chr1	3.285279E	data	7503	95	2107	69	16	14	8000.608	1973.774	2.019153	2.115459	7.71E-09		
chr1:17055.A_76_P01	1	93	6	chr1	1.71E+08	data	1043	89	373	69	14	11	1017.256	253.048	2.0072	2.105955	8.93E-09		
chr1:3419.A_76_P01	1	80	12	chr1	3.419E10	data	1456	89	486	64	13	10	1463.117	364.181	2.007354	2.104838	9.09E-09		


Option of UCSC Track Report

Compatibility with the UCSC Genome Browser view allows you to easily visualize binding data (shown here is a May 2004 human reference sequence based on NCBI Build 35).



QC Report

Shown here is an example of the convenient, at-a-glance visualization for key experimental elements. The QC report shows a plot of the immunoprecipitated extract (IP) vs. the whole cell extract (WCE) and allows identification of enriched vs. unenriched targets.



**Specifications**

Product number (commercial use)	<b>G4477AA, option CML</b>
Product number (academic/non-profit use)	<b>G4477AA, option NFP</b>
License type and duration	<b>One year, single workstation</b>
Operating system	<b>Windows® 2000 Windows® XP SP2 Mac OS X</b>
RAM	<b>Minimum 512 MB (1 GB or more preferred)</b>
Hard disk space	<b>200 MB (software) 2 GB or more (design files and Feature Extraction data files)</b>
Additional software	<b>Internet Explorer 6.0 or later Adobe® Reader 4.0 or later Java 1.4.2_08 or better</b>
CPU	<b>2.0 GHz processor or faster</b>

**About Agilent's****Integrated Biology Solutions**

Agilent Technologies is a leading supplier of life science research systems that enable scientists to understand complex biological processes, determine disease mechanisms, and speed drug discovery. Engineered for sensitivity, reproducibility, and workflow productivity, Agilent's integrated biology solutions include instrumentation, microfluidics, software, microarrays, consumables, and services for genomics, proteomics, and metabolomics applications.

**Agilent online resources  
at your fingertips**

- Ask the Experts — Learn how to get the most out of your Agilent products.
- Special offers — Benefit from savings on Agilent products and services.
- Events calendar — See Agilent at leading industry conferences and trade shows.
- Web links — Go directly to detailed information on the Web.

**For more information**

Learn more:

[www.agilent.com/chem/informatics](http://www.agilent.com/chem/informatics)

Find an Agilent customer center in your country:

[www.agilent.com/chem/contactus](http://www.agilent.com/chem/contactus)

**U.S. and Canada**

1-800-227-9770

[agilent\\_inquiries@agilent.com](mailto:agilent_inquiries@agilent.com)

**Europe**

[info\\_agilent@agilent.com](mailto:info_agilent@agilent.com)

**Asia Pacific**

[adinquiry\\_aplsca@agilent.com](mailto:adinquiry_aplsca@agilent.com)

Research use only. Information, descriptions, and specifications in this publication are subject to change without notice.

Agilent Technologies shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this material.

GenePix® is a registered trademark of Molecular Devices Corporation.

Adobe® is a registered trademark of Adobe Systems Incorporated.

Windows® is a U.S. registered trademark of Microsoft Corporation.

© Agilent Technologies, Inc. 2006

Printed in the U.S.A. April 26, 2006

5989-4360EN

**Agilent Technologies**