


Comparison of Denaturing and Nondenaturing Mobile Phases in SEC: ZORBAX GF-450, +/- GuCl

Application
Biochemical
Robert Ricker

Size-exclusion chromatography is a powerful technique for the size-dependent separation of biomolecules. While standard, non-denaturing mobile phases are often desirable, denaturing mobile phases containing guanidine hydrochloride (or SDS), may be used to insure un-aggregated, monomers for proper size determination. As shown below, denaturation shifts molecules toward larger apparent size and smaller retention volume.



Conditions:
LC: Hewlett Packard HP1050
Column: ZORBAX GF-450 (9.4 x 250), Agilent Part No. 884973-902
UV: 254 nm
Flow: 2.0 mL / min.; ambient
Inj. Vol.: 20 μ L (1 μ g / μ L)

Highlights

- Agilent ZORBAX GF-450 columns separate biomolecules in a size-dependent manner within a molecular-weight range of 10,000-1,000,000 daltons. Note, some molecules (depending upon their characteristics) are retained or excluded by non-SEC mechanisms and are shifted from linearity.
- ZORBAX GF-450 columns are manufactured using extremely hard particles, allowing high flow rates and fast separations.
- When smaller molecules are used in combination with denaturing mobile phases, ZORBAX GF-250 may be used to obtain more-linear separation in the lower molecular-weight range.



Agilent Technologies

*Robert Ricker is an application chemist
based at Agilent Technologies, Wilmington,
Delaware.*

For more information on our products and
services, visit our website at:
www.agilent.com/chem

Copyright© 2002 Agilent Technologies, Inc.
All Rights Reserved. Reproduction,
adaptation or translation without prior
written permission is prohibited, except as
allowed under the copyright laws.

Agilent shall not be liable for errors
contained herein or for incidental or
consequential damages in connection with
the furnishing, performance, or use of this
material.

Information, descriptions, and specifications
in this publication are subject to change
without notice.

Printed in the USA
April 25, 2002
5988-6332EN



Agilent Technologies