

Configure and Maintain Your Agilent ChemStation Computer



For Agilent ChemStation Rev.
B.04.03



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2004, 2005-2010

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G2170-90136

Edition

09/10

Printed in Germany

Agilent Technologies
Hewlett-Packard-Strasse 8
76337 Waldbronn

Software Revision

This guide is valid for revision B.04.03 or higher of the Agilent ChemStation software.

Microsoft® is a U.S. registered trademark of Microsoft Corporation.

Windows 7, Windows Vista, SQL Server, and Visual Studio are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Oracle® is a U.S. registered trademark of Oracle Corporation.

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

In This Guide...

Agilent-delivered ChemStation PC Bundle systems are delivered with the supported Windows operating system pre-installed and are configured for optimum ChemStation performance. Non-Agilent PC Bundle systems require some manual configuration changes in order to provide optimum performance. This manual describes how to configure a non-Agilent PC Bundle system.

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

This chapter identifies and describes all parameters necessary to configure a non-Agilent Windows XP SP 3 system for use with Agilent ChemStation Rev. B.04.03 or higher. It also covers important settings like Internet Security.

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

This chapter identifies and describes all parameters necessary to configure a non-Agilent Windows Vista SP1 system for use with Agilent ChemStation Rev. B.04.03 or higher.

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

This chapter identifies and describes all parameters necessary to configure a non-Agilent Windows 7 system for use with Agilent ChemStation Rev. B.04.03 or higher.

4 Security Settings

This chapter summarizes all security settings for LC/MSD, CE/MSD, LC/MSD Trap, GC 7890A, IO Library Suite 15.5, and Agilent Bootp Service.

5 Additional Modifications

These modifications can only be performed using an Administrator user account

6 Maintaining the PC's File System

To help to optimize the performance of your ChemStation system it is important that you schedule routine maintenance (e.g. deleting unnecessary files). The following instructions should be used as a part of your routine maintenance schedule.

Contents

- 1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher 7**
 - Install Windows XP 8
 - Configure Windows XP Settings 9


- 2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher 25**
 - Install Windows Vista 26
 - Configure Windows Vista Settings 27

- 3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher 45**
 - Install Windows 7 46
 - Configure Windows 7 Settings 47

- 4 Security Settings 67**
 - Security Settings for LC-MSD and CE-MSD Systems 68
 - Security Settings for LC-MSD Trap Systems 72
 - Security Settings for GC 7890A Systems 77
 - Security Settings for IO Library Suite 15.5 80
 - Security Settings for Systems using Agilent Bootp Service 84

- 5 Additional Modifications 85**
 - Install ChemStation 86
 - Profiles and Group Policies to run ChemStation as PowerUser/User 87

- 6 Maintaining the PC's File System 89**
 - Scanning and Repairing NTFS Volumes 90
 - Defragmenting an NTFS Volume 92
 - Scanning for Viruses 94
 - Cleaning up Left-Over Temporary Files 95
 - Performing Backups 97



1

How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Install Windows XP 8

Configure Windows XP Settings 9

This chapter identifies and describes all parameters necessary to configure a non-Agilent Windows XP SP 3 system for use with Agilent ChemStation Rev. B.04.03 or higher. It also covers important settings like Internet Security.



Install Windows XP

NOTE

Keep your PC disconnected from the internet until you have installed the appropriate security patches and hot fixes. It is strongly recommended to install the latest security fixes as supported from Agilent and virus definitions prior to connecting to a network.

- 1 Install Windows XP from recovery CD or Windows XP Professional CD-ROM. During the setup of Windows XP provide the Computer Name, Administrator password and network settings. Select to either join an existing domain or set up the system in workgroup mode.
- 2 Check the Windows version by typing **winver** at the **Start > Run** command prompt: If necessary, install Service Pack 3 for XP. Ignore the option to archive current files and reboot the system.
- 3 To secure your system against viruses please install an AntiVirus program. Running AntiVirus programs might influence the behavior and performance of your computer. Agilent ChemStation was tested to run with Symantecs Norton Antivirus. See “[Scanning for Viruses](#)” on page 94 for more details.

Configure Windows XP Settings

Some changes within this document are mandatory for the ChemStation to work properly on a Windows XP system. Some changes will optimize application performance. Other changes will have a graphical or minor impact.

NOTE

To indicate the relative importance of individual settings, each item is categorized as:

MUST: These changes must be applied.

PERFORMANCE: These changes will improve system performance.

OPTIONAL: Most of these changes will affect the graphical display of the application.

- [OPTIONAL]** 1 General Layout: (**Start Menu**, right-click > **Properties**)
- a **Start Menu Tab:** Select **Classic Start Menu**

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows XP Settings

- b Start Menu > **Customize** button (See [Figure 1](#) on page 10): In **Advanced Start menu options**: Disable the following option:
- **Use Personalized Menus**




Figure 1 Customize Classic Start Menu Dialog

- c Enable in the following options:
- **Display Administrative Tools**
 - **Display Favorites**
 - **Display Log Off** (for Workgroup Systems, this option is only available after configuration step 6)
 - **Display Run**
 - **Enable dragging and dropping**

- Expand Control Panel
 - Expand Network Connections
 - Expand Printers
 - Scroll Programs
 - Show Small Icons in Start menu
- [MUST] 2 Folder Options: (right-click **Start > Explore**)
- a Set **View** to **Details**
 - b Choose **Tools > Folder Options**
 - c **General** tab: for Tasks select **Use Windows classic folders**
 - d **View** tab: Disable the following options:
 - Use simple file sharing
 - Hide extensions for known file types
 - e Enable the following option:
 - Display the full path in the title bar
 - f Apply these settings to all folders by selecting the **Apply to All Folders** button
- [MUST] 3 Internet Explorer: (**Start > Programs > Internet Explorer**) (See [Figure 2](#) on page 11)
- a **Tools > Windows Update**: Apply all critical security patches.




Figure 2 Downloading Critical Security Updates

- [MUST] 4 Register Windows XP Professional with Microsoft.
- [MUST] 5 Windows Logon Options: (For Workgroup Systems, only) (**Start > Settings > Control Panel > User Accounts**)
- a Select **Change the way users log on or off**

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher
Configure Windows XP Settings

- b Disable Use the Welcome Screen
 - c Disable Use Fast User Switching
- [MUST] 6 Power Options: (Start > Settings > Control Panel > Power Options) (See Figure 3 on page 12)
- a Select the Power Schemes tab
 - b From Power Schemes drop down box select one of the power schemes
 - c Set the drop down boxes Turn off hard disks and System standby to Never.
 - d Select the Hibernate tab: disable Enable hibernation .




Figure 3 Power Options Properties Dialog, Power Schemes Tab

- [MUST] 7 Local Security Policy: (Start > Programs > Administrative Tools > Local Security Policy)
- In the Security Settings tree open Local Policies folder, then Security Options folder
 - Double-click the following policy listed in the right hand panel: **Network Access: Sharing and security model for local accounts**
 - In the displayed dialog select the following item from the drop down list: **Classic - local users authenticate as themselves**




Figure 4 Local Security Settings Dialog

- [MUST] 8 Disable Automatic Startup of Windows Messenger: (Start > Programs > Windows Messenger > Tools > Options):
- Select the Preferences tab (See Figure 5 on page 14)
 - In the General section disable Run Windows Messenger when Windows starts.

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows XP Settings

- c. Select **OK** and reboot the system.




Figure 5 Options Dialog, Preferences Tab

- [MUST] 9 Regional and Language Options: (**Start > Settings > Control Panel > Regional and Language Options**) (See [Figure 6](#) on page 15)

- a Set **Regional Options** to **English (United States)** from the drop down list.




Figure 6 Regional and Language Options Dialog, Regional Options Tab

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows XP Settings

- b If regional format other than **English (United States)** is used, the following settings are mandatory. The settings can be defined by clicking the **Customize** button (See [Figure 7](#) on page 16):
- Decimal symbol = . (point)
 - Digit grouping symbol = , (comma)
 - List separator = , (comma)




Figure 7 Customize Regional Options Dialog, Numbers Tab

- c Language for non-Unicode programs: On the **Advanced** Tab select **English (United States)** from the drop down list (See [Figure 8](#) on page 17).




Figure 8 Regional and Language Options Dialog, Advanced Tab

- [PERFORMANCE] 10 Hard disk Properties: (right-click **Start** > **Explore**, right-click on **Local Disk (C:) > Properties**).
- In the **General** tab, disable **Allow Indexing Service to index this disk for fast file searching** (See [Figure 9](#) on page 18)
 - Apply the changes by selecting **OK**.

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows XP Settings

- c In the selection screen, enable the option **Apply changes to C:\, subfolders, and files.**
- d Disable this setting for all non-removable drives.

NOTE

Microsoft applies special sets of permissions to certain folders, such as the \$NTUninstallKB... folders which relate to Windows HotFixes and Patches. These folder permissions result in error messages when applying new settings for the Indexing Service. Select **Ignore**, when these messages appear.




Figure 9 Local Disk (C:) Properties Dialog, General Tab

[PERFORMANCE] 11 Screen layout: (right-click **desktop** > **Properties**):

- a Please change the following items in the **Display Properties** dialog (See [Figure 9](#) on page 18):
- **Themes** tab: Change **Windows XP** to **Windows Classic**
 - **Desktop** tab **Customize Desktop** button:
 - **General** tab: Enable all 4 desktop icons (My Documents, My Computer, My Network Places, Internet Explorer)
 - **General** tab: Disable **Run Desktop Cleanup Wizard every 60 days**
 - **Screen Saver** tab: Change Screen Saver from **Windows XP** to **none**
 - **Appearance** tab > **Effects** button: disable all items, except for
 - **Hide underlined letters for keyboard navigation...**
 - **Settings** tab: 1280 by 1024 resolution, 32 bit color, 75 Hertz refresh or higher. (Please note: some old monitors may only work with 60 Hz.)




Figure 10 Tabs in Display Properties Dialog

[PERFORMANCE] 12 System Properties: **Start** > **Settings** > **Control Panel** > **System**

- a **Advanced** tab > **Performance** > **Settings** button (See [Figure 11](#) on page 20):
- **Visual Effects** tab: Set **Visual Effects** to **Adjust for best performance**.
 - **Advanced** tab > **Virtual Memory**: For optimum performance use the **Change** button to adjust the paging file size to a value of 2 to 3 times of the physical RAM of the PC.
 - **Data Execution Prevention** tab: Enable **Turn on DEP for essential Windows programs and services only**.

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows XP Settings




Figure 11 Performance Options Dialog, Advanced Tab

- b Advanced tab > Startup and Recovery > Settings button (See [Figure 12](#) on page 21):
- **System startup** section: Change both **Time to display...** fields from 30 to 3 sec.
 - **System failure** section: Set **Write debugging information** to none





Figure 12 Startup and Recovery Dialog

- c **Advanced tab > Error Reporting button:** Select **Disable error reporting** and **But notify me when critical error occur**.
- d **Remote tab:** Disable **Allow Remote Assistance invitations to be sent from this computer** and **Allow users to connect remotely to this computer**.
- e **System Restore tab:** Enable **Turn off System Restore on all drives**.
- f **Reboot the system after enabling this setting.**

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows XP Settings

- g Disable hardware acceleration for the video adapter.
- Right-click **desktop** > **Properties** > **Settings Tab** > **Advanced button** > **Troubleshoot tab**
 - Set hardware acceleration to **None**.



- [OPTIONAL] 13 Recycle Bin Properties: (right-click on desktop icon Recycle Bin) Enable the following options (See [Figure 13](#) on page 23):
- a Use one setting for all drives
 - b Do not move files to the Recycle Bin. Remove files immediately when deleted

c Display delete confirmation dialog





Figure 13 Recycle Bin Properties Dialog, Global Tab

[OPTIONAL] 14 Screen layout: (right-click desktop > arrange icons by) Disable **Align to grid**.

1 How to configure Windows XP SP 3 for use with Agilent ChemStation Rev. B.04.03 or higher
Configure Windows XP Settings



2

How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Install Windows Vista 26

Configure Windows Vista Settings 27

This chapter identifies and describes all parameters necessary to configure a non-Agilent Windows Vista SP1 system for use with Agilent ChemStation Rev. B.04.03 or higher.



Install Windows Vista

NOTE

Keep your PC disconnected from the internet until you have installed the appropriate security patches and hot fixes. It's strongly recommended to install the latest security fixes as supported from Agilent Technologies and virus definitions prior to connecting to a network.

- 1 Install Windows Vista from recovery DVD or Windows Vista Business DVD. During the setup of Windows Vista provide the Computer Name, Administrator password and network settings. Select to either join an existing domain or set up the system in workgroup mode.
- 2 Check the Windows version by typing **winver** at the **Start > Start Search** command prompt: If necessary, install Service Pack 1 for Vista. Ignore the option to archive current files and reboot the system.
- 3 To secure your system against viruses please install an AntiVirus program. Running AntiVirus programs might influence the behavior and performance of your computer. Agilent ChemStation was tested to run with Symantecs Norton Antivirus. See [“Scanning for Viruses”](#) on page 94 for more details.

Configure Windows Vista Settings

Some changes within this document are mandatory for the ChemStation to work properly on a Windows Vista system. Some changes will optimize application performance. Other changes will have a graphical or minor impact.

NOTE

To indicate the relative importance of individual settings, each item is categorized as:

MUST: These changes must be applied.

PERFORMANCE: These changes will improve system performance.

OPTIONAL: Most of these changes will affect the graphical display of the application.

NOTE

If User Account Control (UAC) is switched on, some configuration steps will require active confirmation to continue.

- [OPTIONAL] 1 General Layout: (right-click **Start > Properties**)
- a **Start Menu Tab:**
 - Select **Start Menu**
 - In the **Privacy** section select both items
 - b **Start Menu Tab > Customize** button: In **Customize Start Menu** dialog: Disable the following option:
 - **Favorites menu**
 - c Enable in the following options:
 - Computer **Display as a link**
 - **Connect To**
 - Control Panel: **Display as a menu**
 - **Default Programs**
 - Documents: **Display as a link**
 - **Enable context menus and dragging and dropping**
 - Games: **Don't display this item**
 - **Help**

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

- **Highlight newly installed programs**
- **Music: Don't display this item**
- **Network**
- **Open submenu when I pause on them with the mouse pointer**




Figure 14 'Customize Start Menu' Dialog

- **Personal folder: Display as a link**
- **Pictures: Display as a link**

- **Printers**
- **Run command**
- **Search communications**
- **Search favorites and history**
- **Search Files: Search this user's files**
- **Search programs**
- **Sort All Programs menu by name**
- **System administrative tools: Display on the All Programs menu and in the Start menu**
- **Use large icons**

[MUST] 2 Folder Options: (right-click **Start > Explore**)

- a Set **View** to **Details**
- b Choose **Organize > Layout**
 - **Enable Menu Bar**
- c Choose **Tools > Folder Options**
 - **General** tab: for Tasks select **Use Windows classic folders**
 - **View** tab: Disable the option: **Hide extensions for known file types**; Enable the option: **Display the full path in the title bar (classic folders only)**
- d Apply these settings to all folders by selecting the **Apply to Folders** button

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

- [MUST] 3 Got to **Start > All Programs > Internet Explorer**. Perform Windows Update via **Tools > Windows Update** to check for updates and apply all critical security patches (see [Figure 15](#) on page 30).




Figure 15 Downloading Critical Security Updates

- [MUST] 4 Register Windows Vista with Microsoft.
- [MUST] 5 Windows Logon Options: (**Start > search for gpedit.msc**)
- a Navigate to **Local Computer Policy > Computer Configuration > Administrative Templates > System Logon** (see [Figure 16](#) on page 31)
 - Enable **Hide entry points for Fast User Switching** and **Always use classic logon**




Figure 16 Local Group Policy Editor Dialog

- [MUST] 6 Power Options: Start > Control Panel > Power Options (see [Figure 17](#) on page 32)
- a Enable power plan High performance
 - b Click Change Plan settings (see [Figure 18](#) on page 33)
 - c Set the option Put the computer to sleep to Never
 - d Click Change advanced power settings (see [Figure 19](#) on page 34)
Change Turn off hard disk after Settings (Minutes) to Never (0 Minutes)

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

e Save the changes.




Figure 17 Selecting a power Plan




Figure 18 Changing Settings for the Plan: High Performance

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings




Figure 19 Power Options

- [MUST] 7 Local Security Policy: (Start > Programs > Administrative Tools > Local Security Policy)
- Navigate to Security Settings > Local Policies > Security Options
 - Double-click the following policy listed in the right hand panel: **Network Access: Sharing and security model for local accounts**

- c In the displayed dialog select the following item from the drop down list:
Classical - local users authenticate as themselves




Figure 20 Local Security Policy Dialog

- [MUST] 8 Regional and Language Options: (Start > Control Panel > Regional and Language Options) (see Figure 21 on page 36)

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

- a Regional options should be set to **English (United States)** from the drop down list.




Figure 21 Regional and Language Options Dialog

- b If regional format other than **English (United States)** is used, the following settings are mandatory. The settings can be defined by clicking the **Customize this format** button (See [Figure 22](#) on page 37):
- Decimal symbol = . (point)
 - Digit grouping symbol = , (comma)
 - List separator = , (comma)




Figure 22 Customize Regional Options Dialog

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

- c Language for non-Unicode programs: On the **Administrative** Tab click **Change system locale...**. From the drop down list, select **English (United States)** (see [Figure 23](#) on page 38).




Figure 23 Regional and Language Options Dialog, **Administrative** Tab

- [MUST] 9 Power Management for network devices (From **Start** right-click **Network** > **Properties** > **Manage network connections** right-click **Local Area Connection** > **Properties** > **Configure**)

- a On the **Power Management** tab, disable all items (see [Figure 24](#) on page 39).




Figure 24 Network Connection Properties Dialog, Power Management Tab

- [PERFORMANCE] 10 Screen layout: (right-click on desktop > Personalize):
- Click on **Window Color and Appearance**
 - Click on **Open classic appearance properties for more color option.**

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

- c Press **Effects**. Disable all three items.




Figure 25 Effects Dialog

- [PERFORMANCE] 11 System Properties: Start > Control Panel > System: click **Advanced system settings**.
- a On the **Advanced** tab > **Performance** hit the **Settings** button
- **Advanced** tab > **Virtual Memory**: For optimum performance use the **Change** button to adjust the paging file size to a value of 2 to 3 times of the physical RAM on the PC.
 - **Data Execution Prevention** tab: Enable **Turn on DEP for essential Windows programs and services only**.
- b **Advanced** tab > **Startup and Recovery** > **Settings** button (see [Figure 26](#) on page 41):
- **System startup** section: Change both **Time to display ...** fields from **30** to **3** sec.
 - **System failure** section: Set **Write debugging information** to **none**




Figure 26 Startup and Recovery Dialog

c Remote tab

- In the **Remote Assistance** section, disable **Allow Remote Assistance connections to this computer**.
- In the **Remote Desktop** section, enable **Don't allow connections to this computer**.

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings




Figure 27 System Properties Dialog

12 Welcome Center (Start > search for 'gpedit.msc')

- Navigate to **Local Computer Policy > Computer Configuration > Administrative Templates > System > Logon** (see [Figure 28](#) on page 43)
- Enable **Don't display the Getting Started welcome screen at logon**.




Figure 28 Local Group Policy Editor Dialog

[OPTIONAL]

- a Recycle Bin Properties: (right-click on desktop icon Recycle Bin) Enable the following options (See Figure 13 on page 23):
- Custom size; select a size corresponding to approximately 10% of the complete disk space for the drive
 - Check Display delete confirmation dialog

2 How to configure Windows Vista SP1 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows Vista Settings

- b Repeat these steps for all drives of your computer.





Figure 29 Recycle Bin Properties Dialog, General Tab

[OPTIONAL] 13 Screen layout: (right-click on desktop > View) Disable **Align to grid**.



3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Install Windows 7 46

Configure Windows 7 Settings 47

This chapter identifies and describes all parameters necessary to configure a non-Agilent Windows 7 system for use with Agilent ChemStation Rev. B.04.03 or higher.



Install Windows 7

NOTE

Keep your PC disconnected from the internet until you have installed the appropriate security patches and hot fixes. It's strongly recommended to install the latest security fixes as supported from Agilent Technologies and virus definitions prior to connecting to a network.

- 1 Install Windows 7 from recovery DVD or Windows 7 Professional or Enterprise DVD. During the setup of Windows 7 provide the Computer Name, Administrator password and network settings. Select to either join an existing domain or set up the system in workgroup mode.
- 2 To secure your system against viruses please install an AntiVirus program. Running AntiVirus programs might influence the behavior and performance of your computer. Agilent ChemStation was tested to run with Symantecs Norton Antivirus. See [“Scanning for Viruses”](#) on page 94 for more details.

Configure Windows 7 Settings

Some changes within this document are mandatory for the ChemStation to work properly on a Windows 7 system. Some changes will optimize application performance. Other changes will have a graphical or minor impact.

NOTE

To indicate the relative importance of individual settings, each item is categorized as:

MUST: These changes must be applied.

PERFORMANCE: These changes will improve system performance.

OPTIONAL: Most of these changes will affect the graphical display of the application.

NOTE

If User Account Control (UAC) is switched on, some configuration steps will require active confirmation to continue.

- [OPTIONAL] 1 General Layout: (right-click **Start > Properties**)
- a **Start Menu Tab:** In the **Privacy** section select both items
 - b **Start Menu Tab > Customize** button: In **Customize Start Menu** dialog:
 - Disable the following option:
 - **Favorites menu**
 - Enable in the following options:
 - **Computer Display as a link**
 - **Connect To**
 - **Control Panel: Display as a menu**
 - **Default Programs**
 - **Devices and Printers**
 - **Documents: Display as a link**
 - **Enable context menus and dragging and dropping**
 - **Games: Don't display this item**
 - **Help**
 - **Highlight newly installed programs**

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings




Figure 30 Start Menu configuration

- **Music: Don't display this item**
- **Network**
- **Open submenu when I pause on them with the mouse pointer**
- **Personal folder: Display as a link**
- **Pictures: Display as a link**
- **Run command**
- **Search other files and libraries Search public folders**
- **Search programs and Control Panel**

- Sort All Programs menu by name
- System administrative tools: Display on the All Programs menu and in the Start menu
- Use large icons

[MUST] 2 Folder Options: (right-click Start > Open Windows Explorer)

a Choose **Organize > Layout**

- Enable **Menu Bar**
- Enable **Navigation Panel**

b Set **View** to **Details**

c Choose **Tools > Folder Options**


- **View** tab:

Disable the option: **Hide extensions for known file types**;

Enable the option: **Display the full path in the title bar (classic theme only)**

d Apply these settings to all folders by selecting the **Apply to Folders** button

[MUST] 3 Go to **Start > Control Panel > Windows Update**. Check for updates and apply all critical security patches.



[MUST] 4 Register Windows 7 with Microsoft.

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

- [MUST] 5 Select **Change Settings**. Select **Change Settings > Set the Important updates** selection to: **Never check for updates**. Deselect the other update options.




Figure 31 Windows Update Settings

Close the Update screen

- [MUST] 6 Go to **Start > Control Panel > Administrative Tools**. Select **Services**. Right click on the Desktop Window Manager Session Manager. In the context menu select **Properties**. Set the startup type to **Disabled**.




Figure 32 Disable Desktop Service

Confirm by clicking **OK** and close the **Services** window. Close **Administrative Tools**.

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

- [MUST]** 7 Go to **Start > Control Panel > Indexing Options**. Click the **Modify** button. Deselect all drives and locations. Click **OK**. Click **Close**.




Figure 33 Disable file indexing

- [MUST]** 8 Navigate to **Local Computer Policy > Computer Configuration > Administrative Templates > System > Logon**. Enable **Hide entry points for Fast User Switching** and **Always use classic logon**.

- [MUST] 9 Power Options: **Start > Control Panel > Power Options**. See
- Enable power plan **High performance (Show additional plans)**
 - Click **Change Plan settings** (see Figure F on page 28)
 - Set the option **Put the computer to sleep to Never**
 - Click **Change advanced power settings** (see Figure 19 on page 29) Open the nodes for **Hard disk**, **Turn off hard disk after** and set the Minutes to 0 (=Never). Click **OK**
 - Save the changes.




Figure 34 Selecting a power plan

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings




Figure 35 Edit Plan Settings: High Performance



Figure 36 Power Options

- [MUST] 10 Local Security Policy: (Start > Control Panel > Administrative Tools > Local Security Policy)
- Navigate to **Security Settings > Local Policies > Security Options**
 - Double-click the following policy listed in the right hand panel: **Network Access: Sharing and security model for local accounts**

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

- c In the displayed dialog select the following item from the drop down list:
Classical - local users authenticate as themselves




Figure 37 Local Security Policy Dialog: Network

- [MUST] 11 Regional and Language Options: (Start > Control Panel > Regional and Language Options)

- a Regional options should be set to **English (United States)** from the drop down list.




Figure 38 Regional and Language Options Dialog

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

- b If regional format other than **English (United States)** is used, the following settings are mandatory. The settings can be defined by clicking the **Customize this format** button:
- Decimal symbol = . (point)
 - Digit grouping symbol = , (comma)
 - List separator = , (comma)




Figure 39 Customize Regional Dialog

- c Language for non-Unicode programs: On the **Administrative** Tab click **Change system locale...**. From the drop down list, select **English (United States)**





Figure 40 Regional and Language Options Dialog, Administrative Tab

- [MUST] 12 Power Management for network devices (From **Start > Control Panel > Network and Sharing Center**) Select **Change adapter settings**. Right-click **Local Area Connection > Properties > Configure**)

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

- a On the **Power Management** tab, disable all items



- [PERFORMANCE] 13 Screen layout: (right-click on **desktop** > **Personalize**):
- a Click on **Window Color and Appearance**
 - b Click on **Open classic appearance properties for more color option.**

- c Press **Effects**. Disable all three items.




Figure 41 Effects Dialog

- [PERFORMANCE] 14 System Properties: Start > Control Panel > System: click **Advanced system settings**.
- a On the **Advanced** tab > **Performance** hit the **Settings** button
- **Advanced** tab > **Virtual Memory**: For optimum performance use the **Change** button to adjust the paging file size to a value of 2 to 3 times of the physical RAM on the PC.
 - **Data Execution Prevention** tab: Enable **Turn on DEP for essential Windows programs and services only**.

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

b Advanced > Startup and Recovery > Settings button:

- **System startup** section: Change both **Time to display ...** fields from **30** to **3** sec.
- **System failure** section: Set **Write debugging information** to **none**





Figure 42 Startup and Recovery Dialog

c Remote tab

- In the **Remote Assistance** section, disable **Allow Remote Assistance connections to this computer**.
- In the **Remote Desktop** section, enable **Don't allow connections to this computer**.



3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings

15 Welcome Center (Start > search for 'gpedit.msc')

- Navigate to **Local Computer Policy > Computer Configuration > Administrative Templates > System > Logon**
- Enable **Don't display the Getting Started welcome screen at logon**.





Figure 43 Local Group Policy Editor Dialog

[OPTIONAL] 16 Recycle Bin Properties: (right-click on desktop icon Recycle Bin) Enable the following options

- Custom size; select a size corresponding to approximately 10% of the complete disk space for the drive
- Check Display delete confirmation dialog


- a Repeat these steps for all drives of your computer.



[OPTIONAL] 17 Screen layout: (right-click on desktop > View). Disable Align Icons to grid.

3 How to configure Windows 7 for use with Agilent ChemStation Rev. B.04.03 or higher

Configure Windows 7 Settings



4 Security Settings

Security Settings for LC-MSD and CE-MSD Systems	68
Security Settings for LC-MSD Trap Systems	72
Security Settings for GC 7890A Systems	77
Security Settings for IO Library Suite 15.5	80
Security Settings for Systems using Agilent Bootp Service	84

This chapter summarizes all security settings for LC/MSD, CE/MSD, LC/MSD Trap, GC 7890A, IO Library Suite 15.5, and Agilent Bootp Service.



Security Settings for LC-MSD and CE-MSD Systems

Security Settings for LC/MSD and CE/MSD Systems

The MSD part of LC/MSD and CE/MSD systems requires exceptions for the Windows Firewall in case the firewall is enabled. Additionally, the communication with the MSD is sensitive to the order NICs in the **Adapters and Bindings** dialog. These settings are described in the following passages.

NOTE

These settings should only be applied on PCs controlling LC/MSD or CE/MSD systems.

Internet Security Settings:

Start > Settings > Control Panel > Security Center

- 1 Click **Manage security settings for Windows Firewall**
- 2 **General** tab (See [Figure 44](#) on page 69)
 - a Select option **On (recommended)**

b Disable option **Don't allow exceptions**




Figure 44 Windows Firewall Dialog, General Tab

3 Exceptions tab (See Figure 45 on page 70)

- a Select option **Display a notification when Windows Firewall blocks a program.**
- b Add **ChemMain.exe** to the list of **Programs and Services**
- c Enable an exception for this program.

4 Security Settings

Security Settings for LC-MSD and CE-MSD Systems

- d Disable exceptions for all other programs listed by default.




Figure 45 Windows Firewall Dialog, Exceptions Tab

Advanced Network Settings:

Start > Settings > Control Panel > Network Connections

- 1 Select menu **Advanced > Advanced Settings**
- 2 **Adapters and Bindings** tab (see [Figure 46](#) on page 71)
 - a Make sure that the Local Area Connection pertaining to the LC/MS or CE/MS NIC is the first item in the list of connections

NOTE

The names of your LAN cards may differ from those shown in the example. You can tell which LAN card is the Instrument LAN by comparing the IP addresses assigned to the LAN cards.

The LAN pertaining to the company intranet > internet will typically have an IP address assigned by organizations Static IP policy or by a DHCP server. Check with your network administrator.




Figure 46 Advanced Settings Dialog, Adapters and Bindings Tab

Security Settings for LC-MSD Trap Systems

Security Settings for LC/MSD Trap Systems

The MSD part of LC/MSD Trap systems requires exceptions for the Windows Firewall in case the firewall is enabled. These settings are described in the following passages.

NOTE


These settings should only be applied on PCs controlling LC/MSD Trap systems.

Internet Security Settings

Start > Settings > Control Panel > Security Center

1 Click **Manage security settings for Windows Firewall**

- 2 On the **General** tab (See below):
 - Select option **On (recommended)**
 - Disable option **Don't allow exceptions**




4 Security Settings


Security Settings for LC-MSD Trap Systems

3 On the **Exceptions** tab (See below):

- Select option **Display a notification when Windows Firewall blocks a program.**



- 4 On the **Advanced** tab (see below): In the **Network Connection Settings** section, de-select the Local Area Connection pertaining to the LC/MS NIC. Typically, this is **Local Area Connection 2**.



4 Security Settings

Security Settings for LC-MSD Trap Systems

- 5 From the list of Resources select **Change the way Security Center alerts me** (See below)
 - Disable **Automatic Updates**
 - Disable **Firewall**
 - Enable **Virus Protection**



NOTE

With these settings, you will not receive any alerts regarding the state of Windows Firewall. Be sure that you monitor the state of Windows Firewall yourself

Security Settings for GC 7890A Systems

Security Settings for GC 7890A Systems

GC 7890A systems require exceptions for the Windows Firewall in case the firewall is enabled. These settings are described in the following sections.

NOTE

These settings should only be applied on PCs controlling GC 7890A systems.

Internet Security Settings


Start > Settings > Control Panel > Security Center

1 Click **Manage security settings for Windows Firewall**


4 Security Settings

Security Settings for GC 7890A Systems

- 2 On the **General** tab (See below):
 - Select option **On (recommended)**
 - Disable option **Don't allow exceptions**



- 3 On the **Exceptions** tab (See below):
- Select option **Display a notification when Windows Firewall blocks a program**.
 - Add **AgtCSRCAModuleDrv.exe** to the list of **Programs and Services** and enable an exception for this program.
 - Disable exceptions for all other programs listed by default.



Security Settings for IO Library Suite 15.5

Security Settings for IO Library Suite 15.5

Using the IO Library Suite 15.0 requires exceptions for the Windows Firewall in case the firewall is enabled. These settings are described in the following sections.

NOTE


These settings should only be applied on PCs having the IO Library Suite 15.0 installed.

Internet Security Settings

Start > Settings > Control Panel > Security Center

1 Click **Manage security settings for Windows Firewall**

- 2 On the **General** tab (See below):
- Select option **On (recommended)**
 - Disable option **Don't allow exceptions**



4 Security Settings


Security Settings for IO Library Suite 15.5

3 Exceptions tab (See below)

- Select option **Display a notification when Windows Firewall blocks a program.**
- Add **Siclland.exe** to the list of **Programs and Services**
- Enable an exception for this program.
- Add the following ports to the list of **Programs and Services**:
 - LxiAllow: 5044 (UDP)
 - LxiAllow: 5044 (TCP)
 - LxiAllow: 320 (UDP)
 - LxiAllow: 319 (UDP)
 - SunRpcAllow: 111 (UDP)
 - SunRpcAllow: 111 (UPD)
- Enable exceptions for these ports.
- Disable exceptions for all other programs listed by default.

NOTE

The above listed programs and ports are added to the list of exceptions for the firewall during installation of the IO Library Suite 15.5 per default.



4 Security Settings


Security Settings for Systems using Agilent Bootp Service

Security Settings for Systems using Agilent Bootp Service

If you use Agilent Bootp Service for administration of IP addresses for Agilent instruments on a LAN, your system requires exceptions for the Windows Firewall in case the firewall is enabled. These settings are described in the following passages.

NOTE

These settings should only be applied on PCs using Agilent Bootp Service.



5 Additional Modifications

Install ChemStation 86

Profiles and Group Policies to run ChemStation as PowerUser/User 87

These modifications can only be performed using an Administrator user account

Install ChemStation

- 1 Install the ChemStation Software to C:\Chem32 (logged on with **Administrator** user account).


Profiles and Group Policies to run ChemStation as PowerUser/User

The ChemStation system needs to be installed using an operating system Administrator user account. To run the system with the configured options outlined in this document, all Users and PowerUsers should use the same settings.

Some of the outlined settings can be simplified by using the Windows Group Policy or by copying user profiles in order to simplify the User/PowerUser setup. But some changes (e.g. performance settings) can only be set manually using the Administrator account.

5 Additional Modifications

Profiles and Group Policies to run ChemStation as PowerUser/User



6 Maintaining the PC's File System

Scanning and Repairing NTFS Volumes	90
Defragmenting an NTFS Volume	92
Scanning for Viruses	94
Cleaning up Left-Over Temporary Files	95
Performing Backups	97
Windows XP Backup	97
Automating and Scheduling Windows XP Backup	97
Travan Tape Drives	98
Digital Audio Tape (DAT) Drives	98
Magneto-Optical (MO) Disk Drives	98
Recordable CD-ROM	99
Recordable DVD	99
Network Servers	99

To help to optimize the performance of your ChemStation system it is important that you schedule routine maintenance (e.g. deleting unnecessary files). The following instructions should be used as a part of your routine maintenance schedule.


NOTE

You will need Administrator capabilities to perform these tasks.



Scanning and Repairing NTFS Volumes

- 1 Double-click **My Computer** on the desktop.
- 2 Select the C: drive.
- 3 Select **File > Properties**.
- 4 From the **Tools** tab, click **Check Now** in the **Error-checking** section.



- 5 Select both **Check disk** options **Automatically fix file system errors** and **Scan for and attempt recovery of bad sectors**.

- 6 Click the **Start** button.



- 7 The warning message **The disk check could not be performed because the disk check utility needs exclusive access to some Windows files on the disk. These files can be accessed only by restarting Windows. Do you want to schedule this disk check to occur the next time you restart the computer?** is displayed. Press **Yes**.



- 8 Restart the computer. During startup, **Check Disk** will be performed.
- 9 Repeat the steps outlined above for the other drives on your computer.


Defragmenting an NTFS Volume

Compared to the FAT file system, fragmentation of files on an NTFS volume is greatly reduced. By design, fragmentation only occurs on an NTFS volume if a file increases in size after it is on the drive. This means that backing up a volume to tape and restoring it from tape will result in an unfragmented volume.


Windows XP Professional includes a defragmentation utility.

Perform the following steps to use the defragmentation utility provided by Microsoft Windows:

- 1 Select **Start > Programs > Accessories > System Tools > Disk Defragmenter**.
- 2 Select the drive (e.g. C:) you want to analyze and click the **Analyze** button.



- 3 If the result of the analysis indicates that the drive needs defragmentation, click **Defragment**. Otherwise, defragmentation is not necessary at this time, and you can click **Close**.



CAUTION

Defragmentation of your harddisk while ChemStation is acquiring data can cause loss of data.

- Never defragment your hard disk drive during data acquisition

- 4 Repeat these steps for the other drives on your system.

Scanning for Viruses

If the Agilent ChemStation computer is connected to the network in order to communicate with other systems (LAN, modem, BBS, etc.) or if you share files with other computer users, we recommend scanning the drives for computer viruses. If you download programs from the internet or exchange documents or programs through e-mail or flexible disks, a virus-check utility is mandatory.

Many different utilities are available for virus scanning. Many of these utilities also come with periodic updates as new viruses are discovered. You have to be aware that these utilities will degrade the performance of the hard disk when scanning for viruses. Additionally, other side effects may influence the behavior of your system.

Agilent ChemStation was tested to work together with Symantecs Norton Antivirus. We recommend that you scan all disks of your ChemStation computer at least once per week for viruses. It is essential that you perform the LiveUpdate to obtain the latest virus definitions as provided by Symantec as reaction on new viruses.

NOTE

Never perform a virus scan during data acquisition, as this will degrade the performance of the hard disk and the acquired data might not be written to the hard disk fast enough, causing a power fail or buffer overrun message in the instrument log book and possible loss of data.

NOTE

With Norton AntiVirus versions prior to Norton Antivirus 2005, Windows Security Center will be aware that you have Norton AntiVirus installed, but it will not be able to determine its working status. Please either refer to Windows for help on how to disable the Security Alert for Virus Protection and monitor Norton AntiVirus yourself or refer to Symantec support to obtain an update which adds compatibility to Windows Security Center.

Cleaning up Left-Over Temporary Files

Occasionally, temporary files may accumulate in the directory specified by the TEMP environment variable. These files generally are left open when Windows is abnormally terminated, for example, by powering your computer down without closing Windows first. The Windows temporary files are named -XXXXXXX.TMP, where XXXXXXXX are characters and numbers produced by the program that created the temporary file. To recover the temporary space, you should delete these files after closing all currently running applications.

To determine which directory is currently being used for temporary files, type **SET** at the command prompt. This will echo back the current settings for all user and system environmental variables.

- 1 Select **Start > Programs > Accessories > System Tools > Disk Cleanup**.
- 2 Select the drive (e.g. C:) you want to clean up.




- 3 Select **Temporary Files** and **Temporary Internet Files**. You may select other types of files to delete, depending on your understanding of these file types and the way you use your PC. It is possible to delete vital files with the Disk Cleanup utility, so do not delete file types that you do not understand. Do

6 Maintaining the PC's File System

Cleaning up Left-Over Temporary Files

not select the option **Compress old files** on a ChemStation computer. Press **OK**.



4 Repeat these steps for the other drives on your system.

Performing Backups

In order to protect your valuable data, an appropriate backup strategy is mandatory. We recommend to archive your data files in regular intervals. If your data is no longer required to be available online, you should also archive them and remove them from your hard disk. This will free up valuable storage space on your computers hard disk as data acquired from a diode array or mass spectrometer detector may easily exceed several megabytes in size.

Windows XP Backup

Windows XP Professional includes compatible backup programs that support many standard tape devices. The backup utility also allows for backing up to a single file, which can then be placed on an offline medium. Both backup utilities allow the Windows Registry to be backed up as well as data.

Automating and Scheduling Windows XP Backup

Windows XP allows automation and scheduling of Backup jobs through the use of batch files and the Scheduler service. The Scheduler is available by default for Windows XP.

CAUTION

Do not perform a backup while ChemStation is acquiring data.

This will degrade the performance of the hard disk and the acquired data might not be written to the hard disk fast enough, causing a power fail or buffer overrun message in the instrument log book and possible loss of data. It may also result in inconsistencies in the file structure stored on the backup medium.

→ Be sure to perform any backup with a Windows-based backup utility during times where no data acquisition is scheduled.

Travan Tape Drives

Tape drives are the most common backup devices available today. The IDE tape drives combine lowest initial hardware cost with acceptable performance (up to 4 GB > hour). Travan tapes are a reasonable backup solution for most laboratories that run their Agilent ChemStation on single workstation PCs. The minimum storage capacity per tape is currently around 2.5 GB of uncompressed data.

You will need an appropriate backup utility that supports the drive.

NOTE

Never run an MS-DOS-based tape backup utility from Windows. This may result in corrupt file information on the tape.

Digital Audio Tape (DAT) Drives

Digital audio tape (DAT) drives are suited to handle large amounts of data (storage capacity is typically 2 to 6 GB or more), but require a SCSI interface installed in the PC.

Magneto-Optical (MO) Disk Drives

MO disk drives are high storage capacity devices (at least 1.3 GB) connected externally to a SCSI interface that must be installed in your PC. MO disk drives are capable of reading and writing *write-once* media and rewritable disks.

No specific backup software is required. The supporting drivers are shipped with the operating system or the device and allow for accessing the additional MO disk as if it was a floppy disk.

NOTE

If the rewritable disk is formatted with NTFS, it cannot be hot-swapped, i.e. ejected and exchanged while the system is running. To change disks, you need to shut down Windows.

Recordable CD-ROM

CD-R (CD-Recordable) systems allow you to record up to 700 > 800 MB of data on a single CD-R disk. The CD reader > recorder is usually connected to your PC the IDE interface. The CD-R media have a very long life, making them ideally suited to archiving. The CD-R disks can be read by other computers fitted with CD-ROM drives and appropriate software.

Recordable DVD

Recordable DVDs are available in several formats, e.g. DVD+R and DVD-R. They allow you to record up to 4.7 GB (Single Layer) > 8.5 GB (Double Layer) of data on a single DVD. The DVD reader > recorder is usually connected to your PC the IDE interface. DVD±R disks have life times comparable to CD-ROM disks. In contrast to some other DVD formats, DVD±R disks are not rewritable and thus are especially suited for long-term archiving. The DVD±R disks can be read by other computers fitted with DVD drives and appropriate software.

Network Servers

Network servers can also be used for data backup purposes, provided that the network system administrator has implemented an appropriate backup strategy for data stored on the servers disk. One strategy is to allocate a portion of the servers hard disk space for your Agilent ChemStation. This portion on the server is called a *share* which you can use for your standard file operations (access, copy, evaluate files) within Windows.

The advantage of this approach is that backup will be performed centrally by a system administrator for all users that have accounts (logon capabilities) on the system.

NOTE

Using a 3.5-inch floppy disk drive is not a reasonable backup solution for the Agilent ChemStation.

Index

A

anti virus 8, 26, 46, 94

B

backup
strategy 97

C

CD-ROM
recordable 99
configure settings
Windows 7 47
Windows Vista 27
Windows XP 9

D

defragmentation 92
digital audio tape drives 98
disc clean-up 95
drives
for backup 97

F

folder options
Windows 7 49
Windows Vista 29
Windows XP 11

M

magneto-optical disk drives 98

N

network

server 99

NTFS volumes 90

P

Power Options
Windows XP 12
power user 87

R

recordable CD-ROM 99
repairing
NTSF 90

S

security settings
Bootp service 84
GC 7890A 77
LC/MSD and CE/MSD 68
security settings
LC/MSD trap 72
security 8
server
network 99

T

temporary files 95
travan tape
drive 98

U

user profiles 87

V

virus
scan 94

W

Windows 7
configure 47
installation 46
Windows Vista
configure 27
installation 26
Windows XP
configure 9
installation 8
write-once media 98

In This Book

This manual describes how to configure a non-Agilent PC Bundle system.

The first three chapters of this document, identify and describe all parameters necessary to configure a non-Agilent PC Bundle system as well as important settings such as Internet Security.

The fourth and fifth chapter, describe instrument specific security settings and additional configuration tasks, respectively.

The last chapter, describes how to maintain your computer to keep its performance in an optimal status.

© Agilent Technologies 2004, 2005-2010

Printed in Germany
09/10

G2170-90136
G2170-90136

G2170-90136